

Restek's PLOT Column Family —

The Benchmark For Performance!

- Innovative bonding process minimizes particle release.
- More consistent flow means stable retention times.
- Outstanding peak symmetry improves impurity analysis.

See what makes Restek's new **Rt®-Silica BOND** column the best on the market! **p. 11**

PLOT Columns Available In:

Fused Silica

Rt[®] Column Phases:

Rt[®]-Alumina BOND/MAPD

Rt[®]-Alumina BOND/Na₂SO₄

Rt[®]-Alumina BOND/KCl

Rt[®]-Alumina BOND/CFC

Rt[®]-MSieve 5A

Rt[®]-Silica BOND **NEW!**

Rt[®]-Q-BOND

Rt[®]-QS-BOND

Rt[®]-S-BOND

Rt[®]-U-BOND

Metal

MXT[®] Column Phases:

MXT[®]-Alumina BOND/MAPD

MXT[®]-Alumina BOND/Na₂SO₄

MXT[®]-MSieve 5A

MXT[®]-Q-BOND

MXT[®]-S-BOND

Next Generation of Porous Layer Open Tubular (PLOT) Columns


- Stabilized particle layers improve robustness and reproducibility of retention and flow.
- Fully compatible with valve switching and Deans switching systems.
- Highly efficient, reproducible analyses; ideal for permanent gases, solvents, and hydrocarbons.
- Innovative manufacturing procedure reduces particle generation and improves performance of porous polymer and molecular sieve PLOT columns.
- Wound on a 7"-diameter, 11-pin cage unless otherwise noted.

Porous layer open tubular (PLOT) columns are very beneficial for solving application problems, especially for the analysis of volatile compounds. PLOT columns have a unique selectivity, allowing for the separation of gaseous compounds at room temperature. Due to the adsorption mechanism of the supports used in PLOT columns, permanent gases and light hydrocarbons can be resolved at room temperature; columns can then be programmed to higher temperatures to elute higher boiling compounds.

Traditional PLOT Columns Offer Poor Stability

The traditional PLOT column is built with a 5–50 μm layer of particles adhered to the tubing walls. Because this layer of particles generally lacks stability, PLOT columns must be used very carefully, as particle release is common and can cause unpredictable changes in retention time and flow behavior. Traditional PLOT columns also must generally be used in conjunction with particle traps to prevent the contamination of valves, injectors, and GC detectors. Detectors contaminated with particles typically generate electronic noise, which shows up chromatographically as a spike in the baseline. In extreme cases, detector flow can be obstructed by particle buildup. Particles can also affect valves by becoming lodged in the valve and causing leaks or restricting flow. Figure 1 shows an example of blockage caused by particle accumulation inside a Press-Tight[®] connector.

Figure 1: Particles released from traditional PLOT columns can cause blockages.



Restek[®] PLOT Columns Offer Improved Stability to Minimize Particle Release

Restek has developed technology and procedures to manufacture PLOT columns with concentric stabilized adsorption layers. These next generation PLOT columns show a constant flow behavior (permeability) and have significantly improved mechanical stability, resulting in easier operation, better chromatography, and reduced particle release. Greater particle stability means more reproducible retention times, virtually no spiking, and longer column lifetimes. This innovative Restek[®] stabilization chemistry is currently applied to all fused silica and metal PLOT columns featured in this brochure.

Consistent Flow Restriction Factor (F) Guarantees Reproducible Flow

Thick layers of particles are difficult to deposit in a homogeneous layer, and in traditionally manufactured PLOT columns, this results in variable coating thicknesses. The positions where the layer is thicker act as restrictions and affect flow (Figure 2). Depending on the number and intensity of these restrictions, traditional PLOT columns often show greater variation in flow restriction than wall coated open tubular (WCOT) columns. In practice, conventional PLOT columns with the same dimensions can differ in flow by a factor of 4 to 6 when operated at the same nominal pressure. For applications where flow is important, such as with Deans switching, the nonreproducible flow behavior of most commercially available PLOT columns is a problem.


Figure 2: Inconsistent coating thicknesses result in restrictions that cause significant variation in flow.



In order to measure flow restriction reproducibility, Restek introduced a new factor: the flow restriction factor (F). This factor is based on the retention time of an unretained marker compound, as measured on both coated and uncoated tubing using the same backpressure setting (Equation 1). For quality control purposes, methane is used as the marker when evaluating porous polymer columns, and helium is used for testing molecular sieve 5A columns.

Flow restriction factor determination can be used to assess both the degree of column restriction and the reproducibility of the column coating process. Flow restriction can also be calculated (Equation 2). Figure 3 shows typical results for PLOT columns manufactured using a conventional process. Because of the difference in flow restriction, individual columns have very different flow characteristics. In contrast, Figure 4 shows results for columns made using our Rt[®]-QS-BOND (bonded porous polymer) PLOT column process. Clearly, our manufacturing process results in greater consistency in both column coating thickness and flow restriction, which results in more stable retention times and better performance in Deans and related flow switching techniques. Flow restriction factors are specified on the certificate of analysis (CofA) included with every Restek[®] PLOT column, and the values are listed on the report.

Figure 3: Traditional PLOT columns show significant flow variability, indicating inconsistent column coating thicknesses.




Equation 1: Flow restriction factor (F) is used to demonstrate coating consistency.

$$F = \frac{t_{R1} \text{ of unretained component (uncoated tubing)}}{t_{R2} \text{ of unretained component (coated column)}}$$

t_R = retention time

Note: F values will always be <1 as the coated column always has more restriction than the uncoated column.

Figure 4: PLOT columns from Restek offer consistent flow restriction, giving more reproducible results column-to-column.



Equation 2: Percent flow restriction of coated column.

$$\% \text{ restriction} = (1 - F) \times 100$$

See what makes Restek's new **Rt[®]-Silica BOND** column the best on the market! **p. 11**

Restek's PLOT columns are exceptionally robust, featuring concentric stabilized coating layers. They allow for more consistent gas flows and are recommended for applications sensitive to variation in retention time or flow. These PLOT columns are a significant advance in technology and are ideal for efficient, reproducible analyses of permanent gases, solvents, and hydrocarbons.

PLOT Column Phase Cross-Reference: Similar Selectivity

Restek [®] Rt [®] and MXT [®] Columns	Porous Layer	Supelco	Alltech	Agilent (J&W, Varian, Chrompack)	Quadrex
Alumina BOND/Na ₂ SO ₄	Aluminum oxide	Alumina-Sulfate	AT-Alumina	GS-Alumina, CP-Al ₂ O ₃ /NA ₂ SO ₄	—
Alumina BOND/KCl	Aluminum oxide	Alumina-Chloride	—	GS-Alumina KCl, HP PLOT Al ₂ O ₃ , CP-Al ₂ O ₃ /KCl	—
Alumina BOND/CFC	Aluminum oxide	—	—	unique product	—
Alumina BOND/MAPD	Aluminum oxide	—	—	Select Al ₂ O ₃ , MAPD	—
M sieve 5A	Molecular sieve 5A	Molsieve 5A	AT-Molesieve	HP PLOT Molesieve, CP-Molesieve 5A	PLT-5A
Q-BOND	100% Divinylbenzene	Supel-Q-PLOT	AT-Q	HP PLOT Q, CP-PoraPlot Q, PoraBond Q	—
QS-BOND	Intermediate polarity porous polymer	—	—	GS-Q	—
S-BOND	DVB vinylpyridine polymer	—	—	CP-PoraPlot S	—
U-BOND	DVB ethylene glycol-dimethylacrylate polymer	—	—	HP PLOT U, CP-PoraPlot U, CP-PoraBond U	—
Silica BOND	Bonded silica	—	—	CP Silica PLOT, GS-GasPro	—



i tech tip

Traces of water in the carrier gas and samples will affect the retention and selectivity of alumina. If exposed to water, the retention times will shorten. The column can be regenerated by conditioning for 15–30 minutes at 200 °C under normal carrier gas flow. Periodic conditioning ensures excellent run-to-run retention time reproducibility.

Unless noted, the maximum programmable temperature for an Rt®-Alumina BOND column is 200 °C. Temperatures higher than the stated maximum temperature can cause irreversible changes to the porous layer adsorption properties.

Rt®-Alumina BOND Columns

Restek® Rt®-Alumina BOND columns are highly selective for C1–C5 hydrocarbons and separate all saturated and unsaturated hydrocarbon isomers above ambient temperatures. The reactivity of the aluminum oxide stationary phase is minimized to improve column response for polar unsaturates, such as dienes, and the column's sensitivity (or response) ensures linear and quantitative chromatographic analysis for these compounds. Strong bonding prevents particle generation and release, which allows valve switching without harming the injection or detection systems. And because they are stable up to at least 200 °C, Rt®-Alumina BOND columns can be regenerated to restore full efficiency and selectivity by conditioning at their maximum temperature if water is adsorbed. High capacity and loadability give you exceptionally symmetrical peaks, making these columns ideal for volatile hydrocarbon separations at percent levels, as well as impurity analyses at ppm concentrations. Restek® Rt®-Alumina BOND PLOT columns are manufactured on fused silica tubing; select phases are also available on metal MXT® tubing.

Guaranteed Reproducibility

To ensure reproducible retention times and predictable flow behavior column-to-column, each Rt®-Alumina BOND column is extensively tested. A hydrocarbon test mix confirms proper phase retention and selectivity. To calculate *k* (retention or capacity factor), which is a measure of phase retention, 1,3-butadiene is used, while selectivity is measured using retention indices for propadiene and methyl acetylene. The resolution of *trans*-2-butene and 1-butene is also verified, and to measure efficiency, plates per meter are checked using 1,3-butadiene.

Rt®-Alumina BOND/Na₂SO₄ Columns (fused silica PLOT) (Na₂SO₄ deactivation)

- Acetylene and propadiene elute after butanes.
- Best separation for butene isomers (impurities in butene streams).
- Methyl acetylene elutes after 1,3-butadiene.
- Cyclopropane (impurity in propylene) elutes well before propylene.
- Also available on metal MXT® tubing! (See page 12.)

ID	ID	temp. limits	30-Meter cat.#	50-Meter cat.#
0.25 mm	4 µm	to 200 °C	19775	—
0.32 mm	5 µm	to 200 °C	19757	19758
0.53 mm	10 µm	to 200 °C	19755	19756


similar phases

GS-Alumina, CP-Al₂O₃/
Na₂SO₄, Alumina-Sulfate,
AT-Alumina

did you know?

Restek draws our own fused silica tubing and applies our own proprietary stationary phases. By fully managing our production streams, we are able to ensure unparalleled reliability and stability.

Refinery Gas on Rt®-Alumina BOND (Na₂SO₄)



Column Rt®-Alumina BOND/Na₂SO₄, 50 m, 0.53 mm ID, 10 µm (cat.# 19756)
Sample Refinery gas
Injection 10 µL split
Inj. Vol.: Taper (2 mm) (cat.# 20795)
Liner: 200 °C
Inj. Temp.:
Split Vent
Flow Rate: 80 mL/min

Oven
Oven Temp: 45 °C (hold 1 min) to 200 °C at 10 °C/min (hold 3.5 min)
Carrier Gas H₂, constant pressure (8.0 psi, 55.2 kPa)
Linear Velocity: 74 cm/sec @ 45 °C
Detector FID @ 200 °C



Rt®-Alumina BOND/KCl Columns (fused silica PLOT)
(KCl deactivation)


- Restek's lowest polarity alumina column.
- Low moisture sensitivity reduces the need for frequent regeneration.
- Acetylene elutes before *n*-butane.
- Methyl acetylene (impurity in 1,3-butadiene) elutes before 1,3-butadiene.

similar **phases**

GC-Alumina KCl, HP-PLOT
Al₂O₃/KCl, CP-Al₂O₃/KCl,
Alumina-Chloride

ID	df	temp. limits	30-Meter cat.#	50-Meter cat.#
0.25 mm	4 µm	to 200 °C	19776	—
0.32 mm	5 µm	to 200 °C	19761	19762
0.53 mm	10 µm	to 200 °C	19759	19760

Refinery Gas on Rt®-Alumina BOND/KCl



GC_PC01088

Column Rt®-Alumina BOND/KCl, 50 m, 0.53 mm ID, 10 µm (cat.# 19760)
Sample Refinery gas standard
Injection
 Inj. Vol.: 10 µL split
 Liner: Taper (2 mm) (cat.# 20795)
 Inj. Temp.: 200 °C
 Split Vent
 Flow Rate: 80 mL/min
Oven
 Oven Temp: 45 °C (hold 1 min) to 200 °C at 10 °C/min (hold 3.5 min)
Carrier Gas H₂, constant pressure (8.0 psi, 55.2 kPa)
 Linear Velocity: 74 cm/sec @ 45 °C
Detector FID @ 200 °C
 Make-up Gas
 Type: N₂
 Data Rate: 20 Hz
Instrument HP5890 GC

did you **know?**

All Restek® PLOT columns come standard on a 7"-diameter, 11-pin cage. Metal MXT® columns are also available coiled to 3.5" diameter by adding the suffix -273 to the part number (See page 12). If you need more information, please call your local Restek® representative.



For more chromatograms, search our extensive library at
www.restek.com/chromatograms



Rt®-Alumina BOND/CFC Columns (fused silica PLOT)

- Improved inertness for chlorofluorocarbon (CFC) compounds.
- Highly selective alumina-based column separates most CFCs.
- High retention and capacity for CFCs.

Restek
exclusive!

i tech tip

Especially when valve switching or backflushing is used, Restek recommends using particle traps to help prevent detector spikes and valve rotor scratches.

Visit www.restek.com/plot for specialized PLOT column particle traps.

The Alumina BOND/CFC adsorbent is ideal for retaining halogenated compounds, especially CFCs (chlorinated fluorocarbons) like Freon® products. It offers high selectivity, allowing a wide range of CFC isomers to be resolved at above ambient temperatures. The Rt®-Alumina BOND/CFC column is thoroughly deactivated to reduce the reactivity of alumina. Even though there is still some residual reactivity for some mono- or di-substituted CFCs, the majority of these compounds can be accurately quantified from main stream processes or in impurity analyses.


ID	df	temp. limits	30-Meter cat.#
0.53 mm	10 µm	to 200 °C	19763

Impurity Analysis of 1,1,1,2-Tetrafluoroethane (CFC-134a) on Rt®-Alumina BOND/CFC

Rt®-Alumina BOND/CFC columns
provide the best selectivity for
impurity analysis of CFCs.


Peaks

1. Chloropentafluoroethane (CFC-115)
2. Dichlorodifluoromethane (CFC-12)
3. Chlorodifluoromethane (CFC-22)
4. 1,1,1,2-Tetrafluoroethane (CFC-134a)



Column Rt®-Alumina BOND/CFC, 30 m, 0.53 mm ID (cat.# 19763)
Sample 1,1,1,2-Tetrafluoroethane
Conc.: Neat
Injection
Inj. Vol.: 500 µL split
Oven
Oven Temp.: 80 °C (hold 6 min) to 140 °C at 10 °C/min (hold 2 min)
Carrier Gas He
Detector FID
Notes Gas sampling, purity analysis

Note that tailing peaks are common in CFC analyses due to overloading normally employed for this type of work.



For more chromatograms, search our extensive library at
www.restek.com/chromatograms




Rt®-Alumina BOND/MAPD Columns (fused silica PLOT)

- Optimized deactivation produces maximum response when analyzing trace levels of acetylene, methyl acetylene, and propadiene.
- Stable response factors make this column ideal for process-type applications where recalibration must be minimized.
- High loadability reduces peak tailing and improves separations.
- Extended temperature range up to 250 °C for fast elution of high molecular weight (HMW) hydrocarbons and accelerated column regeneration following exposure to water.
- Also available on metal MXT® tubing! (See page 12.)

Restek's R&D chemists have optimized the deactivation technology applied to our Rt®-Alumina BOND/MAPD column for improved analysis of trace concentrations of polar hydrocarbons like acetylene, methyl acetylene, and propadiene in hydrocarbon streams containing higher levels of C1–C5 hydrocarbons. Our Alumina PLOT deactivation produces an incredibly inert column that offers superior reproducibility and stable response factors to maximize the number of analyses before recalibration is required. Its high sample capacity reduces peak tailing, thereby improving the separation of target compounds. In addition, a 250 °C maximum operating temperature lets you more quickly elute hydrocarbons up to dodecane and reduces regeneration time when the column is exposed to water from samples or carrier gases.

ID	df	temp. limits	30-Meter cat.#	50-Meter cat.#
0.25 mm	4 µm	to 250 °C	19781	—
0.32 mm	5 µm	to 250 °C	19779	19780
0.53 mm	10 µm	to 250 °C	19777	19778

1,3-Butadiene on Rt®-Alumina BOND/MAPD (Purity Analysis)



Column	Rt®-Alumina BOND/MAPD, 50 m, 0.53 mm ID, 10.0 µm (cat.# 19778)	Oven	70 °C (hold 5 min) to 200 °C at 10 °C/min (hold 0 min)
Sample Injection	Crude 1,3-butadiene	Carrier Gas	He, constant pressure (20 psi, 137.9 kPa)
Inj. Vol.:	10 µL split	Detector	FID @ 250 °C
Liner:	2.0 mm ID straight inlet liner (cat.# 20712)	Make-up Gas Flow Rate:	30 mL/min
Inj. Temp.:	200 °C	Make-up Gas Type:	N ₂
Split Vent Flow Rate:	100 mL/min	Data Rate:	20 Hz
		Instrument	HP5890 GC

did you know?

All Restek® PLOT columns come standard on a 7"-diameter, 11-pin cage. Metal MXT® columns are also available coiled to 3.5" diameter by adding the suffix -273 to the part number. If you need more information, please call your local Restek® representative.



also available!
Metal MXT®
PLOT Columns
 See page 12 for more information.



did you know?

Rt[®]-Msieve 5A PLOT columns are designed for efficient separation of Ar/O₂ and other permanent gases, including CO.

itech tip

Because molecular sieve materials are very hydrophilic, they will adsorb water from the sample or carrier gas. Water contamination can have a detrimental effect on peak symmetry and can reduce the resolution of all compounds. If water contamination occurs, simply regenerate your Rt[®]-Msieve 5A PLOT column by conditioning at 300 °C with dry carrier gas flow for 3 hours.



also available! Metal MXT[®] PLOT Columns

See page 12 for more information.

Molecular Sieve 5A PLOT Columns


Restek's molecular sieve 5A PLOT columns are designed for efficient separation of Ar/O₂ and other permanent gases, including CO. Special coating and deactivation procedures ensure chromatographic efficiency and the integrity of the porous layer coating. Molecular sieves have very high retention, allowing separations of permanent gases at temperatures above ambient. Our deactivation technology also allows CO to elute as a sharp peak. Additionally, our unique immobilization process guarantees that the uniform particles remain adhered to the tubing—even after continuous valve-cycling.

Our revolutionary molecular sieve 5A PLOT columns separate Ar/O₂ and H₂/He at ambient temperature or above (see chromatogram). These columns also are an excellent choice for rapid separation of permanent gases in refinery or natural gas. Available on fused silica and metal MXT[®] tubing. (See page 12.)

Rt[®]-Msieve 5A Columns (fused silica PLOT)

ID	df	temp. limits	15-Meter cat.#	30-Meter cat.#
0.25 mm	20 µm	to 300 °C	19773	—
0.32 mm	30 µm	to 300 °C	19720	19722
0.53 mm	50 µm	to 300 °C	19721	19723

Separation of Argon/Oxygen and Other Permanent Gases on Rt[®]-Msieve 5A



Column Rt[®]-Msieve 5A, 30 m, 0.53 mm ID, 50 µm (cat.# 19723)
Sample Permanent gases
Injection Sample valve
Sample Loop Vol.: 5 µL
Valve Name: 6-port Valco® valve
Inj. Temp.: 200 °C
Valve Temp.: Ambient
Oven
Oven Temp: 27 °C (hold 5 min) to 100 °C at 10 °C/min (hold 5 min)
Carrier Gas He, constant flow
Flow Rate: 5.0 mL/min
Detector Valco® helium ionization detector @ 150 °C



Porous Polymer Columns

Porous polymers are unique, highly retentive stationary phases with a wide application range that are able to elute both polar and nonpolar compounds. They are very hydrophobic, so water has no impact on retention times and even elutes as a good chromatographic peak. The Q-BOND is our most nonpolar and widely used porous polymer column; functional groups can be added to increase polarity (i.e., QS-, S-, and U-BOND). The process used to manufacture porous polymer PLOT columns causes the particles to adhere strongly to the walls of the tubing, so there is virtually no particle generation. You get reproducible performance from column to column, including selectivity and flow.

Rt®-Q-BOND Columns (fused silica PLOT)

100% divinylbenzene

- Nonpolar PLOT column incorporating 100% divinylbenzene.
- Excellent for analysis of C1 to C3 isomers and alkanes up to C12.
- High retention for CO₂ simplifies gas analysis; CO₂ and methane separated from O₂/N₂/CO (Note: O₂/N₂/CO not separated at room temperature).
- Use for analysis of oxygenated compounds and solvents.
- Maximum temperature of 300 °C.
- Also available on metal MXT® tubing! (See page 12.)


ID	df	temp. limits	15-Meter cat.#	30-Meter cat.#
0.25 mm	8 µm	to 280/300 °C	19764	19765
0.32 mm	10 µm	to 280/300 °C	19743	19744
0.53 mm	20 µm	to 280/300 °C	19741	19742

Our porous polymer PLOT columns are not moisture sensitive, making them ideal for applications where moisture is of major concern.

similar phases

HP PLOT Q, CP-PoraPLOT Q, CP-PoraBond Q, Supel-Q-PLOT, AT-Q


Permanent Gases Rt®-Q-BOND (PLOT)



Column Rt®-Q-BOND, 30 m, 0.32 mm ID, 10 µm (cat.# 19744)
Sample 2-5 mole %
Injection
Inj. Vol.: 30 µL split (split ratio 40:1)
Inj. Temp.: 30 °C
Oven
Oven Temp: 30 °C
Carrier Gas H₂, constant flow
Linear Velocity: 38 cm/sec.
Detector HP-µTCD @ 200 °C

Restek® porous polymer PLOT columns cover a wide range of polarities

least polar



Q-BOND

QS-BOND

S-BOND

U-BOND

most polar



Rt®-QS-BOND Columns (fused silica PLOT)

porous divinylbenzene homopolymer


- Intermediate polarity PLOT column incorporating low 4-vinylpyridine.
- Separates ethane, ethylene, and acetylene to baseline.
- Designed for the best possible separation between all C2 isomers.

similar phases

GS-Q

ID	df	temp. limits	15-Meter cat.#	30-Meter cat.#
0.25 mm	8 µm	to 250 °C	19767	19768
0.32 mm	10 µm	to 250 °C	19739	19740
0.53 mm	20 µm	to 250 °C	19737	19738

Refinery Gas Mixture on Rt®-QS-BOND



Peaks	8. Propane	16. <i>trans</i> -2-Butene
1. Air	9. Propadiene	17. Isopentane
2. Methane	10. Isobutane	18. 1-Pentene
3. Carbon dioxide	11. Isobutylene	19. 2-Methyl-2-butene
4. Ethylene	12. 1-Butene	20. <i>n</i> -Pentane
5. Acetylene	13. 1,3-Butadiene	21. <i>cis</i> -2-Pentene
6. Ethane	14. <i>n</i> -Butane	22. <i>n</i> -Hexane
7. Propylene	15. <i>cis</i> -2-Butene	

Column Rt®-QS-BOND, 30 m, 0.53 mm ID, 20 µm (cat.# 19738)
Sample Refinery gas standard
Injection
Inj. Vol.: 20 µL split
Liner: 2 mm (cat.# 20712)
Inj. Temp.: 200 °C
Split Vent
Flow Rate: 35 mL/min
Oven
Oven Temp: 40 °C (hold 2 min) to 225 °C at 15 °C/min (hold 5 min)
Carrier Gas He, constant pressure (11.5 psi, 79.3 kPa)
Linear Velocity: 68 cm/sec @ 40 °C
Detector TCD @ 225 °C
Make-up Gas
Type: He
Data Rate: 20 Hz
Sensitivity Mode: He/H₂
Instrument HP5890 GC

Rt®-S-BOND Columns (fused silica PLOT)

porous divinylbenzene homopolymer


- Midpolarity PLOT column, incorporating high 4-vinylpyridine.
- Use for the analysis of nonpolar and polar compounds.
- Also available on metal MXT® tubing! (See page 12.)

similar phases

CP-PoraPlot S

ID	df	temp. limits	15-Meter cat.#	30-Meter cat.#
0.25 mm	8 µm	to 250 °C	19769	19770
0.32 mm	10 µm	to 250 °C	19747	19748
0.53 mm	20 µm	to 250 °C	19745	19746

Solvent Mixture on Rt®-S-BOND



Peaks	15. Benzene
1. Methanol	16. 1,2-Dimethoxyethane
2. Ethanol	17. Trichloroethylene
3. Acetonitrile	18. 1,4-Dioxane
4. Acetone	19. Pyridine
5. Dichloromethane	20. Dimethylformamide
6. 1,1-Dichloroethene	21. Methylcyclohexane
7. Nitromethane	22. Toluene
8. <i>trans</i> -1,2-Dichloroethylene	23. 2-Hexanone
9. <i>cis</i> -1,2-Dichloroethylene	24. Chlorobenzene
10. Tetrahydrofuran	25. Ethylbenzene
11. Chloroform	26. <i>m</i> -Xylene
12. Ethyl acetate	27. <i>p</i> -Xylene
13. 1,2-Dichloroethane	28. <i>o</i> -Xylene
14. 1,1,1-Trichloroethane	

Column Rt®-S-BOND, 30 m, 0.53 mm ID, 20 µm (cat.# 19746)
Sample Solvent mixture
Injection
Inj. Vol.: 1.0 µL split
Liner: Taper (4 mm) (cat.# 20798)
Inj. Temp.: 200 °C
Split Vent
Flow Rate: 100 mL/min
Oven
Oven Temp: 120 °C to 220 °C at 5 °C/min (hold 5.0 min)
Carrier Gas He, constant pressure (4.2 psi, 29.0 kPa)
Linear Velocity: 40 cm/sec @ 120 °C
Detector FID @ 220 °C



Rt®-U-BOND Columns (fused silica PLOT)

divinylbenzene ethylene glycol/dimethylacrylate

- Restek's highest polarity porous polymer column.
- Polar PLOT column, incorporating divinylbenzene ethylene glycol/dimethylacrylate.
- Highly inert for the analysis of polar and nonpolar compounds.

similar **phases**

HP-PLOT U, CP-PoraPLOT U, CP-PoraBond U


also **available!**

**Metal MXT®
PLOT Columns**

See page 12 for more information.

ID	df	temp. limits	15-Meter cat.#	30-Meter cat.#
0.25 mm	8 µm	to 190 °C	19771	19772
0.25 mm	12 µm	to 190 °C	19782	—
0.32 mm	10 µm	to 190 °C	19751	19752
0.53 mm	20 µm	to 190 °C	19749	19750

Formaldehyde on Rt®-U-BOND



Peaks

1. Air
2. Carbon dioxide
3. Formaldehyde
4. Water
5. Methanol

Column Rt®-U-BOND, 30 m, 0.53 mm ID, 20 µm (cat.# 19750)
Sample Formaldehyde (manual headspace)
Injection
 Inj. Vol.: 10 µL split (split ratio 10:1)
 Liner: 2 mm split Precision® liner w/wool (cat.# 20823)
 Inj. Temp.: 200 °C
 Split Vent
 Flow Rate: 40 mL/min
Oven
 Oven Temp: 100 °C (hold 1 min) to 150 °C at 25 °C/min (hold 3 min)
Carrier Gas He, constant pressure (7.7 psi, 53.1 kPa)
 Linear Velocity: 39 cm/sec @ 100 °C
Detector TCD @ 200 °C
Make-up Gas
 Type: He
 Data Rate: 20 Hz
 Sensitivity Mode: He/H₂
Instrument HP5890 GC

Rt®-Silica BOND Columns (fused silica PLOT)

- Versatile column ideal for analysis of light hydrocarbons, sulfur gases, halocarbons, and carbon dioxide.
- Individually QC tested with sensitive C4 probes to ensure consistent selectivity.
- Proprietary manufacturing process practically eliminates particle release, reducing downtime due to obstructed FID jets.
- Bonded silica stationary phase minimizes impact of water, resulting in reproducible retention times for water-containing samples.


similar **phases**

GS-GasPro, CP-SilicaPLOT

NEW!

ID	temp. limits	15-Meter cat.#	30-Meter cat.#	60-Meter cat.#
0.32 mm	-80 to 260 °C	19784	19785	19786

Saturated and Unsaturated Hydrocarbons on Rt®-Silica BOND PLOT Column



Peaks

1. Methane
2. Ethane
3. Ethylene
4. Acetylene
5. Propane
6. Cyclopropane
7. Propylene
8. Propadiene
9. n-Butane
10. 1-Butene
11. Methyl acetylene
12. 1,3-Butadiene
13. trans-2-Butene
14. Isobutylene
15. cis-2-Butene
16. Isopentane
17. n-Pentane

Column Rt®-Silica BOND, 30 m, 0.32 mm ID (cat.# 19785)
Sample Custom DCG gas standard
Diluent: Nitrogen
Conc.: 1 mole percent
Injection
 Inj. Vol.: 15 µL split (split ratio 35:1)
 Liner: 2 mm straight Sky® inlet liner (cat.# 23313.1)
 Inj. Temp.: 250 °C
Oven
 Oven Temp.: 60 °C (hold 2 min) to 175 °C at 2 °C/min
Carrier Gas He, constant flow
Flow Rate: 3.3 mL/min
Detector FID @ 260 °C
Make-up Gas
 Flow Rate: 50 mL/min
 Make-up Gas Type: N₂
Hydrogen flow: 40 mL/min
Air flow: 400 mL/min
Data Rate: 10 Hz
Instrument Agilent 7890A GC





7" diameter,
11-pin cage

3.5" coil

Metal MXT® PLOT Columns

Restek chemists have developed technology that allows many of our popular PLOT columns to be made on Siltek®-treated stainless steel. These columns have the same characteristics and performance as fused silica PLOT columns (pages 4–11), but offer several additional benefits for process GCs and field applications.

Advantages of metal MXT® PLOT columns include:

- Can be made in small coil diameters—perfect for tight spaces.
- Rugged material withstands rough handling and shock.
- Designed for robust performance in process GCs and field instruments.
- Available in 3.5"-coil diameter or 7"-diameter, 11-pin cage.

ID	df	temp. limits	3.5" coil 15-Meter cat.#	7" diameter 11-pin cage 15-Meter cat.#	3.5" coil 30-Meter cat.#	7" diameter 11-pin cage 30-Meter cat.#
MXT-Msieve 5A						
0.25 mm	20 µm	to 300 °C	79717-273	79717	—	—
0.53 mm	50 µm	to 300 °C	—	—	79723-273	79723
MXT-Alumina BOND/Na₂SO₄						
0.53 mm	10 µm	to 200 °C	—	—	79714-273	79714
MXT-Alumina BOND/MAPD						
0.53 mm	10 µm	to 250 °C	—	—	79728-273	79728
MXT-Q-BOND						
0.25 mm	8 µm	to 300 °C	79718-273	79718	—	—
0.53 mm	20 µm	to 280/300 °C	—	—	79716-273	79716
MXT-S-BOND						
0.53 mm	20 µm	to 250 °C	—	—	79712-273	79712

MXT® Low Dead Volume Connector Kits for Metal Columns

- Connect a guard column/transfer line to an MXT® stainless steel column.
- Low thermal mass tracks rapid oven temperature programming.
- Stainless steel ferrules and nuts.
- Available in “Y” and union configurations.
- Siltek® treatment ensures ultimate inertness.



Each kit contains the MXT® union, two 1/32-inch ferrules and nuts.

Description	qty.	cat.#
For 0.28 mm ID MXT Columns	kit	20397
For 0.32 mm ID MXT Columns	kit	20536
For 0.53 mm ID MXT Columns	kit	20394

PATENTS & TRADEMARKS

Restek® patents and trademarks are the property of Restek Corporation. (See www.restek.com/Patents-Trademarks for full list.) Other trademarks appearing in Restek® literature or on its website are the property of their respective owners. The Restek® registered trademarks used here are registered in the United States and may also be registered in other countries.