

the

# RESTEK


## Advantage

Innovators of High Resolution Chromatography Products

## Optimized Analysis of Brominated Flame Retardants Using an Rtx<sup>®</sup>-500 GC Capillary Column

by Frank Dorman, Ph.D., Innovations Team, Director of Technical Development

- ✓ Elutes decabromodiphenyl ether in 30 minutes.
- ✓ Low bleed for sensitive ECD and MS analyses.
- ✓ Separates other higher molecular weight compounds.



Brominated flame retardants are an emerging environmental concern that present a unique challenge to analysts. One of the most heavily used types of brominated flame retardants are the polybrominated diphenyl ethers (PBDEs). These compounds have the structure shown in Figure 1, with 1 to 10 bromines substituted on the two rings. This makes these compounds fairly heavy (up to approximately 1000amu), thus placing difficult requirements on the gas chromatographic (GC) analytical system. While the analysis may be performed using either electron capture detection (ECD) or mass spectrometric detection (MS), the compounds require a high oven temperature to elute in a reasonable amount of time. This requires an analytical GC column featuring high-temperature fused silica or metal tubing and a high-temperature stationary phase that has both low bleed and the selectivity necessary to separate the PBDE congeners.


Restek chemists, working in conjunction with Karen MacPherson and Eric Reiner at the Ontario Ministry of the Environment, have developed a new fused silica capillary GC column and analytical procedure for separating PBDE congeners in a reasonable

time. The new Rtx<sup>®</sup>-500 column incorporates a carborane-stabilized polydimethylsiloxane polymer in special high-temperature fused silica tubing. The column can be heated to 380°C, and exhibits very low bleed at this extreme temperature. The column combines the stability required for separating higher molecular weight compounds with the sensitivity required for ECD or MS analysis.

PBDE congeners up to decabromodiphenyl ether (PBDE 209) are separated in less than 45 minutes

**Figure 1**

PBDEs can have high molecular weight, due to bromine substitution around the two rings.



using the new Rtx<sup>®</sup>-500 column in a GC-high resolution MS analysis (Figure 2). In a GC-ECD analysis of the same PBDE congeners, decabromodiphenyl ether is eluted in 30 minutes. This same system has been used for baseline separation of toxic PCB congeners and is currently under investigation for analyses of the brominated and chlorinated dioxins and furans. If you must analyze PBDEs or other high molecular weight compounds, we highly recommend the Rtx<sup>®</sup>-500 column for fast separations and reliable quantitation.

(contd. on page 2)



in this issue

**Rtx<sup>®</sup>-500 GC Capillary Column**  
...pg. 1-2  
**Vespel<sup>®</sup> Ring Inlet Seal**  
...pg. 3  
**Leak Detective<sup>™</sup> II Electronic Leak Detector**  
...pg. 3  
**Bringing Packed Column Technology Back to Life**  
...pg. 4

**Rt-Msie<sup>™</sup> 5A & MXT<sup>®</sup>-5A PLOT Columns**

...pg. 5

**HPLC Analysis of Narcotic/Acetaminophen Mixtures**

...pg. 6-7

**More Reliable Results from Semivolatiles Analysis**

...pg. 8-9

**New Analytical Reference Materials**

...pg. 9

**Scotty<sup>®</sup> Transportables**

...pg. 10

**European Pharmacopoeia Reference Materials**

...pg. 11

**Polar Analytes by HPLC/Replacement Instrument Parts**

...pg. 12-13

**Super-Clean<sup>™</sup> Gas-Trapping System for LC/MS**

...pg. 14

**New Reference Books**

...pg. 14

**Cool Tools**

...pg. 15

**Behind the Scenes**

...pg. 16

**RESTEK**

[www.restekcorp.com](http://www.restekcorp.com)

2002 vol. 3

IONAL

**CHROMALYTIC** +61(0)3 9762 2034  
**ECHnology** Pty Ltd


Website NEW : [www.chromalytic.net.au](http://www.chromalytic.net.au) E-mail : [info@chromtech.net.au](mailto:info@chromtech.net.au) Tel: 03 9762 2034 . . . in AUSTRALIA

**Australian Distributors**  
Importers & Manufacturers  
[www.chromtech.net.au](http://www.chromtech.net.au)



**Figure 2**

The Rtx®-500 column completely separates PBDE congeners in less than 45 minutes in a GC/MS analysis.



**NB DecaBDE (last elutor)**  
elutes at ~43 min.


Chromatogram courtesy of Ontario  
Ministry of the Environment

Reference materials courtesy of  
Wellington Laboratories, Guelph,  
Ontario, Canada  
www.well-labs.com  
US Distributor:  
TerraChem Inc. U.S.A.  
8600 Shawnee Mission Pkwy.,  
Suite 305  
Shawnee Mission, KS 66202  
Phone: 913-722-4919  
Toll-free: 877-809-7039  
Fax: 913-722-4669  
Website: <http://www.terrachem.com>  
Email: [info@terrachem.com](mailto:info@terrachem.com)

**Column:** Rtx®-500 30m, 0.25mm, 0.15µm (cat.# 10750); **GC:** Agilent 6890+; **Sample:** Wellington Laboratories BDE Mix C (300pg/µL each in nonane); **Inj.:** 1µL splitless injection, drilled Uniliner® (cat.# 21054); **Oven temp.:** 100°C (hold 0.64 min.) to 110°C @ 10°C/min. (hold 0 min.) to 180°C @ 80°C/min. (hold 23 min.) to 350°C @ 5°C/min.; **Flow rate:** constant @ 1.5mL/min; **Injector temp.:** 300°C; **Instrument configuration:** Micromass Autospec-UltimaNT (High Resolution Mass Spectrometer); **Source Temperature:** 300°C

**Figure 3**

When used with less expensive ECD methodology, an Rtx®-500 column resolves PBDEs in 30 minutes.



### Ordering Information | Rtx®-500 Columns (Fused Silica)

(Crossbond® carborane/dimethyl polysiloxane) Stable to 380°C

ID	df (µm)	temp. limits	30-Meter	60-Meter
0.25mm	0.15	-60°C to 380°C	10750	10751
0.53mm	0.15	-60°C to 380°C	10752	

[www.restekcorp.com](http://www.restekcorp.com)

**HROMalytic** +61(0)3 9762 2034  
**ECHnology** Pty Ltd  
Website NEW : [www.chromalytic.net.au](http://www.chromalytic.net.au) E-mail : [info@chromtech.net.au](mailto:info@chromtech.net.au) Tel: 03 9762 2034 . . . in AUSTRALIA

**for more info**  
request Applications Note 59389A.

**Australian Distributors**  
Importers & Manufacturers  
[www.chromtech.net.au](http://www.chromtech.net.au)



# Vespel® Ring Inlet Seal

## Seals the First Time, Every Time

by Donna Lidgett, GC Accessories Product Marketing Manager

- ✓ Easy-to-use, patent-pending design makes a better seal, easily.
- ✓ Prevents oxygen from damaging your columns.
- ✓ Reduces wear on the injection port body.


In Agilent split/splitless injection ports, it can be difficult to make and maintain a good seal with a conventional metal inlet disk. The metal-to-metal seal dictates that the analyst apply considerable torque to the reducing nut, and, based on our testing, this does not ensure a leak-tight seal. Over the course of oven temperature cycling, metal seals are prone to leaks, which ultimately can degrade the capillary column, and cause other analytical difficulties.

Our Vespel® Ring Inlet Seal greatly improves injection port performance—it seals even after repeated temperature cycles and without retightening the reducing nut! This seal features a Vespel® ring


embedded into its face. This soft Vespel® ring will not harm the critical seal on the injector body, and is outside the sample flow path. Tests using a high sensitivity helium leak detector indicate the Vespel® Ring Inlet Seal seals equally effectively at torques of 5lb. or 60lb. (Figure 1).

Why trust a metal-to-metal seal when you can make leak-tight seals quickly and easily—and more reliably—with the Restek Vespel® Ring Inlet Seal? Use the stainless steel seal for analysis of unreactive compounds. To reduce breakdown and adsorption of active compounds, use the gold-plated or

**Figure 1**  
The Vespel® Ring Inlet Seal achieves leak-tight seals even at low torque, reducing the chance of leaks.



Silcosteel®-treated seals. The gold surface offers better inertness than standard stainless steel; Silcosteel® treatment provides inertness similar to that of fused silica capillary columns.



### Ordering Information | Vespel® Ring Inlet Seals for Agilent 5890/6890 and 6850 GCs

0.8mm ID Vespel Ring Inlet Seal (washers included)		2-pk.	10-pk.
Gold-Plated		21562	21563
Silcosteel®		21564	21565
Stainless Steel		21560	21561
1.2mm ID Vespel Ring Inlet Seal (washers included)		2-pk.	10-pk.
Gold-Plated		21568	21569
Silcosteel®		21570	21571
Stainless Steel		21566	21567

## A Compact, Sensitive Leak Detector For Every GC Analyst



### The Restek Leak Detective™ II

by Donna Lidgett, GC Accessories Product Marketing Manager

- ✓ Fast results—responds to trace leaks in less than 2 seconds.
- ✓ Sensitive—detects trace leaks at  $1 \times 10^{-10}$  cc/sec.; as low as 100ppm.
- ✓ Micro-chip design improves sensitivity and response time over previous models.
- ✓ Compact, ergonomic design is easy to hold and operate with one hand.
- ✓ Battery-operated for portability (one 9 volt) instant auto-zeroing.

Gas leaks in your GC system can increase detector noise, cause baseline instability, waste carrier gas, and damage valuable analytical columns. Leak checks should be a regular part of your GC maintenance program. The new Leak Detective™ II electronic leak detector is the affordable solution for detecting gas leaks. It will identify minute gas leaks that might go undetected by liquid leak detectors.\*

The Leak Detective™ II electronic leak detector incorporates micro-chip technology and a new

design, to give you better sensitivity and faster response time in a more compact unit. An auto-zero feature allows you to instantly zero the leak detector with a push of a button, and the ergonomic design brings all the controls to your fingertips for easy use. The unit responds in less than two seconds to trace leaks of gases with thermal conductivities different than air. Leaks are indicated by an audible alarm, as well as by an LED readout. For easy, sensitive, and reliable leak detection, order a new Leak Detective™ II electronic leak detector today.



### Ordering Information | Leak Detective II™

Description	qty.	cat.#
Leak Detective™ II Leak Detector (9 volt, Battery-Operated)	ea.	20413

\*Never use liquid leak detectors on a capillary system because liquids can be drawn into the column.

**Caution:** NOT designed for determining leaks of combustible gases. A combustible gas detector should be used for determining combustible gas leaks in possibly hazardous conditions.



# Packed Column Technology Has New Life

## With Restek Innovations

by Neil Mosesman, GC Columns Product Marketing Manager

- ✓ Optimized tubing and support ensure unsurpassed inertness for trace analyses.
- ✓ Bonded phases shorten conditioning times, greatly reduce bleed, and prolong column lifetimes.

For years column manufacturers have proclaimed that packed GC columns have gone the way of the dinosaurs. While packed columns do not fill every need, in some situations their use is dictated by methodology, and, for some analysts, they are the column of preference. Restek innovations have brought new life to the "expired" packed column in all three aspects of column technology: the support, the phase, and the tubing. Silcoport™ support provides unsurpassed inertness for trace analyses. Our bonded phase packings are a revolution in packed column technology. They significantly shorten conditioning times, greatly reduce column bleed, and prolong column lifetimes. SilcoSmooth™ tubing combines the inertness of glass with the durability of stainless steel.


The great sensitivity of modern detection systems and a progressing need to reduce detection limits place challenging demands on a chromatography column. Silcoport™ diatomaceous earth support is the modern solid support that we developed to meet these demands. Unlike conventional DMDCS

deactivation, we use a proprietary fused silica deactivation technology and a special mixture of deactinants to ensure the greatest inertness (Figure 1) without changing the polarity of the stationary phase. Each batch of support is carefully tested to confirm a uniform particle size distribution that ensures columns with maximum efficiency.

By applying our experience in stationary phase synthesis in conjunction with our unique Silcoport™ packing deactivation process we create completely bonded packing materials. To encompass a wide range of applications, we offer Rtx®-1 and Rtx®-5, bonded methyl silicone phases, and Stabilwax®, a bonded Carbowax® phase. Each phase is completely cross-linked on Silcoport™ support. In side-by-side comparisons with conventional nonbonded methyl silicone phase columns, Rtx®-1 and Rtx®-5 columns have lower bleed, improved peak shape, and longer useful lives (Figure 2). Evaluations with an Rtx®-1 column show retention times are repeatable after only 30 minutes of conditioning.

**Figure 1**

A trace level analysis of pesticides, including labile endrin and DDT, demonstrates the inertness of Silcoport™ support.



### Ordering Information | Bonded Packed Column Stationary Phases:

Liquid Phases Bonded on 100/120 Silcoport™ P	Stainless Steel Tubing				SilcoSmooth™ Tubing			
	L (ft.)	OD (in.)	ID (mm)	cat.#*	L (m)	OD (in.)	ID (mm)	cat.#*
3% Rtx®-1	6	1/8	2.1	80441-	2	1/8	2	80401-
10% Rtx®-1	6	1/8	2.1	80442-	2	1/8	2	80405-
20% Rtx®-1	6	1/8	2.1	80443-	2	1/8	2	80409-
3% Rtx®-5	6	1/8	2.1	80444-	2	1/8	2	80477-
10% Rtx®-5	6	1/8	2.1	80445-	2	1/8	2	80478-
20% Rtx®-5	6	1/8	2.1	80446-	2	1/8	2	80479-
5% Rtx®-Stabilwax®	6	1/8	2.1	80447-	2	1/8	2	80415-
10% Rtx®-Stabilwax®	6	1/8	2.1	80448-	2	1/8	2	80416-
20% Rtx®-Stabilwax®	6	1/8	2.1	80449-	2	1/8	2	80417-
Rtx®-1 SimDist 2887	25*	1/8	2.1	80450	25*	1/8	2	80000



\*Please include configuration suffix number (refer to our catalog, lit. cat.# 59662).

If your analysis involves reactive compounds, you probably have used fragile, inflexible glass columns, but now you can do better. Made from ultra-smooth, seamless 304 stainless steel, and treated with our innovative SilcoSteel™ deactivation process, SilcoSmooth™ tubing combines the inertness of glass with the strength and flexibility of stainless steel. SilcoSmooth™ tubing can replace glass tubing in virtually any application. For analyses of ppb levels of sulfur-containing compounds, use Sulfinert™ tubing packed with Rt-XLSulfur™ packing. For undemanding applications we can make columns from conventional tubing: stainless steel, Hastelloy®, nickel, copper, or Teflon®.

In combination, Silcoport™ support, our bonded phase packings, and SilcoSmooth™ tubing make packed column GC a viable alternative in many applications in which the technique had been endangered. If you use packed columns, and think that you have to live with limitations, call us. We can provide the column that will give you the separation you need, but with convenience, inertness, and column lifetimes you never expected from a packed column.

**Figure 2**

Bonded packed columns exhibit longer lifetimes than nonbonded columns.



www.restekcorp.com



# Rt-Msieve™ 5A & MXT®-5A

## PLOT Columns

### Superior Analyses of Permanent Gases

by Neil Mosesman, GC Columns Product Marketing Manager

- ✓ Fast, efficient separations at above ambient temperatures.
- ✓ 100% bonding process eliminates the need for particle traps.
- ✓ Stainless steel columns for durability.

Gas-liquid chromatography (GLC), the most common mode of gas chromatography, has limited application in analyses of gases. Subambient temperatures often are required to achieve a separation, and cryogenic cooling systems are costly and inconvenient. Gas-solid chromatography (GSC), in which gaseous analytes are absorbed into the packing particles, rather than into a surface coating, is far more effective for separating gases. Difficult-to-separate small molecules, such as argon and oxygen, butene isomers, and many others, can be separated by GSC at above ambient temperatures.

Just as capillary columns offer important advantages over packed GLC columns, porous layer open tubular columns-PLOT columns-offer significant advantages over packed GSC columns. Their open design gives PLOT columns greater permeability, and their narrow diameter ensures sharper peaks. The open construction also affords a smaller pressure drop per unit length, so longer columns can be used. This means much higher column efficiency and, therefore, superior resolution. In brief, PLOT columns provide faster and more sensitive analyses than packed GSC columns.

Restek PLOT columns are especially effective for separating mixtures of gaseous analytes. Rt-Msieve™ 5A and MXT®-Msieve 5A PLOT columns contain molecular sieve 5A particles that are bonded to the inner surface of the tubing, using a proprietary process that prevents particle dislocation that could damage valves and detection systems. They are designed for fast, efficient separation of argon and oxygen, hydrogen and helium, and other permanent gases, including permanent gases admixed in refinery or natural gas. Special coating and deactivation procedures ensure chromatographic efficiency and the integrity of the porous layer bonding. Finely

controlled pore size allows selective adsorption of specific target compounds, ensuring difficult separations can be made without subambient temperatures. Figure 1 shows an Rt-Msieve™ 5A column can separate oxygen from argon to baseline, at above ambient temperature, in approximately 2 min. Figure 2 shows the permanent gases resolved from methane in 4 minutes.\* Stainless steel MXT®-Msieve 5A PLOT columns offer the same powerful separating capabilities as fused silica Rt-Msieve™ 5A PLOT columns, plus high resistance to physical damage and ability to be coiled to diameters as small as 3.5" (<9cm), making MXT® columns ideal for portable GCs, process control applications, and other demanding situations.

In addition to Rt-Msieve™ 5A and MXT®-Msieve 5A columns, we manufacture PLOT columns for a wide range of other applications. Rt-Alumina™ PLOT columns (Al<sub>2</sub>O<sub>3</sub> solid phase) offer fast, reproducible

performance for determining hydrocarbon purity or monitoring hydrocarbon streams. Porous polymer Rt-QPLOT™ and MXT®-QPLOT columns (nonpolar), Rt-SPLIT™ and MXT®-SPLIT columns (intermediate polarity), and Rt-UPLOT™ and MXT®-UPLOT columns (highly polar) are particularly useful for situations in which water is likely to be encountered. Applications for these columns include permanent gases at subambient temperatures, carbon dioxide and other inorganic gases, hydrocarbon mixtures, and many nonpolar, intermediate polarity, and polar solvents. For more information and example analyses on Restek PLOT columns, refer to our current chromatography products catalog or our website, or request our new PLOT column flyer (lit. cat. #59456).



If your analyses call for difficult separations of gaseous analytes, and neither conventional packed GC columns nor WCOT capillary columns are providing the separations you want, or your analyses depend on costly or time-consuming conditions, a Restek PLOT column can make your work simpler.

#### Ordering Information | Rt-Msieve™ 5A (Fused Silica PLOT) Temp. limit to 300°C

ID	df (µm)	15-Meter	30-Meter
0.32mm	30	19720	19722
0.53mm	50	19721	19723

#### Ordering Information | MXT®-Msieve 5A (Metal PLOT) Temp. limit to 300°C

ID	df (µm)	15-Meter	30-Meter
0.53mm	50	79721	79723



\*Carbon dioxide is difficult to elute from a molecular sieve column, but is isolated easily on an Rt-QPLOT™ porous polymer column at 30°C. Request lit. cat.# 59540 for details.

#### Restek PLOT columns are superior for 5 reasons:

1. Most efficient and consistent analyses.
2. No need for particle traps.
3. Reproducible quality at affordable prices.
4. Most effective phase for your separation: alumina, molecular sieve 5A, or porous polymer.
5. Fused silica columns for most applications, metal columns for exceptional durability.

• 5 •



# HPLC Analysis of Narcotic/ Acetaminophen Admixtures

## What to Do If a Compendium Method Doesn't Work

by Vernon Bartlett, HPLC Innovations Manager

- ✓ Make changes or modifications stepwise, with defined purpose in mind.
- ✓ When possible, create and validate a single method for a range of similar analytes.

Sometimes methods described in the United States Pharmacopoeia (USP), the European Pharmacopoeia (EP), the British Pharmacopoeia (BP), or other compendia do not provide the desired robustness in separation or reproducibility, or results barely pass system suitability requirements. Modifications can be made to improve the methodology, and the results compared statistically to the original. To improve analysis efficiency and reduce costs associated with revalidating and testing, it may be desirable to create and validate a sin-



gle analytical method for a range of similar drug products.

Many narcotics are very similar in structure, often varying by only a single substitution. Morphine, codeine, hydrocodone, and oxycodone are quite similar, for example (Figure 1). Some of these closely related compounds—all but morphine, in fact—might be blended with other analgesics, such as acetaminophen (APAP). USP 25 describes more than 7 different methods to test these raw materials

and admixtures; some of the older methods do not use HPLC as a primary test for purity.

One of the chromatographic applications in USP 25 is for the analysis of oxycodone raw material. After reading the mobile phase section, we saw some potential problems with the method, including:

- 1) The use of methanol in this analysis could lead to high background absorption and loss of linear range, because the analytical wavelength is 206nm, and the UV cutoff for methanol is 235nm. In extreme cases this also can reduce sensitivity—the more energy the background absorbs, the less is available to the analyte.
- 2) An ion-pairing agent (hexane sulfonic acid) is introduced into the mobile phase without a buffer to maintain pH. This could lead to widened peaks, tailing peaks, and retention time drift.
- 3) Triethylamine (TEA) modifier is included in the method. When basic compounds are analyzed on older-type HPLC columns, TEA often is added as competing base, to reduce the tailing caused by acidic silanol activity. If the analytical species are neutral, or have been “neutralized” by an ion-pairing agent, TEA should have no beneficial effect. Adding TEA, a base, to a mobile phase containing sulfonic acids will cause acid/base neutralization, producing a salt and water and reducing the effective concentration of the acidic ion-pairing agent. This could lead to the formation of undesirable side products in the mobile phase that also will absorb in the low UV range, creating noisy baselines. Furthermore, TEA is volatile, and its composition might change over time if the mobile phase is sparged.



Thus, some aspects of the method appear redundant and some might actually compromise the separation. In addition, some of the reagents, such as TEA, might not be necessary for modern columns. After performing the USP 25 method as written, we made some tests to determine actual needs to achieve the system suitability requirements as specified.

With peak shape, separation, and proper analytical technique in mind, we attempted to eliminate some of the perceived problems. We realized that by using 284nm as the detection wavelength, rather than 206nm as used in USP 25, we might not see some impurities, but in real life the material should be tested against some known source for potency. (Note that with the additional reagents removed, both Ultra C8 and Pinnacle II™ C8 columns provided good results at the 206nm wavelength.)

Next we removed the ion pairing agent and the TEA. We elected to keep a 20 mM phosphate buffer system to maintain a pH of 2.5. Then we reduced the temperature from 35°C to 27°C, to determine whether the greater mass transfer and analyte solubility in the mobile phase at 35°C had been masking other potential problems.



These changes led to a slight increase in tailing for all compounds on both Ultra C8 and Pinnacle II™ C8 columns, but this was acceptable, especially because the run time for the analysis was reduced by a factor of 3 and resolution was improved by 59% to 79%. The system passed the system suitability requirements in the USP monograph.

In the next experiment, we re-introduced the ion pair reagent hexane sulfonic acid into the system under the control of the pH 2.5 phosphate buffer system. The run time doubled, relative to the original procedure, demonstrating that TEA did affect the concentration of the ion-pairing agent. Reducing the concentration of ion pairing agent, or

using a shorter chain length ion-pairing agent, might have been a better alternative to adding TEA. The system still passed the system suitability requirements listed by the USP, but the chromatogram was much noisier—and equilibration problems seen in the USP 25 analysis returned.


After reviewing the monographs for admixtures containing structurally related narcotics and acetaminophen, we created a single separation for morphine sulfate, acetaminophen, codeine phosphate, oxycodone HCl, and hydrocodone bitartrate. The goal was to create an adequate separation while keeping the method as simple as possible. We chose an Ultra C18 column and set detection to 235nm. All components, including a small unknown peak, were separated to baseline (Figure 2).

Next, we increased the amount of buffer to 90% (a 5% increase). This simple increase doubled the analysis time. Resolution doubled between most components, with the greatest change between acetaminophen and codeine. The unknown peak disappeared and probably co-eluted with morphine.

We adjusted the mobile phase ratio to 85:15, buffer:organic solvent, using a 90:10 mixture of acetonitrile and methanol as the organic solvent. Resolution improved, relative to the original mobile phase composition, analysis again was under 10 minutes, and the unknown peak returned (Figure 3). For this analysis, these conditions provided the most desirable results.

In summary, the goal of any method should be to achieve the most stable and robust separation. Sometimes methods are more complicated than they need to be, and this can make analysis unnecessarily difficult. Even troubleshooting such methods adds to production costs. When preparing to follow a method always attempt to determine the reason a reagent would be included in a mobile phase. Any change or modification should have an established scientific purpose. By creating more universal methods for analyses of structurally related compounds, it should be possible to reduce costs for supplies, increase laboratory analysis efficiency, and reduce personnel training time.

For chromatograms illustrating the changes in separation that occur with each change in the mobile phase, please request Applications Note #59453. If you encounter problems when analyzing your samples according to an established method, our experienced Technical Service chemists will be glad to help. Contact them at 800-356-1688, ext. 4 or 814-353-1300, ext. 4, or contact your Restek representative.



#### Ordering Information | Ultra C18 5µm Columns

Length	1.0mm ID cat.#	2.1mm ID cat.#	3.2mm ID cat.#	4.0mm ID cat.#	4.6mm ID cat.#
30mm	9174531	9174532	9174533	—	9174535
50mm	9174551	9174552	9174553	—	9174555
100mm	9174511	9174512	9174513	9174514	9174515
150mm	9174561	9174562	9174563	9174564	9174565
200mm	9174521	9174522	9174523	—	9174525
250mm	9174571	9174572	9174573	—	9174575

#### Trident™ Direct HPLC Guard Column System

Choose from three levels of protection!



##### Trident™ Direct high-pressure filter

✓ Protection against particulate matter.



##### Trident™ Direct 1cm guard cartridge holder with filter

✓ Protection against particulate matter.  
✓ Moderate protection against irreversibly adsorbed compounds.



##### Trident™ Direct 2cm guard cartridge holder with filter

✓ Protection against particulate matter.  
✓ Maximum protection against irreversibly adsorbed compounds.

formoreinfo

For information about Trident™ guard columns, request the Trident™ Fast Facts (lit. cat.# 59314 and 59896).



# More Reliable Results From Semivolatiles Analysis

## Using Restek Columns and Standards

by Christopher English, Environmental Innovations Chemist

- ✓ Rtx®-5Sil MS columns resolve critical pairs and minimize bleed.
- ✓ Integral guard column available.
- ✓ 8270MegaMix™ minimizes mixtures needed, has maximum stability.
- ✓ Monitor all relevant semivolatiles at one detector sensitivity—8270 MegaMix™ includes 3 and 4 methylphenol at 0.5x concentrations of other components.

Complex mixtures of semivolatile organic compounds are extracted from water, soil, or solid waste samples, concentrated, and analyzed by gas chromatography. The current compound list for US EPA Method 8270D, for example, includes basic, neutral, and acidic compounds with boiling points from 150°C to 500°C. Other semivolatiles methods are similarly complex. Because these analyses encompass a broad range of compound classes and require low detection limits, and because sample extracts can include non-target contaminants, significant demand is placed on the efficiency, inertness, thermal stability, and sample capacity of the analytical column. These parameters must be optimized to provide good resolution, fast analysis times, and high sample throughput. The column must have adequate sample capacity to handle the high concentrations of contaminants sometimes found in these extracts, while exhibiting the high inertness needed for accurate quantification of target analytes down to low ng/μL levels.


Restek has designed Rtx®-5Sil MS capillary columns to address the demands of semivolatile by GC/MS. Silarylene polymer technology stiffens the siloxane chain, preventing its thermal breakdown (column bleed). The content of this aryl functionality has been adjusted to give excellent efficiency and lower bleed, compared to conventional 5% diphenyl/95%dimethyl phases; Rtx®-5Sil MS columns exhibit excellent inertness and low bleed, even at 330°C. The optimized stationary phase, proprietary deactivations, and inherently low bleed of the Rtx®-5Sil MS phase, combined with the integral guard column, overcome the problems presented by the compounds and conditions inherent to semivolatiles analysis. High column efficiency ensures the resolution needed to quantify critical pairs and structural isomers, as shown by the separation of benzo(b)- and benzo(k)fluoranthene (peaks 83/84) in Figure 1.

Rtx®-5Sil MS columns are available with an integral, deactivated 5- or 10-meter Integra-Guard™ guard column that prevents non-volatile residues from collecting in the analytical column, where they could interfere with the analytes. Made from a continuous length of tubing, innovative Integra-Guard™ columns offer the column-protecting advantages of a guard column without the potential for leaks at the interface.\*

\*For more information about Integra-Guard™ columns, request lit. cat. # 59441.

Our 8270 MegaMix™ eliminates mixing and minimizes preparation time for calibration and laboratory control samples—it combines all current target analytes in EPA Method 8270D. 8270 MegaMix™ components are indicated in bold in the list of analytes in Figure 1. A unique feature of this mix is the inclusion of 3-methyl- and 4-methylphenol at 0.5x the concentration of the other components, so you won't have to adjust reporting limits when analyzing for these compounds. A long shelf life for unopened ampuls of 8270 MegaMix™ minimizes ordering and inventory problems.

If you are monitoring semivolatile analytes according to US EPA Method 8270D, or similar methods, trust Restek Rtx®-5Sil MS columns and 8270 MegaMix™ to help you obtain reliable, consistent results.



[www.restekcorp.com](http://www.restekcorp.com)

• 8 •



1. N-nitrosodimethylamine
2. pyridine
3. 2-fluorophenol
4. phenol-d6
5. phenol
6. aniline
7. bis(2-chloroethyl)ether
8. 2-chlorophenol
9. 1,3-dichlorobenzene
10. 1,4-dichlorobenzene-d4
11. 1,4-dichlorobenzene
12. benzyl alcohol
13. 1,2-dichlorobenzene
14. 2-methylphenol
15. bis(2-chloroisopropyl)ether
16. 4-methylphenol/  
3-methylphenol\*
17. N-nitroso-di-n-propylamine
18. hexachloroethane
19. nitrobenzene-d5
20. nitrobenzene
21. isophorone
22. 2-nitrophenol
23. 2,4-dimethylphenol
24. bis(2-chloroethoxy)methane
25. benzoic acid
26. 2,4-dichlorophenol
27. 1,2,4-trichlorobenzene
28. naphthalene-d8
29. naphthalene
30. 4-chloroaniline
31. hexachlorobutadiene
32. 4-chloro-3-methylphenol
33. 2-methylnaphthalene
34. 1-methylnaphthalene
35. hexachlorocyclopentadiene
36. 2,4,6-trichlorophenol
37. 2,4,5-trichlorophenol
38. 2-fluorobiphenyl
39. 2-chloronaphthalene
40. 2-nitroaniline
41. 1,4-dinitrobenzene
42. dimethylphthalate
43. 1,3-dinitrobenzene
44. 2,6-dinitrotoluene
45. acenaphthylene
46. 1,2-dinitrobenzene
47. 3-nitroaniline
48. acenaphthene-d10
49. acenaphthene
50. 2,4-dinitrophenol
51. 4-nitrophenol
52. dibenzofuran
53. 2,4-dinitrotoluene
54. 2,3,4,6-tetrachlorophenol
55. 2,3,5,6-tetrachlorophenol
56. diethyl phthalate
57. 4-chlorophenyl phenyl ether
58. fluorene
59. 4-nitroaniline
60. 4,6-dinitro-2-methylphenol
61. diphenylamine\*\*
62. azobenzene\*\*\*
63. 2,4,6-tribromophenol
64. 4-bromophenyl phenyl ether
65. hexachlorobenzene
66. pentachlorophenol
67. phenanthrene-d10
68. phenanthrene
69. anthracene
70. carbazole
71. di-n-butylphthalate
72. fluoranthene
73. benzidine
74. pyrene
75. p-terphenyl-d14
76. butyl benzyl phthalate
77. bis(2-ethylhexyl)adipate
78. benzo(a)anthracene
79. chrysene-d12
80. chrysene
81. bis(2-ethylhexyl)phthalate
82. di-n-octyl phthalate
83. benzo(b)fluoranthene
84. benzo(k)fluoranthene
85. benzo(a)pyrene
86. perylene-d12
87. indeno(1,2,3-cd)pyrene
88. dibenzo(a,h)anthracene
89. benzo(ghi)perylene

Bold indicates a component of the 8270D MegaMix™

**Column:** Rtx®-5Sil MS w/ 5-meter Integra-Guard™  
30m, 0.25mm ID, 0.25µm (cat.# 12723-124)

**Instrument:** Agilent 5973 GC/MS

**Sample:** US EPA Method 8270D Mix 1µL, 16 ppm each component (16ng on column) 8270 MegaMix™ (cat.# 31686) Benzoic Acid Standard (cat.# 31415) Benzidine Standard (cat.# 31441) Acid Surrogate Mix (4/89 SOW)(cat.# 31063) B/N Surrogate Standard Mix (4/89 SOW) (cat.# 31062) SV Internal Standard Mix (cat.# 31006)

**Solvent:** dichloromethane

**Inj.:** 1.0µL splitless (hold 0.3 min.), 4mm Drilled Uniliner™ (cat.# 21055)

**Inj. temp.:** 300°C

**Carrier gas:** helium, constant flow

**Flow rate:** 1.1mL/min.

**Dead Time:** 1.8 minutes @ 35°C

**Oven temp.:** 35°C (hold 4 min.) to 245°C @ 25°C/min. (no hold) to 330°C @ 6°C/min. (hold 3 min.)

**Det:** GC/MS

**Transfer line temp.:** 280°C

**Scan range:** 35–550 amu

**Solvent Delay:** 5 min.

**Tune:** DFTPP

**Ionization:** EI

\*Each at 0.5x concentration of other components.

\*\*N-nitrosodiphenylamine (8270-listed analyte) decomposes to diphenylamine (mix component).

\*\*\*1,2-diphenylhydrazine (8270-listed analyte) decomposes to azobenzene (mix component).

## Ordering Information | Rtx®-5Sil MS Columns (Fused Silica)

(Equivalent selectivity of Crossbond® 5% diphenyl/95% dimethyl polysiloxane) Stable to 360°C

ID	df (µm)	temp. limits	15-Meter	30-Meter
0.25mm	0.10	-60 to 330/350°C	12705	12708
	0.25	-60 to 330/350°C	12720	12723
	0.50	-60 to 330/350°C	12735	12738
	1.00	-60 to 325/350°C	12750	12753

### 8270 Matrix Spike Mix (76 components)

See **bold** compounds in Figure 1 peak list.

200µg/mL each (except noted) in methanol:methylene chloride:benzene (80:15:5), 5mL/ampul

Each	5-pk.	10-pk.
31687	31687-510	—
	with data pack	
31687-500	31687-520	31787

### 8270 MegaMix™ (76 components)

See **bold** compounds in Figure 1 peak list.

1,000µg/mL each (except noted) in methylene chloride:benzene (75:25), 1mL/ampul

Each	5-pk.	10-pk.
31686	31686-510	—
	with data pack	
31686-500	31686-520	31786

# New Analytical Reference Materials

## FAMEs, Acetates, BTEX, Glycols

### Food Industry FAME Mix (37 components)

• Includes *trans* FAMEs.

Chain	% by Weight
C4:0	4.0
C6:0	4.0
C8:0	4.0
C10:0	4.0
C11:0	2.0
C12:0	4.0
C13:0	2.0
C14:0	4.0
C14:1( <i>cis</i> -9)	2.0
C15:0	2.0
C15:1( <i>cis</i> -10)	2.0
C16:0	6.0
C16:1( <i>cis</i> -9)	2.0
C17:0	2.0
C17:1( <i>cis</i> -10)	2.0
C18:0	4.0
C18:1( <i>trans</i> -9)	2.0
C18:1( <i>cis</i> -9)	4.0
C18:2(all- <i>trans</i> -9,12)	2.0
C18:2(all- <i>cis</i> -9,12)	2.0
C18:3(all- <i>cis</i> 6,9,12)	2.0
C18:3(all- <i>cis</i> 9,12,15)	2.0
C20:0	4.0
C20:1( <i>cis</i> -11)	2.0
C20:2(all- <i>cis</i> 11,14)	2.0
C20:3 (all- <i>cis</i> 8,11,14)	2.0
C20:3(all- <i>cis</i> 11,14,17)	2.0
C20:4(all- <i>cis</i> 5,8,11,14)	2.0
C20:5(all- <i>cis</i> 5,8,11,14,17)	2.0
C21:0	2.0
C22:0	4.0
C22:1( <i>cis</i> 13)	2.0
C22:2(all- <i>cis</i> 13,16)	2.0
22:6(all- <i>cis</i> 4,7,10,13,16,19)	2.0
C23:0	2.0
C24:0	4.0
C24:1( <i>cis</i> -15)	2.0

30mg/mL in methylene chloride, 1mL/ampul ea.

35077

### 8260B Acetate Mix (7 components)

• Includes methyl acetate and n-amyl acetate.

n-amyl acetate	methyl acetate
butyl acetate	propyl acetate
ethyl acetate	vinyl acetate
isopropyl acetate	

2,000µg/mL each in P&T methanol, 1mL/ampul

Each	5-pk.	10-pk.
30489	30489-510	
w/data pack	30489-500	30589

### BTEX Standard (6 components)

• m- and p-xylene at 1/2 concentration.  
• Contact Restek for future formulations.

benzene	m-xylene*
ethylbenzene	o-xylene
toluene	p-xylene*

2,000µg/mL each in P&T methanol, 1mL/ampul

Each	5-pk.	10-pk.
30488	30488-510	
w/data pack	30488-500	30588

\*1,000µg/mL

### Glycols Standard

• Assay for de-icing compounds.

ethylene glycol	propylene glycol
-----------------	------------------

50,000µg/mL each in DI water, 1mL/ampul

Each	5-pk.	10-pk.
30471	30471-510	
w/data pack	30471-500	30571



# Scotty® Transportables

## For Laboratory or Field Application

by Donna Lidgett, Air Monitoring Product Marketing Manager



- ✓ Portability makes your job easier.
- ✓ 4-Liter, 14-liter, and 48-liter sizes.
- ✓ Long shelf life.

Restek now offers a broad selection of Scotty® Transportables, ranging from pure gases to multi-component mixes. These standards have found

many uses in the analytical lab, including purging, calibrating, and peak identification.

We offer three sizes to choose from. The 4-liter container has a delivery tube for purge gas or connection to sample loop or bag. The 14-liter container has a CGA 160 fitting for connection to an analytical system. The 48-liter cylinder has a CGA 165 connection, and can deliver large volumes of sample.



### Regulators for Scott Transportable Gases



#### Specifications

Maximum Inlet Pressure: 300psig  
Outlet Pressure Range: 2–10psig  
Operating Temperature Range: 35°F to 150°F (2°C to 65°C)  
Outlet Connection: 1/4" Female NPT

#### Materials of Construction

Body: Brass  
Diaphragm: Viton®  
Seat: Acetal  
Seal: Viton®

#### Description qty. cat.#

Regulator with CGA 160 inlet connection for 14L Scott container ea. 22690  
Regulator with CGA 165 inlet connection for 48L Scott container ea. 22691

#### Description

##### Pure Gases

Air, zero (THC <1ppm)  
Argon, 99.995%  
Carbon dioxide, 99.80%  
Hydrogen, 99.99%  
Methane, 99.00%  
Oxygen, 99.60%

##### Two-Component Mixtures

Benzene in air (1ppm)  
Benzene in air (100ppm)  
1,3-Butadiene in nitrogen (10ppm)  
Carbon dioxide in helium (100ppm)  
Carbon dioxide in nitrogen (100ppm)  
Carbon dioxide in nitrogen (1000ppm)  
Ethylene in air (8-10ppm)  
Ethylene in helium (100ppm)  
Hydrogen in helium (100ppm)  
Hydrogen in nitrogen (1%)  
Hydrogen in nitrogen (100ppm)  
Methane in helium (100ppm)  
Methane in nitrogen (100ppm)  
Methane in nitrogen (1%)  
Nitrogen in helium (100ppm)  
Nitrous oxide in nitrogen (1ppm)  
Oxygen in helium (100ppm)  
Oxygen in nitrogen (2%)  
Oxygen in nitrogen (6%)  
1,1,1-Trichloroethane in nitrogen (10ppm)  
Trichloroethylene in nitrogen (10ppm)  
Vinyl chloride in nitrogen (1ppm)  
Vinyl chloride in nitrogen (10ppm)  
Vinyl chloride in nitrogen (50ppm)  
Vinyl chloride in nitrogen (100ppm)  
Vinyl chloride in nitrogen (1000ppm)

##### Multi-Component Mixtures

Carbon monoxide, carbon dioxide, hydrogen and oxygen in nitrogen (0.5% each)  
Carbon monoxide, carbon dioxide, hydrogen and oxygen in nitrogen (1% each)  
Carbon monoxide, carbon dioxide, methane, ethane, ethylene and acetylene in nitrogen (1% each)  
Carbon monoxide, carbon dioxide, nitrogen, and oxygen, (5% each) and methane and hydrogen (4% each) in helium  
Carbon monoxide (7%), carbon dioxide (15%) and oxygen (5%) in nitrogen  
Carbon monoxide (7%), oxygen (7%), carbon dioxide (15%) and methane (4.5%) in nitrogen  
C1-C6 *n*-Paraffins: methane, ethane, propane, butane, pentane, hexane in nitrogen (15ppm each)  
C1-C6 *n*-Paraffins: methane, ethane, propane, butane, pentane, hexane in helium (100ppm each)  
C1-C6 *n*-Paraffins: methane, ethane, propane, butane, pentane, hexane in helium (1000ppm each)  
C1-C6 *n*-Paraffins: methane, ethane, propane, butane, pentane, hexane in nitrogen (100ppm each)  
C2-C4 Alkynes: acetylene, propylene, 1-butylene, 2-butylene in nitrogen (15ppm each)  
C2-C6 Olefins: ethylene, propylene, 1-butene, 1-pentene, 1-hexene in helium (100ppm each)  
C2-C6 Olefins: ethylene, propylene, 1-butene, 1-pentene, 1-hexene in nitrogen (100ppm each)  
Branched Paraffins: 2,2-dimethylbutane, 2,2-dimethylpropane, iso-butane, 2-methylbutane, 2-methylpentane, 3-methylpentane in nitrogen (15ppm each)  
Methane, ethane, ethylene, acetylene, propane, propylene, *n*-butane in nitrogen (15ppm each)  
*n*-butane, iso-butane, *cis*-2-butene, *trans*-2-butene, 1-butene, iso-butylene, 1,3-butadiene, ethyl acetylene in nitrogen (15ppm each)

Shelf Life

Scotty 4 (4 Liter) cat.#

Scotty 14 (14 Liter) cat.#

Scotty 48 (48 Liter) cat.#

2 yrs.	34447	34448	34449
2 yrs.	34456	34457	—
2 yrs.	34450	34451	34452
2 yrs.	—	34453	—
2 yrs.	—	34454	—
2 yrs.	—	34455	—
1 yr.	—	—	34458
1 yr.	—	—	34459
2 yrs.	—	34460	34461
2 yrs.	—	34462	—
2 yrs.	—	34463	34464
2 yrs.	—	34465	34466
2 yrs.	—	34467	34468
2 yrs.	—	34489	—
2 yrs.	—	34469	—
2 yrs.	34470	34471	34472
2 yrs.	—	34473	34474
2 yrs.	34475	34476	34477
2 yrs.	—	34478	—
2 yrs.	34481	34482	34483
2 yrs.	—	34479	—
2 yrs.	—	34484	34485
2 yrs.	—	34480	—
2 yrs.	34486	34487	34488
2 yrs.	34490	34491	34492
2 yrs.	—	—	34493
2 yrs.	—	34494	34495
2 yrs.	—	34496	34497
2 yrs.	—	34498	34499
2 yrs.	—	34500	—
2 yrs.	—	34501	—
2 yrs.	—	34502	—
2 yrs.	34503	34504	34505
2 yrs.	34506	34507	34508
1 yr.	34509	34510	34511
2 yrs.	—	34512	—
2 yrs.	34513	34514	—
2 yrs.	—	34515	34516
2 yrs.	34517	34518	34519
2 yrs.	34520	34521	34522
2 yrs.	34523	34524	34525
2 yrs.	34526	34527	34528
2 yrs.	34535	—	—
2 yrs.	—	34529	34530
2 yrs.	—	34531	34532
2 yrs.	34533	34534	—
1 yr.	34536	—	34537
1 yr.	34538	—	34539



# European Pharmacopoeia Reference Materials

## Formulated for ICH Solvent Guidelines

by Ken Herwehe, Analytical Reference Materials Marketing Manager

- ✓ All ICH Class 1 and Class 2 solvents included.
- ✓ Class 2 solvents separated into 3 mixtures for less complex chromatograms.

The analysis of residual solvents in pharmaceutical products has changed, particularly for products being sold into Europe. The International Conference on Harmonization (ICH) Guidelines for Residual Solvents is becoming the international standard and is being adopted by more Pharmacopoeias, including the United States Pharmacopoeia, every year.

The ICH method and compound list is more extensive than any method previously used and poses new challenges. Compounds in Class 1 are solvents considered to be of highest risk and to be avoided

in pharmaceutical manufacturing. Use of Class 2 compounds is to be limited, as they pose a lower, but present, threat to health. Compounds in Class 3 pose the lowest toxic potential and may be used routinely in manufacturing.


Restek manufactures Rtx®-1301 (cat.# 16085) and Stabilwax® (cat.# 10640) capillary columns and formulated these analytical reference materials to meet the new requirements.\* Please request cat.# 59107, *European Pharmacopoeia Analysis of Residual Solvents* for more information.

**Figure 1**

Excellent resolution of most European Pharmacopoeia Class 1 and Class 2 compounds at the regulation limit concentration, using an Rtx®-1301 column.

- |   |  |
|---|--|
| 1. methanol                                   | 16. methylcyclohexane                  |
| 2. 1,1-dichloroethene                         | 17. 1,4-dioxane                        |
| 3. acetonitrile                               | 18. pyridine                           |
| 4. dichloromethane<br>(methylene chloride)    | 19. toluene                            |
| 5. hexane (C6)                                | 20. 2-hexanone                         |
| 6. cis-1,2-dichloroethene                     | 21. chlorobenzene                      |
| 7. nitromethane                               | 22. DMF                                |
| 8. chloroform                                 | 23. ethylbenzene                       |
| 9. cyclohexane                                | 24. m-xylene                           |
| 10. 1,1,1-trichloroethane                     | 25. p-xylene                           |
| 11. carbon tetrachloride                      | 26. o-xylene                           |
| 12. benzene                                   | 27. N,N-dimethylacetamide              |
| 13. 1,2-dimethoxyethane                       | 28. 1,2,3,4-tetrahydro-<br>naphthalene |
| 14. 1,2-dichloroethane<br>(trichloroethylene) |  |

**Column:** 30m, 0.53mm ID x 3.0µm Rtx®-1301 (cat.# 16085)  
**Inj:** Headspace injection of 28 Class 1 and 2 residual solvents for pharmaceutical processing. Prepared at the regulatory limit concentration, using samples shaken and heated at 80°C for 15 minutes. 1mL headspace injection.  
**Oven temp.:** 40°C (hold 20 min.) to 240°C @ 10°C/min. (hold 20 min.)  
**Inj./det. temp.:** 200°C/250°C  
**FID sensitivity:** 1.1 x 10<sup>-11</sup> AFS  
**Carrier gas:** hydrogen @ 35cm/sec.  
**Split ratio:** 2:1



Restek reference materials are formulated to correspond to the permissible daily exposure.

### European Pharmacopoeia/ ICH Class 1 Mix

benzene	2µg/mL
carbon tetrachloride	4
1,2-dichloroethane	5
1,1-dichloroethene	8
1,1,1-trichloroethane	1500

In water:dimethylsulfoxide (90:10), 1mL/ampul

Each	5-pk.	10-pk.
36228	36228-510	36328

### European Pharmacopoeia/ ICH Class 2 Mix A

(13 components)

chlorobenzene	360µg/mL
cyclohexane	3,880
cis-1,2-dichloroethene	1,870
dichloromethane	600
N,N-dimethylformamide	880
ethylbenzene	369
hexane	290
methylcyclohexane	1,180
toluene	890
1,1,2-trichloroethene	80
m-xylene	1,302
o-xylene	195
p-xylene	304

In dimethylsulfoxide, 1mL/ampul

Each	5-pk.	10-pk.
36229	36229-510	36329

### European Pharmacopoeia/ ICH Class 2 Mix B

(10 components)

acetonitrile	410µg/mL
chloroform	60
1,2-dimethoxyethane	100
N,N-dimethylacetamide	1,090
1,4-dioxane	380
2-hexanone	50
methanol	3,000
nitromethane	50
pyridine	200
1,2,3,4-tetrahydronaphthalene (tetraline)	100

In water:dimethylsulfoxide (90:10), 1mL/ampul

Each	5-pk.	10-pk.
36230	36230-510	36330

### European Pharmacopoeia/ ICH Class 2 Mix C

2-ethoxyethanol	160µg/mL
ethylene glycol	620
formamide	220
2-methoxyethanol	50
N-methylpyrrolidone	4,840
sulfolane	160

In water, 1mL/ampul

Each	5-pk.	10-pk.
36231	36231-510	36331





# Enhanced Retention of Polar Analytes by HPLC

## Using the New Ultra Aqueous C18 Column

by Terrence Reid, HPLC Applications Team Chemist

- ✓ Excellent peak shape for basic analytes.
- ✓ Compatible with 100% aqueous to 100% organic mobile phases.
- ✓ Compatible with MS detection.

The newest addition to our selection of HPLC columns, the Ultra Aqueous C18 column, is designed to enhance the retention of polar compounds by reversed phase HPLC. The Ultra Aqueous C18 stationary phase is a true C18 chain (USP L1), but it is immobilized on the silica surface through a

unique chemistry that creates polar groups on the silica surface, between the C18 chains (Figure 1). This secondary polar character has several benefits. First, polar analytes that are insufficiently retained on a conventional C18 column interact with the polar groups in an Ultra Aqueous C18 column, producing enhanced retention. Second, the polar groups aid the retention of polar compounds by keeping the stationary phase completely wetted, even in 100% aqueous mobile phases. In theory, eliminating organic solvent from the mobile phase should maximize retention in reversed phase HPLC, and this is true for Ultra Aqueous C18 columns. In contrast, many conventional C18 columns lose ability to retain analytes in highly aqueous mobile phases because the C18 chains self-associate or fold down on the silica (Figure 1), a phenomenon sometimes referred to as chain folding. Third, the polar groups on the Ultra Aqueous C18 stationary phase shield analytes from active silanol sites on the silica surface, ensuring excellent peak shape for basic analytes (Figure 2).

Although they were designed to be used with highly aqueous mobile phases, Ultra Aqueous C18 columns also are completely compatible with highly organic mobile phases. The ability to cover the full range of mobile phase composition, from 100% aqueous to 100% organic, is useful for developing gradient methods for analyzing samples containing

both highly polar and nonpolar analytes. For example, Figure 3 shows water-soluble vitamins are eluted from an Ultra Aqueous C18 column with excellent resolution and as sharp, symmetric peaks.

Ultra Aqueous C18 columns also are compatible with MS detection—the minimal noise generated by an Ultra Aqueous C18 column is comparable to background in a blank analysis with no column in line (Figure 4). LC/MS is the analytical approach for a steadily increasing variety of analytes, and the versatile Ultra Aqueous C18 column is an obvious column choice.

If you are analyzing samples containing polar analytes, or mixtures of polar and nonpolar analytes, and are contending with unsatisfactory resolution, insufficient retention of polar compounds, poorly shaped peaks, and/or complicated mobile phases, an Ultra Aqueous C18 column can be the solution to your problems.









Figure 4

Minimal background makes Ultra Aqueous C18 columns ideal for LC/MS.



## Ordering Information | Ultra Aqueous C18 5µm Columns

	1.0mm ID	2.1mm ID	3.2mm ID	4.6mm ID
Length	cat.#	cat.#	cat.#	cat.#
30mm	9178531	9178532	9178533	9178535
50mm	9178551	9178552	9178553	9178555
100mm	9178511	9178512	9178513	9178515
150mm	9178561	9178562	9178563	9178565
200mm	9178521	9178522	9178523	9178525
250mm	9178571	9178572	9178573	9178575

## FREE HPLC poster!

This helpful poster features step-by-step instructions on HPLC setup and storage and many useful reference charts. Request lit. cat.# 59894.

## Replacement Parts for Waters HPLC Instruments

## Designed to Meet Original Specifications

by Greg France, HPLC Product Marketing Manager



Our replacement parts for Waters, Agilent, and Shimadzu instruments are built to meet or exceed the original equipment manufacturers'

(OEM) specifications. To ensure instrument compatibility and the quality you expect from an OEM product, we obtain many of the components from the same suppliers the OEMs use. Currently we offer detector, autosampler, and pump parts. As we learn what additional needs you have we will continue to expand this line. Listed here is a partial selection of Waters parts. For other Waters parts, or for Agilent or Shimadzu parts, please visit our website or call customer service or HPLC technical service. If we don't stock what you need, just ask—we'll do our best to get it for you.

## HPLC Piston Seal Insertion Tool

- Simplify your pump maintenance
- One end removes old piston seal, other easily and securely installs new seal



Description	qty.	cat.#
HPLC Piston Seal Insertion Tool	ea.	21356

Description	Waters Instrument Model #	Waters™ part #	Equivalent Restek Part qty.	cat.#
Inlet Check Valve Assembly	M6KA, 501, 510, 515, 590, 600E	33679, 25214	ea.	25360
Inlet Check Valve Housing	M6KA, 501, 510, 515, 590, 600E	25203	ea.	25361
Inlet Check Valve Rebuild Kit	M6KA, 501, 510, 515, 590, 600E	60495	2-pk.	25362
Outlet Check Valve Assembly (Actuator Style)	M6KA, 501, 510, 515, 590, 600E	25030	ea.	25363
Outlet Check Valve Housing (Actuator Style)	M6KA, 501, 510, 515, 590, 600E	25212	ea.	25364
Outlet Check Valve Rebuild Kit (Actuator Style)	M6KA, 501, 510, 515, 590, 600E	26016	2-pk.	25365
Outlet Check Valve Assembly (Ball & Seat Style)	M6KA, 501, 510, 515, 590, 600E	25216	ea.	25366
Outlet Check Valve Housing (Ball & Seat Style)	M6KA, 501, 510, 515, 590, 600E	25207	ea.	25367
Outlet Check Valve Rebuild Kit (Ball & Seat Style)	M6KA, 501, 510, 515, 590, 600E	26014	2-pk.	25368
Inlet Check Valve Assembly, 225µL (Extended Flow)	M6KA, 501, 510, 515, 590, 600E	60307	ea.	25369
PerformancePLUS™ Check Valve Housing	M6KA, 501, 510, 515, 590, 600E	70000254	2-pk.	25370
Check Valve Rebuild Kit (Extended Flow)	M6KA, 501, 510, 515, 590, 600E	88223	2-pk.	25371
PerformancePLUS™ Check Valve Housing	M6KA, 501, 510, 515, 590, 600E	—	ea.	25372
Round Pump Head w/ Actuator Outlets	M6KA, 510, 590, 600	60058	ea.	25413
Round Pump Head, Ball & Seat Check Valves	M6KA, M45, 501	—	ea.	25414
Check Valve Cartridges	Alliance™	WAT270941	2-pk.	25373
Super Seal™ for Analytical Heads	M6KA, 501, 510, 515, 590, 600E	22946, 22934	ea.	25374
Plunger Seal, Gold (Analytical Heads)	M6KA, 501, 510, 515, 590, 600E	22934	ea.	25375
Plunger Seal, Tan	M6KA, 501, 510, 515, 590, 600E	25384	ea.	25376
Plunger Seal, Black	M6KA, 501, 510, 515, 590, 600E	26613	ea.	25378
Plunger Seal, Gold (EF Heads)	510, 590, 600E	26644	ea.	25380
Seal Wash Plunger Seal	Alliance™	WAT271018	2-pk.	25386
Head Plunger Seal Kit	Alliance™	WAT270938	2-pk.	25387
Insert Seal Parts Kit	M6KA, 501, 510, 515, 590, 600E	60012	kit	25389
Sapphire Plunger	M6KA, 510, 590, 600	25656	ea.	25381
Sapphire Plunger	M45, M501	26524	ea.	25383
Sapphire Plunger	M515	WAT207069	ea.	25384
Sapphire Plunger	616, 625, 626	31788	ea.	25420
Sapphire Plunger	Alliance™	WAT270959	ea.	25385
Single Solvent Inlet Manifold	600E	60034, 60042	ea.	25390
Gradient Proportioning Valve, 12Volt	600E	62037	ea.	25419
Wash Face Seal	Alliance™ 2690	WAT271017	ea.	25428
Wash Tube Seal Kit	Alliance™ 2690	WAT270940	4-pk.	25429
Proportioning Valve	Alliance™ 2690	WAT270927	ea.	25430
Xenon Lamp	474	—	ea.	25405
Deuterium Lamp (UV/Vis)	484	80357	ea.	25406
Deuterium Lamp	996, 2996	WAT052586	ea.	25408
Deuterium Lamp	2487	WAS081142	ea.	25409



# Super-Clean™ Gas-Trapping System for LC/MS

## New Quick-Change System Removes Hydrocarbon Impurities from Nitrogen

by Donna Lidgett, GC Accessories Product Marketing Manager

- ✓ Fast, easy cartridge changes—no tools, no purge time, no potential contamination.
- ✓ Durable, patented full glass/metal design.

The Super-Clean™ Gas-Trapping System is the latest technology in cartridge systems for gas purification, and is ideal for purifying nitrogen for LC/MS systems. Changing cartridges is quick and easy. A two-position base plate, installed in the gas line, allows you to exchange cartridges without introducing oxygen into the system: spring-loaded check valves seal when a cartridge is removed and open only when a new cartridge has been locked in place. There is no need for loosening and tightening fittings every time you change cartridges, and your system cannot become contaminated during the changing process.

To meet the high flow needs of an LC/MS system, the two cartridges are positioned and connected in parallel. The incoming gas stream is split equally between the two cartridges, and the two streams are rejoined after purification but before the gas exits the base plate. This approach allows longer contact between the nitrogen and the charcoal adsorbent, ensuring higher gas purity. The system can deliver up to 20L of 99.9999% pure nitrogen per minute, at a maximum pressure of 11 bar (160psi). Estimated cartridge lifetime is 3 to 6 months.



### Special Offer

Two Super-Clean™ filtration systems:  
(2) 2-Position baseplates & (4) charcoal filters  
cat.# 22063 Offer ends 12/31/02.

### Ordering Information | Super-Clean™ Gas-Trapping for LC/MS

Description	qty.	cat.#
LC/MS 2-Position Base Plate with 1/4" Fittings	ea.	22060
Charcoal Replacement Filters	2-pk.	22061
Super-Clean™ Gas-Trapping System (includes 2-position base plate and 2 charcoal replacement filters)	ea.	22062

## New Reference Books

## Higher Education Doesn't have to be Expensive

by Jack Crissman, Seminar and Educational Products Manager

**For review of these books, visit our website.**

**The Merck Index. 13<sup>th</sup> Edition:** S. Budvari. Merck. 2001. 2560pp., ISBN 0-91-191013-1 cat.# 21383 (ea.)

**A Century of Separation Science:** Haleem J. Issaq. Marcel Dekker, Inc., 2002. 755pp., ISBN 0-8247-0576-9 cat.# 20473 (ea.)

**Chromatography Theory:** Jack Cazes and Raymond P.W. Scott. Marcel Dekker, Inc., 2002. 475pp., ISBN 0-8247-0778-8 cat.# 21573 (ea.)

**HPLC of Biological Macromolecules: Second Edition, Revised and Expanded:** Karen M. Gooding and Fred E. Regnier. Marcel Dekker, Inc., 2002. 777pp., ISBN 0-8247-0665-X cat.# 21574 (ea.)

**Advances in Chromatography, Volume 41:** Phyllis R. Brown and Eli Grushka. Marcel Dekker, Inc., 2001. 425pp., ISBN 0-8247-0509-2 cat.# 21575 (ea.)

**Milestones in the Evolution of Chromatography:** Leslie S. Ettre. ChromSource, Inc., 2002. 220pp., ISBN 0-9717144-0-1 cat.# 20472 (ea.)

**Chromatography in Food Science and Technology:** Tibor Cserháti and Esther Forgács. CRC Press, LLC, 1999. 552pp., ISBN 1-56676-749-0 cat.# 21492 (ea.)

**Gas Chromatographic Techniques and Applications:** Alan J. Handley and Edward R. Ardland. CRC Press, LLC, 2001. 320pp., ISBN 0-8493-0521-7 cat.# 21491 (ea.)

**Handbook of Chemistry and Physics, 83<sup>rd</sup> Edition:** D. R. Lide. CRC Press, LLC, 2002. 2,672pp., ISBN 0-8493-0483-0 cat.# 21442 (ea.)

**Multidimensional Chromatography:** L. Mondello, A. C. Lewis and K. D. Bartle. John Wiley, 2002. 436pp., ISBN 0-471-98869-3 cat.# 21443 (ea.)

**Modern Derivatization Methods for Separation Sciences:** T. Toyo'oka. John Wiley, 1999. 298pp., ISBN 0-471-98364-0 cat.# 21444 (ea.)

### A Wealth of Practical Chromatography Experience, at a Location Near You

Again this year the chromatography wizards from Restek are presenting comprehensive seminars designed to help you minimize downtime and obtain the results you want. We keep these seminars to one day, and hold them at sites all around the country, to minimize the time you spend away from your lab and conserve your travel budget. The low cost of the seminar is an investment that can be quickly returned, because you will improve your lab throughput and spend less time dealing with problems. Choose the topic that suits your work:

**Comprehensive Capillary GC  
Environmental GC Analysis  
Comprehensive HPLC**

**Food, Flavor, and Fragrance Analysis**

Our brochure 2002 Seminars (lit. cat.# 59282A) provides details about these seminars and lists dates and locations. Call, fax, or e-mail your request for the seminars brochure today, or view it on our website. We look forward to meeting you.



### Special Offer! Capillary Column Installation Video (CD-ROM)

Covers the critical points in installing a capillary GC column: instrument preparation, setting gas flows, leak checks, etc., produced by the technical wizards of Restek.  
cat.# 20499, (ea.)

www.restekcorp.com



**HROMalytic** +61(0)3 9762 2034  
**ECHnology** Pty Ltd

Website NEW : [www.chromalytic.net.au](http://www.chromalytic.net.au) E-mail : [info@chromtech.net.au](mailto:info@chromtech.net.au) Tel: 03 9762 2034 . . . in AUSTRALIA

**Australian Distributors**  
Importers & Manufacturers  
[www.chromtech.net.au](http://www.chromtech.net.au)





# cool tools

by Brad Rightnour and Michael Goss, Instrument Innovations Team


For Easier GC Maintenance Try These New Tools from Restek

## Mini Wool Puller/Inserter

A wool plug that is incorrectly positioned or contaminated with finger oils can be more hindrance than help. This inexpensive little tool greatly simplifies the chore of consistently placing contaminant-free wool plugs in an inlet liner, and retrieving a plug when its time to replace it. We suggest you order several packages-or be ready to spend time trying to find out who's borrowed yours. cat.# 20114, (2-pk.)



Place a 1cm plug of loosely bound wool in the liner. Adjust its position with the puller/insertor tool.



Use the hooked end of the puller/insertor to retrieve the plug when it's time to replace it.

Use with conventional 2mm ID or 4mm ID liners and most other liner configurations, but not with double gooseneck liners.

## MS Installation Gauge

Easily pre-seat ferrules for consistent installations in Agilent 5973 MS!

- ✓ Prevents damage to the column end
  - ✓ Ensures leak-free connection
- cat.# 21894, (ea.)



**new**




Install the nut and ferrule onto the column, then insert the installation tool, exposing several centimeters at the exit end.



Tighten the nut.



Score and remove the exposed end of the column making sure of a clean, square cut, then loosen the nut.



The ferrule will be properly seated and should remain in place when light force is applied. Install the column into the GC/MS interface.

## Pre-Cleaned Sample Vials

Ready to Use for Volatiles Analyses

by Donna Lidgett, GC Accessories Product Marketing Manager

- ✓ Container, liner, and closure cleaned, assembled, and ready to use.
- ✓ Open-top caps.
- ✓ Teflon<sup>®</sup>-faced 0.125" silicone septa.
- ✓ Each case lot numbered.
- ✓ Available in clear or amber.

**new**



Description	qty.	cat.#
20ml CLEAR Pre-Cleaned VOA Vials	72-pk.	21798
20ml AMBER Pre-Cleaned VOA Vials	72-pk.	21799
10ml CLEAR Pre-Cleaned VOA Vials	72-pk.	21796
40ml AMBER Pre-Cleaned VOA Vials	72-pk.	21797



# RESTEK

## Behind the Scenes

### Happy Birthday, USA

Restek continues to be a major sponsor of one of the world's largest fireworks displays—Centre County, PA's 4th Fest. More than 100,000 people participated in part or all of the day-long festivities, which culminated in an aerial display of 12,000 shells. Restek also supplied the birthday cake—all 6 feet by 12½ feet of it. Approximately 1100 people had a slice, in celebration of our country's 226th birthday. Restek's founder and head coach, Paul Silvis, is a long-time supporter and co-chair of the event.

### Restek Wins Healthy Workplace Award

Restek Corporation has been awarded the Healthy Workplace Award for Small Sized Companies by the Pennsylvania Psychological Association. This award is given annually to companies that demonstrate a commitment to family support, employee development, employee involvement, community involvement, and health and safety in the workplace. Last year, Restek was awarded honorable mention. The award was presented on June 21 at an official award ceremony in Lancaster.

Restek offers numerous employee-friendly benefits, including contributions toward child care costs, reimbursement for continuing education and development, on-site fitness and recreational facilities with subsidized personal trainers, 401k and employee stock ownership programs, open-book management, and a casual dress code. For the last 2 years, the company has won awards for being among the 100 Best Companies to Work for in Pennsylvania.

### Restek Employees Prove They're Healthy and Involved

Proving our Healthy Workplace Award was no accident, two Restek employees, Becky Wittrig (Innovations Team chemist) and Matt Reilly (ATM technician) took first place in their divisions in a 5K run benefiting a local library, in June. Becky's son, James, Matt's daughter, Leena, and Alex Reid, son of HPLC Applications Team chemist Terry Reid, ran in the children's divisions. Our Employee Action Group made a \$100 sponsorship donation.

Also in June, Restek's Relay for Life team won the Bronze Award at the local Relay event by raising more than \$1500 for the American Cancer Society.

### New Literature

- ✓ Narcotics / Acetaminophen by HPLC Applications Note (lit. cat.# 59453)
- ✓ Organophosphorus Pesticides by Capillary GC Applications Note (lit. cat.# 59359)
- ✓ Calibration Standards for ASTM Method D2887-01 Fast Facts (lit. cat.# 59383A)
- ✓ Certified PAHs in Diesel Fuel #2 Fast Facts (lit. cat.# 59384A)
- ✓ Environmental Gas Standards Fast Facts (lit. cat.# 59276)
- ✓ Pesticide Reference Materials Fast Facts (lit. cat.# 59446)
- ✓ Pinnacle II™ Amino HPLC Columns Fast Facts (lit. cat.# 59385A)
- ✓ Rtx®-5Sil MS Capillary Columns Fast Facts (lit. cat.# 59323)
- ✓ UST Products for Massachusetts Fast Facts (lit. cat.# 59391)
- ✓ UST Products for the Northwest Region Fast Facts (lit. cat.# 59396)
- ✓ UST Products for Wisconsin Fast Facts (lit. cat.# 59392)
- ✓ Rtx®-VMS Capillary Columns New Product Flyer (lit. cat.# 59209A)
- ✓ Bonded PLOT Columns—Flyer (lit. cat.# 59456)
- ✓ Integra-Guard™ Capillary Columns Flyer (lit. cat.# 59441)
- ✓ Review of Restek Literature for Pharmaceuticals Analyses—Flyer (lit. cat.# 59450)
- ✓ Genuine Restek Replacement Parts for Agilent GCs—Catalog (lit. cat.# 59627C)

**Restek USA:** Phone: 800-356-1688 or 814-353-1300 • Fax: 814-353-1309 • [www.restekcorp.com](http://www.restekcorp.com)  
**Restek France:** Phone: 01 60 78 32 10 • Fax: 01 60 78 70 90 • [restekfr@club-internet.fr](mailto:restekfr@club-internet.fr)  
**Restek GmbH:** Phone: 49 06172 2797 0 • Fax: 49 06172 2797 77 • [RESTEK-GMBH@t-online.de](mailto:RESTEK-GMBH@t-online.de)  
**Restek Ireland:** Phone: 44 28 9081 4576 • Fax: 44 28 9081 4576 • [restekeurope@aol.com](mailto:restekeurope@aol.com)  
**Thames Restek UK LTD:** Phone: 01753 624111 • Fax: 01753 624666 • [Sales@Thamesrestek.co.uk](mailto:Sales@Thamesrestek.co.uk)

**Restek Trademarks:** Integra-Guard, Leak Detective, MegaMix, MXT, Pinnacle II, Press-Tight, Resprep, Res-Sil, Rtx, RT-MSieve, RI-XI, Sulfur, Silcoport, SilcoSmooth, Silcosteel, Siltek, Stabilwax, Stx, Sulfinert, Thermolite, Trident, Uniliner, and the Restek logo.

**Other trademarks:** Agilent (Agilent Technologies, Inc.); Aroclor (Monsanto Co.); Carbowax (Union Carbide Corp.); Hastelloy (Hayes International); Microseal (Merlin Instrument Co.); Pittcon (The Pittsburgh Conference); PEEK (Victrex plc); Scotty (Scott Specialty Gases, Inc.); Super-Clean (SGT Middleburg BV); Teflon, Vespel, and Viton (E.I. du Pont de Nemours & Co., Inc.); Alliance, Durapack, Porasil, and Waters (Waters Associates, Inc.)



Lit. Cat.# 59461-INT

©2002 Restek Corp.

Please direct your comments on this publication to Carrie Sprout, Graphic Designer, at [carrie@restekcorp.com](mailto:carrie@restekcorp.com) or call 814-353-1300, ext. 2151



"some" promos / Products / Offers in the ADVNews  
have been since been progressively superceded  
/ UPDATED OR Since Discontinued

CHECK THE latest Restek ADVantage Newsletter, Restek ESSENTIALS  
... Or The Restek Catalog ... Or other Restek publications for updates  
[www.chromtech.net.au](http://www.chromtech.net.au) or NEW site 2015 > [www.chromalytic.net.au](http://www.chromalytic.net.au)



Website NEW : [www.chromalytic.net.au](http://www.chromalytic.net.au) E-mail : [info@chromtech.net.au](mailto:info@chromtech.net.au) Tel: 03 9762 2034 ... in AUSTRALIA

**Australian Distributors**  
Importers & Manufacturers  
[www.chromtech.net.au](http://www.chromtech.net.au)