FLUIDIC PRODUCTS AND INFORMATION








Intelligent Solutions for Life[™]

HRONalytic+61(0)3 9762 2034

ECH 1000 Y Pty Ltd

Website NEW: www.chromalytic.net.au E-mail: info@chromtech.net.au Tel: 03 9762 2034 . . . in AUSTRALIA

Fluidics | Optics | Consumables | Assemblies



For your information, our U.S. Federal Taxpayer Identification number is 01-0736657; European Taxpayer Identification number is DE 146581180.

ORDER PROCESSING/ SHIPPING

Orders placed on our website process and ship as specified during the checkout process. Orders that are placed via fax/email will be subject to a \$100 minimum order quantity and typically process within one business day and ship within two business days. Customers who have established contract pricing terms must place catalog orders directly with their Customer Service representative.

All shipments are FCA unless otherwise noted. For orders not placed via our website, shipping charges will be prepaid and added to the invoice. Or, if requested, we will charge your carrier account number.

TERMS

Credit cards are accepted for orders placed within the United States. Purchase orders may also be used by qualified customers. No international credit card payment is possible. No international shipments are possible. For terms and conditions of sale please visit www.idex-hs.com/ legal-notices. Terms and prices are subject to change without notice.

RETURN OF MATERIALS

Returns must be authorized in advance. Please contact us within 30 days of purchase for your Returns Material Authorization (RMA) Number. A restocking fee may apply. To request a Return Materials Authorization for IDEX Health & Science products visit www.idex-hs.com/return.

Product purchased through a Distributor must be returned directly to that Distributor. Please contact them for Return Material Authorization procedures.

DIRECT FROM OUR GLOBAL **DISTRIBUTION NETWORK**

IDEX Health & Science has established relationships with a global network of Distributors that offer outstanding product support and assistance. In order to meet our customer's needs, and the demands of being a global supplier, we have come to rely on our authorized Distributors to provide value-added service and support of our products. While we are happy to take your order directly, we encourage you to contact your local IDEX Health & Science Distributor.

HOW TO FIND A LOCAL **DISTRIBUTOR FOR FLUIDICS PRODUCTS**

For a complete distributor listing please visit www.idex-hs.com/ distributors or email:

US:

IHSOrders@idexcorp.com

Europe:

CustomerService.hsEurope@ idexcorp.com

Asia:

CustomerService.hsAsia@ idexcorp.com


Our customer service team will help you contact a distributor in your area.

HOW TO FIND A LOCAL **DISTRIBUTOR FOR SEMROCK PRODUCTS**

For a complete distributor listing please visit www.semrock.com/ international-distributors.aspx

The products in this catalog are intended for use with analytical, biotechnology and diagnostic equipment and accessories. THE MANUFACTURERS OF THE PRODUCTS IN THIS CATALOG SHALL HAVE NO LIABILITY WHATSOEVER DUE TO ANY MISUSE OF THEIR PRODUCTS. Of course, the safe use of our products depends on our customers, since it is you who select and control the protective gear and safety procedures used, as well as the pressures, temperatures, solvents, samples, ventilation, and other variables. Product and material performance ratings are provided as guides only. Individual field tests should be performed by customers to dete safe operating parameters given your particular procedures and use. IDEX Health & Science facilities are certified ISO 9001.





WE HAVE A NEW VISION FOR THE WORLD OF FLUIDICS & OPTICS

As your trusted partner in life science instrument development, we continue to invest in leading technologies and capabilities to solve your most demanding challenges.

For this reason, we have strategically integrated innovative optical companies, like Semrock and Melles Griot, to bring you a newly enhanced portfolio of state-of-the-art components and capabilities that are unrivaled in breadth, performance, quality and design.

Now, as the global authority in fluidics and optics, we deliver a complete spectrum of highly engineered solutions, and more powerful platforms than ever before.



Our vision of the complete path goes far beyond just meeting your needs — it anticipates them, with intelligent solutions for life.



GLOBAL PREMIUM DISTRIBUTOR

Biotech AB in Sweden/Europe (Headquarter), Biotech USA LLC and Bionik Inc Japan are Global Premium Distributors of IDEX Health & Science.

FAST DELIVERY

With several weekly shipments from IDEX manufacturing sites in US, together with our extensive Swedish warehouse, we support the European & Global market with short delivery times as well as fast answers on technical questions and a high level of flexibility.

LOCAL DISTRIBUTORS

In almost every European country we work in close co-operation with local distributors that supports you with excellent daily service, long experience and highest level of technical assistance.

FULL IDEX ASSORTMENT

With the full IDEX assortment of components for fluidic systems along with our long experience and know-how in the field we can assure you to always do our very best in supplying you with excellent solutions and products. Customized OEM solutions, especially in the field of degassing, is something we offer to instrument manufacturers globally.



BIOTECH + IDEX =EVERYTHING FOR THE LAB

KNOW-HOW AND A WIDE ASSORTMENT

As Global Premium Distributor for IDEX we have special know-how about their assortment and can provide you with the right mix of products for a successful lab.

We have cutting edge knowledge to tailor the system for laboratories' specific needs for fluidic handling. You will find a wide range of products to suit every possible flow rate and type of fluid. In this field, we offer a number of well-known brands, such as Systec, Upchurch, Rheodyne and Isolation Technologies - all to give you the best solutions.

INNOVATIVE PRODUCTS FOR FLUIDIC SYSTEMS



www.biotech.se




12 **FLUIDIC CONNECTIONS**



- **TUBING** 13 PEEK Tubing 16 19 Stainless Steel Tubing 22 PEEKsil™ Tubing
- PFA Tubing 24 26 **FEP Tubing**
- 27 ETFE Tubing
- 28 Tools

- 30 **FITTINGS** 32 One-Piece Fingertight Fittings 33 Stainless Steel Fittings 34 MicroTight® Fittings 35 Two-Piece Fingertight Fittings Two-Piece SealTight™ Fingertight Fittings 36 NanoTight[™] Fittings & Sleeves 37 LiteTouch® Fittings 38 39 Super Flangeless™ Fittings 42 VacuTight™ Fittings 43 Metric Flangeless Fittings 44 Flanged Fittings 45 Flangeless Fittings MarvelX™ 48
- 50 Fittings Tools 52 **Tubing Sleeves** 55 Plugs & Caps 56 Large Bore Fittings 57 VHP Micro Fittings 58 VHP PK Fittings 59 Stainless Steel VHP Fittings Reusable VHP Fittings 60 61 Assorted Fittings Kits



| CONNECTORS | 94 | FILTERS & FRITS |
|------------------------------|---|---|
| Threaded Adapters | 95 | Stainless Steel Frits |
| High Pressure Mixing Tees | 97 | PEEK Frits |
| High Pressure Unions | 99 | Frit-in-a-Ferrule™ |
| High Pressure PEEK | 100 | Bottom-of-the-Bottle™ F |
| Tees & Crosses | 102 | Inlet Solvent Filters |
| High Pressure | 103 | Inline Filters |
| _ | 105 | Precolumn Filters |
| Tees & Crosses | 107 | Semi-Prep Filters |
| VHP Stainless Steel | 108 | Inline MicroFilters |
| ZDV Unions | 109 | Mini MicroFilters |
| VHP Tees & Crosses | 110 | Precolumn MicroFilters |
| VHP MicroTight® Unions | 111 | Bottle Caps |
| MicroTight® Adapters | 112 | Bottle Cap Plugs |
| Accessories | | 1 3 |
| NanoPort Assemblies | | |
| Low Pressure Manifolds | | |
| Low Pressure Bulkhead Unions | | |
| | Threaded Adapters High Pressure Mixing Tees High Pressure Unions High Pressure PEEK Tees & Crosses High Pressure MicroTight® Unions High Pressure Stainless Steel Tees & Crosses VHP Stainless Steel ZDV Unions VHP Tees & Crosses VHP MicroTight® Unions MicroTight® Adapters Accessories NanoPort Assemblies Low Pressure Manifolds | Threaded Adapters 95 High Pressure Mixing Tees 97 High Pressure Unions 99 High Pressure PEEK 100 Tees & Crosses 102 High Pressure MicroTight® Unions 105 High Pressure Stainless Steel 107 VHP Stainless Steel 108 ZDV Unions 109 VHP Tees & Crosses 110 VHP MicroTight® Unions 111 MicroTight® Adapters 112 Accessories NanoPort Assemblies Low Pressure Manifolds |

Low Pressure Y Connectors

Quick Connect Luer Adapters

Low Pressure Unions Low Pressure Tees & Crosses

 $Luer Tight^{\scriptscriptstyle{\mathsf{TM}}} \ Fittings$

Barbed Connectors

83

84

86 87

89

90

114 VALVES

| 119 | Stand Alone Valve Products |
|-----|----------------------------|
| 120 | Rapid Replacement Pods™ |
| 121 | Manual Valves |
| 404 | C D : |

124 Spare Parts 126 Valve accessories

127 Stainless Steel Sample Loops

128 PEEK Sample Loops

131 Suction Needle Adapter 132 Injection Port Adapters

133 Wrenches, Brackets, & Replacement Fittings

FLOW REGULATING 134 **VALVES**

135 Inline Cartridge Check Valves

Standard 1/4-28 Inline 136 Check Valves

137 Non-Metallic Check Valves

Quick-Stop Luer Inline 139 Check Valve

140 Micro-Splitter Valves 141 Micro-Metering Valves

142 Shut-Off Valves

BACK PRESSURE 143 **REGULATORS**

144 Ultra-Low Volume Back Pressure Regulators (BPRs)

145 **Back Pressure** Regulator Assemblies

Back Pressure 146 Regulator Holders

Back Pressure 147 Regulator Cartridges

148 Pressure Relief Valves

150 DEGASSERS

151 Debubblers

154 Full Stand Alone Degassing Systems

156 COLUMN **HARDWARE**

157 Guard Hardware

158 **Guard Columns**

HROMalytic +61(0)3 9762 2034

ECH nology Pty Ltd

Australian Distributors Importers & Manufacurers www.chromtech.net.au

Website NEW: www.chromalytic.net.au E-mail: info@chromtech.net.au Tel: 03 9762 2034 . . . in AUSTRALIA





MarvelX[™]

Advance Your Fluidic Systems with Next Generation UHPLC Connection Technology

Page 48

DEGASSER TUBING KIT

Ready-to-send kits with the tubing in desired lengths and the fittings you need; attached, labeled and packed! We help you choose the best fittings and tubing for every application!

WE HELP YOU EVERY STEP OF THE WAY:

- Wide selection of fittings and tubing
- Filters, tools, tubing markers and accessories
- Labeling with your logo and company details
- Packaging in boxes or plastic bags
- Long experience in this business
- We help you to find solutions for your needs





NEW VERSION!

- White color
- Better ventilation

Page Full motal housing for

HROMalytic +61(0)3 9762 2034 ECH no logy Pty Ltd Website NEW: www.chromalytic.net.au E-mail: info@chromtech.net.au Tel: 03 9762 2034 . . . in AUSTRALIA

Australian Distributors



DEGAS

GO WITH A HIGH FLOW INTRODUCING THE ALL NEW DEGASI PREP+

The World's first High Flow Rate, in-line Degasser, for Organic Solvents

Page 155





ORGANIC SOLVENTS





Biocompatible Products

Look for components with the blue "Bio" icon, which designates products that use materials anticipated to maintain the integrity and improve the analysis of biological samples in your intended applications.



Navigation Breadcrumbs

Quickly pinpoint any product or product family within this catalog to the same navigation path each product is categorized on our website. Look to the bottom of every product page for easy-to-follow breadcrumbs to locate specific categories and sub-categories in sequential order, then visit: www.idex-hs.com.



Website NEW: www.chromalytic.net.au E-mail: info@chromtech.net.au Tel: 03 9762 2034... in AUSTRALIA

Australian Distributors Importers & Manufacurers www.chromtech.net.au



Welcome to Fluidics



We are the premier provider of intelligently engineered fluidic components, assemblies, and integrated solutions for a wide range of life science applications that require precise control and measurement.

Perform precise sample preparation and analysis with our premium fluidic solutions.

FLUIDICS








Fluidic Connections

IDEX Health & Science has developed a comprehensive line of standard and custom tubing, connectors, fittings, and flow control devices that meet the increasingly demanding requirements of today's high performance analytical fluidic systems. We feature specialty, high-performance polymers and distinct materials designed to work with your system needs. We offer unique products such as biocompatible PEEK-lined stainless steel tubing as well as an assortment of high pressure and fluoropolymer tubing. All of our fittings, filters and frits and connectors come in a variety of materials and styles. We can provide micro and nano-scale dimensions and well as custom forming, assembly and kitting. We also offer our RI detector that provide high resolution and low dispersion detection for HPLC applications.



FLUIDICS





TUBING

Our high quality, versatile tubing is offered in a variety of materials and styles to meet your system requirements. Our high pressure tubing includes biocompatible PEEK selections and well as seamless, pre-cut stainless steel. Our flouropolymer tubing is constructed with genuine Teflon™ FEP and PFA resin, and our unique High Purity PFA. Many of our tubing options are color coded for easy detection and some are translucent making it easy to view the fluid pathway. Our tubing is manufactured to precise tight tolerances to ensure dependable product consistency.

- 16 HIGH PRESSURE TUBING
- 24 FLUOROPOLYMER TUBING
- 28 TOOLS

TUBING OVERVIEW & FUNCTIONS

HIGH PRESSURE TUBING











| TUBING: | PEEK | Capillary PEEK | Fused Silica | Stainless Steel | PEEKsil™ |
|-------------|---|---|---|---|---|
| Page | 16 | 16 | 16 | 19 | 22 |
| Description | Biocompatible, chemically inert to most commonly used solvents, PEEK tubing is flexible, offers a very smooth internal surface, and can be easily cut to desired lengths. • Great alternative for stainless steel tubing in high pressure applications • Many sizes available in color scheme to help identify ID | All the benefits of larger sized PEEK tubing, while serving as an excellent alternative to more traditional fused silica and stainless steel capillary tubing. Capillary PEEK tubing is available in a wide range of micro and nanoscale inner diameters. • Available in common capillary tubing sizes with tight tolerances on OD and ID • Tubing sleeves available for capillary tubing connections | Because of the tight tolerances of fused silica's inner diameters, this tubing is used for micro-scale analyses such as micro and nano-HPLC and capillary electrophoresis. • Most commonly used OD and ID sizes available • High quality, polyimide-clad fused silica • Offered in convenient, two meter lengths | Seamless, pre-cut 316 stainless steel tubing meets the exacting requirements of today's analyses. Thorough preparation guarantees that the tubing is truly ready-to-use, with flat-bur-free ends and a clean finish. • Wide selection of outside and inside diameters and lengths • Pre-cut to ensure burr-free, flat connections • Many sizes feature a color-coded band for easy ID identification | PEEKsil is mechanically strong and has ideal characteristics for sealing with metal or polymer fittings. • Comprised of high quality fused silica sheathed by PEEK tubing • Excellent chemical compatibility • Very tight manufacturing tolerances • Good replacement for stainless steel, PEEK, or standard fused silica |

| PRESSURE T | | |
|------------|--|--|
| | | |
| | | |

OD (outside diameter)

1/32" (785 µm), 1/16" (1.55 mm), 1/8" (3.2 mm)

0.0145" (360 µm), 1/32" (785 µm), 0.020" (0.5 mm)

0.0145" (360 µm)

1/32" (785 µm), 1/16" (1.55 mm), 1/8" (3.2 mm)

0.0145" (360 µm), 1/32" (785 µm), 1/16" (1.55 mm)

| ID (inside diameter) | 0.001" (25 µm)- 0.080" (2.0 mm) | 0.001" (25 µm)– 0.020" (0.50 mm) | 0.0008" (20 μm)–0.006" (150 μm) | 0.004" (100 μm)– 0.080" (2.0 mm) | 0.001" (25 μm)– 0.012" (300 μm) |
|----------------------------|--|---|------------------------------------|--|------------------------------------|
| Operating Temp | -51 to 100 °C | -51 to 100 °C | -51 to 100 °C | -51 to 289 °C | -51 to 100 °C |
| Pressure Rating | 500–10,000 psi (34–690 bar) | 2,000–5,000 psi (138–345 bar) | N/A* | N/A* | 10,000 psi (690 bar) |
| Typical Tolerances | ±0.001" (25 μm) for 1/16" OD tubing; ±0.003" (75 μm) for 1/8" OD tubing | ±0.0005" (12.5 μm) | ±0.0004" (10 μm) | ±0.001" (25 µm) for 1/16" OD tubing; ±0.003" (75 µm) for 1/8" OD tubing | ±0.0004" (10 μm) |
| Refractive Index (Clarity) | Opaque | Opaque | 1.78 | Opaque | Opaque |
| pH Range | 0–14 | 0–14 | 0–10 | 1–14 | 0–10 |
| Sterilization Techniques | Gamma irradiation; ethylene oxide; thermal | Gamma irradiation; ethylene oxide; thermal | Ethylene oxide; thermal | Gamma irradiation; ethylene oxide; thermal | Ethylene oxide; thermal |
| Autoclavable? | Υ | Υ | Υ | Υ | Υ |
| *The manufacturer of this | tubing or material does not pu | hlish this specification | | | |

TUBING OD SIZES

this chapter the appropriate OD size tubing for your application.

1 mm OD

Page(s) 26

2 mm OD

Page(s)

26

3 mm

OD

Page(s) 26



Australian Distributors Importers & Manufacurers www.chromtech.net.au

FLUIDICS

FLUOROPOLYMER TUBING



| TUBING: | PFA | High Purity PFA | 360 μm High Purity PFA | FEP | ETFE |
|-------------|---|--|--|--|--|
| Page | 24 | 24 | 24 | 26 | 27 |
| Description | Offers excellent chemical compatibility, plus due to its inner surface smoothness, PFA tubing tends to be more translucent than PTFE tubing. • Offers higher purity and enhanced translucence when compared with other fluoropolymer tubes • Great for more critical, low pressure applications | This polymer tubing is manufactured from a premium grade of PFA — one of the most contaminant-free polymers on the market. • Offers chemical stability, mechanical strength, and purity for applications such as medical, diagnostic, pharmaceutical, biotechnology, and semiconductor • Excellent replacement for PTFE where gas permeability and surface texture are issues • Clarity of tubing makes PFA an excellent choice for monitoring fluid movement | This tubing offers excellent chemical compatibility, transparency, very low contaminant levels and is available in the most commonly-used outside diameter for capillary tubing applications. Replacement for capillary tubing in low pressure applications where excellent chemical compatibility is required Tubing sleeves available for capillary tubing connections | FEP tubing is a great alternative to traditional PTFE tubing, desirable for use because it is chemically inert to most solvents, easy to cut, and translucent for easy monitoring of solutions passing through. • Great for general, low pressure applications • Many sizes available in multiple colors for easy identification • Tight manufacturing tolerances to ensure product consistency | ETFE is chemically inert and more suitable for higher pressure applications (when using aqueous mobile phases) than PTFE, FEP, and PFA. Additionally, because ETFE is more rigid than PTFE, FEP, and PFA, this tubing better resists inner diameter collapse. • Excellent solvent resistance • More durable and less gas permeable than PTFE, FEP, and PFA • Operating temperatures up to 80 °C |

| FLUOROPOLYMER | TUBING SPECIFICATIONS | S | | | |
|-------------------------------|--|---|---|---|--|
| OD (outside diameter) | 1/16" (1.55 mm), 1/8" (3.2 mm) | 1/16" (1.55 mm), 1/8" (3.2 mm), 3/16" (4.8 mm), 1/4" (6.35 mm) | 0.0145* (360 μm) | 1/16" (1.55 mm), 0.080" (2.0 mm), 0.118" (3.0 mm), 1/8" (3.2 mm), 0.157" (4.0 mm), 3/16" (4.8 mm), 1/4" (6.35 mm), 5/16" (7.94 mm) | 1/16" (1.6 mm), 1/8" (3.2 mm), 1/4" (6.35 mm) |
| ID (inside diameter) | 0.020" (0.50 mm)– 0.062" (1.55 mm) | 0.020" (0.50 mm)- 0.188" (4.80 mm) | 0.002" (50 μm)– 0.006" (150 μm) | 0.003" (0.075 mm) – 0.250" (6.35 mm) | 0.010" (0.25 mm)- 0.188" (4.80 mm) |
| Operating Temp | -51 to 80 °C | -51 to 80 °C | -51 to 80 °C | -51 to 50 °C | -51 to 80 °C |
| Pressure Rating | 500–2,000 psi (34–138 bar) | 250–2,000 psi (17–138 bar) | 1,750–3,500 psi (121–241 bar) | 2,500–4,000 psi (172 - 276 bar) | 250–4,000 psi (17–276 bar) |
| Typical Tolerances | ±0.001" (25 μm) for 1/16" OD tubing; ±0.003" (75 μm) for 1/8" OD tubing | ±0.001" (25 μm) or 1/16" OD tubing | ±0.0005" (12.5 μm) | ±0.001" (25 µm) for 1/16" OD tubing; ±0.003" (75 µm) for 1/8" OD tubing | ±0.001" (25 µm) for 1/16" OD tubing; ±0.003" (75 µm) for 1/8" OD tubing |
| Refractive Index (Clarity) | 1.34 | 1.34 | 1.34 | 1.338 | 1.4 |
| pH Range | 0–14 | 0–14 | 0–14 | 0–14 | 0–14 |
| Sterilization Techniques | Ethylene oxide; thermal | Gamma irradiation; ethylene oxide; thermal | Gamma irradiation; ethylene oxide; thermal | Ethylene oxide; thermal | Ethylene oxide |
| Autoclavable? | Υ | Υ | Υ | Υ | Υ |







- ▶ 1/16" or 1/8" outside diameter available
- > Biocompatible, inert, and easily cut
- Great for high pressure applications
- Maximum continuous use temperature: 100 °C

Our PEEK (polyetheretherketone) polymer tubing is biocompatible, chemically inert to most solvents, and can be used to replace stainless steel tubing in most liquid analytical systems. Unlike stainless steel tubing, PEEK tubing is flexible and can be easily cut to desired lengths.

PEEK tubing has a very smooth internal surface, which causes less turbulence than similarly sized metal tubing, contributing to improved resolution of sample bands. Of all our polymer tubing materials, PEEK is the least permeable to gas (see material properties on our website: www.idex-hs.com).

In addition, much of our 1/16" OD tubing is color-coded so different IDs are easily identified. Our proprietary extrusion process ensures color permanence in our tubing.

Our 5' length tubing is rough cut to approximately 5'1". A trim cut should be made before use, especially for smaller ID tubing. PEEK tubing can be cut easily with a razor blade. However for an improved cut, try our Tubing Cutters on page 28.



Capillary PEEK Tubing

- ▶ 360 µm or 1/32" outside diameter available
- IDs as small as 25 μm (0.001")

Capillary PEEK tubing offers all the benefits of larger sized PEEK tubing, while serving as an excellent alternative to more traditional fused silica and stainless steel capillary tubing (see Application Note, right). The capillary tubing can be coupled to many of the products in the Connectors chapter (starting on page 62) and to some of the valves in the Valves chapter (starting on page 114).



Fused Silica Tubing

-) Five inner diameters with most common capillary outside diameter, 360 μm
- > Cut in convenient lengths, up to 2 m

These products are manufactured from synthetic fused silica with a polyimide coating.

FLUIDICS



Because the thru-hole of our 25 µm ID PEEK tubing is very small, it is possible for some fittings to cause the ID to become occluded. Please use caution, especially with wrench-tightened fittings. For more information, please contact IDEX Health & Science or your local Distributor directly.



- What Size PEEK Tubing Should I Use?
- lt is usually safe to use 1/16" OD x 0.010" ID tubing throughout an analytical HPLC system. With a 0.010" ID, the pressure drop across most tubing lengths is negligible, and the ID is small enough to minimize band broadening.
- High pressure semi-prep LC systems will most likely use 1/8" OD tubing.
- Use our 1/32" OD tubing for the high pressure flow path of some microbore HPLC systems.
- Choose 360 µm OD tubing for most capillary systems.
- > PEEK tubing is also available by the inch. Contact your local Distributor or IDEX Health & Science directly for pricing information.



APPLICATION NOTE

- An independent study conducted by a major pharmaceutical company indicated LC-MS chromatographic performance could be improved in some cases by switching the postcolumn transfer line from fused silica to PEEK polymer tubing. The switch dramatically reduced peak tailing and eliminated the degradation of peak symmetry as injection volume was reduced. For more information, please contact us or order the "Improved LC-MS Results Study" from the "Literature Request" section of our website at www.idex-hs.com.
- To straighten PEEK polymer tubing, first choose a piece of stainless steel tubing with an inner diameter slightly larger than the OD of your tubing and with an appropriate length for the PEEK tubing you wish to straighten. For instance, for 1/16" OD PEEK tubing with a length of 10", choose our U-825 tubing (stainless steel, 1/8" OD x 0.080" ID x 25 cm long, page 19. Slip your PEEK tubing into the stainless steel tubing. Place this "sleeved" tubing into an oven and bake at 425 °F (218 °C) for 30 minutes or 350 °F (177 °C) for 60 minutes. Allow the sleeved tubing to return to room temperature naturally (i.e., do not quench it with water). Once cooled, remove the PEEK tubing from the stainless steel sleeve and inspect it for straightness. If needed, repeat the process until the desired straightness is achieved.



SPECIFICATIONS & DETAILS

| Tubing OD | Tubing ID | OD Tolerance | ID Tolerance |
|--------------|---|--------------------|--------------------|
| PEEK TUBING | G SPECIFICATIONS | | |
| 1/16" | 25 μm | ±0.001" (25 µm) | ±0.0005" (12.5 μm) |
| 1/8" | All | ±0.003" (75 µm) | ±0.003" (75 μm) |
| CAPILLARY P | PEEK TUBING SPECIFICATIONS | | |
| 360 µm | All | ±0.0005" (12.5 μm) | ±0.0005" (12.5 μm) |
| 1/32" | All | ±0.0005" (12.5 μm) | ±0.0005" (12.5 μm) |
| FUSED SILICA | A TUBING, 360 μm OD | | |
| 360 µm | 20 μm (0.0008") | ±0.0004" (10 μm) | ±0.00008" (2 μm) |
| 360 µm | 50 μm (0.002") and 75 μm (0.003") | ±0.0004" (10 μm) | ±0.00012" (3 μm) |
| 360 µm | 100 μm (0.004") and 150 μm (0.006") | ±0.0004" (10 μm) | ±0.00016" (4 µm) |

PEEK Tubing (Cont.)

PEEK Tubing

| | ID | Color | Max. Pressure | |
|---|---|---|---|-----------------------|
| PEEK TUBING | | | | |
| 560 | 0.0025" (65 μm) ID x 5' (1.5 m) | Natural | 7,000 psi (483 bar) | е |
| 560L | 0.0025" (65 µm) ID x 50' (15 m) | Natural | 7,000 psi (483 bar) | e |
| 560XL | 0.0025" (65 µm) ID x 100' (30 m) | Natural | 7,000 psi (483 bar) | e |
| 560M | 0.0025" (65 µm) ID x 1,000" (304 m) | Natural | 7,000 psi (483 bar) | e |
| | | | | |
| 561 | 0.004" (0.10 mm) ID x 5' (1.5 m) | Black | 7,000 psi (483 bar) | е |
| 561L | 0.004" (0.10 mm) ID x 50' (15 m) | Black | 7,000 psi (483 bar) | е |
| 561XL | 0.004" (0.10 mm) ID x 100' (30 m) | Black | 7,000 psi (483 bar) | е |
| 561M | 0.004" (0.10 mm) ID x 1,000' (304 m) | Black | 7,000 psi (483 bar) | e |
| 535 | 0.005" (0.125 mm) ID x 5' (1.5 m) | Red | 7,000 psi (483 bar) | е |
| 535L | 0.005" (0.125 mm) ID x 50' (15 m) | Red | 7,000 psi (483 bar) | е |
| 535XL | 0.005" (0.125 mm) ID x 100' (30 m) | Red | 7,000 psi (483 bar) | e |
| 535M | | | | |
| | 0.005" (0.125 mm) ID x 1,000' (304 m) | Red | 7,000 psi (483 bar) | е |
| 536 | 0.007" (0.175 mm) ID x 5' (1.5 m) | Yellow | 7,000 psi (483 bar) | е |
| 536L | 0.007" (0.175 mm) ID x 50' (15 m) | Yellow | 7,000 psi (483 bar) | е |
| 536XL | 0.007" (0.175 mm) ID x 100' (30 m) | Yellow | 7,000 psi (483 bar) | е |
| 536M | 0.007" (0.175 mm) ID x 1,000' (304 m) | Yellow | 7,000 psi (483 bar) | е |
| 531 | 0.010" (0.25 mm) ID x 5' (1.5 m) | Natural | 7,000 psi (483 bar) | е |
| 531L | 0.010" (0.25 mm) ID x 50' (15 m) | Natural | 7,000 psi (483 bar) | e |
| | | | | |
| 531XL | 0.010" (0.25 mm) ID x 100' (30 m) | Natural | 7,000 psi (483 bar) | е |
| 531M | 0.010" (0.25 mm) ID x ID x 1,000' (304 m) | Natural | 7,000 psi (483 bar) | е |
| 31B | 0.010" (0.25 mm) ID x 5' (1.5 m) | Blue | 7,000 psi (483 bar) | е |
| 31BL | 0.010" (0.25 mm) ID x 50' (15 m) | Blue | 7,000 psi (483 bar) | е |
| 31BXL | 0.010" (0.25 mm) ID x 100' (30 m) | Blue | 7,000 psi (483 bar) | е |
| 31BM | 0.010" (0.25 mm) ID x 1,000' (304 m) | Blue | 7,000 psi (483 bar) | e |
| | | | | |
| 32 | 0.020" (0.50 mm) ID x 5' (1.5 m) | Orange | 6,000 psi (414 bar) | E |
| 32L | 0.020" (0.50 mm) ID x 50' (15 m) | Orange | 6,000 psi (414 bar) | e |
| 32XL | 0.020" (0.50 mm) ID x 100' (30 m) | Orange | 6,000 psi (414 bar) | е |
| 32M | 0.020" (0.50 mm) ID x 1,000' (304 m) | Orange | 6,000 psi (414 bar) | е |
| 33 | 0.030" (0.75 mm) ID x 5' (1.5 m) | Green | 4,000 psi (276 bar) | e |
| 33L | 0.030" (0.75 mm) ID x 50' (15 m) | Green | 4,000 psi (276 bar) | € |
| 33XL | 0.030" (0.75 mm) ID x 100' (30 m) | Green | 4,000 psi (276 bar) | |
| | | | | € |
| 33M | 0.030" (0.75 mm) ID x 1,000' (304 m) | Green | 4,000 psi (276 bar) | e |
| 38 | 0.040" (1.00 mm) ID x 5' (1.5 m) | Natural | 3,000 psi (207 bar) | € |
| 38L | 0.040" (1.00 mm) ID x 50' (15 m) | Natural | 3,000 psi (207 bar) | e |
| 38XL | 0.040" (1.00 mm) ID x 100' (30 m) | Natural | 3,000 psi (207 bar) | e |
| 38M | 0.040" (1.00 mm) ID x 1,000' (304 m) | Natural | 3,000 psi (207 bar) | e |
| EEK TUBING | | | 5,555 [25: (25: 32:) | - |
| | | N I | 4.000 : (07.1) | |
| 534 | 0.062" (1.55 mm) ID x 5' (1.5 m) | Natural | 4,000 psi (276 bar) | е |
| 44 | 0.080" (2.00 mm) ID x 5' (1.5 m) | Natural | 3,000 psi (207 bar) | 6 |
| APILLARY P | EEK TUBING, 360 μm OD | | | |
| 574 | 25 μm (0.001") ID x 5' (1.5 m) | Natural | 5,000 psi (345 bar) | е |
| 70 | 50 μm (0.002") ID x 5' (1.5 m) | Natural | 2,000 psi (138 bar) | е |
| 71 | 100 µm (0.004") ID x 5' (1.5 m) | Red | 2,000 psi (138 bar) | |
| | | | | 6 |
| 72 | 150 μm (0.006") ID x 5' (1.5 m) | Yellow | 2,000 psi (138 bar) | E |
| APILLARY P | EEK TUBING, 1/32" OD | | | |
| 76 | 0.005" (0.125 mm) ID x 5' (1.5 m) | Red | 5,000 psi (345 bar) | e |
| 76L | 0.005" (0.125 mm) ID x 50' (15 m) | Red | 5,000 psi (345 bar) | 6 |
| 76XL | 0.005" (0.125 mm) ID x 100' (30 m) | Red | 5,000 psi (345 bar) | 6 |
| 76M | 0.005" (0.125 mm) ID x 1,000' (304 m) | Red | 5,000 psi (345 bar) | |
| | | | | 6 |
| 77 | 0.007" (0.175 mm) ID x 5' (1.5 m) | Yellow | 5,000 psi (345 bar) | 6 |
| 77L | 0.007" (0.175 mm) ID x 50' (15 m) | Yellow | 5,000 psi (345 bar) | € |
| 77XL | 0.007" (0.175 mm) ID x 100' (30 m) | Yellow | 5,000 psi (345 bar) | 6 |
| 77M | 0.007" (0.175 mm) ID x 1,000' (304 m) | Yellow | 5,000 psi (345 bar) | 6 |
| 81 | 0.010" (0.25 mm) ID x 5' (1.5 m) | Blue | 5,000 psi (345 bar) | 6 |
| 81L | 0.010" (0.25 mm) ID x 50' (1.5 m) | Blue | 5,000 psi (345 bar) | |
| | | | | 6 |
| 81XL | 0.010" (0.25 mm) ID x 100' (30 m) | Blue | 5,000 psi (345 bar) | € |
| 81M | 0.010" (0.25 mm) ID x 1,000' (304 m) | Blue | 5,000 psi (345 bar) | € |
| 68 | 0.015" (0.40 mm) ID x 5' (1.5 m) | Natural | 4,000 psi (276 bar) | e |
| | 0.015" (0.40 mm) ID x 50' (15 m) | Natural | 4,000 psi (276 bar) | e |
| | | Natural | 4,000 psi (276 bar) | |
| 68L | 0.015" (0.40 mm) ID x 100' (30 m) | ivaturar | | 6 |
| 68KL | | | | |
| 68L 68XL 68M | 0.015" (0.40 mm) ID x 1,000' (304 m) | Natural | 4,000 psi (276 bar) | e |
| 668L 668XL 668M 669 | 0.015" (0.40 mm) ID x 1,000' (304 m) 0.020" (0.50 mm) ID x 5' (1.5 m) | Natural Orange | 4,000 psi (276 bar) 3,000 psi (207 bar) | 6 |
| 668L 668XL 668M 669 669L | 0.015" (0.40 mm) ID x 1,000' (304 m) 0.020" (0.50 mm) ID x 5' (1.5 m) 0.020" (0.50 mm) ID x 50' (15 m) | Natural Orange Orange | 4,000 psi (276 bar) 3,000 psi (207 bar) 3,000 psi (207 bar) | e e |
| 668L 668XL 668M 669 669L 669XL | 0.015" (0.40 mm) ID x 1,000' (304 m) 0.020" (0.50 mm) ID x 5' (1.5 m) | Natural Orange | 4,000 psi (276 bar) 3,000 psi (207 bar) | € € |
| 668L 668XL 668M 669 669L 669XL | 0.015" (0.40 mm) ID x 1,000' (304 m) 0.020" (0.50 mm) ID x 5' (1.5 m) 0.020" (0.50 mm) ID x 50' (15 m) | Natural Orange Orange | 4,000 psi (276 bar) 3,000 psi (207 bar) 3,000 psi (207 bar) | e e e |
| 668L 668XL 668M 669 669L 669XL 669M | 0.015" (0.40 mm) ID x 1,000' (304 m) 0.020" (0.50 mm) ID x 5' (1.5 m) 0.020" (0.50 mm) ID x 50' (15 m) 0.020" (0.50 mm) ID x 100' (30 m) 0.020" (0.50 mm) ID x 1,000' (304 m) | Natural Orange Orange Orange | 4,000 psi (276 bar) 3,000 psi (207 bar) 3,000 psi (207 bar) 3,000 psi (207 bar) | e e e |
| 668L 668XL 668M 669 669L 669XL 669M JSED SILICA | 0.015" (0.40 mm) ID x 1,000' (304 m) 0.020" (0.50 mm) ID x 5' (1.5 m) 0.020" (0.50 mm) ID x 50' (15 m) 0.020" (0.50 mm) ID x 100' (30 m) 0.020" (0.50 mm) ID x 1,000' (304 m) A TUBING, 360 µm OD | Natural Orange Orange Orange Orange | 4,000 psi (276 bar) 3,000 psi (207 bar) 3,000 psi (207 bar) 3,000 psi (207 bar) 3,000 psi (207 bar) | e e e |
| 568L 568XL 568M 569 569L 569XL 569M USED SILICA 5-120 | 0.015" (0.40 mm) ID x 1,000' (304 m) 0.020" (0.50 mm) ID x 5' (1.5 m) 0.020" (0.50 mm) ID x 50' (15 m) 0.020" (0.50 mm) ID x 100' (30 m) 0.020" (0.50 mm) ID x 1,000' (304 m) A TUBING, 360 µm OD 20 µm (0.0008") ID x 6.4' (2 m) | Natural Orange Orange Orange Orange Natural | 4,000 psi (276 bar) 3,000 psi (207 bar) 3,000 psi (207 bar) 3,000 psi (207 bar) 3,000 psi (207 bar) 10,000 psi (690 bar) | e e e e e |
| 568L 568XL 568M 569 569L 569XL 569M USED SILICA 5-120 | 0.015" (0.40 mm) ID x 1,000' (304 m) 0.020" (0.50 mm) ID x 5' (1.5 m) 0.020" (0.50 mm) ID x 50' (15 m) 0.020" (0.50 mm) ID x 100' (30 m) 0.020" (0.50 mm) ID x 1,000' (304 m) A TUBING, 360 µm OD 20 µm (0.0008") ID x 6.4' (2 m) 50 µm (0.002") ID x 6.4' (2 m) | Natural Orange Orange Orange Orange Natural Natural | 4,000 psi (276 bar) 3,000 psi (207 bar) 10,000 psi (690 bar) | 6 6 6 6 6 |
| 568L 568XL 568M 569 569XL 569XL 569M USED SILICA 5-120 5-150 | 0.015" (0.40 mm) ID x 1,000' (304 m) 0.020" (0.50 mm) ID x 5' (1.5 m) 0.020" (0.50 mm) ID x 50' (15 m) 0.020" (0.50 mm) ID x 100' (30 m) 0.020" (0.50 mm) ID x 1,000' (304 m) A TUBING, 360 µm OD 20 µm (0.0008") ID x 6.4' (2 m) 50 µm (0.002") ID x 6.4' (2 m) 75 µm (0.003") ID x 6.4' (2 m) | Natural Orange Orange Orange Orange Natural Natural | 4,000 psi (276 bar) 3,000 psi (207 bar) 10,000 psi (690 bar) 10,000 psi (690 bar) 10,000 psi (690 bar) | e e e e |
| 668L 668XL 668M 669 669L 669XL 669M USED SILICA 6-120 | 0.015" (0.40 mm) ID x 1,000' (304 m) 0.020" (0.50 mm) ID x 5' (1.5 m) 0.020" (0.50 mm) ID x 50' (15 m) 0.020" (0.50 mm) ID x 100' (30 m) 0.020" (0.50 mm) ID x 1,000' (304 m) A TUBING, 360 µm OD 20 µm (0.0008") ID x 6.4' (2 m) 50 µm (0.002") ID x 6.4' (2 m) | Natural Orange Orange Orange Orange Natural Natural | 4,000 psi (276 bar) 3,000 psi (207 bar) 10,000 psi (690 bar) | 6 6 6 6 6 |



Stainless Steel Tubing

- > Precut 316 stainless steel
- Available ODs include 1/32", 1/16", and 1/8"
- > Color-coded banding for easy identification of the inner diameter

IDEX Health & Science seamless, precut stainless steel tubing is designed to meet the exacting requirements of today's analyses. We machine cut and polish each end, deburr the inside and outside edges, and passivate the tubing (please see the passivation information on this page). Finally, we flush reagent-grade isopropanol through each piece.

Our thorough preparation and cleaning procedure guarantees tubing that is truly ready-to-use, with flat, burr-free ends and a clean finish. This care is important in achieving zero-dead-volume connections and good chromatographic results.

We offer a variety of precut lengths as well as longer lengths (5' and 25') of some sizes. Cutting the tubing disturbs and roughens the tubing's end surface, so we recommend using our precut tubing whenever possible. If you need to cut tubing to custom lengths, we suggest you then passivate the tubing.



- Our 1/32" OD tubing is designed for enhanced flexibility in high pressure applications.
- Standard 1/16" and 1/8" OD stainless steel tubing is suited for most analytical applications.



NOTE

PEEK polymer tubing can be used to replace stainless steel tubing in most liquid analytical systems. Unlike stainless steel tubing, PEEK tubing is biocompatible, flexible, and can easily be cut to desired lengths. See page 16. All Stainless Steel tubing longer than 1 m is coiled.

The Beauty of Precut Tubing









Tubing cut by a commercially available tubing cutter

File cut tubing

Stainless Steel Tubing Passivation

Stainless steel is naturally self-passivating, forming an oxidized layer on newly created surfaces. IDEX Health & Science takes extra steps to ensure the chemical resistance of our stainless steel tubing by manually passivating before and after the tubing is cut into specified lengths (except in a few cases where size is prohibitive). In the precut stage, the internal wall is acid passivated and flushed. After the tubing is cut, deburred and polished, it is completely submerged in an acid passivation bath and again flushed clean. The table below summarizes the manual passivation steps performed for each size of our stainless steel tubing:

| Tubing OD | Precut Passivation | Postcut Passivation |
|-----------|--------------------|--------------------------|
| 1/32" | All | All |
| 1/16" | All | All, ex. 25' lengths |
| 1/8" | None | All, ex. 3 & 5 m lengths |

Stainless Steel Tubing (Cont.)



SPECIFICATIONS & DETAILS

- Maximum Recommended Operating Temperature: 750 °F (399 °C).
- Rockwell Hardness (B): Maximum of 95.
- Meets ASTM A269 and A213.

| Tubing OD | OD Tolerance | Tubing ID | ID Tolerance |
|-----------|--------------------------------|------------------------------|--------------------------------|
| 1/32" | +0.002"/-0.000" (+50 μm/-0 μm) | All, except 0.004" (0.10 mm) | +0.000"/-0.002" (+0 μm/-50 μm) |
| 1/32" | +0.002"/-0.000" (+50 μm/-0 μm) | 0.004" (0.10 mm) | +0.002"/-0.000" (+50 μm/-0 μm) |
| 1/16" | +0.002"/-0.000" (+50 μm/-0 μm) | All | ±0.001" (25 μm) |
| 1/8" | ±0.003" (75 μm) | All | ±0.003" (75 μm) |





> PEEK polymer tubing is available in all of these sizes, starting on page 16.



Understanding the Maximum Pressure Value of Stainless Steel Tubing

Stainless steel is unique as a material. The Maximum Pressure value listed for each part number is the safe, continuous working pressure limit that IDEX Health & Science has assigned for the tubing. It reflects a safety margin before the tubing begins to "yield" which is well below the tubing's "burst" pressure. For more information, contact IDEX Health & Science or your authorized Distributor.

Stainless Steel Tubing

| Part No. | ID | Length | Color | Maximum Pressure | Qty. |
|-----------------|------------------|-------------|--------|------------------------|------|
| STAINLESS STEEL | ., 1/32" OD | | | | |
| U-1114 | 0.004" (0.10 mm) | 2" (5 cm) | Red | 19,300 psi (1,331 bar) | ea. |
| U-1115 | 0.004" (0.10 mm) | 4" (10 cm) | Red | 19,300 psi (1,331 bar) | ea. |
| U-1116 | 0.004" (0.10 mm) | 8" (20 cm) | Red | 19,300 psi (1,331 bar) | ea. |
| U-1117 | 0.004" (0.10 mm) | 12" (30 cm) | Red | 19,300 psi (1,331 bar) | ea. |
| U-1120 | 0.006" (0.15 mm) | 2" (5 cm) | Yellow | 19,300 psi (1,331 bar) | ea. |
| U-1122 | 0.006" (0.15 mm) | 8" (20 cm) | Yellow | 19,300 psi (1,331 bar) | ea. |
| U-1125 | 0.008" (0.20 mm) | 2" (5 cm) | Clear | 17,800 psi (1,227 bar) | ea. |
| U-1126 | 0.008" (0.20 mm) | 4" (10 cm) | Clear | 17,800 psi (1,227 bar) | ea. |
| U-1128 | 0.008" (0.20 mm) | 12" (30 cm) | Clear | 17,800 psi (1,227 bar) | ea. |
| U-1130 | 0.010" (0.25 mm) | 2" (5 cm) | Blue | 16,200 psi (1,117 bar) | ea. |
| U-1131 | 0.010" (0.25 mm) | 4" (10 cm) | Blue | 16,200 psi (1,117 bar) | ea. |
| U-1132 | 0.010" (0.25 mm) | 8" (20 cm) | Blue | 16,200 psi (1,117 bar) | ea. |
| U-1133 | 0.010" (0.25 mm) | 12" (30 cm) | Blue | 16,200 psi (1,117 bar) | ea. |
| U-1140 | 0.015" (0.40 mm) | 2" (5 cm) | Green | 12,300 psi (848 bar) | ea. |
| U-1141 | 0.015" (0.40 mm) | 4" (10 cm) | Green | 12,300 psi (848 bar) | ea. |
| U-1142 | 0.015" (0.40 mm) | 8" (20 cm) | Green | 12,300 psi (848 bar) | ea. |
| U-1143 | 0.015" (0.40 mm) | 12" (30 cm) | Green | 12,300 psi (848 bar) | ea. |
| U-1145 | 0.018" (0.45 mm) | 2" (5 cm) | Black | 10,000 psi (689 bar) | ea. |
| U-1146 | 0.018" (0.45 mm) | 4" (10 cm) | Black | 10,000 psi (689 bar) | ea. |
| U-1148 | 0.018" (0.45 mm) | 12" (30 cm) | Black | 10,000 psi (689 bar) | ea. |

Stainless Steel Tubing (Cont.)

| t No. | ID | Length | Color | Maximum Pressure | Qty. |
|------------------|-------------------|-------------------------|--------------|------------------------|------|
| INLESS STEEL, 1 | I/16" OD | | | | |
| 2 | 0.005" (0.125 mm) | 2" (5 cm) | Red | 21,600 psi (1,489 bar) | ea. |
| 53 | 0.005" (0.125 mm) | 4" (10 cm) | Red | 21,600 psi (1,489 bar) | ea. |
| 54 | 0.005" (0.125 mm) | 8" (20 cm) | Red | 21,600 psi (1,489 bar) | ea. |
| 55 | 0.005" (0.125 mm) | 12" (30 cm) | Red | 21,600 psi (1,489 bar) | ea. |
| 56 | 0.005" (0.125 mm) | 1.6' (0.5 m) | Red | 21,600 psi (1,489 bar) | ea. |
| 57 | 0.005" (0.125 mm) | 3.2' (1 m) | Red | 21,600 psi (1,489 bar) | ea. |
| 158 | 0.005" (0.125 mm) | 5' (1.5 m) | Red | 21,600 psi (1,489 bar) | ea. |
| 160 | 0.005" (0.125 mm) | 25' (7.6 m) | Red | 21,600 psi (1,489 bar) | ea. |
| 126 | 0.007" (0.175 mm) | 2" (5 cm) | Black | 20,900 psi (1,441 bar) | ea. |
| 127 | 0.007" (0.175 mm) | 4" (10 cm) | Black | 20,900 psi (1,441 bar) | ea. |
| 128 | 0.007" (0.175 mm) | 8" (20 cm) | Black | 20,900 psi (1,441 bar) | ea. |
| 129 | 0.007" (0.175 mm) | 12" (30 cm) | Black | 20,900 psi (1,441 bar) | ea. |
| 130 | 0.007" (0.175 mm) | 1.6' (0.5 m) | Black | 20,900 psi (1,441 bar) | ea. |
| 131 | 0.007" (0.175 mm) | 3.2' (1 m) | Black | 20,900 psi (1,441 bar) | ea. |
| 108 | 0.007" (0.175 mm) | 5' (1.5 m) | Black | 20,900 psi (1,441 bar) | ea. |
| 161 | 0.007" (0.175 mm) | 25' (7.6 m) | Black | 20,900 psi (1,441 bar) | ea. |
| 111 | 0.010" (0.25 mm) | 2" (5 cm) | Blue | 19,700 psi (1,358 bar) | ea. |
| 112 | 0.010" (0.25 mm) | 4" (10 cm) | Blue | 19,700 psi (1,358 bar) | ea. |
| 113 | 0.010" (0.25 mm) | 8" (20 cm) | Blue | 19,700 psi (1,358 bar) | ea. |
| 114 | 0.010" (0.25 mm) | 12" (30 cm) | Blue | 19,700 psi (1,358 bar) | ea. |
| 132 | 0.010" (0.25 mm) | 1.6' (0.5 m) | Blue | 19,700 psi (1,358 bar) | ea. |
| 133 | 0.010" (0.25 mm) | 3.2' (1 m) | Blue | 19,700 psi (1,358 bar) | ea. |
| 106 | 0.010" (0.25 mm) | 5' (1.5 m) | Blue | 19,700 psi (1,358 bar) | ea. |
| 162 | 0.010" (0.25 mm) | 25' (7.6 m) | Blue | 19,700 psi (1,358 bar) | ea. |
| 101 | 0.020" (0.5 mm) | 2" (5 cm) | Yellow | 15,800 psi (1,089 bar) | ea. |
| 102 | 0.020" (0.5 mm) | 4" (10 cm) | Yellow | 15,800 psi (1,089 bar) | ea. |
| 103 | 0.020" (0.5 mm) | 8" (20 cm) | Yellow | 15,800 psi (1,089 bar) | ea. |
| 104 | 0.020" (0.5 mm) | 12" (30 cm) | Yellow | 15,800 psi (1,089 bar) | ea. |
| 134 | 0.020" (0.5 mm) | 1.6' (0.5 m) | Yellow | 15,800 psi (1,089 bar) | ea. |
| 135 | 0.020" (0.5 mm) | 3.2' (1 m) | Yellow | 15,800 psi (1,089 bar) | ea. |
| 105 | 0.020" (0.5 mm) | 5′ (1.5 m) | Yellow | 15,800 psi (1,089 bar) | ea. |
| 163 | 0.020" (0.5 mm) | 25' (7.6 m) | Yellow | 15,800 psi (1,089 bar) | ea. |
| 115 | 0.030" (0.75 mm) | 2" (5 cm) | White | 12,000 psi (827 bar) | ea. |
| 116 | 0.030" (0.75 mm) | 4" (10 cm) | White | 12,000 psi (827 bar) | ea. |
| 117 | 0.030" (0.75 mm) | 8" (20 cm) | White | 12,000 psi (827 bar) | ea. |
| 118 | 0.030" (0.75 mm) | 12" (30 cm) | White | 12,000 psi (827 bar) | ea. |
| 136 | 0.030" (0.75 mm) | 1.6' (0.5 m) | White | 12,000 psi (827 bar) | ea. |
| 137 | 0.030" (0.75 mm) | 3.2' (1 m) | White | 12,000 psi (827 bar) | ea. |
| 107 | 0.030" (0.75 mm) | 5′ (1.5 m) | White | 12,000 psi (827 bar) | ea. |
| 164 | 0.030" (0.75 mm) | 25′ (7.6 m) | White | 12,000 psi (827 bar) | ea. |
| 138 | 0.040" (1.0 mm) | 2" (5 cm) | N/A | 8,100 psi (558 bar) | ea. |
| 139 | 0.040" (1.0 mm) | 4" (10 cm) | N/A N/A | 8,100 psi (558 bar) | |
| 140 | | | N/A N/A | | ea. |
| | 0.040" (1.0 mm) | 8" (20 cm) | | 8,100 psi (558 bar) | ea. |
| 141 | 0.040" (1.0 mm) | 12" (30 cm) | N/A | 8,100 psi (558 bar) | ea. |
| 142 | 0.040" (1.0 mm) | 1.6' (0.5 m) | N/A | 8,100 psi (558 bar) | ea. |
| 143 | 0.040" (1.0 mm) | 3.2' (1 m) | N/A | 8,100 psi (558 bar) | ea. |
| 144 | 0.040" (1.0 mm) | 5' (1.5 m) | N/A | 8,100 psi (558 bar) | ea. |
| 165 | 0.040" (1.0 mm) | 25' (7.6 m) | N/A | 8,100 psi (558 bar) | ea. |
| 145 | 0.046" (1.15 mm) | 2" (5 cm) | N/A | 5,800 psi (400 bar) | ea. |
| 146 | 0.046" (1.15 mm) | 4" (10 cm) | N/A | 5,800 psi (400 bar) | ea. |
| 147 | 0.046" (1.15 mm) | 8" (20 cm) | N/A | 5,800 psi (400 bar) | ea. |
| 148 | 0.046" (1.15 mm) | 12" (30 cm) | N/A | 5,800 psi (400 bar) | ea. |
| 149 | 0.046" (1.15 mm) | 1.6' (0.5 m) | N/A | 5,800 psi (400 bar) | ea. |
| 150 | 0.046" (1.15 mm) | 3.2′ (1 m) | N/A | 5,800 psi (400 bar) | ea. |
| 51 | 0.046" (1.15 mm) | 5′ (1.5 m) | N/A | 5,800 psi (400 bar) | ea. |
| AINLESS STEEL, 1 | /8" OD | | | | |
| 325 | 0.080" (2.0 mm) | 10" (25 cm) | N/A | 7,600 psi (524 bar) | ea. |
| 800 | 0.080" (2.0 mm) | 3.2' (1 m) | N/A | 7,600 psi (524 bar) | ea. |
| 803 | 0.080" (2.0 mm) | 9.8' (3 m) | N/A | 7,600 psi (524 bar) | ea. |
| 805 | 0.080" (2.0 mm) | 16' (5 m) | N/A | 7,600 psi (524 bar) | ea. |
| HI | ROMalytic +61(0 | 13 9762 2034 Australian | Distributors | | |
| | | , | Manufacurers | | |
| | ECHnois | | mtech.net.au | | |






PEEKsil Tubing

- > PEEK covered fused silica
- ▶ 1/32" and 1/16" outside diameters with a wide variety of inside diameters
- > Precut to numerous standard lengths

PEEKsil's sheathing is mechanically strong and has ideal characteristics for sealing with many styles of fittings. The fused silica core provides a consistent and rigid fluid pathway with very tight tolerances and industry-accepted chemical properties. Together, this makes PEEKsil tubing ideal for numerous applications. In fact, PEEKsil can be used as a direct replacement for conventional stainless steel or PEEK tubing in many analytical systems.

Like traditional fused silica tubing, PEEKsil has excellent chemical compatibility and extremely low adsorption characteristics, especially when compared with stainless steel.

Please Note: Do not cut this tubing. It should be used at its precut lengths because of permanent damage caused by conventional cutters.





SPECIFICATIONS & DETAILS

| Tubing OD | OD Tolerance | Tubing ID | ID Tolerance |
|-----------|------------------|--------------|------------------|
| | | 25 μm | ±0.00004" (1 μm) |
| 1/32" | ±0.0008" (20 μm) | 50–100 μm | ±0.00012" (3 μm) |
| 1/16" | ±0.0012" (30 μm) | 0.15-0.30 mm | ±0.0002" (5 μm) |



SPECIFICATIONS & DETAILS

Because PEEKsil tubing has fused silica tubing at its core, the pressure rating for this tubing is determined by the inner diameter of the tubing. The following chart highlights the Maximum Pressure values for this tubing, as determined by SGE International Pty., Ltd., the manufacturer of this tubing:

| Tubing ID | Maximum Pressure |
|------------|------------------------|
| 25 μm | 25,000 psi (1,723 bar) |
| 50 μm | 20,000 psi (1,379 bar) |
| 75–100 μm | 15,000 psi (1,034 bar) |
| 150–175 μm | 8,500 psi (586 bar) |
| 200–300 μm | 6,000 psi (414 bar) |

The pressure ratings provided are indicative of the performance capabilities of the tubing. The actual pressure limits achievable will depend upon the fittings used, the quality of the receiving port, and other factors. Contact IDEX Health & Science or your authorized Distributor for more information

PEEKsil[™] Tubing

| (10 cm) EKSIL TUBING, 1/32" OD | ID | Length | Color | Qty. |
|-----------------------------------|-----------------------------------|--------------------------|---------|--------------|
| 5 | 0.001" (25 μm) | 2" (5 cm) | Orange | 2-pk |
| 510 | 0.001" (25 μm) | 4" (10 cm) | Orange | 2-pk |
| 515 | 0.001" (25 μm) | 6" (15 cm) | Orange | 2-pk |
| 520 | 0.001" (25 μm) | 8" (20 cm) | Orange | 2-pk |
| 550 | 0.001" (25 μm) | 1.6' (50 cm) | Orange | 2-pk |
| 05 | 0.002" (50 μm) | 2" (5 cm) | Natural | 2-pk |
| 010 | 0.002" (50 μm) | 4" (10 cm) | Natural | 2-pk |
| 015 | 0.002" (50 μm) | 6" (15 cm) | Natural | 2-pk |
| 5020 | 0.002" (50 µm) | 8" (20 cm) | Natural | 2-pk |
| EKSIL TUBING, 1/32" OD | 0.002 (30 μπ) | 0 (20 cm) | Natarai | 2 μκ |
| 050 | 0.002" (50 µm) | 1.6' (50 cm) | Natural | 2-pk |
| 55 | 0.002' (30 μm) 0.003" (75 μm) | 2" (5 cm) | Black | 2-pk |
| 510 | | 4" (10 cm) | Black | |
| 515 | 0.003" (75 μm) | | | 2-pk |
| | 0.003" (75 μm) | 6" (15 cm) | Black | 2-pk |
| 520 | 0.003" (75 μm) | 8" (20 cm) | Black | 2-pk |
| 550 | 0.003" (75 μm) | 1.6' (50 cm) | Black | 2-pk |
| 005 | 0.004" (100 μm) | 2" (5 cm) | Red | 2-pk |
| 0010 | 0.004" (100 μm) | 4" (10 cm) | Red | 2-pk |
| 0015 | 0.004" (100 μm) | 6" (15 cm) | Red | 2-pk |
| 0020 | 0.004" (100 μm) | 8" (20 cm) | Red | 2-pk |
| 10050 | 0.004" (100 μm) | 1.6' (50 cm) | Red | 2-pk |
| 505 | 0.006" (150 µm) | 2" (5 cm) | Purple | 2-pk |
| 15010 | 0.006" (150 μm) | 4" (10 cm) | Purple | 2-pk |
| 15015 | 0.006" (150 μm) | 6" (15 cm) | Purple | 2-pk |
| 5020 | 0.006" (150 μm) | 8" (20 cm) | Purple | 2-pk |
| 5050 | 0.006" (150 μm) | 1.6' (50 cm) | Purple | 2-pk |
| EEKSIL TUBING, 1/16" OD | | | | |
| 255 | 0.001" (25 µm) | 2" (5 cm) | Orange | 5-pk |
| 2510 | 0.001" (25 μm) | 4" (10 cm) | Orange | 5-pk |
| 515 | 0.001" (25 μm) | 6" (15 cm) | Orange | 5-pk |
| 2520 | 0.001" (25 μm) | 8" (20 cm) | Orange | 5-pk |
| 2550 | 0.001" (25 µm) | 1.6' (50 cm) | Orange | 2-pk |
| 05 | 0.002" (50 μm) | 2" (5 cm) | Natural | 5-pk |
| 5010 | 0.002' (50 μm) | 4" (10 cm) | Natural | 5-pk |
| 5015 | 0.002' (50 μm) | 6" (15 cm) | Natural | 5-pk |
| 5020 | | 8" (20 cm) | Natural | |
| | 0.002" (50 μm) | | | 5-pk |
| 050 | 0.002" (50 μm) | 1.6' (50 cm) | Natural | 2-pk |
| ⁷⁵⁵ | 0.003" (75 μm) | 2" (5 cm) | Black | 5-pk |
| 7510 | 0.003" (75 μm) | 4" (10 cm) | Black | 5-pk |
| 7515 | 0.003" (75 μm) | 6" (15 cm) | Black | 5-pk |
| 7520 | 0.003" (75 μm) | 8" (20 cm) | Black | 5-pk |
| 7550 | 0.003" (75 μm) | 1.6' (50 cm) | Black | 2-pk |
| 1005 | 0.004" (100 μm) | 2" (5 cm) | Red | 5-pk |
| 10010 | 0.004" (100 μm) | 4" (10 cm) | Red | 5-pk |
| 10015 | 0.004" (100 µm) | 6" (15 cm) | Red | 5-pk |
| 10020 | 0.004" (100 μm) | 8" (20 cm) | Red | 5-pk |
| 10050 | 0.004" (100 μm) | 1.6' (50 cm) | Red | 2-pk |
| 505 | 0.006" (150 μm) | 2" (5 cm) | Purple | 5-pk |
| 5010 | 0.006" (150 μm) | 4" (10 cm) | Purple | 5-pk |
| 5015 | 0.006" (150 μm) | 6" (15 cm) | Purple | 5-pk |
| 5020 | 0.006" (150 μm) | 8" (20 cm) | Purple | 5-pk |
| 5050 | 0.006" (150 µm) | 1.6' (50 cm) | Purple | 2-pk |
| 17515 | 0.007" (175 µm) | 6" (15 cm) | Yellow | 5-pk |
| 17520 | 0.007" (175 µm) | 8" (20 cm) | Yellow | 5-pk |
| 17550 | 0.007" (175 µm) | 1.6' (50 cm) | Yellow | 2-pk |
| 2005 | 0.007 (173 µm) 0.008" (200 µm) | 2" (5 cm) | Blue | 2-рк 5-рк |
| 0015 | 0.008" (200 µm) | 6" (15 cm) | Blue | |
| | | | | 5-pk |
| 20020 | 0.008" (200 µm) | 8" (20 cm) | Blue | 5-pk |
| 20050 | 0.008" (200 μm) | 1.6' (50 cm) | Blue | 2-pk |
| 3005 | 0.012" (300 μm) | 2" (5 cm) | Gray | 5-pk |
| 30010 | 0.012" (300 μm) | 4" (10 cm) | Gray | 5-pk |
| 30015 | 0.012" (300 μm) | 6" (15 cm) | Gray | 5-pk |
| 30020 | 0.012" (300 μm) | 8" (20 cm) | Gray | 5-pk |
| 30050 | 0.012" (300 um) | 1.6′ (50 cm) | Gray | 2-pk |
| HPON | Alytic +61(0)3 9762 2034 | Australian Distributors | | |
| | TO 11013 3702 20.14 | Importers & Manufacurers | | |

PFA Tubing



- ▶ 1/16" and 1/8" ODs available
- Excellent solvent resistance and low gas permeability
- ➤ Constructed with genuine Teflon[™] PFA resin

PFA (perfluoroalkoxyalkane) tubing offers excellent solvent resistance (virtually identical to FEP and PTFE) while adding several advantages. These include smoother surface texture, higher continuous service temperature and superior polymer purity. The recommended maximum operating temperature for our PFA tubing is 80 °C.



-) 360 μ m, 1/16", 1/8", 3/16", and 1/4" outside diameters available
- > PFA HP and PFA HP Plus Grades available
- > Virtually contaminant free
- Constructed with genuine Teflon™ High Purity PFA resin

PFA High Purity (HP) tubing offers all of the benefits of standard PFA tubing, with the additional benefit of being manufactured from a premium grade of PFA that is one of the most contaminant-free polymers available. In PFA HP, we offer tubing with the following outer diameters: 1/16", 1/8", 3/16", and 1/4".

PFA High Purity (HP) Plus tubing carries all of the benefits of PFA HP tubing, with the additional benefits of increased ability to withstand repeated flexing; improved resistance to stress cracking when exposed to aggressive fluorosurfactants; and smoother, clearer walls. In PFA HP Plus, we offer tubing with the following outer diameters: 360 μ m, 1/16", and 1/8".

(Please Note: Due to the physical nature of the 360 µm OD tubing, we recommend using our A-350 Polymer Tubing Cutter from page 28 when cutting this tubing to length. Additionally, extra care should be taken to ensure fittings are not overtightened and to ensure the tubing is not stretched once secured in place, to ensure the dimensional stability of the tubing.)



SPECIFICATIONS & DETAILS

| Tubing OD | OD Tolerances | Tubing ID | ID Tolerance | | | |
|-----------------|--|-----------|--------------------|--|--|--|
| PFA TUBING SPE | PFA TUBING SPECIFICATIONS | | | | | |
| 1/16" | ±0.001" (25 μm) | All | ±0.001" (25 μm) | | | |
| 1/8" | ±0.003" (75 μm) | All | ±0.003" (75 μm) | | | |
| HIGH PURITY PF | HIGH PURITY PFA TUBING SPECIFICATIONS | | | | | |
| 1/16" | ±0.001" (25 μm) | All | ±0.001" (25 μm) | | | |
| 1/8" | ±0.003" (75 μm) | All | ±0.003" (75 μm) | | | |
| 3/16" | ±0.003" (75 μm) | All | ±0.003" (75 μm) | | | |
| 1/4" | ±0.004" (100 μm) | All | ±0.004" (100 μm) | | | |
| 360 μM OD PFA I | 360 μM OD PFA HP TUBING SPECIFICATIONS | | | | | |
| 360 µm | ±0.0005" (12.5 μm) | All | ±0.0005" (12.5 μm) | | | |

FLUIDICS



| Part No. | ID | Length | Color | Max. Pressure | Qty. |
|-------------------|------------------|------------------|---|---------------------|--------------|
| PFA TUBING, 1/16' | " OD | | | | |
| 1500 | 0.020"(0.50 mm) | 5' (1.5 m) | Natural | 2,000 psi (138 bar) | ea. |
| 1512L | 0.020"(0.50 mm) | 50' (15 m) | Natural | 2,000 psi (138 bar) | ea. |
| 1512M | 0.020"(0.50 mm) | 1,000' (304 m) | Natural | 2,000 psi (138 bar) | ea. |
| 1502 | 0.030" (0.75 mm) | 5' (1.5 m) | Natural | 1,000 psi (69 bar) | ea. |
| 1514L | 0.030" (0.75 mm) | 50' (15 m) | Natural | 1,000 psi (69 bar) | ea. |
| 1514M | 0.030" (0.75 mm) | 1,000' (304 m) | Natural | 1,000 psi (69 bar) | ea. |
| 1503 | 0.040" (1.0 mm) | 5' (1.5 m) | Natural | 500 psi (34 bar) | ea. |
| 1507L | 0.040" (1.0 mm) | 50' (15 m) | Natural | 500 psi (34 bar) | ea. |
| 1507M | 0.040" (1.0 mm) | 1,000' (304 m) | Natural | 500 psi (34 bar) | ea. |
| PFA TUBING, 1/8" | OD | | | | |
| 1509-5 | 0.062" (1.55 mm) | 5' (1.5 m) | Natural | 500 psi (34 bar) | ea. |
| 1509L | 0.062" (1.55 mm) | 50' (15 m) | Natural | 500 psi (34 bar) | ea. |
| PFA HP TUBING, 1 | | | | | |
| 1622-5 | 0.020"(0.50 mm) | 5' (1.5 m) | Natural | 2,000 psi (138 bar) | ea. |
| 1622L | 0.020"(0.50 mm) | 50' (15 m) | Natural | 2,000 psi (138 bar) | ea. |
| 1622M | 0.020"(0.50 mm) | 1,000' (304 m) | Natural | 2,000 psi (138 bar) | ea. |
| 1632-5 | 0.030" (0.75 mm) | 5' (1.5 m) | Natural | 1,000 psi (69 bar) | ea. |
| 1632L | 0.030" (0.75 mm) | 50' (15 m) | Natural | 1,000 psi (69 bar) | ea. |
| 1632M | 0.030" (0.75 mm) | 1,000' (304 m) | Natural | 1,000 psi (69 bar) | ea. |
| PFA HP TUBING, 1 | | | | | |
| 1641-5 | 0.062" (1.55 mm) | 5' (1.5 m) | Natural | 500 psi (34 bar) | ea. |
| 1641L | 0.062" (1.55 mm) | 50' (15 m) | Natural | 500 psi (34 bar) | ea. |
| PFA HP PLUS TUBI | NG, 1/16" OD | | | | |
| 1902-5 | 0.010 (0.25 mm) | 5' (1.5 m) | Natural | 3,000 psi (207 bar) | ea. |
| 1902L | 0.010 (0.25 mm) | 50' (15 m) | Natural | 3,000 psi (207 bar) | ea. |
| 1902M | 0.010 (0.25 mm) | 1,000' (304 m) | Natural | 3,000 psi (207 bar) | ea. |
| 1907-5 | 0.020"(0.50 mm) | 5' (1.5 m) | Natural | 2,000 psi (138 bar) | ea. |
| 1907L | 0.020"(0.50 mm) | 50' (15 m) | Natural | 2,000 psi (138 bar) | ea. |
| 1907M | 0.020"(0.50 mm) | 1,000' (304 m) | Natural | 2,000 psi (138 bar) | ea. |
| 1912-5 | 0.030" (0.75 mm) | 5' (1.5 m) | Natural | 1,000 psi (69 bar) | ea. |
| 1912L | 0.030" (0.75 mm) | 50' (15 m) | Natural | 1,000 psi (69 bar) | ea. |
| 1912M | 0.030" (0.75 mm) | 1,000' (304 m) | Natural | 1,000 psi (69 bar) | ea. |
| PFA HP PLUS TUBI | NG, 1/8" OD | | | | |
| 1921-5 | 0.062" (1.55 mm) | 5' (1.5 m) | Natural | 500 psi (34 bar) | ea. |
| 1921L | 0.062" (1.55 mm) | 50' (15 m) | Natural | 500 psi (34 bar) | ea. |
| H | ROMalytic | 61(0)3 9762 2034 | Australian Distributors Importers & Manufacurers www.chromtech.net.au | | |
| | | | 03 9762 2034 in AUSTRALIA | 2 | 5 (of 172) 2 |







FEP Tubing

- > Great for moderate-to-low pressure applicationse
- **1/16"**, 1/8", 3/16", 1/4", or 5/16" outside diameter available
- ▶ 1 mm, 2 mm, or 3 mm outside diameter available
- Maximum continuous use temperature: 50 °C
- > Constructed with genuine Teflon™ FEP resin

With virtually identical chemical resistance to PFA at a lower price, FEP tubing is great for general, low pressure applications. Compared to PTFE, FEP (fluorinated ethylene propylene) tubing is held to tighter tolerances and has lower gas permeability (see material properties on our website: www.idex-hs.com).

 $\label{eq:much of our FEP Tubing - even the color-tinted options - is translucent, making it$ possible to watch fluid flow. Using different colored tubing can help identify transfer lines in multisolvent systems. Color coding also allows easy identification of the tubing thru-hole size. Black FEP tubing is available for light-sensitive applications (such as enzymatic and chemi-luminescent reactions) and entering/exiting flow cells.

SPECIFICATIONS & DETAILS

| Tubing Size | OD Tolerances | ID Tolerances |
|-------------|-------------------|-------------------|
| 1/16" OD | ±0.001" (25 μm) | ±0.001" (25 μm) |
| 1/8" OD | ±0.003" (75 μm) | ±0.003" (75 μm) |
| 3/16" OD | ±0.004" (0.10 mm) | ±0.004" (0.10 mm) |
| 5/16" OD | ±0.004" (0.10 mm) | ±0.004" (0.10 mm) |
| 1 mm OD | ±0.001" (25 μm) | ±0.001" (25 μm) |
| 2 mm OD | ±0.003" (75 μm) | ±0.003" (75 μm) |
| 3 mm OD | ±0.003" (75 μm) | ±0.003" (75 μm) |

| Part No. | ID | Length | Color | Max. Pressure | Qty |
|---------------|------------------|----------------|---------|---------------------|-----|
| FEP TUBING, 1 | /16" OD | | | | |
| 1527-5 | 0.010" (0.25 mm) | 5' (1.5 m) | Natural | 3,000 psi (207 bar) | ea. |
| 1527L | 0.010" (0.25 mm) | 50' (15 m) | Natural | 3,000 psi (207 bar) | ea. |
| 1527XL | 0.010" (0.25 mm) | 100' (30 m) | Natural | 3,000 psi (207 bar) | ea. |
| 1527M | 0.010" (0.25 mm) | 1,000' (304 m) | Natural | 3,000 psi (207 bar) | ea. |
| 1548-5 | 0.020" (0.50 mm) | 5' (1.5 m) | Natural | 2,000 psi (138 bar) | ea. |
| 1548L | 0.020" (0.50 mm) | 50' (15 m) | Natural | 2,000 psi (138 bar) | ea. |
| 1548XL | 0.020" (0.50 mm) | 100' (30 m) | Natural | 2,000 psi (138 bar) | ea. |
| 1548M | 0.020" (0.50 mm) | 1,000' (304 m) | Natural | 2,000 psi (138 bar) | ea. |
| 1520-5 | 0.030" (0.75 mm) | 5' (1.5 m) | Natural | 1,000 psi (69 bar) | ea. |
| 1520L | 0.030" (0.75 mm) | 50' (15 m) | Natural | 1,000 psi (69 bar) | ea. |
| 1520XL | 0.030" (0.75 mm) | 100' (30 m) | Natural | 1,000 psi (69 bar) | ea. |
| 1520M | 0.030" (0.75 mm) | 1,000' (304 m) | Natural | 1,000 psi (69 bar) | ea. |
| FEP TUBING, 1 | /8" OD | | | | |
| 1521-5 | 0.062" (1.55 mm) | 5' (1.5 m) | Natural | 500 psi (34 bar) | ea. |
| 1521L | 0.062" (1.55 mm) | 50' (15 m) | Natural | 500 psi (34 bar) | ea. |
| 1521XL | 0.062" (1.55 mm) | 100' (30 m) | Natural | 500 psi (34 bar) | ea. |
| FEP TUBING, 3 | /16" OD | | | | |
| 1524L | 0.125" (3.2 mm) | 50' (15 m) | Natural | 500 psi (34 bar) | ea. |
| 1524XL | 0.125" (3.2 mm) | 100' (30 m) | Natural | 500 psi (34 bar) | ea. |
| FEP TUBING, 1 | /4" OD | | | | |
| 1650L | 0.188" (4.8 mm) | 50' (15 m) | Natural | 250 psi (17 bar) | ea. |
| 1650XL | 0.188" (4.8 mm) | 100' (30 m) | Natural | 250 psi (17 bar) | ea. |
| FEP TUBING, 1 | .0 MM OD | | | | |
| 1671L | 0.020" (0.50 mm) | 50' (15 m) | Natural | 500 psi (34 bar) | ea. |
| 1671XL | 0.020" (0.50 mm) | 100' (30 m) | Natural | 500 psi (34 bar) | ea. |
| FEP TUBING, 2 | .0 MM OD | | | | |
| 1673L | 0.40" (1.0 mm) | 50' (15 m) | Natural | 500 psi (34 bar) | ea. |
| 1673XL | 0.40" (1.0 mm) | 100' (30 m) | Natural | 500 psi (34 bar) | ea. |
| FEP TUBING, 3 | .0 MM OD | | | | |
| 1677L | 0.080" (2.0 mm) | 50' (15 m) | Natural | 500 psi (34 bar) | ea. |
| 1677XL | 0.080" (2.0 mm) | 100' (30 m) | Natural | 500 psi (34 bar) | ea. |

FLUIDICS



ETFE Tubing

- > Excellent chemical resistance
- Constructed with genuine Tefzel[™] resin
- > Holds pressure up to 4,000 psi (276 bar)
- ▶ 1/16" or 1/8" outside diameter available
- Maximum continuous operating temperature: 80 °C

ETFE (ethylene-tetrafluoroethylene) tubing is an excellent fluoropolymer product that offers several benefits over tubing manufactured from PTFE, FEP, or PFA. These benefits include enhanced pressure holding capabilities, increased mechanical stability and lower gas permeability.



APPLICATION NOTE

ETFE tubing is an ideal choice for the fluid pathway between the vacuum degasser and the system's pump. Its low gas permeability will help ensure the mobile phase solvents do not regas while in transit.



RELATED PRODUCTS

Other tubing materials and dimensions may be available. Please contact IDEX Health & Science or your local representative directly.



SPECIFICATIONS & DETAILS

| Tubing OD | Tubing ID | OD/ID Tolerances |
|-----------|--|------------------|
| 1/16" OD | 0.010" (0.25 mm), 0.020" (0.50 mm), 0.030" (0.75 mm) | ±0.001" (25 μm) |
| 1/16" OD | 0.040" (1.0 mm) | ±0.002" (50 μm) |
| 1/8" OD | All | ±0.003" (75 μm) |

| Part No. | ID | Length | Color | Max. Pressure | Qty. |
|----------------|------------------|----------------|---------|---------------------|------|
| ETFE TUBING, 1 | /16" OD | | | | |
| 1529 | 0.010 (0.25 mm) | 5' (1.5 m) | Natural | 4,000 psi (276 bar) | ea. |
| 1529L | 0.010 (0.25 mm) | 50' (15 m) | Natural | 4,000 psi (276 bar) | ea. |
| 1529XL | 0.010 (0.25 mm) | 100' (30 m) | Natural | 4,000 psi (276 bar) | ea. |
| 1529M | 0.010 (0.25 mm) | 1,000' (304 m) | Natural | 4,000 psi (276 bar) | ea. |
| 1516 | 0.020"(0.50 mm) | 5' (1.5 m) | Natural | 3,000 psi (207 bar) | ea. |
| 1516L | 0.020"(0.50 mm) | 50' (15 m) | Natural | 3,000 psi (207 bar) | ea. |
| 1516XL | 0.020"(0.50 mm) | 100' (30 m) | Natural | 3,000 psi (207 bar) | ea. |
| 1516M | 0.020"(0.50 mm) | 1,000' (304 m) | Natural | 3,000 psi (207 bar) | ea. |
| 1528 | 0.030" (0.75 mm) | 5' (1.5 m) | Natural | 2,000 psi (138 bar) | ea. |
| 1528L | 0.030" (0.75 mm) | 50' (15 m) | Natural | 2,000 psi (138 bar) | ea. |
| 1528XL | 0.030" (0.75 mm) | 100' (30 m) | Natural | 2,000 psi (138 bar) | ea. |
| 1528M | 0.030" (0.75 mm) | 1,000' (304 m) | Natural | 2,000 psi (138 bar) | ea. |
| 1517 | 0.040" (1.0 mm) | 5' (1.5 m) | Natural | 500 psi (34 bar) | ea. |
| 1517L | 0.040" (1.0 mm) | 50' (15 m) | Natural | 500 psi (34 bar) | ea. |
| 1517XL | 0.040" (1.0 mm) | 100' (30 m) | Natural | 500 psi (34 bar) | ea. |
| 1517M | 0.040" (1.0 mm) | 1,000' (304 m) | Natural | 500 psi (34 bar) | ea. |
| ETFE TUBING, 1 | /8" OD | | | | |
| 1530 | 0.062" (1.55 mm) | 5' (1.5 m) | Natural | 1,000 psi (69 bar) | ea. |
| 1530L | 0.062" (1.55 mm) | 50' (15 m) | Natural | 1,000 psi (69 bar) | ea. |
| 1530XL | 0.062" (1.55 mm) | 100' (30 m) | Natural | 1,000 psi (69 bar) | ea. |
| 1648 | 0.093" (2.4 mm) | 5' (1.5 m) | Natural | 500 psi (34 bar) | ea. |
| 1648L | 0.093" (2.4 mm) | 50' (15 m) | Natural | 500 psi (34 bar) | ea. |
| 1648XL | 0.093" (2.4 mm) | 100' (30 m) | Natural | 500 psi (34 bar) | ea. |



Tools



Fused Silica Tubing Cutters

We offer a precision cutter for fused silica tubing — SGT's Shortix™ Cutter (FS-315). This cutter ensures clean, trouble-free cutting of fused silica tubing, providing better cuts than any other product on the market. It also includes a built-in magnifying glass to examine the cut tubing ends. Order the FS-315-02 Maintenance Kit, as needed, to replace a worn or damaged cutting wheel.

When using traditional fused silica tubing cutters, only a small part of the tubing wall is scratched, then the tubing is snapped or pulled in two, often resulting in a jagged, uneven cut. With a Shortix Cutter, a clean cut is made every time, regardless of skill or experience, as the cut is made by rotating a diamond blade around the entire circumference of the tubing.

Please Note: The FS-315 Fused Silica Tubing Cutters are designed to cut only tubing with ODs of 350 μ m-780 μ m and IDs of 100 μ m-350 μ m.



Polymer Tubing Cutters

For 1/16", 1/8", 3/16", 1/4", and 5/16" OD tubing

A flat, 90°, burr-free cut is difficult to obtain with most commercial polymer tubing cutters. Our experts have designed several tubing cutters specifically to cut polymer tubing. This line of tubing cutters includes a standard cutter for 1/16" and 1/8" OD tubing (A-327), and another for large bore tubing (A-329). Each has guide holes to ensure precise cutting. These cutters are durable, reliable, and easy to operate. Five replacement blades are included with each tool.



Capillary Polymer Tubing Cutters

Our A-350 Cutter is designed to cut capillary-sized polymer tubing. The cutter makes clean, perpendicular cuts without collapsing thin capillary walls. A set of ten tubing sleeves, required for cutting, are included with each cutter, along with five replacement blades. The included tubing sleeves are for cutting 360 µm OD polymer capillary tubing. Alternative sleeves are available for cutting 1/32" OD tubing. All tubing sleeves are 2" long and constructed with genuine Teflon™ FEP resin.



NOTE

- The A-350 Capillary Polymer Tubing Cutter can be used to cut tubing OD sizes other than 360 µm and 1/32". Simply use the proper NanoTight™ Tubing Sleeve found on page 52. Please note, however, that these sleeves are shorter than those listed on this page, and therefore will last through fewer cuts.
- Our tubing cutters are material specific: the A-327, A-329, and A-350 should only be used to cut <u>polymer</u> tubing, where as the FS-315 should only be used to cut fused silica tubing.
- Qty. Part No. Description FS-315 Shortix Fused Silica Tubing Cutter ea Capillary Polymer Tubing Cutter* for 360 μ m-1/32" OD tubing Includes (1) F-262x 10-pack of sleeves and (1) M-438-03 wrench A-350 F-262x Replacement Sleeves for A-350, 0.0155" ID, Green, for cutting 360 μm OD tubing 10-pk A-327 Standard Polymer Tubing Cutter* for 1/16" and 1/8" OD tubing ea Large Bore Polymer Tubing Cutter* for 3/16" – 5/16" OD tubing A-329 ea. Replacement Blades for A-350, A-327, and A-329 5-pk * Includes (1) A-328 5-pack of replacement blades







FITTINGS

We offer a wide and diverse selection of fittings to meet your system requirements. A "fitting" refers to a complete product ready to assemble and connect tubing into a part. This could be a one-piece connector or a nut and ferrule packaged together. A "Nut" indicates the male or female threaded product sold separately, and a "Ferrule" is sold separately when indicated in the description. For your convenience we ship most Fittings and Ferrules in 10-packs. We are dedicated to providing the most reliable, proven products on the market. We have implemented more stringent testing protocols and a generous safety margin to our ratings to ensure your safety.

- 32 CONED FITTINGS
- 39 FLAT-BOTTOM FITTINGS
- 48 MARVELXTM

- 50 FITTINGS TOOLS
- 56 LARGE BORE FITTINGS
- 57 VHP FITTINGS
- 61 FITTINGS KITS

FLUIDICS



All testing is performed with water at room temperature unless otherwise specified. Please contact IDEX Health & Science directly for further details. Results may vary depending on the material of the receiving port and tubing, actual tubing diameters (with stated tolerances), temperature and solvents used. If a pressure range is listed for a product's specification, the pressure rating depends on the tubing material used. The lower end of the range will represent testing performed on softer tubing such as FEP, and the higher end of the range will represent testing performed on harder tubing such as Stainless Steel. For more detail, please see the product specification sheets on our website, www.idex-hs.com, or contact us directly.



What Threads Do I Have?

Hold your fitting over the thread silhouettes below to identify the threads.

| U.S. CUSTOMARY THREADS | |
|------------------------|---|
| 6-40 | *************************************** |
| 6-32 | *************************************** |
| 10-32 | |
| 1/4-28 | |
| 5/16-24 | |
| 1/2-20 | |
| METRIC THREADS | |
| M4 x 0.7 | |
| M6 x 1 | |



One-Piece Fingertight Fittings

- > The original One-Piece Fingertight Fitting
- All polymer construction
- Versions available for 1/16", 1/32" or 1/8" OD tubing
- > Convenient and easy to use
- Reusable




NOTE

- Some of the fittings on this page are available in additional colors. Please contact your distributor or us for more information.
- ➤ Fingertight is generally equal to 3–4 in-lbs (0.34–0.45 N·m).

Our One-Piece Fingertight Fittings provide convenience and ease of use because the ferrule will not stick in a receiving port and the fitting is more easily found if dropped. The fittings for 1/16" OD tubing and 10-32 coned ports are available in a variety of colors, materials and lengths to suit virtually every application.

Beyond the standard 10-32 fittings, also featured in this product family are specialty fittings for specific applications. Our M-645 Fitting is a direct replacement for the 6-40 threaded VICI® (Valco) fitting. The P-100 can be used in 1/4-28 coned ports for 1/8" OD tubing including some of the inlet filters starting on page 102.

RheFlex® One-Piece Fittings are included in many of the manual valves, starting on page 121. The One-Piece RheFlex M4 Fittings, for use with our MX Nano-Scale Modules, are listed on page 61.



| Part No. | Description | Port | Pressure Rating | Head Style | Material | Qty. |
|-----------|---|--------------|--|----------------|------------------------------|-------|
| ONE-PIECE | FINGERTIGHT FITTINGS | | | | | |
| 6000-282 | Fingertight Fitting for 1/16" OD Tubing | 10-32 Coned | 5,000 psi (345 bar) | ChromTRAC knob | PEEK, Natural | 10-pk |
| F-100x | Fingertight Fitting for 1/16" OD Tubing | 10-32 Coned | 4,000 psi (276 bar) | Diamond Knurl | PCTFE, Red | 10-pk |
| F-120x | Fingertight Fitting for 1/16" OD Tubing | 10-32 Coned | 5,000 psi (345 bar) | Standard Knurl | PEEK, Natural | 10-pk |
| F-130x | Fingertight Fitting for 1/16" OD Tubing, Long | 10-32 Coned | 5,000 psi (345 bar) | Standard Knurl | PEEK, Natural | 10-pk |
| M-645x | Fingertight Fitting for 1/32" OD Tubing | 6-40 Coned | 1,750-3,250 psi (121-224 bar) | Headless Knurl | PEEK, Natural/PCTFE, Natural | 10-pk |
| P-100 | Finaertiaht Fittina for 1/8" OD Tubina | 1/4-28 Coned | 1.000 psi (69 bar) | Diamond Knurl | PCTFE, Natural | ea. |
| | HROMalytic +61(0)39 | /bZ ZU34 | Australian Distributors Importers & Manufacurers | | | |



Stainless Steel Fittings

1 NOTE

- Do not use metal fittings in plastic ports, as this can damage the port.
- The recommended torque to tighten these fittings is 20 in-lbs (2.25 N·m).

These 316 Stainless Steel Fittings are rated to 20,000 psi (1,380 bar) when wrench tightened. Choose IDEX Health & Science standard fittings, or select from the Common Valve Fittings or other manufacturer-compatible offerings.

Standard Stainless Steel Fittings



Common Valve Fittings



Beckman® Compatible Fittings



VICI® (Valco) Compatible Fittings



Website NEW: www.chromalytic.net.au E-mail: info@chromtech.net.au Tel: 03 9762 2034 . . . in AUSTRALIA

| Part No. | Description | Port | Pressure Rating | Head Style | Material | Qty. |
|----------|--|----------------------|------------------------|------------|----------|-----------|
| STANDAR | D STAINLESS STEEL FITTINGS | | | | | |
| C-235x | Nut for 1/8" OD Tubing | 1/4-28 Coned | 20,000 psi (1,380 bar) | 5/16" Hex | SST | 10-pk |
| C-236x | Ferrule for 1/8" OD Tubing | 1/4-28 Coned | 20,000 psi (1,380 bar) | _ | SST | 10-pk |
| U-400x | Nut for 1/16" OD Tubing | 10-32 Coned | 20,000 psi (1,380 bar) | 1/4" Hex | SST | 10-pk |
| U-401x | Ferrule for 1/16" OD Tubing | 10-32 or M6 Coned | 20,000 psi (1,380 bar) | _ | SST | 10-pk |
| U-450x | Nut for 1/16" OD Tubing | M6 Coned | 20,000 psi (1,380 bar) | 5/16" Hex | SST | 10-pk |
| COMMON | VALVE FITTINGS | | | | | |
| 6000-082 | Fitting for 1/8" OD Tubing | 5/16-24 Coned | 20,000 psi (1,380 bar) | 5/16" Hex | SST | ea. |
| 6000-083 | Ferrule for 1/8" OD Tubing | 5/16-24 Coned | 20,000 psi (1,380 bar) | _ | SST | 5-pk |
| 6000-209 | Fitting for 1/16" OD Tubing | 10-32 Coned | 20,000 psi (1,380 bar) | 1/4" Hex | SST | 10-pk |
| 6000-210 | Ferrule for 1/16" OD Tubing | 10-32 Coned | 20,000 psi (1,380 bar) | _ | SST | 10-pk |
| 6000-211 | Fitting for 1/16" OD Tubing, Long | 10-32 Coned | 20,000 psi (1,380 bar) | 1/4" Hex | SST | 10-pk |
| 6000-262 | Fitting for 1/16" OD Tubing, Extra Long | 10-32 Coned | 20,000 psi (1,380 bar) | 1/4" Hex | SST | 10-pk |
| MANUFA | CTURER COMPATIBLE FITTINGS | | | | | |
| U-320x | Nuts for 1/16" OD Tubing, Valco/VICI Compatible | 10-32 Coned | 20,000 psi (1,380 bar) | 1/4" Hex | SST | 10-pk |
| U-321x | Ferrule for 1/16" OD Tubing, Valco/VICI Compatible | 10-32 Coned | 20,000 psi (1,380 bar) | _ | SST | 10-pk |
| U-410x | Nuts for 1/16" OD Tubina. Waters Compatible | 10-32 Coned | 20.000 psi (1,380 bar) | 5/16" Hex | SST | 10-pk |
| | HROMalytic +61(0)3 9762 2034 | importers & manuracu | irers | | 22 /- | £ 472\ 2(|



MicroTight Fittings

- Comprehensive Fitting System for Connecting Capillary Tubing
- Made from PEEK Polymer

MicroTight® One-Piece Fittings are designed for use with the NanoPort™ and MicroTight Unions, Adapters and Inline MicroFilters. Specifically made for 360 µm OD tubing, 1/32" OD tubing, or our MicroTight Tubing Sleeves (see page 52), these fittings make superior fingertight connections with capillary tubing. MicroTight Fittings withstand temperatures up to 125 °C.

The MicroTight family also includes a female nut matched with one of five dedicated ferrules for connecting specific tubing ODs.

Use the P-277 Extender Tool to tighten standard micro knurl 6-32 fittings in hard-to-reach places. Tighten micro headless 6-32 fittings with our N-290 Tool. See page 50 for more information.



MicroTight fittings and MicroFerrules

While the MicroTight Female Nuts may be used with any of the separate MicroFerrules, the MicroFerrules themselves are port-specific and are thus not interchangeable. Additionally, the one-piece MicroTight fittings are also port-specific and should not be exchanged.




Standard Head Fitting for use with 360 µm OD tubing



Standard Head Fitting for use with 1/32" OD tubing



F-152 MicroFerrule for 360 µm OD tubing



P-416 Female Nut 5/16-24 internal threads



F-124Hx Headless Fitting for use with 360 μm OD tubing



Standard Head Plug



Standard Head Fitting for use with MicroTight Sleeves



MicroFerrule for 1/32" OD tubing





MicroFerrule for 1/16" OD tubing



MicroFerrule for 0.025" OD tubing

P-416BLK Female Nut 5/16-24 internal threads



P-416G Female Nut 5/16-24 internal threads

RELATED PRODUCTS

Connectors for Capillary Tubing can be found on page 75. Very High Pressure fittings for capillary tubing can be found on page 58. Capillary tubing is featured on page 16.

Part No. Description Port Pressure Rating **Head Style** Material Qty. F-124Hx MicroTight Fitting for 360 µm OD Tubing 6-32 Coned 5,000 psi (345 bar) Micro Headless Knurl PEEK Blue 10-pk F-124Sx PEEK Blue MicroTight Fitting for 360 µm OD Tubing 6-32 Coned 5.000 psi (345 bar) Standard Micro Knurl 10-pk F-125Hx MicroTight Fitting for MicroTight Tubing Sleeves 6-32 Coned 4,000 psi (276 bar) Micro Headless Knurl PEEK Natural 10-pk F-125x MicroTight Fitting for MicroTight Tubing Sleeves 4,000 psi (276 bar) PEEK Natural 6-32 Coned Standard Micro Knurl 10-pk MicroTight Fitting for 1/32" OD Tubing F-1265x 5.000 psi (345 bar) Standard Micro Knurl PEEK Red 10-pk P-555x MicroTight Plug 6-32 Coned 5.000 psi (345 bar) Standard Micro Knurl PEEK Natural 10-pk MICROFER F-112 MicroFerrule for 1/32" OD Tubing 5/16-24 Coned 5,000 psi (345 bar) PEEK Natural F-132 MicroFerrule for 1/16" OD Tubing 5/16-24 Coned 5,000 psi (345 bar) PEEK Natural ea F-152 MicroFerrule for 360 µm OD Tubing 5/16-24 Coned 5.000 psi (345 bar) PEEK Natural ea. F-152BLK MicroFerrule for 360 µm OD Tubing 5/16-24 Coned 5,000 psi (345 bar) PEEK Black ea F-172 PEEK Black MicroFerrule for MicroTight Tubing Sleeves 5/16-24 Coned 4,000 psi (276 bar) ea. P-116 5/16-24 Coned 5.000 psi (345 bar) PEEK Black ea P-416 MicroTight Female Nut 5/16-24 Coned 4.000-5.000 psi (276-345 bar) Female Knurl PEEK Natural ea. P-416BLK MicroTight Female Nut 5/16-24 Coned 4,000-5,000 psi (276-345 bar) Female Knurl PEEK Black ea P-416G MicroTight Female Nut 5/16-24 Coned 4 000-5 000 psi (276-345 bar) Female Knurl PEEK Green ea



Two-Piece Fingertight Fittings

- Designed to connect tubing to 10-32 coned ports
- Ferrules available for directly connecting 1/16", 1/32", 360 μm, or 190 μm OD tubing
- > Economical, replace only the ferrule



APPLICATION NOTE

Some Fingertight Nuts feature wings in addition to a knurled head, which provide more leverage when tightening the fitting into a receiving port. Choose our single or double-winged design.

Please Note: customers can use the standard knurl head fittings with our tightening tools found on page 50.



Two-Piece Fingertight Fittings feature a separate ferrule. Use a standard knurled head fitting for traditional fingertight applications, or use a fitting with wings built into the head for extra tightening leverage. A stainless steel hex headed fitting can be used for applications where a wrench may be needed for added tightening torque.

The M-215 Conductive Perfluoroelastomer Ferrule is designed for mass spectrometer electrospray applications. Unlike most graphite ferrules, the elastomeric properties of this ferrule let you use it through many tightening/retightening cycles. It also eliminates any possibility of graphite contamination in your system. Like graphite ferrules, you can apply voltage through a metallic port block or metallic nut, allowing voltage to translate to the flow path through the ferrule.







F-151PCTFE Ferrule for 360 μm OD tubing



F-113 PEEK Ferrule for 1/32" OD tubing



F-142N ETFE Ferrule for 1/16" OD tubing



Conductive Perfluoroelastometer Ferrule for 360 µm OD tubing



10-32 PEEK Double-Winged Nut with F-142x PEEK Ferrule





| Part No. | Description | Port | Pressure Rating | Head Style | Material | Qty. | | |
|--|---|-------------|---------------------|----------------|-------------------------------|-------|--|--|
| TWO-PIECE FINGERTIGHT FITTINGS (INCLUDES F-142 FERRULES) | | | | | | | | |
| F-300x | Fingertight Fitting for 1/16" OD Tubing | 10-32 Coned | 6,000 psi (414 bar) | Double Wing | PEEK Natural | 10-pk | | |
| F-330x | Fingertight Fitting for 1/16" OD Tubing, Long | 10-32 Coned | 6,000 psi (414 bar) | Standard Knurl | PEEK Natural | 10-pk | | |
| F-331x | Fingertight Fitting for 1/16" OD Tubing | 10-32 Coned | 6,000 psi (414 bar) | Standard Knurl | PEEK Natural | 10-pk | | |
| REPLACEN | REPLACEMENT FERRULES | | | | | | | |
| F-113 | Ferrule for 1/32" OD Tubing | 10-32 Coned | 6,000 psi (414 bar) | _ | PEEK Natural | ea. | | |
| F-142 | Ferrule for 1/16" OD Tubing | 10-32 Coned | 6,000 psi (414 bar) | _ | PEEK Natural | ea. | | |
| F-142N | Ferrule for 1/16" OD Tubing | 10-32 Coned | 4,000 psi (276 bar) | _ | ETFE Natural | ea. | | |
| F-148 | Ferrule for 190 µm OD tubing | 10-32 Coned | 6,000 psi (414 bar) | _ | PCTFE Natural | ea. | | |
| F-151 | Ferrule for 360 µm OD Tubing | 10-32 Coned | 6,000 psi (414 bar) | _ | PCTFE Natural | ea. | | |
| M-215 | Conductive Ferrule for 360 um OD tubina | 10-32 Coned | 1.500 psi (103 bar) | _ | Conductive Perfluoroelastomer | ea. | | |
| | HROMalytic +61(0)3 97 | 62 2034 | alian Distributors | | | | | |

ECH no logy Pty Ltd

Website NEW: www.chromalytic.net.au E-mail: info@chromtech.net.au Tel: 03 9762 2034 . . . in AUSTRALIA



Two-Piece SealTight Fingertight Fittings

- Several nut lengths and head styles to fit into a variety of applications
- Designed to connect 1/16" OD tubing to 10-32 coned ports
- **>** Hold up to 9,000 psi (620 bar)





Overtightening these fittings on fluoropolymer (e.g., FEP, PFA, and ETFE) tubing can cause the ID of your tubing to collapse.



- ➤ Find tightening tools for these fittings on page 50.
- Try the F-350x FlushNut[™] for the ultimate streamline design.

The dual compression created by the specially designed nut and ferrule enables our SealTight™ Fittings system to outperform standard finger tightened fittings. The forward cone of the SealTight Ferrule provides gripping power and a leak-free seal via conventional compression by the receiving port. The slotted end creates the second compression zone in conjunction with a SealTight Nut. All SealTight Nuts are for use with 1/16″ OD tubing and are designed to be used with the F-192x Ferrule. A wide variety of fitting head styles are available for various space constraints.



F-193x 10-32 Short PEEK Hex Head Nut, with F-192x Ferrule







F-192x SealTight Ferrule, for 1/16" OD tubing



| Part No. | Description | Port | Pressure Rating | Head Style | Material | Qty. | | |
|--|---|-------------------|-------------------------------|----------------|------------|-------|--|--|
| SEALTIGHT TWO-PIECE FITTINGS (INCLUDES F-192 FERRULES) | | | | | | | | |
| F-193x | SealTight Fitting for 1/16" OD Tubing, Short | 10-32 Coned | 7,000-9,000 psi (483-620 bar) | 1/4" Hex | PEEK Black | 10-pk | | |
| F-195x | SealTight Fitting for 1/16" OD Tubing, Short | 10-32 Coned | 7,000-9,000 psi (483-620 bar) | Standard Knurl | PEEK Black | 10-pk | | |
| F-196x | SealTight Fitting for 1/16" OD Tubing, Long | 10-32 Coned | 7,000–9,000 psi (483–620 bar) | Standard Knurl | PEEK Black | 10-pk | | |
| F-287x | SealTight Fitting for 1/16" OD Tubing, Long | 10-32 Coned | 7,000-9,000 psi (483-620 bar) | Knurl-1/4" Hex | PEEK Black | 10-pk | | |
| F-350x | SealTight Fitting for 1/16" OD Tubing, FlushNut | 10-32 Coned | 7,000-9,000 psi (483-620 bar) | FlushNut | SST | 10-pk | | |
| REPLACEMENT FERRULES | | | | | | | | |
| F-192x | SealTight Ferrule for 1/16" OD Tubing | 10-32 or M6 Coned | 7.000-9.000 (483-620 bar) | _ | PEEK/Black | 10-pk | | |

FLUIDICS



NanoTight Fittings & Sleeves

- For connecting 1/16" OD or capillary tubing using tubing sleeves to standard 10-32 coned ports
- Multiple nut styles available
- Nuts manufactured from PEEK polymer, ferrules manufactured from ETFE




- Find tightening tools for these head styles on page 50.
- NanoTight Tubing sleeves start on page 52.

NanoTight™ Fittings and Sleeves are designed to connect 70 µm–1 mm OD capillary tubing to any standard 10-32 coned port normally intended for 1/16″ OD tubing using the NanoTight Tubing Sleeves on page 52. The fittings can also be used to connect any 1/16″ OD tubing. The ETFE ferrule material is softer than PEEK, making it a good candidate for connecting thin walled semi-rigid tubing such as FEP and ETFE into 10-32 ports with minimal constricting to the inner diameter.

Select from our expansive line of PEEK NanoTight Fittings, featuring several head style and length options. Each 10-pack of nuts includes ten ETFE F-142Nx ferrules.



F-333Nx Short Headless Nut with F-142Nx Ferrule




Website NEW: www.chromalytic.net.au E-mail: info@chromtech.net.au Tel: 03 9762 2034... in AUSTRALIA



LiteTouch Fittings

- Helps prevent twisting of polymer tubing
- > High pressure with fingertight convenience
- Doptions available for 1/32", 1/16", or 1/8" OD tubing







The stainless steel nuts on page 59 can also be used with the LiteTouch ferrules on this page.

The LiteTouch® Fittings System grips tubing at two compression points (see diagram), holding to high pressures with Fingertight convenience. It also prevents polymer tubing from twisting, a potential problem when using standard Fingertight fittings. LiteTouch Fittings are available for use with 1/32", 1/16", or 1/8" OD tubing sizes, and for 10-32 or 1/4-28 coned ports.


For those space-limited applications where nut heads interfere with each other, try the FlushNut™ Fittings. (FlushNut Fittings require a tightening tool. Please see page 50 for more information about these products.)

To avoid collapsing the ID of your tubing, the LiteTouch system can be used on hard tubing only, such as stainless steel and PEEK polymer tubing. The LiteTouch Ferrule System is not recommended for repeated use in plastic ports.



F-354x





LT-210x 1/4-28 PEEK Double-Winged Nut for 1/8" OD tubing



LT-110x 10-32 PEEK Nut for 1/32" and 1/16" OD tubing



C-235x 1/4-28 Stainless Steel Nut for 1/8" OD tubing



F-364x 1/4-28 Stainless Steel FlushNut for 1/8" OD tubing



LT-132x PEEK Ferrule with Stainless Steel Lock Ring for 1/32" OD tubing



LT-100x PEEK Ferrule with Stainless Steel Lock Ring for 1/16" OD tubing



LT-200x PEEK Ferrule with Stainless Steel Lock Ring for 1/8" OD tubing

| Part No. | Description | Port | Pressure Rating | Head Style | Material | Qty. |
|------------|--|--------------|---------------------|----------------|------------------|-------|
| LITETOU | CH NUTS | | | | | |
| C-235x | LiteTouch Nut for 1/8" OD Tubing | 1/4-28 Coned | 4,500 psi (310 bar) | 5/16" Hex | SST | 10-pk |
| F-354x | LiteTouch Nut for 1/16" or 1/32" OD Tubing, FlushNut | 10-32 Coned | 5,000 psi (345 bar) | FlushNut | SST | 10-pk |
| F-364x | LiteTouch Nut for 1/8" OD Tubing, FlushNut | 1/4-28 Coned | 4,500 psi (310 bar) | FlushNut | SST | 10-pk |
| LT-110x | LiteTouch Nut for 1/16" or 1/32" OD Tubing | 10-32 Coned | 5,000 psi (345 bar) | Standard Knurl | PEEK Natural | 10-pk |
| LT-210x | LiteTouch Nut for 1/8" OD Tubing | 1/4-28 Coned | 4,500 psi (310 bar) | Double Wing | PEEK Natural | 10-pk |
| LITETOU | CH FERRULES | | | | | |
| LT-100x | LiteTouch Ferrule for 1/16" OD Tubing | 10-32 Coned | 5,000 psi (345 bar) | _ | PEEK Natural/SST | 10-pk |
| LT-132x | LiteTouch Ferrule for 1/32" OD Tubing | 10-32 Coned | 5,000 psi (345 bar) | _ | PEEK Natural/SST | 10-pk |
| LT-200x | LiteTouch Ferrule for 1/8" OD Tubing | 1/4-28 Coned | 4,500 psi (310 bar) | _ | PEEK Natural/SST | 10-pk |
| * When use | ed with a stainless steel 10-32 nut from page 59. | | | | | |





Super Flangeless Fittings

Our Super Flangeless™ Fittings provide the highest pressure holding capability in a flat-bottom fitting system that we offer. Our unique design eliminates loosening of fittings due to tubing twist and holds tight even through vibration. Our high pressure fittings are excellent for tubing assemblies and those times when connections need to be broken frequently.

Highest pressure holding flat-bottom fitting system we offer

- Eliminates loosening of fittings due to tubing twist
- > Excellent for Tubing Assemblies
- > Holds tight even through vibration



ASSEMBLY HINT

Make sure the locking ring is oriented correctly! The flattened end of the ring should face towards the nut with the narrow end of the ferrule towards the ring.



Enlarged to show detail

6-40 & 6-32 options (for 1-16" OD tubing)




Super Flangeless Ferrule for 1/16" OD tubing 6-40 Nut shown with M-650x Ferrule (not included)

6-32 Nut shown with M-650x Ferrule (not included)

10-32 options (for 1-16" OD tubing)



SUPER FLANGELESS FITTINGS SYSTEM



Ferrules for M6x1, 1/4-28, 5/16-24





(0.48 cm)

P-360x

for 1/8" OD tubing



for 1/8" OD tubing

Super Flangeless[™] Fittings (Cont.)

M6x1, 1/4-28, 5/16-24 Options for 1/32"-3/16" OD Tubing



(1/4-28, ≤ 1/16" OD Tubing)

(1/4-28, > 1/16" OD Tubing)

F-356x

F-364x



LT-105x (1/4-28, ≤ 1/16" OD Tubing) C-235x (1/4-28, > 1/16" OD Tubing)



P-213x (M6X1, ≤ 1/16" OD Tubing) P-337x (M6X1, > 1/16" OD Tubing) P-232x

(1/4-28, ≤ 1/16" OD Tubing) P-336x (1/4-28, > 1/16" OD Tubing)



P-287x (1/4-28, ≤ 1/16" OD Tubing) P-387x (1/4-28, > 1/16" OD Tubing)



P-219x (M6X1, ≤ 1/16" OD Tubing) P-319x (M6X1, > 1/16" OD Tubing) LT-115x (1/4-28, ≤ 1/16" OD Tubing)




P-217x (M6X1, ≤ 1/16" OD Tubing) P-317x

(M6X1, > 1/16" OD Tubing) P-246x, P-255x, P-281x P-331x, P-381x (1/4-28, > 1/16" OD Tubing)



LT-210x (1/4-28, > 1/16" OD Tubing)



P-420x (1/4-28, ≤ 1/16" OD Tubing) F-156x (1/4-28, > 1/16" OD Tubing)



(5/16-24, ≤ 1/16" OD Tubing) P-137x (5/16-24, 3/16" OD Tubing)

One-Piece Super Flangeless™ Fittings

- ➤ All-PEEK construction
- > For 1/16" OD and 1/8" OD tubing
- ▶ M6x1 and 1/4-28 options
- Finger tight (2–3 in-lbs / 0.23–0.34 N⋅m)
- > Extremely easy to use
- > Reusable one piece design that requires no swaging



P-329x M6X1, for 1/8" OD Tubing P-249x 1/4-28, for 1/16" OD Tubing P-349x 1/4-28, for 1/8" OD Tubing



P-229x M6X1, for 1/16" OD Tubing

ou Flancologa Tubing OD / Throad Commonica.

| Super Flangeless Tubir | ig OD / Thread Comparison | | | | | |
|------------------------|---------------------------|-------|--------|----------|----------|-------|
| | 1/32" | 1/16" | 1.8 MM | 2.5 MM | 1/8" | 3/16" |
| 6-40 | | ✓ | | | | |
| 6-32 | | ✓ | | | | |
| 10-32 | | ✓ | | | | |
| M6x1 | ✓ | ✓ | ✓ | ✓ | ✓ | |
| 1/4-28 | ✓ | ✓ | ✓ | ✓ | ✓ | |
| 5/16-24 | | | | | | ✓ |

FLUIDICS

Super Flangeless[™] and One-Piece Super Flangeless Fittings

| Part No. | Description | Port | Pressure Rating | Head Style | Material | Qty. |
|-----------|---|--------------------------|------------------------------|----------------|------------------|-------|
| | LANGELESS FERRULES FOR 1/32", 1/16", 1/8", 3/16", | | | | | |
| 1-250x | Super Flangeless Ferrule for 1/16" OD Tubing | 10-32 Flat-Bottom | 1,000–5,000 psi (69–345 bar) | _ | PEEK Natural/SST | 10-р |
| 1-650x | Super Flangeless Ferrule for 1/16" OD Tubing | 6-32 or 6-40 Flat Bottom | 750–3,750 psi (52–259 bar) | _ | PEEK Natural/SST | 10-p |
| 248x | Super Flangeless Ferrule for 1/32" OD Tubing | 10-32 Flat-Bottom | 2,500 psi (172 bar) | _ | ETFE Green/SST | 10-p |
| -250x | Super Flangeless Ferrule for 1/16" OD Tubing | 1/4-28 or M6 Flat Bottom | 2,500 psi (172 bar) | _ | PEEK Natural/SST | 10-p |
| -259x | Super Flangeless Ferrule for 1/16" OD Tubing | 1/4-28 or M6 Flat Bottom | 1,350 psi (93 bar) | _ | ETFE Yellow/SST | 10-p |
| -260x | Super Flangeless Ferrule for 1/16" OD Tubing | 1/4-28 or M6 Flat Bottom | 1,850 psi (128 bar) | _ | PEEK Natural/SST | 10-p |
| -350x | Super Flangeless Ferrule for 1/8" OD Tubing | 1/4-28 Flat Bottom | 2,500 psi (172 bar) | _ | PEEK Natural/SST | 10-p |
| -352x | Super Flangeless Ferrule for 1/16" OD Tubing | 1/4-28 or M6 Flat Bottom | 2,500 psi (172 bar) | _ | PEEK Black/SST | 10-p |
| -355x | Super Flangeless Ferrule for 1.8 mm OD Tubing | 1/4-28 or M6 Flat Bottom | 2,500 psi (172 bar) | _ | PCTFE Green/SST | 10-p |
| -357-2x | Super Flangeless Ferrule for 2.0 mm OD Tubing | M6 Flat Bottom | 5,000 psi (345 bar) | _ | PEEK Natural/SST | 10-p |
| P-359x | Super Flangeless Ferrule for 1/8" OD Tubing | 1/4-28 Flat Bottom | 1,000 psi (69 bar) | _ | ETFE Yellow/SST | 10-p |
| P-360x | Super Flangeless Ferrule for 1/8" OD Tubing | 1/4-28 Flat Bottom | 1,500 psi (102 bar) | _ | PEEK Natural/SST | 10-p |
| -366x | Super Flangeless Ferrule for 2.5" OD Tubing | 1/4-28 Flat Bottom | 1,000 psi (69 bar) | _ | PEEK Natural/SST | 10-p |
| P-140x | Super Flangeless Ferrule for 3/16" OD Tubing | 5/16-24 Flat Bottom | 500 psi (34 bar) | _ | ETFE Green/SST | 10-р |
| | D 6-32 FITTINGS FOR 1/16" OD TUBING | | | | | |
| M-660x | Super Flangeless Nut for 1/16" OD Tubing | 6-32 Flat Bottom | 750–3,750 psi (52–259 bar) | Micro Headless | PEEK Natural | 10-p |
| | x Super Flangeless Nut for 1/16" OD Tubing | 6-40 Flat Bottom | 750–3,750 psi (52–259 bar) | Micro Headless | PEEK Green | 10-p |
| 0-32 FIT | FTINGS FOR 1/16" OD TUBING | | | | | |
| Л-652x | Super Flangeless Nut for 1/16" OD Tubing | 10-32 Flat Bottom | 1,000–5,000 psi (69–345 bar) | 1/4" Hex | PEEK Green | 10-p |
| И-653x | Super Flangeless Nut for 1/16" OD Tubing | 10-32 Flat Bottom | 1,000–5,000 psi (69–345 bar) | Headless Knurl | PEEK Green | 10-p |
| M-655x | Super Flangeless Nut for 1/16" OD Tubing, Long | 10-32 Flat Bottom | 1,000–5,000 psi (69–345 bar) | 1/4" Hex | PEEK Green | 10-p |
| M6X1 FI | TTINGS FOR 1/16" AND 1/32" OD TUBING | | | | | |
| P-213x | Super Flangeless Nut for 1/16" or 1/32" OD Tubing, Short | M6 Flat Bottom | * | Headless Knurl | PEEK Black | 10-p |
| P-217x | Super Flangeless Nut for 1/16" or 1/32" OD Tubing | M6 Flat Bottom | * | Standard Knurl | PPS Black | 10-p |
| P-219x | Super Flangeless Nut for 1/16" or 1/32" OD Tubing, Short | M6 Flat Bottom | * | Standard Knurl | PEEK Black | 10-p |
| M6X1 FI | TTINGS FOR 1.8 MM, 20. MM, 2.5 MM, 1/8" OD TUBI | NG | | | | |
| -317x | Super Flangeless For >1/16"-≤ 1/8" OD Tubing | M6 Flat Bottom | * | Standard Knurl | PPS Black | 10-р |
| -319x | Super Flangeless Nut for 1/8" OD Tubing, Short | M6 Flat Bottom | * | Standard Knurl | PEEK Black | 10-р |
| P-337x | Super Flangeless For >1/16"-≤ 1/8" OD Tubing, Short | M6 Flat Bottom | * | Headless Knurl | PEEK Black | 10-р |
| P-357x | Super Flangeless Fitting for 2.0 mm OD Tubing | M6 Flat Bottom | * | Standard Knurl | PEEK Black, | 10-р |
| | | | | | Natural/SST | |
| | TTINGS FOR 1/16" AND 1/32" OD TUBING | | | | | |
| -356x | Super Flangeless Nut for 1/16" or 1/32" OD Tubing, FlushNut | 1/4-28 Flat Bottom | * | FlushNut | SST | 10-p |
| T-105x | Super Flangeless Nut for 1/16" or 1/32" OD Tubing, Short | 1/4-28 Flat Bottom | * | 1/4" Hex | SST | 10-p |
| T-115x | Super Flangeless Nut for 1/16" or 1/32" OD Tubing, Short | 1/4-28 Flat Bottom | | Standard Knurl | PEEK Natural | 10-p |
| P-232x | Super Flangeless Nut for 1/16" or 1/32" OD Tubing, Short | 1/4-28 Flat Bottom | * | Headless Knurl | PEEK Natural | 10-p |
| P-246x | Super Flangeless Nut for 1/16" or 1/32" OD Tubing | 1/4-28 Flat Bottom | * | Standard Knurl | PFA Natural | 10-p |
| P-255x | Super Flangeless Nut for 1/16" or 1/32" OD Tubing | 1/4-28 Flat Bottom | * | Standard Knurl | PEEK Natural | 10-р |
| P-281x | Super Flangeless Nut for 1/16" or 1/32" OD Tubing | 1/4-28 Flat Bottom | * | Standard Knurl | PPS Natural | 10-p |
| P-287x | Super Flangeless Nut for 1/16" or 1/32" OD Tubing | 1/4-28 Flat Bottom | * | Headless Knurl | PPS Natural | 10-p |
| P-420x | Super Flangeless Nut for 1/16" or 1/32" OD Tubing, Female | 1/4-28 Flat Bottom | * | Female Knurl | PEEK Natural | 10-p |
| I/4-28 FI | TTINGS FOR 1.8MM, 2.5 MM, 1/8" OD TUBING | | | | | |
| :-235x | Super Flangeless Nut for 1/8" OD Tubing | 1/4-28 Flat Bottom | * | 1/4" Hex | SST | 10-p |
| -156x | Super Flangeless Nut for 1/8" OD Tubing, Female | 1/4-28 Flat Bottom | * | Female Knurl | PEEK Black | 10-p |
| -364x | Super Flangeless Nut for 1/8" OD Tubing, FlushNut™ | 1/4-28 Flat Bottom | * | FlushNut | SST | 10-p |
| T-210x | Super Flangeless Nut for 1/8" OD Tubing | 1/4-28 Flat Bottom | * | Double Wings | PEEK Natural | 10-p |
| P-331x | Super Flangeless Nut for 1/8" OD Tubing | 1/4-28 Flat Bottom | * | Standard Knurl | PEEK Natural | 10-p |
| P-336x | Super Flangeless Nut for 1/8" OD Tubing, Short | 1/4-28 Flat Bottom | * | Headless Knurl | PEEK Natural | 10-р |
| P-381x | Super Flangeless Nut for 1/8" OD Tubing | 1/4-28 Flat Bottom | * | Standard Knurl | PPS Natural | 10-p |
| -387x | Super Flangeless Nut for 1/8" OD Tubing | 1/4-28 Flat Bottom | * | Standard Knurl | PPS Natural | 10-p |
| 5/16-24 F | FITTINGS FOR 1/16", 1/8", 3/16" OD TUBING | | | | | |
| P-137x | Super Flangeless Fitting for 3/16" OD Tubing | 5/16-24 Flat Bottom | * | Standard Knurl | PEEK Black | 10-p |
| P-141x | Super Flangeless Fitting for 1/16" OD Tubing | 5/16-24 Flat Bottom | * | Standard Knurl | PEEK Natural | 10-p |
| | CE SUPER FLANGELESS FITTINGS FOR 1/16" AND 1/ | | | | | |
| 2-229x | One Piece Super Flangeless Fitting for 1/16" OD Tubing | M6 Flat Bottom | 1,000 psi (69 bar) | Standard Knurl | PEEK | 10-p |
| 2-249x | One Piece Super Flangeless Fitting for 1/16" OD Tubing | 1/4-28 Flat Bottom | 1,000 psi (69 bar) | Standard Knurl | PEEK | 10-p |
| -329x | One Piece Super Flangeless Fitting for 1/16" OD Tubing | M6 Flat Bottom | 1,000 psi (69 bar) | Standard Knurl | PEEK | 10-p |
| 2-349x | One Piece Super Flangeless Fitting for 1/16" OD Tubing | 1/4-28 Flat Bottom | 1,000 psi (69 bar) | Standard Knurl | PEEK | 10-p |
| | rating of nut depends on the ferrule used. | | , | | | . о р |





VacuTight Fittings


- For 1/16" or 1/8" OD tubing connections into 10-32, 1/4-28, or M6 flat-bottom ports
- > Vacuum Rated to 25 in-Hg (84 kPa)
- Improve transfer volume consistency

VacuTight Fittings are designed to provide airtight, dependable connections under vacuum and low pressure conditions. Many of the VacuTight Nuts have streamlined profiles for use in systems requiring a large number of connections in a small area. Furthermore, the VacuTight Ferrule's small size ensures sufficient nut/thread engagement, even in shallow ports. These features make VacuTight Fittings ideal for "combichem" high throughput screening, clinical diagnostic, and other automated liquid handling applications.

The configuration of the VacuTight flat-bottom ferrules prevents overcompression and tubing ID reduction that can occur with many coned fittings. The result is more consistent aspirating and dispensing volumes across all system connections.

The VacuTight fittings can also work well in some positive pressure applications. The pressure range for each fitting is listed below and depends upon the tubing used for the connection. Please contact your distributor or IDEX Health & Science for more information. Additionally, please note that some of the VacuTight fittings have changed in color from red to black; however, this color change does not affect product performance.

All VacuTight Nuts must be used exclusively with VacuTight Ferrules.



- ¹ The dimensions shown apply to P-930x, P-931x, P-938x, P-942x, and P-948x
- ² The dimensions shown apply to P-945x. ³ The dimensions shown apply to P-946x.

| Part No. | Description | Port | Pressure Rating | Head Style | Material | Qty. |
|----------|--|--------------------------|---------------------------|----------------|--------------|-------|
| VACUTIGH | HT FITTINGS (INCLUDES P-840 OR P-940 FER | RRULES) | | | | |
| P-842x | VacuTight Fitting for 1/16" OD Tubing, Short | 10-32 Flat-Bottom | 400-800 psi (27-55 bar) | 1/4" Hex | PEEK Red | 10-pk |
| P-844x | VacuTight Fitting for 1/16" OD Tubing, Short | 10-32 Flat-Bottom | 400–800 psi (27–55 bar) | Headless Knurl | PEEK Red | 10-pk |
| P-846x | VacuTight Fitting for 1/16" OD Tubing, Long | 10-32 Flat-Bottom | 400-800 psi (27-55 bar) | Headless Knurl | PEEK Red | 10-pk |
| P-930x | VacuTight Fitting for 1/16" OD Tubing | 1/4-28 Flat-Bottom | 400–800 psi (27–55 bar) | Standard Knurl | Delrin Red | 10-pk |
| P-931x | VacuTight Fitting for 1/16" OD Tubing | M6 Flat-Bottom | 400-800 psi (27-55 bar) | Standard Knurl | Delrin Red | 10-pk |
| P-938x | VacuTight Fitting for 1/16" OD Tubing | 1/4-28 Flat-Bottom | 400–800 psi (27–55 bar) | Standard Knurl | PEEK Natural | 10-pk |
| P-942x | VacuTight Fitting for 1/8" OD Tubing | 1/4-28 Flat-Bottom | 500-1,000 psi (34-69 bar) | Standard Knurl | Delrin Red | 10-pk |
| P-945x | VacuTight Fitting for 1/8" OD Tubing, Short | M6 Flat-Bottom | 500-1,000 psi (34-69 bar) | Standard Knurl | Delrin Black | 10-pk |
| P-946x | VacuTight Fitting for 1/8" OD Tubing | 1/4-28 Flat-Bottom | 500-1,000 psi (34-69 bar) | Headless Knurl | Delrin Red | 10-pk |
| P-948x | VacuTight Fitting for 1/8" OD Tubing | 1/4-28 Flat-Bottom | 500-1,000 psi (34-69 bar) | Standard Knurl | PEEK Natural | 10-pk |
| REPLACE | MENT FERRULES | | | | | |
| P-840 | VacuTight Ferrule for 1/16" OD Tubing | M6 or 1/4-28 Flat-Bottom | 400-800 psi (27-55 bar) | _ | ETFE Red | ea. |
| P-940x | VacuTiaht Ferrule for 1/8" OD Tubina | M6 or 1/4-28 Flat-Bottom | 500-1.000 psi (34-69 bar) | _ | ETFE Red | 10-pk |



Metric Flangeless Fittings

> For 1/16", 1.8 mm, 2.0 mm, 2.5 mm, 3.0 mm, 4.0 mm, or 1/8" OD tubing

Convenience of flangeless fittings for metric tubing sizes and M6

flat-bottom ports

0.19" (0.48 cm) 0.19" (0.48 cm)

P-342x ETFE Ferrule for 1.8 OD tubing

0.19" (0.48 cm)

0.19" (0.48 cm)

ETFE Ferrule for 2.0 mm OD tubing (0.48 cm)

P-363Rx




0.19" (0.48 cm) P-343x

P-353x ETFE Ferrule for 2.5 mm OD tubing

ETFE Ferrule for 3.0 mm OD tubing





P-344x ETFE Ferrule for 4.0 mm OD tubing P-347x PEEK Nut for 1.8, 2.0, 3.0 mm or 1/8" OD tubing



Delrin® Nut for 1/16" OD tubing

tubing to flat-bottom ports when paired with the appropriate M6, 1/4-28, or 5/16-24 Flangeless Nuts. We also offer M6-threaded nuts to connect 1/16" or 1/8" OD tubing, plus a tubing sleeve to facilitate 1.0 mm OD tubing connections. Please refer to the "Metric Connections" chart on this page for information regarding which nuts and ferrules to use with your tubing.

Metric Flangeless Ferrules are designed to connect 1.8, 2.0, 2.5, 3.0, or 4.0 mm OD

METRIC CONNECTIONS

Use this chart to determine the low pressure fittings needed to connect metric and English-sized tubing into the indicated ports.

| Tubing Size | Port | Ferrules | Nuts |
|-------------|----------|---------------------------------------|---|
| 1.0 mm | M6 | P-200x (w/F-252 sleeve, not included) | P-207x, P-207Sx, P-247x |
| | 1/4-28 | P-200x (w/F-252 sleeve, not included) | Any 1/4-28 nut for 1/16" OD tubing from page 47 |
| 1.8 mm | M6 | P-342x | P-307x, P-307Sx, P-347x |
| | 1/4-28 | P-342x | Any 1/4-28 nut for 1/8" OD tubing from page 47 |
| 2.0 mm | M6 | P-363Rx | P-307x, P-307Sx, P-347x |
| | 1/4-28 | P-363Rx | Any 1/4-28 nut for 1/8" OD tubing from page 47 |
| 2.5 mm | M6 | P-353x | P-307x, P-307Sx, P-347x |
| | 1/4-28 | P-353x | Any 1/4-28 nut for 1/8" OD tubing from page 47 |
| 3.0 mm | M6 | P-343x | P-307x, P-307Sx, P-347x |
| | 1/4-28 | P-343x | Any 1/4-28 nut for 1/8" OD tubing from page 47 |
| 4.0 mm | 5/16-24 | P-344x | XP-132x from page 56 |
| 1/16" | M6 M6 | P-200x P-840 | P-207x, P-207Sx, P-247x, P-931, page 42 |
| 1/8" | M6 M6 | P-300x P-940x | P-307x, P-307Sx, P-347x, P-945x, page 42 |

RELATED PRODUCTS

| | Page | | Page |
|--|------|---|------|
| MORE METRIC-SIZED PRODUCTS | | | |
| High Pressure Polymer Fittings | 58 | Low Pressure Unions | 84 |
| High Pressure Stainless Steel Fittings | 59 | Bulkhead Unions | 82 |
| Luer Adapters | 87 | PEEK (1.8 mm OD and Capillary) and Fused Silica Tubing | 16 |
| Metric Threaded Adapters | 65 | PEEKsil™ Tubing | 22 |
| External NPT Adapters | 66 | FEP Tubing (1.0–4.0 mm OD) and PFA Capillary Tubing | 26 |
| VacuTight™ Fittings | 42 | Frit-In-A-Ferrule™ | 99 |
| Plugs and Cans | 55 | | |

In addition, many of our 1/4-28 threaded Filters, Valves and Flow Control Accessories can be converted to accept 1.8, 2.0, 2.5 and 3.0 mm tubing, using the ferrules listed for 1/4-28 ports in the "Metric Connections" table, this page.


| Part No. | Description | Port | Pressure Rating | Head Style | Material | Qty. |
|----------|---|--|---------------------|----------------|--------------|-------|
| METRIC F | LANGELESS NUTS | | | | | |
| P-207x | Flangeless Nut for 1/16" OD Tubing | M6 Flat-Bottom | 2,000 psi (138 bar) | Standard Knurl | Delrin Black | 10-pk |
| P-2075x | Flangeless Nut for 1/16" OD Tubing, Short | M6 Flat-Bottom | 2,000 psi (138 bar) | Standard Knurl | Delrin Black | 10-pk |
| P-247x | Flangeless Nut for 1/16" OD Tubing, Short | M6 Flat-Bottom | 2,000 psi (138 bar) | Standard Knurl | PEEK Black | 10-pk |
| P-307x | Flangeless Nut for 1.8 mm, 2.0 mm, 3.0 mm, 1/8" OD Tubing | M6 Flat-Bottom | 500 psi (34 bar) | Standard Knurl | Delrin Black | 10-pk |
| P-307Sx | Flangeless Nut for 1.8 mm, 2.0 mm, 3.0 mm, 1/8" OD Tubing | M6 Flat-Bottom | 500 psi (34 bar) | Standard Knurl | Delrin Black | 10-pk |
| P-347x | Flangeless Nut for 1.8 mm, 2.0 mm, 3.0 mm, 1/8" OD Tubing | M6 Flat-Bottom | 500 psi (34 bar) | Standard Knurl | PEEK Black | 10-pk |
| FLANGE | LESS FERRULES | | | | | |
| F-252x | 1/16" OD Tubing Sleeve for 1.0 mm ID Tubing | M6 or 1/4-28 Flat-Bottom | 500 psi (34 bar) | _ | FEP Purple | 10-pk |
| P-200x | Flangeless Ferrule for 1/16" OD Tubing | M6 or 1/4-28 Flat-Bottom | 2,000 psi (138 bar) | _ | ETFE Blue | 10-pk |
| P-300x | Flangeless Ferrule for 1/8" OD Tubing | M6 or 1/4-28 Flat-Bottom | 500 psi (34 bar) | _ | ETFE Yellow | 10-pk |
| P-342x | Flangeless Ferrule for 1.8 mm OD Tubing | M6 or 1/4-28 Flat-Bottom | 500 psi (34 bar) | _ | ETFE Green | 10-pk |
| P-343x | Flangeless Ferrule for 3.0 mm OD Tubing | M6 or 1/4-28 Flat-Bottom | 500 psi (34 bar) | _ | ETFE Orange | 10-pk |
| P-344x | Flangeless Ferrule for 4.0 mm OD Tubing | 5/16-24 | 250 psi (17 bar) | _ | ETFE Natural | 10-pk |
| P-353x | Flangeless Ferrule for 2.5 mm OD Tubing | M6 or 1/4-28 Flat-Bottom | 500 psi (34 bar) | _ | ETFE Natural | 10-pk |
| P-343B4 | Flangeless Ferrula for 2 0 mm OD Tubing | M6 or 1/4-28 Flat-Rottom | 500 nei (34 bar) | _ | ETFE Red | 10-pk |
| | HROMalytic +61(0)3 9762 2034 | Australian Distributors Importers & Manufacurers | | | | |



Flanged Fittings

- > Fittings for 1/16" or 1/8" OD tubing, supplied with nut and 316 stainless steel washer
- Multiple head styles and materials available; contact IDEX Health & Science for more information
- For 1/4-28 and M6 flat-bottom ports
- Some color options available; call for more information


Flanged Fittings are compatible with most standard 1/4-28 or M6 Flat-Bottom flanged fittings. The Delrin® (acetal resin) nut resists cross threading or loosening during use.



P-401x1

Flanged Fitting for 1/16" OD tubing*

- ¹ The dimensions shown apply to all square-head Flanged Fittings
- * Flanged tubing not included



P-482G3

Flanged Fitting for 1/16" OD tubing*

- ³ The dimensions shown apply to all knurled-head Flanged Fittings
- * Flanged tubing not included

RELATED PRODUCTS

> For an alternative to flanging tubing, we highly recommend the Flangeless Fittings found on page 45, the Super Flangeless Fittings found on page 39, or the VacuTight™ Fittings on page 42.

| Part No. | Description | Port Geometry | Head Style | Material (Nut/Washer) | Qty. |
|---------------------|--|---|----------------|-----------------------|------------|
| FLANGED FITT | INGS (INCLUDES STAINLESS STEEL WASHERS) | | | | |
| P-401x | Flanged Fitting for 1/16" OD Tubing | 1/4-28 Flat-Bottom | 5/16" Square | Delrin Black/SST | 10-pk |
| P-482BLK | Flanged Fitting for 1/16" OD Tubing | 1/4-28 Flat-Bottom | Standard Knurl | Delrin Black/SST | ea. |
| P-501x | Flanged Fitting for 1/8" OD Tubing | 1/4-28 Flat-Bottom | 5/16" Square | Delrin Black/SST | 10-pk |
| P-982BLKx | Flanged Fitting for 1/16" OD Tubing | M6 Flat-Bottom | Standard Knurl | Delrin Black/SST | 10-pk |
| P-1082BLKx | Flanged Fitting for 1/8" OD Tubing | M6 Flat-Bottom | Standard Knurl | Delrin Black/SST | 10-pk |
| REPLACEMENT | Γ WASHERS | | | | |
| P-407x | Washer for 1/16" OD Tubing | 1/4-28 Flat-Bottom | _ | SST | 10-pk |
| P-507x | Washer for 1/8" OD Tubing | 1/4-28 Flat-Bottom | _ | SST | 10-pk |
| P-1087x | Washer for 1/8" OD Tubina | M6 Flat-Bottom | _ | SST | 10-pk |
| | HROMalytic +61(0)3 9762 2034 | Australian Distributors Importers & Manufacurers www.chromtech.net.au | | | |
| Website NEW : w | www.chromalytic.net.au E-mail : info@chromtech.net.au To | el: 03 9762 2034 in AUSTRA | LIA | 44 (| of 172) 20 |



THE CONVENIENCE OF FLANGELESS FITTINGS



Flangeless Fittings

Flangeless Fittings eliminate the need to flange tubing. This removable and reusable system provides several benefits:

Convenience: Flangeless Fittings are easy to use. Just slip the nut and ferrule over the tubing and finger tighten the assembly into your receiving port. In tests, it is shown that the ideal amount of torque to achieve expected part performance should be approximately 3–4 in-lbs (0.34–0.45 N·m). Check out the line of special tightening tools designed to adapt to many standard torque wrenches, on page 50 and the adjustable torque driver, VHP-4000 on page 51.

Minimal Down-Time: Component replacement is quick, taking only a few seconds — unlike the significant time required to flange tubing.

Cost-Effectiveness: Repairing a flanged tubing assembly requires a costly flanging tool or the purchase of a complete replacement assembly, including a new length of tubing and a set of fittings. The Flangeless Fittings system typically requires only one new ferrule at minimal cost when repairing a connection.

The 1/4-28 and M6 Flangeless Fittings for 1/16", 1/8", and metric sized OD tubing are summarized on the following page and listed on page 47.

1/4-28 Flangeless Fittings – Nuts



OLUMN HARDWARE

Flangeless Fittings (Cont.)



NOTE

The XP-340x ferrule is designed for use with shallow receiving ports, such as those used on some low pressure valves.



XP-340x

ETFE Small Valve Ferrule

The XF-368x FlushNut is an excellent choice for applications where port-to-port spacing is limited; see page 31 for more information on this innovative product line. As an alternative, consider one of the "headless" fittings shown on this page.

Ferrules



Standard 1/16" P-200x XP-200Nx



Small Valve 1/8" XP-340x



Small Valve 1/16" XP-240x



Standard 4.0 mm



Standard 1/8" and Metric Ferrules XP-300x



RELATED PRODUCTS

- ➤ For the Large Bore Flangeless Fittings, please refer to page 56.
- Nuts for M6 threaded ports are on page 43; nuts for 5/16-24 threaded ports are on page 56.

Flangeless Fittings for 1/16" and 1/8"OD Tubing

- > Wide variety of materials and geometries to fit most applications
- > Excellent replacement for flanged fittings
- Convenient and easy to use
- > Fittings and ferrules packaged together for easy ordering convenience

Flangeless Fittings for 1/16" OD Tubing, and for 1/8" OD Tubing are available in a variety of materials. The replacement ferrules are manufactured from inert ETFE, and are sold in a colored version or ETFE's natural color as the N option. The smaller ferrules, XP-240x and XP-340x are designed for shallow ports.

The fittings shown on this page are packaged together with the ferrules in convenient 10-packs. The 1/16" version include the XP-200x ferrule, and the 1/8" version include the XP-300x ferrule. Please visit our website, www.idex-hs.com, for single packaging options.

For higher pressure and temperature applications consider our Super Flangeless™ found on page 39.

Lock Nut

The XP-312x Lock Nut is for use with any 1/4-28 male Flangeless Fitting. Use this product in applications where vibrations can loosen fittings.

To Use: Thread the lock nut onto the male fitting. When the male fitting is firmly seated into the receiving port, tighten the lock nut down against the receiving port to securely hold the male fitting in place.



XP-312x Lock Nut White Delrin

FLUIDICS

Flangeless Fittings

| Part No. | Description | Port | | Pressure Rating | Head Style | Material | Qty. |
|----------------|---|------------------|--|---------------------|----------------|---------------|----------|
| FLANGEL | ESS FITTINGS FOR 1/16" OD TUBING | | | | | | |
| XF-358x | Flangeless Fitting for 1/16" OD Tubing, FlushNut | 1/4-28 Flat- | Bottom | 2,000 psi (138 bar) | FlushNut | SST | 10-pk |
| XLT-111x | Flangeless Fitting for 1/16" OD Tubing | 10-32 Flat-B | ottom | 2,500 psi (172 bar) | Standard Knurl | PEEK Natural | 10-pk |
| XP-201x | Flangeless Fitting for 1/16" OD Tubing | 1/4-28 Flat- | Bottom | 2,000 psi (138 bar) | Standard Knurl | Delrin Black | 10-pk |
| XP-202x | Flangeless Fitting for 1/16" OD Tubing | 1/4-28 Flat- | Bottom | 2,000 psi (138 bar) | Standard Knurl | Delrin Red | 10-pk |
| XP-218x | Flangeless Fitting for 1/16" OD Tubing | 1/4-28 Flat- | Bottom | 2,000 psi (138 bar) | Standard Knurl | ETFE Natural | 10-pk |
| XP-230x | Flangeless Fitting for 1/16" OD Tubing | 1/4-28 Flat- | Bottom | 2,000 psi (138 bar) | Standard Knurl | PEEK Natural | 10-pk |
| XP-235x | Flangeless Fitting for 1/16" OD Tubing, Short | 1/4-28 Flat- | Bottom | 2,000 psi (138 bar) | Standard Knurl | PEEK Natural | 10-pk |
| XP-238x | Flangeless Fitting for 1/16" OD Tubing | 1/4-28 Flat- | Bottom | 2,000 psi (138 bar) | Standard Knurl | Delrin Purple | 10-pk |
| XP-245x | Flangeless Fitting for 1/16" OD Tubing | 1/4-28 Flat- | Bottom | 2,000 psi (138 bar) | Standard Knurl | PFA Natural | 10-pk |
| XP-286x | Flangeless Fitting for 1/16" OD Tubing | 1/4-28 Flat- | Bottom | 2,000 psi (138 bar) | Headless Knurl | PPS Natural | 10-pk |
| REPLACE | MENT FERRULES FOR 1/16" OD TUBING | | | | | | |
| XP-200x | Flangeless Ferrule for 1/16" OD Tubing | 1/4-28 Flat- | Bottom | 2,000 psi (138 bar) | _ | ETFE Blue | 10-pk |
| XP-200Nx | Flangeless Ferrule for 1/16" OD Tubing | 1/4-28 Flat- | Bottom | 2,000 psi (138 bar) | _ | ETFE Natural | 10-pk |
| XP-240x | Flangeless Ferrule for 1/16" OD Tubing, Small Valve | 1/4-28 or 10 | -32 Flat-Bottom | 2,500 psi (172 bar) | _ | ETFE Natural | 10-pk |
| LANGEL | ESS FITTINGS FOR 1/8" OD TUBING (INCLUDE | S P-300 FERR | ULES) | | | | |
| XF-368x | Flangeless Fitting for 1/8" OD Tubing, FlushNut | 1/4-28 Flat- | Bottom | 500 psi (34 bar) | FlushNut | SST | 10-pl |
| KP-301x | Flangeless Fitting for 1/8" OD Tubing | 1/4-28 Flat- | Bottom | 500 psi (34 bar) | Standard Knurl | Delrin Black | 10-pl |
| (P-302x | Flangeless Fitting for 1/8" OD Tubing | 1/4-28 Flat- | Bottom | 500 psi (34 bar) | Standard Knurl | Delrin Red | 10-pl |
| (P-305x | Flangeless Fitting for 1/8" OD Tubing | 1/4-28 Flat- | Bottom | 500 psi (34 bar) | Standard Knurl | Delrin Green | 10-pl |
| (P-308x | Flangeless Fitting for 1/8" OD Tubing, Short | 1/4-28 Flat- | Bottom | 500 psi (34 bar) | Standard Knurl | Delrin Black | 10-pl |
| XP-315x | Flangeless Fitting for 1/8" OD Tubing | 1/4-28 Flat- | Bottom | 500 psi (34 bar) | Standard Knurl | ETFE Natural | 10-pl |
| XP-330x | Flangeless Fitting for 1/8" OD Tubing | 1/4-28 Flat- | Bottom | 500 psi (34 bar) | Standard Knurl | PEEK Natural | 10-pl |
| XP-335x | Flangeless Fitting for 1/8" OD Tubing, Short | 1/4-28 Flat- | Bottom | 500 psi (34 bar) | Standard Knurl | PEEK Natural | 10-pl |
| REPLACE | MENT FERRULES FOR 1/8" OD TUBING | | | | | | |
| XP-300x | Flangeless Ferrule for 1/8" OD Tubing | 1/4-28 Flat- | Bottom | 500 psi (34 bar) | _ | ETFE Yellow | 10-pk |
| XP-300Nx | Flangeless Ferrule for 1/8" OD Tubing | 1/4-28 Flat- | Bottom | 500 psi (34 bar) | _ | ETFE Natural | 10-pl |
| XP-340x | Flangeless Ferrule for 1/8" OD Tubing, Small Valve | 1/4-28 Flat- | Bottom | 500 psi (34 bar) | _ | ETFE Natural | 10-pl |
| KP-312x | Lock Nut for Flangeless Nuts | 1/4-28 Flat- | Bottom | _ | _ | Delrin White | 10-pl |
| | HROMalytic +61(0)3 976 | 2 2034 Im | stralian Distribut porters & Manufacur www.chromtech.net.a | rers | | | |
| Website N | EW: www.chromalytic.net.au E-mail: info@chromtech | net.au Tel: 03 9 | 762 2034 in Al | USTRALIA | | 47 (of | f 172) : |



MarvelX

- > Finger-tight to 19,000 psi
- > Reusable up to 200 times
- > Zero Dead Volume
- Biocompatible
- Limit Wear
- Increase Product Life
- Many IDs and Lengths Available



Learn more about MarvelX Next Generation UPHLC Fittings at www.idex-hs.com/MarvelX

1 NOTE


MarvelX tubing includes a sleeve that assists in product identification, with ID, length and part number information:



Minimum recommended bend-radius with MarvelX tubing is 1/4" (~6.35 mm).

MarvelX[™] Next Generation UHPLC Fittings

MarvelX UHPLC Connection Systems have been expertly designed for easy routing throughout your instrument, while providing consistent performance and superior re-usability. Built with convenient, removable stainless steel fittings and changeable, precision-cut flexible tubing, MarvelX can be used up to 200 times! The connection system is compatible with 10-32 coned receiving ports and is absolutely finger-tight — no tool required. MarvelX utilizes our unique next-generation patent-pending technology to auto-adjust to various port depths. This ensures zero dead volume and delivers better chromatography results with sleek, simple, and reliable functionality. In addition to our powerful Stainless Steel version, MarvelX offers a truly biocompatible option in PEEK-Lined Stainless Steel.





SPECIFICATIONS & DETAILS

Pressure Capability

19,000 psi (~1,310 bar) for routine use; up to 23,000 psi max over pressure for PEEK-Lined versions; up to 29,000 psi max over pressure for Stainless Steel versions.

Installation Method Finger-tight, 1/8 – 1/4 turn after initial resistance (~2 in-lbs)

Tubing Type 1/32" OD flexible 316 Stainless Steel with 1/16" OD rigid tube ends

Fitting Type 10-32 threaded, removable 316 Stainless Steel

Wetted Materials PEEK-Lined versions: PEEK | Stainless Steel versions: PEEK and 316 Stainless Steel


Maximum Use Temperature 120 °C

NOTE: The above performance specifications apply to use with appropriately-designed receiving ports under optimal conditions, using water at up to 120 °C for the testing process. If different conditions are used, the expected pressure threshold will be different.





- Route tubing to the target port.
- Slide fitting onto the tubing end via slot.
- Slowly finger-tighten to first resistance; continue tightening 1/8turn minimum, to 1/4-turn maximum.





Conventional coned fittings require a ferrule in conjunction with a fitting for proper sealing. They depend on complex techniques, including tools, to improve sealing performance, which significantly increases probability of extra internal volume and poor chromatography results. The excessive force needed for tightening increases wear of expensive components and the likelihood of replacement, adding to overall costs.



MarvelX UHPLC fittings do not depend on ferrules. They seal at the bottom of the port, without complex techniques, which significantly reduces required torque and enables many more connects and disconnects. Furthermore, they are virtually impossible to over-tighten by hand, limiting wear and increasing product life. An enhanced proprietary tip design also ensures zero dead volume (ZDV) and better chromatography results.



MarvelX

| | 070 MM | 150 MM | 250 MM | 350 MM | 500 MM | 600 MM |
|--------------|----------------------|--------------|--------------|--------------|--------------|--------------|
| Length | Part No. | Part No. | Part No. | Part No. | Part No. | Part No. |
| | TAINLESS STEEL ASSEM | BLIES* | | | | 1 20 2 1 2 2 |
| 25 μm ID | UPFP-6025070 | UPFP-6025150 | UPFP-6025250 | UPFP-6025350 | UPFP-6025500 | UPFP-6025600 |
| 50 μm ID | UPFP-6050070 | UPFP-6050150 | UPFP-6050250 | UPFP-6050350 | UPFP-6050500 | UPFP-6050600 |
| 75 µm ID | UPFP-6075070 | UPFP-6075150 | UPFP-6075250 | UPFP-6075350 | UPFP-6075500 | UPFP-6075600 |
| 100 μm ID | UPFP-6100070 | UPFP-6100150 | UPFP-6100250 | UPFP-6100350 | UPFP-6100500 | UPFP-6100600 |
| 150 µm ID | UPFP-6150070 | UPFP-6150150 | UPFP-6150250 | UPFP-6150350 | UPFP-6150500 | UPFP-6150600 |
| 300 μm ID | UPFP-6300070 | UPFP-6300150 | UPFP-6300250 | UPFP-6300350 | UPFP-6300500 | UPFP-6300600 |
| STAINLESS ST | EEL ASSEMBLIES* | | | | | |
| 100 μm ID | UPFS-6100070 | UPFS-6100150 | UPFS-6100250 | UPFS-6100350 | UPFS-6100500 | UPFS-6100600 |
| 125 µm ID | UPFS-6125070 | UPFS-6125150 | UPFS-6125250 | UPFS-6125350 | UPFS-6125500 | UPFS-6125600 |
| 254 μm ID | UPFS-6254070 | UPFS-6254150 | UPFS-6254250 | UPFS-6254350 | UPFS-6254500 | UPFS-6254600 |
| DEDLACEMEN | IT DADTC+ | | | | | |

To order Replacement Tubing, simply add the letter "T" to the end of any of the part numbers listed above. Example: UPFP-6020570T is the replacement tubing for UPFP-6020570.

Replacement Tubing Replacement Fittings UPN-61032 - Includes 3 replacement fittings.

*Product availability and lead times may vary depending on the configuration. Contact Customer Service at +1 800 426 0191 or email Customer Service.hs@idexcorp.com for details.





Fittings Tools

Tightening Tools for VHP & Other Fittings


- > Configured for the optimum torque to provide assurance of a strong connection
- > Prolongs the lifetime of reusable fittings by not overtightening
- Available for multiple fitting head styles

This new line of tightening tools is designed for the VHP fittings and can also be used with any fitting in this chapter described to have a corresponding head style to the tool listed below. There are three styles of tightening tools available for various applications. The Torque Tools (VHP-1000 and VHP-2000) are breakaway torque wrenches designed to deliver a precise amount of torque to the fitting system. These torque wrenches come calibrated according to ISO 6789:2003 (\pm 6% of setting) and have been tested extensively with the reusable VHP fittings on page 60. Choose the appropriate torque delivered and the proper head style to work with the VHP fittings, increasing the ease of use with these fittings.

The VHP-4000 Torque Driver couples with the specially designed Extender Tools listed below and provides an externally adjustable torque setting. This tool along with the appropriate Extender Tools will tighten any IDEX Health & Science knurled polymer fitting in your system. Reference the head style found in the tables at the bottom of each page for information on the proper Extender Tool to select.

Because of the small hex-head on the M4 fittings (VHP-900 and VHP-920), a custom wrench, the VHP-9000, is available in the table at the bottom of the page 51.


VHP-1000 VHP Wrench 1/4 in Hex 10 in-lbs (1.13 N·m)




Extender Tools

These tools can be used to tighten most of our knurled nuts in hard to reach places. See the application note on this page for knurl size and corresponding extender tool.


For precise tightening, the extender tools listed with 1/4" hex drives are designed to adapt to any torque wrench with a female 1/4" socket, such as the VHP-4000 Torque Driver on page 51. The tools featured on this page also include the FlushNut™ wrenches, used to tighten the FlushNuts found throughout this chapter and described in detail on page 51.





APPLICATION NOTE

The drawings represent actual size of the various knurled head designs of the IDEX Health & Science nuts featured in this chapter. Select the appropriate extender tool for the knurl pattern of the nut you've selected.







FLUIDICS





Removal Tool

Use the LT-300 Removal Tool to detach LiteTouch® and Super Flangeless™ Ferrules from tubing. Simply slide the appropriate tool blade slot between the lock ring and the ferrule body. With a slight twist, the ring will pop off, releasing the ferrule from the tubing. Please Note: This Removal Tool will not work with the LT-135x Ferrule System.

Wrenches

For your convenience, we offer wrenches in three standard sizes. You will need two A-304 wrenches to tighten most nuts into unions found on page 36 (for union 1593, you need one A-304 and one A-320 wrench).

The IDEX Wrench is slotted to fit over 1/16" and 1/8" OD tubing, and has 1/4" and 5/16" internal hex ends, to engage with the heads of the hex-head fittings most commonly used with IDEX Health & Science valves and the stainless steel fittings listed on page 59.

Fittings Tools

| Description | Use With Head Style | Torque Delivered | Qty. |
|--|--|--|--|
| NG TOOLS | | | |
| Extender Tool to Torque Driver | FlushNut (1/4-28) | _ | ea. |
| Extender Tool to Torque Driver | Micro Headless | _ | ea. |
| Extender Tool to Torque Driver | 1/4" Hex | _ | ea. |
| Extender Tool to Torque Driver | Female Nut Knurl | _ | ea. |
| Extender Tool to Torque Driver | Micro Nut Knurl | _ | ea. |
| Extender Tool to Torque Driver | Standard Nut Knurl | _ | ea. |
| Extender Tool to Torque Driver | Headless Nut Knurl | _ | ea. |
| Standard Knurl Torque Tool | Standard Knurl | 4 in-lbs (0.45 N·m) | ea. |
| VHP Torque Tool | 1/4" Hex | 10 in-lbs (1.13 N·m) | ea. |
| VHP Torque Tool | 1/4" Hex | 14 in-lbs (1.58 N·m) | ea. |
| VHP Torque Driver | Extender Tool 1/4" Drive | Adjustable between 2–12 in-lbs (0.23–1.35 N·m) | ea. |
| 4 mm Wrench | 4 mm Hex | _ | ea. |
| DLS | | | |
| Description | Material | | Qty. |
| Extender Tool for Standard Head Nuts, with 1/4" Hex Drive | Aluminum | | ea. |
| Extender Tool for Standard Head Nuts | Delrin® | | ea. |
| Extender Tool for Standard Head Nuts | Aluminum | | ea. |
| Extender Tool for Standard Head Nuts, Short | Aluminum | | ea. |
| Extender Tool for Headless Nuts | Aluminum | | ea. |
| Extender Tool for Headless Nuts, with 1/4" Hex Drive | Aluminum | | ea. |
| Extender Tool for Standard Micro Nuts | Aluminum | | ea. |
| Extender Tool for Micro Headless Nuts | Aluminum | | ea. |
| Extender Tool for Female Nuts, with 1/4" Hex Drive | Aluminum | | ea. |
| US TOOLS | | | |
| Wrench, 1/4" x 5/16" | Steel | | ea. |
| Wrench, 1/2" x 9/16" | Steel | | ea. |
| Wrench, 3/8" x 7/16" | Steel | | ea. |
| IDEX Wrench, 1/4" x 5/16" | Steel | | ea. |
| FlushNut Wrench for 10-32 Threaded Fittings | Steel/Plastic Handle | | ea. |
| FlushNut Wrench for 1/4-28 Threaded Fittings | Steel/Plastic Handle | | ea. |
| Removal Tool for LiteTouch and Super Flangeless Ferrules | Steel/Plastic Handle | | ea. |
| Swaging Tool for TinyTight Fittings, for 6–40 MINSTAC Port | SST | | ea. |
| | Extender Tool to Torque Driver Standard Knurl Torque Tool VHP Torque Tool VHP Torque Tool VHP Torque Tool VHP Torque Driver 4 mm Wrench OLS Description Extender Tool for Standard Head Nuts, with 1/4" Hex Drive Extender Tool for Standard Head Nuts Extender Tool for Standard Head Nuts Extender Tool for Standard Head Nuts, Short Extender Tool for Standard Head Nuts, Short Extender Tool for Headless Nuts Extender Tool for Headless Nuts, with 1/4" Hex Drive Extender Tool for Standard Micro Nuts Extender Tool for Standard Micro Nuts Extender Tool for Female Nuts, with 1/4" Hex Drive US TOOLS Wrench, 1/4" x 5/16" Wrench, 3/8" x 7/16" IDEX Wrench, 1/4" x 5/16" FlushNut Wrench for 10-32 Threaded Fittings FlushNut Wrench for 1/4-28 Threaded Fittings Removal Tool for LiteTouch and Super Flangeless Ferrules | Extender Tool to Torque Driver Extender Tool Torque Tool VHP Torque Tool VHP Torque Tool VHP Torque Tool VHP Torque Driver Extender Tool 1/4" Hex VHP Torque Driver Extender Tool 1/4" Drive Extender Tool for Standard Head Nuts, with 1/4" Hex Drive Extender Tool for Standard Head Nuts Extender Tool for Standard Head Nuts Extender Tool for Standard Head Nuts, Short Extender Tool for Standard Head Nuts, Short Extender Tool for Headless Nuts Extender Tool for Headless Nuts Extender Tool for Headless Nuts Extender Tool for Standard Micro Nuts Extender Tool for Standard Micro Nuts Extender Tool for Standard Micro Nuts Extender Tool for Female Nuts, with 1/4" Hex Drive Aluminum Extender Tool for Standard Micro Nuts Extender Tool for Female Nuts, with 1/4" Hex Drive Aluminum Extender Tool for Female Nuts, with 1/4" Hex Drive Aluminum Extender Tool for Standard Micro Nuts Extender Tool for Standard Shirt Steel Wrench, 1/4" x 5/16" Steel Wrench, 1/4" x 5/16" Steel Wrench, 1/4" x 5/16" Steel FlushNut Wrench for 10-32 Threaded Fittings Steel/Plastic Handle FlushNut Wrench for 1/4-28 Threaded Fittings Steel/Plastic Handle | Extender Tool to Torque Driver Extender Tool Torque Driver Extender Tool Torque Driver Extender Tool Torque Tool VHP Torque Tool VHP Torque Tool VHP Torque Driver Arm Wench Arm Hex Description Extender Tool for Standard Head Nuts, with 1/4" Hex Drive Extender Tool for Standard Head Nuts Extender Tool for Female Nuts, with 1/4" Hex Drive Extender Tool for Female Nuts, with 1/4" Hex Drive Extender Tool for Female Nuts, with 1/4" Hex Drive Extender Tool for Female Nuts, with 1/4" Hex Drive Extender Tool for Female Nuts, with 1/4" Hex Drive Extender Tool for Standard Micro Nuts Extender Tool for Standard Micro Nuts Extender Tool for Standard Head Nuts Extender Tool for Headless Nuts Extender Tool for Headles |



Tubing Sleeves



MicroTight® Tubing Sleeves

- Manufactured from PEEK polymer
- > Pressure rated to 4,000 psi (276 bar)
- > Color-coded for easy inner diameter identification

IDEX Health & Science MicroTight Tubing Sleeves feature an outer diameter of 0.025" and offer a wide assortment of inner diameters to help facilitate capillary tubing connections with our MicroTight accessories. Because the sleeves are manufactured from PEEK polymer, they carry an upper temperature threshold of 125 °C.

To use these sleeves properly, choose a sleeve with an inner diameter 0.001"-0.002" (25–50 µm) larger than the outer diameter of your capillary tubing. Then, slip the sleeve over your flow path tubing, such that your tubing extends all the way through the sleeve, but not beyond the end of the sleeve. Choose the correct fitting that corresponds with your receiving port, slide it over the sleeved flow path tubing and connect as normal.



NanoTight[™] Tubing Sleeves

- Manufactured from FEP fluoropolymer
- Pressure rated to 4,000 psi (276 bar)
- ▶ Outer diameter of 1/16" the most popular size used on most instrumentation

NanoTight Tubing Sleeves are manufactured using FEP fluoropolymer and precisely cut to a 1.6" length. A wide assortment of sleeves is available, ensuring the availability of a NanoTight sleeve for most applications. Many of the sleeves feature a light color tint that can help more easily identify the inner diameter for future orders. Because FEP is the base polymer for these sleeves, there is a maximum recommended continuous operating temperature of 50 $^{\circ}$ C.

Our NanoTight sleeves were designed primarily for use with the NanoTight fittings, found on page 37 and also work well with the Super Flangeless™ fittings for 1/16" OD tubing on page 39. For tubing sleeves that can be used effectively with stainless steel fittings and at higher temperatures, consider using the PEEK Tubing Sleeves, found below.



1/16" OD PEEK Tubing Sleeves

- For connecting capillary tubing to standard 10-32 ports
- > Require the use of wrench tightened stainless steel nuts
- Pressure rated to 6,000 psi (414 bar)

Like the NanoTight™ FEP Sleeves on the previous page, these PEEK Tubing Sleeves are designed to be used with 1/16" OD, 10-32 threaded fittings to adapt capillary tubing to standard coned ports. Made of PEEK polymer, these 1.3" long sleeves can be used up to 125 °C.

These sleeves require a wrench tightened nut to achieve proper sealing. We recommend our SealTight™ fittings on page 36. Many researchers also use a stainless steel nut and ferrule with these sleeves, such as our U-400 and U-401 (page 33).



FLUIDICS



1/32" OD PEEK Tubing Sleeves

These 1.6" long 1/32" OD PEEK Tubing Sleeves can be used with any fitting designed for 1/32" OD tubing when smaller tubing must be connected. Select the appropriate sleeve from the product listing for your capillary tubing OD size. The 1/32" OD PEEK Tubing Sleeves have

a maximum recommended temperature of 125 °C and have a pressure rating of 5,000 psi (345 bar).



These 1.6" long sleeves facilitate connecting capillary tubing into ports designed for 1/32" OD tubing. Please refer to the product listing below to select the appropriate sleeve for your capillary OD size. These sleeves can be used at up to 50 °C and have a pressure rating of 1,750 psi (121 bar).

Clockwise, starting at top:



- ▶ 1/32" OD FEP Tubing Sleeves, shown with F-126Sx Fitting
- Fittings and tubing only shown to highlight how sleeves are designed to be used; they are not included with the sleeves



Why use Sleeves?

Because most capillary tubing connections are made into coned receiving ports, where the port is not designed to be used with capillary tubing directly, special care must be used to ensure a good connection. While custom ferrules can help make these connections, they only offer a fixed-length nose — and because most tubing pockets will vary slightly in length, this can lead to leaking or dead volume.

To help save overall expense while maintaining a concentric connection with minimal dead volume, IDEX Health & Science recommends the use of sleeves. Because sleeves are not permanently attached to a ferrule, they can easily adapt to varying tubing pocket depths. Additionally, because they are manufactured using extruded polymer tubing, you are assured of the concentricity of the resultant connection.



Tubing Sleeves (Cont.)



Use 1/32" OD PEEK or FEP Sleeves to connect capillary tubing with the following:

- The F-113 Ferrule and Two-Piece Fingertight Fittings for 10-32 ports (page 35).
- The F-112 and P-416BLK MicroTight® Fittings (page 34) 1/32" OD PEEK Tubing Sleeves only.
- ▶ The 1/32" OD MicroTight Fittings on page 34.
- The RheFlex M4 Fitting (page 61) for MX Module applications; the M-645 Valco®-Compatible Fitting (page 33) for Valco Nanovolume® valve applications.

Tubing Sleeves

| Part No. | ID | | For Tubing OD Size | Color | Qty |
|--------------|---|---|--------------------|---------|-----|
| MICROTIG | HT PEEK TUBING SLEEVES AND KITS, 0.025" OD |) | | | |
| F-180 | 125 μm (0.005") | | 70–110 μm | Red | ea. |
| -181 | 180 μm (0.007") | | 125–165 μm | Yellow | ea. |
| -182 | 230 μm (0.009") | | 175–215 μm | Natural | ea. |
| -183 | 280 μm (0.011") | | 225–265 μm | Blue | ea. |
| -184 | 330 μm (0.013") | | 275–315 μm | Orange | ea. |
| -185 | 395 μm (0.0155") | | 340–380 µm | Green | ea. |
| -186 | 455 μm (0.018") | | 400–440 μm | Black | ea. |
| -187 | 535 μm (0.021") | | 480–520 μm | Natural | ea. |
| -188 | 152 μm (0.006") | | 95–135 μm | Purple | ea. |
| 328 | MicroTight Tubing Sleeve Kit, contains (6) each of the sle | eeve sizes listed above | | _ | ea. |
| 1356 | MicroTight Connector Kit, contains: a 10-pack of each N MicroTight Adapters; and (2) MicroTight P-720 Unions | ficroTight Tubing Sleeve (F-180–F-187); (2) P-770 | _ | _ | ea. |
| NANOTIG | HT FEP TUBING SLEEVES, 1/16" OD | | | | |
| -237 | 125 µm (0.005") | | 70–110 μm | Red | ea. |
| -238 | 180 μm (0.007") | | 125–165 μm | Yellow | ea. |
| -239 | 215 µm (0.0085") | | 160–200 μm | Natural | ea. |
| 240 | 280 μm (0.011") | | 225–265 μm | Blue | ea. |
| 241 | 330 μm (0.013") | | 275–315 μm | Orange | ea. |
| 242 | 395 μm (0.0155") | | 340–380 µm | Green | ea. |
| 243 | 455 µm (0.018") | | 400–440 μm | Black | ea. |
| 244 | 535 µm (0.021") | | 480–520 μm | Natural | ea. |
| 245 | 610 µm (0.024") | | 555–595 μm | Red | ea. |
| 246 | 685 µm (0.027") | | 630–670 µm | Yellow | ea. |
| 247 | 840 µm (0.033") | | 785–825 µm | Green | ea. |
| 252 | 1.07 mm (0.042") | | 1 mm | Purple | ea. |
| EEK TUBI | NG SLEEVES FOR 1/16" OD FITTINGS | | | | |
| -225 | 125 µm (0.005") | | 70–110 μm | Red | ea. |
| 226 | 180 µm (0.007") | | 125–165 µm | Yellow | ea. |
| 227 | 230 µm (0.009") | | 175–215 μm | Yellow | ea. |
| 228 | 250 μm (0.011") | | 225–265 μm | Blue | ea. |
| -229 | 330 µm (0.013") | | 275–315 μm | Natural | ea. |
| -230 | 405 μm (0.016") | | 350–390 μm | Orange | ea. |
| 231 | 560 µm (0.022") | | 505–545 μm | Natural | ea. |
| 232 | 785 μm (0.031") | | 730–770 µm | Natural | ea. |
| 233 | 865 µm (0.034") | | 785–825 µm | Blue | ea. |
| 234 | 685 µm (0.027") | | 630–670 µm | Yellow | ea. |
| | NG SLEEVES FOR 1/32" OD FITTINGS | | 030 070 µ111 | TCHOW | cu. |
| -381 | 180 μm (0.007") | | 125–165 µm | Yellow | ea. |
| -382 | 205 μm (0.008") | | 150–190 μm | Natural | |
| | | | | | ea. |
| -384 -385 | 255 μm (0.010") | | 200–240 µm | Blue | ea. |
| | 380 μm (0.015") | | 325–365 μm | Natural | ea. |
| 386 | 510 μm (0.020") | | 455–495 μm | Orange | ea. |
| 387 | 250 µm (0.011") | | 225–265 µm | Red | ea. |
| -388 | 330 μm (0.013") | | 275–315 μm | Black | ea. |
| | IG SLEEVES FOR 1/32" OD FITTINGS | | | | |
| -374 | 280 μm (0.011") | | 225–265 μm | Blue | ea. |
| -375 | 330 µm (0.013") | | 275–315 μm | Orange | ea. |
| -376 | 395 um (0.0155") | | 340–380 μm | Green | ea. |


FLUIDICS



Plugs & Caps

Seal 6-32, 6-40, 10-32, 1/4-28, M6, or 5/16-24 threaded ports or fittings Use our plugs to close off unused ports in valves and multi-port connectors. Our color-coded 10-32 threaded plugs are perfect for identifying stored columns that have different packing materials, or in which different mobile phases have been utilized. Cap off tubing with one of the PEEK or ETFE caps presented on this page and the appropriate fittings from this chapter.


To help determine which plug or cap is best suited for your application, please visit www.idex-hs.com for detailed chemical compatibility data.



P-316 PFA Plug for 1/4-28 flat-bottom ports



U-467RDelrin® Column Plug
for 10-32 coned ports



P-755 ETFE Cap for 1/4-28 flat-bottom fittings



| Part No. | Description | Head Style | Material | Qty. |
|----------|--|--|---------------|------|
| PLUGS | | | | |
| P-120 | Plug for 1/4-28 Coned Ports for 1/8" OD Tubing | Standard Knurl | PCTFE Natural | ea. |
| P-123 | Plug for 1/4-28 Flat-Bottom Ports | 5/16" Hex | ETFE Natural | ea. |
| P-309 | Plug for 1/4-28 Flat-Bottom Ports | Standard Knurl | Delrin Black | ea. |
| P-311 | Plug for 1/4-28 Flat-Bottom Ports | Standard Knurl | ETFE Natural | ea. |
| P-314 | Plug for M6 Flat-Bottom Ports | Standard Knurl | ETFE Black | ea. |
| P-316 | Plug for 1/4-28 Flat-Bottom Ports | Standard Knurl | PFA Natural | ea. |
| P-321 | Plug for 1/4-28 Flat-Bottom Ports, FlushNut™ | FlushNut | PEEK Natural | ea. |
| P-520 | Plug for 10-32 Coned Ports | 5/16" Hex | SST | ea. |
| P-550 | Plug for 10-32 Coned Ports, Extra Long | Standard Knurl | PEEK Natural | ea. |
| P-551 | Plug for 10-32 Coned Ports | Standard Knurl | PEEK Natural | ea. |
| P-555 | Plug for 6-32 Coned Ports | Standard Micro Knurl | PEEK Natural | ea. |
| P-556 | Plug for 5/16-24 Flat-Bottom Ports | Standard Knurl | PEEK Natural | ea. |
| P-558 | Plug for 6-40 Flat-Botton Ports | Micro Headless Knurl | PEEK Green | ea. |
| P-849 | Plug for 10-32 Flat-Bottom Ports | Standard Knurl | Delrin Black | ea. |
| U-467R | Plug for 10-32 Coned Ports | Standard Knurl | Delrin Red | ea. |
| VHP-600 | VHP Plug for 10-32 Coned Ports | 3/8" Hex | PK-SST | ea. |
| CAPS | | | | |
| P-754 | Cap for 10-32 Coned Ports | Standard Knurl | ETFE Yellow | ea. |
| P-755 | Cap for 1/4-28 Flat-Bottom Ports | Standard Knurl | ETFE Black | ea. |
| P-756 | Cap for M6 Flat-Bottom Ports | Standard Knurl | ETFE Blue | ea. |
| | HROMalytic +61(0)3 9762 2034 | Australian Distributors Importers & Manufacurers | | |

ECH no logy Pty Ltd

Website NEW: www.chromalytic.net.au E-mail: info@chromtech.net.au Tel: 03 9762 2034 . . . in AUSTRALIA

Large Bore Fittings

- > 5/16-24 or 1/2-20 threads
- For use with 1/16", 1/8", 3/16", 1/4", 5/16", 3.0 mm, or 4.0 mm OD tubing

NOTE

Each of the Large Bore Fittings shown on this page comes in a convenient 10-pack and is packaged with the most popularly chosen Ferrule option. The Fittings can be ordered separately by removing the preceding letter "X" from the part number. Additionally, to connect metric-sized tubing with outer diameters less than 4.0 mm to 5/16-24 threaded ports, reference the chart on page 43 to choose the correct nut/ ferrule combination.



RELATED PRODUCTS

| | Page |
|--------------------------|---------|
| MORE LARGE BORE PRODUCTS | |
| 5/16-24 Coned Fittings | 41 |
| Barbed Adapters | 90 |
| Threaded Adapters | 63 |
| Plugs | 55, 112 |
| Y Connector | 90 |
| ETFE Tubing | 27 |
| FEP, PFA Tubing | 26, 24 |
| Tubing Cutter | 28 |
| Pressure Relief Valve | 148 |
| Solvent Inlet Filters | 102 |
| Bottle Caps | 55 |
| Semi-Prep Inline Filters | 107 |
| | |



XP-130x

PEEK Nut, for 1/8" OD tubing shown with P-300 Flangeless Ferrule (page 43)



XP-136x

0.634

(1.60 cm)

PEEK Nut, for 1/16" OD tubing shown with P-200 Flangeless Ferrule (included and found on page 43)





XU-620x

PEEK Nut, for 1/4" OD tubing shown with inverted U-650 Flangeless Ferrule for coned ports (included and found on this page)

0.88

XP-131x

PEEK Nut, for 1/8" OD tubing shown with P-359 Super Flangeless™ Ferrule (included and found on page 41)



PEEK Nut, for 1/16" OD tubing shown with P-259 Super Flangeless Ferrule (included and found on page 41) PEEK Nut, for 3/16" OD tubing shown with P-140 Super Flangeless Ferrule



0.88" (2.24 cm)

0.50" (1.27 cm)

0.88

(2.24 cm)

0.16" /(0.41 cm)

XP-132x

(0.94 cm)

XP-141x

PEEK Nut, for 3/16" and 4.0 mm OD tubing shown with P-133 Flangeless Ferrule (included and found on this page)

5/16-24

XU-655x

0.91

PEEK Nut, for 1/4" OD tubing shown with U-650 Flangeless Ferrule (included and found on this page)

XU-662x PEEK Nut, for 5/16" OD tubing shown with U-660 Flangeless Ferrule (included and found on this page)

| Part No. | Description | Port | Pressure Rating | Head Style | Material (Nut/Washer) | Qty. |
|----------|--|---------------------|---------------------|----------------|------------------------------|-------|
| LARGE BO | ORE FITTINGS | | | | | |
| XP-130x | Flangeless Fitting for 1/8" OD tubing | 5/16-24 Flat-Bottom | 500 psi (34 bar) | Standard Knurl | PEEK Natural/ETFE Yellow | 10-pk |
| XP-131x | Super Flangeless Fitting for 1/8" OD tubing | 5/16-24 Flat-Bottom | 1,000 psi (69 bar) | Standard Knurl | PEEK Natural/ETFE Yellow/SST | 10-pk |
| XP-132x | Flangeless Fitting for 3/16" OD tubing | 5/16-24 Flat-Bottom | 500 psi (34 bar) | Standard Knurl | PEEK Natural/ETFE Blue | 10-pk |
| XP-136x | Flangeless Fitting for 1/16" OD tubing | 5/16-24 Flat-Bottom | 2,000 psi (138 bar) | Standard Knurl | PEEK Natural/ETFE Blue | 10-pk |
| XP-137x | Super Flangeless Fitting for 3/16" OD tubing | 5/16-24 Flat-Bottom | 500 psi (34 bar) | Standard Knurl | PEEK Black/ETFE Green/SST | 10-pk |
| XP-141x | Super Flangeless Fitting for 1/16" OD tubing | 5/16-24 Flat-Bottom | 1,350 psi (93 bar) | Standard Knurl | PEEK Natural/ETFE Yellow/SST | 10-pk |
| XU-620x | Flangeless Fitting for 1/4" OD tubing | 1/2-20 Coned | 250 psi (17 bar) | Large Knurl | PEEK Red/ETFE Natural | 10-pk |
| XU-655x | Flangeless Fitting for 1/4" OD tubing | 1/2-20 Flat-Bottom | 250 psi (17 bar) | Large Knurl | PEEK Black/ETFE Natural | 10-pk |
| XU-662x | Flangeless Fitting for 5/16" OD tubing | 1/2-20 Flat-Bottom | 250 psi (17 bar) | Large Knurl | PEEK Black/ETFE Natural | 10-pk |
| REPLACE | MENT FERRULES | | | | | |
| P-133x | Flangeless Ferrule for 3/16" OD tubing | 5/16-24 Flat-Bottom | 500 psi (34 bar) | _ | ETFE Blue | 10-pk |
| P-133Nx | Flangeless Ferrule for 3/16" OD tubing | 5/16-24 Flat-Bottom | 500 psi (34 bar) | _ | ETFE Natural | 10-pk |
| P-140x | Super Flangeless Ferrule for 3/16" OD tubing | 5/16-24 Flat-Bottom | 500 psi (34 bar) | _ | ETFE Green | 10-pk |
| U-650x | Flangeless Ferrule for 1/4" OD tubing | 1/2-20 Flat-Bottom | 250 psi (17 bar) | _ | ETFE Natural | 10-pk |
| U-660x | Flangeless Ferrule for 5/16" OD tubing | 1/2-20 Flat-Bottom | 250 psi (17 bar) | _ | ETFE Natural | 10-pk |

FLUIDICS



VHP Micro Fittings

Micro Fittings are specifically designed for use with microferrules. They are manufactured from a proprietary PEEK blend (PK) which allow them to be used at higher temperatures (up to 200° C) and higher pressures ideal for UHPLC applications.

VHP MicroFerrules and Fittings are made from a proprietary high performance PEEK polymer blend, a material which is unique in its ability to enable the use of capillary tubing in UHPLC environments. The new high pressure MicroFerrules are available for use with 1/32" or 360 μm OD tubing, and they are incorporated into several of our VHP products for capillary tubing.





Caution: While the proprietary blend of the PK fittings will allow a fitting to attain a higher pressure and minimal cold flow properties relative to pure PEEK, some fittings molded of PK are known to be conductive. Use caution when employing PK fittings in high voltage applications.



MicroTight fittings and MicroFerrules

While the MicroTight Female Nuts may be used with any of the separate MicroFerrules, the MicroFerrules themselves are port-specific and are thus not interchangeable. Additionally, the one-piece MicroTight fittings are also port-specific and should not be exchanged.

| art No. | Description | Port | Pressure Rating | Required Torque | Head Style | Material | Qty. |
|-----------|---------------------------------------|-----------------|---|-----------------------|--------------|---------------|-------------|
| K MICRO I | FERRULES AND FEMALE NUTS | | | | | | |
| -416 | Female Nut for Microferrule | 5/16-24 Coned | 15,000 psi (1,035 bar) | 4.0 in-lbs (0.45 N·m) | Female Knurl | PEEK, Natural | ea. |
| -416BLK | Female Nut for Microferrule | 5/16-24 Coned | 15,000 psi (1,035 bar) | 4.0 in-lbs (0.45 N·m) | Female Knurl | PEEK, Black | ea. |
| K-112 | VHP MicroFerrule for 1/32" OD Tubing | 5/16-24 Coned | 15,000 psi (1,035 bar) | _ | _ | PK | ea. |
| K-152 | VHP MicroFerrule for 360 um OD Tubina | 5/16-24 Coned | 15.000 psi (1.035 bar) | _ | _ | PK | ea. |
| | HROMalytic +61(0)3 | 5 9/62 2034 Imj | stralian Distributors porters & Manufacurers www.chromtech.net.au | | | | |
| | ECHnolog | | | | | | 7 (- 5 4 7 |

Website NEW: www.chromalytic.net.au E-mail: info@chromtech.net.au Tel: 03 9762 2034... in AUSTRALIA



VHP PK Fittings



> Find unions, tees and crosses for VHP applications on page 75, and 74.

Ultra High Performance fittings are manufactured from a proprietary PEEK blend (PK) which allow them to be used at higher temperatures (up to 200 °C) and higher pressures.

The VHP PK One-Piece fittings are available for 10-32 coned, 6-32 coned, or M4 coned ports, and Two-Piece fittings are available to connect either 1/16" or 1/32" OD tubing into 10-32 coned ports in multiple styles.









PK-110x 10-32 PK Fitting shown with PK-132x for 1/32" OD tubing



PK-126Hx 6-32 PK MicroTight® Fitting for 1/32" OD tubing



PK-126x 6-32 PK MicroTight® Fitting for 1/32" OD tubing



UH-904x M4, 1/32" Fitting for IDEX Health & Science

| Part No. | Description Port | | Pressure Rating | Required Torque | Head Style | Material | Qty. |
|------------|---|-------------|------------------------|-----------------------|----------------------|----------|-------|
| PK VHP ON | E-PIECE FITTINGS | | | | | | |
| PK-120BLKx | PK One-Piece Fitting for 1/16" OD Tubing | 10-32 Coned | 12,000 psi (827 bar) | 8.0 in-lbs (0.90 N·m) | Standard Knurl | PK | 10-pk |
| PK-126Hx | PK One-Piece Headless Fitting for 1/32" OD Tubing | 6-32 Coned | 15,000 psi (1,035 bar) | 3.0 in-lbs (0.34 N·m) | Headless Micro Knurl | PK | 10-pk |
| PK-126x | PK One-Piece Fitting for 1/32" OD Tubing | 6-32 Coned | 15,000 psi (1,035 bar) | 3.0 in-lbs (0.34 N·m) | Standard Micro Knurl | PK | 10-pk |
| UH-904x | PK One-Piece Fitting for 1/32" OD Tubing | M4 Coned | 15,000 psi (1,035 bar) | 4.0 in-lbs (0.45 N·m) | Headless Knurl | PK | 10-pk |
| PK VHP FIT | TINGS (LITETOUCH® STYLE, NUTS AND FER | RULES SOLD | SEPARATELY) | | | | |
| PK-100x | PK Ferrule for 1/16" OD Tubing | 10-32 Coned | 16,500 psi (1,140 bar) | _ | _ | PK | 10-pk |
| PK-110x | PK Nut for 1/16" OD Tubing | 10-32 Coned | 16,500 psi (1,140 bar) | 8.0 in-lbs (0.90 N·m) | Standard Knurl | PK | 10-pk |
| PK-132x | PK Ferrule for 1/32" OD Tubina | 10-32 Coned | 16.500 psi (1.140 bar) | _ | _ | PK | 10-pk |
| | HROMalytic +61(0)3 97 | b2 2U34 | ustralian Distributors | | | | |

FLUIDICS



Stainless Steel VHP Fittings

- Pressure rated to 30,000 psi (2,070 bar)
- Double compression ferrule design
- Available with 10-32 threads for 1/16" OD tubing and M4 threads for 1/32" OD tubing



APPLICATION NOTE

In order to seal up to the stated pressure rating, the VHP-200-01 ferrule requires 20 in-lbs (2.25 N·m) of torque. Similar ferrules on the market require tightening torque of at least 30 in-lbs (3.3 N·m), which can result in a restricted tubing passage, as shown in the picture below. This restriction can increase turbulence and add a 'throttling' effect to the fluid pathway, resulting in mixing and other potential chromatographic problems.

IDEX Health & Science VHP-200



Uniform Tubing Passage

Constricted Tubing Passage

The all Stainless-Steel VHP Fittings include a unique ferrule system with two compression points to provide twice the grip of a standard ferrule. This design also allows the bite on the tubing to be less concentrated and does not restrict the inner diameter, as discussed in the Application Note. The ferrules for 1/16" OD tubing and 10-32 coned ports are two pieces, while the grooved ferrule for 1/32" OD tubing and M4 coned ports is a one-piece design for easier handling, but it will act as two pieces with double compression on the tubing as it is tightened down.



VHP-200x VHP 10-32 Fitting for 1/16" OD tubing



VHP-700x VHP 6-40 Fitting for 1/32" OD tubing

| Part No. | Description | Port | Pressure Rating | Required Torque | Head Style | Material | Qty. | | | | |
|---|--------------------------|-------------|-------------------------|----------------------|------------|----------|-------|--|--|--|--|
| STAINLESS STEEL VHP FITTINGS (INCLUDES NUT AND FERRULE) | | | | | | | | | | | |
| VHP-200x | VHP Fitting for 1/16" OD | 10-32 Coned | 30,000 psi (2,070 bar) | 20 in-lbs (2.25 N·m) | 1/4" Hex | SST | 10-pk | | | | |
| VHP-700x | VHP Fitting for 1/32" OD | 6-40 Coned | 30,000 psi (2,070 bar) | 20 in-lbs (2.25 N·m) | 4 mm Hex | SST | 10-pk | | | | |
| STAINLESS S | TEEL VHP FERRULES | | | | | | | | | | |
| VHP-200-01x | VHP Ferrule for 1/16" OD | 10-32 Coned | 30.000 psi (2.070 bar) | 20 in-lbs (2.25 N·m) | _ | SST | 10-pk | | | | |
| | TIPONE - | | Australian Distributors | | | | | | | | |





Reusable VHP Fittings

- > Pressure rated up to 25,000 psi (1,720 bar)
- > Patented innovative design
- Capable of up to ten repeat assembly cycles with no impact on pressure holding ability or carry-over
- > Available in 10-32 threads for 1/16" OD tubing and M4 threads for 1/32" OD tubing
- > Materials of construction: stainless steel and proprietary PEEK polymer blend (PK)
- Quick component replacement, minimal downtime



> Find tightening tools on page 50 designed to deliver the torque necessary for these fittings.

IDEX Health & Science introduces an innovative line of Very High Pressure (VHP) fittings, designed to withstand extreme pressures. This patented line of ground-breaking fitting systems is perfect for use within the increasingly demanding requirements of today's high performance analytical systems.

The Reusable VHP fittings can be reused when following the tightening torque specification listed below. With a polymer front ferrule, there is no damage to the tubing or receiving port, also increasing the life of these components.



10-32 VHP Fitting for 1/16" OD tubing




10-32 VHP Fitting, Long for 1/16" OD tubing



for 1/32" OD tubing









Assorted Fittings Kits

0.25" (0.64 cm) 0.89" (2.26 cm) 6000-360 M4 Fitting M4 threads for 1/32" OD tubing

0.45" (1.14 cm) 0.25" (0.64 cm) 0.95" (2.41 cm) 6000-254 10-32 PEEK Nut with 6000-251 PEEK Ferrule (RD ChromTRAC suffix)



RELATED PRODUCTS

- For PEEK tubing sleeves that can be used with these M4 RheFlex fittings, see page 52.
- For reusable fittings that both work in UHPLC applications and can help ensure the tubing is fully inserted into the receiving port, see page 60.

RheFlex® M4 Fittings

- Incorporates M4 coned threads for 1/32" OD tubing
- Pressure rated to 5,000 psi (345 bar)

Our RheFlex M4 Fitting is designed to connect 1/32" OD tubing in MX Series II™ valves (see Actuated Valves, starting on page 119). This PEEK fitting has a one piece design, which eliminates the need for a separate nut and ferrule. The M4 Fitting design provides dependable zero dead volume connections for micro and nano applications. Due to the unique RheFlex gripping design, the M4 Fitting will hold to 5,000 psi (345 bar) on PEEK or with a PEEK tubing sleeve on fused silica tubing. A PEEK M4 Plug is also available.

Use ChromTRAC™ knobs with the RheFlex M4 Fitting for fingertight convenience and to color-code connections.

Two-Piece RheFlex Fingertight Fittings

The RheFlex Precision Two-Piece PEEK Fittings sets provide inert, biocompatible connections for instrumentation. These fittings have a reliable, time-tested design. Each 1/16" fittings set contains a 10-32 threaded nut and a specially-designed PEEK ferrule. Three lengths of the 1/16" nut are available: Standard, Short, and Extra Long. RheFlex Fingertight Fittings are rated for use up to 7,000 psi (483 bar). Also offered in this product line is the 6000-078 fitting, designed to connect 1/8" OD tubing into our manual preparative-scale injection valves. (See page 115 for more information on these valves.) View the online product bulletin at: www.idex-hs.com.

ChromTRAC™

▶ Brightly colored knobs to easily track inlets and outlets of valves, columns, and detectors

All ChromTRAC-compatible RheFlex fittings offer the ChromTRAC knob option. Specify the ChromTRAC two letter suffix for the color choice when ordering. Please see the ChromTRAC Suffix Codes table below. For example, to order red ChromTRAC knobs with the RheFlex One-Piece Fitting on this page, specify 6000-282RD. No suffix indicates black knobs. View the online product bulletin for RheFlex fittings at: www.idex-hs.com.

ThromTRAC Suffice Codes

Website NEW: www.chromalytic.net.au E-mail: info@chromtech.net.au Tel: 03 9762 2034 . . . in AUSTRALIA

| CIIIOIIIIIAC | Julia Codes | | |
|--------------|-------------|------|--|
| CODE | COLOR | CODE | COLOR |
| BL | Blue | WH | White |
| GN | Green | YL | Yellow |
| GY | Gray | MC | Multi-color (two each of blue, green, gray, red, and yellow) |
| RD | Red | | |

Add these letter suffixes to the end of the seven-digit part numbers of the 10-32 and M4 threaded RheFlex Fittings listed below.

| э. | Description | Port | Pressure Rating | Head Style | Material | Qty. |
|---------------------------|--|--|---|----------------|---------------|---------|
| LEX ONE-PIECE | FITTINGS | | | | | |
| 60 | RheFlex Fitting for 1/32" OD Tubing | M4 Coned | 5,000 (345 bar) | 1/4" Hex | PEEK, Natural | 10-pk |
| EX TWO-PIECE | FITTINGS (INCLUDES FERRULES) | | | | | |
| 78 | RheFlex Fitting for 1/8" OD Tubing | 5/16-24 Coned | 5,000 psi (345 bar) | 5/16" Hex | PEEK, Natural | ea. |
| 54 | RheFlex Fitting for 1/16" OD Tubing | 10-32 Coned | 7,000 psi (483 bar) | ChromTRAC knob | PEEK, Natural | 10-pk |
| 55 | RheFlex Fitting for 1/16" OD Tubing, Short | 10-32 Coned | 7,000 psi (483 bar) | 1/4" Hex | PEEK, Natural | 10-pk |
| ACEMENT FERR | ULES | | | | | |
| 79 | RheFlex Ferrule for 1/8" OD Tubing | 5/16-24 Coned | 7,000 psi (483 bar) | ChromTRAC knob | PEEK, Natural | 5-pk |
| 51 | RheFlex Ferrule for 1/16" OD Tubina | 10-32 Coned | 7.000 psi (483 bar) | ChromTRAC knob | PEEK, Natural | 10-pk |
| HR | ECH no logy Pty Ltd | Importers & I | Manufacurers | | | |
| ACEMENT FERRI 79 51 | ULES RheFlex Ferrule for 1/8" OD Tubing RheFlex Ferrule for 1/14" OD Tubina **Political Service** +61(0)3 9762 2034 | 5/16-24 Coned 10-32 Coned Australian I Importers & II | 7,000 psi (483 bar) 7,000 psi (483 bar) Distributors Manufacurers | ChromTRAC knob | PEEK, N | Natural |





CONNECTORS

Connectors are designed to securely join tubing together or to facilitate the joining of tubing to other fluid pathway components. We offer multiport connectors with different thread and port configurations to meet your system requirements and connection needs. Some of our connectors feature a True ZDV (Zero Dead Volume) internal configuration that helps minimize the formation of dead volume in your fluidic pathway. Our versatile adapters help bring two connectors with different configurations together. Connectors are manufactured from 316 stainless steel or from inert polymers to ensure chemical compatibility with the fluid passing through. Peristaltic tube connectors are ideal for making connections with soft-walled, peristaltic tubing. Our extensive line of connectors includes tees, crosses, Luer Adapters, barbed and threaded adapters, and a variety of other options.

| 63 | THREADED ADAPTERS | 79 | ACCESSORIES |
|----|----------------------------------|-----------|-------------------------|
| 67 | HIGH PRESSURE MULTIPORT | 80 | NANOPORT ASSEMBLIES |
| | CONNECTORS | 81 | LOW PRESSURE |
| 72 | ULTRA HIGH PRESSURE MULTIPORT | | MULTIPORT CONNECTORS |
| | CONNECTORS | 90 | LUER ADAPTERS |
| 77 | MICROTIGHT® ADAPTERS | 91 | PERISTALTIC TUBE |





Threaded Adapters

Threaded adapters in a variety of configurations

- English, Metric, and NPT threaded adapters offered
- Bring together connectors with different threads
- Manufactured from inert polymers PEEK, PCTFE, ETFE, and PTFE


Two of the many challenges researchers face regularly, are trying to use one style of fitting for all connections, or trying to join two different sizes of tubing. To assist in overcoming these challenges we have engineered one of the most extensive threaded adapter lines available.

Threaded Adapters come in a wide variety of configurations to meet your system requirements. They are designed to effectively bring together connectors with different threads. We offer them in English, Metric, and NPT versions. Manufactured from inert polymers and stainless steel they deliver excellent chemical resistance.



English Threaded Adapters

Our versatile English Threaded Adapters are used specifically to securely attach connectors with different threads. We designed these adapters to work with English to English threaded geometries. Manufactured from Stainless Steel, PEEK, or Tefzel™ (ETFE), they deliver excellent solvent resistance.



P-627

0.020" (0.50 mm) thru-hole Includes (1) F-300 for 1/16" OD tubing

Threaded Adapters (Cont.)




NOTE

When using an adapter with male (external) threads, we recommend you first attach the adapter body into the receiving port, and then connect your tubing and fitting into the head of the adapter body.



- Use the 6000-076 Adapter to connect 1/16" OD tubing to the Preparative-Scale Injector Valve (page 66).
- You may not need an adapter to connect 1/16" OD tubing into your flat-bottom port. A less expensive alternative is to use a Flangeless Nut and Ferrule starting on page 45 or a Super Flangeless™ Nut and Ferrule starting on page 39.






Here are application ideas using two of our popular adapters:

- Many injection valves used in HPLC systems have 10-32 coned ports designed to accept 1/16" OD tubing. However, this may be a problem if large injection volumes are required (in excess of 10 mL). The most popular loops for large volume samples are made from 1/8" OD tubing, making it impossible to connect these larger volume loops to your injection valve. The solution: use our P-654 Adapter and the appropriate fittings for your sample loop. This set-up allows connection of 1/8" OD sample loop leads to your injection valve.
- Another potential application is connecting tubing to low-pressure solenoid valves with 1/4-28 flat-ottom ports. Most low-pressure valves of this type have very shallow threaded ports, which typically preclude the use of our Flangeless Fittings. However, by first threading our P-671 Adapter into the valve port(s), you can effectively use standard 1/4-28 fittings to connect your tubing into the backside of the adapter body. This also saves "wear and tear" on the threads in the valve ports.



Metric Threaded Adapters

Our versatile Metric Threaded Adapters are used specifically to effectively attach connectors with different threads. We designed these adapters to work with English to Metric threaded geometries. Manufactured from Stainless Steel, PEEK or Tefzel™ (ETFE) they deliver excellent solvent resistance.





- For an alternative to the Female M6 Adapters presented above, try a P-602 or P-622 Low Pressure Metric Union from page 84, along with the appropriate Metric Flangeless Fittings on page 43.
- To direct connect your tubing into a flat-bottom port, find the appropriate Flangeless or Super Flangeless™ Fittings on page 45 and page 39 respectively.
- Need metric fittings for your connections? See page 43.

Threaded Adapters (Cont.)



National Pipe Thread Adapters

These adapters make connections to female 1/8" and 1/4" National Pipe Thread (NPT) ports.

Manufactured from PEEK polymer, our NPT Adapters are durable and chemically resistant. We provide versions with either 1/4-28 or 5/16-24 flat-bottom threads, suitable for most low pressure applications.

Please Note: Wrap the threads on the NPT side of these adapters with thread seal tape (plumber's tape) to ensure a leak-free seal.



Replacement fittings for these adapters are located on the pages indicated below:

| | Page(s) |
|-----------------------------|---------|
| 1/4-28 for 1/8" OD tubing | 41 |
| 5/16-24 for 1/8" OD tubing | 41, 56 |
| 5/16-24 for 3/16" OD tubing | 41 |

Other tubing/fitting combinations are available. For more information, please contact your local Distributor or IDEX Health & Science directly.



U-5101/8" NPT to 1/4-28 Flat-Bottom Female Adapter for 1/8" OD tubing Includes (1) XP-308 Fitting





Our U-500 and U-510 NPT Adapters are great for attaching 1/8" OD fluoropolymer sparging lines to sparging gas tank regulating valves. Simply thread the appropriately-sized NPT Adapter into the valve's receiving port and then attach your sparging tubing to the adapter body using the fittings provided.

Threaded Adapters

| Part No. | Description | | | Includes | Thru-hole | Swept Volume | Pressure Rating | Qty. |
|----------|---|--------------|-----------|------------|------------------|--------------|---------------------|------|
| ENGLISH | THREADED ADAPTERS | | | | | | | |
| 6000-076 | PEEK Adapter, 5/16-24 C, M to 10-32 C, F | : | | N/A | 0.066" (1.70 mm) | 49.8 µL | 3,000 psi (207 bar) | ea. |
| P-135 | PEEK Adapter, 5/16-24 FB, F to 1/4-28 F | | | N/A | 0.080" (2.05 mm) | 4.1 µL | 1,000 psi (69 bar) | ea. |
| P-627 | PEEK Adapter, 10-32 C, F to 1/4-28 FB, F | | | (1) F-300 | 0.020" (0.50 mm) | 0.30 μL | 1,000 psi (69 bar) | ea. |
| P-681 | PCTFE Adapter, 5/16-24 FB, F to 1/4-28 F | В, М | | N/A | 0.125" (3.20 mm) | 96.6 µL | 1,000 psi (69 bar) | ea. |
| P-684 | PCTFE Adapter, 1/2-20 FB, F to 1/4-28 FE | 3, M | | N/A | 0.130" (3.30 mm) | 121.7 μL | 250 psi (17 bar) | ea. |
| P-718 | PCTFE Adapter, 5/16-24 FB, M to 1/4-28 | FB, F | | N/A | 0.040" (1.00 mm) | 10.3 μL | 1,000 psi (69 bar) | ea. |
| U-659 | PEEK Adapter, 5/16-24 FB, F to 1/2-20 FB | , F | | (1) XU-655 | Tapered** | 42.0 µL | 250 psi (17 bar) | ea. |
| U-665 | PEEK Adapter, 1/2-20 FB, F to 1/4-28 FB, | F | | (1) XU-655 | 0.063" (1.60 mm) | 6.6 µL | 250 psi (17 bar) | ea. |
| P-652 | PEEK Adapter, 1/4-28 FB, F to 10-32 C, N | | | N/A | 0.030" (0.75 mm) | 6.7 µL | 1,000 psi (69 bar) | ea. |
| P-654 | PEEK Adapter, 1/4-28 FB, F to 10-32 C, N | , Extra Long | | N/A | 0.030" (0.75 mm) | 9.5 µL | 1,000 psi (69 bar) | ea. |
| P-669-01 | PEEK Adapter, 10-32 C, F to 1/4-28 FB, M | | | N/A | 0.040" (1.00 mm) | 6.6 µL | 1,000 psi (69 bar) | ea. |
| P-671 | PTFE Adapter, 1/4-28 FB, F to 1/4-28 FB, | М | | N/A | 0.040" (1.00 mm) | 8.0 µL | 1,000 psi (69 bar) | ea. |
| P-672 | PEEK Adapter, 1/4-28 FB, F to 10-32 FB, I | Μ | | N/A | 0.050" (1.25 mm) | 11.4 µL | 1,000 psi (69 bar) | ea. |
| METRIC I | M6 THREADED ADAPTERS | | | | | | | |
| P-626 | PEEK Adapter, 10-32 C, F to M6 FB, F | | | (1) F-300 | 0.020" (0.50 mm) | 0.3 µL | 1,000 psi (69 bar) | ea. |
| P-650 | PEEK Adapter, M6 FB, F to 10-32 C, M St | andard | | N/A | 0.030" (0.75 mm) | 6.7 µL | 1,000 psi (69 bar) | ea. |
| P-670 | PCTFE Adapter, M6 FB, F to 1/4-28 FB, N | 1 | | N/A | 0.030" (0.75 mm) | 2.6 µL | 1,000 psi (69 bar) | ea. |
| P-673 | PCTFE Adapter, 5/16-24 FB, F to M6 FB, | М | | N/A | 0.040" (1.00 mm) | 9.9 µL | 1,000 psi (69 bar) | ea. |
| P-694 | PCTFE Adapter, 1/4-28 FB, F to M6 FB, N | 1 | | N/A | 0.040" (1.00 mm) | 11.3 µL | 1,000 psi (69 bar) | ea. |
| P-920-01 | PEEK Adapter, 10-32 C, F to M6 FB, M | | | N/A | 0.040" (1.00 mm) | 8.0 µL | 1,000 psi (69 bar) | ea. |
| 1/8" MAL | E NATIONAL PIPE THREAD ADAPTERS | | | | | | | |
| Part No. | Description | Color | Tubing OD | Includes | Thru-hole | Swept Volume | Pressure Rating | Qty. |
| U-510 | PEEK 1/8" NPT, M to 1/4-28 FB, F Adapter | Red | 1/8" | (1) XP-308 | 0.062" (1.60 mm) | 17.3 µL | 500 psi (34 bar) | ea. |
| U-514 | PEEK 1/8" NPT, M to 5/16-24 FB, F Adapter | Natural | 3/16" | (1) XP-132 | 0.125" (3.2 mm) | 70.4 µL | 500 psi (34 bar) | ea. |
| 1/4" MAL | E NATIONAL PIPE THREAD ADAPTERS | | | | | | | |
| U-500 | PEEK 1/4" NPT, M to 1/4-28 FB, F Adapter | Red | 1/8" | (1) XP-308 | 0.062" (1.60 mm) | 17.3 µL | 500 psi (34 bar) | ea. |
| U-504 | PEEK 1/4" NPT, M to 5/16-24 FB, F Adapter | Natural | 3/16" | (1) XP-132 | 0.125" (3.2 mm) | 70.4 µL | 500 psi (34 bar) | ea. |

- F = Female (internal) threads; M = Male (external) threads; XL = extra long; C = Coned; FB = Flat-Bottom
- * The pressure rating of this adaptive exceeds the pressure holding ability of the fittings and tubing used with it.

 ** Thru-hole tapers from 0.188" (4.80 mm) to 0.125" (3.20 mm).





High Pressure Mixing Tees

Mixing Tees utilize a specifically engineered internal geometry to efficiently mix two fluid streams into one combined stream. Mixing Tees are ideal for microbore or analytical gradient HPLC. These mixing tees are specifically designed for high pressure applications.

2 0.67" (1.70 cm)

U-466 and U-466SStatic Mixing Tees with F-300 Fingertight Fittings for 1/16" OD tubing

Static Mixing Tees

- > PEEK body with two-piece fingertight fittings
- Low swept volume

Static Mixing Tees are ideal for microbore or analytical gradient HPLC. They have a low swept volume of 2.2 μ L (includes frit volume) and are designed for flow rates of 0.5 to 3 mL/min and a maximum pressure of 5,000 psi (345 bar). The back pressure caused by the tee is typically only 10 to 20 psi (0.7 to 1.4 bar) at these flow rates. The thru-holes are 0.020" (0.50 mm) and the center port features a 10 μ m UHMWPE or stainless steel frit that aids mixing.



NOTE

- Turbulent mixing of solvents often increases outgassing. To maintain a bubble-free fluid pathway, we recommend solvent degassing when using this product.
- The frit incorporated into our U-466 and U-466S Static Mixing Tees is not replaceable. If it becomes clogged, the Mixing Tee must be replaced.

PEEK Frit Outlet

M-540
Micro Static Mixing Tee 0.010" thru-hole with fittings included (tubing and tubing sleeves not included)

Micro Static Mixing Tee

- Constructed of inert PEEK and PCTFE
- Low swept volume of 0.95 μL
- Designed for flow rates of 20–250 µL/min

Our Micro Static Mixing Tee utilizes a specifically engineered internal geometry to efficiently mix two fluid streams into one combined stream. The center port also features a 0.5 µm porosity PEEK polymer frit to aid in mixing. This frit adds a maximum of 20 psi (1.4 bar) back pressure to most systems (within the stated flow rate range). The Mixing Tee handles a maximum pressure of 5,000 psi (345 bar) when directly connecting 1/16" OD tubing, or up to 4,000 psi (276 bar) with capillary tubing when using our NanoTight[™] Fittings and Tubing Sleeves (page 37).



- ➤ See our Vacuum Degassing Systems on page 154.
- Our standard Static Mixing Tees are designed for flow rates from 0.5 mL/min to 3 mL/min.



High Pressure Mixing Tees (Cont.)



APPLICATION NOTE

Several researchers use our PEEK MicroTee to introduce ionizing voltage to their fluid stream just prior to a Mass Spectrometer¹. MicroTees are well suited for this application due to advantageous internal geometry and PEEK polymer's electrical resistance. The materials required for this setup are as follows: one gold or platinum conducting wire, one P-775 or P-875 MicroTee (this page), one MicroTight Tubing Sleeve (page 52) for the conducting wire (as needed to accommodate wire diameter), and at least two more MicroTight Tubing Sleeves (page 52) to connect your capillary tubing.

To set up a similar connection, first thread your wire through the appropriate tubing sleeve, if necessary, with the wire extending beyond both ends of the sleeve. Slip the female nut included with the MicroTee over the wire or sleeved wire, followed by the ferrule — ensuring the wire (and its sleeve) extends well past the end of the ferrule tip. Align the tip of the wire with the thru-hole of the MicroTee and gently insert the wire until it bottoms out. Now finger tighten the female nut into place. Attach your flow path tubing to the MicroTee's two other available ports, following the instructions provided with the MicroTee.

Begin fluid flow through the tee and apply voltage to the conducting wire lead. This setup typically provides effective electrospray ionization in applications having a flow rate of 100 µL/min or greater.

¹One such paper describing pioneering electrospray work: Protein Identification at the Low Femtomole Level from Silver-Stained Gels Using a New Fritless Electrospray Interface for Liquid Chromatography-Microspray and Nanospray Mass Spectrometry. Christine L. Gatlin, Gerd R. Kleemann, Lara G. Hays, Andrew J. Link, John R. Yates III (1998) Analytical Biochemistry 263, 93-101.

MicroTee & Cross for Capillary Tubing


- Direct connect 1/16", 1/32", 360 μm OD tubing, plus other capillary tubing
- Low swept volume

Use our MicroTees and MicroCrosses to join capillary tubing. All of these products are made entirely of PEEK and have 0.006" (0.150 mm) thru-holes, with resulting swept volumes ranging from 29 to 81 nL.





NOTE

▶ Use only the ferrules supplied with each connector — they are not interchangeable. Replacement ferrules and female nuts are available on page 35. For MicroUnions, MicroTees, and MicroCrosses for UHPLC applications, see page 72.









High Pressure Mixing Tees

| Part No. | Description | Threads | Includes | Thru-hole | Swept Volume | Pressure Rating | Qty. |
|----------|--|---------------|-------------------------|-------------------|--------------|---------------------|------|
| STATIC | MIXING TEE | | | | | | |
| U-466 | PEEK Static Mixing Tee for 1/16" OD Tubing, 10 μm UHMWPE Frit | 10-32 Coned | (3) F-300 | 0.020" (0.50 mm) | 2.2 μL | 5,000 psi (345 bar) | ea. |
| U-466S | PEEK Static Mixing Tee for 1/16" OD Tubing, 10 µm SST Frit | 10-32 Coned | (3) F-300 | 0.020" (0.50 mm) | 2.2 µL | 5,000 psi (345 bar) | ea. |
| MICRO | STATIC MIXING TEE | | | | | | |
| M-540 | PEEK Micro Static Mixing Tee, for 1/16" OD Tubing | 5/16-24 Coned | (3) F-132/P-416 | 0.010" (0.250 mm) | 0.95 µL | 5,000 psi (345 bar) | ea, |
| MICROT | TEE, MICROCROSS AND MICROELBOW | | | | | | |
| P-775 | PEEK MicroTee for MicroTight Sleeves | 5/16-24 Coned | (3) F-172, (3) P-416 | 0.006" (0.150 mm) | 29 nL | 4,000 psi (276 bar) | ea. |
| P-777 | PEEK MicroCross for MicroTight Sleeves | 5/16-24 Coned | (4) F-172, (4) P-416 | 0.006" (0.150 mm) | 38 nL | 4,000 psi (276 bar) | ea. |
| P-875 | PEEK MicroTee with Mounting Hole, for MicroTight Sleeves | 5/16-24 Coned | (3) F-172, (3) P-416 | 0.006" (0.150 mm) | 29 nL | 4,000 psi (276 bar) | ea. |
| P-885 | PEEK MicroTee for 1/32" OD Tubing | 5/16-24 Coned | (3) F-112, (3) P-416 | 0.006" (0.150 mm) | 29 nL | 5,000 psi (345 bar) | ea. |
| P-887 | PEEK MicroCross for 1/32" OD Tubing | 5/16-24 Coned | (4) F-112, (4) P-416 | 0.006" (0.150 mm) | 38 nL | 5,000 psi (345 bar) | ea. |
| P-888 | PEEK MicroTee for 360 μm OD Tubing | 5/16-24 Coned | (3) F-152, (3) P-416BLK | 0.006" (0.150 mm) | 29 nL | 5,000 psi (345 bar) | ea. |
| P-889 | PEEK MicroCross for 360 µm OD Tubing | 5/16-24 Coned | (4) F-152, (4) P-416BLK | 0.006" (0.150 mm) | 38 nL | 5,000 psi (345 bar) | ea. |
| P-890 | PEEK MicroTee for 1/16" OD Tubing | 5/16-24 Coned | (3) F-132, (3) P-416 | 0.006" (0.150 mm) | 58 nL | 5,000 psi (345 bar) | ea. |
| P-891 | PFFK MicroCross for 1/16" OD Tubina | 5/16-24 Coned | (4) F-132. (4) P-416 | 0.006" (0.150 mm) | 81 nL | 5,000 psi (345 bar) | ea. |

Importers & Manufacurers



High Pressure

PEEK ZDV Unions

Our PEEK zero-dead-volume (ZDV) Unions come complete with two F-300 Fingertight Fittings for 1/16" OD tubing and are pressure rated to 5,000 psi (344 bar).







Bio-Inert UHPLC Unions

- Unique, Patent-Pending Process allows a fully-PEEK fluid contact area combined with the strength of stainless steel
- Pressure rated to 17,400 psi (1,200 bar)
- Two inner diameters available: 0.008" and 0.016"


These unions are specifically engineered for Bio-Inert UHPLC applications. Combining the physical strength of 316 stainless steel with the inertness and biocompatibility of an all-PEEK fluid pathway, these unions will work well in applications where pressures reach up to 17,400 psi (1,200 bar) — without allowing metal contact by the fluid.

Neither union comes with fittings, but can be paired successfully with any 10-32 coned fitting that uses a polymer nose or ferrule.


Note: All-stainless steel fittings should NOT be used with these unions, as they will damage the internal conical seat.

NanoTight[™] Union

NanoTight Unions improve capillary tubing connections in several ways. The internal design of the union greatly reduces the incidence of tubing misalignment. When using 1/16" OD tubing sleeves (found on page 52) to connect capillary tubing, the webbed thru-hole minimizes breaking of fused silica while adding only miniscule swept volume. The results are fewer blockages, fewer flow rate reductions and fewer back pressure problems.



Bio-Inert UHPLC Union 0.008" thru-hole for 1/16" OD tubing



| Part No. | Description | Threads | Includes | Thru-hole | Swept Volume | Pressure Rating | Qty. | | | |
|-----------------|---|-------------|--------------|------------------|--------------|------------------------|------|--|--|--|
| BIO-INER | BIO-INERT UHPLC UNIONS | | | | | | | | | |
| UP-700 | Bio-Inert UHPLC Union for 1/16" OD Tubing, Natural (Tan) | 10-32 Coned | N/A | 0.008" (0.20 mm) | 0.05 μL | 17,400 psi (1,200 bar) | ea. | | | |
| PEEK ZD' | PEEK ZDV UNIONS | | | | | | | | | |
| P-704 | PEEK Union for 1/16" OD Tubing | 10-32 Coned | (2) F-300 | 0.020" (0.50 mm) | 0.28 μL | 5,000 psi (344 bar) | ea. | | | |
| P-742 | PEEK Union for 1/16" OD Tubing | 10-32 Coned | (2) F-300 | 0.010" (0.25 mm) | 0.07 μL | 5,000 psi (344 bar) | ea. | | | |
| P-760 | PEEK Union for 1/16" OD Tubing | 10-32 Coned | (2) F-300 | 0.050" (1.25 mm) | 1.2 µL | 5,000 psi (344 bar) | ea. | | | |
| NANOTIGHT UNION | | | | | | | | | | |
| P-779 | PFFK NanoTight Union for 1/16" OD Tubing and Tubing Sleeves | 10-32 Coned | (2) F-331N | 0.005" (125 µm) | 8 nL | 5,000 psi (344 bar) | ea. | | | |
| | HROMalytic +61(0)3 9762 2034 | | Distributors | | | | | | | |

Importers & Manufacurers

ECH no logy Pty Ltd


Website NEW: www.chromalytic.net.au E-mail: info@chromtech.net.au Tel: 03 9762 2034 . . . in AUSTRALIA



High Pressure PEEK Tees & Crosses

- > Highest pressure holding flat-bottom fitting system we offer
- Eliminates loosening of fittings due to tubing twist
- > Excellent for Tubing Assemblies
- > Holds tight even through vibration


Our PEEK Tees and Crosses include high pressure F-300 PEEK Fingertight Fittings allowing maximum operating pressures to 3,500 psi (241 bar) when used with 1/16" OD PEEK or stainless steel tubing.



0.020" thru-hole with F-300 Fittings

PEEK 7-Port Manifold

Combine several streams into one or split one fluid stream into several. This PEEK 7-Port Manifold comes complete with F-331 Fingertight Fittings for 1/16" OD tubing and offers a pressure rating of 5,000 psi (345 bar). Seal unused ports with any of our polymer 10-32 coned plugs on page 55.





| Part No. | Description | Threads | Includes | Thru-hole | Swept Volume | Pressure Rating | Qty. |
|----------|--|-------------|---|------------------|--------------|---------------------|------|
| PEEK TEI | ES AND CROSSES | | | | | | |
| P-727 | PEEK Tee for 1/16" OD Tubing | 10-32 Coned | (3) F-300 | 0.020" (0.50 mm) | 0.57 μL | 3,500 psi (241 bar) | ea. |
| P-728 | PEEK Tee for 1/16" OD Tubing | 10-32 Coned | (3) F-300 | 0.050" (1.25 mm) | 3.0 µL | 3,500 psi (241 bar) | ea. |
| P-729 | PEEK Cross for 1/16" OD Tubing | 10-32 Coned | (4) F-300 | 0.020" (0.50 mm) | 0.72 μL | 3,500 psi (241 bar) | ea. |
| PEEK MA | ANIFOLD | | | | | | |
| P-170 | PFFK 7-Port Manifold for 1/16" OD Tubing | 10-32 Coned | (7) F-331 | 0.020" (0.50 mm) | 2.2 μL | 5,000 psi (345 bar) | ea. |
| | HROMalytic +61(| | Australian Dis Importers & Man www.chromted | ufacurers | | | |
| | LUII | | | 70 (-6470) | | | |

Website NEW: www.chromalytic.net.au E-mail: info@chromtech.net.au Tel: 03 9762 2034 . . . in AUSTRALIA



High Pressure MicroTight® Unions

Conductive MicroTight Union

The Conductive MicroTight Union manufactured by IDEX Health & Science provides an excellent opportunity to introduce voltage into an electrospray or capillary electrophoresis system. With an extremely low internal volume of 16 nL, this union can be placed inline with 360 µm OD capillary tubing. Mount and apply voltage to these unions using our Insulating Mounting Bracket below.



Conductive MicroTight Union for 360 µm OD tubing with fittings and Capsule Union included

APPLICATION NOTE

For an example of using a Conductive MicroTight Union in a pressure driven ion preconcentration application see: "Self-Sealed Vertical Polymeric Nanoporous Junctions for High Throughput Nanofluidic Applications."


Sun Jae Kim and Jong Yoon Han. Analytical Chem. 2008, 80: 3507-3511.



Easily integrate the Conductive MicroTight Union into your system with our Insulating Mounting Bracket, shown on page 79.

MicroTight® Connectors for Capillary Tubing

Connect two pieces of capillary tubing with our PEEK MicroTight Connectors. The True ZDV Unions allow two pieces of tubing to connect directly to each other — using the included gauge plug to ensure proper alignment. The standard union and elbow both feature a 0.006" (0.150 mm) thru-hole, adding only a small amount of additional flow-path volume to help ensure proper chromatographic results.



| Part No. | Description | Threads | Includes | Thru-hole | Swept Volume | Pressure Rating | Qty. | | |
|--------------------|---|---------------|---------------------------------------|-------------------|--------------|---------------------|------|--|--|
| MICRO [*] | MICROTIGHT UNIONS | | | | | | | | |
| P-720 | PEEK True ZDV Union for MicroTight Sleeves | 6-32 Coned | (2) F-125, (1) P-553 | N/A | N/A | 4,000 psi (276 bar) | ea. | | |
| P-771 | PEEK True ZDV Union for 1/32" OD Tubing | 6-32 Coned | (2) F-126S, (1) P-553 | N/A | N/A | 5,000 psi (345 bar) | ea. | | |
| P-772 | PEEK Union for 360 µm OD Tubing | 5/16-24 Coned | (2) F-152, (2) P-416BLK | 0.006" (0.150 mm) | 5 nL | 5,000 psi (345 bar) | ea. | | |
| P-874 | PEEK MicroElbow for MicroTight Sleeves | 5/16-24 Coned | (2) F-172, (2) P-416 | 0.006" (0.150 mm) | 20 nL | 4,000 psi (276 bar) | ea. | | |
| REPLAC | REPLACEMENT GAUGE PLUGS (TO ACHIEVE TRUE ZDV CONNECTIONS WITH OUR P-720 AND P-771 UNIONS) | | | | | | | | |
| P-553 | Gauge Plug, Delrin® | 6-32 Coned | N/A | N/A | N/A | N/A | ea. | | |
| CONDL | CONDUCTIVE MICROTIGHT UNIONS | | | | | | | | |
| M-572 | Conductive Union for 360 um OD Tubina. PEFK/SST | 5/16-24 Coned | (2) F-152. (2) P-416BI K. (1) M-128NF | 0.011" (0.279 mm) | 16 nL | 5,000 psi (345 bar) | ea. | | |

HROM +61(0)3 9762 2034

ECH 1010 97 Pty Ltd


Website NEW: www.chromalytic.net.au E-mail: info@chromtech.net.au Tel: 03 9762 2034 . . . in AUSTRALIA

IICROTIGHT UNIONS



High Pressure Stainless Steel Tees & Crosses

These 316 stainless steel connectors come complete with 10-32 stainless steel fittings for use with 1/16" OD tubing and are rated to 20,000 psi (1,380 bar). They are compatible with any 10-32 coned threaded fittings.



U-428 Stainless Steel Tee 0.020" thru-hole with U-400 and U-401 Fittings



U-430 Stainless Steel Cross 0.020" thru-hole with U-400 and U-401 Fittings

| Part No. | Description | Threads | Includes | Thru-hole | Swept Volume | Pressure Rating | Qty. | | | |
|----------|---|---------------------|----------------------|------------------|--------------|------------------------|------|--|--|--|
| VHP TEE | VHP TEE FOR 1/16" OD TUBING | | | | | | | | | |
| U-428 | Stainless Steel Tee for 1/16" OD Tubing | 10-32 Coned | (3) U-400, (3) U-401 | 0.020" (0.50 mm) | 0.57 μL | 20,000 psi (1,380 bar) | ea. | | | |
| U-429 | Stainless Steel Tee for 1/16" OD Tubing | 10-32 Coned | (3) U-400, (3) U-401 | 0.040" (1.00 mm) | 2.1 µL | 20,000 psi (1,380 bar) | ea. | | | |
| U-430 | Stainless Steel Cross for 1/16" OD Tubing | 10-32 Coned | (4) U-400, (4) U-401 | 0.020" (0.50 mm) | 0.72 μL | 20,000 psi (1,380 bar) | ea. | | | |
| U-431 | Stainless Steel Cross for 1/16" OD tubina | 10-32 Coned | (4) U-400. (4) U-401 | 0.040" (1.00 mm) | 2.5 μL | 20,000 psi (1,380 bar) | ea. | | | |
| | HPOM = 1545 +61(0) | Australian Distribu | tors | | | | | | | |



VHP Stainless Steel ZDV Unions


Supplied with fittings for 1/16" OD or 1/8" OD tubing

- > Manufactured from 316 stainless steel
- All union assemblies rated to 20,000 psi (1,380 bar) or higher



It is possible to order the products on this page without the fittings. Simply use a -01 at the end of the product number to order the union body without fittings.

Our high pressure, zero-dead-volume (ZDV) unions are precision machined from 316 stainless steel, carefully passivated, then thoroughly rinsed. Each comes complete with stainless steel nuts and ferrules.





Website NEW: www.chromalytic.net.au E-mail: info@chromtech.net.au Tel: 03 9762 2034... in AUSTRALIA





VHP Tees & Crosses







VHP MicroCross

for 360 µm OD tubing with fittings included



Our VHP Stainless Steel Tees and Crosses are precision machined from durable stainless steel. It is mechanically designed for bringing together three or four pieces of tubing. Our VHP Tees & Crosses have an extremely high pressure rating of 30,000 psi (2070 bar).

VHP Tees & Crosses for Capillary Tubing

- Direct-connect either 360 μm or 1/32" OD tubing no sleeves required!
- Available in both tee and cross configurations
- Pressure rated to 15,000 psi (1,034 bar)

To help facilitate multi-port connections in UHPLC applications, our experts have developed a line of MicroTees and MicroCrosses, manufactured from stainless steel and featuring small thru-holes and very low internal volume. Additionally, the stainless steel construction allows these products to be used in applications where electrical conductivity is desired.

Included with the MicroTees and MicroCrosses are the VHP MicroFerrules found on page 59. The P-278 Extender Tool on page 33 can be used to tighten the female nuts that are included with these connectors.

APPLICATION NOTE

Why 1/32" OD Tubing and 360 µm OD Tubing?

IDEX Health & Science has focused strongly on the development of a variety of connectors and accessories for 1/32" OD tubing and 360 µm OD tubing. We have focused on these specific sizes due to their overwhelming popularity in analytical instruments, especially where micro and nano-scale analyses are being performed. By creating products designed for these popular sizes, the overall connection is easier to make and generally holds to increased pressures over connections where tubing sleeves are involved.

VHP Tee for 1/16" OD Tubing

IDEX Health & Science offers this Very High Pressure (VHP) Tee Connector, designed to bring three pieces of tubing together. The all-316 stainless steel connector is designed for 1/16" OD tubing and is pressure rated to 30,000 psi (2,070 bar).

| Part No. | Description | Threads | Includes | Thru-hole | Swept Volume | Pressure Rating | Qty. |
|----------|---|---------------|---|------------------|--------------|------------------------|------|
| VHP TEE | FOR 1/16" OD TUBING | | | | | | |
| UH-427 | VHP Tee for 1/16" OD Tubing, SST | 10-32 Coned | (3) VHP-200 | 0.020" (0.50 mm) | 0.57 μL | 30,000 psi (2,070 bar) | ea. |
| VHP TEES | S & CROSSES FOR CAPILLARY TUBING | | | | | | |
| UH-700 | VHP MicroTee for 1/32" OD Tubing, PEEK/SST | 5/16-24 Cone | d (3) PK-112, (3) P-416 | 0.010" (0.25 mm) | 84 nL | 15,000 psi (1,034 bar) | ea. |
| UH-750 | VHP MicroTee for 360 µm OD Tubing, PEEK/SST | 5/16-24 Coned | d (3) PK-152, (3) P-416BLK | 0.010" (0.25 mm) | 84 nL | 15,000 psi (1,034 bar) | ea. |
| UH-752 | VHP MicroCross for 360 um OD Tubina. PFFK/SST | 5/16-24 Coned | d (4) PK-152. (4) P-416BI K | 0.010" (0.25 mm) | 101 nL | 15,000 psi (1,034 bar) | ea. |
| | HROMalytic +61(0)39 | 762 2034 | Australian Distributors Importers & Manufacurers www.chromtech.net.au | > | I | | |

Website NEW: www.chromalytic.net.au E-mail: info@chromtech.net.au Tel: 03 9762 2034 . . . in AUSTRALIA



VHP MicroTight® Unions

- Featuring stainless steel bodies and PK/PEEK fittings
- > Pressure rated up to 15,000 psi (1,034 bar)
- Options to direct-connect both 1/32" OD tubing and 360 µm OD tubing

VHP Unions for Capillary Tubing

IDEX Health & Science has expanded its line of specialized fittings and connectors for UHPLC applications to include several innovative unions and adapters.


Two of these products — the UH-432 and UH-436 — follow the design of our popular Mini MicroFilters (see page 109) and allow a convenient union between either 1/32" OD tubing or $360~\mu m$ OD tubing. Each features a stainless steel union body and a unique stainless steel union capsule, enabling both excellent chemical compatibility as well as conductivity, making these a great choice for electrical interfacing in certain LC-MS applications. Each is also coupled with direct-connect ferrules made from our proprietary PEEK polymer blend (PK), allowing tubing connections up to 15,000 psi (1,034 bar). (*Please Note: While these connectors can be used at elevated pressures, they are not recommended for applications above 100 °C.*)

The UH-632 is a more traditionally designed connector, incorporating internally threaded ports. The union (UH-632) features a true ZDV (zero dead volume) connection between both tubes. This unique product is coupled with our one-piece Ultra-High Performance Fingertight fittings manufactured from our proprietary PEEK polymer blend, allowing them to be used in high temperature applications (up to 200 °C) at pressures up to 6,000 psi (414 bar) — or use these connectors at room temperature up to 15,000 psi (1,034 bar)!

The 1959-01 is a new VHP union designed to accept the popular M4x0.7 threaded fittings for 1/32" OD tubing. These unions will work nicely with both the VHP-900 fittings (found on page 59) as well as the reusable VHP-920 (found on page 60).



VHP MicroTight® Unions (Cont.)





What is a True ZDV Union?

True zero dead volume (ZDV) unions are designed so that the two joined pieces of tubing butt perfectly together as shown in the image to the right. These products have no swept volume contained within the union body. The fluid moves directly from one tube into another in this type of connector.

When using true ZDV unions, it is important to take care to ensure connecting tubing has burr-free 90 degree ends. Find tubing cutters on page 28 to assist with cleanly cutting polymer and fused silica tubing. Gauge plugs are supplied with True ZDV Unions to assist with assembly. With the gauge plug inserted into one side of the union, a hard stop is created for the tubing to bottom out against as it is connected to the opposite port. The gauge plug is removed and then the second piece of tubing is connected, using the first piece of tubing to bottom out against resulting in the two tubes joined together in the center of the union.



- > Find replacement VHP fittings on page 59.
- > Find Fused Silica tubing on page 16.
- Find 1/32" OD Stainless Steel tubing on page 19.
- To achieve 15,000 psi (1,034 bar) with the female threaded fittings used with some of these products, use the P-278 extender tool found on page 50.

VHP MicroTight Unions

| Part No. | Description | Threads | Includes | Thru-hole | Swept Volume | Pressure Rating | Qty. |
|----------|--|--|---|-------------------|-----------------|------------------------|------|
| VHP UN | IIONS FOR CAPILLARY TUBING | | | | | | |
| UH-432 | VHP Union for 1/32" OD Tubing, PEEK/SST | 5/16-24 Coned | (2) PK-112, (2) P-416 | 0.006" (0.150 mm) | 5 nL | 15,000 psi (1,034 bar) | ea. |
| UH-436 | VHP Union for 360 µm OD Tubing, PEEK/SST | 5/16-24 Coned | (2) PK-152, (2) P-416BLK | 0.006" (0.150 mm) | 5 nL | 15,000 psi (1,034 bar) | ea. |
| UH-632 | VHP True ZDV Union for 1/32" OD Tubing, PEEK/SST | 6-32 Coned | (2) PK-126, (1) P-553 Gauge Plug | N/A | N/A | 15,000 psi (1,034 bar) | ea. |
| 1959-01 | VHP Union for 1/32" OD Tubina. SST | M4x0.7 | N/A (Fittings must be ordered separately) | 0.007" (178 µm) | 16 nL | 30,000 psi (2,070 bar) | ea. |
| | HROMalytic +61(0 | · Control of the cont | Australian Distributors Importers & Manufacurers www.chromtech.net.au | | | | |

Website NEW: www.chromalytic.net.au E-mail: info@chromtech.net.au Tel: 03 9762 2034... in AUSTRALIA



MicroTight® Adapters

Convenient adapters for common 1/16" OD to capillary tubing

- Direct connect to 1/32" OD or 360 µm OD tubing options available
- VHP adapters pressure rated to 12,000 psi (828 bar)



While many 10-32 coned fittings are interchangeable, coned fittings using different threads are generally not interchangeable. As such, IDEX Health & Science recommends that only the style of coned fittings that accompanies these connectors be used for replacements.

Create a true zero dead volume (ZDV) connection between 1/16" OD tubing and capillary tubing with our MicroTight Adapters.

For Very High Pressure applications the UH-630 will connect 1/16" OD to 1/32" OD tubing in an inline true ZDV connection with the ability to withstand 12,000 psi (828 bar)! The materials of construction also allow this product to be used up to 200 °C, which reduces the pressure rating to 8,000 psi (552 bar). For more information on the fittings used with the VHP adapter, please see page 59.



VHP MicroTight Adapter for 1/16" and 1/32" OD tubing with fittings included



UH-906VHP MicroTight Adapting Cross
10-32 Coned for 1/16" OD tubing and 5/16-24 Coned for 360 μm OD tubing



1958-01 VHP MicroTight Adapter 10-32 Coned for 1/16" OD tubing and M4x0.7 for 1/32" OD tubing



MicroTight ZDV Adapter for 1/16" to 1/32" OD tubing





VHP MicroTight Adapter for 1/16" and 360 µm OD tubing with fittings included



UH-753VHP MicroTight Adapting Tee 360 μm (2 ports) to 10-32 C for 1/16" OD tubing (1 port)



UH-631-01 VHP MicroTight Adapter 10-32 Coned for 1/16" OD tubing and 6-40 Coned for 1/32" OD tubing fittings not included



MicroTight ZDV Adapter for 1/16" OD to MicroTight Sleeves with fittings included



P-882
MicroTight ZDV Adapter
for 1/16" to 360 µm OD tubing
with fittings included

MicroTight® Adapters (Cont.)



- > Replacement 6-32 fittings are on page 34.
- Replacement F-120 style nuts are on page 32 (when ordering, replace the "x" with an "R" or "B" to order either red or blue fittings).
- **)** Use this list to find micro flow products outside this chapter.

| | Page |
|---|------|
| $360~\mu m, 510~\mu m$ (0.020"), and 1/32" OD PEEK Tubing | 16 |
| 360 μm OD Fused Silica Tubing | 16 |
| 1/16" and 1/32" OD PEEKsil™ Tubing | 22 |
| 1/32" OD FEP Tubing | 26 |
| 360 μm OD High Purity PFA Tubing | 24 |
| 510 μm (0.020") and 1/32" OD Stainless Steel Tubing | 19 |
| Polymer Capillary and Fused Silica Tubing Cutters | 28 |
| MX Series II [™] Injection and Switching Valves | 119 |
| Manual Injection Valves | 121 |
| Micro Injection Port Adapters | 132 |
| Micro-Splitter Valves | 140 |
| Micro-Metering Valves | 141 |
| Microbore Guard Column | 158 |
| Ultra-Low Volume Back-Pressure Regulators | 144 |
| Nonmetallic 10-32 Micro-Volume Inline Check Valve | 137 |
| | |

MicroTight® Adapters

| | • | | | | | | |
|-----------|--|----------------------|--|-----------|--------------|-------------------------|------|
| Part No. | Description | Threads | Includes | Color | Swept Volume | Pressure Rating | Qty. |
| MICROTI | GHT ADAPTERS | | | | | | |
| P-770 | PEEK Micro Adapter, True ZDV, for 1/16" OD Tubing to MicroTight Tubing Sleeve | 10-32 C to 6-32 C | (1) F-120, (1) F-125, (1) P-554 | Natural | N/A | 4,000 psi (276 bar) | ea. |
| P-881 | PEEK Micro Adapter, True ZDV, for 1/16" to 1/32" OD Tubing | 10-32 C to 6-32 C | (1) F-120R, (1) F-126S, (1) P-554 | Red | N/A | 5,000 psi (345 bar) | ea. |
| P-882 | PEEK Micro Adapter, True ZDV, for 1/16" to 360 µm OD Tubing | 10-32 C to 6-32 C | (1) F-120B, (1) F-124S, (1) P-554 | Blue | N/A | 5,000 psi (345 bar) | ea. |
| UH-630 | Stainless Steel VHP Micro Adapter, for 1/16" to 1/32" OD Tubing | 10-32 C to 6-32 C | (1) PK-120BLK, (1) PK-126, (1) P-554 | SST/Black | N/A | 12,000 psi (827 bar) | ea. |
| UH-634 | Stainless Steel VHP Micro Adapter, for 1/16" to 360 µm OD Tubing | 10-32 C to 6-32 C | (1) PK-120BLK, (1) PK-124, (1) P-554 | SST/Black | N/A | 12,000 psi (827 bar) | ea. |
| UH-753 | Stainless Steel VHP Micro Adapting Tee, for 1/16" to 360 µm OD Tubing | 10-32 C to 5/16-24 C | (2) P-416BLK, (2) PK-152 | SST/Black | 152 nL | 15,000 psi (1,035 bar)* | ea. |
| 1958-01 | Stainless Steel VHP Micro Adapter, for 1/16" to 1/32" OD Tubing | 10-32 C to M4x0.7 C | N/A | SST | 16 nL | 30,000 psi (2,070 bar)* | ea. |
| UH-631-01 | Stainless Steel VHP Micro Adapter, for 1/16" to 1/32" OD Tubing | 10-32 C to 6-40 C | N/A | SST | 13 nL | 30,000 psi (2,070 bar)* | ea. |
| UH-906 | Stainless Steel VHP Micro Adapting Cross, for 1/16" to 360 µm OD Tubing | 10-32 C to 5/16-24 C | (2) PK-120BLK, (2) P-416BLK, (2) PK-152 | SST/Black | 0.11 μL | 15,000 psi (1,035 bar)* | ea. |
| REPLACE | MENT GAUGE PLUGS (TO ACHIEVE | TRUE ZDV CONNE | CTIONS WITH THE ABOVE AD | APTERS) | | | |
| P-554 | Delrin® Gauge Plug | 10-32 C | | White | N/A | N/A | ea. |
| C = Coned | satis and a sanda was an the Cation and | | | | | | |

* Pressure rating depends upon the fitting used.





cessories



Insulating Mounting Bracket, shown with lead wire and Conductive MicroTight Union, not included.

Insulating Mounting Bracket

Use our Insulating Mounting Bracket to easily integrate the Conductive MicroTight Union (shown on page 71) or our Conductive Mini MicroFilters (on page 109) into your system or lab.

The product snaps into place. Voltage from your lead wire is conducted through the attaching stainless steel nut and screw (included), then onto the mounted product via the stainless steel clip.

The bracket's base includes two holes (#2 screw clearance) for easy mounting onto any lab surface. Dimensions are 1.25" L x 0.45" W x 0.63" H.




> For lab-on-a-chip applications

- Doptions to connect 1/16" OD Tubing directly, or 360µm and 1/32" OD Tubing with tubing sleeves
- Wetted materials: PEEK and perfluoroelastomer

NanoPort Assemblies

NanoPort Assemblies provide consistent fluid connections for chip-based analyses. NanoPort connections will bond to a variety of substrate materials with the use of Loctite.™

All NanoPort components are made of inert, biocompatible PEEK polymer (nuts and ports), Perlast® perfluoroelastomer (gaskets), and ETFE (ferrules). Their unique design also prevents adhesive contamination of the fluid path. And NanoPort connections add no additional volume to the fluid path, virtually eliminating dead volume traditionally associated with chip-based fluid connections.





NOTE

Our NanoPort Assembly will readily connect 1/16" OD tubing with the included fittings. To connect 1/32" OD or 360µm OD, tubing sleeves for each size are included in each assembly.

Adhesive is not included in the N-333 NanoPort Assembly. Please contact IDEX Health & Science for bonding information or use common bonding adhesives such as Loctite.

| Part No. | Description | Threads | For Chip Hole | Tubing OD | Qty. |
|-------------|---------------------------|-------------------------------|---|-----------|-------|
| NANOPORT A | | | | | , |
| 10-32 Coned | NanoPort Assembly | | | | |
| N-333 | F-333N | F-142N | Up to 0.063" (1.6 mm) | 1/16" | ea. |
| NANOPORT I | REPLACEMENT PARTS | | | | |
| F-333Nx | Headless Fittings | 10-32 C | Up to 0.063" (1.6 mm) | 1/16″ | 10-pk |
| F-142Nx | Ferrules | 10-32 C | Up to 0.063" (1.6 mm) | 1/16" | 10-pk |
| Gaskets | | | | | |
| N-123-02 | Gasket. For all assemblie | s except 6-32 Coned Assemblie | s N/A | N/A | ea. |
| | HROMaly | +61(0)3 9762 2034 | Australian Distributors Importers & Manufacurers www.chromtech.net.au | | |

Website NEW: www.chromalytic.net.au E-mail: info@chromtech.net.au Tel: 03 9762 2034 . . . in AUSTRALIA




Low Pressure Manifolds


Choose a 5, 7, or 9 Port Manifold to combine several streams into one, or split one fluid stream into several. Each PEEK manifold comes complete with 1/4-28 Super Flangeless™ Fittings for either 1/16" or 1/8" OD tubing, with pressure ratings of 2,000 psi (138 bar) and 500 psi (34 bar), respectively.

A few useful applications include:

- Multiport mixing chamber
- Gas sparging splitting union
- Sample injection onto multi-well plates or a multiple direction flow path union







| Part No. | Description | Threads | Includes | Thru-hole | Swept Volume | Pressure Rating | Qty. |
|-----------|--|-----------|----------------------|------------------|--------------|--------------------|------|
| MANIFO | LDS | | | | | | |
| Standard | 4 | | | | | | |
| P-150 | PEEK 7-Port Manifold for 1/16" OD Tubing | 1/4-28 FB | (7) P-255, (7) P-250 | 0.040" (1.00 mm) | 42.0 µL | 1,000 psi (69 bar) | ea. |
| P-154 | PEEK 5-Port Manifold for 1/16" OD Tubing | 1/4-28 FB | (5) P-255, (5) P-250 | 0.040" (1.00 mm) | 22.3 µL | 1,000 psi (69 bar) | ea. |
| P-155 | PEEK 5-Port Manifold for 1/8" OD Tubing | 1/4-28 FB | (5) P-331, (5) P-359 | 0.062" (1.60 mm) | 53.8 μL | 500 psi (34 bar) | ea. |
| P-190 | PEEK 9-Port Manifold for 1/8" OD Tubing | 1/4-28 FB | (9) P-331, (9) P-359 | 0.062" (1.60 mm) | 160 μL | 500 psi (34 bar) | ea. |
| P-191 | PEEK 9-Port Manifold for 1/16" OD Tubing | 1/4-28 FB | (9) P-255, (9) P-250 | 0.040" (1.00 mm) | 139 µL | 1,000 psi (69 bar) | ea. |
| ED E1-+ F | 7 - 44 | | | | | | |



Low Pressure Bulkhead Unions

- Designed for plumbing tubing through equipment housing
- > For use with standard 10-32 coned or


1/4-28 flat-bottom threaded fittings




Thread PEEK Bulkhead Unions directly through your equipment housing to connect internal tubing to the outside. Each union has unique 3/8-24 external threads and comes complete with a stainless steel nut and lock washer to hold it in place. Requires a 3/8" hole to mount. The recommended torque limit for these unions is 15 in. – lbs (1.7 N·m).



Bulkhead Union includes stainless steel nut/lock washer





1/4-28 threads 0.67" (1.70 cm) 0.38" (0.97cm)

P-430 PEEK Elbow comes with Flangeless Fittings



Elbow Connectors

Use these Elbow Connectors to easily navigate tight corners. One Elbow is designed for use with 1/16" OD tubing and has a 0.020" (0.50 mm) thru-hole. Use 1/8" OD tubing with the other Elbow, which has a 0.062" (1.6 mm) thru-hole. Both come complete with 1/4-28 PEEK nuts and ETFE ferrules, and are pressure rated to 1,000 psi (69 bar).

Large Bore Union

> 5/16-24 flat-bottom threads

Use any of the 5/16-24 fittings on page 55 and the appropriate ferrule to create a true zero dead volume (ZDV) connection with the P-134 Union.



- Stainless Steel Bulkhead Unions are also available. Please contact us for more information.
- To use Elbows in higher pressure applications, simply replace the provided fittings with Super Flangeless™ Nuts and Ferrules, found on page 39.


| Part No. | Description | Threads | Color | Includes | Thru-hole | Swept Volume | Qty. |
|----------|--------------------------------|---------------------|---------|--------------------|------------------|--------------|------|
| BULKHEA | D UNIONS | | | | | | |
| P-440 | PEEK Bulkhead Union | 10-32 Coned | Natural | (1) SST Nut/Washer | 0.020" (0.50 mm) | 1.9 µL | ea. |
| P-441 | PEEK Bulkhead Union | 1/4-28 Flat-Bottom | Red | (1) SST Nut/Washer | 0.040" (1.00 mm) | 2.9 µL | ea. |
| P-441N | PEEK Bulkhead Union | 1/4-28 Flat-Bottom | Natural | (1) SST Nut/Washer | 0.040" (1.00 mm) | 2.9 µL | ea. |
| ELBOW C | ONNECTORS | | | | | | |
| P-430 | PEEK Elbow for 1/16" OD Tubing | 1/4-28 Flat-Bottom | Natural | (2) XP-235 | 0.020" (0.50 mm) | 1.4 µL | ea. |
| P-432 | PEEK Elbow for 1/8" OD Tubing | 1/4-28 Flat-Bottom | Natural | (2) XP-335 | 0.062" (1.60 mm) | 13.6 μL | ea. |
| LARGE BO | ORE UNION | | | | | | |
| P-134 | PFFK True 7DV Union | 5/16-24 Flat-Bottom | Natural | N/A | N/A | N/A | ea. |
| - | HPOM | | Austra | lian Distributors | | | |



Low Pressure Y Connectors

PEEK Y Connectors are designed to split a stream or join two streams together, just like a tee. However, the configuration of a tee can lead to turbulent flow and solvent outgassing, which increases baseline noise and reduces sensitivity. The geometry of a Y connector creates less turbulence and thus can improve analytical results.

All of these Y Connectors use 1/4-28 Flangeless fittings, except P-515 which uses 5/16-24 fittings (to accommodate larger tubing).



PEEK Y comes with Flangeless Fittings

| Part No. | Description | Threads | Includes | Thru-hole | Swept Volume | Pressure Rating | Qty. |
|----------|----------------------------|------------|------------|------------------|--------------|--------------------|------|
| Y CONN | NECTORS | | | | | | |
| P-512 | PEEK Y for 1/16" OD Tubing | 1/4-28 FB | (3) XP-235 | 0.020" (0.50 mm) | 1.7 µL | 1,000 psi (69 bar) | ea. |
| P-513 | PEEK Y for 1/8" OD Tubing | 1/4-28 FB | (3) XP-335 | 0.040" (1.00 mm) | 6.0 µL | 500 psi (34 bar) | ea. |
| P-514 | PEEK Y for 1/8" OD Tubing | 1/4-28 FB | (3) XP-335 | 0.060" (1.50 mm) | 13.6 μL | 500 psi (34 bar) | ea. |
| P-515 | PEEK Y for 3/16" OD Tubing | 5/16-24 FB | (3) XP-132 | 0.125" (3.20 mm) | 47.7 μL | 500 psi (34 bar) | ea. |




> Manufactured from PEEK, ETFE, Delrin®, polypropylene, or PCTFE

Available with 1/4-28, M6, or 10-32 flat-bottom threads

Low Pressure Unions


Our Low Pressure Unions are available in a variety of polymers, providing several lowcost and chemically-resistant options. The union assemblies below include fittings as shown in the table. The unions in the right column do not include fittings, allowing for customizing the fitting selection. In some cases, a union can be configured to connect two different tubing sizes—for example, if 1/4-28 Flangeless fittings for 1/16" and 1/8" OD tubing were selected from page 45 they can be used with the P-603 union to connect the two different tubing sizes.

Low Pressure PEEK Union Assemblies




PEEK Unions with Flangeless Fittings (1/4-28 internal threads)


Low Pressure ETFE Union Assemblies



Low Pressure Standard Unions



Low Pressure Metric Unions




P-602 & P-622 Metric Unions (M6 internal threads)





VacuTight[™] Union





- lacksquare To use connectors in higher pressure applications, simply replace the provided fittings with Super Flangeless™ Nuts and Ferrules, found on page 39.
- Use any of the 10-32 flat-bottom fittings on 39 and 42 to make an inline connection with our VacuTight Union. This product is designed for use with 1/16" OD tubing.

Low Pressure Unions

| Part No. | Description | Color | Threads | Includes | Thru-hole | Swept Volume | Pressure Rating | Qty. |
|------------|---|-----------------------|-------------------|------------------|------------------|--------------|--------------------|------|
| PEEK UN | IION ASSEMBLIES | | | | | | | |
| P-702 | PEEK Union for 1/16" OD Tubing | Natural | 1/4-28 FB | (2) XP-235 | 0.020" (0.50 mm) | 0.41 μL | 1,000 psi (69 bar) | ea. |
| P-703 | PEEK Union for 1/8" OD Tubing | Natural | 1/4-28 FB | (2) XP-335 | 0.050" (1.25 mm) | 2.57 µL | 1,000 psi (69 bar) | ea. |
| ETFE UN | ION ASSEMBLIES | | | | | | | |
| P-630 | ETFE True ZDV Union for 1/16" OD Tubing | Natural | 1/4-28 FB | (2) P-200N/P-245 | N/A | N/A | 1,000 psi (69 bar) | ea. |
| P-631 | ETFE True ZDV Union for 1/8" OD Tubing | Natural | 1/4-28 FB | (2) P-300N/P-345 | N/A | N/A | 1,000 psi (69 bar) | ea. |
| P-710 | ETFE Union for 1/16" OD Tubing | Natural | 1/4-28 FB | (2) XP-245 | 0.030" (0.75 mm) | 0.93 μL | 1,000 psi (69 bar) | ea. |
| STANDA | RD UNIONS | | | | | | | |
| P-603 | Delrin True ZDV Standard Union | Natural | 1/4-28 FB | N/A | N/A | N/A | N/A* | ea. |
| P-620 | Polypropylene True ZDV Standard Union | Natural | 1/4-28 FB | N/A | N/A | N/A | N/A* | ea. |
| P-623 | ETFE True ZDV Standard Union | Natural | 1/4-28 FB | N/A | N/A | N/A | N/A* | ea. |
| METRIC | UNIONS | | | | | | | |
| P-602 | Delrin Metric Union | Black | M6 FB | N/A | 0.020" (0.50 mm) | 0.41 μL | N/A* | ea. |
| P-622 | ETFE Metric Union | Blue | M6 FB | N/A | 0.020" (0.50 mm) | 0.41 µL | N/A* | ea. |
| MALE UI | NION | | | | | | | |
| P-645 | PCTFE Male Union | Natural | 1/4-28 FB | N/A | 0.062" (1.60 mm) | 61.3 µL | 500 psi (34 bar) | ea. |
| VACUTIO | SHT UNION | | | | | | | |
| P-845-01 | PEEK Union for 1/16" OD Tubing | Red | 10-32 FB | N/A | 0.020" (0.50 mm) | 0.20 μL | N/A* | ea. |
| * Pressure | Rating depends on Fittings selected. See pressu | ire rating for fittir | ngs on appropriat | e page. | | | | |

FB = Flat-Bottom

HROMalytic +61(0)3 9762 2034 ECH no logy Pty Ltd Website NEW: www.chromalytic.net.au E-mail: info@chromtech.net.au Tel: 03 9762 2034 . . . in AUSTRALIA

Australian Distributors Importers & Manufacurers www.chromtech.net.au



Low Pressure Tees & Crosses

0.49" (1.24 cm) Mounting Holes 0.13" diameter (0.33 cm) 1.08" (2.74 cm) Bio P-632 & P-633 0.37" (0.94 cm) ETFE Tees come with Flangeless Fittings 1/4-28 threads 1.08 0.13" diameter

Our Low Pressure Tees and Crosses are available in two inert polymers and can handle pressures to 500 psi (34 bar) or 1,000 psi (69 bar), depending upon the configuration of the products. Each is designed with handy mounting holes. All ETFE Tees and Crosses ship complete with 1/4-28 PFA Flangeless nuts and ETFE ferrules, while their PEEK polymer counterparts ship with 1/4-28 PEEK nuts and ETFE ferrules. Replacement fittings are located on page 47.



NOTE

comes with Flangeless Fittings

P-722 PEEK Cross

To order just the body of one of our tees and crosses without fittings, simply add a '-01' to the part number — e.g., P-632-01.

RELATED PRODUCTS

- > Seal off unused ports with any of our 1/4–28 flat-bottom plugs found on page 55.
- To use the PEEK polymer versions of our Tees and Crosses in higher pressure applications, simply replace the provided fittings with Super Flangeless™ Nuts and Ferrules, found on page 39.
- ▶ High Pressure Tees, Crosses, and a 7-Port Manifold (all with 10-32 threaded ports) are on page 74.

| Part No. | Description | Threads | Includes | Thru-hole | Swept Volume | Pressure Rating | Qty. |
|-----------|--------------------------------|--------------------|--|------------------|--------------|--------------------|------|
| LOW PRESS | SURE TEES AND CROSSES | | | | | | |
| P-632 | ETFE Tee for 1/16" OD Tubing | 1/4-28 Flat-Bottom | (3) P-245, (3) P-200N | 0.020" (0.50 mm) | 2.9 µL | 1,000 psi (69 bar) | ea. |
| P-633 | ETFE Tee for 1/8" OD Tubing | 1/4-28 Flat-Bottom | (3) P-345, (3) P-300N | 0.050" (1.25 mm) | 17.5 μL | 500 psi (34 bar) | ea. |
| P-634 | ETFE Cross for 1/16" OD Tubing | 1/4-28 Flat-Bottom | (4) P-245, (4) P-200N | 0.020" (0.50 mm) | 3.8 µL | 1,000 psi (69 bar) | ea. |
| P-635 | ETFE Cross for 1/8" OD Tubing | 1/4-28 Flat-Bottom | (4) P-345, (4) P-300N | 0.050" (1.25 mm) | 22.8 µL | 500 psi (34 bar) | ea. |
| P-712 | PEEK Tee for 1/16" OD Tubing | 1/4-28 Flat-Bottom | (3) XP-235 | 0.020" (0.50 mm) | 2.9 µL | 1,000 psi (69 bar) | ea. |
| P-713 | PEEK Tee for 1/8" OD Tubing | 1/4-28 Flat-Bottom | (3) XP-335 | 0.050" (1.25 mm) | 17.5 μL | 500 psi (34 bar) | ea. |
| P-714 | PEEK Tee for 1/16" OD Tubing | 1/4-28 Flat-Bottom | (3) XP-235 | 0.040" (1.00 mm) | 11.4 µL | 1,000 psi (69 bar) | ea. |
| P-722 | PEEK Cross for 1/16" OD Tubing | 1/4-28 Flat-Bottom | (4) XP-235 | 0.020" (0.50 mm) | 3.8 µL | 1,000 psi (69 bar) | ea. |
| P-723 | PFFK Cross for 1/8" OD Tubina | 1/4-28 Flat-Bottom | (4) XP-335 | 0.050" (1.25 mm) | 22.8 µL | 500 psi (34 bar) | ea. |
| | HROMalytic | +61(0)3 9762 2034 | Australian Distribu Importers & Manufacus www.chromtech.net. | rers | | | |

Website NEW: www.chromalytic.net.au E-mail: info@chromtech.net.au Tel: 03 9762 2034 . . . in AUSTRALIA



Quick Connect Luer Adapters

- Delrin®, polypropylene, ETFE, or PEEK Versions
- Adapts luers to 1/4-28, 10-32, 5/16-24, or M6 threaded ports

These luer adapters were designed to work in a variety of applications. By connecting any male luer to any female luer, you can create your own quick connect union or adapter. Each Quick Connect Luer Adapter conforms to ISO requirements for medical luer taper configuration and performance (45 psi/3.1 bar).

Find fittings to connect tubing to the threaded ports of these adapters in the Fittings chapter, starting on page 30.

Please Note: Our Female Quick Connect Luer Adapters can be used with any of the Male Luers on this page, i.e., those with and without lock hubs.




- Our A-626 Bottle Cap Plug (page 55) can be used to plug any of the female luer adapters on this page.
- > To prevent a chemical spill when disconnecting your solvent reservoir tubing from the pump, try our Quick-Stop Luer Check Valve on page 139.
- To economically prime an HPLC pump, simply remove the 10-32 fitting on the outlet check valve (standard on most models), insert a P-642 luer adapter, attach a syringe (such as our B-310) and draw the mobile phase through the pump head.




P-604, P-618, P-624 Female Luer to 1/4-28 Male (luer end of P-604 different than shown)



P-619, P-625 Male Luer to 1/4-28 Male



P-628 Female Luer to 1/4-28 Female



P-629 Female Luer to 10-32 Female



P-642 Female Luer to 10-32 Male




P-719 Female Luer to 10-32 Male





Male Luer to 10-32 Female



0.75" (1.91 cm)

8**8**

Quick Connect Luer Adapters (Cont.)


Luer-To-MicroTight® Adapter

> Easily connect 360 μm OD tubing to a syringe

The Luer-to-MicroTight Adapter is ideal for infusing sample into lab-on-a-chip devices. This product is made entirely of biocompatible PEEK polymer and introduces only 14 nL of additional volume to the flow path. Use it to directly connect a luer-tip syringe or other product that terminates with a standard male luer to $360 \, \mu m \, OD$ capillary tubing without tubing sleeves (see photo). MicroTight Fittings are included.



P-662 Luer-to-MicroTight Adapter, shown with a B-310 Syringe (see table below) and PEEK capillary tubing (page 16), not included.



P-662 Luer-To-MicroTight Adapter for Luer to 360 µm OD tubing with fittings included

Quick Connect Luar Adapters

| Part No. | Description | Body Material | Lock Hub Material | Thru-hole | Qty. |
|-------------|---------------------------------------|-------------------------------------|-------------------|------------------|------|
| QUICK CONN | IECT LUER ADAPTERS | | | | |
| P-604 | F Luer to 1/4-28 FB, M | Nat. Delrin | N/A | 0.05" (1.3 mm) | ea. |
| P-618 | F Luer to 1/4-28 FB, M | Nat. PP | N/A | 0.05" (1.3 mm) | ea. |
| P-619 | M Luer to 1/4-28 FB, M | Nat. PP | None * | 0.05" (1.3 mm) | ea. |
| P-624 | F Luer to 1/4-28 FB, M | Nat. ETFE | N/A | 0.05" (1.3 mm) | ea. |
| P-625 | M Luer to 1/4-28 FB, M | Nat. ETFE | None * | 0.04" (1.0 mm) | ea. |
| P-628 | F Luer to 1/4-28 FB, F | Nat. ETFE | N/A | 0.04" (1.0 mm) | ea. |
| P-629 | F Luer to 10-32 C, F | Nat. ETFE | N/A | 0.04" (1.0 mm) | ea. |
| P-642 | F Luer to 10-32 C, M | Nat. ETFE | N/A | 0.05" (1.3 mm) | ea. |
| P-655 | M Luer to 1/4-28 FB, F | Red PEEK | Black PEEK | 0.04" (1.3 mm) | ea. |
| P-656 | M Luer to 10-32 C, F | Nat. PEEK | Black PEEK | 0.05" (1.3 mm) | ea. |
| P-657 | M Luer to M6 FB, F | Black PEEK | Black PEEK | 0.05" (1.3 mm) | ea. |
| P-658 | F Luer to 1/4-28 FB, F | Red PEEK | N/A | 0.05" (1.3 mm) | ea. |
| P-659 | F Luer to 10-32 C, F | Nat. PEEK | N/A | 0.05" (1.3 mm) | ea. |
| P-660 | F Luer to M6 FB, F | Black PEEK | N/A | 0.05" (1.3 mm) | ea. |
| P-661 | F Luer to 5/16-24 FB, M | Nat. ETFE | N/A | 0.05" (1.3 mm) | ea. |
| P-675 | M Luer to 1/4-28 FB, F | Red ETFE | Natural PP | 0.05" (1.3 mm) | ea. |
| P-677 | M Luer to M6 FB, F | Black ETFE | Natural PP | 0.05" (1.3 mm) | ea. |
| P-678 | F Luer to 1/4-28 FB, F | Red ETFE | N/A | 0.05" (1.3 mm) | ea. |
| P-680 | F Luer to M6 FB, F | Black ETFE | N/A | 0.05" (1.3 mm) | ea. |
| P-683 | M Luer to 1/4-28 FB, M | Nat. PEEK | Black PEEK | 0.04" (1.0 mm) | ea. |
| P-686 | F Luer to M6 FB, M | Black ETFE | N/A | 0.05" (1.3 mm) | ea. |
| P-719 | F Luer to 10-32 C, M | Nat. PEEK | N/A | 0.05" (1.3 mm) | ea. |
| SYRINGE WIT | H MALE LUER LOCK | | | | |
| B-310 | 10 cc Disposable Luer-Lock Syringe. F | or use with any Female Luer Adapter | | 0.05" (1.3 mm) | ea. |
| LUER-TO-MIC | ROTIGHT ADAPTER | | | | |
| P-662 | Luer-to-MicroTight Adapter | (1) F-152, (1) P-416 | 0.006" (0.150 mm) | 45 psi (2.4 bar) | ea. |

F = Female (internal) threads; M = Male (external) t PP = Polypropylene; FB = Flat-Bottom; C = Coned * Slip-type male luer. HROMalytic +61(0)3 9762 2034 ECH no logy Pty Ltd

Australian Distributors Importers & Manufacurers www.chromtech.net.au Website NEW: www.chromalytic.net.au E-mail: info@chromtech.net.au Tel: 03 9762 2034 . . . in AUSTRALIA

89




LuerTight Fittings


- Luer fittings for fluoropolymer tubing
- Quick disconnect and barbless
- For 1/16" and 1/8" OD tubing

Our LuerTight fittings are specifically designed to connect fluoropolymer tubing without barbs or nuts! By integrating ferrules into the luer bodies, LuerTights will reliably hold your tubing in place while giving you the convenience of a luer connection. An inline set of LuerTight fittings provides a quick and easy disconnection option. LuerTight connections are also less bulky and more economical than nut-to-luer style fittings.

The bodies of these products are manufactured from polypropylene and the ferrules, where used, are made of ETFE.



P-837 LuerTight Fittings System for 1/16" OD tubing



LuerTight Fittings System for 1/8" OD tubing



LuerTight fittings are designed to be used exclusively within the LuerTight family. Combining LuerTight fittings with non-LuerTight luer products may result in a poor connection.

| Part No. | Description | Includes | Thru-hole | Pressure Rating | Qty. |
|---------------|--|---|-----------|-----------------|------|
| LUERTIGH | T FITTINGS SYSTEMS | | | | |
| P-837 | LuerTight System for 1/16" OD Tubing | (1) P-835, (1) P-836, (1) P-830T | N/A | 100 psi (7 bar) | ea. |
| P-838 | LuerTight System for 1/8" OD Tubing | (1) P-830, (1) P-831, (1) P-830T | N/A | 100 psi (7 bar) | ea. |
| LUERTIGH | T FITTING COMPONENTS | | | | |
| P-830 | Female Fitting for 1/8" OD Tubing | (1) Ferrule | N/A | N/A | ea. |
| P-830T | Set Plug to swage Ferrules into P-835 and P-830 | N/A | N/A | N/A | ea. |
| P-831 | Male Fitting for 1/8" OD Tubing | No Ferrule Required | N/A | N/A | ea. |
| P-835 | Female Fitting for 1/16" OD Tubing | (1) Ferrule | N/A | N/A | ea. |
| P-836 | Male Fitting for 1/16" OD Tubing | (1) Ferrule | N/A | N/A | ea. |
| Female = inte | ernal receiving luer pocket: Male = external luer nose (surr | ounded by internally-threaded locking ring) | | | |



Barbed Connectors

Our Type 1 Barbed Unions have been engineered to effectively join two pieces of soft-walled tubing together. This type of connector is typically the connector of choice for joining two peristaltic tubes with similar inner diameters together. Our unions are manufactured from either polypropylene or nylon.



Barbed to Barbed Adapters

- Adapters on this page feature various luer to barb adaptations
- Adapters on the next page feature a variety of barb-to-barb connectors

Use these barbed adapters to connect peristaltic-type flexible tubing for general, low pressure applications. The polypropylene used to manufacture the majority of these products is a Class VI material. Due to the low melt point of polypropylene (PP), these adapters are not autoclavable, however, they can be sterilized via gamma radiation. There are also Barb to Female Luer-Lock connectors available from ETFE, which has superior solvent resistance and a higher temperature rating (80 °C).



Barbed "Y" Adapters


Our Barbed "Y" Adapters, manufactured from polypropylene, are engineered to effectively join three pieces of soft-walled tubing together in a Y configuration, offering less turbulence and gentler mixing of fluids than a traditional Tee Connector. This type of connector works well for joining three peristaltic tubes with similar inner diameters together.



Thread to Barbed Adapters

- \blacktriangleright Three barb sizes, for 1/16", 1/8", and 3/16" ID flexible tubing
- Adapt to 1/4-28 flat-bottom, 5/16-24 flat-bottom, or 10-32 coned receiving ports

These adapters make it easy to connect flexible tubing to any standard 1/4-28 flat-bottom or 10-32 coned receiving port. Simply thread the adapter into a receiving port and slip tubing over the barbed stem to create a reliable low pressure connection.





- To connect low pressure fluoropolymer tubing, try the LuerTight™ Adapters on page 89.
- To connect peristaltic tubing to low pressure fluoropolymer tubing, see page 92.
- ▶ For peristaltic tubing, see page 15.



Swivel Barb Adapters


- ▶ Barb connection spins freely from the nut to prevent twist during installation
- Manufactured from polypropylene

The Swivel Barb Adapters from IDEX Health & Science are made up of two captive pieces acting as a one-piece fitting for ease of use. Manufactured from polypropylene and available in three barb sizes, the Swivel Barb will facilitate connection between flexible tubing to a 1/4-28 flat-bottom port. The barbed insert spins freely from the threaded nut in order to prevent the tubing from twisting during installation.



Luer to Barbed Adapters

Our Luer to Barbed Adapters are an excellent choice when connecting between soft-walled tubing and luer-based products, such as a syringe or a low-pressure filter, for example. We offer several different configurations, allowing you to connect various sizes of soft-walled tubing to receiving ports that accept a male luer "slip" style connection; a male luer "lock" style connection; and a female-luer style of connector.



Barbed Connectors (Cont.)



Conical Adapters


- Direct connect 1/16" and 1/8" OD rigid and semi-rigid tubing to peristaltic tubing
- ▶ Accept 0.020"-1/8" (0.50-3.2 mm) ID peristaltic tubing
- Biocompatible flow path with excellent chemical compatibility

Conical Adapters provide a reliable connection between rigid/semi-rigid tubing and peristaltic-type flexible tubing, such as Tygon® and PharMed®. These adapters are composed of a PEEK polymer female nut, our Super Flangeless™ ferrule system and an ETFE or PEEK conical adapter body. The narrow coned end of the adapter body allows peristaltic tubing to slide on more easily than it does onto conventional barbed adapters. Peristaltic tubing is also easier to remove from our Conical Adapters, since no cutting or excessive pulling is required.



APPLICATION NOTE

To help secure peristaltic tubing more firmly to the Conical Adapters, simply attach a cable tie to the outside of the peristaltic tubing once it has been placed onto the Adapter body.








Peristaltic Tubing Adapters

These unique adapters connect peristaltic tubing to standard 1/16" or 1/8" OD tubing. A specially-designed nose allows the peristaltic tubing to simply press fit over the nose and then be held tightly in place by the retainer sleeve. Your 1/16" OD tubing may then be connected with the Flangeless Fittings supplied with the adapter. To connect your peristaltic tubing to tubing with a different OD, simply replace the supplied fittings with your choice of Flangeless Fittings from page 45.

One popular application for these adapters is to use them as "stops" for your peristaltic pump. By doing so, you can reduce the amount of peristaltic tubing required for your flow path, thus reducing cost.



Barbed Connectors

| P-801 P-802 BARBED P-860 P-861 P-862 P-863 | TO BARBED ADAPTERS 0.06" (1.5 mm) 0.12" (3.0 mm) "Y" CONNECTORS 0.06" (1.5 mm) 0.10" (2.5 mm) | | | Polypropylene Polypropylene | | | ea. |
|--|---|------------------------|--------------------|--------------------------------|------------------------|---------------------------------|----------|
| P-802 BARBED P-860 P-861 P-862 P-863 | 0.12" (3.0 mm) "Y" CONNECTORS 0.06" (1.5 mm) | | | | | | |
| BARBED P-860 P-861 P-862 P-863 | "Y" CONNECTORS 0.06" (1.5 mm) | | | | | | |
| P-860 P-861 P-862 P-863 | 0.06" (1.5 mm) | | | | | | ea. |
| P-860 P-861 P-862 P-863 | 0.06" (1.5 mm) | | | | | | |
| P-862 P-863 | | | | Polypropylene | | | ea. |
| P-862 P-863 | | | | Polypropylene | | | ea. |
| P-863 | 0.12" (3.0 mm) | | | Polypropylene | | | ea. |
| | 0.18" (4.8 mm) | | | Polypropylene | | | ea. |
| P-864 | 0.25" (6.4 mm) | | | Polypropylene | | | ea. |
| | TO BARBED ADAPTERS | | | Готургоругене | | | ea. |
| Part No. | | | | Material | Thusada | Thur, hala | 04. |
| | Description | | | ETFE | Threads 10-32 Coned | Thru-hole 0.04" (1.0 mm) | Qty. |
| P-663 | Barb Adapter, 1/16" (1.55 mm) ID Tubing | | | | | | ea. |
| P-646 | Barb Adapter, 1/16" (1.55 mm) ID Tubing | | | ETFE | 1/4-28 Flat-Bottom | 0.04" (1.0 mm) | ea. |
| P-647 | Barb Adapter, 1/8" (3.20 mm) ID Tubing | | | ETFE | 1/4-28 Flat-Bottom | 0.08" (2.0 mm) | ea. |
| P-648 | Barb Adapter, 3/16" (4.75 mm) ID Tubing | | | ETFE | 1/4-28 Flat-Bottom | 0.10" (2.5 mm) | ea. |
| P-668 | Barb Adapter, 1/16" (1.55 mm) ID Tubing | | | PEEK | 1/4-28 Flat-Bottom | 0.04" (1.0 mm) | ea. |
| P-689 | Barb Adapter, 3/16" (4.75 mm) ID Tubing | | | ETFE | 5/16-24 Flat-Bottom | 0.10" (2.5 mm) | ea. |
| P-692 | Barb Adapter, 0.020" to 1/32" (0.50 to 0.8 | 0 mm) ID Tubing | | PEEK | 1/4-28 Flat-Bottom | 0.02" (0.5 mm) | ea. |
| SWIVEL | BARB ADAPTERS | | | | | | |
| D-646 | Swivel Barb Adapter, 1/16" (1.55 mm) ID | Tubing | | Polypropylene | 1/4-28 Flat-Bottom | 0.03" (0.75 mm) | ea. |
| D-647 | Swivel Barb Adapter, 3/32" (2.40 mm) ID | Tubing | | Polypropylene | 1/4-28 Flat-Bottom | 0.056" (1.5 mm) | ea. |
| D-648 | Swivel Barb Adapter, 1/8" (3.20 mm) ID To | ubing | | Polypropylene | 1/4-28 Flat-Bottom | 0.08" (2.0 mm) | ea. |
| BARB TC | SLIP-TYPE MALE LUER ADAPTERS | | | | | | |
| Part No. | Description | | | Material | | | Qty. |
| P-854 | Male Luers (Slip-type) for use with 1/16" I | D (1.55 mm) Tubing | | PP | | | ea. |
| | A=0.046" B=0.064" C=0.090" D=0.129" I | | | | | | Cu. |
| | type male luer fittings are for use in systems | | are not desired. | | | | |
| BARB TO | MALE LUER WITH LOCK RING ADA | PTERS | | | | | |
| P-850 | Male Luers with Lock Ring for use with 1/ A=0.049" B=0.065" C=0.090" E=0.583" F | | | PP | | | ea. |
| P-851 | Male Luers with Lock Ring for use with 3/ | | | PP | | | |
| | A=0.071" B=0.100" C=0.139" E=0.681" F Male Luers with Lock Ring for use with 1/ | | | | | | ea. |
| P-852 | A=0.099" B=0.132" C=0.184" E=0.777" F | =0.436" | | PP | | | ea. |
| BARB TC | FEMALE LUER-LOCK CONNECTOR | | | | | | |
| P-857 | Female Luer Connectors for use with 1/10 A=0.030" B=0.063" C=0.106" D=0.100" I | | | PP | | | ea. |
| P-858 | Female Luer Connectors for use with 3/32 A=0.056" B=0.102" C=0.145" D=0.100" I | | | PP | | | ea. |
| P-859 | Female Luer Connectors for use with 1/8 | ' ID (3.20 mm) Tubing | | PP | | | ea. |
| | A=0.080" B=0.135" C=0.187" D=0.100" I For use with 1/16" (1.55 mm) ID Tubing | E=0.733" F=0.253" | | | | | cu. |
| P-870 | A=0.030" B=0.063" C=0.106" D=0.100" I | E=0.598" F=0.253" | | ETFE | | | ea. |
| P-872 | For use with 1/8" (3.20 mm) ID Tubing A=0.080" B=0.137" C=0.187" D=0.100" I | E=0.733" F=0.253" | | ETFE | | | ea. |
| PERISTAI | TIC TUBING ADAPTERS | | | | | | |
| Part No. | Description | Tubing OD | | Peristaltic Tubing ID |) | Thru-Hole | Qty. |
| P-757 | Standard Adapter | up to 0.180" (4.5 | 5 mm) | 0.048"- 0.110" (1.20 | 1–2.80 mm) | 0.030" (0.75 mm) | ea. |
| P-767 | Large Bore Adapter | up to 0.250" (6.3 | 5 mm) | 0.100"- 0.150" (2.55 | i–3.80 mm) | 0.070" (1.78 mm) | ea. |
| CONICA | ADAPTER ASSEMBLIES | | | | | | |
| Part No. | Description | Rigid or Semi-Rig | aid Tubina OD | Peristaltic Tubing ID |) | Thru-Hole | |
| P-794 | Conical Adapter | 1/16" | J | 0.020"-0.030" (0.50 | | 0.020" (0.50 mm) | ea. |
| P-797 | Conical Adapter | 1/16" | | 1/16"-3/32" (1.55 m | | 0.040" (1.0 mm) | ea. |
| P-798 | Conical Adapter | 1/8" | | 1/16"-3/32" (1.55 m | | 0.040" (1.0 mm) | ea. |
| P-799 | Conical Adapter | 1/8" | | 3/32"-1/8" (2.40 mr | | 0.060" (1.5 mm) | ea. |
| | ADAPTER REPLACEMENT PARTS | 170 | | 3/32 1/0 (2.401111 | 11 3.20 mmy | 0.000 (1.5 11111) | Cu. |
| | | Earling With | | Material | | | |
| Part No. | Description | For Use With | | Material | | | |
| F-156 | Female Nut, 1/8", 1/4-28 | P-798, P-799 | | Black PEEK | | | ea. |
| P-420 | Female Nut, 1/16", 1/4-28 | P-794, P-797 | | Natural PEEK | | | ea. |
| P-259 | Super Flangeless Ferrule, 1/16" | P-794, P-797 | | Yellow ETFE/SST | | | ea. |
| P-359 | Super Flangeless Ferrule, 1/8" | P-798, P-799 | | Yellow ETFE/SST | | | ea. |
| P-691 | Conical Adapter Body | P-799 | | Natural ETFE | | | ea. |
| P-692 | Conical Adapter Body | P-794 | Augustian Dist | Natural PFFK | | | ea. |
| | HROMalytic | +61(0)3 9762 2034 | Australian Distr | | | | |
| | ECH | Ology Pty Ltd | www.chromtech. | | | 00 (| . 470\ 0 |
| Website I | IEW : www.chromalytic.net.au E-mail : in | fo@chromtech.net.au Te | el: 03 9762 2034 i | in AUSTRALIA | | 93 (of | f 172) 2 |





FILTERS & FRITS

Our Filters offer an optimal way to filter your solvents, preventing pump cavitation and system damage. We offer different style filters for specific system specifications. Our filters protect your system from particulate matter from the solvent that may otherwise damage expensive hardware.

We offer a complete line of Frits manufactured from two different materials: PEEK and stainless steel. Both materials offer a variety of sizes of frit discs, as well as being available in numerous porosities. All our frits are designed with exceptional uniform porosity and a long filtration life.

95 FRITS

100 FILTERS

111 BOTTLE CAPS & PLUGS



Stainless Steel Frits

Our Analytical-scale 316 Stainless Steel Frits are available in 0.5 µm or 2 µm porosity—the most common HPLC filtration ratings. Each frit includes a PCTFE or PEEK polymer sealing ring.


Many of the frits shown have the common 0.250" (0.64 cm) and 0.254" (0.64 cm) ODs, which allow them to be used in many of the Precolumn and Inline Filters found starting on page 103. Choose the larger diameter faces and/or larger porosity frits for faster flow rates. Choose frits with a smaller diameter face and/or smaller porosity for applications sensitive to extra flow path volume.




To Clean Or Not To Clean?

It is rarely worth the time and effort to clean frits, given the relatively low cost of replacements. Furthermore, cleaning may leave some debris embedded in the frit pores. If the washed frit is accidently returned to your instrument in a reverse orientation, any remaining debris could be flushed out and deposited further down the fluid path. If this frit is being used as a column head frit, the debris may be washed directly onto the column bed.

0.5 µm Stainless Steel Frits



2 µm Stainless Steel Frits





Frits without the polymer rings cannot be used with our standard Precolumn and Inline Filter assemblies.

Semi-Prep Stainless Steel Frits

Many of these frits come complete with a PCTFE, ETFE, or PTFE sealing ring. Choose from 2 μ m, 5 μ m, 10 μ m, and 20 μ m filtration porosities and a range of diameters to match your intended flow rate and filtration requirements.

2 µm Semi-Prep Stainless Steel Frits





Stainless Steel Frits (Cont.)

5 µm Semi-Prep Stainless Steel Frits



C-417

10 µm Semi-Prep Stainless Steel Frits









20 µm Semi-Prep Stainless Steel Frits









Stainless Steel Frits

| Part No. | Porosity | Disc Diameter | Disc Thickness | Ring OD | Ring Material | Frit Volume | Qty. |
|--------------|--------------------|------------------|------------------|------------------|---------------|-------------|------|
| STAINLESS ST | EEL FRITS | | | | | | |
| A-100 | 2 µm | 0.094" (0.24 cm) | 0.062" (0.16 cm) | 0.250" (0.64 cm) | PEEK | 1.7 µL | ea. |
| A-101 | 2 µm | 0.062" (0.16 cm) | 0.062" (0.16 cm) | 0.250" (0.64 cm) | PEEK | 0.7 μL | ea. |
| A-102 | 0.5 µm | 0.062" (0.16 cm) | 0.062" (0.16 cm) | 0.250" (0.64 cm) | PEEK | 0.6 μL | ea. |
| A-103 | 0.5 µm | 0.094" (0.24 cm) | 0.062" (0.16 cm) | 0.250" (0.64 cm) | PEEK | 1.4 µL | ea. |
| A-420 | 2 µm | 0.125" (0.32 cm) | 0.062" (0.16 cm) | 0.200" (0.51 cm) | PCTFE | 3.0 µL | ea. |
| C-128-31 | 0.5 µm | 0.038" (0.10 cm) | 0.028" (0.07 cm) | 0.125" (0.32 cm) | PEEK | 0.1 μL | ea. |
| C-140-30 | 0.5 µm | 0.188" (0.48 cm) | 0.062" (0.16 cm) | 0.254" (0.65 cm) | PCTFE | 6.5 µL | ea. |
| C-401 | 2 µm | 0.125" (0.32 cm) | 0.062" (0.16 cm) | 0.250" (0.64 cm) | PEEK | 3.0 µL | ea. |
| C-402 | 2 μm | 0.188" (0.48 cm) | 0.062" (0.16 cm) | 0.254" (0.65 cm) | PEEK | 7.8 µL | ea. |
| C-407 | 2 µm | 0.038" (0.10 cm) | 0.028" (0.07 cm) | 0.062" (0.16 cm) | PCTFE | 0.1 μL | ea. |
| C-408 | 2 µm | 0.038" (0.10 cm) | 0.028" (0.07 cm) | 0.125" (0.32 cm) | PEEK | 0.1 µL | ea. |
| C-409 | 0.5 µm | 0.038" (0.10 cm) | 0.028" (0.07 cm) | 0.062" (0.16 cm) | PCTFE | 0.1 μL | ea. |
| C-420 | 2 μm | 0.038" (0.10 cm) | 0.028" (0.07 cm) | 0.192" (0.49 cm) | PCTFE | 0.1 µL | ea. |
| C-425 | 0.5 μm | 0.038" (0.10 cm) | 0.028" (0.07 cm) | 0.192" (0.49 cm) | PCTFE | 0.1 μL | ea. |
| SEMI-PREP ST | AINLESS STEEL FRIT | rs | | | | | |
| A-105 | 10 µm | 0.062" (0.16 cm) | 0.062" (0.16 cm) | 0.250" (0.64 cm) | PCTFE | 0.9 μL | ea. |
| A-106 | 10 µm | 0.094" (0.24 cm) | 0.062" (0.16 cm) | 0.250" (0.64 cm) | PCTFE | 2.0 µL | ea. |
| A-107 | 10 μm | 0.189" (0.48 cm) | 0.062" (0.16 cm) | 0.254" (0.65 cm) | PCTFE | 9.1 µL | ea. |
| A-120 | 20 μm | 0.125" (0.32 cm) | 0.062" (0.16 cm) | 0.250" (0.64 cm) | PCTFE | 3.7 µL | ea. |
| A-122 | 20 μm | 0.188" (0.48 cm) | 0.062" (0.16 cm) | 0.254" (0.65 cm) | PCTFE | 9.7 µL | ea. |
| A-224 | 20 μm | 0.188" (0.48 cm) | 0.062" (0.16 cm) | 0.254" (0.65 cm) | PTFE | 9.7 µL | ea. |
| A-331 | 10 μm | 0.750" (1.91 cm) | 0.062" (0.16 cm) | 0.880" (2.24 cm) | ETFE | 141.9 µL | ea. |
| A-332 | 2 μm | 0.750" (1.91 cm) | 0.062" (0.16 cm) | 0.880" (2.24 cm) | ETFE | 141.9 µL | ea. |
| A-337 | 20 μm | 0.750" (1.91 cm) | 0.062" (0.16 cm) | 0.880" (2.24 cm) | ETFE | 152 μL | ea. |
| A-343 | 2 um | 0.625" (1.59 cm) | 0.062" (0.16 cm) | 0.750" (1.91 cm) | PCTFE | 112.6 µL | ea. |

Australian Distributors Importers & Manufacurers www.chromtech.net.au



PEEK Frits

- Inert, biocompatible, and metal-free
- > Uniform porosity, longer filtration life
- > Sealing rings manufactured from PCTFE

Patented IDEX Health & Science PEEK Frits offer exceptionally uniform porosity. This property ensures longer filtration life and consistent frit-to-frit swept volumes. The PEEK polymer frit discs are biocompatible and inert to most solvents, making them well-suited for bioanalytical applications. PEEK's robust properties make these products suitable for low and high pressure applications.

Disc rings, included on all PEEK frits, are made of PCTFE and are slightly thicker than the frit disc, providing enhanced sealing and excellent chemical resistance. PCTFE surrounded PEEK frits can be used up to 80 °C.

0.5 µm PEEK Frits



_0.062" (0.16 cm) diameter 0.062" (0.16 cm) thick



0.195" (0.5 cm) diameter 0.062" (0.64 cm) (0.16 cm) thick

A-701

A-703

A-707

2 µm PEEK Frits



_0.062" (0.16 cm) diameter 0.062" (0.16 cm) thick



0.091" (0.23 cm) diameter 0.062" (0.16 cm) thick



0.125" (0.32 cm) diameter 0.062" (0.16 cm) thick

A-700



.0.188" (0.48 cm) diameter 0.062" (0.16 cm) thick 0

A-702

___ 0.062" (0.16 cm) diameter 0.062" (0.16 cm) thick 0.200" 0.1 (0.51 cm) (0

0.125" (0.32 cm) diameter 0.062" (0.16 cm) thick

A-706


A-708

0.200" / (0.51 cm)


A-710

A-704

2 µm Semi-Prep PEEK Frits



$5\,\mu m$ and $10\,\mu m$ PEEK Frits



HROMalytic +61(0)3 9762 2034

ECH 100 9 Y Pty Ltd

Website NEW: www.chromalytic.net.au E-mail: info@chromtech.net.au Tel: 03 9762 2034 . . . in AUSTRALIA

97 (of 172) 2017Sept

PEEK Frits (Cont.)



NOTE

- > The thickness dimension in the part drawings and the pricing tables represents the thickness of the frit disc not the frit ring. Frit rings are often slightly thicker to ensure a proper seal. When tightened into a filter holder the ring compresses to nearly match the thickness of the frit disc.
- The manufacturing process may cause some slight color variance in our PEEK frits. This does not affect their quality or performance. Frit dimensions are approximate. Actual batch-to-batch frit dimensions may vary slightly.




RELATED PRODUCTS

Any 0.247" to 0.254" diameter frit (including polymer ring) can be used with the Standard HPLC Inline Solvent Filters on page 102 and the Standard Precolumn Filters on page 105.



The term "frit volume" refers to the volume of the various fluid pathways that comprise the matrix of a frit. A standard frit is a mass of small particles fused together through a controlled process of compression and heat. Because of their shape, there are gaps between the fused particles. Fluid makes its way through these gaps, creating a pathway from one side of the frit to the other (see the diagram, below, where the white circles represent frit particles, and the black area represents the void between the particles.)

Generally, when the frit particles increase in size, the frit's porosity increases as well. The larger the particles, the larger the gaps between particles. Cumulatively, these gaps comprise what is known as "frit volume." Using gravimetric determination, it has been experimentally shown that the total volume of any given frit may range from 18%-30%, depending upon the porosity of the frit.



Frit volume is calculated by determining what the mass of the frit would be if it were a solid block of material of equal size. Then the solid mass of the frit is multiplied by the percentage assigned to the porosity to determine the theoretical frit volume.

20% for 0.5 µm frits 24% for 2 μm frits

26% for 5 µm frits

30% for 20 µm frits

28% for 10 µm frits

From a chromatographic perspective, it's important to know the volume of the frit used in your system. It is possible for a frit to negatively impact your chromatography if the total frit volume is too large and if it is placed in an area through which the sample will pass. To avoid frit-related problems like band broadening and loss of resolution, most inline filters placed after the sample introduction point (e.g., between the injection valve and the column) are smaller in size and porosity than inline filters that are placed in areas before the sample is introduced into the flow path (e.g., between the pump and the injection valve).

PEEK Frits

| Part No. | Porosity | Disc Diameter | Disc Thickness | Ring OD | Ring Material | Frit Volume | Qty. |
|--------------|----------|------------------|------------------|------------------|---------------|-------------|------|
| PEEK FRITS | | | | | | | |
| A-700 | 2 µm | 0.062" (0.16 cm) | 0.062" (0.16 cm) | 0.250" (0.64 cm) | PCTFE | 0.7 μL | ea. |
| A-701 | 0.5 µm | 0.062" (0.16 cm) | 0.062" (0.16 cm) | 0.250" (0.64 cm) | PCTFE | 0.6 μL | ea. |
| A-702 | 2 µm | 0.091" (0.23 cm) | 0.062" (0.16 cm) | 0.250" (0.64 cm) | PCTFE | 1.7 µL | ea. |
| A-703 | 0.5 µm | 0.092" (0.23 cm) | 0.062" (0.16 cm) | 0.250" (0.64 cm) | PCTFE | 1.4 µL | ea. |
| A-704 | 2 µm | 0.125" (0.32 cm) | 0.062" (0.16 cm) | 0.250" (0.64 cm) | PCTFE | 3.0 µL | ea. |
| A-706 | 2 µm | 0.188" (0.48 cm) | 0.062" (0.16 cm) | 0.250" (0.64 cm) | PCTFE | 7.1 µL | ea. |
| A-707 | 0.5 µm | 0.195" (0.5 cm) | 0.062" (0.16 cm) | 0.250" (0.64 cm) | PCTFE | 6.1 µL | ea. |
| A-708 | 2 µm | 0.062" (0.16 cm) | 0.062" (0.16 cm) | 0.200" (0.51 cm) | PCTFE | 0.7 µL | ea. |
| A-710 | 2 µm | 0.125" (0.32 cm) | 0.062" (0.16 cm) | 0.200" (0.51 cm) | PCTFE | 3.0 µL | ea. |
| SEMI-PREP PE | EK FRITS | | | | | | |
| A-720 | 10 µm | 0.125" (0.32 cm) | 0.070" (0.18 cm) | 0.250" (0.64 cm) | PCTFE | 4.2 µL | ea. |
| A-722 | 10 µm | 0.197" (0.5 cm) | 0.070" (0.18 cm) | 0.250" (0.64 cm) | PCTFE | 9.9 µL | ea. |
| OC-802 | 2 µm | 0.460" (1.17 cm) | 0.070" (0.18 cm) | 0.560" (1.42 cm) | PCTFE | 46.4 µL | ea. |
| OC-803 | 10 µm | 0.460" (1.17 cm) | 0.072" (0.18 cm) | 0.560" (1.42 cm) | PCTFE | 57.2 μL | ea. |
| OC-805 | 5 µm | 0.460" (1.17 cm) | 0.058" (0.15 cm) | 0.560" (1.42 cm) | PCTFE | 41.1 µL | ea. |
| OC-815 | 5 um | 0.293" (0.74 cm) | 0.062" (0.16 cm) | 0.375" (0.95 cm) | PCTFE | 17.8 µL | ea. |

HROM = 1 y tic +61(0)3 9762 2034 ECH no logy Pty Ltd nail : info@chromtech.net.au Tel: 03 9762 2034 . . . in AUSTRALIA



Frit-in-a-Ferrule™

- > Seals and filters simultaneously
- Less expensive and more convenient than traditional inline filter systems
- Available in both Flangeless and Super Flangeless™ versions

Now you can filter at any point in your system where 1/16" or 1/8" OD tubing is used in a flat-bottom 1/4-28, M6 or 5/16-24 connection.

Our Frit-In-A-Ferrule product line is designed to seal and filter simultaneously by incorporating a frit into the body of a flat-bottom ferrule. This simple design allows you to eliminate traditional inline filters and reduce the number of additional connections in your system.





P-372 Flangeless Frit-In-A-Ferrule for 1/8" OD tubing

P-276Super Flangeless Frit-In-A-Ferrule for 1/16" OD tubing



| Part No. | Description | Porosity | Frit Material | Frit Diameter | Frit Thickness | Swept Volume | Maximum Pressure | Qty. |
|-----------|--|----------|---------------|---------------|----------------|--------------|---------------------|------|
| FRIT-IN- | A-FERRULE FOR 1/16" OD TUBING | | | | | | | |
| P-270 | Super Flangeless, Natural PEEK, SST lock ring | 2 µm | SST | 0.062" | 0.062" | 0.74 μL | 2,500 psi (172 bar) | ea. |
| P-272 | Flangeless, Green PCTFE | 2 µm | SST | 0.062" | 0.062" | 0.74 μL | 2,000 psi (138 bar) | ea. |
| P-273 | Flangeless, Blue PCTFE | 0.5 µm | SST | 0.062" | 0.062" | 0.61 μL | 2,000 psi (138 bar) | ea. |
| P-274 | Super Flangeless, Natural PEEK, SST lock ring | 2 µm | PEEK | 0.046" | 0.030" | 0.20 μL | 2,500 psi (172 bar) | ea. |
| P-275 | Super Flangeless, Black PEEK, SST lock ring | 0.5 µm | PEEK | 0.046" | 0.030" | 0.16 μL | 2,500 psi (172 bar) | ea. |
| P-276 | Super Flangeless, Stainless Steel, Natural ETFE, SST lock ring | 10 µm | SST | 0.062" | 0.062" | 0.90 μL | 2,500 psi (172 bar) | ea. |
| FRIT-IN- | A-FERRULE FOR 1/8" OD TUBING | | | | | | | |
| P-372 | Flangeless, Green PCTFE | 2 µm | SST | 0.094" | 0.062" | 1.69 µL | 500 psi (34 bar) | ea. |
| P-373 | Flangeless, Blue PCTFE | 0.5 µm | SST | 0.094" | 0.062" | 1.41 µL | 500 psi (34 bar) | ea. |
| P-374 | Super Flangeless**, Natural PEEK, SST lock ring | 2 µm | PEEK | 0.094" | 0.042" | 1.15 μL | 2,500 psi (172 bar) | ea. |
| * Swept v | olumes include/reflect theoretical frit volume values. | | | | | | | |

** The 1/8" Super Flangeless versions cannot be used in M6 ports.







Bottom-of-the-Bottle[®] Filters

Our uniquely designed Bottom-of-the-Bottle™ Filters effectively protect your system by filtering out particulate matter that my otherwise damage expensive hardware.

1.70" (4.32 cm) 0.50" (1.27 cm)

sparging gas (optional use)

Stainless Steel Bottom-of-the-Bottle Solvent Filters

- Draws solvent from within 1/8" of the bottom of the bottle
- Replaceable stainless steel filter cups
- > Versions for 1/8" and 3/16" OD tubing
- Materials of construction: PEEK, ETFE, and 316 Stainless Steel

Patented Stainless Steel Bottom-of-the-Bottle Solvent Filter Assemblies feature a 2 µm or 10 µm replaceable stainless steel filter cup and a design that allows solvent to be drawn from within 1/8" of the bottom of your solvent bottle. The filter cups are inexpensive and easy to replace, making this an economical, trouble-free choice.

Port for solvent intake line Port for your sparging line (optional use) Filters and Filters solvent disperses incoming from the

Bottom-of-the-Bottle Inlet Solvent Filter

Maximum Flow Rates: 2 μm-up to 10 mL/min. 10 μm-up to 40 mL/min


bottom of your solvent bottle

All-PEEK Bottom-of-the-Bottle Solvent Filters

- Most recommended filtering unit
- ▶ 100% PEEK polymer construction
- > Easy operation no fittings required

These biocompatible filters are made from 100% PEEK polymer, including the two built-in PEEK frits. The bottom frit (2 µm or 10 µm) will draw solvents from within 0.080'' (2.0 mm) of the bottom of the solvent bottle. The 2 μ m frit on the side may be used for a 1/8" OD helium sparging line.

To use, simply press fit your appropriately sized fluoropolymer tubing firmly into the top holes. That's it!



Maximum Flow Rate: up to 30 mL/min




UHMWPE Bottom-of-the-Bottle™ Solvent Filters

- > Replaceable filter cup
- Economical
- Materials of construction: UHMWPE, ETFE
- > Versions for 1/16" and 1/8" OD tubing

The design of the UHMWPE solvent filters allows tubing to pass through to the bottom of the filter cup, enabling the filter to draw solvent from within 0.10" (2.5 mm) of the bottom of your solvent bottle.

Please Note: UHMWPE is a hydrophobic material. To establish proper surface wetting, you may need to prime the filter with methanol or acetonitrile.



Bottom-of the-Bottle Filters

| Part No. | Description | Porosity For Tubing Size | | Includes | Qty. |
|-------------|--|--------------------------|----------|------------|------|
| STAINLESS S | TEEL BOTTOM-OF-THE-BOTTLE SOLVENT | FILTERS | | | |
| A-550 | SST Filter Assembly, with A-520 filter cup | 10 µm | 1/8" OD | (1) XP-130 | ea. |
| A-551 | SST Filter Assembly, with A-522 filter cup | 2 µm | 1/8" OD | (1) XP-130 | ea. |
| A-520x | SST Replacement Solvent Filter Cups, 10-pk | 10 µm | _ | _ | ea. |
| A-522x | SST Replacement Solvent Filter Cups, 10-pk | 2 µm — | | _ | ea. |
| ALL-PEEK BI | OCOMPATIBLE BOTTOM-OF-THE-BOTTLE S | SOLVENT | FILTERS | | |
| A-435 | PEEK Filter | 2 µm | 1/8" OD | _ | ea. |
| A-437 | PEEK Filter, for small-neck (GL-38) bottles | 2 µm | 1/8" OD | _ | ea. |
| A-438 | PEEK Filter, for small-neck (GL-38) bottles | 10 µm | 1/8" OD | _ | ea. |
| A-440 | PEEK Filter | 10 µm | 1/8" OD | _ | ea. |
| A-441 | PEEK Filter | 10 µm | 3/16" OD | _ | ea. |
| A-451 | PEEK Filter | 10 µm | 1/16" OD | _ | ea. |
| UHMWPE BIO | OCOMPATIBLE BOTTOM-OF-THE-BOTTLE S | SOLVENT | FILTERS | | |
| A-445 | UHMWPE Filter Assembly | 10 µm | 1/16" OD | (1) XP-245 | ea. |
| A-446 | UHMWPE Filter Assembly | 10 µm | 1/8" OD | (1) XP-345 | ea. |
| A-427 | UHMWPF Replacement Solvent Filter Cups. 5-pk | 10 um | _ | _ | ea. |



Inlet Solvent Filters

- Large surface areas prevent pump cavitation
- Disposable
- 2 μm, 10 μm, and 20 μm pore sizes available
- General use and prep filters for higher flow applications



APPLICATION NOTE

Why Use An Inlet Solvent Filter?

- To filter out particulate matter from the solvent that may otherwise damage expensive hardware. (Use a 10 µm or 20 µm version for this purpose. The A-309 and A-230A filters have an added "Bottom of the Bottle™" feature to help draw solvent to within 1/8" of the bottom of your solvent bottle.)
- To prevent particulates originating from the sparging system from entering the mobile phase reservoir and to help disperse the sparging gas efficiently. (Use a 2 µm filter for this purpose.)
- To hold your tubing in place at the bottom of the bottle. (Most stainless steel filter options work best for this purpose.)

Note: It is usually a good idea to change the inlet filter as part of your semi-annual or annual preventative maintenance program. It is good practice to filter your solvents to prevent pump damage. Our 316 stainless steel filters provide that protection.

Because filters should be changed periodically, we make it easy to replace them without tools. For those filters using a plastic fitting, the tubing can be reconnected by finger tightening the fitting into the new filter. The filters with stems allow easy insertion into the inlet tubing.



| Part No. | Description | Porosity | Material | For Tubing Size | Includes | Max. Suggested Flow Rate* | Qty. |
|----------|---|-----------|------------|-----------------|-----------------------|---------------------------|------|
| INLET SO | LVENT FILTERS FOR ANALYTICAL HPLO | | | | | | |
| A-242 | Inlet Solvent Filter with One-Piece Fitting | 2 µm | PCTFE, SST | 1/8" OD | (1) P-100 | 10 mL/min | ea. |
| A-243 | A-242, 5-pack | 2 µm | PCTFE, SST | 1/8" OD | (5) P-100 | 10 mL/min | ea. |
| A-228 | Inlet Solvent Filter with stem | 2 µm | SST | 1/8" ID | _ | 80 mL/min | ea. |
| A-302 | Inlet Solvent Filter with stem | 10 µm | SST | 1/16" ID | _ | 40 mL/min | ea. |
| A-302A | Inlet Solvent Filter with Flangeless Fittings | 10 µm | PCTFE, SST | 1/8" OD | (1) XP-315 | 40 mL/min | ea. |
| A-309 | Inlet Solvent Filter with stem | 10 µm | SST | 1/16" ID | _ | 40 mL/min | ea. |
| A-231A | Inlet Solvent Filter with Flangeless Fittings | 20 µm | PCTFE, SST | 3/16" OD | (1) XP-132 | 100 mL/min | ea. |
| A-310 | Inlet Solvent Filter with stem | 10 µm | SST | 1/8" ID | _ | 40 mL/min | ea. |
| INLET SO | LVENT FILTERS FOR PREPARATIVE HPL | C SYSTEMS | | | | | |
| A-225 | Inlet Solvent Filter with stem | 20 µm | SST | 1/16" ID | _ | 100 mL/min | ea. |
| A-225A | Inlet Solvent Filter with Flangeless Fittings | 20 µm | PCTFE, SST | 1/8" OD | (1) P-315, (1) P-300N | 100 mL/min | ea. |
| A-227A | Inlet Solvent Filter with Flangeless Fittings | 10 µm | PCTFE, SST | 1/4" OD | (1) XU-655 | 100 mL/min | ea. |
| A-230A | Inlet Solvent Filter with Flangeless Fittings | 20 µm | PCTFE, SST | 1/4" OD | (1) XU-655 | 100 mL/min | ea. |
| A-311 | Inlet Solvent Filter with stem | 10 µm | SST | 1/16" ID | _ | 100 mL/min | ea. |
| A-311A | Inlet Solvent Filter with Flangeless Fittings | 10 µm | PCTFE, SST | 1/8" OD | (1) XP-315 | 100 mL/min | ea. |
| | . 10 | , , | | | | | |





Inline Filters

- > Specially engineered for inline filtration
- Versions include Micro, Standard, and Semi-Preparative
- Bio-inert and stainless steel options offered
- Variety of porosities, application appropriate




Fittings

All Standard Inline Solvent Filters have 10-32 threads for 1/16" OD tubing, allowing the use of most standard chromatography high pressure fittings. Our Inline Filters are specially engineered for inline filtration. It is specifically designed to help prevent particulate contamination from clogging sensitive equipment. It is ideally suited for placement along the flow path line between the pump and injection valve/autosampler. We offer a variety of porosities for your application.

Standard Inline Solvent Filters

- For 1/16" OD tubing
- Versions for Standard HPLC (6,000 psi/414 bar) and UHPLC (25,000 psi/1,725 bar)
- Replacement frits availableVersions for Standard HPLC (6,000 psi/414 bar) and UHPLC (25,000 psi/1,725 bar)
- > Help prevent particulate contamination from clogging sensitive equipment
- Ideally suited for placement along the flow path line between the pump and injection valve/autosampler

Inline filter assemblies that begin with the letter "A" are engineered for standard HPLC applications (up to 6,000 psi/414 bar). Inline Filter Assemblies that begin with the "VHP" prefix are suitable for use in UHPLC systems, where pressures can reach 25,000 psi (1,725 bar).






104

Inline Filters (Cont.)

Biocompatible Standard Inline Filters

- ▶ 0.5 µm and 2 µm versions available
- > Features 100% PEEK flow path

Our A-430 and A-431 Inline Filters consist of a stainless steel body and two PEEK end fittings. Maximum recommended flow rate is 25 mL/min for the A-430 Filter and 10 mL/min for the A-431 Filter. And, you get the added benefit of biocompatibility since all wetted surfaces are PEEK. When you need to replace the frit, simply dispose of the end fitting that contains the frit and replace it with a new one.



Inline Filters

| Part No. | Description | Porosity | For Tubing Size | Threads | Includes | Swept Volume | Pressure Rating | Qty. |
|-----------|---|----------|-----------------|-------------|-------------|--------------|------------------------|-------|
| STANDA | RD INLINE SOLVENT FILTERS | | | | | | | |
| A-313 | Solvent Filter Assembly | 20 µm | 1/16" OD | 10-32 Coned | (1) A-224 | 12.3 µL | 6,000 psi (414 bar) | ea. |
| A-314 | Solvent Filter Assembly | 2 μm | 1/16" OD | 10-32 Coned | (1) A-100 | 4 μL | 6,000 psi (414 bar) | ea. |
| A-100 | Replacement Frits, Stainless Steel, ea. | 2 μm | N/A | _ | _ | 1.4 µL | N/A | ea. |
| A-224 | Replacement Frits, Stainless Steel, ea. | 20 µm | N/A | _ | _ | 9.7 μL | N/A | ea. |
| VHP-500 | Inline VHP Filter | 0.5 µm | 1/16" OD | 10-32 Coned | (5) VHP-501 | 1.2 µL | 25,000 psi (1,725 bar) | ea. |
| VHP-505 | Inline VHP Filter | 0.2 µm | 1/16" OD | 10-32 Coned | (5) VHP-506 | 1.1 µL | 25,000 psi (1,725 bar) | ea. |
| VHP-501 | Replacement Inline VHP Frit | 0.5 µm | N/A | N/A | N/A | 0.60 μL | N/A | ea. |
| VHP-506 | Replacement Inline VHP Frit | 0.2 µm | N/A | N/A | N/A | 0.54 μL | N/A | ea. |
| BIOCOM | IPATIBLE INLINE FILTERS | | | | | | | |
| A-430 | Biocompatible Filter Assembly | 2 µm | | 10-32 Coned | (1) A-429 | 7.1 µL | 6,000 psi (414 bar) | ea. |
| A-431 | Biocompatible Filter Assembly | 0.5 µm | | 10-32 Coned | (1) A-428 | 5.9 µL | 6,000 psi (414 bar) | ea. |
| A-428x | PEEK Filter End Fittings, Black PEEK body, 10-pk | 0.5 µm | | 10-32 Coned | _ | 5.7 µL | N/A | 10-pk |
| A-429x | PEEK Filter End Fittings, Natural PEEK body, 10-pk | 2 µm | | 10-32 Coned | _ | 6.9 µL | N/A | 10-pk |
| *Swept vo | lumes include/reflect theoretical frit volume values. | | | | | | | |

*Swept volumes include/reflect theoretical frit volume values SST = Stainless Steel

ECH no logy Pty Ltd nail : info@chromtech.net.au Tel: 03 9762 2034 . . . in AUSTRALIA

Australian Distributors Importers & Manufacurers www.chromtech.net.au



Precolumn Filters


Our economical Precolumn Filters offer secure protection for analytical columns in HPLC and UHPLC. We offer traditional versions that can successfully connect tubing on both sides and our direct-connect versions attach to the inlet port of most standard columns. All versions feature a 10-32 coned ports for 1/16" OD tubing.

Standard Precolumn Filters

- > Economical protection for larger columns and injections
- Traditional versions connect tubing on both sides
- Direct-connect versions attach to the inlet port of most standard columns
- All versions feature 10-32 coned ports for 1/16" OD tubing

These are designed to protect columns by filtering out particulate matter originating from the sample or from rotor seal wear.

- Assemblies that begin with the letter "A" are traditional versions for standard HPLC
- Assemblies that begin with "VHP" are direct-connect versions for UHPLC applications
- Versions that begin with "9085" are direct-connect for standard HPLC and must be used with polymer fittings




106

Precolumn Filters (Cont.)

Biocompatible Precolumn Filters

- > Pre-assembled with either 0.5 μm or 2 μm porosity frits
- Great column protection
- > Feature PEEK bodies and PCTFE-surrounded PEEK frits

Biocompatible Precolumn Filters have 0.020" (0.50 mm) diameter thru-holes and 8° distribution cones for minimal band spreading and mixing. The bodies of these filters are manufactured from biocompatible PEEK polymer and are pressure rated to 5,000 psi (345 bar). These filters are designed for use with 1/16" OD tubing, which can be connected to these filters using standard Fingertight fittings.



Precolumn Filters

| Part No. | Description | Porosity | For Tubing Size | Threads | Includes | Swept Volume* | Pressure Rating | Qty. |
|------------|---|----------|-----------------|-------------|-------------|---------------|------------------------|-------|
| STANDAR | D PRECOLUMN FILTERS | | | | | | | |
| A-315 | Solvent Filter Assembly | 2 µm | 1/16" OD | 10-32 Coned | (1) A-101 | 1.4 µL | 6,000 psi (414 bar) | ea. |
| A-318 | Solvent Filter Assembly | 0.5 µm | 1/16" OD | 10-32 Coned | (1) A-102 | 0.84 µL | 6,000 psi (414 bar) | ea. |
| A-101 | Replacement Frits, Stainless Steel, ea. | 2 µm | N/A | _ | _ | 0.74 µL | N/A | ea. |
| A-102 | Replacement Frits, Stainless Steel, ea. | 0.5 µm | N/A | _ | _ | 0.61 μL | N/A | ea. |
| VHP-550 | Precolumn VHP Filter | 0.5 µm | 1/16" OD | 10-32 Coned | (5) VHP-551 | 1.9 µL | 20,000 psi (1,380 bar) | ea. |
| VHP-555 | Precolumn VHP Filter | 0.2 µm | 1/16" OD | 10-32 Coned | (5) VHP-556 | 1.8 µL | 20,000 psi (1,380 bar) | ea. |
| VHP-551 | Replacement Precolumn VHP Frit Assembly | 0.5 µm | N/A | N/A | N/A | 1.9 µL | N/A | ea. |
| VHP-556 | Replacement Precolumn VHP Frit Assembly | 0.2 µm | N/A | N/A | N/A | 1.8 µL | N/A | ea. |
| 9085-05-10 | ColumnSaver Precolumn Filter, with SST frit | 0.5 µm | 1/16" OD | 10-32 Coned | N/A | 3.1 µL | 6,000 psi (414 bar) | 10-pk |
| 9085-20-10 | ColumnSaver Precolumn Filter, with SST frit | 2 µm | 1/16" OD | 10-32 Coned | N/A | 3.1 µL | 6,000 psi (414 bar) | 10-pk |
| BIOCOMP | ATIBLE PRECOLUMN FILTERS | | | | | | | |
| A-355 | Solvent Filter Assembly, Biocompatible | 2 µm | | 10-32 Coned | (1) A-700 | 1.4 µL | 5,000 psi (345 bar) | ea. |
| A-356 | Solvent Filter Assembly, Biocompatible | 0.5 µm | | 10-32 Coned | (1) A-701 | 1.3 µL | 5,000 psi (345 bar) | ea. |
| A-700 | Replacement Frit, PEEK Polymer | 2 µm | | _ | _ | 0.74 µL | N/A | ea. |
| A-701 | Replacement Frit, PEEK Polymer | 0.5 µm | | _ | _ | 0.61 µL | N/A | ea. |

*Swept volumes include/reflect theoretical frit volume values.






Semi-Prep Filters

Biocompatible Semi-Prep Filters consist of a stainless steel body, two PEEK end fittings, and a separate PEEK frit. These filters are ideal for many higher flow analytical, semi-prep and preparative applications. Best of all, if the filter becomes clogged, simply unscrew the assembly, remove the frit and replace it. The frits are interchangeable.

Semi-Prep Inline Filters

- Designed for high-flow applications
- Economical protection for larger columns and injections
- > SFC and HPLC compatible





Iso-Prep Filter Shown with standard 10-32 stainless steel nuts and ferrules (not included)



Biocompatible Semi-Prep Inline Filters

- Versions for 1/16", 1/8", 3/16", 1/4", and 5/16" OD tubing
- ▶ 100% PEEK flow path

Biocompatible Semi-Prep Filters consist of a stainless steel body, two PEEK end fittings, and a separate PEEK frit. These filters are ideal for many higher flow analytical, semi-prep and preparative applications. Best of all, if the filter becomes clogged, simply unscrew the assembly, remove the frit and replace it. The frits are interchangeable.

| Part No. | Description | Porosity | Threads | Includes | Swept Volume* | Pressure Rating | Qty. |
|-----------|---|----------|---------------------|-------------|---------------|---------------------|------|
| SEMI-PR | REP INLINE FILTERS | | | | | | |
| A-330 | Semi-Prep Filter Assembly | 10 µm | 10-32 Coned | (1) A-331 | 223 μL | 7,500 psi (517 bar) | ea. |
| A-360 | Semi-Prep Filter Assembly | 10 µm | 5/16-24 Flat Bottom | (1) A-331 | 235 μL | 3,500 psi (207 bar) | ea. |
| A-331 | Stainless Steel Frits, Natural ETFE ring | 10 µm | N/A | N/A | 142 µL | N/A | ea. |
| A-332 | Stainless Steel Frits, Natural ETFE ring | 2 µm | N/A | N/A | 122 µL | N/A | ea. |
| A-337 | Stainless Steel Frits, Natural ETFE ring | 20 µm | N/A N/A 152 µL N/A | | ea. | | |
| ISO-PRE | P FILTERS | | | | | | |
| 1020-05 | 21.2 mm Filter Holder | 0.5 µm | 10-32 Coned | (1) 7031-05 | 203 uL | 8,000 psi (552 bar) | ea. |
| 1020-20 | 21.2 mm Filter Holder | 2 µm | 10-32 Coned | (1) 7031-20 | 196 uL | 8,000 psi (552 bar) | ea. |
| 7031-05 | 21.2 mm Replacement Filter | 0.5 µm | N/A | N/A | 122 uL | 8,000 psi (552 bar) | ea. |
| 7031-20 | 21.2 mm Replacement Filter | 2 µm | N/A | N/A | 115 uL | 8,000 psi (552 bar) | ea. |
| BIOCON | MPATIBLE SEMI-PREP INLINE FILTERS | | | | | | |
| A-410 | Biocompatible Filter Assembly | 2 µm | 10-32 Coned | (1) OC-802 | 89 µL | 6,000 psi (414 bar) | ea. |
| A-411 | Biocompatible Filter Assembly | 10 µm | 10-32 Coned | (1) OC-803 | 103 μL | 6,000 psi (414 bar) | ea. |
| A-510 | Biocompatible Filter Assembly | 5 µm | 5/16-24 Flat Bottom | (1) OC-805 | 89 µL | 500 psi (34 bar) | ea. |
| OC-802 | PEEK Frit, Green PCTFE ring | 2 µm | N/A | N/A | 46 µL | N/A | ea. |
| OC-803 | PEEK Frit, Natural PCTFE ring | 10 µm | N/A | N/A | 57 μL | N/A | ea. |
| OC-805 | PEEK Frit, Natural PCTFE ring | 5 µm | N/A | N/A | 50 μL | N/A | ea. |
| *Swept vo | lumes include/reflect theoretical frit volume values. | | | | | | |

Swept volumes include/reflect theoretical frit volume values.








Inline MicroFilters

- > 100% biocompatible PEEK polymer option available
- Miniscule 240 nL void volume
- Two versions: direct connect 1/32" OD tubing or use MicroTight® tubing sleeves for 70-520 µm OD capillary tubing

Our Inline MicroFilters protect your column from particles originating in the mobile phase or sample, or from pump seal and sample injection valve wear. These filters have a 0.006" (150 μ m) thru-hole. Choose the M-520 with a $0.5~\mu$ m 100% PEEK frit to connect to capillary tubing using the MicroTight tubing sleeves (page 52). You may also directly connect 1/32" OD tubing using the M-525 which contains a 0.5 µm PEEK frit.







| Part No. | Description | Porosity | For Tubing Size | | Threads | Includes | Swept Volume | Pressure Rating | Qty. |
|----------|---|-----------|-------------------|------------|--------------------------|----------------------|--------------|---------------------|-------|
| INLINE I | MICROFILTERS | | | | | | | | |
| M-520 | Inline MicroFilter Assembly, PEEK Frit | 0.5 µm | MicroTight Tubing | Sleeve | MicroTight Tubing Sleeve | (5) M-120, (2) F-125 | 240 nL | 4,000 psi (276 bar) | ea. |
| M-525 | Inline MicroFilter Assembly, PEEK Frit | 0.5 µm | 1/32" OD | | 1/32" OD | (5) M-140, (2) F-126 | 240 nL | 4,000 psi (276 bar) | ea. |
| REPLAC | EMENT INLINE MICROFILTER EN | ID-FITTIN | GS | | | | | | |
| M-120x | End-Fittings, Black, with PEEK Frit | 0.5 µm | MicroTight Tubing | Sleeve | MicroTight Tubing Sleeve | N/A | 216 nL | N/A | 10-pk |
| M-140x | End-Fittings. Natural, with PEFK Frit | 0.5 um | 1/32" OD | | 1/32" OD | N/A | 216 nL | N/A | 10-pk |
| | HROMalytic +61(0)3 9762 2034 ECH nology Pty Ltd Australian Distributors Importers & Manufacurers www.chromtech.net.au | | | | | | | | |
| Website | NEW: www.chromalytic.net.au E-ma | | | 108 (of 17 | 72) 2 ⁽ | | | | |

FLUIDICS



Mini MicroFilters

- Total volume as low as 10 nL
- Conductive version for CEC and mass spectrometry applications



APPLICATION NOTE

The Mini MicroFilters can be used to pack capillary tubing. Simply place one of these filters on the effluent side of the capillary tubing, then slurry pack. Once packed, place a filter at the head of the tubing. This creates a reliable capillary column without fusing the silica to make frits or pressing filter paper inside the capillary tubing.

Increase the Life of Your Column

Why use a Precolumn Filter when there is a frit at the head of the column itself? Changing the column frit is extremely difficult to do without disturbing the column packing. A Precolumn Filter provides relatively inexpensive insurance against column damage, and changing its frit is easy. A Precolumn Filter placed between the sample injection valve and the HPLC column protects the column from particles originating in the sample and from pump and valve seal wear.

Our Inline Mini MicroFilter Assemblies filter effectively with internal volumes low enough to ensure reliable chromatographic results — even at nanoliter per minute flow rates! Internal volumes of these encapsulated filters are as low as 85 nL with the micro-screen



and 10 nL to 22 nL with the frit disc option.



SPECIFICATIONS & DETAILS

Because of the size-specific nature of the ferrules included with each Mini MicroFilter assembly, please note that these ferrules are not interchangeable with other MicroFerrules for different tubing sizes.

Filter Capsule Color Identification











What's the Difference Between Precolumn & Inline Filters?

You may have noticed that the bodies of Precolumn and Inline Filters look similar, and as such, you may have wondered what the differences are. Because Precolumn Filters, by definition, are typically placed in a volume-sensitive area immediately preceding the column, these filters usually feature smaller thru-holes and smaller frit diameters. In contrast, Inline Filters are often placed where the internal volume is not as critical and where longer life and less fluid restriction is more important.

| Part No. | Description | Porosity | Frit Type | For use with Tubing | Includes | Swept Volume | Pressure Rating | Qty. | |
|----------|---------------------------------------|----------|------------|---------------------|---------------------------------|--------------|---------------------|------|--|
| MINI MI | MINI MICROFILTER ASSEMBLY | | | | | | | | |
| M-547 | Mini MicroFilter Assembly | 1 µm | SST Frit | 1/32" (790 µm) OD | (5) M-133, (2) F-112, (2) P-416 | 22 nL | 4,000 psi (276 bar) | ea. | |
| M-548 | Mini MicroFilter Assembly | 1 µm | Ti Frit | 1/32" (790 µm) OD | (5) M-134, (2) F-112, (2) P-416 | 22 nL | 4,000 psi (276 bar) | ea. | |
| REPLAC | REPLACEMENT MINI MICROFILTER CAPSULES | | | | | | | | |
| Part No. | Description | Porosity | Frit Type | For Use With | Material | Swept Volume | | Qty. | |
| M-121 | Filter Capsule | 1 µm | SST Screen | M-530 and M-531 | PEEK | 85 nL | | 2-pk | |
| M-125 | NanoFilter Capsule | 1 µm | SST Frit | M-537 and M-538 | PEEK | 10 nL | | 2-pk | |
| M-126 | NanoFilter Capsule | 1 µm | Ti Frit | M-537 and M-538 | PEEK | 10 nL | | 2-pk | |
| M-131 | Filter Capsule | 1 µm | SST Screen | M-543 | PEEK | 85 nL | | 2-pk | |
| M-133 | NanoFilter Capsule | 1 µm | SST Frit | M-547 and M-548 | PEEK | 10 nL | | 2-pk | |
| M-134 | NanoFilter Capsule | 1 µm | Ti Frit | M-547 and M-548 | PEEK | 10 nL | | 2-pk | |
| M-128 | Conductive NanoFilter Capsule | 1 µm | SST Frit | M-534 | SST/PEEK | 10 nL | | 2-pk | |
| | | | | | | | | | |






Precolumn MicroFilters

- Direct connects to columns with 10-32 threads
- Total void volume of 0.5 µL
- > Two versions: direct connect 1/16" OD tubing or use MicroTight® tubing sleeves for 70-520 µm OD capillary tubing

The Precolumn MicroFilters directly connect into your microbore or analytical column. Total theoretical void volume is only $0.5\,\mu L$ (includes frit volume) and the PEEK tubing used in the assembly of these units has a 0.005" (125 µm) ID, virtually eliminating any mixing of the sample with the mobile phase.



Precolumn MicroFilter Assembly (Includes indicated products)

| Part No. | Description | Porosity | For Tubing Size | Threads | Includes | Swept Volume* | Pressure Rating | Qty. |
|----------|---|----------|--------------------------|-------------|---------------------------------|---------------|---------------------|------|
| PRECO | LUMN MICROFILTER ASSEMBLIES | | | | | | | |
| M-500 | Precolumn MicroFilter Assembly, SST Frit | 0.5 µm | MicroTight Tubing Sleeve | 10-32 Coned | (5) C-425, (1) F-172, (1) P-416 | 0.5 μL | 4,000 psi (276 bar) | ea. |
| M-510 | Precolumn MicroFilter Assembly, PEEK Frit | 0.5 µm | MicroTight Tubing Sleeve | 10-32 Coned | (5) A-735, (1) F-172, (1) P-416 | 0.5 μL | 4,000 psi (276 bar) | ea. |
| M-550 | Precolumn MicroFilter Assembly, SST Frit | 0.5 µm | 1/16" OD | 10-32 Coned | (5) C-425, (1) F-132, (1) P-416 | 0.5 μL | 4,000 psi (276 bar) | ea. |
| M-560 | Precolumn MicroFilter Assembly, PEEK Frit | 0.5 µm | 1/16" OD | 10-32 Coned | (5) A-735, (1) F-132, (1) P-416 | 0.5 μL | 4,000 psi (276 bar) | ea. |
| REPLAC | REPLACEMENT PRECOLUMN MICROFILTER FRITS (FRIT DIAMETER X FRIT THICKNESS X OVERALL DIAMETER) | | | | | | | |
| A-735 | PEEK Frits, 0.045" x 0.031" x 0.192" | 0.5 µm | N/A | N/A | N/A | 216 nL | N/A | ea. |
| C-420 | SST Frits, 0.038" x 0.028" x 0.192" | 2 µm | N/A | N/A | N/A | 101 nL | N/A | ea. |
| C-425 | SST Frits. 0.038" x 0.028" x 0.192" | 0.5 um | N/A | N/A | N/A | 101 nL | N/A | ea. |



Bottle Caps

- Extremely simple no threaded ports or fittings
- Manufactured from ETFE and Polypropylene




- A self-regulating sparging system can help reduce helium consumption and improve pump performance. Set this up by pressing your tubing through the appropriate holes in your bottle cap and attaching each line to a filter. Sparge your mobile phase with an inert gas (preferably helium) for 15-20 minutes. Then reduce the outlet pressure of the sparging gas to a maximum of 5 psi (0.34 bar) and insert a plug (A-626 or A-628) into the remaining port of the cap. The sparging gas will shut off once the incoming pressure equals the pressure inside the reservoir. As the mobile phase is consumed and the internal pressure lowers, sparging gas will enter to keep the system pressurized and degassed. Please Note: If gas leaks while pressurizing the bottle, try removing the sealing ring from the bottle, as it sometimes interferes with the sealing of these bottle caps.
- One concern with sparging systems is the possibility of solvent backing up the sparging inlet line. This can occur if the gas tank completely evacuates with the regulating valves open, creating a vacuum in the tubing. Solvent backup may damage sparging system components and cause cross-contamination of mobile phase reservoirs. To help prevent solvent backup, install the CV-3010 Inline Check Valve (page 135) along the tubing line that runs between the gas supply and the solvent bottle.
- For a more efficient degassing system, please see the HPLC Vacuum Degassing Systems on page 154.
- Please see the Quick-Stop Luer Check Valve on page 139 for another solvent inlet Application Note.

If you are looking for a bottle cap that is quick and easy to use, but still allows many connect ion options, we have just what you need! The Bottle Caps fit standard GL-45 (1 L) or smaller-neck GL-38 (4 L) glass bottles.

Each cap has three holes. With two of the holes you simply push your tubing straight through. The third hole, with a luer taper, can be used for a number of options. Any male luer (such as a luer-lock syringe) will fit snugly in this hole, or you can use the A-626 or A-627 Plug. Exceptions are the A-610 Bottle Caps. Please see the note below.







A-622 GL-38 Bottle Cap for 1/8" tubing



NOTE

The A-610 Bottle Cap has a slightly different configuration than other caps. One hole accepts 3/16" OD tubing, the typical size used with some Waters® systems. The remaining two holes accept 1/8" OD tubing. Unlike the other caps, the A-610 does not have a tapered luer hole. If desired, use our A-628 Plug or A-629 Filter Plug for one of the 1/8" holes.



To ensure a tight seal, use fluoropolymer tubing with these bottle caps (page 55).



A-622 for 1/8" OD tubina. Black

Australian Distributors
Importers & Manufacurers
Importers & Manufacurers
Website NEW: www.chromalytic.net.au E-mail: info@chromtech.net.au Tel: 03 9762 2034 . . . in AUSTRALIA

www.idex-hs.com
111 (of 172) 2017Sept



Bottle Cap Plugs

Use the A-626 Bottle Cap Plug to seal the third "tapered" luer hole found in most IDEX Health & Science Bottle Caps. Or, use the A-628 Plug to seal any unused 1/16" or 1/8" bottle cap holes.

Alternatively, try the A-627 or A-629 Filter Bottle Cap Plug to cap an unused hole in your bottle cap. The 20 µm stainless steel frit in these products prevents foreign matter from contaminating your solvent while leaving the bottle open to the atmosphere, thus allowing fluid to be pulled out without creating a vacuum (generally not used with sparging applications). All plug bodies are manufactured from ultra-high molecular weight polyethylene (UHMWPE).



Bottle Cap Plug




A-629 Filter Bottle Cap Plug

| Part No. | Description | | | Qty. |
|-----------|---|--|----|-----------------|
| BOTTLE CA | P PLUGS | | | |
| A-626 | Bottle Cap Plug for luer hole, UHMWPE | | | ea. |
| A-627 | Filter Bottle Cap Plug for luer hole, UHMWPE with 20 μm | n stainless steel frit | | ea. |
| A-628 | Bottle Cap Plug for 1/16", 1/8" or 3/16" hole, UHMWPE | | | ea. |
| A-629 | Filter Bottle Cap Plug for 1/16". 1/8" or 3/16" hole. UHN | IWPF with 20 um stainless steel frit | r | ea. |
| | HROM=1541C +61(0)3 9762 2034 | Australian Distributors Importers & Manufacurers | GS | www.iday.hs.com |

Website NEW: www.chromalytic.net.au E-mail: info@chromtech.net.au Tel: 03 9762 2034 . . . in AUSTRALIA

112 (of 172) 2017Sept



Valves

Our valves are an integral part of advanced fluid-handling solutions for a wide range of analytical instrumentation and clinical diagnostic systems. Our valve options include manual valves for lower frequency use and rotary shear valves that meet the high duty cycle requirements of UHPLC and also come in high and low pressure versions to meet your system requirements. We also offer check valves when there is a need to limit the fluid flow to one direction. Our Back Pressure Regulators products are designed to enhance system performance through outgassing prevention. All of our valve products, components, tools, and accessories are designed keeping our customer's system needs first.

115 VALVE OVERVIEW & FUNCTIONS

115 ROTARY SHEAR VALVES

119 STAND ALONE VALVE PRODUCTS

121 MANUAL VALVES

124 SPARE PARTS

126 VALVE ACCESSORIES

134 FLOW REGULATING VALVES

143 BACK PRESSURE REGULATORS

HRON a y ti c +61(0)3 9762 2034

ECH 100 9 Y Pty Ltd

Website NEW: www.chromalytic.net.au E-mail: info@chromtech.net.au Tel: 03 9762 2034 . . . in AUSTRALIA

VALVE OVERVIEW & FUNCTIONS

| Valve Module | Flow Configurations | Page |
|------------------|---|------|
| ACTUATED VALVES | | |
| UP TO 15,000 PSI | | 119 |
| Switching | 2-Position, 6-Port 2-Position, 10-Port | |
| Injection | For Injection, add the appropriately sized Sample Loop to the Switching valves above | |
| Selection | 6-Position, 7-Port | |
| UP TO 6,000 PSI | | 119 |
| Switching | 2-Position, 6-Port (Analytical and Nano Scale) 2-Position, 10-Port (Analytical and Nano Scale) | |
| Injection | For Injection, add the appropriately sized Sample Loop to the Switching valves above 2-Position, 6-Port (vertical port) | |
| Selection | 6-Position, 7-Port | |
| UP TO 125 PSI | | 119 |
| Switching | 2-Position, 6-Port 2-Position, 6-Port (Double 3-Way) | |
| Selection | 6-Position, 7-Port 10-Position, 11-Port | |

| Valve Module | Flow Configurations | Page |
|-----------------|---|------|
| MANUAL VALVES | | |
| UP TO 9,000 PSI | | 123 |
| Injection | 2-Position, 6-Port (Front-Loading, 9,000 psi) | |
| UP TO 6,000 PSI | | 123 |
| Switching | 2-Position, 6-Port (Analytical and Micro Scale) | |
| Injection | 2-Position, 6-Port | |
| Selection | 6-Position, 7-Port | |
| UP TO 1,000 PSI | | 123 |
| Switching | 2-Way, Right Angle 4-Position, 4-Port 3-Way, T-Shape 4-Position, 4-Port 4-Way, Diagonal Flow 4-Position, 4 Port | |
| Injection | 2-Position, 6-Port | |
| Selection | 6-Position, 7-Port | |

Rotary Shear Valves

Our Rotary Shear Valves were developed in tandem with the evolution of liquid chromatography, where combinations of elevated system pressures, aggressive chemicals, and ever-diminishing fluid volumes continually challenged system manufacturers who required highly precise fluid control and delivery. Today, many other disciplines utilize Rotary Shear Valves for their versatility, reliability, repeatability, long system uptime, and easy preventive maintenance.



Valve Overview and Functions (Cont.)

Choosing a Rotary Shear Valve

Evaluating some simple variables will assist you in choosing the best valve for your needs.

Identify the Operating Pressure of Your Instrument or Application

Valves are designed to repeatedly deliver specific fluids to different locations in a fluidic circuit. Achieving fluidic precision at 15,000 psi requires different valve-design features than those required to achieve fluidic precision at 100 psi. A wide variety of variables such as valve architectures, metals, polymers, coatings, actuation speeds, and manufacturing techniques have been tested to achieve the fluidic accuracy and precision required for the full array of pressure conditions in life science applications. In this catalog, we define four separate pressure groupings:

| Up to 15,000 psi (1,035 bar) | UHPLC/Fast Chromatography |
|------------------------------|---------------------------------------|
| Up to 6,000 psi (410 bar) | HPLC |
| Up to 1,000 psi (69 bar) | Medium Pressure Applications |
| Up to 125 psi (8.5 bar) | Low Pressure/ Atmospheric Pressure |

Identify the Range of Flow Rates in Your System

Because Rotary Shear Valves have been used most often in chromatography systems, certain flow rate ranges have evolved functionally. However, these ranges can apply to any system, not just chromatography:

- Micro/Nano Scale flow rates less than 100 µL per minute
- Analytical Scale flow rates from $100 \, \mu L$ to $10 \, mL$ per minute
- > Prep (or Semi-Prep) Scale flow rates greater than 10 mL per minute

Decide What You Want the Valve to Do

In this chapter Rotary Shear Valves perform three functions:

- > Switching one or more flow paths to a different destination under pressure
- > Injection into a flowing stream under pressure
- > Selection/distribution of a variety of system liquids by means of a common port
- > Read more about valve functions on page 114.

Identify Whether You Want Automated or Manual Control

An automated valve offers more sophisticated functionality. Choose an automated valve if the application requires fast, consistent flow-stream switching. Some other advantages of automated valves include control options (PC- or instrument-triggered), higher torque operation, valve-position feedback, or very small flow paths.

Choose a manual valve if your application involves low frequency of use, demands operator control, or involves injection of smaller sample volumes. (See page 122 for more on Single Mode vs. Dual Mode operation.)

Identify the Chemical Compatibility Requirements Related to Your Fluids

Consulting the chemical compatibility chart in the Technical Resources section at the back of The IDEX Health & Science Laboratory Products catalog helps identify what valve materials to use and avoid — in your application. You can also find Chemical Compatibility information at www.idex-hs.com under Materials and Tools.

Identify Fluidic Connection Requirements in Your System

The rotary shear valves in this catalog accommodate one or more of the following tubing outer diameters: 1/8", 1/16", or 1/32".

Effects of Valves & Tubing on Resolution

The effect of tubing on analytical and microscale analyses can be significant. Since dispersion caused by tubing is proportional to the fourth power of diameter, large bore tubing should be avoided when performing analytical scale or microscale analyses. Tubing ID size ≤ 0.25 mm (0.010") is recommended.

Consider a system with injection and column switching valves and analytical columns with small-bore connecting tubing. The chromatograms below, made using a typical analytical chromatograph, show these effects. Scheme A is the control (injection valve → column → detector) with no valve in the system. In Schemes B and C, two model 7060 Six-Position Switching Valves were placed side by side (injection valve → valve #1 → column -> valve #? -> detector)


Australian Distributors HROMalytic +61(0)3 9762 2034 Importers & Manufacurers www.chromtech.net.au ECH no logy Pty Ltd Website NEW: www.chromalytic.net.au E-mail: info@chromtech.net.au Tel: 03 9762 2034... in AUSTRALIA The injection valve and detector were connected to these valves by the same tubing used in the control. The extra tubing pieces required to connect the valves to the column were a 10 cm length for valve #1-to-column, and a 35 cm length for column-to-valve #2. The diameters of these tubes are indicated in the experimental details, to the right.

Comparison of Observed Column Plates of Analytical and MicroScale Injection Valves

| | 7725 | 8125 | Δ | |
|----------|------|------|------|--|
| k' = 0.6 | 2930 | 5054 | 72% | |
| k' = 1.5 | 4653 | 6904 | 48% | |
| k' = 7.9 | 7875 | 8305 | 5.0% | |

UV detector: 1 μ L volume, 4 mm path. Sample volume: 2 μ L, partial-filling method. Column: 2 mm ID x 100 mm long, 4 μ m C-18. True plates of column = 11,570.

Effects of Valves and Tubing on Resolution



Conclusion: These sequential chromatograms show the effect of adding volume to the flow path through the addition of components.

- (A) Establishes a baseline quality of separation with the minimum volume of liquid in the flow path.


 (B) Adding a valve plus smaller-ID tubing, and thereby increasing the liquid volume only marginally, barely
- affects the separation. However in


 (C) Adding a valve plus larger-ID tubing, thereby increasing the liquid volume in the flow path to a greater degree, distinctly impairs the quality of the separation and the detectable sample.

WHAT IS MAKE-BEFORE-BREAK™, AND WHEN DOES IT MATTER?


Make-Before-Break is a unique design feature of certain dual-mode manual injection valves.

Flow paths of model 7725(i) and 9725(i) with MBB design






Position B




Position C


To maintain a constant, desired high-pressure flow, our Make-Before-Break (MBB®) design creates continuous flow between the LOAD and INJECT positions that virtually eliminates pressure transient shock to the system. A passage in the stator face makes a new connection before old connections break. The MBB design — an improvement over bypass-style injectors — does not dilute the sample and is easy to maintain and troubleshoot.



Switching Valves

Switching valves dynamically alternate between two fluid paths without manually disconnecting plumbing. In Chromatography, these valves can be used for column switching, backflushing, sample enrichment, and other techniques. In Diagnostic or Sequencing applications, the switching valve may alternate flow paths to enable back flushing or other fluidic tasks within the instrument.






Flow path of Two-Position, Ten-Port Switching Valve

Our switching valves operate between two positions, and may have 6 or 10 ports on the face of the stator (2/6 or 2/10). The flow paths connect ports around the circumference of the stator. The manual switching valves (to 1,000 psi) described on page 123 have different flow path geometry as noted.

Injection Valves

Our injection valves are a form of switching valve. Injection valves can be automated or manual, and they are generally utilized in the two-position, six-port (2/6) configuration and have a sample loop attached.



Flow path of Two-Position, Six-Port Injection Valve

The purpose of an Injection valve is to introduce a sample into a flowing stream of liquid. Some Switching valves become Injection valves by the addition of a Sample Loop (a defined length of tubing and fittings configured to match the angle of the valve ports). Sample is loaded and held in the loop until injection is triggered, either manually or automatically.

Injection valves are classified as either Single or Dual Mode based on how the Sample Loop can be filled. A Single Mode Injection valve requires complete filling of the sample loop and is configured for Rear loading, generally in an auto-sample configuration. A Dual Mode Injection valve allows either partial or complete filling of the loop, and introduces sample by syringe through the needle port built into the valve shaft. Complete filling of the sample loop in both the Dual and Single Mode Injection valves provides greater repeatability injection to injection. (See the Application Note, page 131 for greater detail on partial vs. complete loop filling.)

Valve Overview and Functions (Cont.)

Selection Valves

Selection valves enable discrete connections among multiple system liquids (mobile phase, reagents, buffers) by means of a common port (inlet or outlet) connected to a number of different reciprocal ports. In Diagnostic or Sequencing applications, the selection valve alternates between different reagents or sample streams.


Numerous configurations exist among selection valves (e.g., 6-position 7-port, or 10-position 11-port), but these valves typically operate between more than two positions. The ports are usually spaced radially, or outward in some manner around the center port of the stator.



Flow path of Six-Position, Seven-Port Selector Valve



Six column selection using two selection valves.







Stand Alone Valve Products

An automated valve offers more sophisticated functionality. Choose an automated valve if the application requires fast, consistent flow-stream switching. Some other advantages of automated valves include control options (PC- or instrument-triggered), higher torque operation, valve-position feedback, or very small flow paths.

MX Series II

MXT to 15,000 psi (1,035 bar)





Add our MX Series II™ actuated valves to your existing instrument or use in stand-alone lab configurations. MX valves can be controlled remotely or operated manually using the push-button front panel with LED position indicator. MX valves connect to your instrument or PC through contact closure, BCD, serial port, or USB. Commands can be sent to the MX valves using your chromatography software or the included proprietary software for timed-events programmability.


Available flow rates include options for Analytical, Micro/Nano, or Semi-Prep in a range of pressure capabilities. Valve liquid ends are available in materials chosen to be chemically inert and biocompatible. Routine maintenance using authorized RheBuild® kits (page 124) or — for the higher-pressure MXP and MXP valves — the Rapid Replacement Pods™ (page 120) assures optimal performance.

| Part No. | Description | Ports, Connections | Wetted Material | Rapid Replacement Pod | Qty. |
|-------------------|--|--|-----------------|-----------------------|------|
| ACTUATED V | ALVES UP TO 15,000 PSI (1,035 BAR) | | | | |
| SWITCHING | | | | | |
| MXT715-000 | 2-Position, 6-Port | 10-32 Ports for 1/16" OD Tubing | UltraLife | PD715-000 | ea. |
| MXT715-102 | 2-Position, 10-Port | 10-32 Ports for 1/16" OD Tubing | UltraLife | PD715-102 | ea. |
| INJECTION | For Injection, add the appropriately sized Sar | mple Loop to the Switching valve | s above | | |
| SELECTION | | | | | |
| MXT715-105 | 6-Position, 7-Port | 10-32 Ports for 1/16" OD Tubing | UltraLife | PD715-105 | ea. |
| All of these MX | Valves include a set of 1/16" fittings. Replacement Fitting | gs for MXT valves can be located on page | e 133. | | |
| ACTUATED V | ALVES UP TO 6,000 PSI (410 BAR) | | | | |
| SWITCHING | | | | | |
| MXP7900-000 | 2-Position, 6-Port | 10-32 Ports for 1/16" OD Tubing | DuraLife®* | PD7900 | ea. |
| MXP7960-000 | 2-Position, 10-Port | 10-32 Ports for 1/16" OD Tubing | DuraLife | PD7960 | ea. |
| MXP7980-000 | 2-Position, 6-Port, Nano, 5,000 psi (345 bar) | M4 Ports for 1/32" OD Tubing | DuraLife II | PD7980 | ea. |
| MXP7986-000 | 2-Position, 10-Port, Nano, 5,000 psi (345 bar) | M4 Ports for 1/32" OD Tubing | DuraLife II | PD7986 | ea. |
| MXP9900-000 | 2-Position, 6-Port, Biocompatible, 5,000 psi (345 bar) | 10-32 Ports for 1/16" OD Tubing | PEEK | PD9900 | ea. |
| MXP9960-000 | 2-Position, 10-Port, Biocompatible, 5,000 psi (345 bar) | 10-32 Ports for 1/16" OD Tubing | PEEK | PD9960 | ea. |
| INJECTION | For Injection, add the appropriately sized Sa | mple Loop to the Switching valve | s above | | |
| MXP7920-000 | 2-Position, 6-Port, Vertical Port | 10-32 Ports for 1/16" OD Tubing | DuraLife | PD7920 | ea. |
| SELECTION | | | | | |
| MXP7970-000 | 6-Position, 7-Port | 10-32 Ports for 1/16" OD Tubing | DuraLife II** | PD7970 | ea. |
| ** DuraLife II is | proprietary material combination of SST and an advance a proprietary material combination consisting of Titaniur valves include a set of 1/16" fittings. Replacement Fitting | n and an advanced polymer. | e 133. | | |
| ACTUATED V | ALVES UP TO 125 PSI (8.5 BAR) | | | | |
| SWITCHING | | | | | |
| MXX777-601 | 2-Position, 6-Port | Accepts Either 1/16" or 1/8" Tubing | RPC-7* | 1/16" and 1/8" | ea. |
| MXX777-603 | 2-Position, Double Three Way | Accepts Either 1/16" or 1/8" Tubing | RPC-7 | 1/16" and 1/8" | ea. |
| MXX777-612 | 2-Position, 6-Port, Large Bore | Accepts Either 1/16" or 1/8" Tubing | RPC-7 | 1/16" and 1/8" | ea. |
| SELECTION | | | | | |
| MXX777-605 | 6-Position, 7-Port | Accepts Either 1/16" or 1/8" Tubing | RPC-7 | 1/16" and 1/8" | ea. |
| MXX777-616 | 6-Position, 7-Port, Large Bore | Accepts Either 1/16" or 1/8" Tubing | RPC-7 | 1/16" and 1/8" | ea. |
| MXX778-605 | 10-Position, 11-Port | Accepts Either 1/16" or 1/8" Tubing | RPC-7 | 1/16" and 1/8" | ea. |
| * RPC-7 Propries | tary Polymer Combination. | | | | |

All of these MXX valves include a set of 1/16" and 1/8" ferrules. Replacement Fittings for MXX valves can be located on page 133.

HROW
461(0)3 9762 2034

Website NEW: www.chromalytic.net.au E-mail: info@chromtech.net.au Tel: 03 9762 2034 . . . in AUSTRALIA







Rapid Replacement Pods

For IDEX Health & Science MX Series II Valves

- Zero downtime maintenance
- Improves lab throughput

To help keep your instrument online and performing at maximum precision, select the exact Rapid Replacement Pod for your higher pressure MX Series II valves. Replacement pods are easily exchanged as part of scheduled preventive maintenance, or in an emergency, a pod can be substituted quickly while the original is examined and maintained at your convenience. The pod kit contains complete instructions for removal and replacement.



| Part No. | Description | For Valve Part No. | Qty. |
|--------------------|---------------------------------|-------------------------|------|
| TO 15,000 PSI (1,0 | 35 BAR) | | |
| SWITCHING | | | |
| PD715-000 | Rapid Replacement Pod | MXT715-000 | ea. |
| PD715-102 | Rapid Replacement Pod | MXT715-102 | ea. |
| SELECTION | | | |
| PD715-105 | Rapid Replacement Pod | MXT715-105 | ea. |
| UP TO 6,000 PSI (4 | 410 BAR) | | |
| SWITCHING | | | |
| PD7900 | Rapid Replacement Pod | MXP7900-000 | ea. |
| PD7960 | Rapid Replacement Pod | MXP7960-000 | ea. |
| PD7980 | Rapid Replacement Pod | MXP7980-000 | ea. |
| PD7986 | Rapid Replacement Pod | MXP7986-000 | ea. |
| PD9900 | Rapid Replacement Pod | MXP9900-000 | ea. |
| PD9960 | Rapid Replacement Pod | MXP9960-000 | ea. |
| INJECTION | | | |
| PD7920 | Rapid Replacement Pod | MXP7920-000 | ea. |
| SELECTION | | | |
| PD7970 | Rapid Replacement Pod | MXP7970-000 | ea. |
| H | ROM = 1 y tic +61(0)3 9762 2034 | Australian Distributors | |

FLUIDICS



Manual Valves

Choose a manual valve if your application involves low frequency of use, demands operator control, or involves injection of smaller sample volumes.







7725i

Manual Injection Valve Up to 9,000 psi (600 bar)

7060

Manual Switching Valve Up to 7,000 psi (483 bar)

3725i-038

Manual Switching Valve Up to 7,000 psi (483 bar)




Manual Switching Valve Options



4-Way Diagonal Flow Switching Valves (V-100D, V-101D)




Right Angle Flow Switching Valves (V-100L, V-101L)



3-Way Flow Switching Valves (V-100T, V-101T)









Manual Valves (Cont.)




Switching Valve Applications

Protect sensitive system components (such as a column) during a cleaning cycle with our Diagonal Flow Switching Valve ("D"). This valve eliminates the need to remove, plug and reconnect a low pressure column (see below).






A typical application for a Right Angle Flow Switching Valve ("L") is column switching, allowing two columns to use one detector. Detector switching is another common application for this valve (see below). Plug off the extra port with the included plug.




Position 1



Position 2

Your detector switching application may require the flexibility of routing the column effluent to both detectors simultaneously while retaining the ability to isolate each detector. Use our 3-Way Flow Switching Valve ("T"), plugging off the fourth port with the included plug.



Characteristics of Manual Sample Injection Valves

| Scale | Partial Filling Volumes (Range) | Sample Loop Sizes (Range) | Wetted Materials | Max. psi (bar)¹ | Max. T (°C) | MBB ² | Model ³ |
|-------------|------------------------------------|--|---|---|---|--|---|
| Analytical | 1 μL–2.5 mL 1 μL–5.0 mL | 2 μL–5.0 mL 2 μL–10 mL | 316 SST, Vespel® PEEK, ETFE, ceramic | 7,000 (483) 5,000 (340) | 80° 50° | Yes Yes | 7725, 7725i 9725, 9725i |
| Micro | 0.1 μL–500 μL | 5 μL–1.0 mL | 316 SST, PEEK, Vespel, ceramic | 7,000 (483) | 80° | No | 8125 |
| Preparative | 100 μL–10 mL | 2.0 mL-20 mL | 316 SST, PEEK PEEK | 5,000 (340) 4,000 (276) | 50° 50° | Yes Yes | 3725(i)-038, 3725i |
| Analytical | Not Applicable | 5 μL–5.0 mL 5 μL–10 mL | 316 SST, Vespel PEEK, ETFE, Ceramic | 7,000 (483) 5,000 (340) | 150° 50° | No No | 7010 9010 |
| | Analytical Micro Preparative | Scale Volumes (Range) Analytical 1 μL-2.5 mL 1 μL-5.0 mL Micro 0.1 μL-500 μL Preparative 100 μL-10 mL | Scale Volumes (Range) Sizes (Range) Analytical 1 μL-2.5 mL 1 μL-5.0 mL 2 μL-5.0 mL 2 μL-10 mL Micro 0.1 μL-500 μL 5 μL-1.0 mL Preparative 100 μL-10 mL 2.0 mL-20 mL Analytical Not Applicable 5 μL-5.0 mL | Scale Volumes (Range) Size's (Range) Materials Analytical 1 μL−2.5 mL 1 μL−5.0 mL 2 μL−5.0 mL 2 μL−10 mL 316 SST, Vespel® PEEK, ETFE, ceramic Micro 0.1 μL−500 μL 5 μL−1.0 mL 316 SST, PEEK, Vespel, ceramic Preparative 100 μL−10 mL 2.0 mL−20 mL 316 SST, PEEK PEEK Applicable 5 μL−5.0 mL 316 SST, Vespel | Scale Volumes (Range) Sizes (Range) Materials psi (bar)¹ Analytical 1 μL-2.5 mL 1 μL-5.0 mL 2 μL-5.0 mL 2 μL-10 mL 316 SST, Vespel® 7,000 (483) 5,000 (340) Micro 0.1 μL-500 μL 5 μL-1.0 mL 316 SST, PEEK, Vespel, ceramic Vespel, ceramic 7,000 (483) Preparative 100 μL-10 mL 2.0 mL-20 mL 316 SST, PEEK PEEK 5,000 (340) 4,000 (276) Apalytical Not Applicable 5 μL-5.0 mL 316 SST, Vespel 7,000 (483) | Scale Volumes (Range) Sizes (Range) Materials psi (bar)¹ T (°C) Analytical 1 μL-2.5 mL 1 μL-5.0 mL 2 μL-5.0 mL 2 μL-10 mL 316 SST, Vespel® PEEK, ETFE, ceramic 7,000 (483) 5,000 (340) 80° 50° Micro 0.1 μL-500 μL 5 μL-1.0 mL 316 SST, PEEK, Vespel, ceramic 7,000 (483) 80° Preparative 100 μL-10 mL 2.0 mL-20 mL 316 SST, PEEK PEEK 5,000 (340) 4,000 (276) 50° 50° Analytical Not Applicable 5 μL-5.0 mL 316 SST, Vespel 7,000 (483) 150° | Scale Volumes (Rănge) Sizes (Range) Materials psi (bar)¹ T (°C) MBB² Analytical 1 μL-2.5 mL 1 μL-5.0 mL 2 μL-5.0 mL 2 μL-10 mL 316 SST, Vespel® PEEK, ETFE, ceramic 7,000 (483) 5,000 (340) Vespel, ceramic 80° 7,000 (340) Vespel, ceramic No Preparative 100 μL-10 mL 2.0 mL-20 mL 2.0 mL-20 mL 316 SST, PEEK PEEK 5,000 (340) 4,000 (276) 50° 50° Yes Yes Yes Apalytical Not Applicable 5 μL-5.0 mL 316 SST, Vespel 7,000 (483) 7,000 (483) 150° No No |

This is the maximum pressure to which the valve can be adjusted. Some models are shipped from the factory set for lower pressures.

MBB (Make-Before-Break**) is a design that provides uninterrupted flow when switching between LOAD and INJECT. MBB also greatly reduces transient pressure shocks.

Models with an "i" suffix have a built-in position sensing switch. Models 8125 and 9010 also have a built-in switch.





SPECIFICATIONS & DETAILS

| Part No. | Stator Passage Diameter | Factory Set Pressure | Maximum Field Set Pressure | Maximum Temperature (°C) |
|-----------------------|-------------------------|----------------------|----------------------------|--------------------------|
| 3000 (PEEK) | 1.0 mm (0.040") | 3,000 psi (207 bar) | 4,000 psi (276 bar) | 50° |
| 7000, 7010 (SST) | 0.6 mm (0.024") | 5,000 psi (340 bar) | 7,000 psi (483 bar) | 150° |
| 7000L (SST) | 1.0 mm (0.040") | 3,000 psi (207 bar) | 5,000 psi (340 bar) | 150° |
| 7030 (SST) | 0.6 mm (0.024") | 5,000 psi (340 bar) | 7,000 psi (483 bar) | 150° |
| 7030L (SST) | 1.0 mm (0.040") | 3,000 psi (207 bar) | 5,000 psi (340 bar) | 150° |
| 7060 (SST) | 0.4 mm (0.016") | 5,000 psi (340 bar) | 7,000 psi (483 bar) | 80° |
| SST = Stainless Steel | | | | |

Manual Valves

| Part No. | Description | Tubing/Fitting Size | Wetted Material | Configuration | Qty. |
|------------------------|--|--|----------------------------------|---------------|------|
| MANUAL VA | LVES UP TO 9,000 PSI (600 BAR) | | | | |
| INJECTION | | | | | |
| 7725i-188 | 2-Position, 6-Port, 9,000 psi (600 bar) | 10-32 Ports for 1/16" OD Tubing | Stainless Steel, PEEK, Ceramic | Front loading | ea. |
| MANUAL VA | LVES UP TO 6,000 PSI (410 BAR) | | | | |
| SWITCHING | | | | | |
| 3000 | 2-Position, 6-Port, Prep Scale | 5/16-24 Ports for 1/8" OD Tubing | PEEK | _ | ea. |
| 7000 | 2-Position, 6-Port, Large Bore | 10-32 Ports for 1/16" OD Tubing | Stainless Steel & Vespel® | _ | ea. |
| 7000L | 2-Position, 6-Port, Large Bore | 10-32 Ports for 1/16" OD Tubing | Stainless Steel & Vespel | _ | ea. |
| 7030 | 2-Position, 6-Port | 10-32 Ports for 1/16" OD Tubing | Stainless Steel & Vespel | Double 3-Way | ea. |
| 7030L | 2-Position, 6-Port, Large Bore | 10-32 Ports for 1/16" OD Tubing | Stainless Steel & Vespel | Double 3-Way | ea. |
| INJECTION ¹ | | | | | |
| 7010 | 2-Position, 6-Port Single Mode | 10-32 Ports for 1/16" OD Tubing | Stainless Steel & Vespel | 20 μL* | ea. |
| 9010 | 2-Position, 6-Port Single Mode (Switching, Injection) | 10-32 Ports for 1/16" OD Tubing | PEEK, ETFE, Ceramic | 20 μL* | ea. |
| 3725-038 | 2-Position, 6-Port, Prep Scale Dual Mode | 5/16-24 Ports for 1/8" Tubing | Stainless Steel & PEEK | 10 mL* | ea. |
| 3725i | 2-Position, 6-Port, Prep Scale Dual Mode with Switch | 5/16-24 Ports for 1/8" Tubing | PEEK | 10 mL* | ea. |
| 3725i-038 | 2-Position, 6-Port, Prep Scale Dual Mode with Switch | 5/16-24 Ports for 1/8" Tubing | Stainless Steel & PEEK | 10 mL* | ea. |
| 7725 | 2-Position, 6-Port, Analytical Scale Dual Mode | 10-32 Ports for 1/16" OD Tubing | Stainless Steel, Ceramic, Vespel | 20 μL* | ea. |
| 7725i | 2-Position, 6-Port, Analytical Scale Dual Mode with Switch | 10-32 Ports for 1/16" OD Tubing | Stainless Steel, Ceramic, Vespel | 20 μL* | ea. |
| 8125** | 2-Position, 6-Port, Micro Scale Dual Mode with Switch | 10-32 Ports for 0.020" (0.5 mm) or 1/16" Tubing | Stainless Steel, Ceramic, Vespel | 5 μL* | ea. |
| 9725 | 2-Position, 6-Port, Analytical Scale Dual Mode | 10-32 Ports for 1/16" OD Tubing | PEEK, ETFE, Ceramic | 20 μL* | ea. |
| 9725i | 2-Position, 6-Port, Analytical Scale Dual Mode with Switch | 10-32 Ports for 1/16" OD Tubing | PEEK, ETFE, Ceramic | 20 μL* | ea. |
| SELECTION | | | | | |
| 7060 | 6-Position, 7-Port | 10-32 Ports for 1/16" OD Tubing | Stainless Steel & Vespel | 6-Way | ea. |
| * Shine with : | sample loop of indicated volume attached to ports 1 and 4 | | | | |

Ships with a sample loop of indicated volume attached to ports 1 and 4.
 The 8125 requires special ferrules for 0.020" (0.5 mm) tubing. 8125-084-0.5 mm ferrule for 8125; 8125-086-0.5 mm ferrule for 8125 — 4-pk.

| | escription escription | Tubing/Fitting Size | Wetted Material | Configuration | Includes | Qty. |
|-------------------|--|-----------------------------------|-----------------|----------------------|------------|------|
| SWITCHING | | | | | | |
| V-100D 4-F | Position, 4-Port, 500 psi (34 bar) | 1/4-28 Ports for 1/16" OD Tubing | PEEK, PTFE | Double Diagonal | * | ea. |
| V-101D 4-F | Position, 4-Port, Bulkhead, 500 psi (34 bar) | 1/4-28 Ports for 1/16" OD Tubing | PEEK, PTFE | Double Diagonal | * | ea. |
| V-100L 4-F | Position, 4-Port, 500 psi (34 bar) | 1/4-28 Ports for 1/16" OD Tubing | PEEK, PTFE | Right-Angle "L" | ** | ea. |
| V-101L 4-F | Position, 4-Port, Bulkhead, 500 psi (34 bar) | 1/4-28 Ports for 1/16" OD Tubing | PEEK, PTFE | Right-Angle "L" | ** | ea. |
| V-100T 4-F | Position, 4-Port, 500 psi (34 bar) | 1/4-28 Ports for 1/16" OD Tubing | PEEK, PTFE | Single "T" | *** | ea. |
| V-101T 4-F | Position, 4-Port, Bulkhead, 500 psi (34 bar) | 1/4-28 Ports for 1/16" OD Tubing | PEEK, PTFE | Single "T" | *** | ea. |
| INJECTION Fo | or Injection, add the appropriately sized Samp | le Loop to the Switching valves a | above | | | |
| V-450 2-F | Position, 6-Port, 1,000 psi (69 bar) | 1/4-28 Ports for 1/16" OD Tubing | Polyimide, PTFE | Injection | (6) XP-235 | ea. |
| V-451 2-F | Position, 6-Port, Bulkhead Version, 1,000 psi (69 bar) | 1/4-28 Ports for 1/16" OD Tubing | Polyimide, PTFE | Injection | (6) XP-235 | ea. |
| V-540 2-F | Position, 6-Port, 1,000 psi (69 bar) | 1/4-28 Ports for 1/8" OD Tubing | Polyimide, PTFE | Injection | (6) XP-335 | ea. |
| V-541 2-F | Position, 6-Port, Bulkhead Version, 1,000 psi (69 bar) | 1/4-28 Ports for 1/8" OD Tubing | Polyimide, PTFE | Injection | (6) XP-335 | ea. |
| SELECTION | | | | | | |
| V-240 6-F | Position, 7-Port, 1,000 psi (69 bar) | 1/4-28 Ports for 1/16" OD Tubing | Polyimide, PTFE | Multi-port Selection | (6) XP-235 | ea. |
| V-241 6-F | Position, 7-Port, Bulkhead Version, 1,000 psi (69 bar) | 1/4-28 Ports for 1/16" OD Tubing | Polyimide, PTFE | Multi-port Selection | (6) XP-235 | ea. |
| V-340 6-F | Position, 7-Port, 1,000 psi (69 bar) | 1/4-28 Ports for 1/8" OD Tubing | Polyimide, PTFE | Multi-port Selection | (6) XP-335 | ea. |
| V-341 6-F | Position, 7-Port, Bulkhead Version, 1,000 psi (69 bar) | 1/4-28 Ports for 1/8" OD Tubing | Polyimide, PTFE | Multi-port Selection | (6) XP-335 | ea. |

* (4) P-218BLK, (4) P-240. ** (4) P-218BLK, (4) P-240, (1) P-309. *** (4) P-218BLK, (4) P-240, (1) P-309.



Spare Parts

We offer a full line of genuine spare parts to assist with your use of our valve products. We offer RheBuild® Kits designed for specific valve models. Rotor Seal and Stators are commonly replaceable parts.



Rotor Seals & Stators

The rotor seal is the polymeric disc that makes a high pressure seal against the stator or stator face seal. The seal wears with use and is one of the only parts that may need routine replacement.

Stators are available in 316 stainless steel, PEEK and proprietary materials. Typically, stators need replacement only if the ports or sealing surfaces become damaged. Avoid damage from use of improper injection needles by referring to the "Using Proper Syringe Needles" Application Note on page 130.

Please Note: Rotor seals for MX Series II™ Modules are available in RheBuild® Kits on this page. Stators for MX Series II Modules are available on this page. MX (Series I) Module rotor seals are available in RheBuild Kits on this page.



RheBuild® Kits

RheBuild Kits are available for most valve products. Included in each individualized RheBuild Kit are all parts, tools, and instructions to maintain precision performance of your particular product. RheBuild Kits eliminate individual part ordering.



APPLICATION NOTE

How to Select the Right Rotor Seal

The standard rotor seal in many of our manual valves is made from a Vespel® blend. This polyimide has low wear and high chemical resistance. Vespel tolerates a pH range of 0 to 10. Solutions more basic than pH 10 attack Vespel which damages the rotor seal. If you use any solutions above pH 10, our experts recommend a PEEK blend rotor seal. PEEK offers a high chemical resistance and versatility, and will tolerate the entire pH range from 0 to 14. ETFE blend rotor seals are appropriate for use in applications where PEEK is not generally acceptable, such as when methylene chloride or DMSO in higher concentrations is being used.





APPLICATION NOTE


How to Avoid Pressure Transients

Air in the sample loop can cause an instantaneous system pressure drop that eventually returns to a normal level. Air causes the pressure to drop when the injector moves from the LOAD to the INJECT position. When large sample loops (≥ 100 µL) are partially loaded, air present in the needle port tube is pushed into the sample loop (see Figure 1). Air can also enter the sample loop from siphoning which occurs when the vent line is higher than the injection port. In either case, upon injection, the system pressure collapses the air bubble, causing pressure to drop momentarily.

A pressure drop in the system caused by air results in changes in retention time, artifact peaks, and affects column performance.

Avoid pressure drops by removing the air in the needle port tube. Do this by flushing about 1 mL of mobile phase with a luer syringe with needle port cleaner. Keep the needle port tube filled with mobile phase by occasional flushing. Adjust the vent line(s) so the outlet is at the same horizontal level as the needle port (see Figure 2).

For additional injection troubleshooting, refer to our Troubleshooting Guide for HPLC Injection Problems. You may download the Guide from the IDEX Health & Science web site: www.idex-hs.com under Education & Tools.






Figure 1 Air present in the needle port tube is pushed by the syringe during loading into the




Figure 2 Pathway of the flushing mobile phase using the Needle Port Cleaner, Part # 7125-054 (see page 131) when the injector is in INJECT

| Part No. | For Valve Model No. | Description | Qty. |
|----------------------|--|---------------------------------|------------|
| VESPEL BLEND | | - Countries | , |
| 7000-016 | 7000L, 7040L | Vespel Rotor Seal | ea. |
| 7010-039 | 7010, 7000, 7040 | Vespel Rotor Seal | ea. |
| 7030-003 | 7030, 9030 | Vespel Rotor Seal | ea. |
| 7030-014 | 7030L | Vespel Rotor Seal | ea. |
| 7060-070 | 7060, 7066 | Vespel Rotor Seal | ea. |
| 7060-064 | 7060L | Vespel Rotor Seal | ea. |
| 7125-047 | 7125, 7725, 9725 | Vespel Rotor Seal | ea. |
| 7410-038 | 7410 | Vespel Rotor Seal | ea. |
| 7413-013 | 7413 | Vespel Rotor Seal | ea. |
| 8125-038 | 8125 | Vespel Rotor Seal | ea. |
| ETFE BLEND RC | TOR SEALS | | |
| 7010-071 | 7010, 7010-087, 7000, 7040 | ETFE Rotor Seal | ea. |
| 7030-015 | 7030, 9030 | ETFE Rotor Seal | ea. |
| 7060-074 | 7060, 7066, 9060 | ETFE Rotor Seal | ea. |
| 7125-079 | 7125, 7125-081, 7725 | ETFE Rotor Seal | ea. |
| 8125-097 | 8125 | ETFE Rotor Seal | ea. |
| 9010-051 | 9010 | ETFE Rotor Seal | ea. |
| 9125-082 | 9125, 9725 | ETFE Rotor Seal | ea. |
| PEEK BLEND RO | | DEEK D. C. I | |
| 3725-018 | 3725, 3725-038 | PEEK Rotor Seal PEEK Rotor Seal | ea. |
| 9010-065 8125-119 | 7000, 7010, 9010 | PEEK Rotor Seal | ea. |
| 9125-095 | 8125 7125, 7725, 9125, 9725 | PEEK Rotor Seal | ea. ea. |
| | MX SERIES II MODULES | TEER ROLOI Seal | ea. |
| 7123-548 | MXT715-000 | Stator | ea. |
| 7123-550 | MXT715-105 | Stator | ea. |
| 7123-568 | MXT715-102 | Stator | ea. |
| 7770-229 | MXP7920-000 | Stator | ea. |
| 7980-004 | MXP7980-000 | Stator | ea. |
| 7986-004 | MXP7986-000 | Stator | ea. |
| 7900-146 | MXP9900-000 | Stator | ea. |
| 7900-179 | MXP7900-000 | Stator | ea. |
| 7900-183 | MXP7970-000 | Stator | ea. |
| 7960-014 | MXP7960-000 | Stator | ea. |
| 9960-002 | MXP9960-000 | Stator | ea. |
| STATORS FOR C | OTHER IDEX HEALTH & S | CIENCE VALVES | |
| 3725-006 | 3725, 3710-038, 3000-038 and 3030-038 | Stator | ea. |
| 7010-069 | 7000L, 7030L, 7040L | Stator | ea. |
| | 7010, 7125, 7000, 7030 | | |
| 7010-040 | and 7040 | Stator | ea. |
| 7060-039 | 7060 and 7066 | Stator | ea. |
| 7060-065 | 7060L, EV501-100 | Stator | ea. |
| 7123-047 | PR/EV500-100 | Stator | ea. |
| 7123-127 | PR/EV750-107 | Stator | ea. |
| 7123-128 | PR/EV700-107 | Stator | ea. |
| 7123-142 | PR/EV500-104, EV501-104 | Stator | ea. |
| 7123-145 7123-147 | PR/EV550-104, EV551-104 PR/EV550-100 | Stator | ea. |
| 7123-147 | PR/EV500-101 | Stator Stator | ea. |
| 7123-149 | PR/EV550-101 | Stator | ea. |
| 7123-180 | PR703-100 and EV700-105 | Stator | ea. |
| 7123-221 | PR753-100 and EV750-105 | Stator | ea. |
| 7123-223 | PR/EV700-112 | Stator | ea. |
| 7410-041 | 7410 and 7413 | Stator | ea. |
| 7520-030 (inlet) | 7520 | Stator | ea. |
| 7520-035 (outlet) | 7520 | Stator | ea. |
| 7650-002 | PR/EV700-102 | Stator | ea. |
| 7725-010 | 7725(i) | Stator | ea. |
| 7750-070 | 7750 | Stator | ea. |
| 7750-038 | PR/EV700-100 | Stator | ea. |
| 8125-098 | 8125 | Stator | ea. |
| 9125-043 | 9125, 9010, 9030 and 9725(i) | Stator | ea. |
| 9650-009 | PR/EV750-102 | Stator | ea. |

9750-021

PR/FV750-100

Stator HROM = 1 y tic +61(0)3 9762 2034

| Part No. | Description | Qty. |
|----------------|--|------|
| RHEBUILD H | CITS FOR MX SERIES II™ VALVES | - |
| 7150-999 | RheBuild Kit for MXT715-000 (includes 2 rotor seals) | ea. |
| 7152-999 | RheBuild Kit for MXT715-102 (includes 2 rotor seals) | ea. |
| 7155-999 | RheBuild Kit for MXT715-105 (includes 2 rotor seals) | ea. |
| 7920-999 | RheBuild Kit for MXP7920-000 and MXP7900-000 | ea. |
| 7960-999 | RheBuild Kit for MXP9960-000 (includes rotor seal and stator face seal) | ea. |
| 7961-999 | RheBuild Kit for MXP7960-000 | ea. |
| 7970-999 | RheBuild Kit for MXP7970-000 | ea. |
| 79801-999 | RheBuild Kit for MXP7980-000 | ea. |
| 79861-999 | RheBuild Kit for MXP7986-000 | ea. |
| 7900-999 | RheBuild Kit for MXP9900-000 (includes rotor seal and stator face seal) | ea. |
| RHEBUILD I | KITS FOR MANUAL VALVES | |
| 3725-999 | RheBuild Kit for models 3725, 3725i, 3725-038, 3735i-038 | ea. |
| 7010-997 | RheBuild Kit including Stator for model 7010 | ea. |
| 7010-998 | RheBuild Kit, pH Upgrade Kit for model 7000 | ea. |
| 7010-999 | RheBuild Kit for model 7010 and 7010-type Valves | ea. |
| 7125-999 | RheBuild Kit for models 7125 and 7126 | ea. |
| 7410-999 | RheBuild Kit for model 7410 | ea. |
| 7520-999 | RheBuild Kit for models 7520 and 7526 (includes inlet stator and seal) | ea. |
| 7725-999 | RheBuild Kit for models 7725 and 7725i | ea. |
| 7788-999 | RheBuild Kit for model 7725i-188 | ea. |
| 8125-999 | RheBuild Kit for models 8125 and 8126 | ea. |
| 9010-999 | RheBuild Kit for model 9010 | ea. |
| 9125-999 | RheBuild Kit for models 9125 and 9126 | ea. |
| 9725-999 | RheBuild Kit for models 9725 and 9725i; 7725(i) pH upgrade kit | ea. |
| RHEBUILD I | CITS FOR MX SERIES I™ VALVES | |
| 7900-999 | RheBuild Kit for models MX7900-000, MX7925-000, MX9900-000, MX9925-000 | ea. |
| 7960-999 | RheBuild Kit for model MX7960-000 | ea. |
| 7980-999 | RheBuild Kit for model MX7980-000 | ea. |
| 7984-999 | RheBuild Kit for model MX7984-000 | ea. |
| 7986-999 | RheBuild Kit for model MX7986-000 | ea. |
| | KITS FOR LABPRO™ & EV AUTOMATED STRUMENTS | |
| 1006-999 | RheBuild Kit for model PR/EV100-106 | ea. |
| 5001-999 | RheBuild Kit for models PR/EV500-101 and PR/EV550-101 | ea. |
| 5100-999 | RheBuild Kit for models PR/EV500-100 and PR/EV550-100 | ea. |
| 5104-999 | RheBuild Kit for models PR/EV500-104 and PR/EV550-104 | ea. |
| 7004-999 | RheBuild Kit for models PR/EV700-104 and PR/EV750-104 | ea. |
| 7112-999 | RheBuild Kit for models PR/EV700-112 and PR/EV750-112 | ea. |
| 7501-999 | RheBuild Kit for models PR/EV700-100 and PR/EV750-100 | ea. |
| 7502-999 | RheBuild Kit for models PR/EV700-102 and PR/EV750-102 | ea. |
| 7507-999 | RheBuild Kit for models PR/EV700-107 and PR/EV750-107 | ea. |
| 7531-999 | RheBuild Kit for models PR703-100 and PR753-100 | ea. |
| n Distributors | | |





VALVE ACCESSORIES

Our valve accessories include a variety of products that work with and are specific to our valve mechanics. From Sample Loops, driver boards, or mounting brackets we offer a wide array of accessories to meet your system requirements. We also include tools that work specifically with our valves and valve components.

- 127 STAINLESS STEEL SAMPLE LOOPS
- **128** PEEK SAMPLE LOOPS
- 131 SUCTION NEEDLE ADAPTER
- **132** INJECTION PORT ADAPTERS
- 133 WRENCHES, BRACKETS, & REPLACEMENT FITTINGS



Stainless Steel Sample Loops

APPLICATION NOTE

How to Properly Install Sample Loops: Stainless Steel

Stainless steel sample loops are supplied with fittings that are not swaged onto the tube. It is important that the loop be completely bottomed in the injector port before the ferrule is swaged onto the tube. The depth of the tubing holes may vary slightly from port to port and from valve to valve. A fitting made up in one port may leave a small cavity in another port. The cavity causes high dispersion and peak distortion such as fronting, tailing, or broadening. It is good practice to label loop ends so they will be replaced in the same, respective ports that were used in swaging the ferrules. Hint: swaging ferrules separately on each side, into each respective valve port makes loop installation easier.

These high quality stainless steel sample loops have burr-free, square-cut ends to ensure a flush connection to valve ports. The size designations of loops are nominal. The actual volumes can differ from the theoretical designations because the ID tolerance varies depending on the tubing tolerance of the metal tubing bore. Accuracy of large metal loops (1.0 mm, 0.040" bore) is about $\pm 5\%$, intermediate loops (0.5 mm, 0.020" bore) $\pm 10\%$, and small loops (0.2 mm, 0.007" bore) $\pm 30\%$.

Since both standards and unknowns are usually analyzed using the same sample loop, knowledge of the actual, accurate volume is rarely needed. If the sample loop volume must be known, it is best to calibrate the loop in place on the valve so the flow passages in the valve are also taken into account. An alternative to calibration is to use a dual mode injector and partial-filling method of loading. See the "Sample Loop Loading" Application Note on page 131.

Model 7725 Injector loops are not interchangeable with loops for the model 7125. The port angle for the 7725 is 30° whereas the port angle for the 7125 is 20° requiring the loops to have a different shape. Model 8125 Micro-Scale Sample Injector requires special loops in the $5.0~\mu L$ to $50~\mu L$ range. The 8125 sample loops are made with 0.5 mm (0.020″) OD tubing.

To install the sample loop:







Figure 1 Cut-away view of stainless steel sample loop installation

- Take one end of the loop and place the nut (1) and ferrule (2) onto the tubing (3) with the threaded portion of the nut and tapered portion of the ferrule toward the end. See Figure A. 2 Insert the tubing into port (4). Confirm that the tubing is bottomed in the valve port as shown in Figure A. While firmly pressing down on the tubing, hand-tighten the nut as tight as possible. With the IDEX Wrench (page 51), designed especially for fittings, tighten one quarter turn past finger tight. Remove the loop to confirm the ferrule is swaged onto the tube. Repeat steps a-d with the other end of the loop while the swaged end remains outside
- 6 Reinstall each end of the loop to their respective ports. See Figure C.

the valve port. See Figure B.

| Part No. | Volume | Tubing | Qty. |
|--------------|--|-------------------------------------|--------------|
| STAINLESS ST | TEEL LOOPS FOR 7125, 7010 INJE | CTION VALVES (DO NOT USE FOR 7725) | |
| 7020 | 5 μL Sample Loop | 0.18 mm (0.007") ID x 1/16" OD | ea. |
| 7021 | 10 μL Sample Loop | 0.30 mm (0.012") ID x 1/16" OD | ea. |
| 7022 | 20 μL Sample Loop | 0.51 mm (0.020") ID x 1/16" OD | ea. |
| 7023 | 50 μL Sample Loop | 0.51 mm (0.020") ID x 1/16" OD | ea. |
| 7024 | 100 μL Sample Loop | 0.51 mm (0.020") ID x 1/16" OD | ea. |
| 7025 | 200 μL Sample Loop | 0.76 mm (0.030") ID x 1/16" OD | ea. |
| 7026 | 500 μL Sample Loop | 0.76 mm (0.030") ID x 1/16" OD | ea. |
| 7027 | 1.0 mL Sample Loop | 0.76 mm (0.030") ID x 1/16" OD | ea. |
| 7028 | 2.0 mL Sample Loop | 1.0 mm (0.040") ID x 1/16" OD | ea. |
| 7029 | 5.0 mL Sample Loop | 1.0 mm (0.040") ID x 1/16" OD | ea. |
| 1876 | 10 mL Sample Loop | 2.0 mm (0.080") ID x 1/8" OD | ea. |
| STAINLESS ST | TEEL LOOPS FOR 3725-038, 37251 | -038 INJECTION VALVES | |
| 3065-018 | 2.0 mL Sample Loop | 2.0 mm (0.080") ID x 1/8" OD | ea. |
| 3065-019 | 5.0 mL Sample Loop | 2.0 mm (0.080") ID x 1/8" OD | ea. |
| 3065-023 | 10 mL Sample Loop | 2.0 mm (0.080") ID x 1/8" OD | ea. |
| 3065-025 | 20 mL Sample Loop | 2.0 mm (0.080") ID x 1/8" OD | ea. |
| | FEEL LOOPS FOR 7725, 7725I, PR/ ALVES (DO NOT USE FOR 7125) | EV700-100, PR/EV703-100, MX MODULE | |
| 7755-020 | 5 μL Sample Loop | 0.18 mm (0.007") ID x 1/16" OD | ea. |
| 7755-021 | 10 μL Sample Loop | 0.30 mm (0.012") ID x 1/16" OD | ea. |
| 7755-022 | 20 μL Sample Loop | 0.30 mm (0.012") ID x 1/16" OD | ea. |
| 7755-023 | 50 μL Sample Loop | 0.51 mm (0.020") ID x 1/16" OD | ea. |
| 7755-024 | 100 μL Sample Loop | 0.51 mm (0.020") ID x 1/16" OD | ea. |
| 7755-025 | 200 μL Sample Loop | 0.76 mm (0.030") ID x 1/16" OD | ea. |
| 7755-026 | 500 μL Sample Loop | 0.76 mm (0.030") ID x 1/16" OD | ea. |
| 7755-027 | 1.0 mL Sample Loop | 0.76 mm (0.030") ID x 1/16" OD | ea. |
| 7755-028 | 2.0 mL Sample Loop | 1.0 mm (0.040") ID x 1/16" OD | ea. |
| 7755-029 | 5.0 mL Sample Loop | 1.0 mm (0.040") ID x 1/16" OD | ea. |
| 1876 | 10 mL Sample Loop | 2.0 mm (0.080") ID x 1/8" OD | ea. |
| STAINLESS ST | FEEL LOOPS FOR 8125 INJECTOR (| (USE 7755-024 TO 7755-029 FOR VOLUM | IES > 50 μL) |
| 8020 | 5 μL Sample Loop | 0.20 mm (0.008") ID x 0.020" OD | ea. |
| 8021 | 10 μL Sample Loop | 0.20 mm (0.008") ID x 0.020" OD | ea. |
| 8022 | 20 μL Sample Loop | 0.25 mm (0.010") ID x 0.020" OD | ea. |
| 8023 | 50 ul. Sample Loop | 0.30 mm (0.012") ID x 0.020" OD | |
| 2 0762 2024 | Australian Distributors | | |



PEEK Sample Loops

Flexible PEEK sample loops are alternatives to stainless steel loops. PEEK loop ends are provided with clean, straight cuts for easy valve installation.

PEEK polymer is inert to almost all organic solvents and is biocompatible, giving PEEK loops added versatility. Natural PEEK is used for these sample loops. Like metal loops, the size designations of PEEK loops are nominal. The actual volumes can differ from the theoretical designations because of the tolerance of the tubing bore. Accuracy of large PEEK loops (0.8 mm, 0.030" bore) is about $\pm 14\%$, intermediate loops (0.5 mm, 0.020") $\pm 21\%$, and small loops (0.2 mm, 0.007") $\pm 65\%$.

PEEK loops are also supplied with unswaged RheFlex® fittings but do not require the same swaging precaution. The fittings can reposition along the loop tubing when the fitting is reinserted in the ports for correct loop installation.

Sample Mobile Phase Tube Wall

Figure 1 Schematic of sample flow through mobile phase between tubing walls




Figure 2 Sample mass (observed peak area) vs. volume of sample dispensed from the syringe, in units of loop volumes, injected onto the column from our dual mode injector such

APPLICATION NOTE

Fluidic Movement in Tubes

Q: "Why can I load only up to half of the volume of the loop in partial-filling method?"

A: Sample occupies $2~\mu L$ of loop for every $1~\mu L$ loaded from the syringe. For example, $10~\mu L$ of sample spreads out over the entire length of a $20~\mu L$ loop. Any additional sample loaded will overflow the end of the loop and exit out to waste. Reproducibility is poor because the volume of sample in the loop is different from the known volume originally loaded by your syringe.

Fluid spreads in a parabolic shape through a tube instead of moving in one plug because the velocity is different at the center of the tube than at the walls. The velocity at the center of the tube is twice the average velocity, and near the wall the velocity is almost zero, creating a parabolic shape. This fluidic movement is called laminar flow. See Figure 1.

In dual mode injection valves (see "Sample Loop Loading" Application Note on page 131) the sample from the syringe needle loads directly into the sample loop. The sample volume is known since there is no sample waste. The laminar flow phenomenon accounts for the shape of the plot as shown in Figure 2. Note that the plot has three regions:

- 1 Partial-Filling Region. When the volume dispensed is less than half the loop volume, the curve is linear. Sample has not reached the end of the loop. Within this region, performance depends on the syringe and operator.
- 2 Nonlinear Region. When the volume dispensed is between half the loop volume and about two loop volumes, the curve is nonlinear. Sample is lost from the loop, so reproducibility is poor. If you dispense a volume equal to the loop size, you are in this region of poor performance.
- 3 Complete-Filling Region. When the volume of sample dispensed is several loop volumes, the loop contains only pure sample, undiluted by residual mobile phase. Within this region, reproducibility is highest.

In the single mode injection valves the sample must pass through a connecting passage before it reaches the sample loop. Since some of the sample dispensed from the syringe remains in the connecting passageway, an unknown amount enters the sample loop. Therefore, single mode injection valves achieve high reproducibility only by using the complete-filling method.

129



PEEK Physical Strength Characteristics

Although PEEK material is compatible with virtually all solvents, there are many factors that affect burst pressure of PEEK tubing. Factors such as increases in inner diameter, temperature, exposure time, and concentration of organic solvents affect the degradation of PEEK. Other solvents such a THF, methylene chloride and DMSO cause PEEK tubing to swell while concentrated nitric acid and sulfuric acid weaken the tubing.



How to Find and Fix Common Sample Injector Leaks

Leaks cause valuable sample loss. Nobody wants that. The key to the valve holding pressure is the integrity of the sealing surfaces. If there is a scratch on the sealing surface, or the needle seal in the rotor seal is damaged, a leak may appear. It is also important to realize what appears to be a leak can instead be a result of siphoning. The following are the three most common situations in which fluid leaks occur.

- 1 If fluid leaks out of the needle port only while loading the loop (i.e., while pushing down on the plunger of the syringe), the problem is most likely that the needle seal or the needle port fitting in the loop filler port is not gripping the syringe needle tightly enough. Tighten the needle seal grip by pushing with the eraser end of a pencil on the needle port (See Figure 1). The tightening reduces the hole diameter of the needle seal and port fitting.
- 2 If fluid leaks continuously from the needle port or vent lines and/or from the stator-tostator ring interface, replace the rotor seal and/or stator face assembly. Scratches on the rotor seal or cracks in the stator face assembly allow mobile phase to escape and cause cross port leakage. Genuine IDEX Health & Science RheBuild® Kits are listed on page 124.
- 3 If fluid leaks from the needle port and/or vent lines but eventually stops, the cause is most likely siphoning and not a leak. Siphoning occurs if the vent lines are lower or higher than the needle port. Adjust the vent line(s) so that the outlet is at the same horizontal level as the needle port to prevent siphoning. (See Figure 2).

For other leakage or injection troubleshooting, refer to our Troubleshooting Guide for HPLC Injection Problems. You may download the Guide from our web site: www.idex-hs.com under Education & Tools.




Figure 1 To reform the needle seal, push the eraser end of a pencil against the needle port







Figure 2 Needle port level compared to the level of vent line outlet:

- (A) siphoning occurs when the vent line outlet is above the needle port level
- (B) siphoning does not occur if the vent line outlet is the same horizontal level as the needle port

130

PEEK Sample Loops (Cont.)






Figure 1 A square cut needle: (A) stops against the stator face assembly; The tip of a pointed needle (B) slips into the stator face and the tip breaks off as the valve rotates



Using Proper Syringe Needles


With front-loading injection valves it is important to use the correct needle when loading the sample loop. An incorrect needle will damage the valve and can cause poor reproducibility. When the needle is too short the tip will not reach the needle seal. When the needle is too small in diameter the seal will not grip tightly enough. Needles with a beveled tip can damage the rotor seal and stator face assembly (see Figure 1). The needle should be #22 gauge (0.028"-0.0285" / 0.72 mm), and 90° point style (square cut end). Model 3725i requires a #16 gauge (0.0645"-0.0655"/ 1.65 mm) needle. Never use a beveled, pointed, or tapered needle.

Needle specifications are not critical when using a Loop Filler Port to load the sample loop. However, it is important to tighten the needle port fitting around the needle if using a syringe needle with a slightly smaller diameter than 0.7 mm (0.028").

If the loading method used is complete-filling, a syringe without a needle can be used. A syringe fitted with a Needle Port Cleaner can be used with a front-loading valve (Figure 2A) or with a Loop Filler Port (Figure 2B).

Needle port accessories are listed on page 132.





(A) Syringe fitted with Needle Port Cleaner (Part # 7125-054) loading a front-loading valve (model 7725); (B) loading a Loop Filler Port (Part # 7012)

PEEK Sample Loops

| Part No. | Volume | Tubing | Valco No. | Qty. |
|-----------|--|--------------------------------|-----------|------|
| PEEK LOOF | PS FOR 3725, 3725I INJECTION VALVES | | | |
| 3055-018 | 2.0 mL Sample Loop | 1.6 mm (0.062") ID x 1/8" OD | N/A | ea. |
| 3055-019 | 5.0 mL Sample Loop | 1.6 mm (0.062") ID x 1/8" OD | N/A | ea. |
| 3055-023 | 10 mL Sample Loop | 2.0 mm (0.080") ID x 1/8" OD | N/A | ea. |
| 3055-025 | 20 mL Sample Loop | 2.0 mm (0.080") ID x 1/8" OD | N/A | ea. |
| PEEK LOOP | S FOR 9725, 9010, PR/EV750-100, PR/EV753-100 INJECTION VALVI | ES | | |
| Part No. | Volume | Bore / Tubing | Valco No. | |
| 9055-020 | 5.0 µL Sample Loop | 0.18 mm (0.007") ID x 1/16" OD | SL5CWPK | ea. |
| 9055-021 | 10 μL Sample Loop | 0.25 mm (0.010") ID x 1/16" OD | SL10WPK | ea. |
| 055-022 | 20 μL Sample Loop | 0.25 mm (0.010") ID x 1/16" OD | SL20WPK | ea. |
| 9055-023 | 50 μL Sample Loop | 0.51 mm (0.020") ID x 1/16" OD | SL50WPK | ea. |
| 055-024 | 100 μL Sample Loop | 0.51 mm (0.020") ID x 1/16" OD | SL100WPK | ea. |
| 9055-025 | 200 μL Sample Loop | 0.51 mm (0.020") ID x 1/16" OD | N/A | ea. |
| 9055-026 | 500 μL Sample Loop | 0.76 mm (0.030") ID x 1/16" OD | SL500WPK | ea. |
| 9055-027 | 1.0 mL Sample Loop | 0.76 mm (0.030") ID x 1/16" OD | SL1KCWPK | ea. |
| 9055-028 | 2.0 mL Sample Loop | 0.76 mm (0.030") ID x 1/16" OD | SL2KCWPK | ea. |
| 055-029 | 5.0 mL Sample Loop | 0.76 mm (0.030") ID x 1/16" OD | N/A | ea. |
| 055-033 | 10 mL Sample Loop | 0.76 mm (0.030") ID x 1/16" OD | N/A | ea. |
| PEEK LOOF | PS FOR 7725, 7725I, PR/EV700-100 | | | |
| 123-227 | 1 μL Sample Loop (models PR/EV700-100 and EV750-100 only) | Internal groove | N/A | ea. |
| 7755-015 | 2 μL Sample Loop (models 7725, 7725i, and 9725(i) only) | Internal groove | N/A | ea. |
| REPLACEM | ENT RHEFLEX FITTINGS FOR PEEK LOOPS | | | |
| 000-078 | Nut/Ferrule Set, Natural PEEK, 5/16-24, for 1/8" OD loops | | | ea. |
| 000-079 | Ferrules, Natural PEEK, for 1/8" OD loops | | | 5-pk |
| 000-251 | Ferrules, Natural PEEK, for 1/16" OD loops | | | 10-р |
| 5000-254 | Nut/Eerrule Sets, Natural PEEK, 10-32, for 1/16", OD loops | | | 10-p |



Suction Needle Adapter

Our adaptable Loop Filler Ports (Part #7012 and 9012) are used to load sample from syringe needles or luer tips. The Needle Port (Part #9013) conserves sample by minimizing the volume between the needle and the valve.



APPLICATION NOTE

Dual Mode Sample Loop Loading: Partial-Filling vs. Complete-Filling

Partial-Filling

Use the partial-filling method if you need to conserve sample, or if you want to vary sample volume frequently.

In partial-filling, the syringe sets the volume injected onto the column. There is no sample waste, and the volume injected onto the column is equal to that dispensed from the syringe. Reproducibility is 1.0% relative standard deviation (RSD). The volume of the sample loaded is limited to half the sample loop volume. For example, the most you can load into a 200 μL sample loop is 100 μL.

Complete-Filling


Use the complete-filling method if you have plenty of sample, if you do not vary sample volume, or if you need high reproducibility.

In complete-filling, the loop sets the volume loaded onto the column. Use excess sample (two to five loop volumes) to replace all the mobile phase in the loop. See Figure 2. Change the loop to vary the sample volume. Reproducibility is typically 0.1% RSD for loop sizes $\geq 5 \,\mu$ L. Accuracy is limited as loop volumes are nominal.


Q: "Which method should I use and which IDEX Health & Science sample injection valves use this method?

A: There are two types of injection valves available: dual mode and single mode. Dual mode injection valves allow both partial- and complete-filling whereas single mode injection valves allow only complete-filling. See manual injection valves, page 123.

If you are collecting experimental data, sample is scarce, and/or you want to use different sample volumes, a dual mode injector with a large volume sample loop is appropriate. Only dual mode injection valves allow the partial-filling method for easily varying your volumes (up to half your sample loop volume) by setting the syringe volume. Once you begin routine analysis, and/or you have an abundance of sample, either a dual mode or single mode injector is appropriate. Both types of injection valves allow the complete-filling method in which you overfill the sample loop. Complete-filling maximizes the reproducibility of your results.



Column LOAD



INJECT

Part No Description Qty. 7012 Stainless Steel Loop Filler Port ea 7125-054 Needle Port Cleaner ea. 9012 PEEK Loop Filler Port ea 9013 PEEK Needle Port ea 9125-076 Suction Needle Adapter (for Models 7725 and 9725) ea Australian Distributors HROMalytic +61(0)3 9762 2034 Importers & Manufacurers www.chromtech.net.au ECH nology Pty Ltd

nail : info@chromtech.net.au Tel: 03 9762 2034 . . . in AUSTRALIA

131 (of 172) 2017Sept




Injection Port Adapters

- For 360 μm OD tubing
- Mount on bracket or bulkhead

To introduce sample, connect 360 μ m OD capillary tubing to an Injection Port Adapter Assembly. This adapter accepts standard 22 gauge Hamilton-style injection syringe needles. No additional swept volume is added to the fluid pathway by this adapter, as the needle butts directly against the connecting tubing during injections. The adapter can be bulkhead mounted or mounted with the V-447 Kits.

To introduce a sample directly into a 10-32 port, purchase a M-432-03 separately.




M-432 Micro Injection Port Adapter Assembly



| • | A | В | С |
|----------------------|----------|-------|----------|
| FOR 360 µm OD TUBING | | | |
| M-432 and V-447 | P-416BLK | F-152 | M-432-03 |

> For use with Injection Valves on page 123.

This simple, biocompatible adapter is designed specifically for the Injection Valves on page 123 and can also convert any 1/4-28 flat-bottom port into a port that can accept a standard 22 gauge HPLC injection needle. This injection port adapter is adjustable, so you can create a snug fit around the needle to prevent any leaking of the analyte. In addition, this product features an internal stop that prevents you from inserting the needle too far, eliminating the possibility of damaging the valve with the needle tip.



| Part No. | Description | Qty. |
|---------------|---|------|
| MICRO INJEC | TION PORT ADAPTER | |
| FOR 360 µm 0 | DD TUBING | ea. |
| F-152 | Replacement MicroFerrule for M-432, Natural PEEK | ea. |
| M-432 | Micro Injection Port Adapter Assembly | ea. |
| M-432-03 | Replacement Tubing/Fitting Assembly for M-432 & M-433 | ea. |
| P-416BLK | Replacement Female Nut for M-432, Black PEEK | ea. |
| V-447 | Micro Injection Port Adapter Assembly Actuator Mounting Kit Includes (1) M-432 with mini-actuator bracket and (2) mounting screws | ea. |
| 1/4-28 FLAT-B | OTTOM INJECTION PORT ADAPTER | |
| P-295 | Adjustable Injection Port Adapter | ea. |
| P-296 | Replacement Tubing/Ferrule Assembly | ea. |
| | Australian Distributors | |



Wrenches, Brackets, & Replacement Fittings



Valve Wrenches

- > For convenient wrench-tightening of fittings on high pressure rotary shear valves
- > For removal of knobs on Manual Valves

The smartly designed IDEX Wrench is a double-ended slotted socket wrench that fits over 1/16" and 1/8" OD tubing. It easily loosens and tightens 1/4" and 5/16" hex head stainless steel or PEEK fittings. The "Z" shape of the IDEX Wrench provides ideal leverage for changing sample loops and fittings, and keeps one end from restricting the use of the other.

The V-103 is an Allen (hex-key) wrench designed to remove the knob from our V-101 valves (page 123). The V-104 is an Allen wrench that can be used to remove the knob from our Medium Pressure Selection and Injection Valves (also found on page 123).



Mounting Brackets

Our mounting brackets and panels of different shapes and sizes organize and provide a sturdy support for IDEX Health & Science valves. The Ring Stand Mounting Bracket now allows the valves to mount onto common laboratory equipment.

MXX Replacement Fittings

Use these replacement Ferrules and O-rings for 1/8" and 1/16" tubing with the MXX Series II valves shown on page 119. Please see the part number chart below for a list of individual part numbers.

| Part No. | Description | Description | | |
|--------------|---|---|-----------------|--|
| VALVE WRENCH | IES | | | |
| 6810 | IDEX Wrench | | ea. | |
| MOUNTING BRA | ACKET ACCESSORIES | | | |
| 7160 | Mounting Panel | | ea. | |
| 7160-010 | Valve Angle Bracket | | ea. | |
| 7160-029 | Ring Stand Mounting Bracket | | ea. | |
| VALVE BRACKE | | | | |
| M-615-1 | Mounting Bracket for IDEX Health & Science Swit | Mounting Bracket for IDEX Health & Science Switching Valves | | |
| M-615-2 | Mounting Bracket for IDEX Health & Science Inje | ction and Selection Valves | ea. | |
| REPLACEMENT | FITTINGS | | | |
| 7770-039 | Ferrules for 1/8" OD Tubing | | 25-pk | |
| 7770-040 | Ferrules for 1/8" Tubing | | 50-pk | |
| 7770-044 | Ferrules for 1/16" OD Tubing | | 25-pk | |
| 7770-124 | O-rinas for 1/16" OD Tubina | | 25-pk | |
| | HROMalytic +61(0)3 9762 2034 | Australian Distributors Importers & Manufacurers www.chromtech.net.au | www.idex-hs.con | |

Website NEW: www.chromalytic.net.au E-mail: info@chromtech.net.au Tel: 03 9762 2034... in AUSTRALIA

134







FLOW REGULATING VALVES

Our Flow Regulating Valves include specifically designed valves that are used to control or stop the flow of a stream and are ideal for use if your application involves low frequency of use or demands operator control. A variety of types and styles of valves allow you to manage directional flow. In addition, we offer replacement cartridges for all of our flow regulating valves.

- 135 **CHECK VALVES**
- **140** MICRO-SPLITTER VALVES
- MICRO-METERING VALVES 141
- SHUT-OFF VALVES 142



Inline Cartridge Check Valves

- Low cracking pressures
- Less than 150 μL internal volume
- Materials of construction: PEEK; perfluoroelastomer (CV-3001); gold-plated stainless steel spring (CV-3001); ethylene propylene (CV-3011); and stainless steel spring (CV-3011)

Our cartridge-style Inline Check Valves are designed to limit flow to one direction. These assemblies withstand system pressures of 1,000 psi (69 bar). The cracking pressures for the Inline Check Valve Cartridges are 1.5 psi (0.1 bar) for the CV-3001 and 3 psi (0.2 bar) for the CV-3011. Tolerance on the cracking pressure for CV-3001 is \pm 0.5 psi (0.03 bar) and \pm 1.5 psi (0.1 bar) on CV-3011.



| Part No. | Description | Includes | Swept Volume | Qty. |
|-------------------|--|---|--------------|------------------|
| INLINE CAR | TRIDGE CHECK VALVES | | | |
| CV-3000 | Inline Check Valve Assembly for 1/16" OD tubing | (1) CV-3001, (2) XP-215 | 96 µL | ea. |
| CV-3001 | Inline Check Valve Cartridge for CV-3000 | | 91 µL | ea. |
| CV-3010 | Inline Check Valve Assembly for 1/8" OD tubina | (1) CV-3011, (2) XP-315 | 100 µL | ea. |
| | HROMalytic +61(0)3 9762 2034 | Australian Distributors Importers & Manufacurers www.chromtech.net.au | | |
| Website NEW | : www.chromalytic.net.au E-mail : info@chromtech.net.au Te | I: 03 9762 2034 in AUSTRALIA | | 135 (of 172) 201 |



Add back-flow protection to any 1/4-28 flat-bottom port

- ▶ 15 psi (1 bar) and 3 psi (0.2 bar) cracking pressure versions
- > Excellent chemical resistance
- Materials of construction: PEEK; PCTFE; perfluoroelastomer; PTFE (CV-3301 and CV-3302); stainless steel (CV-3301 and CV-3302); or gold-plated stainless steel (CV-3315 and CV-3316)

Standard 1/4-28 Inline Check Valves

Connect these Inline Check Valves to any 1/4-28 flat-bottom port. Then thread your 1/4-28 flat-bottom fitting into the check valve to connect the tubing. Once installed, the spring-actuated sealing system eliminates back flow, helping to prevent upstream contamination or damage. In addition, the unique design of this product eliminates the additional tubing cuts and connections required to install conventional inline check valves.









Standard, Outlet 1/4-28 FB Male to 1/4-28 FB Female 3 psi (0.2 bar) cracking pressure



- > 1/4-28 Inline Check Valves and Non-Metallic Check Valves with 1/4-28 flat-bottom ports (next page) can be used with any 1/4-28 Flangeless, Super Flangeless™, and VacuTight™ fitting on pages 45, 39, and 42, respectively, of the Fittings Chapter.
- Micro-Volume Inline Check Valves and Non-Metallic Check Valves with 10-32 coned ports (next page) can be used with any 10-32 polymer Fingertight or SealTight[™] fitting on page 36. Connect capillary tubing using the optional ferrules listed on page 35 or the NanoTight[™] Fittings and Tubing Sleeves on page 37.

| Part No. | Description | Cracking Pressure | Qty. |
|----------------------|---|-------------------|------|
| STANDARD 1/4-2 | 28 INLINE CHECK VALVES | | |
| CV-3301 | Inlet Check Valve, 1/4-28 FB, M to 1/4-28 FB, F* | 15 psi (1 bar) | ea. |
| CV-3302 | Outlet Check Valve, 1/4-28 FB, M to 1/4-28 FB, F* | 15 psi (1 bar) | ea. |
| CV-3315 | Inlet Check Valve, 1/4-28 FB, M to 1/4-28 FB, F* | 3 psi (0.2 bar) | ea. |
| CV-3316 | Outlet Check Valve, 1/4-28 FB, M to 1/4-28 FB, F* | 3 psi (0.2 bar) | ea. |
| * M = Male (external |) threads: F = Female (internal) threads: C = Coned: FB = Flat-Bottom | | |





Non-Metallic Check Valves

- Cracking pressure of 8 psi (0.6 bar)
- Excellent chemical resistance
- Materials of construction: PEEK and perfluoroelastomer, suitable for biological applications

Our Non-Metallic Check Valves are biocompatible and delivers a low cracking pressure. With a swept volume of only 7.4 µL, our Check Valve is perfect for applications where low flow path volume is critical, such as delivery to lab-on-a-chip, single-cell analysis and microor nano-LC post-column derivatization. Once installed, this check valve helps prevent back flow and the potential for contamination or damage to sensitive upstream equipment.

10-32 Micro-Volume Inline Check Valves

With a swept volume of only 7.4 µL, our 10-32 Micro-Volume Inline Check Valves are perfect for applications where low flow path volume is critical, such as delivery to lab-on-achip, single-cell analysis and micro- or nano-LC post-column derivatization. Once installed, this check valve helps prevent back flow and the potential for contamination or damage to sensitive upstream equipment.



Micro-Volume Inline 10-32 C Female to 10-32 C Female



Check valves are specified by:

- **> Cracking Pressure:** the pressure required for the valve to open in the direction of the arrow.
- **Maximum Pressure:** the maximum pressure the valve can experience in the reverse direction without leaking backwards.
- **)** Back Pressure Created: the amount of back pressure generated by the check valve with 50 mL/min room temperature water flowing in the direction of the arrow.

SPECIFICATIONS & DETAILS

| | Swept Volume | Thru-Hole | Max. Pressure Rating | Back Pressure Created | Cracking Pressure Tolerance |
|-------------------|-----------------|------------------|-------------------------|--------------------------|--------------------------------|
| TANDARD 1/4 | -28 FB | | | | |
| V-3301, V-3302 | 20 μL | 0.020" (0.50 mm) | 2,000 psi (138 bar) | 45 psi (3.1 bar) | ± 5 psi (0.34 bar) |
| V-3315, V-3316 | 16 μL | 0.020" (0.50 mm) | 2,000 psi (138 bar) | 10 psi (0.7 bar) | ± 1.5 psi (0.10 bar) |
| ONMETALLIC | : 10-32 CONED N | IICRO-VOLUME | | | |
| V-3500 | 7.4 µL | 0.010" (0.25 mm) | 3,000 psi (207 bar) | 25 psi (1.7 bar) | ± 5 psi (0.34 bar) |

10-32 Micro-Volume Inline Check Valves



Non-Metallic Check Valves (Cont.)

- > Low cracking pressure of 1 psi (0.07 bar)
- Multiple configurations for different applications
- > Excellent chemical resistance
- Materials of construction: PEEK and perfluoroelastomer

APPLICATION NOTE

- The CV-3320 or CV-3321 style can be connected to any 1/4-28 flat-bottom port for trouble-free back flow protection.
- > When using a pump after the analytical column, consider placing a CV-3330 Check Valve after the column to prevent fluid from the post-column pump from flowing backwards through the column. This product also serves as an excellent nonmetallic alternative to our CV-3010 (page 135) in sparging applications where the mobile phase may be corrosive to the stainless steel or ethylene propylene components inside the CV-3010 assembly.
- ▶ The CV-3335 Inlet and CV-3336 Outlet Check Valves allow tubing larger than 1/16" OD (up to 1/8") to be connected into a 10-32 coned internal port. Use both of these check valves when attaching a larger-volume sample loop to an analyticalscale injection valve. This setup limits the flow of the sample into the loop to one direction, minimizing back flow and sample carry-over.
- The CV-3340 is useful in virtually any high pressure fluid pathway using 1/16" or smaller OD tubing, where limiting the direction of flow is desirable.

1/4-28 & 10-32 Inline Check Valves

Our 1/4-28 & 10-32 Non-Metallic Inline Check Valves provide excellent backflow protection for sensitive equipment along with outstanding chemical resistance guaranteed by the PEEK polymer and perfluoroelastomer construction. Metal-free composition makes these check valves perfect for use with corrosive fluids or biological samples.

These check valves function well up to moderately-high pressure applications. Low internal volume also allows them to be used in areas where flow path volume is important; however, higher flow rates can pass through with minimal pressure drop.





SPECIFICATIONS & DETAILS

| | Swept Volume | Max. Pressure Rating | Back Pressure Created | Cracking Pressure Tolerance |
|------------------|--------------|----------------------|-----------------------|-----------------------------|
| CV-3320, CV-3321 | 37 μL | 2,000 psi (138 bar) | 30 psi (2.1 bar) | ± 0.5 psi (0.03 bar) |
| CV-3330 | 34 µL | 2,000 psi (138 bar) | 30 psi (2.1 bar) | ± 0.5 psi (0.03 bar) |
| CV-3335, CV-3336 | 49 µL | 2,000 psi (138 bar) | 30 psi (2.1 bar) | ± 0.5 psi (0.03 bar) |
| CV-3340 | 34 µL | 2,000 psi (138 bar) | 30 psi (2.1 bar) | ± 0.5 psi (0.03 bar) |
| CV-3322, CV-3323 | 49 µL | 2,000 psi (138 bar) | 30 psi (2.1 bar) | ± 0.5 psi (0.03 bar) |
| CV-3324, CV-3325 | 182 μL | 2,000 psi (138 bar) | 30 psi (2.1 bar) | ± 0.5 psi (0.03 bar) |



NOTE

Upon initial use — or following a period of extended inactivity — the cracking pressure for these check valves may be somewhat higher than the stated cracking pressure.

1/4-28 & 10-32 Inline Check Valves

| Part No. | Description | Cracking Pressure | Thru-Hole | Qty. |
|----------|---|-------------------|------------------|------|
| NONMETA | ALLIC 1/4-28 AND 10-32 INLINE CHECK VALVES | | | |
| CV-3320 | Inlet Check Valve, 1/4-28 FB, M to 1/4-28 FB, F* | 1 psi (0.07 bar) | 0.020" (0.50 mm) | ea. |
| CV-3321 | Outlet Check Valve, 1/4-28 FB, M to 1/4-28 FB, F* | 1 psi (0.07 bar) | 0.020" (0.50 mm) | ea. |
| CV-3322 | Inlet Check Valve, 1/4-28 FB, M to 1/4-28 FB, F* | 1 psi (0.07 bar) | 0.040" (1.0 mm) | ea. |
| CV-3323 | Outlet Check Valve, 1/4-28 FB, M to 1/4-28 FB, F* | 1 psi (0.07 bar) | 0.040" (1.0 mm) | ea. |
| CV-3324 | Inlet Check Valve, 1/4-28 FB, M to 1/4-28 FB, F* | 1 psi (0.07 bar) | 0.060" (1.60 mm) | ea. |
| CV-3325 | Outlet Check Valve, 1/4-28 FB, M to 1/4-28 FB, F* | 1 psi (0.07 bar) | 0.060" (1.60 mm) | ea. |
| CV-3330 | Inlet/Outlet Check Valve, 1/4-28 FB, F to 1/4-28 FB, F* | 1 psi (0.07 bar) | 0.020" (0.50 mm) | ea. |
| CV-3335 | Inlet Check Valve, 1/4-28 FB, F to 10-32 C, M* | 1 psi (0.07 bar) | 0.020" (0.50 mm) | ea. |
| CV-3336 | Outlet Check Valve, 1/4-28 FB, F to 10-32 C, M* | 1 psi (0.07 bar) | 0.020" (0.50 mm) | ea. |
| CV-3340 | Inlet/Outlet Check Valve, 10-32 C, F to 10-32 C, F* | 1 psi (0.07 bar) | 0.020" (0.50 mm) | ea. |
| | | | | |

* M = Male (external) threads: F = Female (internal) threads: C = Coned: FB = Flat-Bottom





Quick-Stop Luer Inline Check Valve

- Check valve protection with luer convenience
- > Remains open when engaged
- > Materials of construction: PEEK, perfluoroelastomer, and gold-plated stainless steel spring



- ▶ 1/4-28 Inline Check Valves and Non-Metallic Check Valves with 1/4-28 flatbottom ports (next page) can be used with any 1/4-28 Flangeless, Super Flangeless[™], and VacuTight™ fitting on pages 45, 39, and 42, respectively, of the Fittings Chapter.
- Micro-Volume Inline Check Valves and Non-Metallic Check Valves with 10-32 coned ports (page 137) can be used with any 10-32 polymer Fingertight or SealTight[™] fitting on page 36. Connect capillary tubing using the optional ferrules listed on page 35 or the NanoTight™ Fittings and Tubing Sleeves on page 37.

The Quick-Stop Luer Check Valve is designed to provide inline luer connect/disconnect convenience without the mess and hazard of spills. Just connect the valve assembly to your inline tubing using standard 1/4-28 flat-bottom fittings (see pages 39 – 46). The check valve is automatically opened once the luer connection is engaged, allowing flow in either direction. Disconnecting the luer union causes the check valve to close. Please see the "Application Note" on this page for specific ideas regarding use of this valve.



APPLICATION NOTE

Inlet Solvent Reservoir:

Quickly change your solvent on the low pressure end of an HPLC system, while preventing potentially hazardous spills! Just install a Quick-Stop Luer Check Valve Assembly between your solvent reservoir and the pump, with the valve towards the bottle. The valve will prevent solvent leakage from the line coming from the reservoir, while the check valves in your pump prevent spills from the line leading to the pump. With both lines still full of solvent, this system also helps reduce the need to reprime your pump.

FIA Sample Injection:

The Quick-Stop Luer Check Valve provides a practical means to introduce a sample into FIA and other low pressure systems, when used in conjunction with a P-612 Pressure Relief Valve Tee (page 148). Simply connect the Tee into the appropriate flow path line with the included fittings and thread the P-697 Quick-Stop Luer Valve onto the 1/4-28 male end of the Tee. Sample can then be introduced conveniently by using a standard luer-tipped syringe. The check valve is automatically opened when the syringe is attached and closed when the syringe is removed.

Post Column Derivitization:

For post-column derivitization, place a CV-3000 Inline Check Valve on the effluent side of your column to prevent derivatizing agents from flowing backwards and poisoning the column. Placement on the post-column reagent line will also prevent mobile phase from contaminating the reagent if the auxiliary pump fails.

Helium Sparging Tank Protection:

Try the CV-3010 Assembly, designed specifically for degassing (sparging) lines to prevent solvent backup if the sparging gas runs out. This check valve will help prevent potential solvent cross-contamination and damage to the gas regulating valve.

| Part No. | Description | | Include | s | Swept Volume | Qty. |
|-------------|--|---|----------|-----------------|--------------|------------------|
| QUICK-STO | P LUER CHECK VALVE | | | | | |
| P-696 | Quick-Stop Luer Check Valve Assembly | | (1) P-69 | 7, (1) P-655 | 127 µL | ea. |
| P-697 | Quick-Stop Luer Check Valve | | | | 107 μL | ea. |
| P-699 | Bulkhead Ouick-Stop Luer Valve | | (1) nut/ | lock washer set | 107 μL | ea. |
| | HROMalytic +61(0)3 9762 2034 | Australian Distributors Importers & Manufacurers www.chromtech.net.au | | VALVE | | |
| Website NEV | V : www.chromalytic.net.au E-mail : info@chromtech.net.au Te | el: 03 9762 2034 in AUSTR | RALIA | | | 139 (of 172) 201 |



Micro-Splitter Valves

- > For interfacing LC-MS systems
- Adjustable split stream flow rates
- Versions for up to 800 psi (55 bar) and up to 4,000 psi (276 bar)



APPLICATION NOTE


- With an incoming flow rate of 1 mL/min using room temperature water and equal pressures on both outlet lines, the minimum split flow rate is 2 µL/min for the standard micro-splitter valves and 4.8 µL/min for the high pressure micro-splitter valves.
- All Micro-Splitter Valves have been tested at flow rates to 100 mL/min, with a maximum resulting pressure drop of only 45 psi (3.1 bar) when the valve is fully opened.

Our Micro-Splitter Valves are designed to accurately split and control a low-flow stream off a single incoming supply.

Choose between 1/4-28 flat-bottom and 10-32 coned threaded versions.

The High Pressure Micro-Splitter Valves are designed to operate successfully up to 4,000 psi (276 bar) and the standard Micro-Splitter valves are pressure rated to 800 psi (55 bar).

The Graduated Valve offers many of the benefits and features of Micro-Splitter Valves, plus the ability to adjust and set the split flow to repeatable settings. This allows documentation of settings and the resulting flow rates for easier method development. The graduations also make it easier to employ the valve in a system used to run multiple analyses that require different split flow rates.



P-451 Standard Micro-Splitter Valve



SPECIFICATIONS & DETAILS

| Part No. | Valve Type | Threads | Internal Volume ¹ (closed/fully open) | Max. Operating Pressure |
|-----------|----------------------|---------|--|-------------------------|
| P-450 | Standard | 1/4-28 | 2.1 / 4.1 μL | 800 psi (55 bar) |
| P-451 | Standard | 10-32 | 1.2 / 2.8 µL | 800 psi (55 bar) |
| P-460S, T | High Pressure | 10-32 | 1.2 / 2.8 µL | 4,000 psi (276 bar) |
| P-470 | High Pres. Graduated | 10-32 | 1.2 / 2.8 µL | 4,000 psi (276 bar) |
| 4 mm | | | (a a a a a a a a a a a a a a a a a a a | |


¹ The supply and waste port thru-holes have IDs of 0.020" (0.50 mm). The ID for the split-stream port thru-hole is 0.020" (0.50 mm) in standard versions; in capillary versions it is 0.010" (0.25 mm).

| Part No. | Description | Includes | Qty. | | | |
|------------------|---|------------|------|--|--|--|
| MICRO-SPLI | TTER VALVES | | | | | |
| P-450 | Standard, 1/4-28, Biocompatible | (3) XP-235 | ea. | | | |
| P-451 | Standard, 10-32, Biocompatible | (3) F-120 | ea. | | | |
| P-460S | High Pressure, 10-32, with Stainless Steel Needle | (3) F-120 | ea. | | | |
| P-460T | High Pressure, 10-32, with Titanium Needle | (3) F-120 | ea. | | | |
| GRADUATE | GRADUATED MICRO-SPLITTER VALVES | | | | | |
| P-470 | High Pressure Graduated, 10-32, with Stainless Steel Needle | (3) F-120 | ea. | | | |
| * Use with the I | MicroTight Tubing Sleeves, found on page 52. | | | | | |



Micro-Metering Valves

- > Flow rates as low as 3.5 µL/min*
- → 1/4-28 flat-bottom and 10-32 coned designs available
- **> Materials of construction:** PEEK, PTFE
- * At 1.0 mL/min incoming flow rate with room temperature water.



1/4-28 threaded PEEK Micro-Metering Valve

For fine control of fluid flow rates, Micro-Metering Valves can reduce outgoing flow to as low as 3.5 µL/min*. These needle valves are perfect for use with peristaltic pump fluid-transfer applications, mass spectrometry, and fraction collection.

Our Micro-Metering Valves can also be used to regulate gas flow in helium sparging lines and as a flow-dependent variable back pressure regulator. For flow independent regulation of back pressure, please see page 147.

Flow path materials are PEEK polymer and PTFE. All versions of this valve have 0.020" (0.50 mm) thru-holes.

🌣 APPLI

APPLICATION NOTE

Back Pressure Considerations

The Micro-Splitter Valves are designed to work when both effluent flow path pressures are nearly identical. However, the split flow path will often have higher back pressure than the waste flow path, making it hard to achieve any split flow at all. There are two possible solutions. Place a back pressure regulator (page 143) on the waste flow path that is equal to or slightly greater than the pressure on the split flow path. Or, switch the two effluent pathways such that the split flow pathway is attached to the "waste" port on the valve and the waste flow pathway is attached to the "split" port on the valve. (Please Note: This second method may result in a loss of adjustment sensitivity.)

Multi-Column and Detector Systems

Does your work require analyses with multiple columns and detectors that use the same mobile phase? If so, install one of our High Pressure Micro-Splitter Valves after your injector. A single injection can then be split to two separate columns and detector systems, at two different flow rates. This economical set-up eliminates the need for an additional pump and injector valve, while allowing data to be obtained simultaneously.

Post-Detector Interfacing

Use a Standard Micro-Splitter Valve to route fluid exiting an initial detector to other devices, such as a mass spectrometer and a fraction collector. The valve will split and reduce the flow rate to that required for MS interfacing, while diverting the remainder of the flow to the collector (a back pressure regulator may also be required for this set up, available on page 143).

Other Applications

These valves are also suited for other applications, such as adapting a standard HPLC system to handle microbore analyses. For more information and plumbing diagrams for this application and those listed above, please contact your local distributor or IDEX Health & Science directly.

| Part No. | Material | OD Tubing | Thru-hole | Internal Volume* | Includes | Qty. | | |
|-------------------|---|-----------|------------------|---------------------|------------|------|--|--|
| MICRO-METERI | MICRO-METERING VALVES | | | | | | | |
| P-445 | PEEK, Black | 1/16" | 0.020" (0.50 mm) | 7.7 µL | (2) XP-230 | ea. | | |
| P-446 | PEEK, Black | 1/16" | 0.020" (0.50 mm) | 7.2 µL | (2) F-120 | ea. | | |
| P-447 | PEEK, Black | 1/8" | 0.020" (0.50 mm) | 7.7 µL | (2) XP-330 | ea. | | |
| * Maximum interna | * Maximum internal volume, with valve fully open. | | | | | | | |







Shut-Off Valves

- > Biocompatible, all-polymer flow path
- Available for 1/16" and 1/8" OD tubing
- > Pressure rated to 500 psi (34 bar)

Stop a flow stream quickly with IDEX Health & Science biocompatible Shut-Off Valves. The bodies are manufactured from either PEEK or ETFE, and both versions feature a PCTFE rotor, making them highly resistant to chemical attack. The blue colorant used in some valve configurations has proven not to leach out with common HPLC solvents.

Connect semi-rigid or rigid tubing, such as PEEK, stainless steel or fluoropolymer, with the 1/4-28 Flangeless Fittings provided. Soft tubing, such as PharMed® or Tygon® (see page 14), may be connected to these valves using our 1/4-28 barbed adapters, found on page 90.





| Part No. | Material | OD Tubing | Thru-hole | Internal Volume* | Includes | Qty. |
|--------------------|------------------------------|-----------|-----------------|---------------------|-----------------------|------|
| SHUT-OFF VALV | ES, BIOCOMPATIBLE | | | | | |
| P-721 | ETFE, Natural | 1/8" | 0.040" (1.0 mm) | 10.0 μL | (2) P-335, (2) P-300N | ea. |
| P-732 | PEEK, Natural | 1/16" | 0.020" (0.5 mm) | 2.5 µL | (2) XP-235 | ea. |
| P-733 | PEEK, Natural | 1/8" | 0.040" (1.0 mm) | 10.0 µL | (2) XP-335 | ea. |
| P-782 | ETFE, Blue | 1/16" | 0.020" (0.5 mm) | 2.5 µL | (2) XP-235 | ea. |
| P-783 | ETFE, Blue | 1/8" | 0.040" (1.0 mm) | 10.0 μL | (2) XP-335 | ea. |
| * Maximum internal | volume with valve fully open | | | | | |





BACK PRESSURE REGULATORS

Back Pressure Regulators (BPR) are designed to enhance system performance through outgassing prevention and improved pump check valve efficiency. It includes 5 and 20psi assemblies (replacement cartridges not available), a variety of pressure rated cartridges and assemblies, PEEK and stainless steel BPR holders, high pressure adjustable BPR for pressure between 2000 and 5000psi and ultra-low volume BPRs set to 100 and 500psi.

- 144 ULTRA-LOW VOLUME BACK PRESSURE REGULATOR
- **145** BACK PRESSURE REGULATOR ASSEMBLIES
- 146 BACK PRESSURE REGULATOR HOLDERS
- 147 BACK PRESSURE REGULATOR CARTRIDGES
- **148** PRESSURE RELIEF VALVES






Ultra-Low Volume Back Pressure Regulators (BPRs)

Wetted flow path materials: PEEK, perfluoroelastomer, and ETFE

- Available pressure settings of 100 or 500 psi (7 or 34 bar)
- Low swept volume of only 6 μL

Our Ultra-Low Volume Back Pressure Regulators (BPRs) were developed to minimize swept volume, which is especially important for multi-detector applications. With a maximum swept volume of only 6 μ L*, it is nearly impossible to detect these BPRs as part of your fluid pathway. To minimize the swept volume added to your flow path, we recommend trimming the length of the attached tubing. And because the flow path is completely polymeric, you are assured of biocompatibility.



Please Note: Our Ultra-Low Volume Back Pressure Regulators cannot be used as check valves due to their unique internal design. Try our Micro-Volume Inline Check Valve on page 137.

* The maximum internal swept volume listed above is for the back pressure regulator only and does not include the volume of the attached tubing lines



SPECIFICATIONS & DETAILS

| | Back Pressure Setting psi (bar) | Flow Rate Recommendations | Recommended Pressure Range psi (bar) | 1/16" OD Tubing |
|-------|------------------------------------|--|---|-----------------|
| M-410 | 100 2 (7)2 | Optimal: 100 µL–1 mL/min Max.: 4 mL/min | 40–150 (3–10) | PEEK, 0.010" ID |
| M-412 | 500 ² (34) ² | Optimal: 100 µL–1 mL/min Max.: 4 mL/min | 250–525 (17–36) | PEEK, 0.010" ID |
| M-420 | 100 ³ (7) ³ | Optimal: 3–8 mL/min Max.: 10 mL/min | 40–150 (3–10) | PEEK, 0.020" ID |

- ¹ All data generated using water at room temperature.
 ² Set at a flow rate of 0.5 mL/min.
- 3 Set at a flow rate of 5 mL/min.

| Part No. | Description | Pressure Setting | Tubing OD | Includes | Swept Volume | Qty. |
|-----------------|--------------------------------------|--|-------------|----------|--------------|------|
| ULTRA-LOW V | OLUME BPRs | | | | | |
| M-410 | Low Flow | 100 psi (7 bar) | 1/16" | XP-230 | 6 μL | ea. |
| M-412 | Low Flow | 500 psi (34 bar) | 1/16" | XP-230 | 6 μL | ea. |
| M-420 | Hiah Flow | 100 psi (7 bar) | 1/16" | XP-230 | 6 μL | ea. |
| | HROMalytic | +61(0)3 9762 2034 Australian D. Importers & M. www.chromto | anufacurers | | | |
| Website NEW : v | www.chromalytic.net.au E-mail : info | in AUSTRALIA | | 144 (o | f 172) 20 | |





Back Pressure Regulator Assemblies


Choose from our line of Biocompatible and Stainless Steel BPR Assemblies, each complete with a replaceable, factory preset cartridge (except the 5 and 20 psi versions).

Our BPR Assemblies create incremental back pressures ranging from 5 to 1,000 psi (0.3)

Our BPR Assemblies create incremental back pressures ranging from 5 to 1,000 psi (0.3 to 69 bar). The Biocompatible BPR Assemblies feature a PEEK holder; polymer-based fittings; biocompatible BPR cartridges and wrenches for tightening. Stainless Steel BPR Assemblies feature the same biocompatible BPR cartridges with a 316 stainless steel holder and polymer fittings.







High Pressure Adjustable BPR Includes One-Piece Fingertight Fittings for 1/16" OD tubing

High Pressure Adjustable Back Pressure Regulator

Materials of construction: PEEK, perfluoroelastomer, and PTFE

The biocompatible P-880 High Pressure Adjustable BPR offers the flexibility to adjust your system back pressure between 2,000 and 5,000 psi (138 and 345 bar), independent of the flow. Only 10% fluctuation in pressure generally occurs with flow rates of 0.1–10 mL/min. Lower or higher flow rates will lead to greater fluctuations in pressure. To achieve the desired back pressure setting, simply turn the thumbscrew while monitoring your system pressure. Because this product creates such high back pressure, please check system component specifications prior to using to avoid damaging any sensitive components.

| Part No. | Pressure Setting | Holder Material | Includes | Swept Volume | Qty. |
|----------|--------------------|-----------------|----------------------------------|--------------|------|
| | | | | | |
| P-790 | 5 psi (0.3 bar) | PEEK | (2) XP-215 | 134 µL | ea. |
| P-791 | 20 psi (1.4 bar) | PEEK | (2) XP-215 | 134 µL | ea. |
| P-785 | 40 psi (2.8 bar) | PEEK | (1) P-761, (2) XP-215 | 131 µL | ea. |
| P-786 | 75 psi (5.2 bar) | PEEK | (1) P-762, (2) XP-215 | 131 μL | ea. |
| P-787 | 100 psi (7 bar) | PEEK | (1) P-763, (2) XP-215 | 131 μL | ea. |
| P-788 | 250 psi (17 bar) | PEEK | (1) P-764, (2) XP-235 | 102 μL | ea. |
| P-789 | 500 psi (34 bar) | PEEK | (1) P-765, (2) P-250, (2) LT-115 | 96 μL | ea. |
| P-455 | 1,000 psi (69 bar) | PEEK | (1) P-796, (2) P-250, (2) LT-115 | 89 µL | ea. |
| U-605 | 40 psi (2.8 bar) | SST | (1) P-761, (2) XP-201 | 129 µL | ea. |
| U-606 | 75 psi (5.2 bar) | SST | (1) P-762, (2) XP-201 | 129 μL | ea. |
| U-607 | 100 psi (7 bar) | SST | (1) P-763, (2) XP-201 | 129 μL | ea. |
| U-608 | 250 psi (17 bar) | SST | (1) P-764, (2) XP-201 | 99 µL | ea. |
| U-609 | 500 psi (34 bar) | SST | (1) P-765, (2) XP-201 | 93 µL | ea. |
| U-610 | 750 psi (52 bar) | SST | (1) P-795, (2) P-250, (2) LT-115 | 91 μL | ea. |
| | | | | | |







Back Pressure Regulator Holders

P-465 PEEK and U-469 Stainless Steel BPR Holders work with any of our replacement BPR Cartridges. Each holder comes with fittings for 1/16" OD tubing (see below). The U-469 Holder is surface-treated to prevent galling, a potential problem with large, threaded metal parts.

Please Note: These Back Pressure Regulator Holders are designed to allow each cartridge to operate at its stated pressure setting when tightened to 20 in-lbs. of torque. To approximate this level of torque, first finger tighten the Holder, then tighten an additional 1/8-1/4 turn with the supplied wrenches.





- Using PEEK tubing and supplied fittings. Using stainless steel tubing and supplied fittings

| Part No. | Pressure Setting | Holder Material | Includes | Swept Volume | Qty. |
|--------------------|------------------------------------|------------------------|---|--------------|---------------|
| BPR HOLDERS | | | | | |
| P-465 | Biocompatible BPR | PEEK | (2) P-250, (2) LT-115 | 7 μL | ea. |
| U-469 | High Pressure BPR | SST | (2) F-300 | 4 μL | ea. |
| | HROMalytic | +61(0)3 9762 2034 | Australian Distributors Importers & Manufacurers www.chromtech.net.au | | |
| Website NEW: w | www.chromalytic.net.au E-mail : in | fo@chromtech.net.au Te | l: 03 9762 2034 in AUSTRALIA | 146 | i (of 172) 20 |



Back Pressure Regulator Cartridges



APPLICATION NOTE

Small gas bubbles often form as solvent moves from the high pressure of an HPLC column to the low pressure environment leading to the detector. This outgassing can cause erratic baseline readings and loss of sensitivity. Placing a BPR (usually a 40–100 psi) after the detector provides an excellent, low-cost method for reducing this problem by maintaining enough back pressure on the mobile phase to keep gases dissolved in solution.

A back pressure regulator can also be used as a pump preload for low and fluctuating pressure applications. Many of today's pumps require a steady back pressure to function properly. Install an IDEX Health & Science BPR (usually 500–1,000 psi) between the pump and the injector to enhance pump performance.

Caution: Do not exceed the maximum operating pressure of your system — please refer to the operating manuals for your system components before choosing the appropriate BPR.

- > Proven outgassing protection
- > Flow-independent pump preload for greater pump efficiency
- > 5 to 1,000 psi cartridges and assemblies available

Back Pressure Regulators are designed to enhance system performance through outgassing prevention and improved pump check valve efficiency.

IDEX Health & Science back pressure regulators include:

- ▶ 5 and 20 psi assemblies (replacement cartridges not available)
- **)** 40, 75, 100, 250, 500, 750, and 1,000 psi cartridges and assemblies
- > PEEK and stainless steel BPR holders
- ▶ High pressure adjustable BPR for pressures between 2,000 and 5,000 psi
- Ultra low volume BPRs set to 100 and 500 psi (page 144)
 For flow control options try the Micro-Metering Valves found on page 141.

Back Pressure Regulator Replacement Cartridges

Materials of construction: PEEK, ETFE, perfluoroelastomer, and gold-plated stainless steel

These replacement cartridges will operate in any of the standard BPR holders shown on this page. These cartridges create back pressures from 40 to 1,000 psi (2.8 to 69 bar)—all independent of flow except as noted below.

The recommended operating flow rate range for our BPR Cartridges is 0.1 mL-10 mL/min. Within this range, the amount of back pressure created by the BPR Cartridges and Assemblies will not vary more than $\pm 10\%$. Lower or higher flow rates may result in larger pressure fluctuations.



ECH no logy Pty Ltd

Website NEW: www.chromalytic.net.au E-mail: info@chromtech.net.au Tel: 03 9762 2034 . . . in AUSTRALIA

P-761 40 psi BPR Cartridge



P-796 1,000 psi BPR Cartridge

COLOR CODING

| Part No. | Pressure Setting | Body | End-Cap | Swept Volume | Qty. |
|-----------------------|--------------------|---------------|-------------------------|--------------|------|
| BPR CARTRIDGES | | | | | |
| P-761 | 40 psi (2.8 bar) | Tan | Blue | 125 μL | ea. |
| P-762 | 75 psi (5.2 bar) | Tan | Yellow | 125 µL | ea. |
| P-763 | 100 psi (7 bar) | Tan | Red | 125 μL | ea. |
| P-764 | 250 psi (17 bar) | Tan | White | 95 μL | ea. |
| P-765 | 500 psi (34 bar) | Tan | Green | 89 µL | ea. |
| P-795 | 750 psi (52 bar) | Black | Blue | 87 μL | ea. |
| P-796 | 1.000 psi (69 bar) | Black | Green | 83 µL | ea. |
| H | ROMalytic +61 | 0)3 9762 2034 | Australian Distributors | | |

147 (of 172) 2017Sept




Pressure Relief Valves

Prevent system over-pressurization

Our Pressure Relief Valves are ideal for preventing system over-pressurization. These products protect system components by diverting fluid flow automatically when inline pressure exceeds the set limit. Choose between preset 100 psi (7 bar) and 5 psi (0.3 bar) assemblies, both shipped with Flangeless Fittings. The 100 psi version is a good, general purpose valve, while the 5 psi version is perfect for protecting syringe and peristaltic pump systems. The void volume of both relief valves is low due to the small 0.020" (0.50 mm) thru-holes in the valve tee body.

If you wish to have the Pressure Relief Valve open at a different pressure than 5 or 100 psi, simply combine one of the other replacement Back Pressure Regulator Assemblies listed on page 145 with the P-612 Pressure Relief Valve Tee. Choose the P-612S for larger bore tubing and higher flow applications.



| Part No. | Description | Pressure Settin | g Tubing OD | Includes | Swept Volume | Qty. |
|--------------|--------------------------|-----------------|-------------------------|------------------|--------------|------|
| PRESSURE REL | IEF VALVES | | | | | |
| U-455 | Pressure Relief Assembly | 5 psi (0.3 bar) | 1/16" | XP-201 | 148 µL | ea. |
| U-456 | Pressure Relief Assembly | 100 psi (7 bar) | 1/16" | XP-201, wrenches | 139 μL | ea. |
| P-612 | Pressure Relief Tee Only | | 1/16" | XP-201 | 14 μL | ea. |
| P-612S | Pressure Relief Tee Only | | 3/16" | XP-201 | 348 µL | ea. |
| | HROMalytic +6 | 1/013 9762 2034 | Australian Distributors | | | |

FLUIDICS



egassers

Degassers improve fluidic instrument precision and reliability by removing dissolved gases from fluids before they outgas and form problem causing bubbles. Three main types of bubble removing products are available. AF based degassers offer the widest range of chemical compatibility and are used to eliminate retention shifts and baseline fluctuations. Silicone based degassers offer the highest flow rate capabilities for water based systems such as diagnostic and life science instrumentation to improve dispense accuracy and reliability. Poridex based products provide rapid bubble remove for locations where bubble introduction cannot be avoided.

> 151 **DEBUBBLERS DEGASSING SYSTEMS**



Debubblers



APPLICATION NOTE

- Liquid handling
- IVD
- > HPLC/UHPLC
- O₂ and CO₂ removal

In medical analyzers, bubbles interfere with critical volumetric reagent dispenses and cause sample failures, wasting time and money. Because bubbles adhere to nearly every part of a dispensing system, high velocity or induced turbulent flow is often used to displace and discharge bubbles from the flow stream and into a waste area. These alternative processes waste reagents and are time consuming, unpredictable, and may additionally require designing the system to recognize bubbles are present. Regardless of how the systems are designed, aqueous systems will always be subject to the laws of physics that cause out-gassing during changes in fluid temperature, pressure, or chemicals mixture. In fluid applications like these, debubblers are the optimal solution to capture and remove formed bubbles to prevent sample dispense inaccuracies, and degassing is ideal to prevent downstream bubble formation from recurring.

Remove Bubbles, Dissolved Gas, or Both!

Dissolved gases and bubbles in system liquids cause dispense volume anomalies in many instruments, negatively affecting both dispense precision and analytical accuracy. Now you have a choice of components for actively removing bubbles with or without also removing dissolved system gases. Online Vacuum Degassing offers operating convenience, high efficiency and low operating costs compared to other common degassing technologies.

Debubbler/Degasser

Combines Vacuum Degassing with Active Bubble Removal

- Improves instrument performance reduces downtime due to bubble formation.
- Fewer false positives due to reduction of partial reagent dispenses.
- Easily integrates into any pump, degassing tray, or stand-alone degassing application.
- Designed for use with water based solutions with no surfactants. Active degassers are recommended for other solutions.

Active Debubbler

Remove Bubbles in Fluid Stream Before or After the Pump

- Improves instrument performance reduces downtime due to bubble formation.
- > Fewer false positives due to reduction of partial reagent dispenses.
- Easily integrates into any pump, degassing tray, or stand-alone degassing application.

Transfer-Line Degasser

Removes Dissolved Gases During Fluid Transfer


- ▶ Eliminates baseline fluctuations for improved detector sensitivity.
- Coaxial design reduces number of connections, improves reliability.
- > Single lumen design increases degassing reliability.




| | ACTIVE DEBUBBLER | DEBUBBLER/ DEGASSER | TRANSFER-LINE DEGASSER |
|--|---------------------|------------------------|---------------------------|
| Perfect for applications that require dissolved gas like oxygen for reaction kinetics | ~ | | |
| Improves dispense precision by capturing and removing bubbles | ~ | * | |
| Eliminates false positives and reduces reagent waste by improving instrument performance | ~ | ~ | |
| Easily integrates into fluidic path | ✓ | * | * |
| Creates stable instrument performance across system and environmental fluctuations | ~ | ~ | ~ |
| Prevents the formation of bubbles downstream of the degasser | | * | * |
| Eliminates fluctuations for improved detector sensitivity and accuracy by preventing bubble formation | | ~ | ~ |
| Improves fluidic system reliability because coaxial design reduces the number of connections | | | * |
| Flexible design can be implemented as transfer line in new instruments or existing instruments that don't have space available | | | ~ |
| Minimizes fluidic system internal volumes to reduce reagent cost | | | ~ |


Debubblers (Cont.)

Overall Dimensions Please note: These drawings are not actual size.





9000-1544 Debubbler / Degasser, 2.5 mL




9000-1545 Debubbler/Degasser, 5 mL





9000-1549 Transfer-Line Debubbler, 1.1 meter





Transfer-Line Degasser Implementation







Gas bubbles are actively removed from a flowing liquid stream by vacuum via the PORIDEX membrane.

| | ACTIVE DEBUBBLERS | DEGASSER/DEBUBBLERS | 200 KPA (30 PSI) @ 25 °C |
|--|--|---|--|
| Bubble Removal (volume of air removed/min @ 10 mL/min H₂O) | Up to 30 cc | Up to 30 cc | N/A |
| Degassing Efficiency [†] @ 1 mL/min H₂O | N/A | $2.5 \text{ mL size: } 36\% \text{ O}_2 \text{ removal, } 5.0 \text{ mL size: } 55\% $ $\text{O}_2 \text{ removal}$ | $<$ 4 ppm dissolved O_2 at 5 mL/min |
| Membrane Material | PORIDEX® | PORIDEX | PORIDEX |
| Wetted Materials | PORIDEX, Polyolefin, FEP, ETFE, Ultem® | PORIDEX, Polyolefin, FEP, ETFE, Ultem | PORIDEX, Polyolefin, FEP, ETFE |
| Solvent Compatibility | Solutions > 50% aqueous. Not compatible with detergent concentrations > 0.05%. | Solutions > 50% aqueous. Not compatible with detergent concentrations > 0.05%. | Solutions > 50% aqueous. Not compatible with detergent concentrations > 0.05%. |
| Standard Bubble Trap Volume | 2.5 / 5.0 mL | 2.5 / 5.0 mL | N/A |
| Transfer-Line Volume | N/A | 2.5 / 5.0 mL | < 4 mL |
| Max. Operating Pressure | 200 kPa (30 psi) @ 25 °C | 200 kPa (30 psi) @ 25 °C | 200 kPa (30 psi) @ 25 °C |
| Max. Operating Temperature | 40 °C | 40 °C | 40 °C |
| Recommended Vacuum Level | Minimum 16 kPa absolute | Minimum 16 kPa absolute | Minimum 16 kPa absolute |
| Liquid Connection | 1/4-28 fitting system | 1/4-28 fitting system | 1/4-28 fitting system |
| Vacuum Connection | Tubing vacuum port(s) for 1/8" (3 mm) ID elastomeric tubing | Tubing vacuum port(s) for 1/8" (3 mm) ID elastomeric tubing | Tubing vacuum port(s) for 1/8" (3 mm) ID elastomeric tubing |
| Pressure Drop | 0.8 mm Hg / mL / min (assumes laminar flow and viscosity of 1 cP) | 0.8 mm Hg / mL / min (assumes laminar flow and viscosity of 1 cP) | 0.8 mm Hg / mL / min (assumes laminar flow and viscosity of 1 cP) |
| † Debubblina / deaassina efficiency o | can be optimized based on flow rate. fluid to be d | declassed, and das to be removed. | |







Dissolved gases are actively removed from a flowing liquid stream by vacuum via the IDEX Health & Science AF® membrane.

APPLICATION NOTE

Why Degas Your Mobile Phase?

Dissolved air in HPLC mobile phases can result in flow rate instability and baseline disturbance.

Flow rate instability: Non-degassed mobile phase can outgas in the pump head, causing bubbles to be formed and trapped inside the head or check valves. These bubbles can cause flow disturbances and pressure fluctuations, resulting in flow rate instability.

Baseline disturbance: As the mobile phase passes through the column, it experiences a large pressure drop. Non-degassed mobile phase can outgas due to this pressure differential, causing air bubbles to form. Air bubbles passing through or lodging in the flow cell cause detection disturbances, exhibited as baseline noise.

Why Use a Degassing System?

Helium sparging is a common means of degassing HPLC solvents. This method has its drawbacks, however. Helium tanks are expensive and bulky, and solvent backup and contamination are concerns. In addition, helium sparging can change the composition of a premixed mobile phase over time, due to the difference in the evaporation rates of mobile phase components.

In contrast, the IDEX Health & Science Degassing System has none of these drawbacks, and it is extremely fast and efficient at removing dissolved gases — more efficient than helium sparging or PTFE-based degassing systems.

Tubing Connections

We recommend ETFE tubing (page 27) be used to limit regassing of mobile phase between the degasser and your pump. ETFE is recommended because of its superior impermeability to gases (compared to PTFE, FEP, and PFA tubing). Applicable flangeless fittings for 1/8" OD tubing are found on page 45.

GPC and HFIP Applications

Standard degassing chambers, with PEEK bulkhead unions, are <u>not recommended</u> for GPC applications or for use with HFIP (hexafluoroisopropanol). Special GPC "hardened" versions are available. Please contact us for more information.



Website NEW: www.chromalytic.net.au E-mail: info@chromtech.net.au Tel: 03 9762 2034... in AUSTRALIA

Degassing tubing is flexible and therefore can be coiled to shorten the overall length or used to transfer the fluid within an instrument to the next desired location.

Debubblers







- Analytical and Prep scale models
- Ultra-high degassing efficiency
- > Low volume, easy to prime
- Patented control eliminates baseline fluctuations
- Inert flow path
- > 5+ year lifetime

Full Stand Alone egassing Systems

No Troubles with Bubbles Anymore!

- prevent bubble formation in your Fluidic System with BIOTECH DEGASi line of Degassers.

Dissolved gasses in a fluidic system can often cause troubles. When the pressure or the temperature changes, the dissolved gasses can form bubbles which affect the accuracy, precision and performance of your equipment. On-line degassing is a very efficient way of removing dissolved gasses from the liquid and preventing bubble formation.







GPC

ELUENT



LIMITED

VOLUME



SEMI

PREP

| Systec AF™ Internal Volume | 480 µl | 480 µl | 100 µl | 925 µl |
|--|--------|--------|--------|--------|
| No of Channels | 2-5 | 2-5 | 2-5 | 2-5 |
| Biocompatible Flow Path | Υ | N | Υ | Υ |
| Approximate Max Flow Per Channel (ml/min) | 3 | 3 | 0.5 | 6 |
| Inner Dimension Flow Path (mm) | 1 1 / | 1 1 / | 1 1 / | 1 1/1 |




FLUIDICS



| Part Number | Number of Channels | Internal Volume |
|-----------------------|--------------------|-----------------|
| BIOTECH DEGASI CLASSI | | |
| 0001-6352-A | 2 | 480 µl |
| 0001-6353-A | 3 | 480 µl |
| 0001-6354-A | 4 | 480 µl |
| 0001-6355-A | 5 | 480 µl |
| BIOTECH DEGASI GPC | | |
| 0001-6622 | 2 | 480 µl |
| 0001-6623 | 3 | 480 µl |
| 0001-6624 | 4 | 480 µl |
| 0001-6625 | 5 | 480 µl |
| BIOTECH DEGASI MICRO | | |
| 0001-6352-S | 2 | 100 μΙ |
| 0001-6353-S | 3 | 100 μΙ |
| 0001-6354-S | 4 | 100 μΙ |
| 0001-6355-S | 5 | 100 μΙ |
| | | |

| Part Number | Number of Channels | Internal Volume |
|-----------------------|--------------------|-----------------|
| BIOTECH DEGASI SEMI-P | REP | |
| 0001-6352-L | 2 | 925 µl |
| 0001-6353-L | 3 | 925 µl |
| 0001-6354-L | 4 | 925 µl |
| 0001-6355-L | 5 | 925 µl |
| BIOTECH DEGASI PREP | | |
| 0001-2053 | 2 | 5.3 ml |
| 0001-6482 | 2 | 8.4 ml |
| 0001-6484 | 2 | 13.8 ml |
| BIOTECH DEGASI PREP+ | | |
| 0001-0120 | 1 | 23 ml |
| 0001-0220 | 2 | 23 ml |
| 0001-0420 | 4 | 23 ml |
| BIOTECH DEGASI COMPA | АСТ | |
| 0004-2285 | 2 | 285 μΙ |
| 0004-4285 | 4 | 285 μΙ |
| 0004-6285 | 6 | 285 µl |
| BIOTECH DEGASI HIGH F | LOW | |
| HF.500-S Stand Alone | 1 | 60 ml |
| HF.500-A OEM version | 1 | 60 ml |
| | | |



NEW!

A B B B C











HIGH

FLOW RATE





PREP

ORGANIC SOLVENTS

COMPACT APPLICATION

HIGH FLOW

| Systec AF™ Internal Volume | 5.3 / 8.4 / 13.8 ml | 23 ml | 285 μΙ | 60 ml (silicone) |
|--|---------------------|--------|---------|------------------|
| No of Channels | 2 | 1-4 | 2, 4, 6 | 1 |
| Biocompatible Flow Path | N/Y/Y | Υ | Υ | Υ |
| Approximate Max Flow Per Channel (ml/min) | 15 / 25 / 50 | 75-100 | 2 | 500 |
| Inner Dimension Flow Path (mm) | 1.91 | N/A | 0.89 | N/A |





Column Hardware

We offer an extensive line of HPLC and UHPLC Column Hardware that has been optimized to enable selectivity, efficiency, and high-quality separation performance in your flow paths. To browse the full line of our column hardware portfolio please visit: www.idex-hs.com/column-hardware.html

- Biocompatible materials for LC and UHPLC columns
- > Accessories for column protection and packing

Our line of column hardware includes protective accessories and connection products that enhance column functionality. Our columns come in a variety of sizes and materials to meet your system requirement.

157 GUARD HARDWARE

158 GU

GUARD COLUMNS



Guard Hardware

Our Guard Hardware portfolio offers a variety of guard cartridges and holders to meet your system requirements. These include guards, guard holders, analytical cartridges, and tools designed to protect your system's valuable columns and help ensure high performance and excellent retention time.



Prep Scale Guard Holders

- > 21.2 mm and 30 mm ID column protection
- Improves plate count and symmetry
- New anti-rotation feature aids guard holder assembly

Iso-Prep™ Guard is a guard cartridge system designed to protect valuable prep columns. It offers superior column protection for adsorptive samples and a proven sample distribution mechanism via a precision machined holder. Iso-Prep Guard is ideal for protecting prep columns with no degradation of peak shape or plate count.

The high-performance guard protects columns in two ways. First, it acts as a filter, trapping particles in the frits. Second, when the guard cartridge is packed with the same material as the prep column, it removes compounds that irreversibly adsorb to the packing material.



Tools


Our tools includes the Iso-Prep Guard Scraper that is specifically designed as a tool used to dress a cartridge upon packing in order to set the frit at the correct depth.



- > 60:40 Acetonitrile:Water
- **>** 50 mL/min

- > Kromasil 10 µm C18
- ▶ Backpressure: 100 psi (7 bar)

10 μ C18 100 x 21.2 mm, 60:40 Acetonitrile:Water, 20 mL/min



Without Iso-Prep Guard — 38,150 Plates/M 1.24 As

With Iso-Prep Guard — 41,920 Plates/M 1.20 As

| Part No. | Description | | Qty. |
|---------------|---|-------------------------|------|
| PREP SCALE GU | JARD HOLDERS | | |
| 9197-P | Iso-Prep Guard Holder | | ea. |
| 9197-P-AR | Iso-Prep Guard Holder, 21.2mm, Anti Rotation | | ea. |
| 9197-20 | Iso-Prep Guard Cartridge, 21.2 mm x 1 cm, 1 Frit | | ea. |
| 9196-P | Iso-Prep Guard Holder, 30 mm | | ea. |
| 9196-P-AR | Iso-Prep Guard Holder, 30 mm, Anti Rotation | | ea. |
| 9196-20 | Iso-Prep Guard Cartridge, 30 mm x 1 cm, 2 µm Frit | | ea. |
| 8083-MOD | Iso-Prep Guard, 30 mm Frit, 2 μm | | ea. |
| TOOLS | | | |
| 9197-S | Iso-Prep Guard Finishing Tool | | ea. |
| 9196-S | Iso-Prep Guard Finishing Tool, 30 mm | | ea. |
| | HROM = 1 y & i = +61(0)3 9762 2034 | Australian Distributors | |

Bio



Guard Columns

We offer wide selection of pre-packed, cartridge style guard columns in addition to pack-it-yourself hardware. There are options for micro flow applications as well as analytical-scale applications. These guard columns can be immediately implemented into a system for your convenience.



Guard Column Kits

- ▶ 100% biocompatible flow path
- Pressure rated to 4,000 psi (276 bar)
- > Wetted materials are Titanium and PEEK
- Reusable holder complete with fingertight fittings

Insert one of our analytical guard columns between the injection valve and column of your HPLC system to extend the life of your column and help ensure reproducible results. Convenient, prepacked PEEK polymer cartridges complete the system and are available in a variety of bonded phases to match your column chemistry held in place by Titanium frits.

The C-270 Stainless Steel Guard Column Holder is engineered for high-pressure applications to 4,000 psi (276 bar). Each of these holders is surface treated to prevent galling*, a potential problem with threaded metal parts.

The flow path of the C-270 Guard Column Holder is biocompatible. Each comes complete with fittings for 1/16" OD tubing, and can be used with any of the C-28X or C-7XX guard column cartridges listed on the following pages.

* Galling is a form of "cold welding." When two fittings manufactured from the same metal are wrench-tightened too tightly, they can "weld" together, making it virtually impossible to separate the two components.



Microbore Guard Columns


- Ideal for Microbore HPLC
- Easily dry packed (or slurry packed with adapter)
- Made of PEEK polymer and stainless steel

This ultralow volume guard column (1.0 mm ID x 2 cm length) is ideal for narrow-bore chromatography. The unpacked guard column allows you to exactly match the chemistry of your column, resulting in optimum column protection. The total packing volume of $16.2 \, \mu L$ ensures maximum column efficiency and analytical column protection.

Frits often become plugged before a guard column is contaminated. The two 0.5 μ m frits included with this guard column can be changed in minutes. Optional 2 μ m frits may be purchased separately (C-408).



Packing Material Specifications: The cartridges on this page are packed with 5 µm or 10 µm basedeactivated 80 Å spherical silica.








Analytical Guard Columns

▶ Easy to pack

The C-130B is our most popular guard column. HPLC users find this column easy to pack and extremely economical. This narrow-bore short column (2.0 mm ID x 2 cm length) creates only a slight pressure increase with virtually no detectable theoretical plate loss when used with a 3 mm ID or larger column. The 2 μ m frits are easy to change, prolonging the life of the guard column. With only 62 μ L packing volume per guard column, a 3 g bottle of packing material will pack about 30 guard columns.






Guard Column Holders

- ▶ 10 mm ID column protection
- Convenient cartridge system
- Easy to pack

The internal volume of our semi-prep guard column is just 780 μ L, which only requires approximately 1.50 g of packing material — ideally the same material used in your semi-prep column. The C-1000 Holder will hold to high pressures, and is specially treated to prevent galling*. Use standard 10-32 coned fittings (not included) to connect your 1/16" OD tubing.

* Galling is a form of "cold welding." When two fittings manufactured from the same metal are wrench-tightened too tightly, they can "weld" together, making it virtually impossible to separate the two components.





C-1000 Semi-Prep Guard Column
Fittings, tubing, and column shown are not included.



Guard Columns (Cont.)



Accessories

Our accessories include packing adapters and seal, packing funnels, and other useful products to work with your column systems.



Guard Column Cartridges

Our biocompatible Guard Column Cartridges are conveniently prepackaged and are offered in multiple packs. We offer these cartridges in a variety of bonded phases to match your column chemistry. These cartridges are reusable and economical.



APPLICATION NOTE

Why Use A Guard Column?

A guard column can increase the life of your analytical column significantly. Use a guard column with the same packing as your column — it will act as a chemical filter, removing strongly retained materials in your sample that might otherwise contaminate your analytical column. And, it is more economical to replace a guard column cartridge than to buy a new analytical column.



All Guard Columns featured on this page include 10-32 Coned threads. Use any of the 10-32 coned fittings starting on page 32 to connect tubing to these guard columns.



Signs Indicating the Guard Column Needs to be Changed

- System pressure build-up
- > Faster than usual retention times
- Reduced resolution

Guard Columns

Website NEW: www.chromalytic.net.

| Part No. | Description | Includes | Qty |
|--------------------------|--|------------------------------|-----|
| GUARD COLUMN KITS | | | |
| C-281 | 2.0 mm ID C18 Cartridges (6-pk) with (1 |) C-270 Assembly | ea. |
| MICROBORE GUARD | COLUMNS | | |
| 1.0 mm ID x 2 cm Unpack | red | | |
| C-128 | Guard Column | (2) C-128-31 | ea. |
| C-128-31 | 0.5 µm Stainless Steel Replacement Frit | | ea. |
| C-408 | 2 µm Stainless Steel Replacement Frit | | ea. |
| ANALYTICAL GUARD | COLUMNS | | |
| 2.0 mm ID x 2 cm Unpack | ced | | |
| C-130B | Guard Column | (2) A-100 | ea. |
| C-130-20 | Packing Funnel | | ea. |
| A-100 | 2 µm Stainless Steel Replacement Frit | | ea. |
| A-103 | 0.5 µm Stainless Steel Replacement Frit | t | ea. |
| BIOCOMPATIBLE GUA | ARD COLUMN HOLDERS | | |
| 10 mm ID x 1 cm | | | |
| C-270 | High Pressure, Stainless Steel, with (2) F | -200 Fittings, Biocompatible | ea. |
| GUARD COLUMN HOL | LDERS | | |
| 10 mm ID x 1 cm | | | |
| C-1000 | Semi-Prep Cartridge Guard Column Ho | older | ea. |
| C-1035 | Semi-Prep Cartridge | | ea. |
| ACCESSORIES | | | |
| C-128-40 | Slurry Packing Adapter | | ea. |
| C-130-40 | Slurry Packing Adapter | | ea. |
| C-1030 | Threaded Frit Cap with 2 µm Stainless S | Steel Frit | ea. |
| GUARD COLUMN CAR | RTRIDGES | | |
| 2.0 mm ID x 1 cm, 10 μm | Silica | | |
| C-282 | Reversed Phase C18 | Reversed Phase C18 | |
| 4.3 mm ID x 1 cm, 5 μm S | iilica | | |
| C-751 | 4.3 mm ID C18 Cartridges (6-pk) with (1 |) C-270 Assembly | ea. |
| C-752 | Reversed Phase C18 | | 10- |

| | FLUIDICS |
|--|--|
| | |
| | · |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | Ž |
| | FLUIDIC CONNECTIONS |
| | |
| | 20 |
| | |
| | |
| | |
| | |
| | |
| | VALVES |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | <u> </u> |
| | ASSS |
| | DEGASSERS |
| | |
| | |
| | |
| | |
| | 8 |
| | The second secon |
| | |
| | |
| | COLUMN HARDWARE |
| HROM = 15 +61(0)3 9762 2034 Australian Distributors | |
| ECH nology Pty Ltd www.chromtech.net.au | 404 (-5.470) 2047 |
| Website NEW: www.chromalytic.net.au E-mail: info@chromtech.net.au Tel: 03 9762 2034 in AUSTRALIA | 161 (of 172) 2017 <mark>Sept</mark> |

Educational Materials

Help is only a click away to access a wealth of product and application support information, including:

Chemical Compatibility

Use our compatibility charts to discover what materials are best suited for your application.

www.idex-hs.com/chemical-compatibility

Chromatography Forum

Browse or join this group that is devoted to the analytical and scientific communities.

www.idex-hs.com/chromatography-forum

Conversion Tools

Find useful Conversion Tools that allow you to easily convert temperatures, measurements, pressure ratings, and other functions.

www.idex-hs.com/conversions

FAOs

Explore our FAQs section to assist in finding the answers to issues you may be facing.

www.idex-hs.com/faqs www.semrock.com/technical-faq.aspx

Fittings Primer

Learn the complexities and functions of nutand-ferrule fittings using our Fittings Primer.

www.idex-hs.com/fittings-primer

Glossary of Terms

Our Glossary of Terms provides a convenient list of definitions and/or explanations to language that is applicable in our industry.

www.idex-hs.com/glossary

HPLC Center

Visit our HPLC Center to learn about high- and ultra-high performance liquid chromatography.

www.idex-hs.com/hplc-center

Injectors Troubleshooting

This guide presents common problems when using manual sample injection valves.

www.idex-hs.com/injectors

IP Rating

Use our IP Rating chart to discover the environmental care of our electrical equipment.

www.idex-hs.com/ip-rating

Materials Guide

Our Materials Guide allows you to view material property and solvent compatibility.

www.idex-hs.com/materials-guide

Optical Filter Plotting Tool

Dur SearchLight™ optical filter plotting tool allows you to easily select elements for your system, visualize the spectral properties, and quickly calculate critical parameters such as signal brightness, autofluorescence levels, and signal-to-noise ratio.

searchlight.semrock.com

Polymer Information

➤ Study the Polymer Information chart to learn about analytical solvent compatibility for HPLC.

www.idex-hs.com/polymer-information

Pump School

Jump over to Pump School, an introduction to rotary positive-displacement (PD) pumps.

www.idex-hs.com/pump-school

Standard Port Drawings

We offer downloadable schematics of our Standard Port Drawings that allow you to easily create your own internal ports using our fittings.

www.idex-hs.com/standard-port-drawings

Tech Tips

Examine various topics that will maximize the precision performance of our products.

www.idex-hs.com/tech-tips

Webinars & Videos

IDEX Health & Science hosts ongoing educational webinars on specific topics, and has also created a product- and process-focused video library.

www.idex-hs.com/webinars-and-videos

Literature & Downloads

Our literature and downloads section provides a wide selection of valuable information about our products and our company. You can access and download a variety of support literature, manuals, designed to assist in the use of our products and services, including:

Operation Manuals & Instruction Cards

Product operation manuals and instruction cards are here for your convenience to download and use as reference.

www.idex-hs.com/manuals

Optical Filter Journals and White Paper Articles

Download and share our extensive collection of white papers and journal reprint articles written by our expert PhD staff. Browse varied topics, from coatings for ultrafast lasers, to polarization, super-resolution microscopy, tunable bandpass filters, and more.

www.semrock.com/tech-library.aspx

Optical Filter Technical Notes

- Visit our General Filter information to learn more about the many optical filters we offer, how to find filter orientation, how to measure OD, and more.
- Learn the basics of Fluorescence Filters, including different types of fluorescence filters, applications, uses, and information.
- Our Raman & Laser Optical System Filters information pages give you an in-depth look into what you need to know and how to avoid issues integrating into your Raman & Laser Optical systems.

www.semrock.com/tech-notes.aspx

Optical Systems, Laser, & Shutter Technical Resources

Have a question or need some technical support? Consult our library of documentation, manuals, data sheets and specifications on our entire family of optical products, including:

- > Technical Note
- White papers
- Journal articles

www.mellesgriot.com/store/technical-library

Support Literature

Download informational IDEX Health & Science product and OEM catalogs, brochures, white papers, data sheets and more.

www.idex-hs.com/support-literature



NDEX

Index by Part Number

| 0001-6352-A | 154 | |
|----------------------|-----|-------------|
| 0001-6353-A | 154 | |
| 0001-6354-A | 154 | |
| 0001-6355-A | 154 | 1006-999 12 |
| 0001-6622 | 154 | 1020-05 10 |
| 0001-6623 | 154 | 1020-20 10 |
| 0001-6624 | 154 | 1328 5- |
| 0001-6625 | 154 | 1356 5- |
| 0001-6352-S | 154 | 1531 18 |
| 0001-6353-S | 154 | 1531B 18 |
| 0001-6354-S | 154 | 1532 18 |
| 0001-6355-S | 154 | 1533 18 |
| 0001-6352-L | 154 | 1534 18 |
| 0001-6353-L | 154 | 1535 18 |
| 0001-6354-L | 154 | 1536 18 |
| 0001-6355-L | 154 | 1538 18 |
| 0001-2053 | 155 | 1544 18 |
| 0001-6482 | 155 | 1560 18 |
| 0001-6484 | 155 | 1561 18 |
| 0001-0120 | 155 | 1568 18 |
| 0001-0220 | 155 | 1569 |
| 0001-0420 | 155 | 1570 18 |
| 0004-2285 | 155 | 1571 18 |
| 0004-4285 | 155 | 1572 1 |
| 0004-6285 | 155 | 1574 18 |
| HF.500-S Stand Alone | 155 | 1576 18 |
| HF.500-A OEM version | 155 | 1577 18 |
| | | |

| 1501 | 10 | 7020 | 127 |
|-----------|-----------|----------|-----|
| 1581 | 18 | 7020 | 127 |
| 1593 | 73 127 | 7021 | 127 |
| 1876 | 127 | 7022 | 127 |
| 1958-01 | 78 | 7023 | 127 |
| 1959-01 | 76 | 7024 | 127 |
| 3000 | 123 | 7025 | 127 |
| 3055-018 | 130 | 7026 | 127 |
| 3055-019 | 130 | 7027 | 127 |
| 3055-023 | 130 | 7028 | 127 |
| 3055-025 | 130 | 7029 | 127 |
| 3065-018 | 127 | 7030 | 123 |
| 3065-019 | 127 | 7030-003 | 125 |
| 3065-023 | 127 | 7030-014 | 125 |
| 3065-025 | 127 | 7030-015 | 125 |
| 3255 | 23 | 7030L | 123 |
| 3505 | 23 | 7031-05 | 107 |
| 3725-006 | 125 | 7031-20 | 107 |
| 3725-018 | 125 | 7060 | 123 |
| 3725-038 | 123 | 7060-039 | 125 |
| 3725-999 | 125 | 7060-064 | 125 |
| 3725i | 123 | 7060-065 | 125 |
| 3725i-038 | 123 | 7060-070 | 125 |
| 3755 | 23 | 7060-074 | 125 |
| 5001-999 | 125 | 7112-999 | 125 |
| 5100-999 | 125 | 7123-047 | 125 |
| 5104-999 | 125 | 7123-127 | 125 |
| 6000-076 | 66 | 7123-128 | 125 |
| 6000-078 | 61, 130 | 7123-142 | 125 |
| 6000-079 | 61, 130 | 7123-145 | 125 |
| 6000-082 | 33 | 7123-147 | 125 |
| 6000-083 | 33 | 7123-148 | 125 |
| 6000-209 | | 7123-149 | 125 |
| 6000-210 | 33 | 7123-180 | 125 |
| 6000-211 | | 7123-221 | 125 |
| 6000-251 | 61, 130 | 7123-223 | 125 |
| 6000-254 | 61, 130 | 7123-223 | 130 |
| 6000-255 | 61 | 7123-548 | 125 |
| 6000-262 | 33 | 7123-550 | 125 |
| | 33 32 | | |
| 6000-282 | 32 61 | 7123-568 | 125 |
| 6000-360 | | 7125-047 | 125 |
| 6255 | 23 | 7125-054 | 131 |
| 6505 | 23 | 7125-079 | 125 |
| 6755 | 23 | 7125-999 | 125 |
| 6810 | 51, 133 | 7150-999 | 125 |
| 7000 | 123 | 7152-999 | 125 |
| 7000-016 | 125 | 7155-999 | 125 |
| 7000L | 123 | 7160 | 133 |
| 7004-999 | 125 | 7160-010 | 133 |
| 7010 | 123 | 7160-029 | 133 |
| 7010-039 | 125 | 7410-038 | 125 |
| 7010-040 | 125 | 7410-041 | 125 |
| 7010-069 | 125 | 7410-999 | 125 |
| 7010-071 | 125 | 7413-013 | 125 |
| 7010-997 | 125 | 7501-999 | 125 |
| 7010-998 | 125 | 7502-999 | 125 |
| 7010-999 | 125 | 7507-999 | 125 |
| 7012 | 131 | 7520-030 | 125 |
| | | | |




166

| 7520-035 | 125 | 9000-1549 | 153 | 62515 | 23 |
|-----------|-----|------------|-----|-----------|--------------|
| 7520-999 | 125 | 9010 | 123 | 62520 | 23 |
| 7531-999 | 125 | 9010-051 | 125 | 62550 | 23 |
| 7650-002 | 125 | 9010-065 | 125 | 63005 | 23 |
| 7725 | 123 | 9010-999 | 125 | 65010 | 23 |
| 7725-010 | 125 | 9012 | 131 | 65015 | 23 |
| 7725-999 | 125 | 9013 | 131 | 65020 | |
| 7725i | | 9055-020 | | 65050 | |
| 7725i-188 | | 9055-021 | | 67510 | |
| 7750-038 | | 9055-022 | | 67515 | |
| 7750-070 | | 9055-023 | | 67520 | |
| 7755-015 | | 9055-024 | | 67550 | |
| | | | | | |
| 7755-020 | | 9055-025 | | 79801-999 | |
| 7755-021 | | 9055-026 | | 79861-999 | |
| 7755-022 | | 9055-027 | | 310010 | |
| 7755-023 | | 9055-028 | | 310015 | |
| 7755-024 | | 9055-029 | 130 | 310020 | 23 |
| 7755-025 | 127 | 9055-033 | 130 | 310050 | 23 |
| 7755-026 | 127 | 9085-05-10 | 106 | 315010 | 23 |
| 7755-027 | 127 | 9085-20-10 | 106 | 315015 | 23 |
| 7755-028 | 127 | 9125-043 | 125 | 315020 | 23 |
| 7755-029 | 127 | 9125-076 | 131 | 315050 | 23 |
| 7770-039 | 133 | 9125-082 | 125 | 610010 | 23 |
| 7770-040 | 133 | 9125-095 | 125 | 610015 | 23 |
| 7770-044 | 133 | 9125-999 | 125 | 610020 | 23 |
| 7770-124 | | 9196-20 | | 610050 | |
| 7770-229 | | 9196-P | | 615010 | |
| 7788-999 | | 9196-P-AR | | 615015 | |
| 7900-146 | | 9196-S | | 615020 | |
| 7900-148 | | 9197-20 | | 615050 | |
| | | | | | |
| 7900-183 | | 9197-P | | 617515 | |
| 7900-999 | | 9197-P-AR | | 617520 | |
| 7920-999 | | 9197-S | | 617550 | |
| 7960-014 | | 9650-009 | | 620015 | |
| 7960-999 | | 9725 | | 620020 | |
| 7961-999 | | 9725-999 | | 620050 | |
| 7970-999 | | 9725i | 123 | 630010 | |
| 7980-004 | | 9750-021 | | 630015 | |
| 7980-999 | 125 | 9960-002 | 125 | 630020 | |
| 7984-999 | 125 | 31005 | 23 | 630050 | |
| 7986-004 | 125 | 31505 | 23 | A-100 | 96, 104, 160 |
| 7986-999 | 125 | 32510 | 23 | A-101 | 96, 106 |
| 8020 | 127 | 32515 | 23 | A-102 | 96, 106 |
| 8021 | 127 | 32520 | 23 | A-103 | 96, 160 |
| 8022 | 127 | 32550 | 23 | A-105 | 96 |
| 8023 | 127 | 35010 | 23 | A-106 | 96 |
| 8083-MOD | 157 | 35015 | 23 | A-107 | 96 |
| 8125 | 123 | 35020 | 23 | A-120 | 96 |
| 8125-038 | 125 | 35050 | 23 | A-122 | 96 |
| 8125-097 | 125 | 37510 | | A-224 | 96, 104 |
| 8125-098 | | 37515 | | A-225 | |
| 8125-119 | | 37520 | | A-225A | |
| 8125-999 | | 37550 | | A-227A | |
| 9000-1540 | | 61005 | | A-228 | |
| | | | | A-230A | |
| 9000-1541 | | 61505 | | | |
| 9000-1544 | | 62005 | | A-231A | |
| 9000-1545 | 133 | 62510 | 23 | A-242 | 102 |



| A-243 | 102 | A-702 | 70 | F-IIZ | |
|--------|---------|----------|------------|----------|------------|
| A-302 | 102 | A-703 | 98 | F-113 | 35 |
| A-302A | 102 | A-704 | 98 | F-120x | 32 |
| A-304 | 51 | A-706 | 98 | F-124Hx | 34 |
| A-305 | 51 | A-707 | 98 | F-124Sx | 34 |
| A-309 | 102 | A-708 | 98 | F-125Hx | 34 |
| A-310 | 102 | A-710 | 98 | F-125x | 34 |
| A-311 | 102 | A-720 | 98 | F-126Sx | 34 |
| A-311A | 102 | A-722 | 98 | F-130x | 32 |
| A-313 | 104 | A-735 | 110 | F-132 | 34 |
| A-314 | 104 | B-310 | 88 | F-142Nx | 35, 37, 80 |
| A-315 | 106 | C-128 | 160 | F-142x | 35 |
| A-318 | 106 | C-128-31 | 96, 160 | F-148 | 35 |
| A-320 | 51 | C-128-40 | 160 | F-151 | 35 |
| A-327 | 28 | C-130-20 | 160 | F-152 | 34, 132 |
| A-328 | 28 | C-130-40 | 160 | F-152BLK | 34 |
| A-329 | 28 | C-130B | 160 | F-156 | 41, 93 |
| A-330 | 107 | C-140-30 | 96 | F-172 | 34 |
| A-331 | 96, 107 | C-235x | 33, 38, 41 | F-180 | 54 |
| | 96, 107 | C-236x | 33 | F-181 | 54 |
| A-337 | 96, 107 | C-270 | 160 | F-182 | 54 |
| | 96 | C-281 | 160 | F-183 | 54 |
| A-350 | 28 | C-282 | 160 | F-184 | 54 |
| A-355 | 106 | C-401 | 96 | F-185 | 54 |
| A-356 | 106 | C-402 | 96 | F-186 | 54 |
| A-360 | 107 | C-407 | 96 | F-187 | 54 |
| A-410 | 107 | | 96, 160 | F-188 | 54 |
| A-411 | 107 | | 96 | F-192x | 36 |
| A-420 | 96 | C-420 | 96, 110 | F-193x | 36 |
| | 101 | | 96, 110 | F-195x | 36 |
| A-428x | 104 | | 160 | F-196x | 36 |
| A-429x | 104 | | 160 | | 54 |
| A-430 | 104 | C-1000 | 160 | | 54 |
| A-431 | 104 | C-1030 | 160 | F-227 | 54 |
| | 101 | C-1035 | 160 | F-228 | 54 |
| A-437 | 101 | CV-3000 | 135 | F-229 | 54 |
| A-438 | 101 | CV-3001 | 135 | F-230 | 54 |
| A-440 | 101 | CV-3010 | 135 | F-231 | 54 |
| A-441 | 101 | CV-3301 | 136 | F-232 | 54 |
| A-445 | 101 | CV-3302 | 136 | F-233 | 54 |
| A-446 | 101 | CV-3315 | 136 | F-234 | 54 |
| A-451 | 101 | CV-3316 | 136 | F-237 | 54 |
| A-510 | 107 | CV-3320 | 138 | F-238 | 54 |
| A-520x | 101 | CV-3321 | 138 | F-239 | 54 |
| A-522x | 101 | CV-3322 | 138 | F-240 | 54 |
| A-550 | 101 | CV-3323 | 138 | F-241 | 54 |
| A-551 | 101 | | 138 | F-242 | 54 |
| A-610 | 111 | | 138 | F-243 | 54 |
| A-620 | 111 | | 138 | F-244 | 54 |
| A-622 | 111 | CV-3335 | 138 | F-245 | 54 |
| | 112 | | 138 | | 54 |
| | 112 | | 138 | | 54 |
| | 112 | | 137 | | 54 |
| | 112 | | 93 | | 43 |
| | 111 | | 93 | | 28 |
| | 98, 106 | | 93 | | 36 |
| | 98, 106 | | 32 | | 35 |
| | -, | | | | |





| r-330X | 41 | r-0123 | 140 | P-/U3 | öɔ |
|----------|-------------|----------|-----|-------|-----|
| P-352x | 41 | P-618 | 88 | P-704 | 69 |
| P-353x | 43 | P-619 | 88 | P-710 | 85 |
| P-355x | 41 | P-620 | 85 | P-712 | 86 |
| P-357-2x | 41 | P-622 | 85 | P-713 | 86 |
| P-357x | 41 | P-623 | 85 | P-714 | 86 |
| P-359 | 93 | P-624 | 88 | P-718 | 66 |
| P-359x | 41 | P-625 | 88 | P-719 | 88 |
| P-360x | 41 | P-626 | 66 | P-720 | 71 |
| P-363Rx | 43 | P-627 | 66 | P-721 | 142 |
| P-366x | 41 | P-628 | 88 | P-722 | 86 |
| P-372x | 99 | P-629 | 88 | P-723 | 86 |
| P-373x | 99 | P-630 | 85 | P-727 | 70 |
| P-374x | 99 | P-631 | 85 | P-728 | 70 |
| P-381 | 41 | P-632 | 86 | P-729 | 70 |
| P-387 | 41 | P-633 | 86 | P-732 | 142 |
| P-399 | 51 | P-634 | 86 | P-733 | 142 |
| P-401x | 44 | P-635 | 86 | P-742 | 69 |
| P-407x | 44 | P-642 | 88 | P-754 | 55 |
| P-416 | 34, 57 | P-645 | 85 | P-755 | 55 |
| P-416BLK | 34, 57, 132 | P-646 | 93 | P-756 | 55 |
| P-416G | 34 | P-647 | 93 | P-757 | 93 |
| P-420 | 41, 93 | P-648 | 93 | P-760 | 69 |
| P-430 | 82 | P-650 | 66 | P-761 | 147 |
| P-432 | 82 | P-652 | 66 | P-762 | 147 |
| P-440 | 82 | P-654 | 66 | P-763 | 147 |
| P-441 | 82 | P-655 | 88 | P-764 | 147 |
| P-441N | 82 | P-656 | 88 | P-765 | 147 |
| P-445 | 141 | P-657 | 88 | P-767 | 93 |
| P-446 | 141 | P-658 | | P-770 | 78 |
| P-447 | 141 | P-659 | | P-771 | 71 |
| P-450 | 140 | P-660 | 88 | P-772 | 71 |
| P-451 | 140 | P-661 | 88 | P-775 | 68 |
| P-455 | 145 | P-662 | 88 | P-777 | 68 |
| P-460S | 140 | P-663 | | P-779 | 69 |
| | 140 | P-668 | | | 142 |
| P-465 | 146 | P-669-01 | 66 | | 142 |
| | 140 | P-670 | | P-785 | 145 |
| | 44 | P-671 | | | 145 |
| P-501x | 44 | P-672 | | | 145 |
| P-507x | 44 | P-673 | | | 145 |
| P-512 | 83 | P-675 | | | 145 |
| | 83 | P-677 | | | 145 |
| | 83 | P-678 | | | 145 |
| | 83 | P-680 | | | 93 |
| | 55 | P-681 | | | 147 |
| | 55 | P-683 | | | 147 |
| | 55 | P-684 | | | 93 |
| | 71 | P-686 | | | 93 |
| | 78 | P-689 | | | 93 |
| | 34, 55 | P-691 | | | 93 |
| | 55 | P-692 | | | 93 |
| | 55 | P-694 | | | 89 |
| | 85 | P-696 | | | 89 |
| | 85 | P-697 | | | 89 |
| | 88 | P-699 | | | 89 |
| | 148 | P-702 | | | 89 |
| . 🗸 | 170 | 1 -1 02 | 00 | | 07 |

| P-837 | 89 | PK-100x | 58 | U-160 | 21 |
|------------|---|------------|----|---|-----|
| P-838 | 89 | PK-110x | 58 | U-161 | 21 |
| P-840 | 42 | PK-112 | 57 | U-162 | 21 |
| P-842x | 42 | PK-120BLKx | 58 | | 21 |
| P-844x | 42 | PK-126Hx | 58 | U-164 | 21 |
| P-845-01 | 85 | PK-126x | | | 21 |
| P-846x | | PK-132x | | | 33 |
| P-849 | | PK-152 | | | 33 |
| P-850 | | U-101 | | | 73 |
| P-851 | | U-102 | | | 33 |
| P-852 | | U-103 | | | 33 |
| P-854 | | U-104 | | - | 73 |
| P-857 | • | U-105 | | | |
| P-858 | | | | | 33 |
| | | U-106 | | - | 73 |
| P-859 | | U-107 | | | 73 |
| P-860 | | U-108 | | | 72 |
| P-861 | | U-111 | | _ | 72 |
| P-862 | | U-112 | | 0 100 1111111111111111 | 72 |
| P-863 | | U-113 | | | 72 |
| P-864 | | U-114 | | | 73 |
| P-870 | | U-115 | | | 73 |
| P-872 | | U-116 | 21 | U-450x | 33 |
| P-874 | 71 | U-117 | 21 | U-455 | 148 |
| P-875 | 68 | U-118 | 21 | U-456 | 148 |
| P-880 | 145 | U-126 | 21 | U-466 | 68 |
| P-881 | 78 | U-127 | 21 | U-466S | 68 |
| P-882 | 78 | U-128 | 21 | U-467R | 55 |
| P-885 | 68 | U-129 | 21 | U-469 | 146 |
| P-887 | 68 | U-130 | 21 | U-500 | 66 |
| P-888 | 68 | U-131 | 21 | U-504 | 66 |
| P-889 | 68 | U-132 | 21 | U-510 | 66 |
| P-890 | 68 | U-133 | 21 | | 66 |
| P-891 | | U-134 | | | 145 |
| P-920-01 | | U-135 | | | 145 |
| P-930x | | U-136 | | | 145 |
| P-931x | | U-137 | | • | 145 |
| P-938x | | U-138 | | | 145 |
| P-940x | | U-139 | | | 145 |
| P-942x | ······ | U-140 | | | 56 |
| P-945x | | U-141 | | | 66 |
| P-946x | | - | | | |
| | | U-142 | | | 56 |
| P-948x | | U-143 | | | 66 |
| P-982BLKx | | U-144 | | | 21 |
| P-1000 | | U-145 | | | 21 |
| P-1082BLKx | | U-146 | | | 21 |
| P-1087x | | U-147 | | | 21 |
| PD715-000 | | U-148 | | | 20 |
| PD715-102 | 120 | U-149 | | U-1115 | 20 |
| PD715-105 | 120 | U-150 | 21 | U-1116 | 20 |
| PD7900 | 120 | U-151 | 21 | U-1117 | 20 |
| PD7920 | 120 | U-152 | 21 | U-1120 | 20 |
| PD7960 | 120 | U-153 | 21 | U-1122 | 20 |
| PD7970 | 120 | U-154 | 21 | U-1125 | 20 |
| PD7980 | 120 | U-155 | 21 | U-1126 | 20 |
| PD7986 | | U-156 | | | 20 |
| PD9900 | 120 | U-157 | | | 20 |
| PD9960 | | U-158 | | | 20 |
| . 2 , , 30 | 120 | 5 100 | ∠1 | | 20 |



| U-1132 | . 20 | UPFP-6300350 |
|--------------|------|--------------|
| U-1133 | .20 | UPFP-6300500 |
| U-1140 | .20 | UPFP-6300600 |
| U-1141 | . 20 | UPFS-6100070 |
| U-1142 | . 20 | UPFS-6100150 |
| U-1143 | .20 | UPFS-6100250 |
| U-1145 | .20 | UPFS-6100350 |
| U-1146 | .20 | UPFS-6100500 |
| U-1148 | | UPFS-6100600 |
| UH-402 | | UPFS-6125070 |
| UH-427 | | UPFS-6125150 |
| UH-432 | | UPFS-6125250 |
| UH-436 | | UPFS-6125350 |
| UH-630 | | UPFS-6125500 |
| UH-631-01 | | UPFS-6125600 |
| UH-632 | | UPFS-6254070 |
| UH-634 | | UPFS-6254150 |
| UH-700 | | UPFS-6254250 |
| UH-750. | | UPFS-6254350 |
| | | |
| UH-752 | | UPFS-6254500 |
| UH-753 | | UPFS-6254600 |
| UH-904x | | V-100D |
| UH-906 | | V-100L |
| UP-700 | | V-100T |
| UPFP-6025070 | | V-101D |
| UPFP-6025150 | | V-101L |
| UPFP-6025250 | | V-101T |
| UPFP-6025350 | | V-240 |
| UPFP-6025500 | | V-241 |
| UPFP-6025600 | | V-340 |
| UPFP-6050070 | .49 | V-341 |
| UPFP-6050150 | .49 | V-447 |
| UPFP-6050250 | . 49 | V-450 |
| UPFP-6050350 | .49 | V-451 |
| UPFP-6050500 | . 49 | V-540 |
| UPFP-6050600 | . 49 | V-541 |
| UPFP-6075070 | . 49 | VHP-200-01x |
| UPFP-6075150 | . 49 | VHP-200x |
| UPFP-6075250 | . 49 | VHP-320x |
| UPFP-6075350 | .49 | VHP-325x |
| UPFP-6075500 | .49 | VHP-500 |
| UPFP-6075600 | .49 | VHP-501 |
| UPFP-6100070 | . 49 | |
| UPFP-6100150 | .49 | |
| UPFP-6100250 | . 49 | |
| UPFP-6100350 | . 49 | |
| UPFP-6100500 | . 49 | |
| UPFP-6100600 | . 49 | |
| UPFP-6150070 | | |
| UPFP-6150150 | | |
| UPFP-6150250 | | \mathbb{R} |
| UPFP-6150350 | | |
| UPFP-6150500 | | |
| UPFP-6150600 | | |
| UPFP-6300070 | | |
| UPFP-6300150 | | |
| 0111-0300130 | . т/ | |

| IPFP-630035049 | VHP-50510 | 04 |
|----------------|----------------|----|
| IPFP-630050049 | VHP-5061 | 04 |
| IPFP-630060049 | VHP-5501 | 06 |
| IPFS-610007049 | VHP-55110 | 06 |
| IPFS-610015049 | VHP-55510 | 06 |
| IPFS-610025049 | VHP-55610 | 06 |
| IPFS-610035049 | VHP-600 | 55 |
| IPFS-610050049 | VHP-700x | 59 |
| PFS-610060049 | VHP-920x | 60 |
| PFS-612507049 | VHP-1000 | 51 |
| IPFS-612515049 | VHP-1001 | 60 |
| IPFS-612525049 | VHP-2000 | 51 |
| IPFS-612535049 | VHP-3200 | 60 |
| IPFS-612550049 | VHP-4000 | 51 |
| IPFS-612560049 | VHP-9000 | 51 |
| IPFS-625407049 | XF-358x | 47 |
| IPFS-625415049 | XF-368x | 47 |
| IPFS-625425049 | XLT-111x | 47 |
| IPFS-625435049 | XP-130x | 56 |
| IPFS-625450049 | XP-131x | 56 |
| IPFS-625460049 | XP-132x | 56 |
| -100D123 | XP-136x | 56 |
| -100L123 | XP-137x | 56 |
| -100T123 | XP-141x | 56 |
| -101D123 | XP-201x | 47 |
| -101L123 | XP-202x | 47 |
| -101T123 | XP-218x | 47 |
| -240123 | XP-230x | 47 |
| -241123 | XP-235x | 47 |
| -340123 | XP-238x | 47 |
| -341123 | XP-245x | 47 |
| -447132 | XP-286x | 47 |
| -450123 | XP-301x | 47 |
| -451123 | XP-302x | 47 |
| -540123 | XP-305x | 47 |
| -541123 | XP-308x | 47 |
| HP-200-01x59 | XP-315x | 47 |
| HP-200x59 | XP-330x | 47 |
| HP-320x60 | XP-335x | 47 |
| HP-325x60 | XU-620x | 56 |
| HP-500104 | | |
| | XU-655x | 56 |
| HP-501104 | XU-655xXU-662x | |



For a complete list of IDEX Health & Science registered trademarks and their proper attribution, please refer to our trademarks list at: www.idex-hs.com/trademarks



UPFP-6300250......49



For ordering and technical support, please contact:

EUROPE | BIOTECH AB | +46 (0)300-56 91 80 | info@biotech.se | www.biotech.se USA | BIOTECH USA LLC | 612-703-5718 | sales@biotechusa.us | www.biotech.se JAPAN | BIONIK INC. | +81 545-38-9125 | info@bionikinc.com | www.bionikinc.com



