


2023-6



2023

RESEARCH & INDUSTRIAL  
PRODUCT CATALOG



# TABLE OF CONTENTS

|  |           |
|--|-----------|
| <b>1. INTRODUCTION</b>   | <b>3</b>  |
| Flow Control Technologies  | 4         |
| Company Overview   | 6         |
| History  | 7         |
| Our Distributors   | 8         |
| <b>2. MICROFLUIDIC EQUIPMENT</b>                                     | <b>9</b>  |
| Pressure Systems (LineUp™, MFCST™ series)                            | 10        |
| Sensors (FLOW UNIT, Flowboard, PRESSURE UNIT)                        | 20        |
| Microfluidic Valves (2-SWITCH™, M-SWITCH™, L-SWITCH™)                | 26        |
| Pressure & Vacuum Sources (FLPG+, Compact Pressure or Vacuum Source) | 30        |
| Sample Reservoirs (P-CAP, Fluiwell, Bottle-CAP series)               | 33        |
| Standalone Accessories (Pressure Manifold, Bubble Trap, ...)         | 36        |
| Surfactant (dSurf)   | 38        |
| Microfluidic Chips (RayDrop, EZ Drop, Drop-Seq, ...)                 | 40        |
| Software Solutions (OxyGEN, SDK)                                     | 46        |
| Tubing & Fitting Kits (Pressure System Kits, Sensor Kits, ...)       | 48        |
| <b>3. SYSTEMS &amp; PLATFORMS</b>                                    | <b>55</b> |
| Omi, Automated Organ-on-chip Platform                                | 56        |
| Aria, Automated Sequential Injection System                          | 58        |
| Microfluidic Complex Emulsion Platform                               | 60        |
| Cell Encapsulation Platform  | 62        |
| <b>4. MICROFLUIDIC PACKS</b>   | <b>63</b> |
| Microfluidic application packs (Cell biology, Droplet generation)    | 64        |
| Starter packs  | 78        |
| Tubing & Fitting Kits (Systems, Packs)                               | 81        |
| <b>5. FLUIGENT INDUSTRIAL</b>  | <b>85</b> |
| Customizable OEM modules   | 86        |
| Standard OEM Equipment   | 89        |
| Custom solutions   | 94        |
| <b>6. SUPPORT &amp; TOOLS</b>  | <b>96</b> |
| Tips   | 97        |
| Customer Support & Services  | 98        |
| Reference Guide  | 99        |
| Product Guide  | 104       |



1

# INTRODUCTION


# FLOW CONTROL TECHNOLOGIES

Fluigent was the first company to introduce pressure-driven flow control to research in microfluidics.

Pressure-driven flow control has multiple advantages compared to conventional syringe and peristaltic pumps for many applications.

Depending on the field of application, shear stress-related flow requirements can be different. Some studies exclude this parameter and other researchers are trying to reproduce in-vivo shear stress conditions. In both cases, precise, pulseless flow control is critical for repeatable results. Peristaltic and syringe pumps generate pulsatile and unstable flows. **Pressure-driven pumps have been shown to provide superior performance.**


## Flow rate profile over the time



The Fluigent Flow Controller used in this experiment is our MFCS™-EX, based on the patented FASTAB™ technology. This technology is the best adapted to manipulating fluid volumes at the sub-microliter scale compared to syringe, peristaltic or piston pumps.

Moreover, in microdroplet generation, **droplet size and frequency** are directly linked to the flow rates of the continuous and dispersed phases. Flow rate stability is critical for having repeatable and monodispersed droplets. **Pressure pumps provide more stable flow** leading to **better experimental data.**

## Flow rate response over the time



The Fluigent Flow Controller used in this experiment is our LineUp Flow EZ™. This component allows you to save experimental time, precious samples and expensive reagents with significantly shorter response times compared to syringe pumps.

## Flow controller system comparison

| Features                     | Fluigent LineUp™ series | Fluigent MFCST™ series | Other Pressure-Based Solutions | Syringe Pumps |
|------------------------------|-------------------------|------------------------|--------------------------------|---------------|
| Free standing (no PC needed) | ✓                       | ✗                      | ✗                              | ~             |
| Modular & stackable          | ✓                       | ✗                      | ✗                              | ~             |
| Very short response time     | ✓                       | ✓                      | ~                              | ✗             |
| Pulseless flow               | ✓                       | ✓                      | ✓                              | ✗             |
| Control & monitor display    | ✓                       | ✗                      | ~                              | ~             |
| Compact                      | ✓                       | ✓                      | ~                              | ✗             |
| Integrated pressure source   | ✗                       | ✓                      | ~                              | ✗             |



# COMPANY OVERVIEW

Fluigent develops, manufactures, sells and supports innovative fluid handling solutions for a variety of rapidly growing applications where fluid control is critical.

Since 2005 we have delivered more than **3000 Fluigent systems**. These include the MFCSTM-EZ (Microfluidic Flow Control System), the LineUpSTM series (Pressure-based flow controllers) and the ESSSTM valves (Easy Switch SolutionsSTM).

RESEARCH  
INSTRUMENTS



FLUIGENT  
INDUSTRIAL




We strive to provide “Smart Microfluidic” solutions. Our modular system architecture enables our customers to focus on their application.



In addition to our Paris area headquarters, we have 2 local subsidiaries to be closer to our customers and provide local service and support: **Fluigent Deutschland GmbH** in Jena, Germany serves Northern and Eastern European customers, and **Fluigent Inc.** near Boston, USA for our customers in the Americas.

# HISTORY

Fluigent is the leader in microfluidic control



# OUR DISTRIBUTORS



At Fluigent, we strongly believe that local customer and technical support is essential to maintaining customer satisfaction.

Since 2006, we have delivered our instruments to more than 40 countries worldwide thanks to our direct offices in France (Fluigent SAS), the USA (Fluigent Inc.), Germany (Fluigent Deutschland GmbH), and our network of distributors and partners.

Through regular training and support, we help our distributors serve our customers from around the world every day. No matter where you are in the world, Fluigent will be there.

To contact your local Fluigent representative, or if you are interested in partnering with Fluigent for OEM development, distribution, or services, please send us an email at [sales@fluigent.com](mailto:sales@fluigent.com).



2

# RESEARCH INSTRUMENTS

# LINEUP™ SERIES

The most advanced flow control system



## LineUp™ series modules

Our LineUp™ product range is the next generation of microfluidic systems. With the **Flow EZ™** or **Push-Pull** modules one can precisely regulate and control pressure & vacuum, the **LINK** and **LINK COM** modules provide communication to a computer or any external instrument using TTL ports, USB cable or Serial port communication. The **Adapt** is used to connect Flow EZ™ modules with different pressure ranges without the need of additional pressure sources. The **P-SWITCH** allows to multiply the outlets of the system and the **SWITCH EZ** controls microfluidic valves. The entire system can be controlled without a PC using local control and can also be monitored by **Fluigent Software** to extend its capabilities and benefit from automation. Select modules you need and combine them together.

*Focus on your  
experiment with the  
local control dial*

*Easily adaptable  
to any setup*

*Start within minute:  
Use with or without  
a PC*

*Economical  
& expandable*

# FLOW EZ™

Microfluidic Flow Controller




## Characteristics

### LineUp Flow EZ™

| Range in mbar | Product Number |
|---------------|----------------|
| 0 to 25       | LU-FEZ-0025    |
| 0 to 69       | LU-FEZ-0069    |
| 0 to 345      | LU-FEZ-0345    |
| 0 to 1000     | LU-FEZ-1000    |
| 0 to 2000     | LU-FEZ-2000    |
| 0 to 7000     | LU-FEZ-7000    |
| 0 to -25      | LU-FEZ-N025    |
| 0 to -69      | LU-FEZ-N069    |
| 0 to -345     | LU-FEZ-N345    |
| 0 to -800     | LU-FEZ-N800    |

The **Flow EZ™** is the most advanced flow controller for pressure-based fluid control. It can be combined with a **FLOW UNIT** to control both pressure and flow rate. It can be used without a PC to minimize benchtop space, or connected to the **LINK** to benefit from **Fluigent Software** features such as graph view in real-time, automation and custom integration.



Microfluidic setup



# Characteristics

## LineUp™ Push-Pull

| Range in mbar | Product Number |
|---------------|----------------|
| -800 to 1000  | ELUPPU1000     |


The **Push-Pull** is the newest version of the Flow EZ™ with the ability to **regulate pressure and vacuum from one unit**. The module controls over the range of **-800 mbar to 1000 mbar**.

# PUSH-PULL

Vacuum Pressure-Based Controller



The **Push-Pull** can be combined with a **FLOW UNIT** to control both pressure and flow rate. It can be used without a PC to minimize benchtop space. When connected to a PC with the **LINK** one can benefit from **Fluigent Software** features.



Microfluidic setup

# P-SWITCH

Pneumatic Valve Controller



Supplied by the **Flow EZ** modules or other external pressure system, one P-SWITCH can pressurize **up to 8 reservoirs**. With a **LINK** and **Fluigent Software**, the module allows for automation behavior and actuation timing of any pneumatic or quake valves.

## Characteristics

### LineUp™ P-SWITCH


Range in mbar

-800 to 2000

Product Number

ELUPSW2000

The **P-SWITCH** is a pneumatic valve controller. By providing two regulated pressure / vacuum sources (any vacuum or pressure from **-800 mbar to +2000 mbar**), each module is able to deliver one of the two provided pressure through **8 independent outlets**.



To learn more about the P-SWITCH applications, see the [application note](#).

# Characteristics

## LineUp™ SWITCH EZ

| Valve compatibility             | Product Number |
|---------------------------------|----------------|
| M-SWITCH™, L-SWITCH™, 2-SWITCH™ | ELUSEZ         |


The **SWITCH EZ** is a module allowing one to control **Fluigent's microfluidic valves**. The module has **6 external connections** and can be combined with other **LineUp™** products to have a complete, compact system for benchtop use.

# SWITCH EZ

## Microfluidic Valve Controller



Connected valves can be **controlled or programmed** either by using the local control directly on the device or by creating a software protocol to automate **valve actuation timing**.



# ADAPT

Pressure Reducer



## Characteristics


### LineUp™ Adapt

Product Number

LU-ADP-0001

The **Adapt** is a **pressure reducer**. Placed between two pressure controllers having different input requirements, it allows one to have a system with many **different pressure channel ranges**. It reduces the outgoing pressure supply from the module to its left and provides the required pressure supply to the module on its right.

Each **LineUp™ Flow EZ** or **Push-Pull** requires a specific amount of pressure supplied. To build a single LineUp™ system with various pressure range, use the **Adapt** to reduce the pressure supply and meet exact requirements.



| Pressure requirement   | Product Number |
|------------------------|----------------|
| 150 mbar               | LU-FEZ-0025    |
|                        | LU-FEZ-0069    |
| 11000 mbar             | LU-FEZ-0345    |
|                        | LU-FEZ-1000    |
| 2100 mbar              | LU-FEZ-2000    |
| 71000 mbar             | LU-FEZ-7000    |
|                        | LU-FEZ-N025    |
|                        | LU-FEZ-N069    |
| -800 mbar              | LU-FEZ-N345    |
|                        | LU-FEZ-N800    |
| 1100 mbar<br>-800 mbar | ELUPPU1000     |

# Characteristics

## LineUp™ LINK / LINK COM

| Communication                       | Product Number       |
|-------------------------------------|----------------------|
| USB port, TTL signal                | LU-LNK-0002 (LINK)   |
| Serial port (RS-232),<br>TTL signal | ELULNK232 (LINK COM) |

# LINK LINK COM

Microfluidic Software Control



LINK COM



LINK

The **LINK** provides communication for all connected **LineUp™ series** modules to a PC for **software control**. Place it at the first position of the chain and connect it to the PC to use Fluigent Software.



Three types of connection and communication are available:



## Related products

|  |  |                                |
|--|--|--------------------------------|
| <b>OxyGEN</b><br>SSFT-OXY                  | Fluigent Software extends the devices capabilities with features such as pressure / flow rate graph views, control in real time and protocol automation. | Provided                       |
| <b>LineUp™ Supply Kit</b><br>LU-SPK-0002   | Power cable and pneumatic tubing to supply the LineUp series modules. Allow for local control without a PC.  | Provided in package            |
| <b>Chain-to-Chain kit</b><br>LU-C2C-0001   | Cable to connect data and power transmission between two modules without combining them together. Allows for flexibility and several pressure sources.   | Accessory                      |
| <b>P-CAP series</b><br>P/N on version      | Air-tight metal cap for reservoir pressurization.  | Required for liquids           |
| <b>FLOW UNIT</b><br>P/N on flow rate range | The FLOW UNIT is a high-precision flow sensor used for direct flow control.  | Required for flow rate control |
| <b>SDK</b><br>FLUIGENT-SDK                 | The Software Development Kit provides integration of devices over multiple programming languages   | Provided                       |

## See more about LineUp™ series



See more

Dedicated webpage with additional information about the complete series. Visit <https://www.fluigent.com/research/instruments/pressure-flow-controllers/lineup-series/>



Tutorials

Fluigent tutorials on Flow EZ™ pressure-based flow controller. Watch us on YouTube at [https://www.youtube.com/channel/UC8wJ\\_15yYGjbcl-DnMBPEww](https://www.youtube.com/channel/UC8wJ_15yYGjbcl-DnMBPEww)




Applications

Discover the use of LineUp™ series in a number of applications. Visit <https://www.fluigent.com/resources-support/expertise/application-notes/>



# MFCS™ SERIES

## Microfluidic Flow Control System



## MFCS™

The **MFCS™** is a microfluidic flow controller. Either **4 or 8 channels** are available with different pressure ranges for **precision operations** in microfluidic experiments. By using the FASTAB™ microfluidic patented technology, the MFCS™ generates a **constant pressure-driven flow rate** that allows for reliable and repeatable experiments.



*Easy to use*

*Adaptable:  
Independent  
pressure channels*

*Reliable and  
reproducible results:  
Pulseless flow*

*Compact:  
Save benchtop space*

## Characteristics

|                   | Range<br>in mbar                     | Product<br>Number |                | Range<br>in mbar | Product<br>Number |
|-------------------|--------------------------------------|-------------------|----------------|------------------|-------------------|
| MFCSTM-EZ         | 0 to 345                             | EZ-00345001       | MFCSTM-EX      | 0 to 345         | EX-00345001       |
|                   | 0 to 1000                            | EZ-01000001       |                | 0 to 1000        | EX-01000001       |
|                   | 0 to 2000                            | EZ-01000002       |                | 0 to 2000        | EX-01000002       |
|                   | 0 to 7000                            | EZ-07000001       |                | 0 to 7000        | EX-07000001       |
|                   | 0 to -345                            | EZ-80345001       |                | 0 to -345        | EX-80345001       |
|                   | 0 to -800                            | EZ-80800001       |                | 0 to -800        | EX-80800001       |
| Unit in mbar      |                                      |                   | Product Number |                  |                   |
| MFCSTM-EZ<br>Base | Basic                                |                   | EZ-11000001    |                  |                   |
|                   | Positive Pressure<br>Source Included |                   | EZ-source-pos  |                  |                   |
|                   | Negative Pressure<br>Source Included |                   | EZ-source-neg  |                  |                   |
| MFCSTM-EX<br>Base | Basic                                |                   | EX-11000008    |                  |                   |
|                   | Positive Pressure<br>Source Included |                   | EX-source-pos  |                  |                   |
|                   | Negative Pressure<br>Source Included |                   | EX-source-neg  |                  |                   |

## MFCSTM series

Our MFCSTM series product range is the first generation of microfluidic systems. Along with the **MCFSM-EZ** or **MCFSM-EX**, a **Manifold** can be added to redirect the pressure to multiple fluid reservoirs. The flow generated can be measured with **FLOW UNITS** and the **Flowboard**.

The MFCSTM can also have an integrated pressure source or be coupled with the **FLPG Plus**, an external pressure source (page 30). This system is controlled by **Fluigent Software Solutions**.



More product specifications  
at [www.fluigent.com](http://www.fluigent.com)



More information on  
flow control technologies  
on pages 4-5

## Related products

|                            |  |                                      |
|----------------------------|--|--------------------------------------|
| <b>OxyGEN</b><br>SSFT-OXY  | Fluigent Software extends the devices capabilities with features such as pressure / flow rate graph views, control of real time and protocol automation. | Provided                             |
| <b>SDK</b><br>FLUIGENT-SDK | The Fluigent Software Development Kit includes full integration of devices interfaces within LabVIEW, MATLAB and other IDEs.                             | Provided                             |
| <b>MFCSTM KITS</b>         | The MFCSTM Low or High Pressure Kits are specially designed to be used with any MFCSTM with any low or high pressure channel.                            | Required                             |
| <b>FLOW UNIT</b>           | The FLOW UNIT is a high-precision individual flow sensor used for direct flow control.   | Required<br>for flow<br>rate control |
| <b>FLOWBOARD</b><br>FLB    | The Flowboard is a hub that communicates between Fluigent Software and up to eight FLOW UNITS.   | Required<br>for flow<br>rate control |

# FLOW UNIT

## Bidirectional Microfluidic Flow Sensor

The **FLOW UNIT** is a **bidirectional flow sensor** compatible with our Software Solutions, our **MFCSTM series**, our **LineUp™ series** or other external control system. The sensor extends the capabilities of Fluigent instruments and allows for **direct control of flow rate**.




*High precision:  
get reliable  
results*

*Adaptable:  
Large range of  
flow rates*

*Flexible:  
Usable with  
any flow control  
system*

*Ease of use:  
Plug & play  
and easy  
combination*



PC free pressure-based flow control

# Characteristics

| FLOW UNIT+                       | XS                         | S                           |                             | M+                        |                            | L+                        |                            |
|----------------------------------|----------------------------|-----------------------------|-----------------------------|---------------------------|----------------------------|---------------------------|----------------------------|
| Product Number                   | FLU-XS                     | FLU-S-D                     |                             | FLU-M+                    |                            | FLU-L+                    |                            |
| Sensor Inner Diameter            | 25 μm                      | 150 μm                      |                             | 400 μm                    |                            | 400 μm                    |                            |
| Maximum Pressure                 | 200 bar                    | 200 bar                     |                             | 12 bar                    |                            | 12 bar                    |                            |
| Wetted Materials                 | PEEK & Quartz Glass        |                             |                             | PPS, stainless steel 316L |                            |                           |                            |
| Calibrated Media                 | Water                      | Water                       | IPA                         | Water                     | IPA                        | Water                     | IPA                        |
| Range                            | 0±1.5 μL /min              | 0±7 μL /min                 | 0±70 μL /min                | 0±2 mL /min               | 0±2 mL /min                | 0±40 mL /min              | 0±40 mL /min               |
| Accuracy (measured value)        | 10% mv above<br>75 nL /min | 5% mv above<br>0.42 μL /min | 20% mv above<br>4.2 μL /min | 5% mv above<br>10 μL /min | 10% mv above<br>50 μL /min | 5% mv above<br>10 μL /min | 10% mv above<br>50 μL /min |
| Lowest detectable flow increment | 3.7 nL/min                 | 10 nL/min                   |                             | 0.25 μL/min               |                            | 25 μL/min                 |                            |

| FLOW UNIT                        | XL                          | L                         |                           | M                        |                          |
|----------------------------------|-----------------------------|---------------------------|---------------------------|--------------------------|--------------------------|
| Product Number                   | FLU-XL                      | FLU-L-D                   |                           | FLU-M-D                  |                          |
| Sensor Inner Diameter            | 1.8 mm                      | 1.0 mm                    |                           | 430 μm                   |                          |
| Maximum Pressure                 | 5 bar                       | 12 bar                    |                           | 100 bar                  |                          |
| Wetted Materials                 | PEEK and Borosilicate Glass |                           |                           |                          |                          |
| Calibrated Media                 | Water                       | Water                     | IPA                       | Water                    | IPA                      |
| Range                            | 0±5mL/min                   | 0±1mL/min                 | 0±10mL/min                | 0±80μL/min               | 0±500μL/min              |
| Accuracy (measured value)        | 5% m.v. above 0.2 mL/min    | 5% m.v. above 0.04 mL/min | 20% m.v. above 0.5 mL/min | 5% m.v. above 2.4 μL/min | 20% m.v. above 25 μL/min |
| Lowest detectable flow increment | 3 μL/min                    | 0.7 μL/min                |                           | 0.06 μL/min              |                          |

Measured Values from 5% to 100% of product range in normal conditions

## Related products

### FLOWBOARD FLB

The Flowboard is a hub that manages communication between Fluigent Software and up to eight FLOW UNITS.

Required with MFCS™ or stand alone

### FLOW UNIT KITS

Tubing & fitting elements dedicated for each range of FLOW UNIT, allow for fast handling and optimized performance.

Required

### LineUp™ Push-Pull ELUPPU1000


Pressure & vacuum controller used to pressurize the sample reservoirs and drive fluids to a microfluidic set-up. Allows for control and monitor in flow rate if combined with a FLOW UNIT.

Required to monitor fluids bidirectionally without a PC

# FLOWBOARD

## Flow Sensor Hub

P/N : FLB



The **Flowboard** is a hub that manages communication between **Fluigent Software** and **up to eight FLOW UNITS**. The **Flow Rate Platform** is designed to be used with any flow control system.


This unique flow rate measurement system provides the **best precision for various flow rate ranges**. Used with any Fluigent pressure system, the **Direct Flow Control algorithm** and Fluigent Software allow for **accurate** and **fast regulation** of the fluids by automatically adjusting pressure to maintain the set flow rate.

*Easy to set  
and use*

*Up to  
8 flow sensors*

*Mix range from  
nL/min to  
mL/min*

*Adaptable :  
work with any  
flow controller*



## Related products

### FLOW UNIT

The FLOW UNIT is a bidirectional flow sensor used for direct flow control.

Required

### MFCST™ series

The MFCST™ series products are designed to control pressure. Combined in the set-up with a Flowboard and FLOW UNIT sensors, it allows for direct flow control.

Required for  
pressure-based  
flow control



# PRESSURE UNIT

## Microfluidic In-Line Pressure Sensor



The **PRESSURE UNIT** is a stand-alone sensor for **continuous, accurate measurement** of the pressure in a fluidic path. The sensors can detect values over the range of **-1000 mbar** (-15 psi) to **7000 mbar** (100 psi). They are directly connected to a PC via USB, and display measurement in real time with **Fluigent Software** interface. User's can output this value for custom software applications using the **Software Development Kit**.

*Compact device:  
dedicated for  
benchtop use*


*Wide range of  
detection:  
from -1000 to  
7000 mbar*

*Ease of use:  
Operate within  
a minute*

*Plug & Play:  
No hub required,  
connect it  
directly to the PC*

## Characteristics

| PRESSURE UNIT  | S                | M                  | L                  |
|----------------|------------------|--------------------|--------------------|
| Product Number | EIPS345          | EIPS1000           | EIPS7000           |
| Pressure range | -345 to 345 mbar | -1000 to 1000 mbar | -1000 to 7000 mbar |



*Microfluidic set-up with in-line pressure detection*

## Related products

### PRESSURE UNIT KIT

EIPSKIT

Tubing & fitting elements dedicated for all range of PRESSURE UNIT, allow for fast handling and optimized performance. **Required**

### LineUp™ Push-Pull

ELUPPU1000

Pressure & vacuum controller used to pressurize the sample reservoirs and drive the fluids toward microfluidic set-up. The PRESSURE UNIT allows for pressure detection and control anywhere it is placed in the system.

**Required to monitor fluids and pressure**

### OxyGEN

SSFT-OXY

Fluigent Software allows for graphical display in real-time of the measured value by the sensor.

**Required**

# MICROFLUIDIC VALVES

Easy Switch Solutions™



**M-SWITCH™**

11-port  
10-position



**Switch EZ**

microfluidic valve  
controller



**L-SWITCH™**

6-port  
2-position



**2-SWITCH™**

3-port  
2-way

The **ESS™** is a versatile **fluid handling platform** for directing fluid flow. It can be **automated** using Fluigent software or **controlled locally** without the need of a PC. The product line offers **three different valves** for a wide range of applications. One valve controller can be combined with the pressure system and Fluigent Software Solutions or used in stand alone mode.

*Econominal:  
reduce reagent  
consumption*

*Accuracy:  
fast valve  
actuation and  
low internal  
volume*

*Versatility:  
Control in real-  
time, Automate  
and Integrate  
as needed*

*Compact:  
dedicated for  
benchtop use &  
useable without  
a PC*

# 2-SWITCH™

P/N : 2SW002


3-port/2-way microfluidic sampling valve



The **2-SWITCH™** is a compact **3-port/2-way bidirectional** microfluidic valve. Using standard connectors it can be integrated to any microfluidic setup. Its unique design allows the combination of multiple units together to save space on your benchtop.

## Characteristics

The **2-SWITCH™** may be used as a manually-operated, stand-alone device or controlled by Fluigent Software for long-term experiments. Its versatility makes it ideal for applications where fluid sorting, switching or periodic sampling are required.



## Related products

**SWITCH EZ**  
ELUSEZ

Microfluidic valve controller that can host up to six 2-SWITCH™ and allow for local actuation without a PC or time-based automation if connected to Fluigent Software with a LINK.

Required

**2-SWITCH™**  
**KIT**  
CTQ-KIT-2SW2

Tubing & fitting elements dedicated the 2-SWITCH™, allows for fast handling and optimized performance.

Required

# M-SWITCH™

P/N : ESSMSW003


11-port/10-position microfluidic bidirectional valve



The **M-SWITCH™** is a bidirectional 11-port / 10-way microfluidic valve for injection or selection of **up to 10 different fluids or chips**. The **flow is bidirectional** in the valve, meaning that the device can be used as a **selector** and a **distributor** for either multiplexing or demultiplexing purposes.

## Characteristics

The **M-SWITCH™** is compatible with Fluigent's pressure controllers **MFCSTM series** and **LineUp™ series**. It is controlled with **LineUp™ SWITCH EZ** allowing for **local control** without a PC or **automated valve actuation** using Fluigent Software.



## Related products

### SWITCH EZ ELUSEZ

Microfluidic valve controller that can host up to three M-SWITCH™ and allow for local actuation without a PC or time-based automation if connected to Fluigent Software with a LINK.

Required

### M-SWITCH™ KIT CTQ-KIT-ESSMSW003

Tubing & fitting elements dedicated the M-SWITCH™, allows for fast handling and optimized performance.

Required

### PRESSURE MANIFOLD CTQ-MANI

The Pressure Manifold is designed to redirect pressure flow from a single controller into up to 10 reservoirs.

Optional

# L-SWITCH™

P/N : LSW001


6-port/2-position microfluidic injection / recirculation valve



The **L-SWITCH™** is a bidirectional 6-port/2-position valve for **precise sample injection** or switching between different fluids. Its configuration make it ideal for **recirculation** in cell culture.

## Characteristics

The **L-SWITCH™** can be used as a useful cell culture tool: a small volume of buffer can be recirculated within a closed loop into the chip for several hours or days. Combined with our **MFCST™-EZ** or **LineUp™ series** it can provide stable flow with a minimal impact on shear stress. The valve also enables one to load and inject a precise volume of fluid. By selecting the position of the **L-SWITCH™** one will choose when to load the fluid then inject it (using a sample loop). Several sample loops are available from 5 µL to 100 µL.



## Related products

**SWITCH EZ**  
ELUSEZ

Microfluidic valve controller that can host up to three L-SWITCH™ and allow for local actuation without a PC or time-based automation if connected to Fluigent Software with a LINK.

Required

**L-SWITCH™ KIT**  
CTO-KIT-LSW

Tubing & fitting elements dedicated the L-SWITCH™, allows for fast handling and optimized performance.

Required



# FLPG+

## Microfluidic Low Pressure Generator



The **Fluigent Low Pressure Generator** is the perfect tool for those who need a pressure source with all accessories included and integrated. This tool is adapted to the **LineUp™ series** and **MFCST™ series** or any other microfluidic pressure-based instrument.

### Characteristics

The FLPG+ compressor supplies **up to plus 2 bar**. This item contains the pressure source, a manual regulator, and a pressure sensor and display. **Three FLPG models** are available depending on requirements

| FLPG+ model   | Product Number |
|---|----------------|
| FLPG+ High Pressure Supply - Up to 8 channels (2 bar) or 16 channels (1 bar)  | FLPG003        |
| FLPG+ Silent Pressure Supply - Up to 4 channels (2 bar) or 8 channels (1 bar) | FLPG005        |
| FLPG+ with incubator aspiration option  | FLPG004        |

### Related products

|                       |   |  |
|-----------------------|---|--|
| <b>MFCST™ series</b>  | The MFCST™ series products are designed to control pressure. If the device has no integrated pump, it needs to be supplied by the FLPG+ compressor. | <b>Required for pressure-based flow control</b>  |
| <b>LineUp™ series</b> | The LineUp™ series products are designed to control pressure. The system requires an external pressure source to be supplied.                       | <b>Required for pressure-based flow control.</b> |

# COMPACT PRESSURE SOURCE



Product Number

E-AC-RX1-2500

The **Compact Pressure Source** is the latest addition to the Fluigent Pressure Source family of products. It also exists in an OEM version, called the **RX**. This compact compressor is designed to be used as a **standalone pressure source**. Packaged in a robust steel enclosure, it provides dried and filtered air at up to **2500 mbar** for one or pressure control systems such as the **LineUp™ series**, **MFCST™ series**, **PX-series**, or other instruments which need compressed air to operate.

*Compact  
and versatile*

*Provides  
condensation-proof  
and filtered air*

*Suited for  
OEM integration*

## Related products

### LineUp™ series

The LineUp™ series products are designed to control pressure. The system requires an external pressure source to be supplied. The Compact Pressure Source has been especially designed to supply properly the LineUp™ series and reduce global system size.

**Required for  
pressure-based  
flow control.**

# VACUUM PUMP

## Compact Vacuum Pump



The **Compact Vacuum Pump** generates a vacuum to **supply any negative pressure channel** or perform **aspiration**. The pump is compact and has a **regulation dial** to provide an outstanding ease of use.

Product Number

EACVACPUMP

*High flow rate:  
up to  
7L/min*

*Ease of use  
and  
compactness*

*Low vibration  
and  
noise level*

*High  
compatibility:  
from 0 to  
-1000 mbar*

## Extended products

### MFCS™ series

The MFCS™ series products are designed to control pressure and vacuum. If the device has no integrated vacuum pump and negative pressure channels, it can be supplied by the Compact Vacuum Pump

Recommended

### Push-Pull

The LineUp™ Push-Pull is a pressure and vacuum controller. The system requires an external vacuum source to perform aspiration. The Compact Vacuum Pump can supply properly the LineUp™ series with negative range and Push-Pull.

Recommended

# P-CAP SERIES

## Tube metal caps



The **P-CAP** is an air-tight metal cap that allows for pressurization of standard lab tubes for microfluidic fluid delivery. The reservoirs and caps are available with different **volume sizes** and **pressurization levels**.

*Autoclavable  
& Incubator  
compatibility*


*No contact  
between the  
sample and the  
P-CAP*

*Suitable for  
long term  
experiments*

*Compatible with  
different tube  
sizes*

## Characteristics

| Name                                   | Max. volume of pressurized liquid | Pressure compatibility                             | Product Number |
|--|-----------------------------------|--|----------------|
| P-CAP for Eppendorf plastic reservoirs | 1.5 mL and 2 mL                   | MFCST <sup>TM</sup>                                | P-CAP2-LP      |
|  |                                   | Flow EZ <sup>TM</sup> & MFCST <sup>TM</sup> 7 bars | P-CAP2-HP      |
| P-CAP for Flacon Tubes                 | 15 mL                             | MFCST <sup>TM</sup>                                | P-CAP15-LP     |
|  |                                   | Flow EZ <sup>TM</sup> & MFCST <sup>TM</sup> 7 bars | P-CAP15-HP     |
|  | 50 mL                             | MFCST <sup>TM</sup>                                | P-CAP50-LP     |
|  |                                   | Flow EZ <sup>TM</sup> & MFCST <sup>TM</sup> 7 bars | P-CAP50-HP     |



## Related products

### FLOW EZ<sup>TM</sup>

The Flow EZ<sup>TM</sup> is the most advanced flow controller.

Optional


### MFCST<sup>TM</sup>

The MFCST<sup>TM</sup> series products are designed to control pressure.

Optional

# FLUIWELL SERIES

## Pressurized fluid reservoirs



The **Fluiwell** is a tool for pressurizing samples inside the vials in order for the fluids to flow through the microfluidic system. The reservoirs and caps are available with different **volume sizes**, **channel numbers**, and **pressurization levels**.

*Autoclavable  
compatibility*


*Compatible  
with any tubing*

*Suitable for  
long term  
experiments*

*Compatible with  
different tube  
sizes*

## Characteristics

| Name        | Max. volume of pressurized liquid | Pressure compatibility | Product Number |
|-------------|-----------------------------------|------------------------|----------------|
| Fluiwell-4C | 0.5 mL                            | MFCS™                  | 14000501       |
|             |                                   | Flow EZ™ & MFCS™ 7 bar | 24000501       |
|             | 2 mL                              | MFCS™                  | 14002001       |
| Fluiwell-1C | 15 mL                             | Flow EZ™ & MFCS™ 7 bar | 24002001       |
|             |                                   | MFCS™                  | 11015001       |
|             | 50 mL                             | Flow EZ™ & MFCS™ 7 bar | 21015001       |
|             |                                   | MFCS™                  | 11050001       |
|             |                                   | Flow EZ™ & MFCS™ 7 bar | 21050001       |



Consult our product guide page 92

## Related products

### FLUIWELL KITS

The Fluiwell Kits ensure the air-tightness of the fluid reservoirs.

Required

# BOTTLE-CAP SERIES

## Pressurized bottle cap



The **Bottle-CAP** series are **air-tight microfluidic adapters** that allows the pressurization of **large volumes of liquid** for microfluidic applications. This cap is compatible with bottles with a **GL-45 thread**. It is compatible with Fluigent's pressure controllers: the **MFCST<sup>™</sup> series** and the **LineUp<sup>™</sup> series**. The number of outlets for the cap can be customized, with **up to 3 fluidic outlets**, to refill easily or to perfuse different microfluidic chips at the same time.

*Precaution : Glass laboratory bottles should not be pressurized above 2 bar.*

*Easy to use  
and  
standardized  
fittings*

*Incubator  
compatibility*

*Robust  
and  
autoclavable*

*Possibility  
to refill  
during  
experiment*

## Characteristics

| Bottle-CAP model | Tubing & Fitting kit | Number of fluidic port |
|------------------|----------------------|------------------------|
| RES-CAP          | CTQ-KIT-BC           | 1 fluidic port         |
| RES-CAP-PCK      | included             | 1 fluidic port         |
| RES-CAP-3P       | included             | 2 fluidic ports        |
| RES-CAP-4P       | included             | 3 fluidic ports        |



# STANDALONE ACCESSORIES



## P-CAP SERIES RACK

(P-CAP-RACK)

The **P-CAP rack** is a **compact holder** compatible with any P-CAP size. The versatility of this microfluidic accessory allows to hold **different sizes of pressurized reservoirs**. The rack can be rotated easily if needed.

### Benefits:

- ▶ *Autoclavable*
- ▶ *Incubator compatibility*
- ▶ *Convenient setup saving space on the bench*

### Extended products :

**P-CAP series:** Metal air-tight cap for standard lab reservoirs (*Recommended*)



## EDUCATIONAL MICROSCOPE

(OSEDUCMIC)

A **cost-effective microscope** for microfluidics. The LED display enables **several people to observe** at the same time. LCD microscope 50-500x, 2000 (digital) magnification transmitted and reflected light SD card reader USB port.

### Benefits:

- ▶ *LCD Microscope display - live view, photo, and video*
- ▶ *Magnification: 50x-500x, 2000x*
- ▶ *Optimal LED illumination with color filters*

# PRESSURE MANIFOLD


(CTQ-MANI)

## 10-Position Pressure Manifold

The **Pressure Manifold** is designed to redirect pressure flow into **multiple reservoirs**. It allows splitting a gas line into up to 10 lines.

### Contents:

- 10 way manifold (x1)
- Luer lock HP (x11)
- Red plugs (x10)
- 30 cm high pressure tubings (4 mm OD x 2.5 mm ID) (x10)



# AIR FLOW REGULATION KIT

(10000001)

A **Regulation Kit** for Pressure Sources consisting of a pressure regulator, an air dryer, and pneumatic tubing.

### Contents:

- Pressure regulator (x1)
- Air drier (x1)
- 6 mm OD pneumatic tubing (5 m)
- G1/4" – 6 mm OD tube fitting (x2)
- G1/8" – 6 mm OD fitting (dryer) (x2)



# BUBBLE TRAP

(CTQ-006BT)

The **Bubble Trap** from Fluigent is a device suitable for aqueous stream flows, to **prevent air bubbles** from entering a microfluidic system. Based on a porous membrane inside the trap, the Bubble Trap prevents high instabilities in the flow rate, high shear stress variation, and damage to the microfluidic system.


### Contents:

- Bubble Trap (x1)
- Membranes (x3)
- 1/4-28 Flat-Bottom for 1/16" OD Tubing and ferrules (x2)
- FEP Tubing 1/16" OD x 0.020" (508 µm) ID (1 m)



# DSURF

## Quality Surfactant



**dSURF** is a high-performance fluorosurfactant for microdroplet generation.

It allows for **high-quality droplet formation** and **long-term stability** in conditions such as dPCR and cell culture experiments.


This experiment was realized by using the dSURF and the Droplet Starter Pack that contains our Fluigent microfluidic chip, the EZ Drop.

*Biocompatible  
with mammal  
cell, yeast and  
bacterias*

*High-  
performance  
fluorosurfactant*

*Reliable Results:  
high droplet  
stability*

*Broad range  
of dye  
compatibility*



*dSURF Biological Applications*

**dSURF** is a **high-performance fluorosurfactant** dedicated to microdroplet generation. Be sure to get reliable results! dSURF being a non-ionic fluorosurfactant, it reduces droplet cross-talk and fits properly with any biological application such as dPCR and cell culture.

dSURF is adaptable to most droplet microfluidic applications. It comes in a 2% formula that can be diluted with our fluorinated oil, **dOIL**, to suit all application requirements. For dPCR experiments, dSURF has demonstrated excellent compatibility with FAM™, HEX™, VIC® and EvaGreen® dyes.

## Characteristics

| Name               | Characteristics | Product Number |
|--------------------|-----------------|----------------|
| dSurf<br>3x4 mL    | 2% formula      | DR-RE-SU-12    |
| dSurf<br>30 mL     | 2% formula      | DR-RE-SU-30    |
| Honey<br>dSurf 5 g | 100% formula    | ODSURF5G       |




## dOIL

dOIL is a pure fluorinated oil (3M™ Novec™ 7500 Engineered Fluid) in which our dSURF surfactant is diluted.

Fluorinated oils have shown several advantages compared to other carrier fluids such as mineral oils. They show **better PDMS compatibility** due to minimum swelling. They are also more **adapted to biochemical experiments** due to low organic compound transfer drop to drop, and they have shown better biocompatibility in **long term in droplet cell culture experiments**.

dSURF is a new generation of fluorosurfactant providing highly **reliable droplet production** and stability even under PCR amplification conditions. Combined with the **droplet pack**, our **biocompatible** emulsion stabilizer also enables the generation of monodispersed droplet of any size.



## Reagents

### dOIL

DR-RE-SU-A1

dOIL is pure Novec™ fluorinated oil in which our dSURF emulsion stabilizer is diluted.

Optional

### DROPLET STARTER PACKAGE

DROPPACK-01

The Droplet Starter Pack is designed for microfluidic droplet experiments.

Optional

### EZ DROP

DROPKIT01

The Droplet Kit is designed to be used with the droplet starter pack.

Optional

# MICROFLUIDIC CHIPS

## RAYDROP SINGLE EMULSION DEVICE



### Characteristics

#### RayDrop Single Emulsion

| Capillary size                       | Product Number |
|--------------------------------------|----------------|
| 30 $\mu\text{m}$ - 150 $\mu\text{m}$ | 1DPRD01        |
| 60 $\mu\text{m}$ - 300 $\mu\text{m}$ | ORDRPSE-60-300 |
| 90 $\mu\text{m}$ - 450 $\mu\text{m}$ | ORDRPSE-90-450 |

The **RayDrop** developed and manufactured by Secoya is a microfluidic **droplet generator** composed of three main fully removable parts: two inserts on each side and a center section containing a nozzle and an outlet capillary. There are four standard microfluidic connections, two on the box for the continuous phase, and one on each insert for the dispersed phase entry and the collecting emulsion outlet.

The RayDrop's design allows for **multiple liquid emulsification**, in the same device **without the need for special coatings**.

*One device  
for multiple  
applications*

*Easy  
cleaning*

*No coating  
needed*

*High  
monodispersity*

## Applications

- ▶ PLGA microparticle
- ▶ Polymer microcapsule
- ▶ Alginate microbeads
- ▶ Liposome Nanoparticle
- ▶ Oil-in-Water droplets
- ▶ Water-in-Oil droplets

## RAYDROP DOUBLE EMULSION DEVICE



## Characteristics

### RayDrop Double Emulsion

| Capillary size       | Product Number    |
|----------------------|-------------------|
| 30µm - 70µm - 150µm  | ORDRPDE-30-70-150 |
| 60µm - 120µm - 300µm | O-DE-RDRPC03-EUP  |
| 90µm - 160µm - 450µm | O-DE-RDRPC02-EUP  |

The **RayDrop** droplet generator developed and manufactured by Secoya is a patented device for **droplet generation** with most of the **advantages of glass chips** such as resistance to strong chemicals and **compatibility with high pressures** ( > 2 bar). Glass chips are a very expensive disposable, are prone to leakage issues and have limited lifetimes as they are very difficult to recover once clogged. The RayDrop device uses **standard fittings** leading to **sealed connections** and its design allows for **easy recovery and cleaning** if clogging occurs.

*Perform double emulsions in one single device*

*Droplet size from 70 µm to 150 µm diameter*

*Easy-to-clean exchangeable nozzle*

*High flexibility (w/o/w) o/w/o)*


*Surface coating-free: No need any surfactant*

## Applications

- ▶ **Double emulsions** can be produced in only one step in a single device
- ▶ Produce **water-oil-water** and **oil-water-oil** droplets **without any surface coatings** needed.
- ▶ The system does **not need any surfactant** for droplet formation

# EZ DROP

## Droplet generation chip



## Characteristics

| Microfluidic Chips | Product Number |
|--------------------|----------------|
| EZ Drop            | DROPKIT01      |

The **EZ Drop** is a PDMS chip dedicated for water-in-oil droplet generation. It comes with the tubing and fittings needed for an experiment. It works with any **LineUp™** pressure channel and **FLOW UNITS** models.

|                    |   |  |  |  |
|--------------------|---|--|--|--|
| Up to<br>10 000 Hz | From 20 $\mu\text{m}$<br>to 100 $\mu\text{m}$<br>droplets | 20 $\mu\text{m}$ , 50<br>$\mu\text{m}$ , and 100<br>$\mu\text{m}$ markers<br>on the PDMS<br>Droplet Chip | Microscope<br>slide<br>dimensions.<br>Easy<br>connections. | Quality<br>monitoring<br>with a QR<br>code |
|--------------------|---|--|--|--|

## Contents

- **EZ Drop chip**, with 3 designs each (x3 chips, 9 designs)
- **Tubing** (250 ID ; 1/32"OD) (2 m)
- **Sleeves** (x2)



# DROP-SEQ

Drop-seq chip



## Characteristics

| Microfluidic Chips | Product Number |
|--------------------|----------------|
| FlowJEM Drop-Seq   | ODROPSEQCHIP   |


Fluigent has worked with FlowJEM to provide the best microfluidic chips for **Drop-seq**. Drop-sequencing (Drop-seq) developed by the McCarroll lab at Harvard Medical School, is a method designed for the **parallel analysis of mRNA expression** in thousands of individual cells following their **encapsulation in tiny droplets**. These droplets (nanolitre scale) are formed by **precisely combining** aqueous and oil flows in a specially designed microfluidic device (Drop-seq chip). Expression profiling can then be carried out in **tens of thousands of cells in a matter of hours**.

## Benefits

- ▶ **Latest Design:** Each Droplet Generation Device is based on the design recommended in the latest McCarroll lab Drop-seq protocol.
- ▶ **Precision Engineered:** Robust Devices durable over a wide range of pressures, temperatures, and flow rates.
- ▶ **26 Droplet Generation Devices Per Chip:** Provides value for money in a chip which lasts. When the life of one device is depleted, simply move onto the next one.
- ▶ **Produces Highly Mono-Dispersed Droplets:** Reliable and consistent generation of droplets of optimal size for Drop-sequencing.
- ▶ **Efficient Production Of Transcript Libraries:** Superior design promotes optimal mixing of component fluids, minimizing bead shearing or premature lysis of cells and mRNA release.



# CELL CULTURE CHIPS



## BE-FLOW (OOC-FLOW-01)

Most easy-to-use device dedicated to cell culture under flow

**BE-Flow** is Beonchip's simplest device dedicated to **cell culture under flow**. It allows the performance of long-term **2D or 3D culture** in two independent channels. BE-Flow is compatible with any microfluidic pump system and can be used with a rocker as its fluid reservoirs are situated by the inlets/outlets. This is an optimal device for **vascular research** where shear stress plays a major role in gene expression.

*Easy  
to use*

*Easy  
to connect*


*No unspecific  
absorption*

*Impermeability*

*Cell recovery  
option*

## BE-DOUBLE FLOW (OOC-DBLE-08)

Explore the interaction between different 2D and 3D cultures




**BE-Doubleflow** is Beonchip's most advanced device. It consists of two perfusable channels connected via a **porous membrane**. It allows for the investigation of different 2D and 3D cultures in a **biomimetic environment** and controls the efficiency of the interaction by selecting the optimal pore size for your application. This is the optimal device when a **hypoxic environment** is needed for studying the **effect of circulating particles** (bacteria, immune system, circulating tumor cells) and for **endothelium/epithelium barrier** when no ALI is needed or when flux plays a role in both sides of the coculture.

*Easy  
to use*

*Easy  
to connect*

*No unspecific  
absorption*

*Cell recovery  
option*



## BE-TRANSFLOW (OOC-TRNS-07)

The most versatile cell culture platform

**Be-Transflow** is Beonchip's most **versatile cell culture platform**. It allows the study of **complex culture configurations** by joining a culture well with a microfluidic channel via a porous membrane. This is the optimal device for **Air Liquid Interface (ALI)** culture, **endothelium/epithelium barrier**, and **crosstalk studies**.


*Easy  
to use*

*Easy  
to connect*

*No unspecific  
absorption*

*Impermeability*

*Cell recovery  
option*



## BE-GRADIENT (OOC-GRAD-05)

Designed for the application of electrochemical gradients to 3D cell cultures

**BE-Gradient** is Beonchip's device for the application of **chemical gradients to 3D cell cultures**. BE-Gradient is compatible with any type of optical microscopy (inverted phase contrast, confocal, fluorescence..). BE-Gradient consists of a central chamber for cell culture and two lateral channels connecting to the central chamber through 3 small micro-channels. The lateral channels are meant to **simulate blood vessels**. 2D culture is also possible for adherent cells not only in the central chamber but also in the lateral channels.

*Easy  
to use*

*Easy  
to connect*

*No unspecific  
absorption*

*Cell recovery  
option*

# OXYGEN

P/N : SSFT-OXY

Control in real time, Automate protocols and Record experimental data



The new way to get full control of your microfluidic setup.

OxyGEN is a **single interface**, with **plug and play capabilities**, available for **common desktop OS (windows, MacOS, Linux)**, that allows you to **control, monitor and automate** all Fluigent products. It combines in one program all the functions and capabilities of our traditional software: A-i-O, MAT, ESS control and much more.

Through its **intuitive dashboard**, OxyGEN is our new reference tool for **real-time control** and for developing **time based protocols** focusing on **pressures, flow rates, volumes, and valve control** in microfluidic experiments.

*Control in  
real time  
pressures, flow  
rates, valves  
and volumes*

*Edit and  
automate  
long time  
protocols*

*Record  
and  
export data*

*Plug & Play  
connection  
& Simulated  
instruments*

# SOFTWARE DEVELOPMENT KIT

P/N : FLUIGENT-SDK

## Custom Software Development



Custom your own software application using the SDK to automate and monitor Fluigent devices

**Fluigent Software Development Kit (SDK)** allows one to fully **integrate Fluigent devices** in a customized application. It has been designed in **several languages**, among the most popular ones in the instrumentation field (LabVIEW, C++, C#.NET, Python, and MATLAB).

The SDK controls all Fluigent pressure and sensor instruments as well as microfluidic valves and advanced regulation loops. One can still use independent SDK for **basic hardware set-ups** or for **specific software requirements**.


## Languages

- ▶ **Python** (ver. 3.1+ minimum)
- ▶ **LabVIEW** (ver. 2016+ minimum)
- ▶ **MATLAB** (ver. R2015a+ minimum)
- ▶ **C++** (ver. 11+ minimum)
- ▶ **C#** (ver. .NET Core 3.1 minimum)

## Operating system

- ▶ **Windows 10, 8 and 7** (32 Bits, 64 Bits)
- ▶ **MacOS**
- ▶ **Linux**

# PRESSURE SYSTEM KITS




## MFCSTM HIGH PRESSURE KIT

(CTQ-KIT-HP-MFCS)

Tubing and connecting kit designed to be used with MFCSTM-EZ for the 7 bar pressure range.

**Kit contents:**

- Transparent pneumatic tubing 4 mm (0.2 m)
- Backflow filters (x4)
- Blue pneumatic tubings 4 mm (x4)
- Male luer black tubings (x4)



## MFCSTM LOW PRESSURE KIT

(CTQ-KIT-LP-MFCS)

Tubing & Fitting Kit compatible with low pressure MFCSTM-EZ – channels from 25 mbar to 2000 mbar ranges and from -800 mbar to -25 mbar ranges.

**Kit contents:**

- Male white luer fittings, 1.6 mm (x4)
- Luer caps (x4)
- Backflow filters (x4)
- 1x3 mm tubing (2 m)




## LineUp™ SUPPLY KIT

(LU-SPK-0002)

The **LineUp™ Supply Kit** contains all the components to provide power and pressure supply to any LineUp™ series module.

### Supply kit contents:

- Power supply 24 V 1,75 A (x1)
- Blue pneumatic tubing 6 mm (x1)
- Female Staubli tubing fitting (x1)
- Power cable – EUR – Little (x1)



## LineUp™ CHAIN-TO-CHAIN KIT

(LU-C2C-0001)

LineUp™ series connection cable

The **Chain-2-Chain Kit** allows to connect **LineUp™ series** modules together without combining them. The cable allows for data and power transmission while the tubing allows pressure to flow from one module to the next. Both the cable and the tubing are approximatively 1.8 m. One can use different gases to pressurize the reservoirs and use **negative and positive modules at the same time**.



## LineUp™ P-SWITCH KIT

(ELUPSWKIT)

Tubing & fitting kit for the **LineUp™ P-SWITCH**.

### Kit contents:

- Pneumatic tubing 3 mm OD (2 m)
- Transparent tubing 4 mm OD (1 m)
- Red plugs 4 mm (x2)
- Red plugs 3 mm (x8)

# SENSOR KITS



## PRESSURE UNIT KIT

(EIPSKIT)

Tubing & Fitting Kit dedicated to the **PRESSURE UNIT** product line.

**Kit contents:**

- Flangeless fittings 1/4-28 (x2)
- Ferrules (x4)
- 1 m x FEP tubing 1/16" OD 0.020" ID (1 m)



## FLOW UNIT XS KIT

(CTQ-KIT-LQ-XS)

A kit containing all the tubing and fittings for our **FLOW UNIT XS**.

**Kit contents:**

- Adapters 10/32" (x2)
- Blue PEEK tubing 1/32x0.010x100ft (LQ) (1 m)
- Sleeves 1/16" → 1/32" (x2)
- Flowmeter fittings (x2)
- 1/16" → 1/32" tubing fittings (x1)
- PEEK filter XS (x1)




## FLOW UNIT S AND M 1/32" KIT

(CTQ-KIT-LQ)

A kit containing all the tubing & fitting elements for the **FLOW UNIT S and M** with 1/32" outer diameter (OD) tubing.

**Kit contents:**

- Flowmeter fittings (x2)
- Sleeves 1/16" → 1/32" (x1)
- Blue PEEK Tubing 1/32x0.010x100ft (1 m)
- Red tubing reducer 1/16" → 1/32" (x1)




## FLOW UNIT S AND M 1/16" KIT

(CTQ-KIT-FU2)

The most recent kit containing all the tubing & fittings for the **FLOW UNIT S and M**. These new fittings allow using **1/16" outer diameter (OD) tubing** with fewer connecting elements, making it the best option for 1/16" tubing.

### Kit contents:

- FLOW UNIT S and M adapters (x2)
- FEP tubing 1/16" OD 0.020" ID (1 m)
- Flangeless fittings 1/4-28 (x2)
- Blue ferrules (x2)



## FLOW UNIT L KIT

(CTQ-KIT-HQ)

A kit containing all the tubing and fittings for the **FLOW UNIT L**.

### Kit contents:

- FEP tubing 1/16" OD 0.020" ID (1.00 m)
- Flangeless fittings 1/4-28 (x2)
- Blue ferrules (x4)



## FLOW UNIT XL KIT

(CTQ-KIT-XL)

Kits containing the tubing and fittings to use the **FLOW UNIT XL**.

### Kit contents:

- FEP tubing 1/16" OD 0.020" ID (1 m)
- Flangeless fittings 1/4-28 (x2)
- Low-pressure union (x1)
- Tube PEEK (x0.1 m)
- Blue ferrules (x4)



# MICROFLUIDIC VALVE KITS



## 2-SWITCH™ KIT

(CTQ-KIT-2SW2)

Tubing & fitting kit for the **2-SWITCH™** microfluidic valve.

**Kit contents:**

- FEP tubing 1/16" OD 0.010" ID (3 m)
- PEEK plugs 1/4-28 black (x2)
- Flangeless fittings 1/4-28 (x6)
- Blue ferrules (x12)
- Fittings 4 mm (x1)
- P-Y tube (x1)



## M-SWITCH™ KIT

(CTQ-KIT-ESSMSW003)

Kit consisting of tubing & fitting elements for the **M-SWITCH™** microfluidic valve.

**Kit contents:**

- FEP tubing 1/16" OD 0.10" ID (2.5 m)
- Sleeves 1/16" to 1/32" (x6)
- PEEK plugs 1/4-28 black (x11)
- Flangeless fittings 1/4-28 (x11)
- Blue ferrules (x22)



## L-SWITCH™ KIT


(CTQ-KIT-LSW)

Tubing & fitting kit for the **L-SWITCH™** microfluidic valve.

**Kit contents:**

- 10/32" adapters (x7)
- FEP tubing 1/16" OD 0.020" ID (1 m)
- Acetal blue plugs 10/32" (x2)
- 10-32 female – female luer adapter (x1)

# SAMPLE RESERVOIRS KITS



## P-CAP 2 ML KIT

(CTQ-KIT-PCAP2)

Tubing & Fitting Kit for **2mL P-CAP** reservoirs.

**Kit contents:**

- FEP tubing 1/16" OD 0.020" ID (0.5 m)
- Flangeless fittings 1/4-28 (x2)
- Blue ferrules (x4)
- O-ring P-CAP 2 mL (x1)



## P-CAP 15 ML KIT

(CTQ-KIT-PC15)

Tubing & fitting kit for **15 mL P-CAP** reservoirs.

**Kit contents:**

- Tubing FEP 1/16" OD 0.020" ID (1 m)
- Flangeless fittings 1/4-28 (x4)
- Blue ferrules (x4)
- O-ring P-CAP15 (x1)




## P-CAP 50 ML KIT

(CTQ-KIT-PC50)

Tubing & fitting kit for **50 mL P-CAP** reservoirs.

**Kit contents:**

- Tube FEP 1/16" OD 0.020" ID (1.50 m)
- Flangeless fittings 1/4-28 (x3)
- Blue ferrules (x6)
- O-ring P-CAP50 (x1)
- Black caps PEEK 1/4-28 (x2)




## FLUIWELL 4C KIT

(CTQ-KIT-F4C)

Tubing & fitting kit for **Fluiwell 4-Channels** (0.5 or 2 mL reservoirs).

**Kit contents:**

- Adaptors 10/32" (x4)
- FEP tubing 1/16" OD 0.020" ID (m) (2 m)
- Fluiwell 4-Channels O-rings (x4)



## FLUIWELL 15 ML KIT

(CTQ-KIT-F1C15)

Tubing & fitting kit for **Fluiwell-1C 15 mL** reservoir.

**Kit contents:**

- Adaptors 10/32" (x2)
- Tubing FEP 1/16" OD 0.020" ID (1 m)
- O-ring Fluiwell 1C-15 mL (x1)




## FLUIWELL 50 ML KIT

(CTQ-KIT-F1C50)

Tubing & fitting kit for **Fluiwell 1C-50 mL** reservoirs.

**Kit contents:**

- Adaptors 10/32" (x2)
- Tubing FEP 1/16" OD 0.020" ID (1 m)
- O-ring Fluiwell 1C-50 mL (x1)



## BOTTLE-CAP KIT

(CTQ-KIT-BC)

A kit containing tubing and fittings for the **Bottle-CAP 2 ports**.

**Kit contents:**

- FEP tubing 1/16" OD 0.020" ID (1 m)
- Flangeless fitting 1/4-28 + blue ferrules (x2)
- Male luer locks (white) to barb 1.6 mm (x1)
- Blue ferrules (x2)
- Black Male Luer tubing (x1)

3  
—

# SYSTEMS & PLATFORMS



# Automated Organ-on-chip Platform

Product Number

O-OMISC-PTF

Omi is an automated platform that helps reproduce the microphysiological behavior of organs inside microfluidic chips. It is compatible any type of chips to sustain different cell culture types or organ on chip models (Gut , Skin...)



Omi, automated organ-on-chip platform, can perform any perfusion protocol. It has the ability to customize and automate any protocol that includes simple perfusion, recirculation, sampling and injection. It meets the needs of both beginners in cell culture research and advanced Organ-on-Chip researchers who seek automation and reproducibility.

|  |                                |                                       |   |
|--|--------------------------------|---------------------------------------|---|
| <i>Cell culture under sterile conditions</i> | <i>Versatile Fits any chip</i> | <i>Automated and fully integrated</i> | <i>Intuitive interface and technology</i> |
|--|--------------------------------|---------------------------------------|---|

|           |  |
|-----------|--|
| Protocols | <b>Perfusion</b><br>Perfuse up to 4 mL of many types of liquids such as cell culture media for a long period of time in a controlled and reproducible way. |
|           | <b>Recirculation</b><br>Unidirectional recirculation of medium while maintaining a constant flow rate for a long period of time.                           |
|           | <b>Sampling</b><br>Sample up to 4 mL of soluble factors secreted in the culture medium for analysis or imagery.  |
|           | <b>Injection</b><br>Inject up to 4 mL of fresh medium during the experiments.  |
|           |  |

# Characteristics

## Performance

|                         |                                |
|-------------------------|--------------------------------|
| Flow rate control       | From 1 $\mu$ L/min to 1 mL/min |
| Maximum pressure        | 600 mbar at maximum            |
| Fluid reservoir volume  | 4 mL                           |
| Perfusion volume        | From 1.5 mL to 4 mL            |
| Min volume recirculated | 2.5 mL                         |

## Applications

### Single channel

Blood vessel on chip

### Dual channel

Blood/brain barrier

Gut-on-a-chip

# Features

## Fully integrated platform

Omi can perform user customized protocol including perfusion, drug delivery, sampling and recirculation.

## Compact and transportable

The device is designed to fit in incubators and can be easily transported to a microscope for imaging or cell analysis while maintaining perfusion.

## Remote control

Set up and monitor your protocols on Omi's application's user friendly interface with WIFI.

## Contamination-free

Inject cell or fresh medium, deliver test compounds, collect samples without disconnecting the microfluidic chip and disturbing the shear stress conditions of the experiments.

## Completely autonomous platform

2-hour battery allows for transport from incubator to imaging system

# Platform contents

- Omi device (x1)
- Power supply kit (x1)
- Tablet (x1)
- Cartridges (x3)
- Pneumatic connector kit (x1)
- Fluidic connector kit (x1)
- Low resistance adaptors (x3)
- High resistance adaptors (x3)

# ARIA, AUTOMATED SEQUENTIAL INJECTION SYSTEM

Unique solution for automating perfusion imaging studies.

**You do the science, Aria does the rest**

**ARIA** is a compact instrument to **automate multiple fluid perfusions**. It allows the user to set up a **custom time schedule** for exposing cells, nucleic acids, etc. to antibodies, chromophores, test compounds, or other liquids.



Reagents can be delivered to a perfusion chamber, microfluidic device, or other system for imaging or analysis. **Up to 10 unique solutions** can be delivered at precise times for **faster** and more **reproducible results**.


*Easy to use:  
protocol  
automation*

*Easy to handle:  
intuitive user  
interface*

*Intelligent software  
adapted to any  
application*

*Compatible  
with biological  
applications*

## ARIA Technology



Stop flow, stable perfusion and sequential injections are all functions controlled by ARIA. These functions can be combined using the Aria software to automate protocols such as live cell based assays, micro dosing, cell perfusion, immunolabeling, periodic injections or calcium imaging.

## Characteristics

| Product & Service   | Product Number | Contents   |
|---|----------------|--|
| Aria L (40 $\mu$ L – 1 mL/min) single output  | OAR-L-02-2SW   | Aria 2-SWITCH™ + FLOW UNIT L   |
| Aria L (40 $\mu$ L – 1 mL/min) single output<br>+ Installation + 1 Year extra warranty    | OAR-L-02-2SW-W | Aria 2-SWITCH™ + FLOW UNIT L<br>+ Installation + 1 year Warranty   |
| Aria M (3.2 $\mu$ L – 80 $\mu$ L) single output   | OAR-M-02-2SW   | Aria 2-SWITCH™ + FLOW UNIT M   |
| Aria M (3.2 $\mu$ L – 80 $\mu$ L) single output<br>+ Installation + 1 Year extra warranty | OAR-M-02-2SW-W | Aria 2-SWITCH™ + FLOW UNIT M<br>+ Installation + 1 year Warranty   |
| Aria L (40 $\mu$ L – 1 mL/min) serial output  | OAR-L-02-MSW   | Aria M-SWITCH™ + FLOW UNIT L   |
| Aria L (40 $\mu$ L – 1 mL/min) serial output<br>+ Installation + 1 Year extra warranty    | OAR-L-02-MSW-W | Aria M-SWITCH™ + FLOW UNIT L<br>+ Installation + 1 year Warranty   |
| Aria M (3.2 $\mu$ L – 80 $\mu$ L) serial output   | OAR-M-02-MSW   | Aria M-SWITCH™ + FLOW UNIT M   |
| Aria M (3.2 $\mu$ L – 80 $\mu$ L) serial output<br>+ Installation + 1 Year extra warranty | OAR-M-02-MSW-W | Aria M-SWITCH™ + FLOW UNIT M<br>+ Installation + 1 year Warranty   |
| Aria Single Output Kit  | 2SW-KIT-AR     | Flangeless fittings ¼-28 (x6) + Ferrules<br>blue (x12) + Tube FEP (2 m) + F 120 x2   |
| Aria Serial Output Kit  | MSW-KIT-AR     | Flangeless fittings ¼-28 (x15)<br>+ Ferrules blue (x30) + Tube FEP (4 m)<br>+ Plug black ¼-28 ref P309 (x15) + F<br>120 (x2) |
| Aria Replacement Tubbing Kit  | CTQ-KIT-AR     | Aria internal replacement tubing kit   |

## Focus on the science, not on the setup



**Aria Automation Software** provides quick and easy navigation to control the experiment and walk away while the experiment is running.

Create custom  
injection/  
perfusion  
sequences


Set injection  
parameters  
based on time  
or volume

Set the desired  
flow rate value  
for each  
injection

Program  
stops flow  
and  
incubation



# MICROFLUIDIC COMPLEX EMULSION PLATFORM



**Secoya**  
FLUIDIFY PHARMA

The emulsion platform, developed by Secoya, with our flow control equipment, is a fast and easy screening system to perform emulsification processes such as simple emulsions and double emulsions.

Save time with integrated, organized, ready to use platform and get monodispersed complex emulsions rapidly.

## Applications

- w/o and o/w droplets
- w/o/w and o/w/o emulsions
- UV polymerised resin microparticles
- UV polymerised resin microcapsules
- PLGA microparticles
- PLGA microcapsules
- Alginate microparticles
- Chitosan microcapsules

**Product Number**

O-DE-STD-PTF

## Platform composition

- LineUp Flow EZ™ 7 bar (x3)
- Power supply kit & USB cable (x1)
- LineUp™ LINK (x1)
- FLOW UNIT M (x2)
- FLOW UNIT L (x1)
- 15 mL P-CAP (x4)
- 50 mL P-CAP (x1)
- RayDrop Double Emulsion (x1)
- Fluidic tubing & fitting kit (x1)
- Pneumatic tubing kit (x1)
- Optical system (x1)
- (Optional) protective hood (x1)


|                                     |   |  |
|-------------------------------------|---|--|
| <i>High monodispersity (2%)</i>     | <i>Double emulsion size: 25 µm – 500 µm</i>   | <i>High Frequency: 5 000 Hz</i>            |
| <i>Adjustable shell thickness</i>   | <i>No surface coating</i>                     | <i>w/o/w and o/w/o double emulsion</i>     |
| <i>Ready-to-use robust platform</i> | <i>Single &amp; Double emulsion available</i> | <i>Adaptable for various microcapsules</i> |

## Performance

| Feature                         | Batch method                                      | Fluigent microfluidic method      |
|---------------------------------|---|-----------------------------------|
| Process                         | Multiple step process                             | Direct double emulsion production |
| Homogeneity                     | Random distribution of single and double emulsion | High homogeneity                  |
| Particle size distribution      | More than 50%                                     | ~2%                               |
| Reproducibility                 | Low   | High                              |
| Live particle size control      | No  | Precise                           |
| Continuous / In-line production | No  | Yes                               |

# CELL ENCAPSULATION PLATFORM

**Secoya**  
FLUIDIFY PHARMA



The cell encapsulation platform, developed and manufactured by Secoya and constituted by Fluigent's flow control equipment and Secoya's Emulsion Technology, complete system for high throughput encapsulation of complex and individual cells within highly monodisperse double emulsion droplets small enough ( $<90\ \mu\text{m}$ ) for further analysis.

## Product Number

O-FACS1-PTF

O-FACS2-PTF

O-FACS3-PTF

O-FACS4-PTF

This fast and easy system allows the encapsulation of cells in aqueous core (such as medium, PBS) and oil shell (such as HFE 7500) double emulsions, providing a powerful tool for biochemical and cellular assays, as it enables the isolation of each cell within microreactors, highlighting its characteristics and concentrating the signals to measurable levels to obtain meaningful biological data.

4  
—

# MICROFLUIDIC PACKS

# MICROFLUIDIC APPLICATION PACKS

Packages for cell biology, droplet & particle generation applications

We built **microfluidic application packages** with pressure/flow controllers, microfluidic chips, and valves fit for your application. In this section, you will find packages to create a mechanical stimulation of your cells, generate a unidirectional and stable flow rate in your cell culture chamber, create highly monodisperse simple or double emulsions, and sort your cells and particles according to their size.

*Easy to use  
with simple  
connections*

*Flexible  
packages and  
instruments*

*Contamination-free  
and designed for  
biology*

## MECHANICAL STIMULATION PACK

### A versatile package for mechanical stimulation of cells

This package is the perfect setup to create mechanical stimulation in your experiments with your own microfluidic chip.

#### Package contents:

- LineUp™ LINK
- Flow EZ 1000 mbar
- PCAP 15 mL (x2)
- LineUp™ Push-Pull (x3)
- FLOW UNIT M
- LineUp™ supply kit
- Tubing & fitting kits



# IMMUNOSTAINING PACK

(O-MULTILABEL-PCK)

**A ready-to-use package for sequential injection of multiple fluids for cell or tissue staining.** The Immunostaining Package provides stable and automated delivery of up to 10 different fluids to a single cell culture chamber.

## Package contents:

- LineUp™ supply kit
- LineUp™ LINK
- Flow EZ 1000 mbar
- PCAP 15 mL (x12)
- FLOW UNIT M
- Tubing & fitting kits
- SWITCH EZ
- M-SWITCH™
- 2-SWITCH™
- 10 position pressure manifold



# RECIRCULATION PACK

(O-RECIR-PCK)

**A ready-to-use package for unidirectional and stable fluid recirculation.** The Recirculation Package allows continuous flow between two reservoirs while maintaining a unidirectional flow-rate in the cell culture chamber.


## Package contents:

- LineUp™ supply kit
- LineUp™ LINK
- Flow EZ 345 mbar (x3)
- PCAP 15 mL (x2)
- FLOW UNIT M
- Tubing & fitting kits
- SWITCH EZ
- 2-SWITCH™ (x2)



# CELL PERFUSION PACK

The science of flow control  
for Organ-On-A-Chip applications



Organ on chip study is an emerging field that brings substantial benefits compared to conventional cell culture. In many labs, considerable effort is put in choosing the right chip design, but the impact of flow control is still undetermined. It is our intent to create awareness on **the importance of flow** and its effects on one's studies.

- 1 LineUp™ – Flow Controller
- 2 FLOW UNIT – Flow Rate Sensor
- 3 P-CAPS – Sample Reservoirs
- 4 BE-FLOW – Microfluidic Chip

## Characteristics

| Perfusion Pack       | Product Number |
|----------------------|----------------|
| Starter Pack         | OOC-ST-001     |
| High Throughput Pack | OOC-HT-001     |



# MICROFLUIDIC CELL SORTING PACK

Microfluidic size cell sorting



The microfluidic **cell sorting package** allows sorting cells or particles according to their size using inertial and Dean drag forces. Fluigent has joined with microfluidic ChipShop to provide a solution with **well-designed microfluidic devices** and **excellent flow control** for sorting experiments.

- 1 FLOW UNIT – Flow rate sensor
- 2 LINK – Connect to the PC
- 3 Flow EZ (2 bar) – Flow controller
- 4 P-CAP – Sample reservoir
- 5 ChipShop sorting chip

Product Number

ESORT-PCK01

microfluidic  
**ChipShop**

Label-free

High  
throughput:  
up to 3 mL/min

Passive sorting:  
does not require  
additional  
external source

Up to 8 cells/  
particles sizes  
sorted per  
experiment



# DROP-SEQ PACK



The **Drop-Seq protocol**, is a high throughput method that enables the **sequencing of the mRNA** from a large number of cells. With this method it is possible to create a **gene expression map of the cell**, or even distinguish cell populations within a tissue. For optimized Drop-Seq experiments in Next-Generation Sequencing. The Fluigent **Drop-Seq Package** allows better **reproducibility** and control of both single cell and bead **encapsulation**.

- 1 LINK – PC control
- 2 Flow EZ – Flow controllers
- 3 P-CAP – Sample reservoirs
- 4 FLOW UNITS – flow rate sensors
- 5 Drop-Seq - microfluidic chip

Product Number

ODROPSEQFPCK

## Package contents

- LineUp FLOW EZ™ (x3)
- LineUp™ LINK (x1)
- FLOW UNIT M (x2)
- FLOW UNIT L (x1)
- P-CAP 15 mL (x1)
- P-CAP 2 mL (x2)
- Drop-Seq device FlowJEM PDMS Microfluidic Devices (x1)
- Fluigent Software
- LineUp™ Supply Kit (x1)
- Tubing & Fitting Kits

*Gain time:  
less than 1 minute to  
obtain droplets*

*Avoid losing cells  
or beads  
during transition time*

*Better control:  
avoid problems such  
as the backflow*

## Drop-Seq chip



Chip allowing for parallel analysis of mRNA expression in thousands of individual cells following their encapsulation in tiny droplets. These nanodroplets are formed by precisely combining aqueous and oil flows in a specially designed microfluidic device. Expression profiling can then be carried out in tens of thousands of cells in a matter of hours.

# ALGINATE BEADS GENERATION PACK



The Fluigent **Alginate beads production station** is a complete system for producing **outstanding monodispersed** Alginates beads with the flexibility to change **particle sizes** production in hundreds of milliseconds without interrupting the production. Its performance results from the combination of Fluigent's **LineUp™** microfluidic pumps and the **RayDrop** droplet generator, a breakthrough technology for **high-quality particles and beads production**.

- 1 Flow EZ™ pressure controllers
- 2 SWITCH EZ valve controller
- 3 LINK for software control
- 4 FLOW UNITs flow sensors
- 5 RayDrop droplet generator
- 6 2-SWITCH™ microfluidic valves
- 7 P-CAP sample reservoirs
- 8 Digital High-Speed Microscope

Product Number

O-SE-ALG-PCK

## Package contents

- LineUp Flow EZ™ 2 bar (x2)
- LineUp™ supply kit (x1)
- LineUp™ LINK (x1)
- FLOW UNIT M (x2)
- P-CAP 15 mL (x3)
- 2-SWITCH™ (x2)
- LineUp™ SWITCH EZ (x1)
- RayDrop Single Emulsion (x1)
- Tubing & fitting kit (x1)

|                                  |   |  |                        |                                 |
|----------------------------------|---|--|------------------------|---------------------------------|
| <i>Continuous<br/>production</i> | <i>High<br/>Monodispersity<br/>(2%)</i> | <i>High<br/>flexibility<br/>on alginate<br/>beads size</i> | <i>Ease<br/>of use</i> | <i>High<br/>reproducibility</i> |
|----------------------------------|---|--|------------------------|---------------------------------|

## Performance

| Feature                         | Extrusion Method | Fluigent Microfluidic Method |
|---------------------------------|------------------|------------------------------|
| Particle size distribution      | Up to 50%        | ~2%                          |
| Reproducibility                 | Low              | High                         |
| Live particle size control      | No               | Precise                      |
| Continuous (in-line) production | No               | Yes                          |

# NANOPARTICLE PRODUCTION PACK



The Fluigent **Liposome Nanoparticle Production Station** is a complete system for **precise and long-term production of liposomes** with high flexibility in terms of the liposome size range. Its performance results from the combination of Fluigent's **LineUp™** microfluidic pumps and the **RayDrop** device, a breakthrough technology for high-quality particles production.

- 1 Flow EZ™ pressure controllers
- 2 SWITCH EZ valve controller
- 3 LINK for software control
- 4 FLOW UNITs flow sensors
- 5 RayDrop droplet generator
- 6 2-SWITCH™ microfluidic valves
- 7 P-CAP sample reservoirs
- 8 Digital High-Speed Microscope

Product Number

O-MIX-LIPO-PCK

## Package contents

- LineUp Flow EZ™ 2 bar (x2)
- LineUp™ supply kit (x1)
- LineUp™ LINK (x1)
- FLOW UNIT M (x1)
- FLOW UNIT L (x1)
- P-CAP 15 mL (x3)
- 2-SWITCH™ (x2)
- LineUp™ SWITCH EZ (x1)
- RayDrop Single Emulsion (x1)
- Tubing & fitting kit (x1)

*Continuous  
/ In-line  
production*

*High flexibility  
in liposome  
nanoparticle size*

*High  
reproducibility*

*Full control  
on liposome  
nanoparticle size*

## Performance

| Feature                         | Extrusion Method | Fluigent Microfluidic Method |
|---------------------------------|------------------|------------------------------|
| Particle size distribution      | Up to 50%        | ~2%                          |
| Reproducibility                 | Low              | High                         |
| Live particle size control      | No               | Precise                      |
| Continuous (in-line) production | No               | Yes                          |

# PLGA MICROPARTICLE PRODUCTION PACK



- |                                 |                                 |
|---------------------------------|---------------------------------|
| 1 Flow EZ™ pressure controllers | 5 RayDrop droplet generator     |
| 2 SWITCH EZ valve controller    | 6 2-SWITCH™ microfluidic valves |
| 3 LINK for software control     | 7 P-CAP sample reservoirs       |
| 4 FLOW UNITs flow sensors       | 8 Digital High-Speed Microscope |

The PLGA microparticle production station is a **robust, high-performance** solution to **generate polymer microparticles** in a **homogenous** and **fully controlled manner**. The performance brought by the **RayDrop droplet generator**, with the combination of polylactic-co-glycolic acid as an encapsulation polymer and ethyl acetate as a solvent provides a **biocompatible solution** lowering both hazard risk and precipitation time. Suitable for **biological applications**, the RayDrop, and its station offer a **semi-automated** solution for one of the most successful drug delivery systems in laboratories and clinics. To solve the current problems of the PLGA microparticles synthesis, droplet-based microfluidics appears to be a powerful tool. Droplet control and generation allow **highly monodispersed** and **continuous production** as compared to batch emulsion methods.

## Characteristics

| Package            | Description     |
|--------------------|-----------------|
| Standard package   | 1DPPL01         |
| Automation package | O-SE-PLGAAP-PCK |

## Package contents

| Content                   | Description  | Package                                |
|---------------------------|--|--|
| Flow EZ™                  | Two 2 bar pressure-based flow controllers                                | Standard Package<br>Automation Package |
| LINK                      | Monitor and automate the setup with Fluigent Software                    | Standard Package<br>Automation Package |
| FLOW UNITS                | Two high precision flow sensors (M and L range)                          | Standard Package<br>Automation Package |
| P-CAPS                    | Air-tight caps for sample reservoirs (one of 50 mL, two of 15 mL)        | Standard Package<br>Automation Package |
| RayDrop                   | Single emulsion droplet generator  | Standard Package<br>Automation Package |
| Fluigent Software         | Software to control in real time and automate instruments                | Standard Package<br>Automation Package |
| 2-SWITCH™                 | Two bidirectional valves   | Automation Package                     |
| SWITCH EZ                 | Microfluidic valves controller   | Automation Package                     |
| Tubing & Fitting Kit      | Dedicated tubing & fitting elements to set up the system                 | Standard Package<br>Automation Package |
| Digital High-Speed Camera | Observe and record results with high resolution and rates up to 7092 fps | Optional                               |

|   |  |  |
|---|--|--|
| <i>Particle monodispersity ~2%</i>                | <i>Particle size control and encapsulation performance</i> | <i>System recovery and cleaning</i>    |
| <i>Biocompatibility: Ethyl acetate as solvent</i> | <i>Continuous / In-line production</i>                     | <i>Semi-automated priming protocol</i> |



# DOUBLE EMULSION GENERATION PACK



The **Fluigent Double Emulsion Production Station** is a complete system for producing outstanding **monodispersed double emulsions in one single device**. Its performance results from the combination of Fluigent's **LineUp™** microfluidic pumps and the **RayDrop** device, a breakthrough technology for **high monodispersed and stable double emulsion**.

- 1 Flow EZ™ pressure controllers
- 2 SWITCH EZ valve controller
- 3 LINK for software control
- 4 FLOW UNITS flow sensors
- 5 RayDrop droplet generator
- 6 2-SWITCH™ microfluidic valves
- 7 P-CAP sample reservoirs
- 8 Digital High-Speed Microscope

Product Number

O-DE-STD-PCK

## Package contents

- LineUp Flow EZ™ 2 bar (x3)
- LineUp™ supply kit (x1)
- LineUp™ LINK (x1)
- FLOW UNIT M (x2)
- FLOW UNIT L (x1)
- 15 mL P-CAP (x1)
- 50 mL P-CAP (x2)
- 2-SWITCH™ (x2)
- LineUp™ SWITCH EZ (x1)
- RayDrop Double Emulsion (x1)
- Tubing & Fitting kit (x1)

|                                   |   |  |
|-----------------------------------|---|--|
| <i>High monodispersity (2%)</i>   | <i>Double emulsion size: 70 µm – 150 µm</i> | <i>High Frequency: 5 000 Hz</i>        |
| <i>Adjustable shell thickness</i> | <i>No surface coating</i>                   | <i>w/o/w and o/w/o double emulsion</i> |

## Performance

| Feature                         | Batch method                                      | Fluigent microfluidic method      |
|---------------------------------|---|-----------------------------------|
| Process                         | Multiple step process                             | Direct double emulsion production |
| Homogeneity                     | Random distribution of single and double emulsion | High homogeneity                  |
| Particle size distribution      | More than 50%                                     | ~2%                               |
| Reproducibility                 | Low   | High                              |
| Live particle size control      | No  | Precise                           |
| Continuous / In-line production | No  | Yes                               |

# MICROFLUIDIC DROPLET PACK

Ideal for biologists and chemists



The **Droplet Package** is designed for all who want to start with **microfluidic droplets experiments**. It includes the **EZ Drop** chips and fully adapted liquid handling solution and accessories.

- 1 LineUp™ pressure-based flow controllers
- 2 LINK for Software control
- 3 High precision flow sensors
- 4 Sample reservoirs
- 5 Droplet generation chips

*Highly monodisperse droplets*

*Most simple droplet generation kit*

*Cost-effective solution*

## Package contents

- Flow EZ™ 2 or 1 bar (x2)
- LineUp™ supply kit (x1)
- LINK (x1)
- FLOW UNIT S (x2)
- P-CAP 2 mL (x2)
- EZ Drop chip (x1)
- Tubing & fitting kits


### Product Number

DROPPACK-01

DROPPACK-02

# EDUCATIONAL PACK

Complete microfluidic setup  
for starting experiment



The **Fluigent Educational Packages** provide a **broad introduction to microfluidics** and its applications by familiarizing the user with general **microfluidic principles** and microfluidic systems. These **ready-to-teach packages** are specifically handy for professors and teachers.

*Flexible offer  
with 4 packages  
available*

*Up to 4 hours of  
practical work  
with solutions*

*A handbook  
for an overview on  
microfluidics*

## FIRST LEVEL: CO-FLOW

(SEDUC-COFLOW)

**A beginner package to discover microfluidics by experimenting visually a pillar concept of microfluidics: laminar flows.**

1 hour guided experiments. Suited for: biologists, (bio)engineers, chemists, but also for high schools with scientific programs.

### Package contents:

- LineUp™ supply kit
- Flow EZ 1000 mbar (x2)
- PCAP 15 mL (x2)
- Co-Flow chip (x3)
- Tubing & fitting kit
- Microscope with SD memory card
- Dye solutions
- Printed Handbook
- Exp. Leaflet Co-Flow
- Accessories bag (tubing cutter, pen, notebook, counter, ruler ...)

# FIRST LEVEL: RESISTANCE

(SEDUC-RESIST)

**Master and take advantage of one powerful tool for optimizing your microfluidic experiments: hydrodynamic resistance.** 1 hour guided experiments:

- Perform resistance measurements on the microfluidics system
- Theoretical calculations related to the experiments

## Package contents:

- |                      |                           |  |
|----------------------|---------------------------|--|
| ■ LineUp™ supply kit | ■ Tubing & fitting kit    | ■ Accessories bag (tubing cutter, pen, notebook, counter, ruler ...) |
| ■ Flow EZ 1000 mbar  | ■ Dye solutions           |  |
| ■ PCAP 15 mL         | ■ Printed Handbook        |  |
| ■ FLOW UNIT M        | ■ Exp. Leaflet Resistance |  |

# SECOND LEVEL: CO-FLOW & RESISTANCE

(SEDUC-RESITCOFL)

**Co-flow and Resistance in one package for a first overview on microfluidics.**

2 hours guided experiments. Suited for: (bio)engineers, physicists, chemists.

## Package contents:

- |                          |                                  |  |
|--------------------------|----------------------------------|--|
| ■ LineUp™ supply kit     | ■ Tubing & fitting kit           | ■ Exp. Leaflet Co-Flow   |
| ■ Flow EZ 1000 mbar (x2) | ■ Microscope with SD memory card | ■ Exp. Leaflet Resistance  |
| ■ FLOW UNIT M            | ■ Dye solutions                  | ■ Accessories bag (tubing cutter, pen, notebook, counter, ruler ...) |
| ■ PCAP 15 mL (x2)        | ■ Printed Handbook               |  |
| ■ Co-Flow chip (x3)      |                                  |  |

# FULL COURSE: CO-FLOW, RESISTANCE & DROPLETS

(SEDUC-DROPLET)

**Get the most complete overview, with experiments pushed to real-world applications: droplet generation.**

4 hours guided experiments. Suited for: (bio)engineers, chemical engineering, physicists, biologists and researchers.

## Package contents:

- |                          |                                  |  |
|--------------------------|----------------------------------|--|
| ■ LineUp™ supply kit     | ■ dOil 120mL                     | ■ Exp. Leaflet Co-Flow   |
| ■ Flow EZ 1000 mbar (x2) | ■ Microbeads bottle              | ■ Exp. Leaflet Resistance  |
| ■ LINK                   | ■ Fluigent Software              | ■ Exp. Leaflet Droplet   |
| ■ FLOW UNIT M (x2)       | ■ Tubing & fitting kit           | ■ Accessories bag (tubing cutter, pen, notebook, counter, ruler ...) |
| ■ PCAP 15 mL (x2)        | ■ Microscope with SD memory card |  |
| ■ Co-Flow chip (x3)      | ■ Dye solutions                  |  |
| ■ EZ Drop chip (x3)      | ■ Printed Handbook               |  |
| ■ dSurf 2% 12mL          |                                  |  |

# SYSTEMS & PACKS

# TUBING & FITTING KITS

## ARIA SINGLE OUTPUT KIT

(2SW-KIT-AR)



Tubing Kit for external **2-SWITCH™** addition for the **Aria** device.

### Kit contents:

- Flangeless fittings 1/4-28 (x6)
- FEP tubing 1/16" OD 0.010" ID (2 m)
- Blue ferrules (x12)
- Adapters 10-32 (x2)
- RJ45 female coupler (x1)
- Flat RJ45 flat cable (x1)

## ARIA REPLACEMENT TUBING KIT


(CTQ-KIT-AR)



Pre-cut set of tubings to replace internal tubings of **Aria** device.

## DROP-SEQ KIT

(ODROPSEQCTQ)



Tubing & fitting kit for the **Drop-Seq** microfluidic chip.

### Kit contents:

- Sleeves 1/16" → 1/32" (x5)
- Flangeless fittings 1/4-28 (x5)
- 1/32"x0.010x100ft (LQ) PEEK tubing (2 m)
- Blue ferrules (x10)




## DOUBLE EMULSION GENERATION KIT

(O-DE-STD-CTK)

A kit containing tubing, fitting, plug, and ferule elements required for the **Double Emulsion Generation**.

### Kit contents:

- 1/32" PEEK tubing with 250 µm inner diameter (ID) (1 m)
- 1/32" PEEK tubing with 150 µm inner diameter (ID) (2 m)
- FEP tubing 1/16" OD 0.010" ID (1 m)
- FEP tubing 1/16" OD 0.020" ID (1 m)
- Sleeves 1/16" → 1/32" (x16)
- PEEK plugs 1/4-28 Black (x4)
- Flangeless fittings 1/4-28 (x5)
- Blue ferrules (x10)
- Inline filters (x3)
- 1/16" unions (x4)
- 10-32 adapters (x10)
- Tube cutter (x1)
- RayDrop support (x1)
- Lateral rails (x4)



## NANOPARTICLE PRODUCTION KIT

(O-MIX-LIPO-CTK)

A kit containing tubing, fitting, plug, and ferule elements required for the **Nanoparticle Production pack**.

### Kit contents:

- 1/32" PEEK tubing with 250 µm inner diameter (ID) (1 m)
- 1/32" PEEK tubing with 150 µm inner diameter (ID) (2 m)
- FEP tubing 1/16" OD 0.010" ID (1 m)
- FEP tubing 1/16" OD 0.020" ID (1 m)
- Sleeves 1/16" → 1/32" (x16)
- PEEK plugs 1/4-28 Black (x4)
- Flangeless fittings 1/4-28 (x5)
- Blue ferrules (x10)
- Inline filters (x3)
- 1/16" unions (x4)
- 10-32 adapters (x10)
- Tube cutter (x1)
- RayDrop support (x1)
- Lateral rails (x4)



## PLGA MICROPARTICLE PRODUCTION KITS: STANDARD PACK

(IDPPLC1)

A kit containing tubing, fitting, plug, and ferrule elements required for the **PLGA Microparticle Production pack**.

### Kit contents:

- RayDrop tubing & fitting kit (x1)
- 1/32" PEEK tubing with 250 µm inner diameter (ID) (1 m)
- FEP tubing 1/16" OD 0.010" ID (3 m)
- transparent tube 4 mm (1 m)
- Sleeves 1/16" → 1/32" (x6)
- PEEK plugs 1/4-28 Black (x2)
- Flangeless fittings 1/4-28 (x8)
- Blue ferrules (x12)
- Fluidic T (x1)
- T fitting 4 mm (x1)
- Flowmeter fittings (x 4)
- Tube cutter (x1)
- RayDrop support (x1)
- Lateral rails (x4)



## PLGA MICROPARTICLE PRODUCTION KIT: AUTOMATION PACK

(IDPPLC2)

A kit containing tubing, fitting, plug, and ferrule elements required for the **PLGA Microparticle Production pack**.

### Kit contents:

- RayDrop tubing and fitting kit (x1)
- 1/32" PEEK tubing with 250 µm inner diameter (ID) (1 m)
- FEP tubing 1/16" OD 0.010" ID (3 m)
- Transparent tube 4 mm (1 m)
- Sleeves 1/16" → 1/32" (x6)
- PEEK plugs 1/4-28 black (x2)
- Flangeless fittings 1/4-28 (x15)
- Blue ferrules alone (x20)
- T fitting 4mm (x1)
- Flowmeter fittings (x4)
- Tube cutter (x1)
- RayDrop support (x1)
- Lateral rails (x4)






## ALGINATE BEADS GENERATION KIT

(O-SE-ALG-CTK)

A kit containing tubing, fitting, plug, and ferrule elements required for the **Alginate bead generation pack**.

### Kit contents:

- 1/32" PEEK tubing with 250 µm inner diameter (ID) (0.5 m)
- 1/32" PEEK tubing with 150 µm inner diameter (ID) (1 m)
- FEP tubing 1/16"OD 0.010"ID (0.5 m)
- FEP tubing 1/16"OD 0.020"ID (0.4 m)
- Sleeves 1/16" → 1/32" (x6)
- 1/16" unions (x3)
- PEEK plugs 1/4-28 Black (x3)
- Flangeless fittings 1/4-28 (x4)
- Blue ferrules (x8)
- Tube cutter (x1)
- RayDrop support (x1)
- Lateral rails (x4)



## RAYDROP SINGLE EMULSION DEVICE KIT

(IDPRDC1)

A kit containing tubing, fitting, plug, and ferrule elements required for using the **RayDrop Single Emulsion Device**.

### Kit contents:

- 1/32" PEEK tubing with 150 µm inner diameter (ID) (0.2 m)
- FEP tubing 1/16" OD 0.010" ID (2 m)
- FEP tubing 1/16" OD 0.020" ID (0.1 m)
- Sleeves 1/16" → 1/32" (x2)
- 1/16" unions (x3)
- PEEK plugs 1/4-28 Black (x3)
- Flangeless fittings 1/4-28 (x4)
- Blue ferrules (x8)
- Tube cutter (x1)
- RayDrop support (x1)
- Lateral rails (x4)

4  
—

# FLUIGENT INDUSTRIAL

# INTRODUCTION

For **integration** into **manufactured systems**, Fluigent has developed a wide range of **OEM products & services**. These are dedicated for **flow control** and **fluid handling** in microfluidic and nanofluidic applications. Equipments can easily be integrated into more complex instruments due to their compact design and **Software Development Kit**.

**Fluigent complete OEM product & service line offers :**

## MICROFLUIDIC OEM EQUIPMENTS

(fluid management solutions, pressure control and generation, switching solutions, flow-rate measurement, ...)

## CUSTOMIZABLE OEM FLOW CONTROL MODULES


from Fluigent proprietary technologies (specific development from mechanical to fluidic, electronic and software, component integration, specific settings, ...)

## FULLY CUSTOM MICROFLUIDIC DEVICE

(Any pressure / flow rate configuration, dedicated PCB, software integration)

# F-OEM

## Modular OEM microfluidic flow controller



Our **F-OEM** offers our highest performance, efficiency, and widest pressure and flow rate ranges to support the most demanding industrial applications, including microfluidic and nanofluidic applications (microchannels, nanochannels, capillaries, lab on a chip...). It is a standalone, modular platform that will perform complex fluidic operations.


*Compact  
standalone platform  
for industrial uses*

*Flexible and  
configurable  
with pressure & switch*

*Contamination-free  
Not in contact with  
liquids*

# P-OEM

## Microfluidic flow management unit



The **OEM Microfluidic modules** are a range of products for **flow control** and **fluid handling** in microfluidic and nanofluidic applications. Dedicated for industrial microfluidic applications they can economically replace conventional syringe pumps or peristaltic pumps.

The multi-channel design is **highly compact**, ranging from **1 to 8 channels** for the OEM microfluidic modules. As fluids are not in contact with the instrument, there is no routine cleaning and cross-contamination is drastically reduced. **Positive or Negative pressure ranges** can be applied independently on different channels of the same module.

The optional advanced multi-channel **Flow-Rate Control Module Software** combined with the Flow Sensor Unit device makes it possible to have **full control of flow rates** in any coupled multi-channel configuration, while keeping the benefits of pressure actuation. Switching solutions can also be integrated to extend the **fluid management capabilities**. OEM versions of the **2-SWITCH** (3-port/2-way) and **M-SWITCH** (11-port/10-position) are available for integration.

*Highly  
customizable  
through numerous  
options*

*Adapted to  
industrial  
usage*

*Quality and  
application  
expertise*

*Compact  
solution*

# PX SERIES

## OEM microfluidic pressure controller



The **PX-series**, is the most versatile industrial pressure controller. It is **fast and stable** thanks to our field proven, patented FASTAB™ technology.

The PX-series is **designed for OEM integration in industrial environments**. Its RS232 and USB ports ensure its communication with any system. The series presents **5 versions** regarding pressure range requirements.

*High-Quality  
Pressure  
Control*

*Unmatched  
Price-  
Performance  
Ratio*

*Designed for  
OEM  
Integration*

*Patented  
technology*

## Characteristics

| PX model | Range in mbar | P/N      |
|----------|---------------|----------|
| PX-345   | 0 to 345      | ICPX3451 |
| PX-1     | 0 to 1000     | ICPX11   |
| PX-2     | 0 to 2000     | ICPX21   |
| PX-V1    | 0 to -600     | ICPXV1   |
| PX-V2    | 0 to -750     | ICPXV2   |

# NIFS

P/N : INIFS

## Non-Invasive Flow Sensor for Industrial Integration

Our OEM non-invasive flow sensor (NIFS) is dedicated to flow rate monitoring and control.



When combined with Fluigent pressure controllers, it allows pressure-based flow rate control without fluid contact, no liquid calibration is needed. It allows for accurate and repeatable measurements of liquid flow rates from 100  $\mu\text{L}/\text{min}$  and to 10 mL/min bi-directionally.

*High performance*  
Industry-proven  
technology

*Versatile*  
use for all liquids, wide  
flow rate range

*Reliable*  
No liquid contact

## Key applications

Including cell culture under flow, blood/plasma perfusion, organ-on-a-chip, single/double emulsions, microparticles or cell-encapsulation.

- **Cell biology**

Many biological experiments require disposable or sterilized components in the fluidic path. The use of NIFS is ideal for these applications.

- **Microfluidic droplet generation**

As there is no calibration needed to manipulate several liquids, the NIFS saves time when users want to generate droplets, a process that involves different fluids (aqueous solutions, oil, surfactants, alcohols...).

- **Biochemical and molecular analysis**

A sterile environment and lack of need for a calibration step are beneficial for drug, enzyme or food analysis through screening with droplet microfluidics.

## Key features

**High accuracy and precision:** With Fluigent's regulation control technology, the NIFS can measure flow rate with 5% accuracy.


**Liquid contamination-free:** As the non-invasive flow sensor is placed in the pneumatic path, it does not touch the liquid. It is suitable for applications where sterility is required. No cleaning is necessary.

**Fluid calibration-free:** Unlike traditional flowmeters, the NIFS functions independently of the liquid properties.

## Related products


### F-OEM

The F-OEM offers our highest performance, efficiency, and widest pressure and flow rate ranges to support the most demanding industrial applications



### PX-SERIES

The PX-series, is the most versatile industrial pressure controller. It is fast and stable thanks to our field proven, patented FASTAB™ technology.





# FS+ SERIES

## Microfluidic OEM flow sensor



Our FS series is dedicated to flow-rate control and monitoring. When combined with Fluigent pressure controllers, it permits pressure-based flow control. It enables precise and accurate measurements of dynamic liquid flow rates from 0 – 1.5  $\mu\text{L}/\text{min}$  and up to 40  $\text{mL}/\text{min}$  bi-directionally.

| FS+ series model                 | XS                      |                          | S                        |                        | M+                        |                        | L+                      |  |
|----------------------------------|-------------------------|--------------------------|--------------------------|------------------------|---------------------------|------------------------|-------------------------|--|
| Part number                      | IFSXS1                  |                          | IFSS1                    |                        | FLU-M+-OEM                |                        | FLU-L+-OEM              |  |
| Sensor Inner Diameter            | 25 μm                   |                          | 150 μm                   |                        | 400 μm                    |                        | 400 μm                  |  |
| Maximum Pressure                 | 200 bar                 |                          | 200 bar                  |                        | 12 bar                    |                        | 12 bar                  |  |
| Wetted Materials                 | PEEK & Quartz Glass     |                          |                          |                        | PPS, stainless steel 316L |                        |                         |  |
| Calibrated Media                 | Water                   | Water                    | IPA                      | Water                  | IPA                       | Water                  | IPA                     |  |
| Range                            | 0±1.5 μL /min           | 0±7 μL /min              | 0±70 μL /min             | 0±2 mL /min            | 0±2 mL /min               | 0±40 mL /min           | 0±40 mL /min            |  |
| Accuracy (measured value)        | 10% mv above 75 nL /min | 5% mv above 0.42 μL /min | 20% mv above 4.2 μL /min | 5% mv above 10 μL /min | 10% mv above 50 μL /min   | 5% mv above 10 μL /min | 10% mv above 50 μL /min |  |
| Lowest detectable flow increment | 3.7 nL/min              | 10 nL/min                |                          | 0.25 μL/min            |                           | 25 μL/min              |                         |  |

Measured Values from 5% to 100% of product range in normal conditions

Short response time  
pressure-based  
direct flow control

High reliability  
and long-term  
stability

Reliable  
industry-proven  
technology

Intuitive  
and easy-to-use  
systems

| FS series model                  | XL                          | L                         |                           | M                        |                          |
|----------------------------------|-----------------------------|---------------------------|---------------------------|--------------------------|--------------------------|
| Product Number                   | IFSXL1                      | IFSL1                     |                           | IFSM1                    |                          |
| Sensor Inner Diameter            | 1.8 mm                      | 1.0 mm                    |                           | 430 μm                   |                          |
| Maximum Pressure                 | 5 bar                       | 12 bar                    |                           | 100 bar                  |                          |
| Wetted Materials                 | PEEK and Borosilicate Glass |                           |                           |                          |                          |
| Calibrated Media                 | Water                       | Water                     | IPA                       | Water                    | IPA                      |
| Range                            | 0±5mL/min                   | 0±1mL/min                 | 0±10mL/min                | 0±80μL/min               | 0±500μL/min              |
| Accuracy (measured value)        | 5% m.v. above 0.2 mL/min    | 5% m.v. above 0.04 mL/min | 20% m.v. above 0.5 mL/min | 5% m.v. above 2.4 μL/min | 20% m.v. above 25 μL/min |
| Lowest detectable flow increment | 3 μL/min                    | 0.7 μL/min                |                           | 0.06 μL/min              |                          |

# RX

## OEM microfluidic pressure source

| Model | P/N    |
|-------|--------|
| RX    | ISRX21 |



The **RX** is designed to be seamlessly integrated into OEM systems or used as a **standalone pressure source**. Packaged in a robust steel enclosure, it provides dried and filtered air at **up to 2500 mbar** for one or multiple pressure control module(s) such as the **PX-series** or other instruments which need compressed air to operate.

*Compact  
& versatile*

*Condensation-proof  
and filtered air*

*Designed for  
OEM Integration*

# FULLY CUSTOM MICROFLUIDIC DEVICE




Fluigent can assist with the design and manufacture **customized platforms** for specific applications and needs, thanks to detailed product design and application engineering capabilities.

- ▶ Pioneers in innovative flow control solutions
- ▶ From product customization to specific solution development
- ▶ High Performance and Robustness
- ▶ Flexibility and compactness
- ▶ In house design and assembly
- ▶ Quality standard certification
- ▶ Broad range of applications
- ▶ Customer needs focused
- ▶ Highly experienced customer service

Whether you're a start-up or a larger organization, Fluigent can drive your microfluidic project from conception to design and manufacturing.

By partnering with us, customers benefit from Fluigent's portfolio of products, technologies, patents, and 15 years of experience with microfluidic technology.

We combine our expertise and knowledge to design and manufacture the highest quality custom microfluidic devices. Partner with Fluigent to accelerate your product time to market.



*Fluigent expertise*  
Microfluidic  
knowledge

*Flexible*  
We adapt to your  
specifications needs

*Cost-efficient*  
Device and  
development

# Fluigent's unique portfolio of technologies


- **DFC Flow Rate Control Algorithm**  
Self-learning algorithm allowing live pressure/flow-rate monitoring and control
- **Micropumps algorithm**  
Fluigent algorithm in micropumps for compact pressure supply and control
- **Fluid mixer & Temperature management**  
Biological vortex mixer. Orbital motion, mass compensated vibration, temperature modules for heating and cooling
- **Robotic pipetting**  
x-y-z robotic stage for fluid sampling into microwells
- **Non-intrusive flow sensing**  
Contact-free, no calibration needed, high precision and accuracy

## Custom software development

We provide custom software based on our well-known end-user software OxyGEN as well as on our SDK library for your custom microfluidic device.



## Development Process



5  
—

# SUPPORT & TOOLS

# TIPS


*Working in a microfluidic environment almost automatically means using fittings and tubing to connect your microfluidic device or your Lab-on-a-chip to the various elements of your microfluidic circuit.*

- *Tubing enables you to connect the various elements of your microfluidic circuit.*
- *Fittings enable you to attach, adapt or adjust the tubing to these elements, ensuring tight connections.*

*Tubing and sleeves are defined by their diameter, length and material.*

## Diameter

Inner diameter (often abbreviated as "ID") is diameter through which the fluid flows. The inner section times the length of the tubing gives you the internal volume of the tubing. The inner diameter plays a significant role in the fluidic resistance to flow brought by the tubing. The smaller the diameter is, the more resistant the tubing will be.



## Length

Usually the tubing is made as short as possible to have smaller internal volumes. It is also a parameter that contributes to the resistance of the tubing.


In order to get a clean interface and prevent any clogging or collapsing of the fluidic path, all tubing should be cut with specifically designed cutters.

## Materials

A wide range of materials are available for the same ID/OD combination. The material should be selected according to the nature of the reagents flowing through the tubing. Be careful to check the chemical compatibility of the tubing before installing it in your application. Some of the most common materials for microfluidic tubing include:

- **PEEK (Polyetheretherketone):** very good chemical resistance and biocompatibility, low non-specific adsorption. Rigid and opaque. For low and high pressure applications. Very small internal diameters available.
- **FEP (Fluorinated ethylene-propylene):** does not react with most chemicals and is biocompatible. Flexible and transparent. Mostly for low-pressure applications (no higher than 7 bar).

# Place your tubing on the page for a quick reference



## OD (OUTER DIAMETER) IDENTIFICATION

- 1/16" (1.6 mm)
- 1/32" (0.8 mm)

## UNIT CONVERSIONS

|                                     |                                      |                  |
|-------------------------------------|--------------------------------------|------------------|
| 1 bar = 14.5 psi                    | 1 mbar = 1.45 X 10 <sup>-2</sup> psi | 1 inch = 2.54 cm |
| 1 psi = 6.90 X 10 <sup>-2</sup> bar | 1 psi = 68.95 mbar                   | 1 cm = 0.39 inch |

# CUSTOMER SUPPORT & SERVICE

Service & Technical Support from Fluigent-certified experts

At Fluigent, we understand that a non-functioning system means time lost in the lab. **Our customer support team is dedicated to performing timely, cost effective repairs.** Our application experts are available to advise you on any questions you may have on the use of our products and how to adapt them to different experimental designs.

|  |   |                                      |
|--|---|--------------------------------------|
| Dedicated<br>Fluigent-certified team<br>& trained partners | Diagnosis through<br>remote session<br>or on-site visit | Customer<br>satisfaction<br>oriented |
|--|---|--------------------------------------|

# REFERENCE GUIDE

| Product Number                                 | Product Name                                | Page | Notes |
|--|---|------|-------|
| <b>Research Instruments - Pressure Systems</b> |   |      |       |
| <a href="#">LU-FEZ-0025</a>                    | LineUp Flow EZ™ 25 mbar                     | 11   |       |
| <a href="#">LU-FEZ-0069</a>                    | LineUp Flow EZ™ 69 mbar                     | 11   |       |
| <a href="#">LU-FEZ-0345</a>                    | LineUp Flow EZ™ 345 mbar                    | 11   |       |
| <a href="#">LU-FEZ-1000</a>                    | LineUp Flow EZ™ 1000 mbar                   | 11   |       |
| <a href="#">LU-FEZ-2000</a>                    | LineUp Flow EZ™ 2000 mbar                   | 11   |       |
| <a href="#">LU-FEZ-7000</a>                    | LineUp Flow EZ™ 7000 mbar                   | 11   |       |
| <a href="#">LU-FEZ-N025</a>                    | LineUp Flow EZ™ -25 mbar                    | 11   |       |
| <a href="#">LU-FEZ-N069</a>                    | LineUp Flow EZ™ -69 mbar                    | 11   |       |
| <a href="#">LU-FEZ-N345</a>                    | LineUp Flow EZ™ -345 mbar                   | 11   |       |
| <a href="#">LU-FEZ-N800</a>                    | LineUp Flow EZ™ -800 mbar                   | 11   |       |
| <a href="#">ELUPPU1000</a>                     | LineUp™ Push-Pull                           | 12   |       |
| <a href="#">ELUPSW2000</a>                     | LineUp™ P-SWITCH                            | 13   |       |
| <a href="#">LU-ADP-0001</a>                    | LineUp™ Adapt                               | 15   |       |
| <a href="#">LU-LNK-0002</a>                    | LineUp™ LINK                                | 16   |       |
| <a href="#">ELULNK232</a>                      | LineUp™ LINK COM                            | 16   |       |
| <a href="#">EZ-00345001</a>                    | MFCST™-EZ 345 mbar                          | 18   |       |
| <a href="#">EZ-01000001</a>                    | MFCST™-EZ 1000 mbar                         | 18   |       |
| <a href="#">EZ-01000002</a>                    | MFCST™-EZ 2000 mbar                         | 18   |       |
| <a href="#">EZ-07000001</a>                    | MFCST™-EZ 7000 mbar                         | 18   |       |
| <a href="#">EZ-80345001</a>                    | MFCST™-EZ -345 mbar                         | 18   |       |
| <a href="#">EZ-80800001</a>                    | MFCST™-EZ -800 mbar                         | 18   |       |
| <a href="#">EX-00345001</a>                    | MFCST™-EX 345 mbar                          | 18   |       |
| <a href="#">EX-01000001</a>                    | MFCST™-EX 1000 mbar                         | 18   |       |
| <a href="#">EX-01000002</a>                    | MFCST™-EX 2000 mbar                         | 18   |       |
| <a href="#">EX-07000001</a>                    | MFCST™-EX 7000 mbar                         | 18   |       |
| <a href="#">EX-80345001</a>                    | MFCST™-EX -345 mbar                         | 18   |       |
| <a href="#">EX-80800001</a>                    | MFCST™-EX -800 mbar                         | 18   |       |
| <a href="#">EZ-11000001</a>                    | MFCST™-EZ Basic Base                        | 18   |       |
| <a href="#">EZ-source-pos</a>                  | MFCST™-EZ Positive Pressure Source Included | 18   |       |
| <a href="#">EZ-source-neg</a>                  | MFCST™-EZ Negative Pressure Source Included | 18   |       |
| <a href="#">EX-11000008</a>                    | MFCST™-EX Basic Base                        | 18   |       |
| <a href="#">EX-source-pos</a>                  | MFCST™-EX Positive Pressure Source Included | 18   |       |
| <a href="#">EX-source-neg</a>                  | MFCST™-EX Negative Pressure Source Included | 18   |       |
| <b>Research Instruments - Sensors</b>          |   |      |       |
| <a href="#">FLU-XS</a>                         | FLOW UNIT XS                                | 20   |       |
| <a href="#">FLU-S-D</a>                        | FLOW UNIT S                                 | 20   |       |



**Research Instruments - Sensors**

|          |                  |    |
|----------|------------------|----|
| FLU-M-D  | FLOW UNIT M      | 20 |
| FLU-L-D  | FLOW UNIT L      | 20 |
| FLU-XL   | FLOW UNIT XL     | 20 |
| FLU-M+   | FLOW UNIT M+     | 20 |
| FLU-L+   | FLOW UNIT L+     | 20 |
| FLB      | Flowboard        | 22 |
| EIPS345  | PRESSURE UNIT S  | 24 |
| EIPS1000 | PRESSURE UNIT M  | 24 |
| EIPS7000 | PRESSURE UNIT XL | 24 |

**Research Instruments - Microfluidic Valves**

|           |                   |    |
|-----------|-------------------|----|
| 2SW002    | 2-SWITCH™         | 27 |
| ESSMSW003 | M-SWITCH™         | 28 |
| LSW001    | L-SWITCH™         | 29 |
| ELUSEZ    | LineUp™ SWITCH EZ | 14 |

**Research Instruments - Pressure and Vacuum Sources**

|               |  |    |
|---------------|--|----|
| FLPG003       | FLPG+ High Pressure Supply             | 30 |
| FLPG004       | FLPG+ with Incubator Aspiration Option | 30 |
| FLPG005       | FLPG+ Silent Pressure Supply           | 30 |
| E-AC-RX1-2500 | Compact Pressure Source                | 31 |
| EACVACPUMP    | Compact Vacuum Pump                    | 32 |

**Research Instruments - Sample Reservoirs**

|             |                                  |    |
|-------------|----------------------------------|----|
| P-CAP2-LP   | P-CAP 2 mL Low Pressure          | 33 |
| P-CAP2-HP   | P-CAP 2 mL High Pressure         | 33 |
| P-CAP15-LP  | P-CAP 15 mL Low Pressure         | 33 |
| P-CAP15-HP  | P-CAP 15 mL High Pressure        | 33 |
| P-CAP50-LP  | P-CAP 50 mL Low Pressure         | 33 |
| P-CAP50-HP  | P-CAP 50 mL High Pressure        | 33 |
| 14000501    | Fluidwell 0.5 mL Low Pressure    | 34 |
| 24000501    | Fluidwell 0.5 mL High Pressure   | 34 |
| 14002001    | Fluidwell 2 mL Low Pressure      | 34 |
| 24002001    | Fluidwell 2 mL High Pressure     | 34 |
| 11015001    | Fluidwell 15 mL Low Pressure     | 34 |
| 21015001    | Fluidwell 15 mL High Pressure    | 34 |
| 11050001    | Fluidwell 50 mL Low Pressure     | 34 |
| 21050001    | Fluidwell 50 mL High Pressure    | 34 |
| RES-CAP     | Bottle-CAP 1 fluidic port        | 35 |
| RES-CAP-PCK | Bottle-CAP 1 fluidic port + Kit  | 35 |
| RES-CAP-3P  | Bottle-CAP 2 fluidic ports + Kit | 35 |
| RES-CAP-4P  | Bottle-CAP 3 fluidic ports + Kit | 35 |

**Research Instruments - Standalone Accessories**

|            |                               |    |
|------------|-------------------------------|----|
| P-CAP-RACK | P-CAP series Rack             | 36 |
| CTQ-MANI   | 10-Position Pressure Manifold | 37 |

|   |  |    |
|---|--|----|
| CTQ-006BT   | Bubble Trap + Kit                              | 37 |
| 10000001  | Air Flow Regulation Kit                        | 37 |
| OSEDUCMIC   | Educational Microscope                         | 36 |
| <b>Research Instruments - Surfactant</b>                |  |    |
| DR-RE-SU-12   | dSurf 3 x 4mL                                  | 38 |
| DR-RE-SU-30   | dSurf 30mL                                     | 38 |
| ODSURF5G  | Honey dSurf 5g                                 | 38 |
| DR-RE-SU-A1   | dOil   | 39 |
| <b>Research Instruments - Microfluidic Chips</b>        |  |    |
| 1DPRD01   | RayDrop single emulsion 30 µm - 150 µm         | 40 |
| ORDRPSE-60-300  | RayDrop single emulsion 60 µm - 300 µm         | 40 |
| ORDRPSE-90-450  | RayDrop single emulsion 90 µm - 450 µm         | 40 |
| ORDRPDE-30-70-150                                       | RayDrop double emulsion 30 µm - 70 µm - 150 µm | 41 |
| DROPKIT01   | EZ Drop droplet generation chip                | 42 |
| ODROPSEQCHIP  | Drop-Seq chip                                  | 43 |
| OOC-FLOW-01   | Be-Flow chip                                   | 44 |
| OOC-DBLE-08   | Be-Doubleflow chip                             | 44 |
| OOC-GRAD-05   | Be-Gradient chip                               | 45 |
| OOC-TRNS-07   | Be-Transflow chip                              | 45 |
| <b>Research Instruments - Software Solutions</b>        |  |    |
| SSFT-OXY  | OxyGEN software                                | 46 |
| FLUIGENT-SDK  | Software Development Kit                       | 47 |
| <b>Research Instruments - Tubing &amp; Fitting Kits</b> |  |    |
| CTQ-KIT-HP-MFCS   | MFCST <sup>TM</sup> High Pressure Kit          | 48 |
| CTQ-KIT-LP-MFCS   | MFCST <sup>TM</sup> Low Pressure Kit           | 48 |
| LU-SPK-0002   | LineUp <sup>TM</sup> Supply Kit                | 49 |
| LU-C2C-0001   | LineUp <sup>TM</sup> Chain-To-Chain            | 49 |
| ELUPSWKIT1  | LineUp <sup>TM</sup> P-SWITCH Kit              | 49 |
| EIPSKIT   | PRESSURE UNIT Kit                              | 50 |
| CTQ-KIT-LQ-XS   | FLOW UNIT XS Kit                               | 50 |
| CTQ-KIT-LQ  | FLOW UNIT S and M 1/32"OD Kit                  | 50 |
| CTQ-KIT-FU2   | FLOW UNIT S and M 1/16"OD Kit                  | 51 |
| CTQ-KIT-HQ  | FLOW UNIT L Kit                                | 51 |
| CTQ-KIT-XL  | FLOW UNIT XL Kit                               | 51 |
| CTQ-KIT-2SW2  | 2-SWITCH <sup>TM</sup> Kit                     | 52 |
| CTQ-KIT-ESSMSW003                                       | M-SWITCH <sup>TM</sup> Kit                     | 52 |
| CTQ-KIT-LSW   | L-SWITCH <sup>TM</sup> Kit                     | 52 |
| CTQ-KIT-PCAP2   | P-CAP 2 mL Kit                                 | 53 |
| CTQ-KIT-PCAP15  | P-CAP 15 mL Kit                                | 53 |
| CTQ-KIT-PCAP50  | P-CAP 50 mL Kit                                | 53 |
| CTQ-KIT-F4C   | Fluiwell 4C Kit                                | 54 |
| CTQ-KIT-FIC15   | Fluiwell 15 mL Kit                             | 54 |
| CTQ-KIT-FIC50   | Fluiwell 50 mL Kit                             | 54 |
| CTQ-KIT-BC  | Bottle-CAP Kit                                 | 54 |

**Systems & Platforms - Automated Systems**

|                       |  |    |
|-----------------------|--|----|
| <b>O-OMISC-PTF</b>    | Omi, automated organ-on-chip platform (single channel)     | 56 |
| <b>O-OMIDC-PTF</b>    | Omi, automated organ-on-chip platform (double channel)     | 56 |
| <b>O-OMI-SA</b>       | Omi standalone   | 56 |
| <b>O-OMI-TAB</b>      | Omi tablet   | 56 |
| <b>O-OMI-CART</b>     | Set of 3 Omi Cartridges (sterile package)                  | 56 |
| <b>O-OMI-LRCC</b>     | Set of 3 Omi Low resistance adaptors (sterile package)     | 56 |
| <b>O-OMI-HRCC</b>     | Set of 3 Omi high resistance adaptors (sterile package)    | 56 |
| <b>O-OMI-DEV</b>      | Omi device   | 56 |
| <b>O-OMI-CTK</b>      | OMI tubing and fitting kit                                 | 56 |
| <b>O-OMI-AF-CTK</b>   | OMI air filter and pneumatic kit                           | 56 |
| <b>O-OMI-SK</b>       | Power supply kit (socket depends on the region)            | 56 |
| <b>O-OMI-LU</b>       |  | 56 |
| <b>OAR-L-02-2SW</b>   | Aria L Single Output                                       | 58 |
| <b>OAR-L-02-2SW-W</b> | Aria L Single Output + Installation + 1 year warranty ext. | 58 |
| <b>OAR-M-02-2SW</b>   | Aria M Single Output                                       | 58 |
| <b>OAR-M-02-2SW-W</b> | Aria M Single Output + Installation + 1 year warranty ext. | 58 |
| <b>OAR-L-02-MSW</b>   | Aria L Serial Output                                       | 58 |
| <b>OAR-L-02-MSW-W</b> | Aria L Serial Output + Installation + 1 year warranty ext. | 58 |
| <b>OAR-M-02-MSW</b>   | Aria M Serial Output                                       | 58 |
| <b>OAR-M-02-MSW-W</b> | Aria M Serial Output + Installation + 1 year warranty ext. | 58 |

**Systems & Platforms - Microfluidic Application Platforms**

|                     |  |    |
|---------------------|--|----|
| <b>O-DE-STD-PTF</b> | Complex Emulsion Production Platform           | 60 |
| <b>O-FACS1-PTF</b>  | Cell Encapsulation Platform, RayDrop 30-70-45  | 62 |
| <b>O-FACS2-PTF</b>  | Cell Encapsulation Platform, RayDrop 30-70-60  | 62 |
| <b>O-FACS3-PTF</b>  | Cell Encapsulation Platform, RayDrop 60-120-60 | 62 |
| <b>O-FACS4-PTF</b>  | Cell Encapsulation Platform, RayDrop 60-120-90 | 62 |

**Microfluidic Packs - Microfluidic Application Packs**

|                         |   |    |
|-------------------------|---|----|
| <b>O-MULTILABEL-PCK</b> | Immunostaining Pack                           | 65 |
| <b>O-RECIR-PCK</b>      | Recirculation Pack                            | 65 |
| <b>OOC-ST-001</b>       | Cell Perfusion Pack                           | 66 |
| <b>OOC-HT-001</b>       | High Throughput Cell Perfusion Pack           | 66 |
| <b>ESORT-PCK01</b>      | Microfluidic Cell Sorting Pack                | 67 |
| <b>ODROPSEQFPCK</b>     | Drop-Seq Pack                                 | 68 |
| <b>O-SE-ALG-PCK</b>     | Alginate Beads Generation Pack                | 70 |
| <b>O-MIX-LIPO-PCK</b>   | Nanoparticle Production Pack                  | 72 |
| <b>IDPPL01</b>          | PLGA Microparticle Production Standard Pack   | 74 |
| <b>O-SE-PLGAAP-PCK</b>  | PLGA Microparticle Production Automation Pack | 74 |
| <b>O-DE-STD-PCK</b>     | Double Emulsion Generation Pack               | 76 |

**Microfluidic Packs - Starter Packs**

|                        |   |    |
|------------------------|---|----|
| <b>DROPPACK-01</b>     | Microfluidic Droplet Pack                                     | 78 |
| <b>SEDUC-COFLOW</b>    | Educational Package First Level Co-Flow                       | 79 |
| <b>SEDUC-RESIST</b>    | Educational Package First Level Resistance                    | 80 |
| <b>SEDUC-RESITCOFL</b> | Educational Package Second Level Co-Flow & Resistance         | 80 |
| <b>SEDUC-DROPLET</b>   | Educational Package Full Course Co-Flow, Resistance, Droplets | 80 |

**Systems & Packs - Tubing & Fitting Kits**

|                   |                        |    |
|-------------------|------------------------|----|
| <b>2SW-KIT-AR</b> | Aria Single Output Kit | 91 |
| <b>MSW-KIT-AR</b> | Aria Serial Output Kit | 58 |

|  |  |    |
|--|--|----|
| <b>CTQ-KIT-AR</b>  | Aria Replacement Tubing Kit                          | 81 |
| <b>ODROPSEQCTQ</b>   | Drop-Seq Kit   | 81 |
| <b>O-DE-STD-CTK</b>  | Double Emulsion Generation Kit                       | 82 |
| <b>O-MIX-LIPO-CTK</b>  | Nanoparticle Production Kit                          | 83 |
| <b>IDPPLC1</b>   | PLGA Microparticle Production Station Standard Kit   | 83 |
| <b>IDPPLC2</b>   | PLGA Microparticle Production Station Automation Kit | 83 |
| <b>O-SE-ALG-CTK</b>  | Alginate Beads Generation Kit                        | 84 |
| <b>IDPRDC1</b>   | RayDrop Single Emulsion Device Kit                   | 84 |
| <b>Fluigent Industrial - Microfluidic OEM Equipment</b>            |  |    |
| <b>ICPX3451</b>  | PX 345 mbar  | 89 |
| <b>ICPX11</b>  | PX 1000 mbar   | 89 |
| <b>ICPX21</b>  | PX 2000 mbar   | 89 |
| <b>ICPXV1</b>  | PX -600 mbar   | 89 |
| <b>ICPXV2</b>  | PX -750 mbar   | 89 |
| <b>ISRX21</b>  | RX Pressure Source                                   | 93 |
| <b>IFSXS1</b>  | FS series XS   | 92 |
| <b>IFSS1</b>   | FS series S  | 92 |
| <b>FLU-OEM-M+</b>  | FS series M+   | 92 |
| <b>FLU-OEM-L+</b>  | FS series L+   | 92 |
| <b>INIFS</b>   | Non-Invasive Flow Sensor                             | 90 |
| <b>IFSM1</b>   | FS series M  | 92 |
| <b>IFSL1</b>   | FS series L  | 92 |
| <b>IFSXL1</b>  | FS series XL   | 92 |
| <b>Fluigent Industrial - Customizable OEM Flow Control Modules</b> |  |    |
| <b>INT-FOEM-3</b>  | Integration board FOEM 3 slots                       | 87 |
| <b>INT-FOEM-4</b>  | Integration board FOEM 4 slots                       | 87 |
| <b>INT-FOEM-EXT-1</b>  | Extension board FOEM 1 slot                          | 87 |
| <b>INT-FOEM-EXT-2</b>  | Extension board FOEM 2 slot                          | 87 |
| <b>INT-FOEM-EXT-3</b>  | Extension board FOEM 3 slot                          | 87 |
| <b>INT-FOEM-EXT-4</b>  | Extension board FOEM 4 slot                          | 87 |
| <b>PRM-FOEM-0025</b>   | Pressure Module FOEM 25 mbar                         | 87 |
| <b>PRM-FOEM-0069</b>   | Pressure Module FOEM 69 mbar                         | 87 |
| <b>PRM-FOEM-0345</b>   | Pressure Module FOEM 345 mbar                        | 87 |
| <b>PRM-FOEM-1000</b>   | Pressure Module FOEM 1000 mbar                       | 87 |
| <b>PRM-FOEM-2000</b>   | Pressure Module FOEM 2000 mbar                       | 87 |
| <b>PRM-FOEM-7000</b>   | Pressure Module FOEM 7000 mbar                       | 87 |
| <b>PRM-FOEM-N025</b>   | Pressure Module FOEM -25 mbar                        | 87 |
| <b>PRM-FOEM-N069</b>   | Pressure Module FOEM -69 mbar                        | 87 |
| <b>PRM-FOEM-N345</b>   | Pressure Module FOEM -345 mbar                       | 87 |
| <b>PRM-FOEM-N800</b>   | Pressure Module FOEM -800 mbar                       | 87 |
| <b>PRM-FOEM-PP</b>   | Push Pull FOEM                                       | 87 |
| <b>SWM-FOEM-4</b>  | Switch module FOEM                                   | 87 |

# PRODUCT GUIDE

## Do you prefer...

Automating your protocol

Live monitoring &amp; control

Customize &amp; develop own application

Go PC-Free

## Software

Use OxyGEN software

Use OxyGEN software

Use the SDK

## Valves

Add a SWITCH EZ

Add an M-SWITCH™

Add an L-SWITCH™

Add a 2-SWITCH™

Add a P-SWITCH

## Do you need to switch between multiple solutions over the course of the experiment?

Sequential injection or Multiple distribution

Recirculation or Sample injection

Sorting or Fluid redirection

Pneumatic or Quake valve actuation

## Do you have a pressure source in your laboratory?

No, with MFCST™ series

No, with LineUp™ series

Yes

## Pressure source

Ask for integrated pressure source

Add an independent FLPG+ or Compact Pressure Source

### How many solutions/ liquids/fluids do you want to handle in your system?

5 to 8

4

1 to 3

I want to extend freely

9 to 32

### Pressure system

Go with MFCST<sup>TM</sup>-EXGo with MFCST<sup>TM</sup>-EZGo with LineUp<sup>TM</sup> series(Use LineUp<sup>TM</sup> P-SWITCH)

### Sensors

Add FLOW UNIT

No additional sensor

Add PRESSURE UNIT

### Do you work with:

Flow rate

Pressure only

Extended pressure  
detection

REQUEST A QUOTE  
or  
VISIT OUR WEBSHOP

# NOTES





## CONTACT

WESTERN EUROPE  
& OTHER AREAS

FLUIGENT  
FRANCE

EASTERN EUROPE

FLUIGENT

DEUTSCHLAND GmbH

NORTH AMERICA

FLUIGENT INC.

USA

Customer Support

fluigent.com  
06 77 77 77 77

[www.fluigent.com](http://www.fluigent.com)

Find us on   

FEB 20