


Purus Gas Systems



Top: Brandon Tarr, Product Development Engineer

Bottom: Steve Kozel, U.S. Sales & Distribution Manager





for **more** info

Questions about which carrier gas purifier to use?

Call 800-356-1688 or 814-353-1300, ext. 4, or contact your Restek representative to discuss your application with our technical service chemists.

Why do I need to use traps and where should I install them?

Carrier gas must contain less than 1ppm of oxygen, water vapor, or any other trace contaminant, to prevent column degradation, shortened column lifetime, and increased stationary phase bleed. Contaminants cause ghost peaks to appear during temperature programming and degrade the validity of analytical data. Make-up gas also should be contaminant-free, or baseline fluctuations and excessive detector noise can occur; detector gases should be free of water and hydrocarbons, or excessive baseline noise can result. Gas purifiers remove these contaminants from gas sources, thereby improving system performance.

Moisture Removal

Moisture in carrier gas lines will prematurely degrade oxygen and hydrocarbon traps and increase detector noise (particularly with ECDs). As a precaution, we highly recommend installing a moisture trap before the hydrocarbon and oxygen traps on all carrier gas lines. Our favorite trap is the Super-Clean™ Ultra-High Capacity Moisture Filter (cat.# 22028, pg. 234).

Hydrocarbon Removal

Use a hydrocarbon trap if your gas has a potential source of hydrocarbon contaminants (e.g., an oil pump in an air compressor) or if you suspect you are observing carrier gas ghost peaks. Install the hydrocarbon trap after the moisture trap, to prevent moisture from degrading the hydrocarbon-trapping ability of the activated carbon in the hydrocarbon trap. We recommend the Super-Clean™ Ultra-High Capacity Hydrocarbon Filter (cat.# 22030, pg. 234).

Oxygen Removal

Oxygen is a column killer. It is present even in UHP gases, as minute leaks at fittings allow oxygen to influx against the concentration gradient. There are many choices for oxygen removal—the Super-Clean™ Ultra-High Capacity Oxygen Filter (cat.# 22029, pg. 234) is popular with Restek chemists. Because oxygen can enter a gas line at any fitting, the oxygen trap should be the last connection before the gas line enters the chromatograph.

tech tip

Carrier Gas Purity

Carrier gas should contain less than 1ppm of oxygen, moisture, or other trace contaminants, to prevent column degradation, increase column lifetime, and decrease stationary phase bleed.

The expense of using high-purity gases in combination with carrier gas line purifiers will be offset by longer column lifetime and less GC maintenance.


did you know?

Trap replacement made simple!

Try the Super-Clean™ Triple Filter Carrier Gas Cleaning Kit (cat.# 22019, pg. 234) —it removes moisture, hydrocarbons, and oxygen in one easy-to-change, economical cartridge.

Purus Gas Systems

We know how important high-purity gas is to the success of your analysis. To provide you with the products and services you need to achieve a first-quality gas delivery system, Restek has created Purus Gas Systems. Purus Gas Systems are application-specific assemblies of the finest products and technical expertise for a systematic approach to gas handling. "Systematic" means that every product needed to deliver research purity gas is available from one source, at the level of quality you expect from Restek. We will not only supply products, we will work with you to design the best gas system for your application.



Restek provides the following total gas system solution:

- Restek gas purifiers provide cost-effective gas purity assurance.
- Restek stainless steel and copper tubing, precleaned and ready to use.
- Parker A-Lok® tube fittings consistently deliver high-quality performance.
- Extensive line of hand tools to make your work easier and faster, including Restek's Electronic Leak Detector.
- Gas generators for an uninterrupted supply of gas.
- Gas regulators for optimum line pressure control of all your chromatography gases.

Restek's Technical Service Team, 800-356-1688 or 814-353-1300, ext. 4, or your Restek representative, can answer your questions and provide system-design advice. From the gas source to your point of use, we offer the products and services that ensure the purity of your gas.


Thermal Gas Purifier



Single-Tube Model



Dual-Tube Model



Dimensions: 13" x 1/2"
(33 x 1.3 cm)

Thermal Gas Purifier*

- Removes oxygen, water, carbon monoxide, carbon dioxide, hydrocarbons (except methane) to ppb levels—pure enough for MS.
- Packed with reactor-grade, pure getter material for maximum efficiency and no contamination.
- Welded end-fittings on getter tubes eliminate leaks.
- Each tube has 12L oxygen and 35L water vapor capacity.
- Maximum flow: 1 liter/minute.

The getter material in Restek's re-engineered line of thermal gas purifiers reacts chemically with impurities in the carrier gas stream. Because the reaction is nonreversible, there is no possibility of contaminants breaking through the thermal gas purifier.

Gas purification is very economical when you use a thermal gas purifier. Getter tubes normally require changing only once per year; heavy use or very impure feed gas may require more frequent getter tube replacement.

| Restek Single-Tube Thermal Gas Purifier, 110 Volt (1 tube included)** | qty. | cat.# |
|---|------|-------|
| 1/8" Fittings | ea. | 21496 |
| 1/4" Fittings | ea. | 21497 |
| Restek Dual-Tube Thermal Gas Purifier, 110 Volt (2 tubes included)** | qty. | cat.# |
| 1/8" Fittings | ea. | 21498 |
| 1/4" Fittings | ea. | 21499 |
| Replacement Getter Tube | qty. | cat.# |
| 1/8" Fittings (Similar to Supelco part# 2-2396) | ea. | 21661 |
| 1/4" Fittings (Similar to Supelco part# 2-2398) | ea. | 21660 |

*Use with helium or nitrogen carrier gas only. Do not use with hydrogen, oxygen, or air—a violent reaction and/or fire will occur.

**Single-tube model dimensions: 15" x 7" x 6" (38 x 17.8 x 15.2 cm). Dual-tube model dimensions: 15" x 10" x 6" (38 x 25.4 x 15.2 cm).

did you know?

Getter material is the physical material placed in a vacuuming tube that binds residual gases to the surface, maintaining or increasing the vacuum.



Always know the status of the Thermal Gas Purifier getter tube: change tube when pressure drops.

Gas Pressure Gauge Kit

- Use an in-line pressure gauge to indicate when the Thermal Gas Purifier getter tube should be replaced.
- Includes 1/8" tee and 0–100 psi (0–689kPa) gauge.


| Description | qty. | cat.# |
|--|------|-------|
| In-line Gas Pressure Gauge Kit for Thermal Gas Purifiers | kit | 21657 |

Restek Gas Management System

- Removes moisture, hydrocarbons, and oxygen from carrier gas, extending column lifetime.
- Produces high-purity carrier gas for most applications.
- Includes one each: moisture, hydrocarbon, and indicating oxygen trap.
- Replacing traps is safe and easy.
- Maximum flow: 1 liter/minute.



Dimensions: 12" x 14" x 3"
(30.5 x 35.6 x 7.6 cm)



did you know?

The Restek Gas Management System removes water vapor (to 10ppb), hydrocarbons (to 0.1ppm), and oxygen (to less than 0.1ppm) with three traps housed in one unit.

| Description | qty. | cat.# |
|---|------|-------|
| Restek Gas Management System (includes fittings for 1/8" and 1/4" gas line) | ea. | 21999 |
| Replacement Traps | qty. | cat.# |
| High-Capacity Moisture Trap, 1/8" Nickel-Plated Brass Fittings | ea. | 21997 |
| Capillary-Grade Hydrocarbon Trap, 1/8" Nickel-Plated Brass Fittings | ea. | 21991 |
| Indicating Oxygen Trap, 1/8" Nickel-Plated Brass Fittings | ea. | 22010 |

Super-Clean™ Gas Filters

- High-purity output ensures 99.9999% pure gas (at max. flow of 2L/min.).
- “Quick connect” fittings for easy, leak-tight cartridge changes.
- Glass inside to prevent diffusion; polycarbonate housing outside for safety.

Designed for fast, simple cartridge changing

Cartridge systems make changing gas filters quick and easy, and Super-Clean™ gas filters are the latest in cartridge-style gas filtration. A base plate allows cartridges to be exchanged without introducing oxygen, moisture and hydrocarbons. Spring-loaded check valves seal when a filter is removed and open only when a new filter has been locked in place. There is no longer a need for loosening and tightening fittings every time a trap is changed, and your system will not become contaminated during the process.

With available 2- or 3-position base plates, you can purify all GC gas streams at one location. Figure 1 shows some possible filter cartridge combinations using these base plates. Any combination is possible because any Restek Super-Clean™ gas filter cartridge can be used with any Restek base plate.

High-purity output for improved sensitivity (Table I)

The Triple Gas Filter cartridge (cat.# 22020) is ideal for purifying carrier gas. This trap contains oxygen, moisture, and hydrocarbon scrubbers and indicators for oxygen and moisture in one cartridge. The purity of your carrier gas after flowing through the Triple Gas Filter is better than six-9s (99.9999% pure at max. flow of 2L/min.), which is ideal for sensitive mass spectrometry (MS) or ECD analyses, and for protecting your analytical columns against damage from contaminated carrier gas.

The Fuel Gas Filter cartridge (cat.# 22022) is perfect for purifying flame ionization detector (FID) fuel gases, removing both moisture and hydrocarbons. Using the Fuel Gas Filter for FID hydrogen and air will produce a stable baseline, improving overall reproducibility and sensitivity.

Figure 1 Filter cartridges can be configured for different applications.




Table I Each Super-Clean™ gas filter provides high-purity outlet gas.

| Type of Filter | Outlet Gas Quality (%) | Maximum Pressure/ Maximum Flow Rates | Use for: | Indicator Color Change | H ₂ O (g) | O ₂ (mL) | Capacity Hydrocarbons (g) | Estimated Lifetime (years) |
|--|------------------------|---|--------------------------------------|---|----------------------|---------------------|---------------------------|----------------------------|
| Moisture cat.# 22028 | >99.9999 | 11 bar 159psi/ 7 L/min. | Inert carrier gas Air Hydrogen | Yellow/orange to clear | 7.2 | — | — | >2 |
| Oxygen cat.# 22029 | >99.9999 | 11 bar 159psi/ 7 L/min. | Inert carrier gas | Green to grey | NA | 1000 | — | >2 |
| Hydrocarbons cat.# 22030 | >99.9999 | 11 bar 159psi/ 7L/min. | Inert carrier gas Air Hydrogen | No indicator | NA | — | 12 ³ | >2 |
| Fuel Gas ¹ cat.# 22022 | >99.9999 | 11 bar 159psi/ 7L/min. | Inert carrier gas Air Hydrogen | Yellow/orange to clear | 3.5 | — | 24 ³ | >1.5 |
| Triple ² cat.# 22020 | >99.9999 | 11 bar 159psi/ 7L/min. | Inert carrier gas | Yellow/orange to clear Green to grey | 1.8 | 500 | 4 ³ | >1 |
| Helium Specific ² cat.# 21982 | >99.9999 | 11 bar 159psi/ 7L/min. | Helium | Yellow/orange to clear Green to grey | 1.8 | 500 | — | >1 |

¹Removes hydrocarbons, moisture.

²Removes hydrocarbons, moisture, oxygen.



did you know?


All Super-Clean™ gas filter cartridges (except hydrocarbon filter cat.# 22030) feature easy-to-read indicators. The indicator code is shown on every trap so there is no confusion about when to replace it.

tech tip

Oxygen and Moisture Traps
Restek highly recommends oxygen and moisture traps for make-up gas when operating sensitive detectors such as electron capture detectors (ECD). The hydrogen reaction gas for sensitive electrolytic conductivity detectors (ELCD) also requires a hydrocarbon trap, to remove trace impurities.

Super-Clean™ gas filters are listed on pages 234 and 235.

Super-Clean™ Gas Filters and Base Plates



All traps measure:
10 7/8" x 1 3/4" (27 x 4.4 cm)
Each base plate unit
measures:
4" x 4" x 1 7/8"
(10.2 x 10.2 x 4.8 cm)

Super-Clean™ Gas Filter and Base Plate Kits

| Description | qty. | cat.# |
|--|------|-------|
| Carrier Gas Cleaning Kit (includes mounting base plate, 1/8" inlet/outlet fittings, and oxygen/moisture/hydrocarbon Triple Gas Filter) | kit | 22019 |
| Fuel Gas Purification Kit (includes mounting base plate, 1/8" inlet/outlet fittings, and hydrocarbon/moisture Fuel Gas Filter) | kit | 22021 |

Replacement Gas Filters

| Description | qty. | cat.# |
|---|------|-------|
| Replacement Triple Gas Filter (removes oxygen, moisture and hydrocarbons) | ea. | 22020 |
| Replacement Fuel Gas Filter (removes moisture and hydrocarbons) | ea. | 22022 |

Gas Filter Bundle Kit

- Ideal for use in combination with 3-position base plate—purchase separately.

| Description | qty. | cat.# |
|-----------------------|------|-------|
| Gas Filter Bundle Kit | kit | 22031 |

Super-Clean™ Ultra-High Capacity Gas Filters

| Description | qty. | cat.# |
|--|------|-------|
| Ultra-High Capacity Hydrocarbon Filter | ea. | 22030 |
| Ultra-High Capacity Moisture Filter | ea. | 22028 |
| Ultra-High Capacity Oxygen Filter | ea. | 22029 |

Helium-Specific Super-Clean™ Gas Filter and Kit

- Specifically designed for purification of helium in GC/MS systems.
- Traps are packed and conditioned using helium.
- Uses standard single-position base plate.

| Description | qty. | cat.# |
|---|------|-------|
| Helium-Specific Carrier Gas Cleaning Kit (includes mounting base plate, 1/8" inlet/outlet fittings, and oxygen/moisture/hydrocarbon Helium-Specific Filter) | kit | 21983 |
| Replacement Helium-Specific Gas Filter (removes oxygen, moisture and hydrocarbons) | ea. | 21982 |

Filter Base Plates

- Standard base plate fittings are 1/8". To adapt to 1/4", order 1/8" to 1/4" tube-end unions.
- Base plates fit all Super-Clean™ gas filters listed above.

| Description | qty. | Brass | new! Stainless Steel | |
|-----------------------------------|------|-------|----------------------|-------|
| | | cat.# | qty. | cat.# |
| Single-Position Filter Base Plate | ea. | 22025 | ea. | 22344 |
| 2-Position Filter Base Plate | ea. | 22026 | ea. | 22345 |
| 3-Position Filter Base Plate | ea. | 22027 | ea. | 22346 |

Wall Mounting Bracket

- Base plates can be mounted by using screws and the mounting holes on the base plate, or by using this optional wall mounting bracket.

| Description | qty. | cat.# |
|--|------|-------|
| Wall Mounting Bracket for Super-Clean™ Base Plates | ea. | 21984 |

Replacement O-Rings for Cartridge Base Plates

- Pack includes 10 large O-rings and 10 small O-rings.

| Description | qty. | cat.# |
|---|--------|-------|
| Replacement O-Rings for Cartridge Base Plates | 20-pk. | 22023 |

1/8-Inch to 1/4-Inch Tube-End Unions

- To adapt 1/8" Super-Clean™ base plate fittings to 1/4", use 1/8" to 1/4" tube-end unions.

| Description | qty. | Brass | Stainless Steel | |
|------------------------------|-------|-------|-----------------|-------|
| | | cat.# | qty. | cat.# |
| 1/8" to 1/4" Tube-End Unions | 5-pk. | 21833 | 2-pk. | 21933 |

Super-Clean™ Gas Trapping System for LC/MS**Quick-change cartridge system for removing hydrocarbon impurities from nitrogen**

- Changing filters is quick and easy.
- Up to 20L of hydrocarbon-free nitrogen per minute.

The Super-Clean™ Gas Trapping System is the latest technology in cartridge-style gas filtration for purifying nitrogen, and is ideal for use in LC/MS systems. The two-position base plate (installed in the gas line) allows cartridges to be exchanged without introducing oxygen into the system. Spring-loaded check valves seal when a cartridge is removed and open only when a new cartridge has been locked in place. There is no need for loosening and tightening fittings every time you change cartridges, and your system cannot become contaminated during the changing process.

To meet the high flow needs of the LC/MS system, the activated charcoal-filled cartridges are positioned and connected in parallel. The incoming gas stream is split equally between the cartridges, and the two streams are rejoined after purification but before the gas exits the base plate. This approach allows longer contact between the nitrogen and the adsorbent, ensuring higher gas purity and eliminating a potential source of contaminants to your analyses.

A handy date wheel, included with the system, indicates the cartridge installation date and the recommended replacement date. A special particle filter is included with the high flow base plate to be installed at the gas outlet. This exchangeable particle cup filter has a nominal pore size of 0.5µm and is recommended to be used to help remove potential dust from the charcoal, which could enter into the nebulizer gas stream and possibly damage the LC/MS system.

Table I Super-Clean™ Gas Filters provide high-purity outlet gas

| Type of filter | Max. Flow | Outlet Gas Quality % | Maximum Pressure | Estimated Lifetime |
|------------------------|-----------|----------------------|------------------|--------------------|
| Hydrocarbon (charcoal) | 20L/min. | 99.9999% | 11 bar/159psi | 3 to 6 months |

Super-Clean™ Gas Trapping System for LC/MS

| Description | qty. | cat.# |
|--|-------|-------|
| Super-Clean™ Gas-Trapping System (2-position base plate, 2 charcoal filters) | ea. | 22062 |
| 2-Position Base Plate (1/4" Fittings) for use with hydrocarbon filters (cat.# 22061) | ea. | 22060 |
| Replacement Hydrocarbon (Charcoal) Filters | 2-pk. | 22061 |



20L of purified
nitrogen per minute!



Stephanie Sunner
& Tim Herring

a plus 1 story

"West Coast companies often have to deal with East Coast companies closing by 2:00 pm Pacific time and grumpy people itching to leave. Not so with Restek! When I called at 2:55 pm (Pacific) on a Friday, Stephanie was kind enough to not only answer my questions, but forward me immediately to Tim at technical service so that I could catch him before he left - and then stayed late to help me set up an account. Tim deserves kudos, too, for not only being extremely friendly, but also being able to identify a part with no serial number and a label from a different company! Awesome service!"

Tim S., Research Scientist, a West Coast pharmaceutical company

also available

Looking for a nitrogen generator for your LC/MS?
Restek offers a full line of Parker LC/MS generators. See **page 245**.



Super-Clean™ Click-On Purification Gas Traps

Click-On Inline Super-Clean™ Purification Gas Traps

- High-purity output ensures 99.9999% pure gas.
- Click-On fittings for easy, leak-tight cartridge changes; brass or stainless steel, 1/4" or 1/8".
- Helium-Specific Triple Gas Trap is ideal for GC/MS.

Using the same features and benefits as the Super-Clean™ base-plates and filters (page 234), Click-On adaptor connectors allow cartridges to be exchanged without introducing oxygen, moisture, and hydrocarbons. Spring-loaded check valves seal when a filter is removed and open only when a new filter has been locked in place. There is no need for loosening and tightening fittings every time a gas trap is changed, and your system will not become contaminated during the process.

The Triple Gas Trap is ideal for purifying carrier gas—it contains oxygen, moisture, and hydrocarbon scrubbers in one cartridge.

The Fuel Gas Trap is ideal for purifying flame ionization detector (FID) fuel gases, removing both moisture and hydrocarbons.

The Helium-Specific Triple Trap is ideal for purifying helium in GC/MS systems. This gas trap is packed and purged under helium and contains oxygen, moisture, and hydrocarbon scrubbers in one cartridge.

Trap replacement depends on the quality of the incoming gas. Use the double connector and install an indicating cartridge after a trap to indicate when the trap should be replaced.

Click-On Traps measure:
8 1/2" x 1 1/4" (21.6 x 3.2cm)



| Filter Type | Gas Quality at Outlet | Maximum Pressure | Maximum Flow (L/min.) | Use For | H ₂ O (g) | Capacity O ₂ (mL) | Hydrocarbons ³ (g) | Estimated Lifetime (years) |
|----------------------------------|-----------------------|------------------|-----------------------|--|----------------------|------------------------------|-------------------------------|----------------------------|
| Moisture cat.#22467 | >99.9999 | 11 bar 160psi | 25 | Inert carrier gas, helium, air, H ₂ | 15 | NA | NA | >3 |
| Oxygen cat.#22468 | >99.9999 | 11 bar 160psi | 25 | Inert carrier gas | NA | 2000 | NA | >3 |
| Hydrocarbon cat.#22466 | >99.9999 | 11 bar 160psi | 25 | Inert carrier gas, helium, air, H ₂ | NA | NA | 24 | >3 |
| Fuel Gas ¹ cat.#22465 | >99.9999 | 11 bar 160psi | 25 | Inert carrier gas, helium, air, H ₂ | 7 | NA | 12 | >2 |
| Triple ² cat.#22464 | >99.9999 | 11 bar 160psi | 25 | Inert carrier gas | 4 | 1000 | 8 | >2 |

¹Removes hydrocarbons, moisture.

²Removes hydrocarbons, moisture, oxygen.

³As *n*-butane.

NOTE: Super-Clean™ Gas Filters are recommended for purifying noncorrosive gases with low concentrations of contaminants. The maximum concentration of oxygen in the incoming gas stream for oxygen purifiers is 0.5%.

Click-On Inline Super-Clean™ Gas Traps and Connector Kits

| Description | qty. | cat.# |
|---|------|-------|
| Carrier Gas Purification Kit, 1/8" Stainless Steel Includes (2) 1/8" SS connectors and (1) oxygen/moisture/hydrocarbon triple trap | kit | 22456 |
| Carrier Gas Purification Kit, 1/8" Brass Includes (2) 1/8" brass connectors and (1) oxygen/moisture/hydrocarbon triple trap | kit | 22457 |
| Carrier Gas Purification Kit, 1/4" Stainless Steel Includes (2) 1/4" SS connectors and (1) oxygen/moisture/hydrocarbon triple trap | kit | 22458 |
| Carrier Gas Purification Kit, 1/4" Brass Includes (2) 1/4" brass connectors and (1) oxygen/moisture/hydrocarbon triple trap | kit | 22459 |
| Fuel Gas Purification Kit, 1/8" Stainless Steel Includes (4) 1/8" SS connectors and (2) hydrocarbon/moisture traps | kit | 22460 |
| Fuel Gas Purification Kit, 1/8" Brass Includes (4) 1/8" brass connectors and (2) hydrocarbon/moisture traps | kit | 22461 |
| Fuel Gas Purification Kit, 1/4" Stainless Steel Includes (4) 1/4" SS connectors and (2) hydrocarbon/moisture traps | kit | 22462 |
| Fuel Gas Purification Kit, 1/4" Brass Includes (4) 1/4" brass connectors and (2) hydrocarbon/moisture traps | kit | 22463 |



Click-On Inline Super-Clean™ Replacement Gas Traps

| Description | qty. | cat.# |
|--|------|-------|
| Click-On Super-Clean™ Replacement Triple Trap (removes oxygen, moisture and hydrocarbons) | ea. | 22464 |
| Click-On Super-Clean™ Replacement Fuel Gas Trap (removes moisture and hydrocarbons) | ea. | 22465 |

Click-On Inline Super-Clean™ Ultra-High Capacity Gas Traps

| Description | qty. | cat.# |
|--------------------------------------|------|-------|
| Ultra-High Capacity Hydrocarbon Trap | ea. | 22466 |
| Ultra-High Capacity Moisture Trap | ea. | 22467 |
| Ultra-High Capacity Oxygen Trap | ea. | 22468 |

Helium-Specific Click-On Inline Super-Clean™ Gas Trap and Connector Kits

| Description | qty. | cat.# |
|--|------|-------|
| Kits | | |
| Helium-Specific Carrier Gas Cleaning Kit, 1/8" Stainless Steel Includes (2) 1/8" SS connectors and (1) oxygen/moisture/hydrocarbon helium-specific triple trap | kit | 22469 |
| Helium-Specific Carrier Gas Cleaning Kit, 1/8" Brass Includes (2) 1/8" brass connectors and (1) oxygen/moisture/hydrocarbon helium-specific triple trap | kit | 22470 |
| Helium-Specific Carrier Gas Cleaning Kit, 1/4" Stainless Steel Includes (2) 1/4" SS connectors and (1) oxygen/moisture/hydrocarbon helium-specific triple trap | kit | 22471 |
| Helium-Specific Carrier Gas Cleaning Kit, 1/4" Brass Includes (2) 1/4" brass connectors and (1) oxygen/moisture/hydrocarbon helium-specific triple trap | kit | 22472 |
| Replacement Trap | | |
| Helium-Specific Replacement Triple Trap (removes oxygen, moisture and hydrocarbons) | ea. | 22473 |

Click-On Inline Super-Clean™ Indicator

- Oxygen: Green to Grey
- Moisture: Beige to Clear

| Description | qty. | cat.# |
|--|------|-------|
| Click-On Inline Super-Clean™ Indicator (oxygen, moisture) | ea. | 22474 |

Click-On Inline Super-Clean™ Connectors

- Click-On connectors allow you to change traps quickly, without introducing oxygen into your system.

| Description | qty. | cat.# |
|--|-------|-------|
| 1/8" Brass Click-On Inline Super-Clean™ Connectors | 2-pk. | 22475 |
| 1/8" Stainless Steel Click-On Inline Super-Clean™ Connectors | 2-pk. | 22476 |
| 1/4" Brass Click-On Inline Super-Clean™ Connectors | 2-pk. | 22477 |
| 1/4" Stainless Steel Click-On Inline Super-Clean™ Connectors | 2-pk. | 22478 |

Click-On Inline Super-Clean™ Double Connector

- Connects any Click-On trap to a Click-On indicator.

| Description | qty. | cat.# |
|--|------|-------|
| Click-On Inline Super-Clean™ Double Connector, stainless steel | ea. | 22479 |

Wall-Mounting Clamps for Click-On Inline Super-Clean™ Gas Traps

| Description | qty. | cat.# |
|---|-------|-------|
| Wall-Mounting Clamps for Click-On Inline Super-Clean™ Gas Traps | 4-pk. | 22480 |

Replacement O-Rings for Click-On Inline Super-Clean™ Connectors

- Pack includes 10 large O-rings and 10 small O-rings.

| Description | qty. | cat.# |
|---|--------|-------|
| Replacement O-Rings for Click-On Inline Super-Clean™ Connectors | 20-pk. | 22481 |



did you know?

Helium-Specific Click-On Inline Super-Clean™ Gas Trap and Kits are designed specifically for purification of helium in GC/MS systems!



tech tip

Install an indicator after the Click-On inline gas filter so there is no confusion about when to replace the traps.



Each connector is 2 3/8" (6cm) in length.



Each double connector is 3" (8cm) in length.



Gas-Specific Purifier Modules



Dimensions: 21" x 1 1/2"
(53.3 x 3.8 cm)

VICI® Mat/Sen® Gas-Specific Purifier Modules

- Replace separate oxygen, moisture, and hydrocarbon traps with one multiple-bed purifier, specific for purifying helium, hydrogen, nitrogen, or air.
- Reduce gas impurities from ppm to low ppb levels.
- Decrease baseline noise and increase GC/MS sensitivity.
- Prepurged with the specified gas, to shorten downtime.

These gas-specific purifier modules offer dramatic reductions in contaminant levels—from ppm to levels that are below the limit of standard analytical detection. Performance is optimized by incorporating a multiple-bed format that progressively lowers concentrations of contaminants at each successive bed. VICI® Mat/Sen® purifiers are guaranteed to produce gases that are purer than 99.9999%, when supplied with gas of 99.995% purity, and are prepurged with the specified gas to speed conditioning. Purifier capacity is approximately four tanks of gas at 99.995% (50ppm) purity, and correspondingly longer for purer gases.

Specifications:

| | |
|--------------------------|-------------------|
| Length | 21" (52.3cm) |
| Diameter | 1.5" (3.8cm) |
| Maximum Inlet Pressure | 1000psi (6895kPa) |
| Maximum Recommended Flow | 500mL/min. |

| | |
|---|---|
| Pressure Drop from 120psi (827 kPa) inlet at 0-500mL/min.: | <0.20psi (1.4kPa) |
| End Fittings | compression, 1/8" or 1/4", stainless steel |
| Shipping Weight | 3.04 lb. (1300 g) |

Please Note: We recommend using an indicating oxygen trap (e.g., cat.# 22029, pg. 234) downstream from a VICI® Mat/Sen® purifier to continually ensure gas purity and indicate absolute change-out time.

it's a fact

The nitrogen module is excellent for LC/MS systems.

| Gas-Specific Purifier Module | Max Pressure | Compression Tube Fittings | | | |
|---------------------------------|------------------|---------------------------|-------|----------|-------|
| | | 1/4-inch | | 1/8-inch | |
| | | qty. | cat.# | qty. | cat.# |
| Helium Purifier Module | 1000psi, 6895kPa | ea. | 22600 | ea. | 22601 |
| Hydrogen Purifier Module | 1000psi, 6895kPa | ea. | 22602 | ea. | 22603 |
| Nitrogen Purifier Module* | 1000psi, 6895kPa | ea. | 22604 | ea. | 22605 |
| Air Purifier Module | 1000psi, 6895kPa | ea. | 22606 | ea. | 22607 |

*Warning: Do not use with nitrogen containing more than 500ppm of oxygen. If the oxygen level in the stream exceeds 500ppm, use an air purifier.

a plus 1 story

"We enjoy doing business with Restek. Their technical knowledge and willingness to back their products help us to maximize the performance of our chromatography instruments."

Jean-François Vergelin, Département de Seine et Marne,
Direction de l'Eau et de l'Environnement, Laboratoire
Département d'Analyse des eaux (Melun, France)

High-Capacity Indicating Oxygen Trap

- Indicator changes color from dark blue to black as oxygen & water are trapped.
- Lasts longer than three smaller traps.
- Use with all carrier gases.
- Ambient operating temperature, 100psi (689kPa) operating pressure.
- Built-in frit traps microparticles.
- Outlet gas purity: $O_2 < 0.1\text{ppm}$ when inlet does not exceed 15ppm.
 $H_2O < 0.5\text{ppm}$ when inlet does not exceed 10ppm.
- Includes cartridge housing and one cartridge.
- Maximum operating pressure: 150psi (1034kPa).
- Maximum flow: 16.5L/min.



Dimensions: 9¹/₄" x 2"
(23.5 x 5.1 cm)

| Description | qty. | cat.# |
|--|------|-------|
| High-Capacity Indicating Oxygen Trap, 1/8" Compression Tube Brass Fittings | ea. | 20624 |
| High-Capacity Indicating Oxygen Trap, 1/4" Compression Tube Brass Fittings | ea. | 20623 |
| Replacement Cartridge (fits 1/4" or 1/8" housing) | ea. | 20625 |
| Replacement O-Rings (5 small O-rings and 5 large O-rings) | kit | 22081 |

Indicating Oxygen Trap

- Indicator changes color from light green to grey as oxygen is trapped.
- Heavy-walled glass body prevents oxygen & water infusion.
- Prepurged for fast stabilization.
- 100psi (689kPa) maximum operating pressure.
- Reduces oxygen to 0.1ppm.

| Description | qty. | cat.# |
|---|------|-------|
| Indicating Oxygen Trap, 1/8" Nickel-Plated Brass Fittings | ea. | 22010 |
| Indicating Oxygen Trap, 1/4" Nickel-Plated Brass Fittings | ea. | 22011 |



Dimensions: 10" x 1¹/₄"
(25.4 x 3.2 cm)

High-Capacity Oxygen Trap

- Removes up to 600mg of oxygen or 2g of water.
- Long life—typically purifies more than five 200ft³ cylinders.
- Reduces oxygen to 15ppb.
- Maximum operating pressure: 250psi (1724kPa).
- Flow: 3L/min. @ 32psi (221kPa).

| Description | qty. | cat.# |
|--|------|-------|
| High-Capacity Oxygen Trap, 1/8" Nickel-Plated Brass Fittings | ea. | 20601 |
| High-Capacity Oxygen Trap, 1/4" Nickel-Plated Brass Fittings | ea. | 20600 |



Dimensions: 11" x 1¹/₂"
(27.9 x 3.8 cm)

Rechargeable Molecular Sieve S-Trap

- Traps water vapor; increases column and oxygen trap lifetime.
- Reduces baseline noise from sensitive detectors such as ECDs and mass spectrometers.
- Activated and ready to use.
- Reduces water to less than 1ppm.
- Fits in GC oven for easy thermal recharging.
- Maximum flow: 1L/min.

| Description | qty. | cat.# |
|--|------|-------|
| Rechargeable Molecular Sieve S-Trap, 1/8" Brass Fittings | ea. | 20686 |
| Rechargeable Molecular Sieve S-Trap, 1/4" Brass Fittings | ea. | 20685 |



Dimensions: 6³/₄" x 5⁵/₈"
(17.1 x 14.3 cm)

tech tip**Carrier Gas Purity**

Carrier gas should contain less than 1ppm of oxygen, moisture, or other trace contaminants, to prevent column degradation, increase column lifetime, and decrease stationary phase bleed.

The expense of using high-purity gases in combination with carrier gas line purifiers will be offset by longer column lifetime and less GC maintenance.

Moisture Traps, Hydrocarbon Traps




Dimensions: 11" x 1 1/2"
(27.9 x 3.8 cm)

High-Capacity Moisture Trap

- Purged with ultra-high-purity helium; ready to use.
- Reduces water to less than 15ppb.
- Maximum operating pressure: 250psi (1724kPa).
- Maximum flow: 1.25L/min.

| Description | qty. | cat.# |
|--|------|-------|
| High-Capacity Moisture Trap, 1/8" Nickel-Plated Brass Fittings | ea. | 21997 |
| High-Capacity Moisture Trap, 1/4" Nickel-Plated Brass Fittings | ea. | 20638 |



Dimensions: 13" x 2"
(33 x 5.1 cm)

Indicating Moisture Trap

- Reduces water to less than 10ppb; indicator changes color from blue to pink at 5% relative humidity.
- Prepurged for fast stabilization.
- Reduces noise from high-sensitivity detectors.
- Heavy-walled glass body prevents oxygen & water infusion.
- 40µm frit prevents microparticulate damage to needle valves and flow controllers.
- Maximum operating pressure: 100psi (689kPa).

| Description | qty. | cat.# |
|---|------|-------|
| Indicating Moisture Trap, 1/8" Nickel-Plated Brass Fittings | ea. | 22014 |
| Indicating Moisture Trap, 1/4" Nickel-Plated Brass Fittings | ea. | 22015 |



Dimensions: 11" x 1 1/2"
(27.9 x 3.8 cm)

Capillary-Grade Hydrocarbon Trap

- Packed with an extremely high surface area, baked coconut shell-based activated carbon.
- Purged with ultra-high-purity helium.
- Reduces organics to 0.1ppm (assuming 100ppm input).
- Maximum operating pressure: 250psi (1724kPa).

| Description | qty. | cat.# |
|---|------|-------|
| Capillary-Grade Hydrocarbon Trap, 1/8" Nickel-Plated Brass Fittings | ea. | 21991 |
| Capillary-Grade Hydrocarbon Trap, 1/4" Nickel-Plated Brass Fittings | ea. | 21992 |



Dimensions: 9 1/4" x 2 1/4"
(23.5 x 5.7 cm)

Refillable Hydrocarbon Trap

- Removes trace impurities from carrier gas: reduces organics to 0.1ppm (assuming 100ppm input).
- 20µm frit prevents gas contamination by purifier particles.
- Stops carrier gas contaminants from interfering with purge & trap systems.
- Refillable and rechargeable.
- Maximum operating pressure: 250psi (1724kPa).
- Maximum flow: 1.25L/min.

| Description | qty. | cat.# |
|--|------|-------|
| Refillable Hydrocarbon Trap, 1/8" Nickel-Plated Brass Fittings | ea. | 22012 |
| Refillable Hydrocarbon Trap, 1/4" Nickel-Plated Brass Fittings | ea. | 22013 |
| Carbon Refill (two recharges) | pint | 20626 |



Dimensions: 6 3/4" x 5 5/8"
(17.1 x 14.3 cm)

Hydrocarbon S-Trap

- Removes hydrocarbons and other contaminants from the gas stream.
- Reduces organics to 0.1ppm (assuming 100ppm input).
- Each trap individually activated to ensure maximum efficiency.
- Fits in GC oven for easy thermal recharging.
- Maximum operating pressure: 60psi (414kPa).

| Description | qty. | cat.# |
|---|------|-------|
| Hydrocarbon S-Trap, 1/8" Brass Fittings | ea. | 22016 |



Dimensions: 6" x 1 3/4"
(15.2 x 4.4 cm)

Indicating Hydrocarbon Trap for Air Compressors

- Pass compressed air from an oil-filled air compressor through this trap, to remove oil vapors and mist.
- Indicator changes from pale pink to deep pink.

| Description | qty. | cat.# |
|--|------|-------|
| Indicating Hydrocarbon Trap for Air Compressors, 1/8" Brass Fittings | ea. | 20637 |
| Indicating Hydrocarbon Trap for Air Compressors, 1/4" Brass Fittings | ea. | 20636 |

High-Capacity Split Vent Trap

- Reduces the release of hazardous materials from the capillary split vent into the lab.
- Lasts approximately one month or 1,500 injections.
- Includes connecting lines and mounting kit.

| Description | qty. | cat.# |
|-------------------------------|-------|-------|
| High-Capacity Split Vent Trap | ea. | 20698 |
| High-Capacity Split Vent Trap | 5-pk. | 20699 |



Dimensions: 6" x 1"
(15.2 x 2.5 cm)

ECD Vent Trap

- Reduces the release of hazardous materials from the ECD vent into the lab.
- Includes connecting lines and mounting kit.

| Description | qty. | cat.# |
|---------------|------|-------|
| ECD Vent Trap | ea. | 22017 |




Dimensions: 6" x 1"
(15.2 x 2.5 cm)

Replacement Chemical Traps for Agilent GCs

- Easy to install.
- Attach to same fittings as original equipment.
- Built-in frits retain fine particles; adsorbents remove both moisture and hydrocarbons.

| Description | Similar to Agilent part # | qty. | cat.# |
|---|----------------------------|--------|-------|
| A) Replacement Split Vent Trap for Agilent 6890/6850 GCs | G1544-80550 | ea. | 22820 |
| B) Replacement Chemical Trap for Agilent 5890 GCs | 05890-61260 | ea. | 21610 |
| C) Replacement Chemical Trap for Agilent 5880 GCs | 19362-60500 | ea. | 21158 |
| D) Split Vent Line for Agilent GCs (32-inch) (Includes all installation hardware.) | 19251-80525 G1544-20620 | 2-pk. | 22800 |
| E) O-rings for Agilent Trap Fittings | 5180-4181 | 25-pk. | 22064 |
| F) Optional Split Vent Trap Assembly for Agilent 6890/6850 GCs | G1544-60610 | kit | 23031 |
| G) Replacement Traps (2) and O-rings (4) | G1544-80530 | kit | 23032 |



Restek Electronic Leak Detector

- Reliable thermal conductivity leak detector.
- Responds to leaks in less than 2 seconds.
- Audible alarm plus LED readout.
- Auto zeros with the touch of a button.
- Built-in rechargeable 7.2-volt battery.

| Description | qty. | cat.# |
|---|------|-----------|
| Leak Detector with 110Volt Battery Charger | ea. | 22451 |
| Leak Detector with 220Volt European Battery Charger | ea. | 22451-EUR |
| Leak Detector with 220Volt UK Battery Charger | ea. | 22451-UK |



Small, compact unit—easy to hold and operate.

Leak Detector Accessory Kit

The kit includes an adaptor fitting to detect leaks in hard-to-reach locations, and a mounting bracket that can be affixed to the wall or GC.



Verify hard-to-reach leaks with the adaptor fitting.



Leak Detector Accessory Kit

| Description | qty. | cat.# |
|---|------|-------|
| Leak Detector Accessory Kit (adaptor fitting for probe, mounting bracket) | kit | 22453 |



Leak Detector is easily accessed when stored in the mounting bracket.



Parker Balston® Hydrogen Generators

- Proton Exchange Membrane (PEM) cell eliminates the need for liquid electrolytes.
- Reliably generate 99.9995% pure hydrogen—for better chromatography.
- Eliminates high-pressure cylinders—greater convenience and improved lab safety.
- Compact unit, requiring only one square foot of bench space.
- Quick and easy to service and maintain; unique display lighting changes color for easy status checks and water level indication.
- Comes with a set of universal power adapters for US, European, and Asian plug types.

Fuel-grade high purity hydrogen generators are safer alternatives to high-pressure gas cylinders. The new Proton Exchange Membrane (PEM) cell eliminates the use of liquid electrolytes with hydrogen generators. Deionized water is all that is required to generate hydrogen for weeks of continuous operation. With an output capacity of up to 510cc/minute, one generator can supply 99.9995% pure hydrogen for up to several FIDs. Based on cylinder gas savings alone, a hydrogen generator pays for itself in one or two years.

Produced and supported by an ISO 9001 registered organization, Parker Balston® hydrogen generators are the first built to meet the toughest laboratory standards in the world: CSA, UL, IEC 1010, and CE Mark. A great safety feature is the built-in sensing circuit, which shuts the generator down if a hydrogen leak is detected.

Specifications

| | |
|--------------------------|---|
| Purity: | 99.9995% pure hydrogen |
| Delivery Pressure: | 10-100psig \pm 1psig (69-689kPa \pm 7kPa) |
| Outlet Port: | 1/8" compression |
| Electrical Requirements: | 100-230VAC/50-60Hz |
| Physical Dimensions: | 17.12"h x 13.46"w x 17.95"d (43.48 x 34.19 x 45.6cm) |
| Shipping Weight: | 40 lbs. (18kg) dry |

new and improved!

Hydrogen PEM generators now come with a set of universal power adapters for US, European, and Asian plug types.



CE

| Description | Capacity | qty. | cat.# |
|------------------------------|------------|------|-------|
| Hydrogen Generator H2PEM-100 | 100cc/min. | ea. | 23065 |
| Hydrogen Generator H2PEM-165 | 165cc/min. | ea. | 23066 |
| Hydrogen Generator H2PEM-260 | 260cc/min. | ea. | 23067 |
| Hydrogen Generator H2PEM-510 | 510cc/min. | ea. | 23068 |

Replacement and Maintenance Components for Hydrogen Generators (for all models listed above)

| | | |
|--|-----|-------|
| Replacement Desiccant Cartridge for H2PEM Generators | ea. | 23069 |
| 6-Month Maintenance Kit for H2PEM Generators (Includes: 1 deionizer cartridge, 1 water filter, 3 environmental filters) | kit | 23070 |
| 24-Month Maintenance Kit for H2PEM Generators (Includes: 1 deionizer cartridge, 1 water filter, 3 environmental filters, 1 water level sensor, 1 water pump, and 1 desiccant cartridge) | kit | 23071 |

tech tip

Gas generators are an economical source of pure gases, and eliminate the inconvenience and danger of high-pressure cylinders.


free literature

Parker Balston® Hydrogen Generators

Download your free copy from www.restek.com.

Fast Facts
lit. cat.# 580053A

The combination of high-purity gas and gas purifying traps can save analytical time in the long run.



Parker Balston® Model FID-1000 and FID-2500 Gas Stations

- Single unit produces UHP zero air from house compressed air and 99.9995% pure hydrogen from deionized water.
- Ideal for supplying up to 5-6 FIDs.
- Eliminates inconvenient and dangerous gas cylinders.
- Silent operation, minimal operator attention required.

Parker Balston® Gas Stations provide both UHP grade hydrogen gas and zero grade air for flame ionization detectors. The system is specifically designed to supply gas to FIDs and to support flame thermionic and flame photometric detectors. The units produce zero air by purifying compressed air to a total hydrocarbon concentration of 0.1ppm or less (measured as methane).

The hydrogen generators produce hydrogen gas from deionized water, using the principle of electrolytic dissociation of water and hydrogen proton conduction through a proton exchange membrane cell.

Specifications

| | |
|---------------------------|---|
| Hydrogen Purity: | 99.9995% |
| Zero Air Purity: | FID-1000: < 0.1ppm total hydrocarbons as methane FID-2500: < 0.05ppm total hydrocarbons as methane |
| Max. Hydrogen Flow Rate: | FID-1000: 90cc/min. FID-2500: 250cc/min. |
| Max. Zero Air Flow Rate: | FID-1000: 1000cc/min. FID-2500: 2500cc/min. |
| Power: | 120VAC/amp, 60Hz, 400 watts |
| Hydrogen Outlet Pressure: | 60 psig (414kPa) |
| Zero Air Outlet Pressure: | 40-125 psig* (276-862kPa) |
| Inlet Connection: | 1/4" NPT (female) |
| Outlet: | 1/8" compression |
| Dimensions: | 16.5"h x 10.5"w x 17"d (42cm x 27cm x 43cm) |
| Weight: | 53 lbs. (24kg) |

*Zero air inlet requires minimum of 40psig (276kPa) compressed air pressure.



Produce zero air and pure hydrogen from one unit!

free literature**FID Gas Stations**


Download your free copy from www.restek.com.

Fast Facts
lit. cat.# 580051

ordering note

For international orders, please add the appropriate power cord suffix from the table below.

Compare baselines produced by a Parker Balston® FID Gas Station and bottled fuel air. The baseline produced by the Parker Balston® Generator is flat, with no fluctuations or peaks; the chromatogram from the bottled air fuel supply has many peaks ranging from 0.25 ppm to -0.25 ppm total hydrocarbons.



| Description | qty. | cat. # |
|--|------|--------|
| Model FID-1000 Gas Station (ideal for 1-2 FIDs) | ea. | 20177 |
| Model FID-2500 Gas Station (ideal for 5-6 FIDs) | ea. | 24913 |
| Replacement Components for FID Gas Stations | | |
| Resin Bed Cartridge for Hydrogen Generators in FID-1000 and FID-2500 Gas Stations | ea. | 24914 |
| Replacement Desiccant Cartridge | ea. | 21671 |
| FID Gas Station Maintenance Kit (Includes 1 desiccant cartridge, 1 resin bed cartridge, 1 filter cartridge) | ea. | 24915 |

International Power Cord Sets

Just add the proper suffix to the catalog number for the gas generator you are ordering.

| Location | qty. | cat.# suffix |
|---|------|--------------|
| United Kingdom (230VAC, 50/50Hz) | ea. | -550 |
| European (230VAC, 50/60Hz) | ea. | -551 |
| IEC Connector Only (230VAC, 50/60Hz) | ea. | -552 |
| Japanese (200VAC, 50/60Hz) | ea. | -556 |
| Japanese for Zero Air (100VAC, 50/60Hz) | ea. | -553 |
| Japanese for Hydrogen (100VAC, 50/60Hz) | ea. | -554 |
| Japanese for Nitrogen (100VAC, 50/60Hz) | ea. | -555 |



| Model | Number of FIDs* |
|------------|-----------------|
| 75-83NA | Up to 3 |
| HPZA-3500 | Up to 11 |
| HPZA-7000 | Up to 23 |
| HPZA-18000 | Up to 60 |
| HPZA-30000 | Up to 100 |

*based on a 300 cc/min. fuel air rate

Parker Balston® Zero Air Generators

- Turn in-house compressed air into ultra-pure air (<0.1ppm total hydrocarbons).
- Remove hydrocarbons to less than 0.1ppm by catalytic oxidation.
- Operate at 40 to 125psi (276-862kPa).
- Typical payback is less than one year, based on cylinder costs.
- Install easily and take up little bench space.
- Maintenance kits include a one year supply of prefilters and final filter.

Specifications

| | | |
|---|-------------------------------|---|
| Maximum Zero Air Flow Rate: | 75-83NA | 1 lpm |
| | HPZA-3500 | 3.5 lpm |
| | HPZA-7000 | 7 lpm |
| | HPZA-18000 | 18 lpm |
| | HPZA-30000 | 30 lpm |
| Outlet Hydrocarbon Concentration (as methane): | 75-83NA | < 0.1 ppm |
| | HPZA-30000 | < 0.1 ppm |
| | Other Models | < .05 ppm |
| Minimum/Maximum Inlet Air Pressure: | 40 psig/125 psig (276/862kPa) | |
| Maximum Inlet Hydrocarbon Concentration (as methane): | 100 ppm | |
| Pressure Drop at Maximum Flow Rate: | 4 psi (28kPa) differential | |
| Maximum Inlet Air Temperature: | 78°F (25°C) | |
| Inlet/Outlet Ports: | 1/4" NPT (female) | |
| Start-up Time to Specified Hydrocarbon Concentration: | 45 minutes | |
| Electrical Requirements: | 75-83NA | 120 VAC/60 Hz, 0.5 amps |
| | Other Models | 120 VAC/60 Hz, 3.5 amps |
| Dimensions: | 75-83NA | 12"h x 10"w x 3"d (30cm x 25cm x 8cm) |
| | Other Models | 16"h x 11"w x 13"d (42cm x 27cm x 34cm) |
| Shipping Weight: | 75-83NA | 7 lbs. (3 kg) |
| | Other Models | 41 lbs. (19 kg) |

| Zero Air Generator | Capacity | qty. | cat. # |
|---|-----------------|-------------|---------------|
| Zero Air Generator Model 75-83NA | 1000cc/min. | ea. | 20684 |
| Zero Air Generator Model 75-83NA with United Kingdom Power Cord | 1000cc/min. | ea. | 20684-550 |
| Zero Air Generator Model HPZA-3500 | 3500cc/min. | ea. | 20680 |
| Zero Air Generator Model HPZA-3500 with European Power Cord | 3500cc/min. | ea. | 20680-551 |
| Zero Air Generator Model HPZA-7000 | 7000cc/min. | ea. | 20681 |
| Zero Air Generator Model HPZA-18000 | 18,000cc/min. | ea. | 20682 |
| Zero Air Generator Model HPZA-30000 | 30,000cc/min. | ea. | 20683 |
| Maintenance Kits (includes a one-year supply of prefilters and final filter) | | qty. | cat. # |
| Maintenance Kit for Model 75-83NA | | kit | 21646 |
| Maintenance Kit for Models HPZA-3500, HPZA-7000, HPZA-18000, HPZA-30000 | | kit | 21647 |
| Replacement Catalyst Towers | Capacity | qty. | cat. # |
| Replacement Catalyst Tower for Model 75-83NA | 1000cc/min. | ea. | 22005 |
| Replacement Catalyst Tower for Model HPZA-3500 | 3500cc/min. | ea. | 22004 |
| Replacement Catalyst Tower for Model HPZA-7000 | 7000cc/min. | ea. | 22006 |
| Replacement Catalyst Tower for Model HPZA-18000 | 18,000cc/min. | ea. | 22007 |
| Replacement Catalyst Tower for Model HPZA-30000 | 30,000cc/min. | ea. | 22008 |

ordering note

For international orders, please add the appropriate power cord suffix from the table below.

free literature

Zero Air Generators

Download your free copy from www.restek.com.

Fast Facts

lit. cat.# 580050

International Power Cord Sets

Just add the proper suffix to the catalog number for the gas generator you are ordering.

| Location | qty. | cat.# suffix |
|---|------|--------------|
| United Kingdom (230VAC, 50/50Hz) | ea. | -550 |
| European (230VAC, 50/60Hz) | ea. | -551 |
| IEC Connector Only (230VAC, 50/60Hz) | ea. | -552 |
| Japanese (200VAC, 50/60Hz) | ea. | -556 |
| Japanese for Zero Air (100VAC, 50/60Hz) | ea. | -553 |
| Japanese for Hydrogen (100VAC, 50/60Hz) | ea. | -554 |
| Japanese for Nitrogen (100VAC, 50/60Hz) | ea. | -555 |

Parker Balston® Nitrogen Gas Generators

- Turn compressed air into ultra-pure nitrogen (up to 99.9995%).
- Flows from 1 to 75+ lpm.
- Require only a compressed air source and 110 volt AC power.
- Safe, reliable, low maintenance.
- Maintenance kits include replacement filters.
- N2-14 and N2-14A can be used for LC/MS.

Specifications

| | Model HPN2-1100 or UHPN2-1100 | Model HPN2-2000 | Model N2-14 or N2-14A |
|---|--|--|---|
| Maximum Nitrogen Flow Rate: | See Flow Table | 2 lpm | 78scfh** at 95% purity |
| Nitrogen Purity: | 99.9999% | 99.99% | 95.0%–99.5% |
| Maximum Nitrogen Outlet Pressure: | See Flow Table | 90 psig | |
| CO Concentration: | < 1.0 ppm | NA | |
| CO ₂ Concentration: | < 1 ppm | < 1 ppm | |
| O ₂ Concentration: | < 1 ppm | < 100 ppm | |
| H ₂ O Concentration: | ≤ 2 ppm | ≤ 2 ppm | |
| Hydrocarbon Concentration ¹ : | < 0.1 ppm | NA | |
| Argon Concentration ² : | 0.9% | 0.9% | |
| Atmospheric Dewpoint: | | | -58°F (-50°C) |
| Suspended Liquids: | | | None |
| Particles > 0.01µm: | | | None |
| Oxygen Analyzer: | | | Included with Model 75-720NA |
| Commercially Sterile: | | | Yes |
| Minimum/Maximum Inlet Pressure: | 60 psig/125 psig (414/862kPa) | 75 psig/120 psig (517/827kPa) | 60 psig/145 psig (414/1,000kPa) |
| Maximum Pressure Drop (99% N ₂ Purity, 125 psig): | | | 10 psig (69kPa) |
| Recommended Inlet Temperature: | ≤ 78°F (25°C) | ≤ 78°F (25°C) | ≤ 68°F (25°C) (Max.) |
| Ambient Operating Temperature: | 60°F–100°F (16°C–38°C) | 40°F–100°F (4°C–38°C) | 110°F (43°C) (Max.) |
| Maximum Air Consumption: | 42 lpm (1.5 scfm)* | 42 lpm (1.5 scfm)* | |
| Inlet Connection: | 1/4" NPT (female) | 1/4" NPT (female) | 1/4" NPT |
| Outlet Connection: | 1/4" compression | 1/4" NPT compression | 1/4" NPT |
| Electrical Requirements ³ : | 120 VAC/60 Hz | 120 VAC/60 Hz | N2-14: None N2-14A: 120 VAC/60 Hz/25 Watts |
| Dimensions: | 35"h x 12"w x 16"d (89cm x 30cm x 41cm) | 35"h x 12"w x 16"d (89cm x 30cm x 41cm) | 50"h x 16"w x 16"d (127cm x 41cm x 41cm) |
| Shipping Weight: | 115 lbs. (52 kg) | 115 lbs. (52 kg) | N2-14: 75 lbs. (34 kg) N2-14A: 80 lbs. (36 kg) |

Flow Table for Models HPN2-2000, HPN2-1100, and UHPN2-1100

| Inlet Air Pressure | Maximum Outlet Flow (cc/min.) | Maximum Outlet Pressure |
|--|-------------------------------|-------------------------|
| Models HPN2-1100 and UHPN2-1100 | | |
| 125 psig (862kPa) | 1100 | 85 psig (586kPa) |
| 110 psig (758kPa) | 1000 | 75 psig (517kPa) |
| 100 psig (689kPa) | 900 | 65 psig (448kPa) |
| 90 psig (621kPa) | 800 | 60 psig (414kPa) |
| 80 psig (552kPa) | 700 | 50 psig (345kPa) |
| 70 psig (483kPa) | 600 | 45 psig (310kPa) |
| 60 psig (414kPa) | 500 | 35 psig (241kPa) |
| Model HPN2-2000 | | |
| 75-120 psig (517-827kPa) | 2000 | 90 psig (621kPa) |

Nitrogen Generators for LC/MS or General Purpose

| | qty. | cat.# |
|---|------|-----------|
| Nitrogen Generator N2-14 (general purpose) 78 scfh** max. flow at 95% purity | ea. | 20677 |
| Nitrogen Generator N2-14 with European Power Cord | ea. | 20677-551 |
| Nitrogen Generator N2-14A (general purpose w/oxygen analyzer) 78 scfh** max. flow at 95% purity | ea. | 21652 |

Nitrogen Generators

| | qty. | cat.# |
|---|------|-------|
| Nitrogen Generator HPN2-2000 (high purity—99.99%) 2.0 lpm max. flow | ea. | 21654 |
| Nitrogen Generator HPN2-1100 (ultra-high purity—99.9995%) 1.1 lpm max. flow | ea. | 21653 |
| Nitrogen Generator UHPN2-1100 (ultra-high purity—99.9995%); [HC<0.1ppm] 1.1 lpm max. flow | ea. | 20697 |

Maintenance Kits

| | qty. | cat.# |
|---|------|-------|
| Maintenance Kit for Models N2-14, N2-14A, 75-72, 75-720NA | kit | 21648 |
| Maintenance Kit for Models HPN2-1100, HPN2-2000, 76-96, 76-92 | kit | 21649 |
| Maintenance Kit for Models UHPN2-1100, 76-94 | kit | 21655 |

*Standard cubic feet per minute.



Model: N2-14

¹Models HPN2-1100 and HPN2-2000 do not remove hydrocarbons.

²Purity specification for nitrogen does not include argon concentration.

³Power consumption is:
Model HPN2-1100 = 25 Watts
Model UHPN2-1100 = 700 Watts
Model HPN2-2000 = 25 Watts

free literature**Nitrogen Generators**

Download your free copy from
www.restek.com.

Fast Facts
lit. cat.# 580052

ordering note

For international orders,
please add the appropriate
power cord suffix from the
table on the previous page.

Introduction to Gas Pressure Regulators



Mike Zezzo
Mid-Atlantic States
Sales Representative
9+ years of service!

**General Purpose or Analytical?**

General-purpose gas regulators usually are best suited for applications involving gases that are less than 99.995% pure: pneumatically-actuated valves and autosamplers, blanketing, inert atmospheres, and any other application not directly integrated with analytical data production. General purpose gas regulators have nylon-reinforced neoprene diaphragms that provide very good pressure control but are prone to air and moisture diffusion and hydrocarbon off-gassing.

Analytical regulators are recommended for applications in which maintaining the purity of a gas or mixture is the overriding concern, i.e., for applications requiring gases that are greater than 99.995% pure. They are commonly used in analytical labs. Analytical regulators have stainless steel diaphragms for pressure control. Stainless steel is not subject to the diffusion and off-gassing associated with neoprene diaphragms, and is easily purged of atmospheric contaminants when put into service.

Dual- or Single-Stage?

Dual-stage gas regulators reduce the source pressure to outlet pressure in two steps. The first stage reduces the inlet pressure to about three times the maximum working pressure. Outlet pressure gas regulation is controlled by the second stage and is set through the use of an adjusting knob. This two-step gas regulation is highly recommended for services requiring a near constant delivery pressure as the source pressure decays, including chromatographic analyses.

Single-stage gas regulators perform the same function as dual-stage gas regulators, but in a single step-down from source pressure to outlet pressure. For this reason, the outlet pressure cannot be as accurately maintained as the source pressure decays. We highly recommend that single-stage gas regulators be used only in circumstances in which the operator can monitor and adjust the regulator as needed, when the regulator is supplied with a nearly constant source pressure, or when additional pressure regulation is supplied downstream.

Brass or Stainless Steel?

Analytical gas regulators made from brass bar stock provide optimum performance for most analytical applications. Brass provides excellent strength and cleanliness and the machined bar stock design has less dead volume than forged-body gas regulators, making purging of atmospheric contaminants faster and more assured.

Gas regulators with stainless steel bodies were designed for delivering corrosive gases that would be incompatible with brass. With the advent of semiconductor manufacturing and high sensitivity analytical techniques, stainless steel also has proven to be a better surface for removing “sticky” atmospheric contaminants that interfere with detectors downstream. Unless these gas regulators are used in an all-stainless-steel system that incorporates welded tubing and special fittings, and in which rigorous cleaning and proper gas management are practiced, the extra expense relative to brass is not justified.

Overview of Restek's Chrome-Plated Brass and Stainless Steel Body Ultra-High-Purity Gas Regulators

These regulators feature metal-to-metal seals throughout for long-term leak-tightness, and a metal diaphragm outlet valve ensures gas purity. Each regulator is helium leak-test-certifiable to 1×10^{-8} scc/sec. and is fully assembled and tested for your convenience. 100psig maximum delivery pressure supports EPC operation. Maximum inlet pressure is 3000psig. Chrome-plated brass bar stock construction minimizes dead volume. Stainless steel construction is more easily purged of atmospheric contaminants, and is more resistant to attack from dry corrosive gases.

Ultra-High-Purity Brass Body Gas Regulators

These regulators are the best choice when using ultra-high-purity carrier gas for sensitive GC applications using MS, PID, HID, or ECD detection methods. They feature reduced internal dead-volume, relative to stainless steel bodies. The stainless steel gas regulator diaphragm minimizes nonmetal parts, to help ensure ultra-high gas purity, and the metal valve diaphragm ensures leak-free shut-off. Oxidation-resistant chrome plating maintains a like-new appearance.

Dual-Stage Ultra-High-Purity Chrome-Plated Brass Gas Regulators

- Oxidation-resistant, chrome-plated.
- Most stable outlet pressure control throughout the life of a high-pressure gas cylinder.
- Secondary pressure regulation not needed.
- Most widely used regulator.
- Less internal volume than stainless steel gas regulators.

Outlet pressure: 0 to 100psig (0-689kPa)
 Outlet gauge: 30" – 0 to 200psig (0-1379kPa)
 Inlet gauge: 0 to 4000psig (0-27,579kPa)
 Outlet assembly: diaphragm valve, 1/4" tube fitting

| Fitting | qty. | cat.# |
|--|------|-------|
| CGA 580 (N ₂ , He, Ar) | ea. | 21667 |
| CGA 350 (H ₂ , P ₂) | ea. | 21668 |
| CGA 590 (Air) | ea. | 21669 |

Single-Stage Ultra-High-Purity Chrome-Plated Brass Gas Regulators

- Oxidation-resistant, chrome-plated.
- Use when there is secondary pressure regulation downstream.
- Identical gas purity protection as with our dual-stage gas regulators.

Outlet pressure: 0 to 100psig (0-689kPa)
 Outlet gauge: 30" – 0 to 200psig (0-1379kPa)
 Inlet gauge: 0 to 4000psig (0-27,579kPa)
 Outlet assembly: diaphragm valve, 1/4" tube fitting

| Fitting | qty. | cat.# |
|--|------|-------|
| CGA 580 (N ₂ , He, Ar) | ea. | 20646 |
| CGA 350 (H ₂ , P ₂) | ea. | 20647 |
| CGA 590 (Air) | ea. | 20648 |

Ultra-High-Purity Chrome-Plated Brass Line Gas Regulator

- Oxidation-resistant, chrome-plated.
- Use where you need to reduce the line pressure by 20psig (138kPa) or more.
- Same purity protection as high-pressure cylinder regulators.

Inlet connections: 1/4" FPT
 Outlet assembly: 1/4" FPT port

| Fitting | Outlet Gauge | Outlet Pressure | qty. | cat.# |
|------------------------|--------------------------------|----------------------|------|-------|
| 1/4" female NPT ports* | 30" - 0 to 100psig (0-689kPa) | 0-50psig (0-345kPa) | ea. | 21666 |
| 1/4" female NPT ports* | 30" - 0 to 200psig (0-1379kPa) | 0-100psig (0-689kPa) | ea. | 22452 |

*Order appropriate male connector, pipe-to-tube fittings.

Male Connector, Pipe-to-Tube Fittings

| Fitting Type | Size (inches) | Parker # | Similar to Swagelok® | qty. | Brass cat.# | Stainless Steel qty. | cat.# |
|------------------|-------------------|----------|----------------------|--------|-------------|----------------------|-------|
| Male Connector | 1/4" to 1/4" NPT | 4 MSC 4N | 400-1-4 | 10-pk. | 21842 | 2-pk. | 21942 |
| Male Connector | 1/8" to 1/4" NPT | 2 MSC 4N | 200-1-4 | 10-pk. | 21844 | 2-pk. | 21944 |
| Tube End Reducer | 1/4" tube to 1/8" | 4 TUR 2 | 200-R-4 | 5-pk. | 21834 | 2-pk. | 21934 |



male connector



tube end reducer

Stainless Steel Gas Pressure Regulators, Switchover Systems

Ultra-High-Purity Stainless Steel Body Gas Regulators

These regulators are the standard for ultra-high-purity and corrosion-resistant pressure regulation. They are more easily purged of atmospheric components, compared to brass gas regulators, making them ideal for the most demanding applications. Gas regulation performance is equal to our brass body gas regulators. For use in all-stainless steel systems where welded tubing and special fittings are used, and rigorous cleaning and proper gas management are practiced.



Dual-Stage Ultra-High-Purity Stainless Steel Gas Regulators

- Most stable outlet pressure control throughout the life of a high-pressure gas cylinder.
- Secondary pressure regulation not needed.

Outlet pressure: 0 to 100psig (0-689kPa)
 Outlet gauge: 30" – 0 to 200psig (0-1379kPa)
 Inlet gauge: 0 to 4000psig (0-27,579kPa)
 Outlet assembly: diaphragm valve, 1/4" tube fitting

| Fitting | qty. | cat.# |
|--|------|-------|
| CGA 580 (N ₂ , He, Ar) | ea. | 20662 |
| CGA 350 (H ₂ , P ₂) | ea. | 20663 |
| CGA 590 (Air) | ea. | 20664 |

Single-Stage Ultra-High-Purity Stainless Steel Gas Regulators

- Use when there is secondary pressure requirement downstream.
- Identical gas purity protection as with our dual-stage gas regulators.

Outlet pressure: 0 to 100psig (0-689kPa)
 Outlet gauge: 30" – 0 to 200psig (0-1379kPa)
 Inlet gauge: 0 to 4000psig (0-27,579kPa)
 Outlet assembly: diaphragm valve, 1/4" tube fitting

| Fitting | qty. | cat.# |
|--|------|-------|
| CGA 580 (N ₂ , He, Ar) | ea. | 20665 |
| CGA 350 (H ₂ , P ₂) | ea. | 20666 |
| CGA 590 (Air) | ea. | 20667 |



Critical Purity Automatic Switchover System for Noncorrosive Service

High-purity automatic switchover systems provide a continuous supply of high purity gas to the laboratory, process, or instrument, to allow you to replace a depleted gas source without interruption in the supply of gas. Continuous supply is achieved by setting the two regulators at slightly different pressures, to discharge one side of the system at a time. These models include flexible, all-stainless-steel pigtailed with armor casing. The CGA connection on each pigtail has a check valve in the gland to prevent contamination and minimize purging requirements.



Switching pressure:
 200psig/170psig
 (1379/1172kPa)
Inlet connections:
 flexible SS pigtailed (36")
Line regulator:
 0 to 100psig (0-689kPa)

| Brass Automatic Switchover System with Line Regulator | qty. | cat.# |
|---|------|----------|
| CGA 580 (N ₂ , He, Ar) | ea. | 20668580 |
| CGA 350 (H ₂ , P ₂) | ea. | 20668350 |
| CGA 590 (Air) | ea. | 20668590 |
| Stainless Steel Automatic Switchover System with Line Regulator | qty. | cat.# |
| CGA 580 (N ₂ , He, Ar) | ea. | 21593580 |
| CGA 350 (H ₂ , P ₂) | ea. | 21593350 |
| CGA 590 (Air) | ea. | 21593590 |

Protocol Station

The protocol station is designed for convenient wall mounting of high-purity gas regulators. Wall mounting provides ease of use, prevents gas regulator damage, and improves safety. Either chrome-plated brass or 316 stainless steel option is complete with a 3-foot, flexible, all-stainless-steel pigtail with armor casing. The CGA connection on the pigtail has an integral check valve in the gland to prevent contamination during cylinder changeout.



| Chrome-Plated Brass Protocol Station* | qty. | cat.# |
|--|------|-------|
| CGA 580 (N ₂ , He, Ar) | ea. | 21347 |
| CGA 350 (H ₂ , P ₂) | ea. | 21348 |
| CGA 590 (Air) | ea. | 21349 |
| Stainless Steel Protocol Station* | qty. | cat.# |
| CGA 580 (N ₂ , He, Ar) | ea. | 21327 |
| CGA 350 (H ₂ , P ₂) | ea. | 21328 |
| CGA 590 (Air) | ea. | 21329 |

*Pressure regulator not included. Order separately.

CGA Fittings

CGA-specified nuts and nipples with internal frit, 1/4-inch NPT nickel-plated brass.

| Description | qty. | cat. # |
|---|------|--------|
| CGA 580 Fitting, (N ₂ , He, Ar) | ea. | 21336 |
| CGA 350 Fitting, (H ₂ , P ₂) | ea. | 21337 |
| CGA 590 Fitting, (Air) | ea. | 21338 |

Flexible Stainless Steel Hoses

| Description | qty. | cat. # |
|--|------|--------|
| Flexible Stainless Steel Hose, 36", 1/4" Female NPT Fittings | ea. | 21339 |
| Flexible Stainless Steel Hose, 18", 1/4" Female NPT Fittings | ea. | 21340 |
| Flexible SS Hose, 36", with Stainless Steel CGA 580 | ea. | 21344 |
| Flexible SS Hose, 36", with Stainless Steel CGA 350 | ea. | 21345 |
| Flexible SS Hose, 36", with Stainless Steel CGA 590 | ea. | 21346 |

Flammable Gas Flash Arrestor—Factory Mutual Approved*

- Gas flow shuts off in the event of a flashback.
- Flame extinguished—flame front prevented from reaching the gas supply.
- No gas flow restriction under normal operating conditions.

| Description | qty. | cat. # |
|--|------|--------|
| Flammable Gas Flash Arrestor, Brass Body | ea. | 21334 |

*Approved for brass body servicing hydrogen, acetylene, propane, or natural gas only.

Backpressure Gas Regulator

Capillary GC inlet systems have backpressure regulators to maintain a constant upstream pressure and rapidly respond to catastrophic leaks. The 0–60psig (0–414kPa) operating range is sufficient to operate a 105m, 0.25mm ID column at its optimum flow rate.

| Description | qty. | cat. # |
|----------------------------|------|--------|
| Backpressure Gas Regulator | ea. | 20635 |

MINICYL Regulator

This compact general purpose regulator has many laboratory applications including air-drying glassware, sparging or evaporating solutions, and controlling pneumatic valves. It is constructed of lightweight aluminum with an elastomer diaphragm. Includes a 0–60psig (0–414kPa) gauge and either 1/8- or 1/4-inch tube fittings.

| Description | qty. | cat. # |
|----------------------------------|------|--------|
| MINICYL Regulator, 1/8" Fittings | ea. | 20610 |
| MINICYL Regulator, 1/4" Fittings | ea. | 20611 |

Cylinder Valve Wrench

This specially-designed wrench enables easy opening of cylinder valves that are fitted with a hand wheel. It is also suitable for removing difficult cylinder caps.

| Description | qty. | cat. # |
|-----------------------|------|--------|
| Cylinder Valve Wrench | ea. | 21321 |

Universal Cylinder Wrench

Use this versatile wrench for tightening gauges and gas regulator CGA fittings to cylinder outlets and pipe thread connections.

| Description | qty. | cat. # |
|---------------------------|------|--------|
| Universal Cylinder Wrench | ea. | 21322 |

Wall-Mount Gas Cylinder Holder

Prevent serious injuries! This holder is designed to prevent free-standing gas cylinders from tipping over and injuring personnel. The cast aluminum holder can be secured to a wall or the side of a work bench. The mount will secure a cylinder from 4–15 inches in diameter.

| Description | qty. | cat. # |
|--------------------------------|------|--------|
| Wall-Mount Gas Cylinder Holder | ea. | 21333 |



Instrument-Grade Tubing



ordering note

Required length in meters
x 3.2808 = length in feet.

Instrument-Grade Welded and Drawn 304 Grade Stainless Steel Tubing

Clean tubing is critical to ensure the delivery of pure gas to your instrument. Restek's stainless steel tubing is specially cleaned for inertness by using the procedure for processing our Silcosteel® and Siltek®-treated products, because scrupulously clean parts are a prerequisite for a quality coating.

| Tubing Dimensions | | Length (per-foot pricing on 26 feet or more) | | |
|-------------------|----------|--|-------------|-------------|
| | | 25 Feet | 26–100 Feet | > 100 Feet* |
| ID (in.) | OD (in.) | cat. # | cat. # | cat. # |
| 0.01" ID | 1/16" OD | 21500 | 21501 | 21502 |
| 0.02" ID | 1/16" OD | 21503 | 21504 | 21505 |
| 0.03" ID | 1/16" OD | 21506 | 21507 | 21508 |
| 0.04" ID | 1/16" OD | 21509 | 21510 | 21511 |
| 0.085" ID | 3/8" OD | 21512 | 21513 | 21514 |
| 0.21" ID | 3/4" OD | 21515 | 21516 | 21517 |

*The availability of long lengths (continuous lengths up to 500 feet) is subject to inventory constraints. Please inquire before ordering.

ordering note

An extra charge is applied for cutting and/or straightening stainless steel and/or copper tubing, calculated from the total number of pieces produced for each line item

precleaned

All tubing is precleaned and ready to use.

Cleaned Copper Tubing

- Adheres to ASTM B-280.
- Precleaned and ready to use.
- Use for plumbing GC systems.

| ID | OD | Wall | qty. | cat. # |
|--------|------|--------|--------|--------|
| 0.065" | 1/8" | 0.030" | 50 ft. | 21590 |
| 0.190" | 3/4" | 0.030" | 50 ft. | 21592 |



GC Installation Kit

This kit contains the tubing and fittings needed to add an additional GC to your lab bench. Also included are four 1/8-inch tees, so carrier, fuel, and other GC gases can be routed to the new inlet or detector from existing gas lines. Order additional parts, such as gas purifiers or regulators, separately to customize the GC installation to your specifications. Kit includes: one tubing cutter, one 1/8-inch x 1/4-inch reamer, one 7/16-inch wrench, one 1/2-inch wrench, four 1/8-inch brass tees, ten 1/8-inch brass nuts, ten brass front and back ferrules, and 50 feet (15.2 meters) of our instrument-grade cleaned 1/8-inch copper tubing.

| Description | qty. | cat. # |
|---------------------|------|--------|
| GC Installation Kit | kit | 21325 |

tech tip

Plumbing a GC

It is essential to use clean chromatographic-grade tubing to plumb a GC. Standard-grade tubing contains residual hydrocarbon contaminants from the drawing process. These contaminants migrate into the carrier gas stream, elevating background noise and causing down time.

1/16-Inch Tubing Cutter

- Produces square, smooth cuts in 1/16-inch tubing.
- Eliminates tubing distortion.
- Replaceable cutting wheel.

| Description | qty. | cat.# |
|----------------------------|-------|-------|
| 1/16" Tubing Cutter | ea. | 20192 |
| Replacement Cutting Wheels | 3-pk. | 20185 |



Ridgid® Tubing Cutter

- Excellent for cutting 1/8- or 1/4-inch metal tubing.
- Compact size is ideal for tight spaces.
- Replaceable cutting wheel.

| Description | qty. | cat.# |
|---|-------|-------|
| Ridgid® Tubing Cutter for 1/8" or 1/4" metal tubing | ea. | 23011 |
| Replacement Cutting Wheels | 2-pk. | 23012 |



Tubing Reamer

- Removes burrs from stainless steel tubing.
- For 1/4- or 1/8-inch tubing.
- Nonslip safety design.

| Description | qty. | cat.# |
|---------------|------|-------|
| Tubing Reamer | ea. | 20134 |



SSI TC-20 Tube Cutting Machine

- Cuts 1/16", 1/8", or 1/4" tubing with inside diameter as small as 0.008".
- Electrically operated bench-top model.
- Handy dressing tool on the swing arm removes burrs and reams tubing.
- Voltage selectable 110–120/220–240 volts, 50–60Hz.*

| Description | qty. | cat.# |
|--------------------------------|-------|-------|
| SSI Tubing Cutter Machine | ea. | 23029 |
| SSI Replacement Cutting Wheels | 3-pk. | 23030 |

*Unit shipped set for 110–120 operating voltage. Switch to 220–240 volts by using alternate fuse and power cord (included).



CE

8" x 6 1/4" x 4 1/4"
(20.3 x 15.9 x 10.8 cm)
Weight: 11 lbs. (5.0 kg)

Tubing Dressing Tool

Same tool as included with the SSI tube cutting machine.

| Description | qty. | cat.# |
|---|------|-------|
| 1/16" Tubing Dressing Tool | ea. | 20188 |
| Replacement Insert for 1/16" Tubing Dressing Tool | ea. | 20189 |
| 1/8" Tubing Dressing Tool | ea. | 20190 |
| Replacement Insert for 1/8" Tubing Dressing Tool | ea. | 20191 |



Tubing Bender

- Bends 1/8-inch, 3/16-inch, or 1/4-inch tubing.
- Accurate left-hand, right-hand, or offset bends.

| Description | qty. | cat.# |
|---------------|------|-------|
| Tubing Bender | ea. | 23009 |



1/16-Inch Tubing Cutting Pliers

- Ideal for cutting 1/16-inch tubing.
- Cuts quickly, reducing distortion.
- Cuts clean, eliminating need for deburring.

| Description | qty. | cat.# |
|-----------------------------|------|-------|
| 1/16" Tubing Cutting Pliers | ea. | 20193 |



Tube & Pipe Fittings

also available

For Siltek®/Sulfinert® and Silcosteel®-CR-treated fittings, see pages 392-393.

Tube & Pipe Fittings

The Instrumentation Group of Parker Balston® Corporation designs and manufactures a top-quality line of components and systems for use in process instrumentation, semiconductor manufacturing, and analytical instrumentation. Parker's product quality and delivery have made them a world-wide leader—and this is the level of quality and service Restek wants to deliver to you.

| Fitting Type | Parker # | Similar to Swagelok® # | Size | qty. | Brass cat. # | 316 Grade Stainless Steel qty. | 316 Grade Stainless Steel cat. # |
|---|----------|------------------------|---------------|--------|--------------|--------------------------------|----------------------------------|
| nut  | 1 Nu 1 | 102-1 | 1/16" | 20-pk. | 21800 | 5-pk. | 21900 |
| | 2 Nu 2 | 202-1 | 1/8" | 40-pk. | 21801 | 10-pk. | 21901 |
| | 4 Nu 4 | 402-1 | 1/4" | 40-pk. | 21802 | 10-pk. | 21902 |
| front ferrule  | 1 FF 1 | 103-1 | 1/16" | 20-pk. | 21803 | 10-pk. | 21903 |
| | 2 FF 2 | 203-1 | 1/8" | 40-pk. | 21804 | 20-pk. | 21904 |
| | 4 FF 4 | 403-1 | 1/4" | 40-pk. | 21805 | 20-pk. | 21905 |
| back ferrule  | 1 BF 1 | 104-1 | 1/16" | 20-pk. | 21806 | 10-pk. | 21906 |
| | 2 BF 2 | 204-1 | 1/8" | 40-pk. | 21807 | 20-pk. | 21907 |
| | 4 BF 4 | 404-1 | 1/4" | 40-pk. | 21808 | 20-pk. | 21908 |
| nut & ferrule set  | — | — | 1/16" | 10-pk. | 21809 | 2-pk. | 21909 |
| | — | — | 1/8" | 20-pk. | 21810 | 5-pk. | 21910 |
| | — | — | 1/4" | 20-pk. | 21811 | 5-pk. | 21911 |
| plug  | 1 BLP 1 | 100-P | 1/16" | 5-pk. | 21815 | 2-pk. | 21915 |
| | 2 BLP 2 | 200-P | 1/8" | 10-pk. | 21816 | 4-pk. | 21916 |
| | 4 BLP 4 | 400-P | 1/4" | 10-pk. | 21817 | 4-pk. | 21917 |
| union  | 1 SC 1 | 100-6 | 1/16" | 3-pk. | 21818 | ea. | 21918 |
| | 2 SC 2 | 200-6 | 1/8" | 5-pk. | 21819 | 2-pk. | 21919 |
| | 4 SC 4 | 400-6 | 1/4" | 5-pk. | 21820 | 2-pk. | 21920 |
| reducing union  | 2 RU 1 | 200-6-1 | 1/8" to 1/16" | 5-pk. | 21823 | ea. | 21923 |
| | 4 RU 1 | 400-6-1 | 1/4" to 1/16" | 5-pk. | 21824 | 2-pk. | 21924 |
| | 4 RU 2 | 400-6-2 | 1/4" to 1/8" | 5-pk. | 21825 | 2-pk. | 21925 |

Parker's (A-Lok) two-piece ferrules and NPT fittings are ideal for installing new equipment, modifying existing instrumentation, or replacing worn connections.

| Size | Parker # | Similar to Swagelok® | qty. | Brass cat.# | 316 Grade Stainless Steel qty. | 316 Grade Stainless Steel cat.# | Fitting Type |
|---------------------------------|----------|----------------------|--------|-------------|--------------------------------|---------------------------------|--|
| 1/16" | 1 ET 1 | 100-3 | 2-pk. | 21826 | ea. | 21926 |  tee |
| 1/8" | 2 ET 2 | 200-3 | 2-pk. | 21827 | ea. | 21927 | |
| 1/4" | 4 ET 4 | 400-3 | 2-pk. | 21828 | ea. | 21928 | |
| 1/8" | 2 ECR 2 | 200-4 | 2-pk. | 21829 | ea. | 21929 |  cross |
| 1/4" | 4 ECR 4 | 400-4 | 2-pk. | 21830 | ea. | 21930 | |
| A 1/8" tube to B 1/16" | 2 TUR 1 | 100-R-2 | 5-pk. | 21831 | 2-pk. | 21931 |  A B tube end reducer |
| 1/4" tube to 1/16" | 4 TUR 1 | 100-R-4 | 5-pk. | 21832 | 2-pk. | 21932 | |
| 1/8" tube to 1/4" | 2 TUR 4 | 400-R-2 | 5-pk. | 21833 | 2-pk. | 21933 | |
| 1/4" tube to 1/8" | 4 TUR 2 | 200-R-4 | 5-pk. | 21834 | 2-pk. | 21934 | |
| 1/8" | 2 PC 2 | 201-PC | 5-pk. | 21835 | 2-pk. | 21935 |  port connector |
| 1/4" | 4 PC 4 | 401-PC | 10-pk. | 21836 | 2-pk. | 21936 | |
| 1/8" tube to 1/4" | 2 PC 4 | 401-PC-2 | 5-pk. | 21837 | 2-pk. | 21937 | |
| A 1/8" tube to B 1/8" NPT | 2 MSC 2N | 200-1-2 | 10-pk. | 21841 | 2-pk. | 21941 |  A B male connector |
| 1/4" to 1/4" NPT | 4 MSC 4N | 400-1-4 | 10-pk. | 21842 | 2-pk. | 21942 | |
| 1/16" to 1/8" NPT | 1 MSC 2N | 100-1-2 | 5-pk. | 21843 | 2-pk. | 21943 | |
| 1/8" to 1/4" NPT | 2 MSC 4N | 200-1-4 | 10-pk. | 21844 | 2-pk. | 21944 | |
| 1/4" to 1/8" NPT | 4 MSC 2N | 400-1-2 | 10-pk. | 21845 | 2-pk. | 21945 | |
| A 1/8" tube to B 1/8" NPT | 2 FSC 2N | 200-7-2 | 5-pk. | 21846 | 2-pk. | 21946 |  A B female connector |
| 1/4" to 1/4" NPT | 4 FSC 4N | 400-7-4 | 5-pk. | 21847 | 2-pk. | 21947 | |
| 1/4" to 1/8" NPT | 4 FSC 2N | 400-7-2 | 5-pk. | 21848 | 2-pk. | 21948 | |
| 1/8" male* | 2A-Q4VN | QC4D-200 | — | — | ea. | 21957 |  male & female quick couplings |
| 1/8" male | 2A-Q4P | QC4S-200 | — | — | ea. | 21958 | |
| 1/8" female* | 2A-Q4CN | QC4B-200 | — | — | ea. | 21959 | |
| 1/4" male* | 4A-Q4VN | QC4D-400 | — | — | ea. | 21960 | |
| 1/4" male | 4A-Q4P | QC4S-400 | — | — | ea. | 21961 | |
| 1/4" female* | 4A-Q4CN | QC4B-400 | — | — | ea. | 21962 | |

*Includes self-sealing shut-off valve.

Gas Valves, Leak Detector, GC Installation Kit



Toggle-valve—
instant on/off
control



Ball valve—leak-free
bidirectional sealing



Plug valve—leak-free at
wide temperature
swings



cat.# 22202

cat.# 22200



cat.# 22209
Optional knob
for accurate
control of low
flow



CE

Parker Balston® Shut-Off Gas Valves

Parker toggle gas valves are ideal for applications in which instant on/off gas control is necessary. They are rated to 200psig at 21°C and have a maximum operating temperature of 148°C. Ball valves have a floating ball to assist sealing and to reduce operating torque, and dual seats to provide leak-tight bidirectional sealing. They are rated to 1500psig at 21°C and have a maximum operating temperature of 177°C. Perfect for instrument supply lines, plug valves work well in any application requiring throttling or on/off operations. Parker plug valves remain leak-free even when subjected to wide temperature swings. They are rated to 3000psig at 21°C and have a maximum operating temperature of 205°C.

| Valve Type | 1/8" Brass | | 1/4" Brass | | 1/2" Stainless Steel | | 3/4" Stainless Steel | |
|------------|------------|-----|------------|-----|----------------------|-----|----------------------|-----|
| | cat.# | ea. | cat.# | ea. | cat.# | ea. | cat.# | ea. |
| Toggle | 22188 | | 22189 | | 22190 | | 22191 | |
| Ball | 22192 | | 22193 | | 22194 | | 22195 | |
| Plug | 22196 | | 22197 | | 22198 | | 22199 | |

Parker Balston® Precision Metering Gas Valves

Precision metering gas valves provide accurate, stable control of low gas and liquid flow rates. The valve stem threads do not contact the fluid stream, making these valves ideal for high-purity applications. The O-ring seal and stem threads are coated with a low vapor pressure, silicone lubricant for optimum performance. An optional vernier turn-counter knob allows repeatable flow settings for standardized operating processes.

| Valve Type | 1/8" Nickel-Plated Brass | | 1/4" Nickel-Plated Brass | | 1/2" Stainless Steel | | 3/4" Stainless Steel | |
|-----------------------|--------------------------|-----|--------------------------|-----|----------------------|-----|----------------------|-----|
| | cat.# | ea. | cat.# | ea. | cat.# | ea. | cat.# | ea. |
| Metering, straight | 22200 | | 22201 | | 22204 | | 22205 | |
| Metering, right angle | 22202 | | 22203 | | 22206 | | 22207 | |

Vernier Knob for Metering Valve: cat.# 22209 (ea.)

Restek Electronic Leak Detector

- Reliable thermal conductivity leak detector.
- Responds to leaks in less than 2 seconds.
- Audible alarm plus LED readout.
- Auto zeros with the touch of a button.
- Built-in rechargeable 7.2-volt battery.



Easy-to-clean probe
assembly.

| Description | qty. | cat.# |
|---|------|-----------|
| Leak Detector with 110Volt Battery Charger | ea. | 22451 |
| Leak Detector with 220Volt European Battery Charger | ea. | 22451-EUR |
| Leak Detector with 220Volt UK Battery Charger | ea. | 22451-UK |

Caution: The Restek Electronic Leak Detector is NOT designed for determining leaks of combustible gases. A combustible gas detector should be used for determining combustible gas leaks under any condition. The Restek Electronic Leak Detector may be used for determining trace amounts of hydrogen in a GC environment only.

Leak Detector Accessory Kit

The kit includes an adaptor fitting to detect leaks in hard-to-reach locations, and a mounting bracket that can be affixed to the wall or GC.

| Description | qty. | cat.# |
|---|------|-------|
| Leak Detector Accessory Kit (adaptor fitting for probe, mounting bracket) | kit | 22453 |

GC Installation Kit

This kit contains the tubing and fittings needed to add an additional GC to your lab bench. Also included are four 1/8-inch tees, so carrier, fuel, and other GC gases can be routed to the new inlet or detector from existing gas lines. Order additional parts, such as purifiers or regulators, separately to customize the GC installation to your specifications. Kit includes: one tubing cutter, one 1/8-inch x 1/4-inch reamer, one 7/16-inch wrench, one 1/2-inch wrench, four 1/8-inch brass tees, ten 1/8-inch brass nuts, ten brass front and back ferrules, and 50 feet (15.2 meters) of our instrument-grade cleaned 1/8-inch copper tubing.

| Description | qty. | cat.# |
|---------------------|------|-------|
| GC Installation Kit | kit | 21325 |

54-Piece Tool Kit

A must-have for every lab—all the tools you need in one place! Set comes with screwdrivers, pliers, wrenches, sockets, scissors, clamps, and more. Durable, zippered, notebook-style carrying case for easy transport.

| Description | qty. | cat.# |
|-------------------|------|-------|
| 54-Piece Tool Kit | kit | 23004 |

Plier Set

Includes 6-inch nose/side cutter, 6-inch wire cutter, and 6-inch adjusting pliers.

| Description | qty. | cat.# |
|-------------|------|-------|
| Plier Set | set | 23033 |

Metric 9 Piece Ball-Point Hex Key Set

Includes 9 metric hex keys (Allen wrenches): 1.5, 2, 2.5, 3, 4, 5, 6, 8, 10mm.

| Description | qty. | cat.# |
|---------------------------------------|------|-------|
| Metric 9 Piece Ball-Point Hex Key Set | set | 22999 |

12 Piece Ball-Point Hex Key Set

Includes 12 hex keys (Allen wrenches): .050", $\frac{1}{16}$ ", $\frac{5}{64}$ ", $\frac{3}{32}$ ", $\frac{7}{64}$ ", $\frac{1}{8}$ ", $\frac{9}{64}$ ", $\frac{5}{32}$ ", $\frac{3}{16}$ ", $\frac{7}{32}$ ", $\frac{1}{4}$ ", and $\frac{5}{16}$ ".

| Description | qty. | cat.# |
|---------------------------------|------|-------|
| 12 Piece Ball-Point Hex Key Set | set | 22998 |

Tubing Bender

- Bends $\frac{1}{8}$ -inch, $\frac{3}{16}$ -inch, or $\frac{1}{4}$ -inch tubing.
- Accurate left-hand, right-hand, or offset bends.

| Description | qty. | cat.# |
|---------------|------|-------|
| Tubing Bender | ea. | 23009 |

Torx® Screwdriver Set

- Set includes TR-10, TR-15, and TR-20.
- Ideal for performing routine maintenance on Agilent 6890 GCs.

| Description | qty. | cat.# |
|-----------------------|------|-------|
| Torx® Screwdriver Set | set | 23034 |

5-in-1 Magnetic Screwdriver

Magnetic power tip holds bits and screws securely.

| Description | qty. | cat.# |
|-----------------------------|------|-------|
| 5-in-1 Magnetic Screwdriver | set | 23002 |

Ratchet Wrenches

Easier to use in confined spaces, compared to adjustable wrenches.

| Description | qty. | cat.# |
|---------------------------------|------|-------|
| $\frac{3}{8}$ " Ratchet Wrench | ea. | 23005 |
| $\frac{1}{2}$ " Ratchet Wrench | ea. | 23006 |
| $\frac{7}{16}$ " Ratchet Wrench | ea. | 23007 |
| $\frac{9}{16}$ " Ratchet Wrench | ea. | 23008 |

Wrench Set

Includes 4-inch, 6-inch, and 8-inch adjustable wrenches.

| Description | qty. | cat.# |
|-------------|------|-------|
| Wrench Set | set | 23001 |

