

# Specialty Columns by Application



## Specialty columns

- Designed for specific classes of compounds.
- Specially deactivated columns.

## Many chromatography markets and applications represented

- Foods, Flavors, & Fragrances.
- Petrochemical.
- Clinical/Forensic.
- Pharmaceutical.
- Environmental.

## Unique stationary phases and applications

- Designed to help solve chromatographic challenges.



## Basic Compounds Analysis

**Rtx®-5 Amine** (low polarity phase; Crossbond® 5% diphenyl/95% dimethyl polysiloxane)

- Application-specific columns for amines and other basic compounds, including alkylamines, diamines, triamines, ethanolamines, and nitrogen-containing heterocyclics.
- Stable to 315°C.

Active basic compounds that otherwise require derivatization, or an alternative analytical technique, can be analyzed on an Rtx®-5 Amine column. The tubing surface is chemically altered to reduce tailing of basic compounds, eliminating the need for column priming. An Rtx®-5 Amine column is ideal for analyzing a wide variety of basic compounds, but breakthrough technology also allows the analysis of neutral compounds, adsorptive compounds with oxygen groups susceptible to hydrogen bonding, or even weakly acidic compounds such as phenols. Every Rtx®-5 Amine column is tested to ensure that it exceeds the requirements for analyzing ppm levels of amines, without priming, and to ensure low bleed at maximum operating temperature.

**Rtx®-5 Amine Columns** (fused silica)

(Crossbond® 5% diphenyl/95% dimethyl polysiloxane)

ID	df (µm)	temp. limits	15-Meter	30-Meter
0.25mm	0.25	-60 to 300/315°C	12320	12323
	0.50	-60 to 300/315°C	12335	12338
	1.00	-60 to 300/315°C	12350	12353
0.32mm	1.00	-60 to 300/315°C	12351	12354
	1.50	-60 to 290/305°C	12366	12369
0.53mm	1.00	-60 to 290/305°C	12352	12355
	3.00	-60 to 280/295°C	12382	12385

restek  
innovation!


similar **phase**

PTA-5

also **available**

See **page 58** for Rtx®-35 Amine columns.

### Excellent peak shapes for amines & phenols on an Rtx®-5 Amine column.



Column: Rtx®-5 Amine, 30m, 0.32mm ID, 1.0µm (cat.# 12354)  
 Inj.: 1.0µL split injection of amines and phenols in water  
 On-column conc.: 22ng  
 Oven temp.: 120°C to 220°C @ 10°C/min.  
 Inj./det. temp.: 305°C  
 Carrier gas: hydrogen  
 Linear velocity: 38cm/sec. set @ 120°C  
 FID sensitivity: 6.4 x 10<sup>-11</sup> AFS  
 Split ratio: 25:1

please **note**

We recommend using base-deactivated fused silica guard columns (**page 29**) and base-deactivated liners (Instrument Supplies section of this catalog) with Rtx®-5 Amine columns.

Table of Contents for  
**Applications**

see **pages 518-519**

## Basic Compounds Analysis

restek  
innovation!**Rtx®-35 Amine** (midpolarity phase; Crossbond® 35% diphenyl/65% dimethyl polysiloxane)


- Application-specific columns for amines and other basic compounds, including alkylamines, diamines, triamines, ethanolamines, and nitrogen-containing heterocyclics.
- Stable to 220°C.

Active basic compounds that otherwise require derivatization, or an alternative analytical technique, can be analyzed on an Rtx®-35 Amine column. The tubing surface is chemically altered to reduce tailing of basic compounds, eliminating the need for column priming. An Rtx®-35 Amine column is ideal for analyzing a wide variety of basic compounds, but breakthrough technology also allows the analysis of neutral compounds, adsorptive compounds with oxygen groups susceptible to hydrogen bonding. Every Rtx®-35 Amine column is tested to ensure that it meets the requirements for analyzing ppm levels of amines, without priming, and to ensure low bleed at maximum operating temperature.

**Rtx®-35 Amine Columns** (fused silica)

(Crossbond® 35% diphenyl/65% dimethyl polysiloxane)

ID	df (µm)	temp. limits	15-Meter	30-Meter
0.25mm	0.50	0 to 220°C	11335	11338
	1.00	0 to 220°C	11350	11353
0.32mm	1.00	0 to 220°C	11351	11354
	1.50	0 to 220°C	11366	11369
0.53mm	1.00	0 to 220°C	11352	11355
	3.00	0 to 220°C	11382	11385

**Sharp ethanolamine peaks, low bleed: Rtx®-35 Amine column.**

1. monoethanolamine
2. diethanolamine
3. triethyleneglycol monomethylether
4. triethanolamine

Column: Rtx®-35 Amine, 30m, 0.32mm ID, 1.0µm (cat.# 11354)  
 Sample: 500µg/mL ethanolamine standard in water  
 Inj.: 1.0µL split (split ratio 10:1),  
 cup splitter inlet liner (cat.# 20709)  
 Inj. temp.: 300°C  
 Carrier gas: helium, constant pressure  
 Linear velocity: 40cm/sec. @ 50°C  
 Oven temp.: 50°C (hold 0.50 min.) to 280°C @15°C/min.  
 Det.: FID @ 300°C

Table of Contents for  
**Applications**

see pages 518-519

## Basic Compounds Analysis

**Stabilwax®-DB** (polar phase; Crossbond® base-deactivated Carbowax® polyethylene glycol)

- Application-specific columns for underivatized amines and other basic compounds, including alkylamines, diamines, triamines, nitrogen-containing heterocyclics. No need for column priming.
- Temperature range: 40°C to 220°C.

Stabilwax®-DB columns reduce adsorption and improve responses for many basic compounds, without analyte derivatization or column priming. For different selectivity of basic compounds, or higher oven temperatures, use an Rtx®-5 Amine or Rtx®-35 Amine column.

Stabilwax®-DB is a bonded stationary phase, but avoid rinsing these columns with water or alcohols.

### Stabilwax®-DB Columns (fused silica)


(Crossbond® Carbowax® polyethylene glycol for amines and basic compounds)

ID	df (µm)	temp. limits	15-Meter	30-Meter	60-Meter
0.25mm	0.25	40 to 210/220°C	10820	10823	
	0.50	40 to 210/220°C		10838	
0.32mm	0.25	40 to 210/220°C	10821	10824	
	0.50	40 to 210/220°C		10839	
	1.00	40 to 210/220°C	10851	10854	10857
0.53mm	0.50	40 to 210/220°C		10840	
	1.00	40 to 210/220°C	10852	10855	10858
	1.50	40 to 210/220°C		10869	

similar **phases**

DB-CAM, Carbowax® Amine, CP Wax 51 for amines

### Hexamethylenediamine on a Stabilwax®-DB column.



Column: Stabilwax®-DB, 30m, 0.32mm ID, 0.25µm (cat.# 10824)  
 Inj.: 0.4µL direct injection of a neat hexamethylenediamine sample  
 On-column conc.: 10 to 1000ng/component  
 Oven temp.: 95°C (hold 6 min.) to 235°C @ 7°C/min. (hold 4 min.)  
 Inj./det. temp.: 250°C  
 Carrier gas: hydrogen  
 Linear velocity: 40cm/sec.  
 FID sensitivity: 2 x 10<sup>-11</sup> AFS

## Acidic Compounds Analysis



**Steve Allison**  
R&D Group Leader  
1+ year of service!

### Stabilwax®-DA (polar phase; Crossbond® acid-deactivated Carbowax® polyethylene glycol)

- Application-specific columns for free (underivatized) acids, some inorganic acids.
- Resistant to oxidative damage.
- Temperature range: 40°C to 250°C.
- Equivalent to USP G25, G35 phases.

Stabilwax®-DA bonded polyethylene glycol has an acidic functionality incorporated into the polymer structure. This permits analysis of acidic compounds without derivatization, significantly reduces adsorption of acids, and increases sample capacity for volatile free acids. Stabilwax®-DA columns last longer and give better peak shapes for high molecular weight acids. Some inorganic acids also chromatograph well on a Stabilwax®-DA column; the limitation is the volatility of the acidic compound.

### Stabilwax®-DA Columns (fused silica)

(Crossbond® Carbowax® polyethylene glycol for acidic compounds)

ID	df (µm)	temp. limits	15-Meter	30-Meter	60-Meter
0.25mm	0.10	40 to 250°C	11005	11008	11011
	0.25	40 to 250°C	11020	11023	11026
	0.50	40 to 250°C	11035	11038	11041
0.32mm	0.10	40 to 250°C	11006	11009	11012
	0.25	40 to 250°C	11021	11024	11027
	0.50	40 to 250°C	11036	11039	11042
	1.00	40 to 240/250°C	11051	11054	11057
0.53mm	0.10	40 to 250°C	11007	11010	11013
	0.25	40 to 250°C	11022	11025	11028
	0.50	40 to 250°C	11037	11040	11043
	1.00	40 to 240/250°C	11052	11055	11058
	1.50	40 to 230/240°C	11062	11065	11068

### similar phases

DB-FFAP, HP-FFAP, NUKOL,  
OV-351, CP-Wax 58 CB, FFAP


### crossbond® technology

reduces bleed, prolongs  
column lifetime, and allows  
rejuvenation through solvent  
rinsing.

### Underivatized alcoholic beverage acids and methyl esters on a Stabilwax®-DA column.

Peak List	Conc. (ppm)
1. ethyl octanoate	100
2. acetic acid	100
3. propionic acid	100
4. isobutyric acid	100
5. 3-decanol	50
6. ethyl decanoate	50
7. ethyl laurate	50
8. cis-lactone	100
9. 2-phenylethanol	50
10. trans-lactone	100
11. methyl myristate	50
12. ethyl myristate	50
13. octanoic acid	100
14. ethyl palmitate	50
15. decanoic acid	100
16. dodecanoic acid	100
17. vanillin	100

Column: Stabilwax®-DA, 30m, 0.18mm ID, 0.18µm (cat.# 550752)  
Inj.: 1µL splitless (hold 0.5 min.) at conc. shown in peak list, in ethyl acetate, 4mm ID splitless liner w/wool (cat.# 20814-202.1)  
Inj. temp.: 240°C  
Carrier gas: hydrogen  
Make-up gas: nitrogen  
Linear velocity: 28psi @ 240°C  
Oven temp.: 70°C to 240°C at 12°C/min. (hold 3 min.)  
Det.: FID



## Enantiomers Analysis

### Cyclodextrin Columns for Analyzing Many Chiral Compounds

By adding  $\beta$  or  $\gamma$  cyclodextrin to our bonded Rtx®-1701 stationary phase, we greatly enhance overall utility and column lifetime for our chiral columns, compared to columns that have pure cyclodextrin stationary phases. Separations of more than one hundred chiral compounds have been achieved using our unique DEX columns, and our columns continue to demonstrate stability after hundreds of temperature program cycles. Refer to the applications section of this catalog for examples, or call our Technical Service chemists or your Restek representative for assistance in matching a column to your chiral analysis.

#### Rt™- $\beta$ DEXm Columns (fused silica)

(permethylated beta cyclodextrin doped into 14% cyanopropylphenyl/86% dimethyl polysiloxane)

ID	df ( $\mu$ m)	temp. limits	30-Meter
0.25mm	0.25	40 to 230°C	13100
0.32mm	0.25	40 to 230°C	13101

Uses: General purpose chiral phase with many published applications.

#### Rt™- $\beta$ DEXsm Columns (fused silica)

(2,3-di-O-methyl-6-O-*tert*-butyl dimethylsilyl beta cyclodextrin doped into 14% cyanopropylphenyl/86% dimethyl polysiloxane)

ID	df ( $\mu$ m)	temp. limits	30-Meter
0.25mm	0.25	40 to 230°C	13105
0.32mm	0.25	40 to 230°C	13104

Uses: Excellent column for most chiral compounds in essential oils.

#### Rt™- $\beta$ DEXse Columns (fused silica)

(2,3-di-O-ethyl-6-O-*tert*-butyl dimethylsilyl beta cyclodextrin doped into 14% cyanopropylphenyl/86% dimethyl polysiloxane)

ID	df ( $\mu$ m)	temp. limits	30-Meter
0.25mm	0.25	40 to 230°C	13107
0.32mm	0.25	40 to 230°C	13106

Uses: Similar in performance to Rt™- $\beta$ DEXsm but provides better resolution for limonene, linalool, linalyl acetate, ethyl-2-methylbutyrate, 2,3-butane diol, and styrene oxides.

#### Rt™- $\beta$ DEXsp Columns (fused silica)

(2,3-di-O-propyl-6-O-*tert*-butyl dimethylsilyl beta cyclodextrin doped into 14% cyanopropylphenyl/86% dimethyl polysiloxane)

ID	df ( $\mu$ m)	temp. limits	30-Meter
0.25mm	0.25	40 to 230°C	13111
0.32mm	0.25	40 to 230°C	13110

Uses: Often useful in dual-column configurations, with the Rt™- $\beta$ DEXsm column, for complex enantiomeric separations.

#### Rt™- $\beta$ DEXsa Columns (fused silica)

(2,3-di-acetoxy-6-O-*tert*-butyl dimethylsilyl beta cyclodextrin doped into 14% cyanopropylphenyl/86% dimethyl polysiloxane)

ID	df ( $\mu$ m)	temp. limits	30-Meter
0.25mm	0.25	40 to 230°C	13109
0.32mm	0.25	40 to 230°C	13108

Uses: Unique selectivity for esters, lactones, and other fruit flavor components.

#### Rt™- $\beta$ DEXcst Columns (fused silica)

(Proprietary cyclodextrin material doped into 14% cyanopropylphenyl/86% dimethyl polysiloxane)

ID	df ( $\mu$ m)	temp. limits	30-Meter
0.25mm	0.25	40 to 230°C	13103
0.32mm	0.25	40 to 230°C	13102

Uses: Proprietary stationary phase, developed specifically for the fragrance industry. Also used for pharmaceutical applications.

#### Rt™- $\gamma$ DEXsa Columns (fused silica)

(2,3-di-acetoxy-6-O-*tert*-butyl dimethylsilyl gamma cyclodextrin doped into 14% cyanopropylphenyl/86% dimethyl polysiloxane)

ID	df ( $\mu$ m)	temp. limits	30-Meter
0.25mm	0.25	40 to 230°C	13113
0.32mm	0.25	40 to 230°C	13112

Uses: Larger organic molecules. Also useful for flavor compounds in fruit juices.

### free literature

#### A Guide to the Analysis of Chiral Compounds by GC

Many example chromatograms in our 24-page chiral analysis guide will help you find the best chiral columns for your applications.

Download your free copy from [www.restek.com](http://www.restek.com)

Technical Guide  
lit. cat.# 59889

### please note

Application-specific chiral column kits are available! See [page 62](#).

### free literature

#### Grape Flavor Analysis, Using an Rt™- $\gamma$ DEXsa GC Column

Download your free copy from [www.restek.com](http://www.restek.com).

Applications Note  
lit. cat.# 59553

#### GC Analysis of Chiral Flavor Compounds in Apple Juices, Using Rt™- $\beta$ DEXsm and Rt™- $\beta$ DEXse Columns

Download your free copy from [www.restek.com](http://www.restek.com).

Applications Note  
lit. cat.# 59546



**Doug Smith**  
R&D Technician  
10+ years of service!

## Enantiomers Analysis

**Convenient chiral column kits, based on sample type. Enantiomeric profile and confirmational identification of individual chiral compounds.**

- Broad range of columns and selectivities.
- Rugged, highly stable columns.
- Extended column lifetime.
- Convenience and cost savings.

For phase descriptions, or to order columns separately, see page 61.

### Pharmaceutical Chiral Column Kits

(fused silica)

Dimensions & Columns	cat.#
30m, 0.25mm ID, 0.25 $\mu$ m Rt <sup>TM</sup> - $\beta$ DEXcst & Rt <sup>TM</sup> - $\beta$ DEXsm columns	13190
30m, 0.32mm ID, 0.25 $\mu$ m Rt <sup>TM</sup> - $\beta$ DEXcst & Rt <sup>TM</sup> - $\beta$ DEXsm columns	13191

### Environmental Chiral Column Kits

(fused silica)

Dimensions & Columns	cat.#
30m, 0.25mm ID, 0.25 $\mu$ m Rt <sup>TM</sup> - $\beta$ DEXcst & Rt <sup>TM</sup> - $\beta$ DEXsm columns	13192
30m, 0.32mm ID, 0.25 $\mu$ m Rt <sup>TM</sup> - $\beta$ DEXcst & Rt <sup>TM</sup> - $\beta$ DEXsm columns	13193

### Juices Chiral Column Kits

(fused silica)

Dimensions & Columns	cat.#
30m, 0.25mm ID, 0.25 $\mu$ m Rt <sup>TM</sup> - $\beta$ DEXse, Rt <sup>TM</sup> - $\beta$ DEXsm & Rt <sup>TM</sup> - $\gamma$ DEXsa columns	13194
30m, 0.32mm ID, 0.25 $\mu$ m Rt <sup>TM</sup> - $\beta$ DEXse, Rt <sup>TM</sup> - $\beta$ DEXsm & Rt <sup>TM</sup> - $\gamma$ DEXsa columns	13195

### Essential Oils Chiral Column Kits

(fused silica)

Dimensions & Columns	cat.#
30m, 0.25mm ID, 0.25 $\mu$ m Rt <sup>TM</sup> - $\beta$ DEXsm, Rt <sup>TM</sup> - $\beta$ DEXse, Rt <sup>TM</sup> - $\beta$ DEXsa, & Rt <sup>TM</sup> - $\beta$ DEXsp columns	13196
30m, 0.32mm ID, 0.25 $\mu$ m Rt <sup>TM</sup> - $\beta$ DEXsm, Rt <sup>TM</sup> - $\beta$ DEXse, Rt <sup>TM</sup> - $\beta$ DEXsa, & Rt <sup>TM</sup> - $\beta$ DEXsp columns	13197

## please note

Application-specific chiral column kits offer convenience and cost savings.

## tech tip

**Chiral selectivity improves significantly by realizing lower elution temperatures.**

**This can be achieved by:**


Faster linear velocities (80cm/sec.) with hydrogen carrier gas.

Slower temperature ramp rates (1–2°C/min.).

Appropriate minimum operating temperature (40 or 60°C).

On-column concentrations of 50ng or less.

### Complex $\gamma$ -lactones mix resolved on an Rt<sup>TM</sup>- $\beta$ DEXcst column.



1. (+/-)- $\gamma$ -heptalactone
2. (+/-)- $\gamma$ -octalactone
3. (+/-)- $\gamma$ -nonalactone
4. (+/-)- $\gamma$ -decalactone
5. (+/-)- $\gamma$ -dodecalactone

Column: Rt<sup>TM</sup>- $\beta$ DEXcst, 30m, 0.32mm ID, 0.25 $\mu$ m (cat.# 13102)  
 Oven temp.: 60°C (hold 1 min.) to 200°C @ 1°C/min.  
 Carrier gas: hydrogen  
 Linear velocity: 40cm/sec. set @ 60°C  
 Det.: FID @ 220°C

Table of Contents for  
**Applications**

see pages 518-519

## Polyunsaturated FAME Analysis

### FAMEWAX (polar phase; Crossbond® polyethylene glycol)

- Application-specific columns for FAMEs, specially tested with a FAME mixture.
- Temperature range: 20°C to 250°C.

The elution order of polyunsaturated FAMEs on FAMEWAX columns is comparable to that on other Carbowax® columns, but baseline resolution is achieved in significantly less time.

### FAMEWAX Columns (fused silica)

(Crossbond® polyethylene glycol)

ID	df (μm)	temp. limits	30-Meter
0.25mm	0.25	20 to 250°C	12497
0.32mm	0.25	20 to 250°C	12498
0.53mm	0.50	20 to 250°C	12499


similar **phase**

Omegawax

### FAMEs (marine oil standard) on a FAMEWAX column.

Column: FAMEWAX, 30m, 0.32mm ID, 0.25μm (cat. # 12498)  
 Inj.: 1μL  
 Conc.: 10,000μg/mL in iso-octane (total FAMES; see breakdown in peak list)  
 Oven temp.: 195–240°C at 5°C/min., 1 min. hold  
 Inj./det. temp.: 250°C/275°C  
 Carrier gas: hydrogen  
 Flow rate: 3mL/min. (constant flow)  
 Split ratio: 100:1

Peak List	Conc. (μg/mL)
1. C14:0	600
2. C14:1	100
3. C16:0	1600
4. C16:1	500
5. C18:0	800
6. C18:1 (oleate)	1300
7. C18:1 (vaccenate)	400
8. C18:2	200
9. C18:3	200
10. C20:0	100
11. C20:1	900
12. C20:2	100
13. C20:4	300
14. C20:3	100
15. C20:5	1000
16. C22:0	100
17. C22:1	300
18. C24:0	100
19. C22:6	1200
20. C24:1	100



free literature

### Foods, Flavors, and Fragrances

Includes important analysis tips and chromatograms for analysis of fats and oils, carbohydrates, vitamins, amino acids, organic acids, preservatives, flavors and fragrances, essential oils, and chiral separations. Retention time indices and complete product listings for all relevant GC and HPLC products also are included.

Minicatalog  
 lit. cat. # 59260A

### Monitoring Volatile Compounds in Food Contact Packaging, Using Purge and Trap GC/MS and an Rtx®-5MS Capillary Column

Download your free copy from  
[www.restek.com](http://www.restek.com)

Applications Note  
 lit. cat. # 59348

## Flavor and Fragrance Compounds Analysis

### Rt™-CW20M F&F (polar phase; Carbowax® polyethylene glycol—not bonded)

- Application-specific columns for flavor and fragrance compounds, specially tested.
- True non-bonded Carbowax® 20M polarity.
- Temperature range: 20°C to 220°C.

### Rt™-CW20M F&F Columns (fused silica)

(Carbowax® polyethylene glycol)

ID	df (μm)	temp. limits	30-Meter	50-Meter
0.25mm	0.25	60 to 220°C	12523	
0.32mm	0.33	60 to 220°C		12539

similar **phases**

HP-20M, Carbowax® 20M

## Flavor &amp; Fragrance Compounds Analysis

please **note**

Our Rt™-CW20M F&F column is the perfect confirmation column for flavor and fragrance analysis. See page 63.

**Rtx®-1 F&F** (nonpolar phase; Crossbond® 100% dimethyl polysiloxane)

- Application-specific columns for flavor and fragrance compounds.
- Stable to 350°C.

Retention index libraries in the flavor and fragrance industry have been compiled from years of data and thousands of compounds. Any slight variation in column selectivity could render the column useless. Rtx®-1 F&F columns are tailored to match the selectivity required in the industry, while offering excellent thermal stability. Our stringent quality testing ensures column-to-column reproducibility and extended column lifetimes over conventional 100% dimethyl polysiloxane columns.

**Rtx®-1 F&F Columns** (fused silica)

(Crossbond® 100% dimethyl polysiloxane)

ID	df (µm)	temp. limits	15-Meter	30-Meter	50-Meter	60-Meter
HP-1	0.25mm	0.25	-60 to 330/350°C		18023	18026
		0.50	-60 to 330/350°C		18038	18041
		1.00	-60 to 320/340°C		18053	18056
0.32mm	0.25	0.25	-60 to 330/350°C		18024	18027
		0.50	-60 to 330/350°C		18039	18042
		1.00	-60 to 320°C		18054	18057
0.53mm	0.50	0.50	-60 to 310/330°C	18037	18040	18043
		1.00	-60 to 310/330°C	18052	18055	18058
		1.50	-60 to 310/330°C	18067	18070	18073

similar **phase**

HP-1

**Fused Silica Manufacturing Group**

Jason Fisher, Sheldon McMurtrie, Michele Richner, Tom Barone, Ken Kline, Linda Holden, Tom Gurecki, Jack Haesler, Aaron Decker, Carolyn Williams, Dale Lucas, Kelsea Miller



Raymond Ciampichini, Pete Rose, David Rhodes, Jackie Glasgow, Kim Shaffer, Marcia Fulton, Santina Newlen, David W. Rhodes, Henry Knepp, Russ Stewart, Jessica Andrus, Tim Wilson, Bob Langford

## Triglycerides in Foods Analysis

**Rtx®-65TG** (high polarity phase; Crossbond® 65% diphenyl/35% dimethyl polysiloxane)

- Application-specific columns, specially tested for triglycerides.
- Stable to 370°C.


The Rtx®-65TG phase resolves triglycerides by degree of unsaturation as well as by carbon number. Because of the chemistry required to achieve 370°C thermal stability, an Rtx®-65TG column should not be used for analyses of polar compounds.

### Rtx®-65TG Columns (fused silica)

(Crossbond® 65% diphenyl/35% dimethyl polysiloxane)

ID	df (μm)	temp. limits	15-Meter	30-Meter
0.25mm	0.10	40 to 370°C	17005	17008
0.32mm	0.10	40 to 370°C	17006	17009
0.53mm	0.10	40 to 370°C	17007	17010

### Rtx®-65TG Structure



save **money!**

Get six columns for the price of five. Call **800-356-1688, ext. 4**, or your Restek representative for details!

crossbond®  
**technology**


reduces bleed, prolongs column lifetime, and allows rejuvenation through solvent rinsing.

please **note**

Triglycerides are often injected via on-column injection. Use 0.53mm retention gaps and appropriate connectors.

- Vu2 Union® (page 222)
- MXT™-Union Connector Kits for Fused Silica (page 225)

### Sharp resolution of butter triglycerides on an Rtx®-65TG column.



## Get More!

Food, Flavor & Fragrance  
Related Articles Online

[www.restek.com/FFF](http://www.restek.com/FFF)

## Detailed Hydrocarbon Analysis (DHA)

**Rtx®-1PONA** (nonpolar phase; Crossbond® 100% dimethyl polysiloxane)

- Application-specific columns meet ASTM and CGSB requirements for detailed hydrocarbon analysis.
- Stable to 340°C.

please **note**

To achieve critical resolutions, a 5-meter tuning column is connected to the analytical column and adjusted to the needed length through a series of trial analyses.

The Rtx®-1PONA polymer was designed to offer the exact polarity necessary to resolve hydrocarbons in the specific order requested by petrochemical companies. In order to meet the demanding resolution and retention criteria of the American Society for Testing and Materials (ASTM) and the Canadian General Standards Board (CGSB), Restek has developed unique quality control tests and specifications for the Rtx®-1PONA column. The measured values for retention (*k*), efficiency (*n*), and stationary phase selectivity (RI) are controlled so that each column exceeds the requirements of the ASTM and CGSB methods.

**Rtx®-1PONA Column** (fused silica)

(Crossbond® 100% dimethyl polysiloxane—optimized for hydrocarbon analysis)


ID	df (μm)	temp. limits	100-Meter
0.25mm	0.50	-60 to 300/340°C	10195

similar **phases****Rtx®-5PONA Tuning Column** (fused silica)

(Crossbond® 5% diphenyl/95% dimethyl polysiloxane)

ID	df (μm)	temp. limits	5-Meter
0.25mm	1.0	-60 to 325°C	10196

Petrocol DH, DB-Petro, HP-PONA


**Sharp, symmetric peak for ethanol (gasoline oxygenate), using an Rtx®-1PONA column.\*\***

C5 efficiency: 613,596 total theoretical plates  
*k'* (C5): 0.489  
*tert*-butanol skewness: 1.25  
 Resolution (*tert*-butanol/2-methylbutene-2): 5.60

1. ethanol
2. C5
3. *tert*-butanol
4. 2-methylbutene-2

**Rtx®-1PONA column  
 produces near  
 symmetrical alcohol peaks!**

\*\*P=paraffins; O=olefins; N=naphthenes; A=aromatics.  
 In alternate terminology: paraffins & isoparaffins = alkanes;  
 naphthenes = cyclic alkanes; olefins = alkenes.

**Critical pairs of gasoline components resolved per ASTM specifications, using an Rtx®-1PONA column.**

Circles indicate critical pairs that must be resolved.

1. ethanol
2. C5
3. *tert*-butanol
4. 2-methylbutene-2
5. 2,3-dimethylbutane
6. methyl *tert*-butyl ether (MTBE)
7. C6
8. 1-methylcyclopentane
9. benzene
10. cyclohexane
11. 3-ethylpentane
12. 1-*tert*-2-dimethylcyclopentane
13. C7
14. 2,2,3-trimethylpentane
15. 2,3,3-trimethylpentane
16. toluene
17. C8
18. ethylbenzene
19. *p*-xylene
20. 2,3-dimethylheptane
21. C9
22. 5-methylnonane
23. 1,2-methylethylbenzene
24. C10
25. C11 (undecane)
26. 1,2,3,5-tetramethylbenzene
27. naphthalene
28. C12 (dodecane)
29. 1-methylnaphthalene
30. C13 (tridecane)

Column: Rtx®-1PONA, 100m, 0.25mm ID, 0.5μm (cat.# 10195) plus Rtx®-5PONA tuning column, 2.62m, 0.25mm ID, 1.0μm, connected via Press-Tight® connector (cat.# 20446)  
 Sample: custom detailed hydrocarbon analysis (DHA) mix, neat  
 Inj.: 0.01μL, split (split ratio 150:1), 4mm cup inlet liner (cat.# 20709)  
 Inj. temp.: 200°C  
 Carrier gas: helium, constant flow  
 Linear velocity: 28cm/sec. (2.3mL/min.)  
 Oven temp.: 5°C (hold 15 min.) to 50°C @ 5°C/min. (hold 50 min.) to 200°C @ 8°C/min. (hold 10 min.)  
 Det.: FID @ 250°C

## Simulated Distillation (C5-C44) Analysis

### Rtx®-2887 (nonpolar phase; Crossbond® 100% dimethyl polysiloxane)

- Application-specific columns for simulated distillation.
- Stable to 360°C.

Rtx®-2887 columns' stationary phase, column dimensions, and film thickness have been optimized to exceed the resolution and skewing factor requirements currently specified in ASTM method D2887. Each column is individually tested to guarantee a stable baseline with low bleed and reproducible retention times. The Crossbond® methyl silicone stationary phase has increased stability compared to packed columns, ensuring stable baselines and shorter conditioning times.

### Rtx®-2887 Column (fused silica)

(Crossbond® 100% dimethyl polysiloxane—for simulated distillation)

ID	df (μm)	temp. limits	10-Meter
0.53mm	2.65	-60 to 360°C	10199

### MXT®-2887 Column (Siltek® treated stainless steel)

(Crossbond® 100% dimethyl polysiloxane—for simulated distillation)

ID	df (μm)	temp. limits	10-Meter
0.53mm	2.65	-60 to 400°C	70199




**Jason Fisher**  
GC Column  
Manufacturing Supervisor

similar **phases**

DB-2887, Petrocol EX2887

### Negligible baseline rise for C5 to C44 hydrocarbons on an Rtx®-2887 column.



1. C5
2. C6
3. C7
4. C8
5. C9
6. C10
7. C11
8. C12
9. C14
10. C16
11. C18
12. C20
13. C24
14. C28
15. C32
16. C36
17. C40
18. C44

Column: Rtx®-2887, 10m, 0.53mm ID, 2.65μm (cat.# 10199)  
 Sample: 1μL direct injection of 0.01-0.1 wt. % C5 to C44 hydrocarbon standard in carbon disulfide  
 Inj. temp.: 360°C  
 Det. temp.: 360°C  
 Carrier gas: helium (constant flow)  
 Linear velocity: 15mL/min. (112cm/sec.)  
 Oven temp.: 35°C to 360°C @ 15°C/min. (hold 5 min.)

## Simulated Distillation (C44-C100) Analysis

### MXT®-1HT Sim Dist/MXT®-1 Sim Dist/MXT®-500 Sim Dist (nonpolar phases)

- Application-specific columns in unbreakable Siltek® treated stainless steel tubing meet all resolution criteria for high temperature simulated distillation (e.g., ASTM Method D2887 Extended).
- MXT®-1HT Sim Dist and MXT®-1 Sim Dist phases offer true methyl silicone polarity; MXT®-500 Sim Dist phase is a carborane siloxane polymer.
- Stable to 430°C.

Manufactured from Siltek® treated stainless steel tubing, MXT® columns are the most durable high temperature GC columns available. As outlined in ASTM Method D6352, high temperature simulated distillation requires a column that can withstand temperatures to 430°C. MXT®-1HT Sim Dist and MXT®-500 Sim Dist columns exhibit excellent peak shape and low bleed, even at 430°C! The unique MXT®-1HT Sim Dist methyl silicone polymer gives the correct retention time/boiling point curve. The MXT®-500 Sim Dist carborane siloxane polymer offers a slight shift in the calculated boiling range distribution for petroleum samples containing aromatic hydrocarbons.

### MXT®-1HT Sim Dist Column (Siltek® treated stainless steel)

ID	df (μm)	temp. limits	5-Meter
0.53mm	0.10	-60 to 430°C	70100

### similar phases

DB-1HT, HT-Simdist CB

### MXT®-1 Sim Dist Column (Siltek® treated stainless steel)

ID	df (μm)	temp. limits	6-Meter
0.53mm	0.15	-60 to 430°C	70101


### MXT®-500 Sim Dist Column (Siltek® treated stainless steel)

ID	df (μm)	temp. limits	6-Meter
0.53mm	0.15	-60 to 430°C	70104

### Polywax® Calibration Materials

Description	qty.	cat.#
Polywax® 655 calibration material	1g	36225
Polywax® 1000 calibration material	1g	36227

### C44-C130 hydrocarbons on an MXT®-1HT Sim Dist column.



Column: MXT®-1HT Sim Dist, 5m, 0.53mm ID, 0.10μm (cat.# 70100)  
 Sample: Polywax® 1000 (cat.# 36227)  
 Solvent: carbon disulfide  
 Sample size: 0.2μL  
 Instrument: Agilent 5890 GC w/GC Racer  
 Injector: on-column (track oven)  
 Carrier gas: hydrogen, constant pressure (1.0psi)  
 Det.: FID @ 430°C  
 Oven temp.: 40°C to 430°C @ 60°C/min. (hold 30 min.)

GC\_PC00543

## Aromatics & Oxygenates in Gasoline Analysis

### Rt™-TCEP (highly polar phase; 1,2,3-tris[2-cyanoethoxy]propane—not bonded)

- General purpose columns, ideal for aromatics and oxygenates in gasoline.
- Temperature range: 0°C to 135°C.

Most gasolines contain aliphatic hydrocarbons up to *n*-dodecane (C12). To improve identification of the aromatics and oxygenates, it is desirable to elute benzene after C11 and toluene after C12. The extremely polar Rt™-TCEP stationary phase provides a retention index for benzene greater than 1100 and permits the separation of alcohols and aromatics from the aliphatic constituents in gasoline.

Rt™-TCEP columns have the same high polarity as TCEP packed columns (precolumns in ASTM Method D4815 for the analysis of petroleum oxygenates), with the efficiency of a capillary column. The result is a column that can separate a wide variety of compounds with an elution pattern unattainable using other high polarity siloxanes.

The Rt™-TCEP column incorporates a nonbonded stationary phase coated on a surface specialized for enhanced polymer stability and extended column lifetime. Solvent rinsing should be avoided. Conditioning is necessary only if the column is to be used at temperatures near the maximum operating temperature.

### Rt™-TCEP Columns (fused silica)

(1,2,3-tris[2-cyanoethoxy]propane)

ID	df (μm)	temp. limits	30-Meter	60-Meter
0.25mm	0.40	0 to 135°C	10998	10999




**Brian Salisbury**  
ASG Chemist  
1+ year of service!

### similar phases

SPB-TCEP, CP-TCEP

### Petroleum oxygenates on an Rt™-TCEP column.



Column: Rt™-TCEP 60m, 0.25mm ID, 0.4μm (cat.# 10999)  
Inj.: 1.0μL split injection, components @ 500ppm.  
Oven temp.: 60°C (hold 5 min.) to 100°C @ 5°C/min. (hold 10 min.)  
Inj./det. temp.: 200°C  
Carrier gas: helium  
Linear velocity: 30cm/sec. set @ 80°C  
FID sensitivity: 6.4 x 10<sup>-11</sup> AFS  
Split flow: 46mL/min.

## Biodiesel Fuels Analysis

new!

## Rtx®-Biodiesel TG

- Linearity for all reference compounds exceeds method requirements.
- Alumaseal™ connector provides leak-free connection; guard column extends column life.
- Low column bleed at high temperatures.
- For glycerine and glyceride analysis, according to ASTM D6584 and EN 14105 methods.

## Rtx®-Biodiesel TG Columns (fused silica)

ID	df (μm)	temp. limits	10-Meter
0.32mm	0.10	to 330/380°C	10292
10-Meter w/2m x 0.53mm IP Guard Column attached using Alumaseal™ Connector			
ID	df (μm)	temp. limits	10-Meter
0.32mm	0.10	to 330/380°C	10291

## MXT®-Biodiesel TG

- Fast analysis times and sharp mono-, di-, and triglyceride peaks.
- Stable at 430°C for reliable, consistent performance.
- Integra-Gap™ built-in retention gap on a 0.53mm ID column eliminates column coupling completely.

## MXT®-Biodiesel TG Columns (Siltek® treated stainless steel)

ID	df (μm)	temp. limits	14-Meter w/2m Integra-Gap™**
0.53mm	0.16	-60 to 380/430°C	70289
ID	df (μm)	temp. limits	10-Meter w/2m x 0.53mm retention gap**
0.32mm	0.10	-60 to 380/430°C	70290


\*Total column length=16 meters.

\*\*Connected with low-dead-volume Alumaseal™ connector.

new!

Integra-Gap™ technology.

- Built-in retention gap
- Eliminates connector




**Michelle Long**  
Innovations Chemist

also available

See pages 678 and 680 for more applications on our new Rtx®-Biodiesel TG and MXT®-Biodiesel TG columns.

### Derivatized B100 and internal standards on an MXT®-Biodiesel TG column with 2m x 0.53mm ID retention gap, according to ASTM D6584.

Note perfect peak shape of solvent, with Alumaseal™ coupling.



Column: MXT®-Biodiesel TG, 10m, 0.32mm ID, 0.1μm with 2m x 0.53mm retention gap (cat.# 70290)  
 Sample: B100 + IS Butanetriol & Tricaprin derivatized with MSTFA as per ASTM D-6584  
 Inj.: 1.0μL cool on-column  
 Inj. temp.: oven track  
 Carrier gas: hydrogen, constant flow  
 Flow rate: 4mL/min.  
 Oven temp.: 50°C (hold 1 min.) to 180°C @ 15°C/min., to 230°C @ 7°C/min., to 430°C @ 30°C/min. (hold 5 min.)  
 Det.: FID @ 430°C

## Blood Alcohol Analysis

### Rtx®-BAC1/Rtx®-BAC2 (proprietary Crossbond® phase)

- Application-specific columns for blood alcohol analysis—achieve baseline resolution in less than 3 minutes. Also excellent for abused inhalant anesthetics,  $\gamma$ -hydroxybutyrate (GHB)/ $\gamma$ -butyrolactone (GBL), glycols, and common industrial solvents.
- Rtx®-BAC2 confirmation column provides four elution order changes under the same conditions.
- Stable to 260°C.

These columns separate to baseline all blood alcohol compounds in blood, breath, or urine, in less than 3 minutes, under isothermal conditions. Isothermal analysis increases productivity by eliminating the need for oven cycling. Confirmation is easily achieved with this tandem set because there are four elution order changes between the two columns.

restek  
**innovation!**

Baseline resolution in less than 3 minutes.

similar **phases**

DB-ALC1, DB-ALC2


### Rtx®-BAC1 Columns (fused silica)

ID	df ( $\mu\text{m}$ )	temp. limits	30-Meter
0.32mm	1.80	-20 to 240/260°C	18003
0.53mm	3.00	-20 to 240/260°C	18001

### Rtx®-BAC2 Columns (fused silica)

ID	df ( $\mu\text{m}$ )	temp. limits	30-Meter
0.32mm	1.20	-20 to 240/260°C	18002
0.53mm	2.00	-20 to 240/260°C	18000

### Rapid, reliable blood alcohol testing, using Rtx®-BAC column.



Inj.: 1.0mL headspace sample of a blood alcohol mix  
Sample conc.: 0.1% per compound  
Oven temp.: 40°C  
Inj./det. temp.: 200°C  
Carrier gas: helium  
Linear velocity: 80cm/sec. set @ 40°C  
FID sensitivity: 1.28 x 10<sup>-10</sup> AFS

### ordering **note**

Rtx®-BAC1 and Rtx®-BAC2 columns are available with Integra-Guard™ built-in guard columns. Get the protection without the connection! See **page 30** for descriptions and ordering information.

## Organic Volatile Impurities (OVI) Analysis



**Rick Morehead**  
Marketing Specialist  
17+ years of service!

**Rtx®-G27** (Crossbond® 5% diphenyl/95% dimethyl polysiloxane)  
with Integra-Guard™ Guard Column

**Rtx®-G43** (Crossbond® 6% cyanopropylphenyl/94% dimethyl polysiloxane)  
with Integra-Guard™ Guard Column

- Application-specific columns for residual solvents in pharmaceutical products. Meet all requirements of USP <467>.
- Analytical column with Integra-Guard™ guard column eliminates connecting problems and leaks.
- Rtx®-G27 stable to 290°C; Rtx®-G43 stable to 240°C.

Some USP <467> methods require the use of a guard column. Our Integra-Guard™ integrated guard column system makes it easy to comply.

please **note**

Analytical Reference Materials for USP <467> are available on **page 510**.

**Rtx®-G27 Column** (fused silica with 5-meter Integra-Guard™ guard column)

(Crossbond® 5% diphenyl/95% dimethyl polysiloxane)

ID	df (µm)	temp. limits	30-Meter with 5-Meter, 0.53mm ID Integra-Guard™ Guard Column
0.53mm	5.00	-60 to 270/290°C	10279-126

**Rtx®-G43 Column** (fused silica with 5-meter Integra-Guard™ guard column)

(Crossbond® 6% cyanopropylphenyl/94% dimethyl polysiloxane)

ID	df (µm)	temp. limits	30-Meter with 5-Meter, 0.53mm ID Integra-Guard™ Guard Column
0.53mm	3.00	-20 to 240°C	16085-126

## free literature

**A Technical Guide for Static Headspace Analysis Using GC**


Download your free copy from [www.restek.com](http://www.restek.com).

Technical Guide  
lit. cat.# 59895A

Searching for a chromatogram?

[www.restek.com](http://www.restek.com)

## Specially-designed Rtx®-G27 and Rtx®-G43 columns resolve residual solvents in USP &lt;467&gt;.



Inj.: 1.0µL direct injection of USP <467> Mix #1 (cat.# 36004)  
Oven temp.: 35°C (hold 5 min.) to 175°C @ 8°C/min. to 260°C @ 35°C/min.  
Inj./det. temp.: 200°C/240°C  
Carrier gas: helium  
Linear velocity: 34cm/sec. set @ 35°C  
FID sensitivity: 1 x 10<sup>-11</sup> AFS  
Recommended liner: Uniliner® inlet liner

## Organic Volatile Impurities (OVI) Analysis

**Rtx®-1301/Rtx®-624** (low to midpolarity phase; Crossbond® 6% cyanopropylphenyl/94% dimethyl polysiloxane)

- General purpose columns for residual solvents, alcohols, oxygenates, and volatile organic compounds.
- Temperature range: -20°C to 240°C.
- Equivalent to USP G43 phase.

Many analysts feel the Rtx®-1301/Rtx®-624 column has the best cyanosilicone bonded stationary phase available, with no other column manufacturer providing lower bleed, longer lifetime, or better inertness. Our polymer is fully characterized to ensure long-term reproducibility, column-to-column consistency, and low bleed—even with sensitive detectors such as ECDs and MSDs.

**Rtx®-1301 (G43) Columns** (fused silica)

(Crossbond® 6% cyanopropylphenyl/94% dimethyl polysiloxane)

ID	df (µm)	temp. limits*	15-Meter	30-Meter	60-Meter	75-Meter	105-Meter
0.53mm	3.00	-20 to 240°C	16082	16085	16088	16076	16091


**Rtx®-624 Columns** (fused silica)

(Crossbond® 6% cyanopropylphenyl/94% dimethyl polysiloxane)

ID	df (µm)	temp. limits	30-Meter	60-Meter
0.32mm	1.80	-20 to 240°C	10970	10972

\*Maximum temperatures listed are for 15- and 30-meter lengths. Longer lengths may have a slightly reduced maximum temperature.

**Rtx®-1301/Rtx®-624 Structure**



similar **phases**

DB-1301, DB-624, HP-1301, HP-624, SPB-1301, SPB-624, VF-1301, VF-624ms, CP-1301, CP-Select 624 CB

please **note**

Rtx®-1301 and Rtx®-624 columns are available with Integra-Guard™ built-in guard columns. Get the protection without the connection! See **page 30** for descriptions and ordering information.


also **available**

**MXT® Columns**


Rugged, flexible, Siltek® treated stainless steel tubing; inertness comparable to fused silica tubing. See **pages 105 and 107** for our MXT®-1301 and MXT®-624 columns.

See our OVI Retention Index on page 721.

**USP Residual Solvent Class 2 Mixture A standard solution on an Rtx®-624 column.**



## Organic Volatile Impurities (OVI) Analysis

Rtx®-5  
Structure**Rtx®-5** (low polarity phase; Crossbond® 5% diphenyl/95% dimethyl polysiloxane)

- General purpose columns for drugs, solvent impurities, pesticides, hydrocarbons, PCB congeners or (e.g.) Aroclor mixes, essential oils, semivolatiles.
- Temperature range: -60°C to 350°C.
- Equivalent to USP G27 and G36 phases.

The 5% diphenyl/95% dimethyl polysiloxane stationary phase is the most popular GC stationary phase and is used in a wide variety of applications. All residual catalysts and low molecular weight fragments are removed from the Rtx®-5 polymer, providing a tight mono-modal distribution and extremely low bleed.

## similar phases

DB-5, HP-5, HP-5MS, Ultra-2, SPB-5, Equity-5, MDN-5

**Rtx®-5 Columns** (fused silica)

(Crossbond® 5% diphenyl/95% dimethyl polysiloxane)

ID	df (µm)	temp. limits*	15-Meter	30-Meter	60-Meter
0.53mm	5.00	-60 to 270/290°C	10277	10279	10283

\*Maximum temperatures listed are for 15- and 30-meter lengths. Longer lengths may have a slightly reduced maximum temperature.

## also available


**MXT® Columns**

Rugged, flexible, Siltek® treated stainless steel tubing; inertness comparable to fused silica tubing. See page 102 for our MXT®-5 columns.

## it's a fact

For exceptional inertness, ultra-low bleed, and unsurpassed performance, choose Rxi®-5ms columns! See pages 33-35, 37.

## Organic volatile impurities on an Rtx®-5 (Rtx®-G27) column.



Rtx®-5 (Rtx®-G27) with 5m phenylmethyl Integra-Guard™ guard column, 30m, 0.53mm ID, 5.0µm (cat. # 10279-126)

Inj.: Headspace injection of common solvents for pharmaceutical processing. Prepared to equal about 500ppm in the bulk pharmaceutical. Samples shaken and heated at 90°C for 15 minutes, 1mL headspace injection.

Oven temp.: 35°C (hold 10 min.) to 100°C @ 5°C/min., to 240°C @ 25°C/min. (hold 5 min.)

Inj./det. temp.: 220°C/240°C

FID sensitivity: 1.05 x 10<sup>-11</sup> AFS

Carrier gas: helium, 35cm/sec. set @ 35°C

Split ratio: 2:1

1. ethylene oxide
2. methanol
3. ethanol
4. diethyl ether
5. 1,1-dichloroethene
6. acetone
7. isopropanol
8. acetonitrile
9. methylene chloride
10. *n*-hexane
11. *n*-propanol
12. methyl ethyl ketone
13. ethyl acetate
14. tetrahydrofuran
15. chloroform
16. 1,1,1-trichloroethane
17. carbon tetrachloride
18. benzene
19. 1,2-dichloroethane
20. heptane
21. trichloroethylene
22. *n*-butanol
23. 1,4-dioxane
24. pyridine
25. toluene

See our OVI Retention Index on page 721.

## Organic Volatile Impurities (OVI) Analysis

### Stabilwax® (polar phase; Crossbond® Carbowax® polyethylene glycol)

- General purpose columns for FAMES, flavor compounds, essential oils, amines, solvents, xylene isomers, and US EPA Method 603 (acrolein/acrylonitrile).
- Resistant to oxidative damage.
- Temperature range: 40°C to 250°C.
- Equivalent to USP G14, G15, G16, G20, and G39 phases.


Our polar-deactivated surface tightly binds the Carbowax® polymer and increases thermal stability, relative to competitive columns. The bonding mechanisms produce a column that can be rejuvenated by solvent washing. Compared to silicone stationary phases, PEG phases are more resistant to damage from strongly acidic or basic volatile compounds, including inorganic acids and volatile inorganic bases.

### Stabilwax® Columns (fused silica)

(Crossbond® Carbowax® polyethylene glycol)

ID	df (µm)	temp. limits	15-Meter	30-Meter	60-Meter
0.32mm	0.25	40 to 250°C	10621	10624	10627
0.53mm	0.25	40 to 250°C	10622	10625	10628

### Stabilwax® Structure



### similar phases

DB-WAX, DB-WAXetr,  
HP-Wax, HP-Innowax,  
Supelcowax 10

### ordering note


Stabilwax® columns are available with Integra-Guard™ built-in guard columns. Get the protection without the connection! See [page 30](#) for descriptions and ordering information.

### also available

#### MXT® Columns

Rugged, flexible, Siltek® treated stainless steel tubing; inertness comparable to fused silica tubing. See [page 106](#) for our MXT®-WAX columns.

### Residual solvents class 1 on a Stabilwax® (G16) column.




- 1,1-dichloroethene
- 1,1,1-trichloroethane
- carbon tetrachloride
- benzene
- 1,2-dichloroethane

Column:	Stabilwax®, 30m, 0.32mm ID, 0.25µm (cat. # 10624)	Headspace Conditions
Sample:	USP Stock Mixture USP<467> Residual Solvents Class 1 Mix (cat. # 36279) in 20mL headspace vial (cat. # 24685), water diluent	Instrument: Overbrook Scientific HT200H
Inj.:	headspace injection (split ratio 1:5), 2mm splitless liner IP deactivated (cat. # 20712)	Syringe temp.: 100°C
Inj. temp.:	140°C	Sample temp.: 80°C
Carrier gas:	helium, constant flow	Sample equil. time.: 45 min.
Flow rate:	2.15mL/min., 35.2cm/sec.	Injection vol.: 1.0mL
Oven temp.:	50°C for 20 min., to 165°C @ 6°C/min. (hold for 20 min.)	Injection speed: setting 8
Det.:	FID @ 250°C	Injection dwell: 5 sec.

## Semivolatiles Analysis

new!

**Rxi®-5Sil MS  
Structure**

Rxi®-5Sil MS column is recommended for US EPA Method 8270.

## similar phases

DB-5MS, VF-5ms, CP-Sil 8  
Low-Bleed/MS

**Rxi®-5Sil MS** (low polarity Crossbond® silarylene phase; similar to 5% diphenyl/95% dimethyl polysiloxane)

- Engineered to be a low bleed GC/MS column.
- Excellent inertness for active compounds.
- General purpose columns—ideal for GC/MS analysis of chlorinated hydrocarbons, phthalates, phenols, amines, organochlorine pesticides, organophosphorus pesticides, drugs, solvent impurities, hydrocarbons.
- Temperature range: -60°C to 350°C.

The Rxi®-5Sil MS stationary phase incorporates phenyl rings in the polymer backbone. This improves thermal stability, reduces bleed, and makes the phase less prone to oxidation. Rxi®-5Sil MS columns are ideal for GC/MS applications requiring high sensitivity, including use in ion trap systems.


**Rxi®-5Sil MS Columns** (fused silica)

(Crossbond®, selectivity similar to 5% diphenyl/95% dimethyl polysiloxane)

ID	df (μm)	temp. limits	15-Meter	30-Meter	60-Meter
0.25mm	0.10	-60 to 330/350°C	13605	13608	
	0.25	-60 to 330/350°C	13620	13623	13626
	0.50	-60 to 330/350°C	13635	13638	
	1.00	-60 to 325/350°C	13650	13653	13697
0.32mm	0.25	-60 to 330/350°C	13621	13624	
	0.50	-60 to 330/350°C		13639	
	1.00	-60 to 325/350°C		13654	
0.53mm	1.50	-60 to 310/330°C		13670	

ID	df (μm)	temp. limits	10-Meter	20-Meter
0.10mm	0.10	-60 to 330/350°C	43601	
0.18mm	0.18	-60 to 330/350°C		43602
	0.36	-60 to 330/350°C		43604

## Semivolatile organics by US EPA Method 8270 on an Rxi®-5Sil MS column.



**Cindy Ross**  
Southeast Sales  
Representative  
21+ years of service!

# Organophosphorus Pesticides Analysis

## Rtx®-OPPesticides/Rtx®-OPPesticides2 (proprietary Crossbond® phases)

- Application-specific columns for organophosphorus pesticides; best column combination for US EPA Method 8141A.
- Low bleed—ideal for GC/FPD, GC/NPD, or GC/MS analyses.
- Stable to 330°C.

restek  
**innovation!**

- Better separations
- Faster analysis

Using sophisticated computer modeling software, we created two stationary phases for separating the 55 organophosphorus pesticides (OPP) listed in EPA Method 8141A. Separation is improved, and analysis time is significantly reduced, compared to other columns. The extended upper temperature limit of these phases (330°C) allows analysts to bake out high molecular weight contamination typically associated with pesticide samples. The low bleed columns are a perfect match for sensitive detection systems.

## Rtx®-OPPesticides Columns (fused silica)

ID	df (μm)	temp. limits	30-Meter
0.32mm	0.50	-20 to 310/330°C	11239
0.53mm	0.83	-20 to 310/330°C	11240


## Rtx®-OPPesticides2 Columns (fused silica)

ID	df (μm)	temp. limits	20-Meter	30-Meter
0.18mm	0.20	-20 to 310/330°C	11244	
0.25mm	0.25	-20 to 310/330°C		11243
0.32mm	0.32	-20 to 310/330°C		11241
0.53mm	0.50	-20 to 310/330°C		11242

did you **know?**

These application-specific Restek columns give fast, efficient analysis of the OPPs listed in EPA Method 8141A

## Organophosphorus pesticides by US EPA Method 8141A on an Rtx®-OPPesticides2 column.



- dichlorvos
- hexamethylphosphoramide
- mevinphos
- trichlorfon
- TEPP
- demeton-O
- thionazin
- tributyl phosphate (IS)
- ethoprop
- naled
- sulfotepp
- phorate
- dicrotophos
- demeton-S
- monocrotophos
- terbufos
- dimethoate
- dioxathion
- fonophos
- diazinon
- disulfoton
- phosphamidon isomer
- dichlorofenthion
- chlorpyrifos methyl
- phosphamidon
- parathion-methyl
- ronnel
- fenitrothion
- aspon
- malathion
- chlorpyrifos
- trichloronate
- parathion-ethyl
- fenthion
- merphos
- chlorfenvinphos
- crotoxyphos
- stirofos
- tokuthion
- merphos oxone (breakdown product)
- ethion
- fensulfothion
- bolstar
- carbofenthion
- famphur
- triphenyl phosphate (SS)
- EPN
- phosmet
- leptophos
- tri-o-cresyl phosphate
- azinphos-methyl
- azinphos-ethyl
- coumaphos

Column: Rtx®-OPPesticides2, 30m, 0.25mm ID, 0.25μm (cat.# 11243)  
Sample: US EPA Method 8141A Custom Standard Mix 1μL 100ppm (100ng on column)

Triphenylphosphate Standard (cat.# 32281)  
Tributylphosphate Standard (cat.# 32280)  
8140/8141 OP Pesticides Calibration Mix A (cat.# 32277)  
8141 OP Pesticides Calibration Mix B (cat.# 32278)

Inj.: Custom Mixes: Call Restek for Information  
1.0μL splitless (hold 0.4 min.), 4mm double  
goose-neck inlet liner (cat.# 20785)

Inj. temp.: 250°C  
Carrier gas: helium, constant flow  
Flow rate: 1.0mL/min.  
Oven temp.: 80°C (hold 0.5 min.) to 140°C @ 20°C/min.  
to 210°C @ 4°C/min. (hold 1 min.) to  
280°C @ 30°C (hold 5 min.)

Det: MS  
Transfer line temp.: 280°C  
Scan range: 35-400amu  
Ionization: EI

## Chlorinated Pesticides Analysis

restek  
innovation!

- Very low bleed
- Faster analysis

**Rtx®-CLPesticides/Rtx®-CLPesticides2** (proprietary Crossbond® phases)

- Application-specific columns for organochlorine pesticides and herbicides.
- Low bleed—ideal for GC/ECD or GC/MS analyses.
- Baseline separations in less than 10 minutes.
- Stable to 340°C.

Improved resolution and faster analyses, compared to 1701 or phenyl phases, make these the pesticide columns of choice. Rtx®-CLPesticides columns are specially designed to overcome the coelutions and analyte breakdown typically encountered in chlorinated pesticide analyses for EPA Methods 8081, 608, and CLP. By achieving baseline resolution of the 20 target analytes, more accurate qualitative data can be obtained, providing reliable identification without GC/MS.

Column bleed, measured by ECD, is extremely low at temperatures up to 330°C, which is critical for baking-out the column to remove high-boiling compounds commonly found in pesticide/PCB extracts. An analysis time of less than 10 minutes improves throughput compared to other stationary phases.

## free literature

**Fast GC Analysis of Chlorinated Pesticides**

Download your free copy from  
www.restek.com

Flyer  
lit. cat.# 59547A

**Rtx®-CLPesticides Columns** (fused silica)

ID	df (µm)	temp. limits	10-Meter	15-Meter	20-Meter	30-Meter	60-Meter
0.10mm	0.10	-60 to 310/330°C	43101				
0.18mm	0.18	-60 to 310/330°C	42101		42102		
0.25mm	0.25	-60 to 320/340°C		11120		11123	11126
0.32mm	0.32	-60 to 320/340°C			<b>new!</b>	11141	
	0.50	-60 to 320/340°C		11136		11139	
0.53mm	0.50	-60 to 300/320°C		11137		11140	

**Rtx®-CLPesticides2 Columns** (fused silica)

ID	df (µm)	temp. limits	10-Meter	15-Meter	20-Meter	30-Meter	60-Meter
0.10mm	0.10	-60 to 310/330°C	43301		43302		
0.18mm	0.14	-60 to 310/330°C	42301		42302		
0.25mm	0.20	-60 to 320/340°C		11320		11323	11326
0.32mm	0.25	-60 to 320/340°C		11321		11324	
	0.50	-60 to 320/340°C				11325	
0.53mm	0.42	-60 to 300/320°C		11337		11340	

**Rtx®-CLPesticides Column Kits**

(Note: Columns are not preconnected in these kits.)

**0.25mm ID Rtx®-CLPesticides Kit** cat.# 11199 (kit),**Includes:**

	cat.#
30m, 0.25mm ID, 0.25µm Rtx®-CLPesticides Column	11123
30m, 0.25mm ID, 0.25µm Rtx®-CLPesticides2 Column	11323
Universal Angled "Y" Press-Tight® Connector	20403
5m, 0.25mm ID Siltek® Guard Column	10026

kit

new!

→ New column dimensions now available.

**0.32mm ID Rtx®-CLPesticides Kit** cat.# 11196 (kit),**Includes:**

	cat.#
30m, 0.32mm ID, 0.32µm Rtx®-CLPesticides Column	11141
30m, 0.32mm ID, 0.25µm Rtx®-CLPesticides2 Column	11324
Universal Angled "Y" Press-Tight® Connector	20403
5m, 0.32mm ID Siltek® Guard Column	10027

kit

**0.53mm ID Rtx®-CLPesticides Kit** cat.# 11197 (kit),**Includes:**

	cat.#
30m, 0.53mm ID, 0.50µm Rtx®-CLPesticides Column	11140
30m, 0.53mm ID, 0.42µm Rtx®-CLPesticides2 Column	11340
Universal Angled "Y" Press-Tight® Connector	20403
5m, 0.53mm ID IP Deactivated Guard Column	10045

kit

also available

For column connectors, see pages 220-227.

**Add a reference mix to your kit order and save!**

Description	suffix #
Organochlorine Pesticide Mix AB #1 (cat.# 32291)	-530
Organochlorine Pesticide Mix AB #2 (cat.# 32292)	-535



## Chlorinated Pesticides Analysis

## Fast GC analyses of chlorinated pesticides on Rtx®-CLPesticides and Rtx®-CLPesticides2 columns.

## Rtx®-CLPesticides &amp; Rtx®-CLPesticides2 columns (0.32mm ID)

## Rtx®-CLPesticides

## Rtx®-CLPesticides2



Columns: Rtx®-CLPesticides, 30m, 0.32mm ID, 0.32 $\mu$ m (cat.# 11141) and Rtx®-CLPesticides2, 30m, 0.32mm ID, 0.25 $\mu$ m (cat.# 11324) with 5m x 0.32mm ID Rxi® deactivated guard tubing (cat.# 10039), connected using Deactivated Universal "Y" Press-Tight® connector (cat.# 20405-261)

Sample: Organochlorine Pesticide Mix AB #2, 8-80 $\mu$ g/mL each component in hexane/toluene (cat.# 32292), Pesticide Surrogate Mix, 200 $\mu$ g/mL each component in acetone (cat.# 32000)

Inj.: 1.0 $\mu$ L splitless (hold 0.3 min.), 4mm single gooseneck inlet liner (cat.# 20799)

Inj. temp.: 250°C

Carrier gas: helium, constant flow

Linear velocity: 60cm/sec. @ 120°C

Oven temp.: 120°C to 200°C @ 45°C/min. to 230°C @ 15°C/min. to 330°C (hold 2 min.) @ 30°C/min.



Det.: ECD @ 330°C

1. 2,4,5,6-tetrachloro-*m*-xylene (SS)
2.  $\alpha$ -BHC
3.  $\gamma$ -BHC
4.  $\beta$ -BHC
5.  $\delta$ -BHC
6. heptachlor
7. aldrin
8. heptachlor epoxide (isomer B)
9.  $\gamma$ -chlordane
10.  $\alpha$ -chlordane
11. endosulfan I
12. 4,4'-DDE
13. dieldrin
14. endrin
15. 4,4'-DDD
16. endosulfan II
17. 4,4'-DDT
18. endrin aldehyde
19. endosulfan sulfate
20. methoxychlor
21. endrin ketone
22. decachlorobiphenyl (SS)

## Rtx®-CLPesticides &amp; Rtx®-CLPesticides2 columns (0.53mm ID)

## Rtx®-CLPesticides

## Rtx®-CLPesticides2



Columns: Rtx®-CLPesticides, 30m, 0.53mm ID, 0.50 $\mu$ m (cat.# 11140) and Rtx®-CLPesticides2, 30m, 0.53mm ID, 0.42 $\mu$ m (cat.# 11340) with 5m x 0.53mm ID Rxi® deactivated guard tubing (cat.# 10054), connected using Universal Siltek® "Y" Press-Tight® connector (cat.# 20486)

Sample: Organochlorine Pesticide Mix AB #2, 8-80 $\mu$ g/mL each component in hexane/toluene (cat.# 32292), Pesticide Surrogate Mix, 200 $\mu$ g/mL each component in acetone (cat.# 32000)

Inj.: 1.0 $\mu$ L splitless (hold 0.3 min.), 4mm single gooseneck inlet liner (cat.# 20799)

Inj. temp.: 250°C

Carrier gas: helium, constant flow

Linear velocity: 45cm/sec. @ 120°C

Oven temp.: 120°C to 200°C @ 45°C/min. to 230°C @ 12.5°C/min. to 325°C (hold 2 min.) @ 30°C/min.

Det.: ECD @ 330°C

restek  
innovation!

## Chlorinated Pesticides Analysis

### Stx™-CLPesticides/Stx™-CLPesticides2 (proprietary Crossbond® phases)

- Application-specific columns for organochlorine pesticides and herbicides.
- Baseline separations in less than 10 minutes.
- Siltek® surface deactivation enhances responses for endrin, DDT, methoxychlor.
- Stable to 330°C.

Many laboratories analyzing organochlorine pesticides struggle with breakdown and adsorption of endrin, DDT, and methoxychlor caused by active sites throughout the analytical system. Siltek® passivation technology enables these columns to offer unsurpassed inertness and the highest responses for active pesticides.

### Stx™-CLPesticides Columns (fused silica with Siltek® deactivation)

ID	df (µm)	temp. limits	15-Meter	30-Meter
0.25mm	0.25	-60 to 310/330°C	11540	11543
0.32mm	0.32	-60 to 310/330°C		<b>new!</b> 11546
	0.50	-60 to 310/330°C	11541	11544

### Stx™-CLPesticides2 Columns (fused silica with Siltek® deactivation)

ID	df (µm)	temp. limits	15-Meter	30-Meter
0.25mm	0.20	-60 to 310/330°C	11440	11443
0.32mm	0.25	-60 to 310/330°C	11441	11444

### Stx™-CLPesticides Kits

(Note: Columns are not preconnected in these kits.)

### it's a fact

These columns are treated with Siltek® deactivation, which provides better responses for endrin, DDT, and methoxychlor.

### ordering note

Kits include Siltek® deactivated guard column.

0.25mm ID Stx™-CLPesticides Kit	cat.# 11190 (kit)	kit
<b>Includes:</b>		<b>cat.#</b>
30m, 0.25mm ID, 0.25µm Stx™-CLPesticides Column		11543
30m, 0.25mm ID, 0.20µm Stx™-CLPesticides2 Column		11443
Universal Angled "Y" Press-Tight® Connector		20403
5m, 0.25mm ID Siltek® Guard Column		10026

### new!

New column dimensions now available.

0.32mm ID Stx™-CLPesticides Kit	cat.# 11193 (kit)	kit
<b>Includes:</b>		<b>cat.#</b>
30m, 0.32mm ID, 0.32µm Stx™-CLPesticides Column		11546
30m, 0.32mm ID, 0.25µm Stx™-CLPesticides2 Column		11444
Universal Angled "Y" Press-Tight® Connector		20403
5m, 0.32mm ID Siltek® Guard Column		10027



**Get More!**


Environmental  
Related Articles Online

[www.restek.com/environmental](http://www.restek.com/environmental)

## Chlorinated Pesticides Analysis


Excellent responses for all US EPA Method 8081  
chlorinated pesticides: Stx™-CLPesticides columns.

**Stx™-CLPesticides**  
30m, 0.32mm ID, 0.5µm  
(cat.# 11544)



1. 2,4,5,6 tetrachloro-*m*-xylene (SS)
2.  $\alpha$ -BHC
3.  $\gamma$ -BHC
4.  $\beta$ -BHC
5.  $\delta$ -BHC
6. heptachlor
7. aldrin
8. heptachlor epoxide
9.  $\gamma$ -chlordane
10.  $\alpha$ -chlordane
11. 4,4'-DDE
12. endosulfan I
13. dieldrin
14. endrin
15. 4,4'-DDD
16. endosulfan II
17. 4,4'-DDT
18. endrin aldehyde
19. methoxychlor
20. endosulfan sulfate
21. endrin ketone
22. decachlorobiphenyl (SS)

**Stx™-CLPesticides2**  
30m, 0.32mm ID, 0.25µm  
(cat.# 11444)



## free literature

For detailed information about  
these columns, request  
Applications Note  
lit. cat.# 59351B

GC\_EV00512

Inj.: 1µL direct injection of 20/40/200ng/mL std. concentration in hexane  
Oven temp.: 110°C (hold 1 min.) to 245°C @ 20°C/min. to 300°C @ 6°C/min.  
Inj./det. temp.: 210°C/310°C  
Carrier gas: helium  
Dead time: 0.8min. @ 120°C  
Inlet liner: Silitex® deactivated Drilled Uniliner® inlet liner (cat.# 21055-214.1)  
Make-up gas: nitrogen

Table of Contents for  
Applications

see pages 518-519



**Jason Thomas**  
Innovations Chemist  
6+ years of service!

## it's a fact

The Stx™-CLPesticides and Stx™-CLPesticides2 column pair will provide <10 minute analysis times if you use the same conditions used for the Rtx®-CLPesticides and Rtx®-CLPesticides2 column pair listed on [page 79!](#)

restek  
innovation!

## PCB Congeners Analysis

### Rtx®-PCB (proprietary Crossbond® phase)


- Unique polymer for PCBs analysis by GC/ECD or GC/MS.
- Good results for other semivolatiles.
- Low polarity; inert to active compounds.
- Stable to 340°C.

### Rtx®-PCB Columns (fused silica)

ID	df (µm)	temp. limits*	20-Meter	30-Meter	40-Meter	60-Meter
0.18mm	0.18	30°C to 320/340°C	41302		41303	41304
0.25mm	0.25	30°C to 320/340°C		13223		13226
0.32mm	0.50	30°C to 320/340°C		13239		

\*Maximum temperatures listed are for 15- and 30-meter lengths. Longer lengths may have a slightly reduced maximum temperature.

### Aroclor 1242/1254/1262 PCBs on Rtx®-PCB: best available resolution of individual congeners.



Column: Rtx®-PCB, 60m, 0.25mm ID, 0.25µm (cat.# 13226)  
 Sample: Aroclor 1242 (cat.# 32009), 1254 (cat.# 32011), 1262 (cat.# 32409), 333ppm each  
 Inj.: 1.0µL splitless (hold 0.75 min.), 4mm single gooseneck inlet liner w/wool (cat.# 22405)  
 Inj. temp.: 280°C  
 Carrier gas: helium, constant flow  
 Flow rate: 1.1mL/min.  
 Oven temp.: 100°C (hold 1 min.) to 200°C @ 30°C/min., to 320°C @ 2°C/min. (hold 1 min.)  
 Det.: MS  
 Transfer line temp.: 280°C  
 Scan range: 50 to 550amu  
 Ionization: EI  
 Mode: scan

## PCB Congeners Analysis

### Stx™-500 (Crossbond® carborane/dimethyl polysiloxane)

- Application-specific columns for brominated flame retardants, coplanar PCB congeners, and other analytes with high boiling temperatures.
- Low bleed—ideal for GC/FPD, GC/NPD, or GC/MS analyses.
- Stable to 380°C.
- Stx™ is used for columns that have been deactivated using Restek's Siltek® deactivation.

The Stx™-500 column gives excellent results for neutral or slightly acidic compounds. It is not recommended for analyses of basic compounds.

### Stx™-500 Columns (fused silica)

(Crossbond® carborane/dimethyl polysiloxane)

ID	df (µm)	temp. limits*	30-Meter	60-Meter
0.25mm	0.15	-60°C to 380°C	10750	10751
0.53mm	0.15	-60°C to 380°C	10752	

\*Maximum temperatures listed are for 15- and 30-meter lengths. Longer lengths may have a slightly reduced maximum temperature.

similar **phase**

HT-8

## Dioxin & Furan Congeners Analysis

### Rtx®-Dioxin (proprietary Crossbond® phase)

- Replacement column for 5% diphenyl phases.
- Improved separations of dioxin or furan congeners.
- Greater thermal stability than 5% diphenyl phases or high-cyano confirmation columns.

### Rtx®-Dioxin Columns (fused silica)


ID	df (µm)	temp. limits	60-Meter
0.25mm	0.15	-60°C to 380°C	10755

restek  
**innovation!**

### Rtx®-Dioxin column separates all five components in the TCDD resolution check mixture.

also **available**

Rtx®-Dioxin2 columns.  
See **page 84**.



Column: Rtx®-Dioxin, 40m, 0.18mm ID, 0.11µm  
Initial temp.: 130°C  
Instrument: Micromass Altima high resolution GC/MS

Chromatography courtesy of Karen MacPherson and Eric Reiner, Ontario Ministry of the Environment, Etobicoke, ON, Canada.

restek  
innovation!

Excellent for dioxins or furans.

## Dioxin & Furan Congeners Analysis


### Rtx®-Dioxin2 (proprietary Crossbond® phase)

- Isomer specificity for 2,3,7,8-TCDD and 2,3,7,8-TCDF achieved with one GC column.
- Thermally stable to 340°C for longer lifetime.
- Unique selectivity for toxic dioxin and furan congeners allow use as a primary or confirmation GC column.

### Rtx®-Dioxin2 Columns (fused silica)

ID	df (μm)	temp. limits	40-Meter	60-Meter
0.18mm	0.18	20°C to 340°C	10759	—
0.25mm	0.25	20°C to 340°C	—	10758

### 2,3,7,8-Tetrachlorodibenzodioxin resolved from other TCDD congeners, using an Rtx®-Dioxin2 column.



Other peak identifications available upon request.

GC\_EV00948


Column: Rtx®-Dioxin2, 60m, 0.25mm ID, 0.25μm (cat.# 10758)  
 Sample: WMS-01 Reference Material, Wellington Laboratories  
 Inj.: Splitless  
 Inj. temp.: 250°C  
 Carrier gas: helium, constant flow  
 Flow rate: 1.5mL/min.  
 Oven temp.: 130°C (hold 1.0 min.) to 200°C @ 40°C/min. to 235°C @ 3.0°C/min. to 300°C @ 5°C/min. (hold 10 min.)  
 Det.: Micromass Ultima high-resolution mass spectrometer  
 Ionization: EI  
 Mode: SIR



**Jack Cochran**

Director of New Business & Technology

### Tetrachlorodibenzofuran congeners on an Rtx®-Dioxin2 column.



Other peak identifications available upon request.

GC\_EV00949

Column: Rtx®-Dioxin2, 60m, 0.25mm ID, 0.25μm (cat.# 10758)  
 Sample: WMS-01 Reference Material, Wellington Laboratories  
 Inj.: Splitless  
 Inj. temp.: 250°C  
 Carrier gas: helium, constant flow  
 Flow rate: 1.5mL/min.  
 Oven temp.: 130°C (hold 1.0 min.) to 200°C @ 40°C/min. to 235°C @ 3.0°C/min. to 300°C @ 5°C/min. (hold 10 min.)  
 Det.: Micromass Ultima high-resolution mass spectrometer  
 Ionization: EI  
 Mode: SIR

### free literature

#### Rtx®-Dioxin2 Column: 2,3,7,8-TCDD and 2,3,7,8-TCDF Specificity in One GC Column

New data available—all 128 dioxin and furan tetra thru octa congeners acquired on the Rtx®-Dioxin2 column.

Download your free copy from [www.restek.com](http://www.restek.com)

Flyer

lit. cat.# 580119A

Chromatograms courtesy of Terry Kolic, Karen MacPherson, Eric Reiner, Ontario Ministry of the Environment, Toronto, Ontario, Canada

## Polycyclic Aromatic Hydrocarbon (PAH) Analysis


### Rt™-PAH (polar, proprietary liquid crystalline phase)

- Specially designed for the analysis of polycyclic aromatic hydrocarbons (PAHs) listed in US EPA methods 610 and 8100.
- Low bleed at 285°C.
- Temperature range: 80°C to 285°C.

### Rt™-PAH Columns (fused silica)

ID	df (μm)	temp. limits	12-Meter
0.25mm	0.15	80°C to 285°C	19733

### Separation of 16 regulated polycyclic aromatic hydrocarbons (PAHs) in less than 15 minutes.



Column: Rt™-PAH, 12m, 0.25mm ID, 0.15μm (cat.# 19733)  
 Sample: 16 component EPA Method 610 PAH standard  
 (20ng/μl of each component in dichloromethane)  
 Inj.: 1.0μL split (split ratio 10:1)  
 Inj. temp.: 225°C  
 Carrier gas: helium, 110kPa column head pressure  
 Oven temp.: 80°C to 220°C @ 40°C/min., 220°C to 285°C @ 8°C/min. (hold 5 min.)  
 Detector: FID @ 290°C

Chromatogram courtesy of J&K Scientific.

### Rt™-LC50 (polar, dimethyl (50% liquid crystal) polysiloxane)

- General purpose column with selectivity for dioxin or furan congeners, or PCB congeners.
- Low bleed at 270°C.
- Temperature range: 100°C to 270°C.

The unique liquid crystalline Rt™-LC50 stationary phase resolves compounds of similar structure and boiling point. It has proven effective for resolving many polycyclic aromatic hydrocarbons; other potential applications include dioxin, furan, or PCB congeners.

### Rt™-LC50 Columns (fused silica)

ID	df (μm)	temp. limits	10-Meter	20-Meter
0.10mm	0.10	100°C to 270°C	19736	—
0.18mm	0.10	100°C to 270°C	19735	—
0.25mm	0.10	100°C to 270°C	—	19734

## Explosives Analysis

restek  
innovation!

**Rtx®-TNT/Rtx®-TNT2** (proprietary Crossbond® phase)

- Application-specific columns for explosives in US EPA Method 8095.
- Low bleed—ideal for ECD analysis.
- Complete analysis in less than 20 minutes.
- Rtx®-TNT2 confirmation column provides 8 elution order changes under same conditions.
- Economical 3-packs.
- Stable to 310°C.

## please note

Polymer specially designed for explosives analysis.

We designed Rtx®-TNT and Rtx®-TNT2 columns specifically for analyses of nitroaromatic compounds by GC/ECD, such as the 16 analytes listed in US EPA Method 8095. They provide better resolution and higher thermal stability than any other currently recommended columns. Operate the Rtx®-TNT primary column and Rtx®-TNT2 confirmation column under identical GC oven temperature programs.

**Rtx®-TNT Columns** (fused silica)

ID	df (µm)	temp. limits	6-Meter/3-pk.
0.53mm	1.50	-20 to 300/310°C	12998


**Rtx®-TNT2 Columns** (fused silica)

ID	df (µm)	temp. limits	6-Meter/3-pk.
0.53mm	1.50	-20 to 300/310°C	12999

## US EPA Method 8095 explosives on Rtx®-TNT and Rtx®-TNT2 columns.

**Rtx®-TNT**


6m, 0.53mm ID, 1.50µm (cat.# 12998)



1. nitrobenzene
2. 2-nitrotoluene
3. 3-nitrotoluene
4. 4-nitrotoluene
5. nitroglycerine
6. 1,3-dinitrobenzene
7. 2,6-dinitrotoluene
8. 2,4-dinitrotoluene
9. 3,4-dinitrotoluene (IS)
10. 1,3,5-trinitrobenzene
11. trinitrotoluene
12. PETN
13. RDX
14. 4-amino-2,6-dinitrotoluene
15. 3,5-dinitroaniline
16. 2-amino-4,6-dinitrotoluene
17. tetryl

**Rtx®-TNT2**

6m, 0.53mm ID, 1.50µm (cat.# 12999)



Inj.: Direct injection using a 1mm Siltek® Uniliner® (cat.# 21052-214.1)  
 On-column conc.: est. 200-1000pg/compound. 8095 Calibration Mix A (cat.# 31607),  
 8095 Calibration Mix B (cat.# 31608), and 3,4-dinitrotoluene (cat.# 31452)  
 Oven temp.: 80°C (hold 1 min.) to 180°C @ 10°C/min. to 300°C @ 30°C/min. (hold 3 min.)  
 Inj. temp.: 250°C  
 Det.: ECD @ 330°C with anode purge  
 Dead time: 4.4 sec.  
 Head pressure: helium @ 3psi (20.7 KPa)  
 Flow rate: helium @ 17mL/min. @ 80°C



**Jarl Snider**  
R&D Chemist  
12+ years of service!

## Volatile Organics Analysis

### Rtx®-VMS (proprietary Crossbond® phase)

- Application-specific columns for volatile organic pollutants by GC/MS.
- Complete separation of US EPA Method 8260 compounds in less than 18 minutes.
- Stable to 260°C.
- No known equivalent phases.

Rtx®-VMS columns offer lower bleed, better selectivity, and overall faster analysis for separating volatile organic compounds, such as those listed in US EPA Method 8260. The Rtx®-VMS stationary phase is a highly stable polymer that provides outstanding analysis of volatile compounds, in combination with sensitive ion traps and Agilent 5973 mass spectrometers. 0.18 and 0.25mm ID columns allow sample splitting at the injection port, eliminating the added expense and maintenance of a jet separator. A 0.45mm or 0.53mm ID column can be directly connected to the purge & trap transfer line in a system equipped with a jet separator.

### Rtx®-VMS Columns (fused silica)

ID	df (µm)	temp. limits	30-Meter	60-Meter	75-Meter
0.25mm	1.40	-40 to 240/260°C	19915	19916	
0.32mm	1.80	-40 to 240/260°C	19919	19920	
0.45mm	2.55	-40 to 240/260°C	19908	19909	
0.53mm	3.00	-40 to 240/260°C	19985	19988	19974

ID	df (µm)	temp. limits	20-Meter	40-Meter
0.18mm	1.00	-40 to 240/260°C	49914	49915

### restek innovation!


- First choice for use with dual purge & traps<sup>1</sup>
- EPA recommended surrogate used.

<sup>1</sup>A.L. Hilling and G. Smith, Environmental Testing & Analysis, 10(3), 15-19, 2001.

Need a column for a  
volatiles analysis?

see page 563

### Rapid analysis of volatile organics in US EPA Method 8260B, on an Rtx®-VMS column.



Column: Rtx®-VMS, 20m, 0.18mm ID, 1.00µm (cat.# 49914)  
 Conc.: 10ppb in 5mL of RO water  
 unless otherwise noted; ketones at 2.5X  
 Concentrator: Tekmar LSC-3100 Purge and Trap  
 Trap: Vocabr 3000 (type K)  
 Purge: 11 min. @ 40mL/min. (ambient temperature)  
 Dry purge: 1 min. @ 40mL/min.  
 Desorb preheat: 245°C  
 Desorb: 250°C for 2 min., flow 40mL/min.  
 Bake: 260°C for 8 min.  
 Interface: 0.53mm ID Silcosteel® tubing transfer line  
 1:40 split at injection port. 1mm ID liner.  
 Oven temp.: 50°C (hold 4 min.) to 100°C @ 18°C/min. (hold 0 min.)  
 to 230°C @ 40°C/min. (hold 3 min.)  
 Carrier gas: helium @ ~1.0mL/min. constant flow  
 Adjust dichlorodifluoromethane to a retention time of 1.03 min. @ 50°C.  
 Detector: Agilent 5973 MSD  
 Scan range: 35-300amu

## Volatile Organics Analysis

restek  
**innovation!**

Polymer specially designed for  
volatiles analysis by PID/ELCD.

**Rtx®-VGC** (proprietary Crossbond® phase)

- Application-specific columns for volatile organic pollutants by GC/PID/ELCD. Excellent separation of trihalomethanes.
- Complete US EPA Method 8021A analysis in less than 28 minutes.
- Stable to 260°C.
- No known equivalent phases.

Using computer modeling techniques, we optimized the Rtx®-VGC column for analysis of volatile organic compounds on GC systems equipped with photoionization (PID) and electrolytic conductivity detectors (ELCD). It performs the most difficult separations of volatile organic compounds, such as those listed in US EPA Methods 502.2 and 8021, providing unsurpassed separation in the fastest analysis time, thereby increasing sample throughput. The Rtx®-VGC column provides  $\geq 85\%$  resolution of trihalomethanes (THMs) from other volatile compounds. This unique column also achieves excellent separation of gases and early eluting compounds.

**Rtx®-VGC Columns** (fused silica)

ID	df ( $\mu\text{m}$ )	temp. limits	30-Meter	60-Meter	75-Meter	105-Meter
0.25mm	1.40	-40 to 240/260°C	19415	19416		
0.32mm	1.80	-40 to 240/260°C	19419	19420		
0.45mm	2.55	-40 to 240/260°C	19408		19409	
0.53mm	3.00	-40 to 240/260°C	19485	19488	19474	19489

ID	df ( $\mu\text{m}$ )	temp. limits	20-Meter	40-Meter
0.18mm	1.00	-40 to 240/260°C	49414	49415

Need a column for a  
volatiles analysis?


see page 563

**Fast separation of US EPA Method 8021 volatile organics: Rtx®-VGC column.**




Primary column, dual-column analysis. Confirmation analysis shown on page 566.

**Rtx®-VGC**

75m, 0.45mm ID, 2.55 $\mu\text{m}$  (cat.# 19409)



See page 565 for the peak list  
and run conditions.

**PID****ELCD**

GC\_EV00416

Acknowledgement: Finnigan 9001 GC,  $\mu\text{Gold}$  Tandem Photoionization/HALL® 2000 Electrolytic Conductivity Detector provided courtesy of Thermo Scientific GC & GC/MS Division, 2215 Grand Avenue Pkwy, Austin, Texas 78728

# Volatile Organics Analysis

## Rtx®-VRX (proprietary Crossbond® phase)

- Application-specific columns for volatile organic pollutants.
- Excellent for US EPA Method 8021 compounds.
- Stable to 260°C.

The Rtx®-VRX stationary phase and optimized column dimensions provide low bleed, excellent resolution, and fast analysis times for volatile compounds.

similar **phase**

DB-VRX

## Rtx®-VRX Columns (fused silica)

ID	df (µm)	temp. limits	30-Meter	60-Meter	75-Meter	105-Meter
0.25mm	1.40	-40 to 240/260°C	19315	19316		
0.32mm	1.80	-40 to 240/260°C	19319	19320		
0.45mm	2.55	-40 to 240/260°C	19308		19309	
0.53mm	3.00	-40 to 240/260°C	19385	19388	19374	19389


ID	df (µm)	temp. limits	20-Meter	40-Meter
0.18mm	1.00	-40 to 240/260°C	49314	49315

Need a column for a  
volatiles analysis?

see page 563

GC COLUMNS | SPECIALTY COLUMNS BY APPLICATION

### Excellent resolution of EPA Method 8021 volatile organics on an Rtx®-VRX column.



## Volatile Organics Analysis

**Rtx<sup>®</sup>-502.2** (proprietary Crossbond<sup>®</sup> diphenyl/dimethyl polysiloxane phase)

- Application-specific columns with unique selectivity for volatile organic pollutants. The Rtx<sup>®</sup>-502.2 column is cited in US EPA Method 502.2 and in many gasoline range organics (GRO) methods for monitoring underground storage tanks.
- Excellent separation of trihalomethanes; ideal polarity for light hydrocarbons and aromatics.
- Stable to 270°C.

similar **phase**

DB-502.2

An Rtx<sup>®</sup>-502.2 column will enable you to quantify all compounds listed in US EPA methods 502.2 or 524.2, whether you use a mass spectrometer or a PID in tandem with an ELCD. The diphenyl/dimethyl polysiloxane based Rtx<sup>®</sup>-502.2 stationary phase provides low bleed and thermal stability to 270°C. A 105-meter column can separate the light gases specified in EPA methods without subambient cooling. Narrow bore columns can interface directly in GC/MS systems.

**Rtx<sup>®</sup>-502.2 Columns** (fused silica)(proprietary Crossbond<sup>®</sup> diphenyl/dimethyl polysiloxane phase)

ID	df (μm)	temp. limits	30-Meter	60-Meter	75-Meter	105-Meter
0.25mm	1.40	-20 to 250/270°C	10915	10916		
0.32mm	1.80	-20 to 250/270°C	10919	10920		10921
0.45mm	2.55	-20 to 250/270°C			10986	
0.53mm	3.00	-20 to 250/270°C	10908	10909		10910

ID	df (μm)	temp. limits	20-Meter	40-Meter
0.18mm	1.00	-20 to 250/270°C	40914	40915

also **available****MXT<sup>®</sup> Columns**

Rugged, flexible, Siltek<sup>®</sup> treated stainless steel tubing; inertness comparable to fused silica tubing. See **pages 106 and 107** for our MXT<sup>®</sup>-502.2 and MXT<sup>®</sup> Volatiles columns.

similar **phase**VOCOL<sup>®</sup>**Rtx<sup>®</sup>-Volatiles** (proprietary Crossbond<sup>®</sup> diphenyl/dimethyl polysiloxane phase)

- Application-specific columns for volatile organic pollutants.
- Stable to 280°C.

Rtx<sup>®</sup>-Volatiles columns were the first columns designed specifically for analyses of the 34 volatile organic pollutants listed in US EPA methods 601, 602, and 624. With these columns, you can quantify all compounds listed in these methods, whether you use a mass spectrometer or a PID in tandem with an ELCD. The diphenyl/dimethyl polysiloxane based Rtx<sup>®</sup>-Volatiles stationary phase provides low bleed and thermal stability to 280°C. Narrow bore columns can interface directly in GC/MS systems.

**Rtx<sup>®</sup>-Volatiles Columns** (fused silica)(proprietary Crossbond<sup>®</sup> diphenyl/dimethyl polysiloxane phase)

ID	df (μm)	temp. limits*	30-Meter	60-Meter	105-Meter
0.25mm	1.00	-20 to 270/280°C	10900	10903	
0.32mm	1.50	-20 to 270/280°C	10901	10904	
0.53mm	2.00	-20 to 270/280°C	10902	10905	10906

ordering **note**

Rtx<sup>®</sup>-Volatiles columns are available with Integra-Guard™ built-in guard columns. Get the protection without the connection! See **page 30** for descriptions and ordering information.

it's a **fact**

Quantify all compounds in US EPA method 601, 602, or 624, using an Rtx<sup>®</sup>-Volatiles column.

## Volatile Organics Analysis

**Rtx®-624** (low to midpolarity phase; Crossbond® 6% cyanopropylphenyl/94% dimethyl polysiloxane)

- Application-specific columns for volatile organic pollutants. Recommended in US EPA methods for volatile organic pollutants.
- Temperature range: -20°C to 240°C.
- Equivalent to USP G43 phase.

The unique polarity of the Rtx®-624 column makes it ideal for analyzing volatile organic pollutants. Although the Rtx®-502.2 column is recommended in many methods, the Rtx®-624 column offers better resolution of early eluting compounds. The Rtx®-624 phase produces greater than 90% resolution of the first six gases in EPA Methods 8260 and 524.2. This stationary phase is especially well-suited for EPA Method 524.2 revision IV since it resolves 2-nitropropane from 1,1-dichloropropanone, which share quantification ion *m/z* 43 and must be separated chromatographically.


### Rtx®-624 Columns (fused silica)

(Crossbond® 6% cyanopropylphenyl/94% dimethyl polysiloxane)

ID	df (µm)	temp. limits	30-Meter	60-Meter	75-Meter	105-Meter
0.25mm	1.40	-20 to 240°C	10968	10969		
0.32mm	1.80	-20 to 240°C	10970	10972		
0.45mm	2.55	-20 to 240°C			10982	
0.53mm	3.00	-20 to 240°C	10971	10973	10974	10975

ID	df (µm)	temp. limits	10-Meter	20-Meter	40-Meter
0.18mm	1.00	-20 to 240°C		40924	40925

### Rtx®-624 Structure



### similar phases

DB-1301, DB-624, HP-1301, HP-624, SPB-1301, SPB-624, VF-1301, VF-624ms, CP-1301, CP-Select 624 CB

### also available

#### MXT® Columns

Rugged, flexible, Siltek® treated stainless steel tubing; inertness comparable to fused silica tubing. See [page 107](#) for our MXT®-624 columns.

## custom standards

Restek is your #1 source for custom analytical reference materials!

- Made to your exact specifications.
- Quick quotations.
- Most orders shipped within 5–10 working days.

Call our reference materials department, use the form on [page 427](#) or at [www.restek.com/solutions](http://www.restek.com/solutions), or contact your Restek representative for more information.



**Kristi Sellers**  
Innovations Chemist  
20+ years of service!