



GC Innovations

MGA #1 uTCD Upgrade

2010 by Chromalytic

Keep your gas products in spec! Monitor gas product purity, natural gas, and ambient air quality.

Sounds expensive and complicated to operate?

Not from SRI! The SRI Multiple Gas

Analyzer #1 uses just ONE gas sampling valve and TWO analytical columns to perform

the same separations that require multiple valves and columns in other systems.

Best of all, the Multiple Gas Analyzer #1 can achieve ppm to 100% concentrations with a single injection!

- ★ Separates multiple gases with a single injection
- ★ Very tolerant of user adjustments and timing variations
- ★ Simpler than other multi-gas capable GC systems
- ★ Multiple gas analysis in a compact unit

MUCH more !

Optional Vacuum Pump Interface for continuous sampling or from gas bags multiple injections

Substitute VICI u-TCD (for SRI TCD) for Capillary Column use;

~10ppM Sensitivity

Restek PLOT MS5A, ShinCarbon ST for Permanent Gases AND CO₂

Heated Valve Oven
original SRI dual column
1m (3') Molecular Sieve packed column separates H₂, O₂, N₂, CH₄ & CO

10 Port Gas Sampling Valve and 1mL Sample Loop

Sample In and Out for Valve Injection

On-column Injector

original SRI dual packed columns
2m (6') Silica Gel packed column separates CO₂ & C₂-C₆

XXX

TCD Detector - universal response, 250ppm to 100% detection range

Optional

FID Detector - Hydrocarbon selectivity, 5ppm detection limits Methanizer in FID body- converts CO and CO₂ to Methane for FID detection

Temperature Programmable Column Oven



ALSO add VICI Mass Flow Controller

replace packed columns with microPacked/PLOT column

Specifically Designed for
Separation of Whole Gas Components and Natural Gas Samples

Multiple Gas Analyzer #1



Separating out the hydrocarbon components of natural gas facilitates accurate BTU quantification. This compositional analysis of a natural gas standard by an SRI Multiple Gas Analyzer #1 shows good separation up to the pentanes. Performing compositional analyses of natural gas product before and after refining helps to maximize process efficiency and profit.



The same instrument produced this chromatogram, separating a sample mix of 1% fixed gas standard and ethane. With the built-in PeakSimple data system, the gas sampling valve was programmed to inject the sample loop contents into the carrier gas stream at 5 seconds, then rotate back at 6 minutes, after CO elution.

The basic Multiple Gas Analyzer #1 has a TCD detector only; this model provides analyses in the 250ppm to 100% range for fixed and natural gases. A second option is a TCD, Methanizer, and FID detector combination which adds 5ppm detection limits for CO, CO₂, and all hydrocarbon peaks; this model is useful for air quality monitoring and other applications. A third option is a TCD-HID detector combination, for detection limits in

the 10ppm range for all analytes...the HID even sees hydrogen! Since we build each GC from the boards up, the Multiple Gas Analyzer #1 may be further customized to suit your application needs. With the optional built-in "whisper-quiet" air compressor, the Multiple Gas Analyzer #1 can be used with the SRI H₂-50 hydrogen generator to separate multiple gases anywhere, without using compressed gas cylinders!



- 8610-0070 Multiple Gas Analyzer #1 GC with TCD detector
- 8610-0071 Multiple Gas Analyzer #1 GC with TCD, Methanizer, FID & built-in Air Compressor
- 8610-0072 Multiple Gas Analyzer #1 GC with TCD & HID detectors
- 8690-0070 Built-in Air Compressor, 120 VAC
- 8690-2270 Built-in Air Compressor, 220 VAC

Dual Cell Microvolume Thermal Conductivity Detector

- Stand-alone unit
- Optimized for capillary chromatography
- Thermal stability to $\pm 0.02^\circ\text{C}$
- Dual filaments capable of independent or referenced (differential) operation



The Valco Microvolume Thermal Conductivity Detector (TCD) is useful in a wide variety of capillary and packed column applications. Constant filament temperature control provides a linear dynamic range permitting measurement of a wide range of concentrations without the need for multiple standards or sample dilution.


Since the detector is non-destructive of the sample and contributes virtually no band spreading, it can be used in series with other detectors without affecting the performance characteristics of either.

Description


The detector consists of the cell housing and the electronics controller. The cell design permits mounting in virtually any orientation with no effect on performance. It can be installed easily on virtually any gas chromatograph, comprising a stand-alone unit requiring nothing else for operation but carrier gas flow.

Each of the two cell chambers is independent of the other, except for block temperature. Filaments can be replaced individually. Front panel controls set the temperature for the cell and for each filament. Since each detector cell can be operated separately or simultaneously, two analyses can be run using a single Valco TCD.

To insure compatibility with any system, two outputs are provided: 0-1 mV full scale attenuated output for recorders and 0-10 V full scale output for data acquisition systems.




Unleaded Gasoline
50 m x 320 micron columns, 0.06 μl valve injection



Gas Standard
30 m x 530 micron PLOT column, 100 μl valve injection

Independent vs. referenced operation

Specifications



Helium Blend

Sample size: 250 µl
 Sample concentration: 100 ppm each
 Column: 10' x 1/16" OD x 0.040" ID
 Molesieve 5Å, micropacked
 Column temp: 65°C
 Detector temp: 100°C
 Filament temp setting: 5.0
 Flow rate
 Channel A: 5.5 ml/min
 Channel B: 5.42 ml/min

Overall

Linear range 1 nanogram to 3 micrograms nC₄
 Minimum detectable approx. 50 picograms n-butane quantity
 Time constant < 150 milliseconds
 Cell temperature Automatic proportional control with control
 ±0.02°C stability
 Maximum cell 300°C
 temperature

Detector assembly

Dimensions 3.12" x 6" x 3.75" high
 (8 cm x 15 cm x 9 cm)
 Gas connections Valco 1/16" zero dead volume fittings
 Single multi-pin 5 foot cable supplied
 connector

Control unit

Dimensions 12" x 8" x 5" high
 (30 cm x 20 cm x 13 cm)
 Electrical connections... Single multi-pin connector
 Operator controls Cell temperature control (40-400°C)
 10-turn filament temperature
 potentiometers (A & B)
 10-turn coarse and fine baseline
 adjustment potentiometers (A & B)
 12 position recorder attenuator output
 switch (A, B, or A-B)
 Filament power on/off switch
 Indicator LEDs Detector heater "on"
 Filament power "on"
 Power requirements Universal 100-250 VAC
 50/60 Hz, 100W maximum

Product numbers

	110 VAC	230 VAC
Dual cell microvolume TCD with:		
nickel/iron filaments	TCD2-NIFE	TCD2-NIFE-220
tungsten/rhenium filaments	TCD2-WRE	TCD2-WRE-220

North America, South America, and Australia/Oceania contact:

VICI® Valco Instruments Co. Inc.

Europe, Asia, and Africa contact:

VICI® VICI AG International

Built-in "Whisper Quiet" Air Compressor

- **Built into the GC Chassis**
- **Powerful enough to supply FID air (300mL/minute)**
- **Convenient—Recommended for Field Work**

The Built-in "Whisper Quiet" Air Compressor provides an infinite and nearly silent supply of air for the FID, FID/DELCD, NPD, FPD, TID, or CCD detector. It mounts unobtrusively inside the 8610 or 310 GC chassis, and delivers unfiltered air to the detector.




With the built-in air compressor, no air cylinders are required. This simplifies field operations, and saves the expense of regularly replacing air cylinders.

8690-0070	Built-in "Whisper Quiet" Air Compressor
8690-2270	Same as above but 220 VAC

Vacuum Pump Interface

- **Draw air samples through traps or load the loop of a gas sampling valve**
- **Enables Data System Control of an external vacuum pump (included)**
- **Extremely reproducible flow through traps**

The Vacuum Pump Interface is a data system controlled main power outlet (120 or 220 VAC) on the side of an 8610 or 310 GC for an external vacuum pump. The PeakSimple data system can turn the power to this receptacle ON/OFF, thus controlling the vacuum pump.



Typically, the vacuum pump is used to draw gaseous samples through the traps for ambient air monitoring applications, or to load the loop of a gas sampling valve by pulling sample gas from a remote location.

Because the vacuum pump can be turned ON for a precise length of time, the gas flow through the traps is very reproducible (approximately 100mL/minute).

8690-0073	Vacuum Pump Interface and Pump
------------------	---------------------------------------

Molecular Sieve 5A PLOT Columns

Restek's molecular sieve 5A PLOT columns are designed for efficient separation of Ar/O₂ and other permanent gases, including CH₄, C₂H₆, and CO. Special coating and deactivation procedures ensure chromatographic efficiency and the integrity of the porous layer coating. Molecular sieves have very high retention, allowing separations of permanent gases at temperatures above ambient. Additionally, our unique immobilization process guarantees that the uniform particles remain adhered to the tubing—even after continuous valve-cycling.

Our revolutionary molecular sieve 5A PLOT columns separate Ar/O₂ and H₂/He at ambient temperature or above (see figure). These columns also are an excellent choice for rapid separation of permanent gases in refinery or natural gas.

Rt®-Msieve 5A Columns (fused silica PLOT)

ID	df	temp. limits	15-Meter	30-Meter
0.25mm	20μm	to 300°C	19773	
0.32mm	30μm	to 300°C	19720	19722
0.53mm	50μm	to 300°C	19721	19723

MXT®-Msieve 5A Columns (Siltek®-treated stainless steel PLOT)

Advantages of metal MXT® PLOT columns include:

- Can be made in small coil diameters—perfect for tight spaces.
- Will not spontaneously break, making them ideal for rugged environments.
- Designed for robust performance in process GCs and field instruments.
- Available in 3.5" coil diameter or 7" diameter 11-pin cage.

ID	df	temp. limits	15-Meter	3.5" coil 30-Meter	7" diameter 11-pin cage 30-Meter
0.25mm	20μm	to 300°C	79717		
0.53mm	50μm	to 300°C		79723-273	79723



advanced
technology

Details on **pages 106-107**.



did you know?

Rt®-Msieve 5A PLOT columns are designed for efficient separation of Ar/O₂ and other permanent gases, including CH₄, C₂H₆, and CO.



tech tip


Because molecular sieve materials are very hydrophilic, they will adsorb water from the sample or carrier gas. Water contamination can have a detrimental effect on peak symmetry and can reduce the resolution of all compounds. If water contamination occurs, reactivate your Rt®-Msieve 5A PLOT column by conditioning at 300 °C with dry carrier gas flow for 3 hours.



tech tip

Carbon dioxide will not elute from molecular sieve columns. Rt®-Q-BOND is a good choice for this analysis.

Permanent gases on an Rt®-Msieve 5A PLOT column.



Column: Rt®-Msieve 5A, 30m, 0.53mm ID, 50μm (cat.# 19723)
Sample: permanent gases (ppm)
Inj.: 5μL sample loop, 6-port Valco® valve, valve temp.: ambient
Inj. temp.: 200°C
Carrier gas: helium, constant flow
Flow: 5mL/min.
Oven temp.: 27°C (hold 5 min.) to 100°C @ 10°C/min. (hold 5 min.)
Det.: Valco® helium ionization detector @ 150°C

did you know?

ShinCarbon ST micropacked columns are another alternative for analyzing permanent gases.

See **page 130** for information.

RESTEK


CHROMALYTIC +61(0)3 9762 2034
ECHnology Pty Ltd

Australian Distributors
Importers & Manufacturers
www.chromtech.net.au

11/12

www.restek.com 109

Water and ethanol in acetone on an Rt®-Q-BOND PLOT column.



Column: Rt®-Q-BOND, 30m, 0.53mm ID, 20µm (cat.# 19742)
 Sample: 0.5% water and ethanol in acetone
 Inj.: 3µL split (split ratio 11:1), 4mm single gooseneck liner w/ wool (cat.# 22405)
 Inj. temp.: 250°C
 Carrier gas: helium, constant flow
 Linear velocity: 28.7cm/sec. @ 200°C
 Oven temp.: 200°C, isothermal
 Det.: TCD @ 260°C

PLOT Column Particle Trap

- Includes two Press-Tight® connectors and a 2.5 m column.
- Protects detector and valves; connects between column and detector or valve.
- Eliminates detector spikes and scratches in valve rotors.

The technology used to adhere particles in PLOT columns is excellent; however, there is still a possibility for particles to dislodge when extreme pressure shocks and gas flow changes are anticipated. This sometimes happens when valve backflush or MS detection is used. In those extreme cases, using particle traps is recommended.

Description	qty.	cat.#	price
PLOT Column Particle Trap, 2.5m, 0.32mm ID with 2 Press-Tight Connectors	ea.	19753	
PLOT Column Particle Trap, 2.5m, 0.53mm ID with 2 Press-Tight Connectors	ea.	19754	



Particle Trap

Restek Customer Service

In the U.S.

Call: 800-356-1688 (ext. 3) or 814-353-1300 (ext. 3)


Monday–Friday 8:00 a.m.–6:00 p.m. ET

Fax: 814-353-1309—24-hours a day

Online: www.restek.com—24-hours a day

Outside the U.S.

Contact your Restek representative:
 Refer to our list on pages 4-5 or visit our website at www.restek.com



Melissa Decker, Customer Service



www.restek.com 111



it's a fact

ShinCarbon ST is an ideal packing material for permanent gases, low molecular weight hydrocarbons, sulfur dioxide, and Freon® gases.

also available

For adapter kits for installing packed/micropacked columns, see **page 133**.

Permanent Gases & Hydrocarbon Analysis

ShinCarbon ST Packed/Micropacked Columns

- Separate permanent gases, including CO/CO₂, without cryogenic cooling.
- Rapid separations of permanent gas/light hydrocarbon mixtures.
- Excellent compatibility with most GC detectors—minimal bleed, minimal baseline rise.
- Preconditioned, less than 30 minutes to stabilize.

Analyze oxygen, nitrogen, methane, carbon monoxide, and carbon dioxide with one column and at room temperature. ShinCarbon ST material, a high surface area carbon molecular sieve (~1,500 m²/g), is the ideal medium for separating gases and highly volatile compounds by gas solid chromatography (GSC). The rapid, above-ambient analyses these columns provide will be a great convenience. Excellent thermal stability of the high surface area carbon, combined with careful conditioning during column manufacturing, ensures low-bleed operation and rapid stabilization when installing a new column. Custom-made ShinCarbon ST columns are available on request.

ShinCarbon ST is a highly stable material. Its 330 °C upper temperature limit minimizes bleed and baseline rise during temperature programming, making the material compatible with most detection systems used for gas analysis, including TCD or HID. All ShinCarbon ST columns are fully conditioned in an oxygen/moisture free environment to prevent contamination. This minimizes stabilization time (less than 30 minutes) when installing a new column which, in turn, minimizes downtime.

ShinCarbon ST 80/100 Columns (packed) (SilcoSmooth® Stainless Steel)*

OD	ID	2-Meter
1/8" SilcoSmooth	2.0mm	80486-

ShinCarbon ST 100/120 Columns (micropacked) (SilcoSmooth® Stainless Steel)**

OD	ID	1-Meter	2-Meter
1/16"	1.0mm	19809 \$245	19808
0.95mm	0.75mm	19810 \$245	

*Please add column instrument configuration suffix number to cat.# when ordering. See chart on the next page.

**Does not include column nuts and ferrules. Optional installation kits can be ordered separately—see page 133.


Chromatogram Search Tool

Search by compound name, synonym, CAS # or keyword

www.restek.com/chromatograms



Separate permanent gases in 10 minutes, without cryogenics.



Column: ShinCarbon ST, 100/120 mesh, 2m, 1mm ID micropacked (cat.# 19808)
Sample: 5µL permanent gases mix, approx. 5 mole % each
Inj. temp.: 100°C
Carrier gas: helium
Flow rate: 10mL/min.
Oven temp.: 40°C (hold 3 min.) to 250°C @ 8°C/min. (hold 10 min.)
Det. HID @ 200°C

please note

For additional chromatograms for ShinCarbon ST columns, see **pages 647, 649, and 652**.