

Some Hint on How to Make a Standard UHPLC Column with + 300 000 Theoretical Plates/Meter



BIOTECH AB
WWW.BIOTECH-SWEDEN.COM

Norikazu Nagae¹, Tomoyasu Tsukamoto¹, Henrik Svennberg²
¹ ChromaNik Technologies Inc. Namiyoke, Minato-ku, Osaka Japan
² Biotech AB, Innovation, Onsala, Sweden

ChromaNik
ChromaNik Technologies Inc.
www.chromanik.co.jp




The race to achieve the highest number of theoretical plate counts has been going for many years already, even though it could and should be argued that a good separation is mainly governed by good selectivity, a reasonable retention and then high column efficiency. However, the race has still been going, and this presentation will show how + 300 000 theoretical plates/meter can be achieved in a standard column for UHPLC. The importance of a good bonding and end capping technique will be discussed and hints about how this is made is revealed. The importance of base particle choice with a good heat transfer capacity will be shown and the importance of small extra column peak broadening will be emphasized. The + 300 000 theoretical plates/meter is shown in two different columns lengths with the inner diameter of 2.1 mm, 150 and 50 mm long respectively.

Effect of end-capping

End-capping with hexamethyltrisiloxane and TMS on C18 silica



Another end-capping with heating on C18 silica, reduce of silanol groups



Comparison of 2 kinds of end-capping

Conventional end-capping: Stationary phase with silanol groups.

Both Sunniest end-capping and Silanol Activity Control: Stationary phase with siloxane bond, more hydrophobic.

Carbon loading and retention factor (reaction at two different temperatures)

End-capping reaction temperature	X °C (more than 200 °C)	(X-40) °C
Carbon loading of only C18	7.0%	7.0%
Carbon loading after end-capping	7.3% (cut off some C18 chains by heat)	7.7% (not cut off C18 chain by heat)
Silanol activity control	Yes	No
Retention factor (k)	10.6	9.5

X °C was too high to keep carbon loading of initial C18. (X-40) °C could keep carbon loading. As a result, carbon loading after end-capping at (X-40) °C was higher than that at X °C.

Regarding retention factor, reversed result was obtained. Less silanol group made hydrophobicity of stationary phase high.



Relationship between retention factor and carbon loading

SunShell C18, 2.6 μm	Retention factor ^a	Carbon loading ^b (%)	Specific surface area ^c (m ² /g)
SunShell C18, 2.6 μm	10.4	7.3	125
Ascentis Express C18, 2.7 μm	9.7	8.0	133
PoroShell C18 EC, 2.7 μm	9.0	8.5	135
Cortecs C18, 2.7 μm	7.7	7.3	113
Accucore C18, 2.6 μm	7.4	8.8	130
Kinetex C18, 2.6 μm	5.4	4.9	102

SunShell C18 showed the largest retention factor although carbon loading was not the highest.


It is considered that less residual silanol group made hydrophobicity of a stationary phase high, then made retention factor large.

Comparison of amitriptyline peak



Effect of frictional heating


Frictional heating phenomenon



Column, core shell C18 2.6 μm, 150 x 4.6 mm; mobile phase, methanol/flow rate, 5 mL/min; column back pressure, 70 MPa; temperature, 23 °C (room temperature).

Temperature rose 14 °C at the outlet of the column by effect of frictional heating.

Efficiency loss by thermal friction*




Column dimension, 100 x 2.1 mm; mobile phase, CH₃CN/H₂O=60/40; temperature, 40 °C; sample, acenaphthene; UHPLC instrument, X-LC with 0.1 mm i.d. connecting tube (JASCO).

Efficiency and back pressure

Acenaphthene	Retention time (min)	Tailing factor	Pressure (MPa)	Theoretical plate (per column)
SunShell C18	29.721	1.04	48.5	613
Company S Core Shell C18	25.533	1.19	53.5	477
Company P Core Shell C18	24.700	0.97	53.8	458
Company W totally porous hybrid C18	14.511	1.01	54.1	269
Company S totally porous monodisperse C18	26.592	1.05	43.9	605

Peak shape and efficiency of pyridine and amitriptyline




Influence of frictional heating is most difficult to undergo for a core shell silica. When using a fine particle such as 2 μm or smaller than 2 μm C18 packings under high pressure condition, a core shell C18 with dense end-capping showed the highest efficiency for not only neutral compounds but also basic compounds.

Effect of extra-column volume


Comparison of an inner diameter of tubing in the UHPLC instrument

Connecting tube used	Instrument	Injector to column	Column to flow cell	Tube volume
X-LC with SUS	SUS 0.1 mm i.d. x 300 mm	PeekSil 0.1 mm i.d. x 200 mm		3.93 μL
X-LC with MarvelX	MarvelX 0.075 mm i.d. x 350 mm	MarvelX 0.075 mm i.d. x 150 mm		2.21 μL


SunShell C18 2 μm, 50 x 2.1 mm



SunShell C18 2.6 μm, 150 x 4.6 mm



SunShell C18 2 μm, 150 x 2.1 mm



Conclusion

- End-capping done by both bonding of a silylation agent and conversion from silanol groups to siloxane bond led to not only high efficiency for basic compounds but also increase in retention.
- Serious losses in column efficiency was caused by frictional heating effects when very fine particles were packed into a column and it was operated at high back pressure. A core shell particle showed the least losses in column efficiency, while a totally porous hybrid particle showed large losses.
- Effect of extra-column volume is important for a column packed with very fine particles. When 0.075 mm i.d. connecting tubes were used, SunShell C18 2 μm column showed + 300,000 theoretical plates/meter.
- In case of a low dead volume column sized 4.6 mm i.d. and 150 mm length, uracil peak showed + 300,000 theoretical plates/meter even if a particle size was 2.6 μm.
- In case of 150 mm length column, efficiency was + 300,000 theoretical plates/meter even if 0.1 mm i.d. connecting tubes were used.