

C18, C18-WP, HFC18-16, RP-AQUA, C8, C30, PFP, PFP&C18, Phenyl, Biphenyl, Cyano, C8-30HT, C4-100, HILIC-Amide, HILIC-S and 2-EP


SunShell

2 µm, 2.6 µm, 3.4 µm, 3.5 µm and 5 µm HPLC column



New Product
SunShell Cyano
Diisopropylcyanopropyl

ChromaNik Technologies Inc.



“SunShell “ is a core shell silica column made by ChromaNik Technologies.

The next generation to Core Shell particle



SUNSHELL

Superficially porous silica

Features of SunShell

- * 1.2 μm , 1.6 μm , 2.3 μm , 3.0 μm and 3.4 μm of core and 0.4 μm , 0.5 μm , 0.2 μm and 0.6 μm of superficially porous silica layer
- * Higher efficiency and higher throughput to compare with totally porous silica with same size
- * Same chemistry as Sunniest technology (reference page 6)
- * Good peak shape for all compounds such as basic, acidic and chelating compounds
- * High stability (pH range for SunShell C18, 1.5 to 10)
- * Low bleeding


SUS column hardware is used in parallel by two companies.




✘ Metal free columns and nano- and microcolumns are listed on page 24




Superficially porous particle




✘ SunShell guard cartridge columns are listed on page 25




Particle size: 2.6 μm
Pore diameter: 16 nm



Particle size: 3.4 μm
Pore diameter: 30 nm



Particle size: 2.6 μm
Pore diameter: 100 nm




Electron micrograph of core shell silica

Core shell and obser

HROMalytic +61(0)3 9762 2034
ECHnology Pty Ltd
 Website NEW : www.chromalytic.net.au E-mail : info@chromotech.net.au Tel: 03 9762 2034 . . . in AUSTRALIA

Australian Distributors
 Importers & Manufacturers
www.chromotech.net.au

on treatment,



Schematic diagram of a core shell silica particle, 2.0, 2.6, 3.4 3.5 and 4.6 μm

HROMalytic +61(0)3 9762 2034
ECHnology Pty Ltd
 Website NEW : www.chromalytic.net.au E-mail : info@chromotech.net.au Tel: 03 9762 2034 . . . in AUSTRALIA

Australian Distributors
 Importers & Manufacturers
www.chromotech.net.au





Why does a 2.6 μm core shell particle show the same performance as a sub 2 μm particle?

Van Deemter Equation

$$H = A d_p + B \frac{D_m}{u} + C \frac{d_p^2}{D_m} u$$


- A term : Eddy diffusion (dp is particle diameter)
- B term : Longitudinal diffusion (Dm is diffusion coefficient)
- C term : Mass transfer




SunShell C18 shows same efficiency as a sub 2 μm C18. In comparison between fully porous 2.6 μm and core shell 2.6 μm (SunShell), SunShell shows lower values for A term, B term and C term of Van Deemter equation. The core shell structure leads higher performance to compare with the fully porous structure.

All terms in Van Deemter Equation reduce.


A term



Wide particle distribution (Conventional silica gel D₉₀/D₁₀=1.50)



Narrow particle distribution (core shell silica D₉₀/D₁₀=1.15)




Packing state of core shell and fully porous silica

The size distribution of a core shell (SunShell) particle is much narrower than that of a conventional totally porous particle, so that the space among particles in the column reduces and efficiency increases by reducing Eddy Diffusion (multi-path diffusion) as the A term in Van Deemter Equation.


Diffusion of a solute is blocked by the existence of a core, so that a solute diffuses less in a core shell silica column than in a totally porous silica column. Consequently B term in Van Deemter Equation reduces in the core shell silica column.

B term


Totally porous silica A solute diffuses in a pore as well as outside of particles.




Core shell silica A core without pores blocks diffusion of a solute.




Difference of longitudinal diffusion





As shown in the left figure, a core shell particle has a core so that the diffusion path of samples shortens and mass transfer becomes fast. This means that the C term in Van Deemter Equation reduces. In other words, HETP (theoretical plate) is kept even if flow rate increases. A 2.6 μm core shell particle shows as same column efficiency as a totally porous sub-2 μm particle.


Considering diffusion of solute within pore



The left figure shows that a diffusion width of a sample in a 2.6 μm core shell particle and a 2 μm totally porous particle. Samples or solutes enter into the particle and move by diffusion, then they go out of a particle. In this moment, sample peak width is broadened. This broadening width is statistically same for 2.6 μm core shell particle and 2 μm fully porous particle. The 2.6 μm core shell particle is superficially porous, so that the diffusion width becomes narrower than particle size. Same diffusion means same efficiency.

Comparison of Performance by Plate/Pressure

	Plate	Back press. (MPa)	Plate/back press.
Sunniest C18 –HT 2.0 μm	9,900	16.7	593
Brand A C18 1.9 μm	7,660	16.3	470
Brand B C18 1.8 μm	10,100	19.6	515
Brand C C18 1.7 μm	11,140	32.0	348
SunShell C18 2.6 μm	9,600	9.7	990



Back pressure and theoretical plate were compared for 2 μm and sub 2 μm C18 and 2.6 μm SunShell C18. All columns showed almost the same theoretical plate except for brand A C18 1.9 μm. However back pressure was not same. Especially Brand C C18 1.7 μm showed the highest back pressure. And SunShell C18 2.6 μm showed the lowest back pressure. On the comparison of theoretical plate per back pressure, SunShell indicated the largest value. This is a big advantage.

Column: 50 x 2.1 mm C18, Mobile phase: Acetonitrile/water=(70/30), Temperature: 25 °C




SUNSHELL


STATIONARY PHASE

Reversed phase


C18, C18-WP (7 page, 16 page, 20 page, SunShell Bio)




RP-AQUA, C30 (16 page, 19 page)




C8, C8-30HT (16 page, 20 page, 21 page)




C4-100 (20 page, 21 page, SunShell Bio)




Phenyl (16 page)




Phenyl (SunShell Bio)



PFP (16 page)




PFP&C18 (16 page)




**All reversed phases except for PFP and PFP&C18 were end-capped at high temperature using SunShell Endcapping technique.


HFC18-16 (20 page)



Biphenyl (16 page)




Cyano (16 page)




HILIC and SFC


2EP (SunShell SFC)




Diol (SunShell SFC, SunShell HILIC)



HILIC-Amide (23 page, SunShell HILIC)



HILIC-Silica (23 page, SunShell SFC, SunShell HILIC)




*Stationary phase for both SFC and HILIC was not end-capped.

SunShell Bio (1000Å), SunShell SFC and SunShell HILIC see individual catalogue.




Unique end-capping by new concept

This figure shows comparison of hydrophobicity between two C18 stationary phases. We developed silanol activity control technique which was a reaction at extremely high temperature. This technique makes residual silanol groups change to siloxane bond. The upper one is a C18 phase with conventional end-capping and the lower one is a C18 phase with both Sunniest end-capping and silanol activity control. A residual silanol group contributes as a polar site and makes hydrophobicity of stationary phase decrease. On the other hand siloxane bond in the lower one doesn't make hydrophobicity decrease. Consequently the lower one is more hydrophobic than the upper one.



End-capping method

- 1) Unique end-capping reagent
 <<Hexamethyltrisiloxane>>
- 2) Secondly TMS end-capping



An end-capping of hexamethyltrisiloxane works as an arm. This arm moves like a Geometrid caterpillar, so that a functional group on the tip of the arm can bond with a silanol group which is located anywhere. Finally TMS reagent is bonded to a remaining silanol group.

SunShell C18, 2 μm, 2.6 μm, 3.5 μm, 5 μm




Characteristics of SunShell C18

	Core shell silica			C18 (USP L1)				
	Particle size	Pore diameter	Specific surface area	Carbon content	Bonded phase	End-capping	Maximum operating pressure ^a	Available pH range
SunShell C18	2.0 μm	9 nm	120 m ² /g	6.5%	C18	Sunniest endcapping	100 MPa or 14504 psi	1.5 - 10
SunShell C18	2.6 μm	9 nm	150 m ² /g	7%	C18	Sunniest endcapping	60 MPa or 8,570 psi	1.5 - 10
NEW SunShell C18	3.5 μm	9 nm	120 m ² /g	6.5%	C18	Sunniest endcapping	60 MPa or 8,570 psi	1.5 - 10
SunShell C18	4.6 μm	9 nm	90 m ² /g	5.5%	C18	Sunniest endcapping	50 MPa or 7,141 psi	1.5 - 10


a) Unless otherwise specified in the column test report

Core Shell particle shows 1.4 to 1.5 times higher plate than fully porous particle.



Theoretical plate and tailing factor


Used columns: SunShell C18 2 μm, Ascentis Express C18 2 μm, Kinetex C18 1.7 μm, Acquity BEH C18 1.7 μm, Titan C18 1.9 μm



*Ascentis Express is a registered trade mark of Sigma Aldrich. Titan is a registered trade mark of Sigma Aldrich. Comparative separations may not be representative of all applications.



Comparison of retention and plate using HPLC




When Semi-Micro LC is used, 5 μm shows 24,000plates, 3.5 μm shows 31,000plates, 2.6 μm shows 38,000plates.

Column size: 150 x 4.6 mm
 Mobile phase: CH₃OH/H₂O=75/25
 Flow rate: 1.0 mL/min
 Temperature: 40 °C
 Sample: 1 = Uracil
 2 = Caffeine
 3 = Phenol
 4 = Butylbenzene
 5 = o-Terphenyl
 6 = Amylbenzene
 7 = Triphenylene
 HPLC: Hitachi LaChrom ELITE
 (Tubing, 0.25 mm i.d.)

	Totally porous silica Sunniest C18, 5 μm		Core shell silica SunShell C18, 2.6 μm		Core shell silica SunShell C18, 3.5 μm		Core shell silica SunShell C18, 5 μm	
	Retention time (t _R)	Retention factor (k)	Retention time (t _R)	Retention factor (k)	Retention time (t _R)	Retention factor (k)	Retention time (t _R)	Retention factor (k)
1) Uracil	1.70	0	1.34	0	1.33	0	1.30	0
6) Amylbenzene	19.96	10.74	16.56	11.36	13.90	9.45	12.43	8.56
Relative value of Amylbenzene	100%	100%	83%	106%	70%	88%	63%	80%

Comparison between porous C18 and SunShell C18 3.5 μm column





SunShell C18 3.5 μm showed the highest theoretical plate and the lowest back pressure. Regarding retention factor, SunShell C18 3.5 μm was the middle.

Column size: 150 x 4.6 mm
 Mobile phase: CH₃OH/H₂O=75/25
 Flow rate: 1.0 mL/min
 Temperature: 40 °C
 Sample: 1 = Uracil
 2 = Caffeine
 3 = Phenol
 4 = Butylbenzene
 5 = o-Terphenyl
 6 = Amylbenzene
 7 = Triphenylene
 HPLC: Conventional HPLC instrument
 (Tubing, 0.25 mm i.d.)



Comparison between normal and semi-micro HPLC



Comparison of chromatograms

Column: SunShell C18, 5 μ m 250 x 4.6 mm
 Mobile phase: CH₃CN/H₂O= 70/30
 Flow rate: 1.0 mL/min
 Temperature: 40 °C
 Pressure: 6.7 MPa
 Detection: UV@250 nm
 Sample: 1 = Uracil
 2 = Toluene
 3 = Acenaphthene
 4 = Butylbenzene




Semi-micro HPLC derives near 100% performance of a core shell column. Even if normal HPLC is used, it derives 80% performance except for a narrow peak whose width is less than 5 second

HPLC: Hitachi LaChrom ELITE

Relationship between Peak width and theoretical plate

Effect of inner diameter of tubing




Average of theoretical plate (n=3)

Inner diameter of tubing	0.06mm	0.1mm	0.25mm
Peak (1)	792	785	246
Peak (2)	7790	7652	3535
Peak (3)	10704	10345	7998
Peak (4)	10113	9772	7689

Column: SunShell C18, 2.6 μ m 50 x 2.1 mm
 Mobile phase: CH₃CN/H₂O=60/40
 Flow rate: 0.3 mL/min Temperature: Ambient
 Tube length: 30 cm (Peek, from the column to the flow cell)
 Instrument: X-LC(JASCO) Response time: 0.01 sec

The above theoretical plate was compared changing the inner diameter of tubing between a column and a flow cell of the detector. A tubing with a large inner diameter has a large dead volume, so that it makes the peak width be wide. As a result, theoretical plate decreases. I recommend to use the tubing with 0.1 mm or less than 0.1 mm inner diameter for core shell columns.

Effect of response time of detector



Column: SunShell C18, 2.6 μ m 100 x 4.6 mm
 Mobile phase: CH₃CN/H₂O=60/40
 Flow rate: 1.8 mL/min Temperature: Ambient
 Sample: Toluene Tube: i.d.0.1mm x 20 cm Peeksil
 Instrument: X-LC(JASCO)


The response time of a detector is important. Regarding uracil, the real peak width is less than 0.8 sec. When the peak width is less than 1 sec, 0.03 sec of response time is needed. Furthermore, the sampling rate of an integrator should be set to be 0.1 sec.

SunShell C18 2 μm

Comparison of core shell 2 μm and totally porous sub 2 μm

Used columns: SunShell C18 2 μm, Ascentis Express C18 2 μm, Kinetex C18 1.7 μm, Acquity BEH C18 1.7 μm, Titan C18 1.9 μm


Separation of standard samples



Column: SunShell C18, 2 μm 100 x 2.1 mm
 Company S Core Shell C18, 2 μm 100 x 2.1 mm
 Company P Core Shell C18, 1.7 μm 100 x 2.1 mm
 Company W C18, 1.7 μm 100 x 2.1 mm
 Company S Monodisperse C18, 1.9 μm 100 x 2.1 mm
 Mobile phase: CH₃OH/H₂O=75/25
 Flow rate: 0.2 mL/min
 Temperature: 40 °C
 Sample: 1 = Uracil, 2 = Caffeine, 3 = Phenol, 4 = Butylbenzene
 5 = o-Terphenyl, 6 = Amylbenzene, 7 = Triphenylene


	Hydrogen bonding (Caffeine/Phenol)	Hydrophobicity (Amylbenzene/Butylbenzene)	Steric selectivity (Triphenylene/o-Terphenyl)
SunShell C18	0.43	1.59	1.41
Company S Core Shell C18	0.37	1.59	1.38
Company P Core Shell C18	0.45	1.57	1.17
Company W C18	0.35	1.55	1.30
Company S Monodisperse C18	0.53	1.58	1.16

Comparison of Pyridine (2) as a basic compound




Column dimension: 100 x 2.1 mm
 Mobile phase: CH₃OH/H₂O=30/70
 Flow rate: 0.2 mL/min
 Temperature: 40 °C
 Detection: UV@250nm
 Sample: 1 = Uracil
 2 = Pyridine
 3 = Phenol

Comparison of Oxine (1) as a metal chelating compound



Column dimension: 100 x 2.1 mm
 Mobile phase: CH₃CN/20mM H₃PO₄=10/90
 Flow rate: 0.2 mL/min
 Temperature: 40 °C
 Detection: UV@250nm
 Sample: 1 = 8-Quinolol (Oxine)
 2 = Caffeine


Comparison of Formic acid (1) as an acidic compound



Column dimension: 100 x 2.1 mm
 Mobile phase: CH₃CN/0.1% H₃PO₄=2/98
 Flow rate: 0.2 mL/min
 Temperature: 40 °C
 Detection: UV@210nm
 Sample: 1 = Formic acid
 2 = Acetic acid
 3 = Propionic Acid


SunShell C18 2 μm

Comparison of Amitriptyline (4) as a strong basic compound




Column dimension: 100 x 2.1 mm
 Mobile phase: CH₃CN/20 mM Phosphate buffer pH 7.0=60/40
 Flow rate: 0.3 mL/min
 Temperature: 40 °C
 Detection: UV@250 nm
 Sample: 1 = Uracil, 2 = Propranolol, 3 = Nortriptyline, 4 = Amitriptyline

Column dimension: 100 x 2.1 mm
 Mobile phase: CH₃CN/10 mM ammonium acetate pH 6.8=40/60
 Flow rate: 0.3 mL/min
 Temperature: 40 °C
 Detection: UV@250 nm
 Sample: 1 = Uracil, 2 = Propranolol, 3 = Nortriptyline, 4 = Amitriptyline



Decreasing of theoretical plate due to frictional heating effect



Core shell silica has a solid core (non-porous silica), so that thermal conductivity is high in the column. There is no influence of reducing theoretical plate by frictional heating.

Regarding totally porous hybrid silica, not only totally porous structure but also including ethylene groups make thermal conductivity be low in the column. It is considered that frictional heating deflects thermal distribution in the column and theoretical plate decreases..


Column: 100 x 2.1 mm
 Mobile phase: CH₃CN/H₂O=60/40
 Temperature: 40 °C
 Sample: Acenaphthene,

Comparison of core shell 2.6 μm columns

Comparison of standard samples among core shell C18s


Used columns

1. Kinetex C18, 2.6 μm
2. Accucore C18, 2.6 μm
3. PoroShell C18 EC, 2.7 μm
4. Ascentis Express C18, 2.7 μm
5. Cortecs C18, 2.7 μm
6. SunShell C18, 2.6 μm




Retention of standard samples and back pressure were compared for six kinds of core shell type C18s. Company A C18 showed only a half retention to compare with SunShell C18. Steric selectivity becomes large when ligand density on the surface is high. SunShell C18 has the largest steric selectivity so that it has the highest ligand density. This leads the longest retention time.

Comparison of pyridine




Residual silanol groups make pyridine be tailing under methanol/water mobile phase condition. SunShell C18 shows a sharp peak for pyridine.

Comparison of Oxine




8-Quinolinol (Oxine) is a metal chelating compound. Metal impurities in the core shell particle leads the tailing for oxine peak.

Comparison of formic acid




Formic acid is used as an indicator for a acidic inertness. SunShell and Company A and C C18 show a sharp peak.



Loading capacity of amitriptyline as a basic compound

Amitriptyline overloads much more at acetonitrile/buffer mobile phase than methanol/buffer. Three kinds of core shell C18s were compared loading capacity of amitriptyline at three different mobile phases.

Common condition: Column dimension, 150 x 4.6 mm, flow rate; 1.0 mL/min, temperature; 40 °C



Physical properties

	Carbon loading (%)	Specific surface area ^a (m ² /g)	Pore volume ^a (mL)	Pore diameter ^a (nm)
SunShell C18	7.3 (7) ^b	125 (150) ^b	0.261	8.34 (9) ^b
Ascentis Express C18	8.0	133 (150) ^b	0.278	8.20 (9) ^b
PoroShell C18 EC	8.5 (8) ^b	135 (130) ^b	0.414	12.3 (12) ^b
Accucore C18	8.8 (9) ^b	130 (130) ^b	0.273	8.39 (8) ^b
Cortecs C18	7.3 (6.6) ^b	113	0.264	9.32
Kinetex C18	4.9 (12 effective) ^b	102 (200 effective) ^b	0.237	9.25 (10) ^b

- a. Measured after sintered at 600 degree Celsius for 8 hours.
- b. Value cited in company brochure or literature.


Comparison column

1. Kinetex C18, 2.6 μm
2. Accucore C18, 2.6 μm
3. PoroShell C18 EC, 2.7 μm
4. Ascentis Express C18, 2.7 μm
5. Cortecs C18 2.7 μm
6. SunShell C18, 2.6 μm




All columns are core shell type. All columns sized 150 x 4.6 mm except for company E show 38,000 to 40,000 plates for a neutral compound. However regarding a basic compound like amitriptyline, SunShell C18 and company C C18 showed a good peak, while Company A, B and D C18 showed a poor peak. Company A C18 overloaded at more than 0.01 μg of amitriptyline while SunShell C18 overloaded at more than from 0.3 to 1 μg of amitriptyline. Surprisingly loading capacity of company A C18 was only one hundredth to compare with SunShell C18 under acetonitrile/20mM phosphate buffer pH7.0=(60:40) mobile phase. Company D C18 always showed poor peak of amitriptyline.

◆ Evaluation of Stability



Stability under acidic pH condition was evaluated at 80 °C using acetonitrile/1% trifluoroacetic acid solution (10:90).

★ Sunshell C18 has kept 90% retention for 100 hours under such a severe condition. SunShell C18 is 5 to 10 times more stable than the other core shell C18.




Stability under basic pH condition was evaluated at 50 °C using methanol/Sodium borate buffer pH 10 (30:70) as a mobile phase. Sodium borate is used as a alkaline standard solution for pH meter, so that its buffer capacity is high.


Elevated temperature of 10 °C makes column life be one third. The other company shows stability test at ambient (room temperature). If room temperature is 25 °C, column life at room temperature (25 °C) is sixteen times longer than that at 50 °C.

★ SunShell C18 is enough stable even if it is used under pH 10 condition. Regarding stability under basic pH condition, there is little C18 column like SunShell C18 except for hybrid type C18. It is considered that our end-capping technique leads high stability.

◆ Continuous analysis under pH9.5 condition




◆ Comparison of particle size



*Measured using Beckman Coulter Multisizer 3 after C18 materials were sintered at 600 degree Celsius for 8 hours. The measured value of each sintered core shell silica is considered to be different from that of the original core shell silica a. Median particle size


SunShell

Peptides (using the 1.0 mm i.d. column)



Column: SunShell RP-AQUA, 2.6 μ m 100 x 1.0 mm
 Mobile phase: A) 0.1 % trifluoroacetic acid (TFA) in water
 B) 0.08 % trifluoroacetic acid (TFA) in acetonitrile
 %B 10% to 30% in 25 min
 Flow rate: 0.15 mL / min
 Temperature: 60 $^{\circ}$ C
 Detection: UV@214 nm
 Sample: 1 = Gly-Tyr, 2 = Val-Tyr-Val, 3 = Met enkephalin,
 4 = Leu enkephalin, 5 = Angiotensin II
 (HPLC peptide standard mixture by Sigma-Aldrich)

Amino Acids derivatized with OPA and FMOC




Column: SunShell C18 2.6 μ m, 150 x 2.1 mm
 Mobile phase: A) 10mM Na_2PO_4 + 10mM $\text{Na}_2\text{B}_4\text{O}_7$ + 0.5mM NaN_3 (pH7.8)
 B) Acetonitrile/Methanol/Water (45/45/10 %V)

Time (min)	0	0.4	12.8	13.8
%B	5	5	50	100

Flow rate: 0.61 mL/min
 Temperature: 40 $^{\circ}$ C
 Detection: UV@338 nm
 Sample: 1=Aspartic acid, 2=Glutamic acid, 3=Serine, 4=Histidine, 5=Glycine,
 6=Threonine, 7=Arginine, 8=Alanine, 9=Tyrosine, 10=Valine, 11=Methionine,
 12=Tryptophan, 13=Pheylalanine, 14=Isoleucine, 15=Leucine, 16=Lysine,
 17=Proline


Dansylated estrogen hormones



Column: SunShell C18 2.6 μ m, 100 x 2.1 mm
 Mobile phase:
 A) H_2O with 0.1% formic acid.
 B) CH_3CN with 0.1% formic acid.
 Gradient program:
 0 - 0.5 min: 10% B
 0.51 - 3.0 min: 10 - 72% B
 3.01 - 6.0 min: 72% B
 6.01 - 7.0 min: 72 - 100% B
 7.01 - 10.0 min: 100% B
 Flow rate: 0.45 mL/min.
 Temperature: 40 $^{\circ}$ C
 Detection: MS(sim), m/z, 522.20, 506.20, 504.20
 Samples: 1. Dansylated estriol, 2. Dansylated 17beta-estradiol,
 3. Dansylated 17alpha-estradiol, 4. Dansylated estrone

Courtesy of Department of Chemistry & Biochemistry, The University of Texas at Arlington

Oolong tea



Column: SunShell C18 2.6 μ m, 75 x 4.6 mm
 Mobile phase:
 A) 0.1% Phosphoric acid
 B) CH_3CN
 Gradient program

Flow rate: 1.0 mL/min,
 Temperature: 25 $^{\circ}$ C
 Detection: UV@250 nm
 Sample: Oolong tea

SunShell C18-WP, RP-AQUA, C8, Phenyl, PFP, PFP&C18, Cyano, 2.6 μm

(Pentafluorophenyl)


◆ Characteristics of SunShell

	Core shell silica			Bonding phase					
	Particle size	Pore diameter	Specific surface area	Carbon content	Bonded phase	USP L code	End-capping	Maximum operating pressure ^a	Available pH range
SunShell C18	2.6 μm	9 nm	150 m ² /g	7%	C18	L1	Sunniest endcapping	60 MPa	1.5 - 10
SunShell C18-WP	2.6 μm	16 nm	90 m ² /g	5%	C18	L1	Sunniest endcapping	60 MPa	1.5 - 10
SunShell RP-AQUA	2.6 μm	16 nm	90 m ² /g	4%	C30	L62	Sunniest endcapping	60 MPa	2 - 8 ^b
SunShell C8	2.6 μm	9 nm	150 m ² /g	4.5%	C8	L7	Sunniest endcapping	60 MPa	1.5 - 9
SunShell Phenyl	2.6 μm	9 nm	150 m ² /g	5%	Phenylhexyl	L11	Sunniest endcapping	60 MPa	1.5 - 9
SunShell Biphenyl	2.6 μm	9 nm	150 m ² /g	5%	Biphenyl	L11	Sunniest endcapping	60 MPa	1.5 - 9
SunShell PFP	2.6 μm	9 nm	150 m ² /g	4.5%	Pentafluorophenyl	L43	TMS endcapping	60 MPa	2 - 8
SunShell PFP&C18	2.6 μm	9 nm	150 m ² /g	6%	Pentafluorophenyl + C18	L43	TMS endcapping	60 MPa	2 - 8
SunShell Cyano	2.6 μm	9 nm	150 m ² /g	2.5%	Diisopropylcyanopropyl	L10	No	60 MPa	2 - 8



NEW

a) Unless otherwise specified in the column test report b) Under 100% aqueous condition

◆ Separation of standard samples






Separation of basic compounds




Separation of fluorobenzene

Separation of cresol isomers


Separation of nucleotides




Separation of amino acids (LC/MS)



Separation of methylhippuric acid isomers




Separation of cis, trans-stilbene




Column: SunShell C18 2.6 μ m, 150 x 4.6 mm i.d.
 SunShell C30 2.6 μ m, 150 x 2,1 mm i.d.
 SunShell Phenyl 2.6 μ m, 150 x 4.6 mm i.d.
 SunShell Biphenyl 2.6 μ m, 150 x 4.6 mm i.d.
 SunShell PFP&C18 2.6 μ m, 150 x 4.6 mm i.d.
 SunShell PFP 2.6 μ m, 150 x 4.6 mm i.d.


Mobile phase: Methanol/water = 75/25
 Acetonitrile/water = 60/40
 Flow rate: 1.0 mL/min and 0.2 mL/min for only C30
 Temperature: 40 °C
 Detection: UV@230 nm
 Sample: 1 = cis-Stilbene, 2 = trans-Stilbene




Separation of branched-chain amino acids



Column: SunShell PFP, PFP&C18, RP-AQUA, Biphenyl 2.6 μ m, 150 x 4.6 mm
 Mobile phase: 0.1% formic acid
 Flow rate: 1.0 mL/min
 Temperature: 40 °C
 Detection: UV@205 nm
 Sample: 1 = L-Valine 2 = L-Isoleucine 3 = L-Leucine




Separation of steroids




Column: SunShell Biphenyl, Phenyl and C18 2.6 μ m, 150 x 4.6 mm
 Mobile phase: Acetonitrile/water = 45/55
 Flow rate: 1 mL/min
 Temperature: 25 °C
 Detection: UV@230 nm

- Peak
1. Estriol
 2. Hydrocortisone
 3. Prednisone
 4. Corticosterone
 5. β -Estradiol
 6. Testosterone
 7. Cortisonacetate



SunShell C30, 2.6 μm




Specification of SunShell C30


	Core shell silica				Bonding phase					
	Particle size (μm)	Core size (μm)	Pore size (nm)	Specific surface area (m^2/g)	Carbon loading (%)	Ligand	USP L category	End-capping	Maximum pressure ^{a)}	pH range
SunShell C30	2.6	1.6	12	95	7	C30	L62	TMS	60 MPa	1.5 - 9

a) Unless otherwise specified in the column test report


Problem of C30 column




Separation of tocopherols



Fast separation of vitamin K1 isomers



Comparison of isomers separation of Vitamin k1



SunShell 2.6 μm C18-WP, HFC18-16, C8-30HT, C4-100






Characteristics of SunShell

For separation of peptides and proteins


	Core shell silica			Bonding phase						
	Particle size	Pore diameter	Specific surface area	Stationary phase	Carbon content	Ligand density	End-capping	Maximum operating pressure ^a	Available pH range	USP L line
SunShell C18-WP	2.6 μm	16 nm	90 m ² /g	C18	5 %	2.5 μmol/m ²	Sunniest endcapping	60 MPa or 8,570 psi	1.5 - 10	L1
SunShell HFC18-16	2.6 μm	16 nm	90 m ² /g	C18	2.5%	1.2 μmol/m ²	Sunniest endcapping	60 MPa or 8,570 psi	1.5 – 9	L1
SunShell C8-30HT	3.4 μm	30 nm	15 m ² /g	C8	0.5%	2.5 μmol/m ²	Sunniest endcapping	60 MPa ^a or 8,570 psi ^a	1.5 – 9	L7
SunShell C4-100	2.6 μm	100 nm	22 m ² /g	C4	0.6%	3 μmol/m ²	Sunniest endcapping	60 MPa ^a or 8,570 psi ^a	1.5 – 8	L26

Note: The SunShell HFC18-30, C8-30, and C4-30 columns will be discontinued once the packaging materials are out of stock.

a) Unless otherwise specified in the column test report
b) 50MPa, 7141psi for 4.6 mm i.d. column




Separation of peptides



SunShell 2.6 μm C8-30HT, C4-100

For separation of peptides and proteins

Comparison of column temperature



Column: SunShell C8-30, 2.6 μm (30 nm) 100 x 2.1 mm
 Mobile phase: A) 0.1% TFA in water
 B) 0.08 % TFA in acetonitrile

Gradient program: Time 0 min 15 min
 %B 20% 65%

Flow rate: 0.5 mL/min,
 Temperature: 25 °C 60 °C or 80 °C


Detection: UV@215 nm,

Sample: 1 = Cytochrome C, 2 = Lysozyme, 3 = BSA,
 4 = Myoglobin, 5 = Ovalbumin

A macromolecule compound like a protein diffuses very slowly, so that an elevated temperature makes a peak be shaper and improves separation. BSA peak seemed to be tailing at 25 degree Celsius. BSA, however, was separated several peaks at 80 degree Celsius.



Comparison of SunShell stationary phase



Column dimension: 100 x 2.1 mm,
 Mobile phase: A) 0.1% TFA in water, B) 0.1 % TFA in Acetonitrile
 Gradient program: Time 0 min 60 min
 %B 20% 65%



Flow rate: 0.5 mL/min, Temperature: 80 °C, Detection: UV@215 nm, Injection volume: 1.0 μL
 Sample: 1 = Cytochrome C, 2 = Lysozyme, 3 = BSA, 4 = Myoglobin, 5 = Ovalbumin,
 UHPLC instrument: HITACHI Chromaster

Comparison of peak width (W0.5, min)

	C4-100	C4-30	C8-100	C8-30	C8-30HT	C8-30HT 0.5 μL	Sample concentration
Cytochrome C	0.167	0.177	0.160	0.155	0.212	0.144	0.050%
Lysozyme	0.164	0.180	0.153	0.166	0.196	0.145	0.050%
BSA	0.308	0.410	0.276	0.514	0.422	0.330	0.100%
Myoglobin	0.197	0.221	0.180	0.199	0.238	0.176	0.050%
Ovalbmin	0.391	0.889	0.247	0.428	0.184	0.176	0.050%

The above table indicated that C4-100 with 1000Å of pore showed a sharper peak than the other. C8-30HT has a thin porous layer and low surface area, so that low sample loadnig made a peak sharper.

Separation of monoclonal antibody



Column dimension: 100 x 2.1 mm i.d.
 Mobile phase: A) water/TFA (100/0.1)
 B) acetonitrile/TFA (100/0.085)
 30-45%B (0-30 min)

Flow rate: 0.4 mL/min
 Temperature: 80 °C
 Detection: UV at 215 nm
 Injection: 5 μL

Sample : Monoclonal antibody purified from cell culture using Protein G

Not eluted Core shell C8, 300 Å

Totally porous C4, 300 Å

Core shell C4, 300 Å

SunShell C4-100, 1000 Å

Separation using nanocolumn SunShell C4-100 (100 x 0.15 mm i.d.)
 Flow rate: 4 μL/min

Retention time/min

Regarding reversed phase separation of monoclonal antibody (IgG), not only core shell C4 with 30 nm pore showed the better separation than totally porous C4, but also 100 nm of pore leded the best separation. Nano column showed almost the same separation of IgG as semi-micro column.

SunShell HILIC-Amide, HILIC-S, 2.6 μm


For Hydrophilic Interaction Chromatography

Characteristics of SunShell HILIC-Amide

	Core shell silica				Bonded phase					
	Particle size	Core size	Pore diameter	Specific surface area	Carbon content	Bonded phase	End-capping	Maximum operating pressure ^a	USP category	Available pH range
SunShell HILIC-Amide	2.6 μm	1.6 μm	9 nm	150 m ² /g	3%	Amide	No	60 MPa or 8,570 psi	L68	2 - 8
SunShell HILIC-S	2.6 μm	1.6 μm	9 nm	150 m ² /g	0%	Bare silica	No	60 MPa or 8,570 psi	L3	1 - 5


a) Unless otherwise specified in the column test report

Stationary phase of HILIC-Amide




R: Hydrophilic group

Stationary phase of HILIC-S




Bare silica




Stationary phase of SunShell HILIC-Amide consists of AMIDE and HYDROPHILIC GROUP, so that this stationary phase is more polar than an individual group. High speed separation is led by core shell structure that derives high efficiency and fast equilibration. HILIC-S is recommended for separation using LC/MS.


Separation of Nucleic acid bases: Comparison of the other core shell hilic columns





Separation of Cyanuric acid and Melamine




Separation of water- soluble vitamins




Artificial sweeteners



1. Aspartame





2. Saccharin




3. Acesulfame K

Column: SunShell HILIC-Amide, 2.6 μ m, 100 x 4.6 mm
 Mobile phase: Acetonitrile: 25 mM phosphate buffer (pH2.5) =8:2
 Flow rate: 1.0 mL/min ,
 Temperature: Ambient
 Detection: UV@215 nm
 Sample: 1 = Aspartame, 2 = Saccharin, 3 = Acesulfame K


Glycoside




1. Helicin



2. Salicin




3. Arbutin



4. Rutin


Column: SunShell HILIC-Amide, 2.6 μ m, 100 x 4.6 mm
 Mobile phase: Acetonitrile:25 mM phosphate Ammonium (pH4.9) =8:2
 Flow rate: 1.0 mL/min
 Temperature: Ambient
 Detection: UV@215 nm
 Sample: 1 = Helicin, 2 = Salicin, 3= Arbutin, 4 = Rutin

Nucleic acid base



Column: SunShell HILIC-S, 2.6 μ m 100 x 2.1 mm
 Mobile phase: 100 mM ammonium acetate (pH3.0) /acetonitrile = 1/9
 Flow rate: 0.4 mL/min
 Temperature: 40 $^{\circ}$ C
 Detection: UV@250 nm
 Sample: 1 = Acenaphthene, 2 = Uridine, 3 = Cytosine

Nucleic acid bases



Column: SunShell HILIC-S, 2.6 μ m 100 x 2.1 mm
 Mobile phase: 100 mM ammonium acetate (pH3.0) /acetonitrile = 1/9
 Flow rate: 0.2 mL/min
 Temperature: 40 $^{\circ}$ C
 Detection: UV@250 nm
 Sample: 1 = Thymine, 2 = Uracil, 3 = Uridine, 4 = Cytosine, 5 = Cytidine





Metal Free Column

PEEK clad SUS column

Metal-free column

- Standard metal-free columns with high chemical stability.
- Column hardware has a pressure resistance of 100 MPa.
- Metal-free columns supplied by two companies.
- Can be packed with all SunShell packing materials.



High pressure-resistant PEEK columns manufactured by Tomoe

*I.D. 2.1 mm, length 50 mm, 100 mm, 150 mm

*The catalog number, change the last digit "1" to "MTF".

For example, for SunShell C18 2.6 μm , 100 x 2.1 mm columns, the model number CB6961 for the standard columns becomes CB696MTF.



IDEX Bioinert System Column

*2.1 and 4.6 mm i.d., lengths 50 mm, 100 mm and 150 mm.

*The catalog number, change the last digit "1" to "M".

For example, for a SunShell C18 2.6 μm , 150 x 4.6 mm column, the catalog number CB6371 of the normal column becomes CB637M.



Nanocolumn, Microcolumn

Nano column: 0.075 mm i.d., 0.1 mm i.d., 0.15 mm i.d.


Micro column: 0.3 mm i.d., 0.5 mm i.d.

Column length: 50 mm, 150 mm

Column tube material:

0.075mm i.d., 0.10 mm i.d., 0.15 mm i.d.; PEEKSIL

0.3 mm i.d., 0.5 mm i.d.; Glass Lined SUS Tubing



Example of prices

	Inner diameter (mm)	0.075	0.1	0.15	0.3	0.5	USP L code
	Length (mm)	Catalog number	Catalog number	Catalog number	Catalog number	Catalog number	
SunShell C18, 2 μm	50	CB1J4P	CB1H4P	CB1K4P	CB1G48	CB1F48	L1
	150	CB1J7P	CB1H7P	CB1K7P	CB1G78	CB1F78	
SunShell C18, 2.6 μm	50	CB6J4P	CB6H4P	CB6K4P	CB6G48	CB6F48	
	150	CB6J7P	CB6H7P	CB6K7P	CB6G78	CB6F78	
SunShell C18, 5 μm	50	-----	-----	-----	CB3G4L	CB3F4L	
	150	-----	-----	-----	CB3G7L	CB3F7L	
SunShell Phenyl, 2.6 μm	50	CP6J4P	CP6H4P	CP6K4P	CP6G48	CP6F48	L11
	150	CP6J7P	CP6H7P	CP6K7P	CP6G78	CP6F78	
SunShell C8-30HT, 3.4 μm	50	C56J4P	C56H4P	C56K4P	C56G48	C56F48	L7
	150	C56J7P	C56H7P	C56K7P	C56G78	C56F78	
SunShell C4-100, 2.6 μm	50	C66J4P	C66H4P	C66K4P	C66G48	C66F48	L26
	150	C66J7P	C66H7P	C66K7P	C66G78	C66F78	

※ Packings (stationary phase) and column sizes other than those listed above can also be manufactured. For details, please contact ChromaNik Technologies.

※ The end-fitting of the column is Parker type.

※ The P, L or 8 at the end of the catalog number indicates the material of the column tubing: P is PEEKSIL, L and 8 are glass-lined tubing. Only L has an upper limit of back pressure of 45 MPa, but the others have an upper limit of back pressure of 80 MPa.

SunShell Guard Cartridge Column



RP & S GUARD CARTRIDGE COLUMN




- * The cartridge column is packed with SunShell C18 (RP) and Core shell silica (S) into a cartridge sized 3 x 2 mm i.d.
- * RP guard cartridge is used for all reversed phases and S guard cartridge for hilic phases.
- * Low dead volume structure
- * Upper pressure limit is more than 60 Mpa
- * Available for 2.1 mm i.d. to 4.6 mm i.d. columns

SunShell C18, 2.6 μ m 50 x 2.1 mm

Without Guard Cartridge column

Back pressure: 10.2 MPa



Mobile phase:

CH₃CN/H₂O=60/40 for 2.1 mm i.d.

CH₃CN/H₂O=70/30 for 4.6 mm i.d.

Flow rate:

0.3 mL/min for 2.1 mm i.d.

1.8 mL/min for 4.6 mm i.d.

Temperature: 25 °C

Detection: UV@250nm

Sample: 1 = Uracil


2 = Ethylbenzoate

3 = Acenaphthene

4 = Butylbenzene

With Guard Cartridge Column RP

Back pressure: 10.5 MPa




SunShell C18, 2.6 μ m 150 x 4.6 mm

Without Guard Cartridge column

Back pressure: 28.8 MPa


$t_R(3) = 3.20$ min
 $N(3) = 35,632$



With Guard Cartridge Column RP

Back pressure: 31.4 MPa

$t_R(3) = 3.21$ min
 $N(3) = 35,204$
1% decrease of plate



Ordering Information of SunShell Guard Cartridge Column

Description	Part number
SunShell Guard Cartridge RP Starter Kit (holder, cartridge, tubing)	CB32CK
SunShell Guard Cartridge RP for exchange (2 PCS)	CB32CC
SunShell Guard Cartridge S Starter Kit (holder, cartridge, tubing)	CS32CK
SunShell Guard Cartridge S for exchange (2 PCS)	CS32CC
SunShell Guard Cartridge holder	HOL2CC

Ordering information of SunShell

	Inner diameter (mm)	1.0	2.1	3.0	4.6	USP category
	Length (mm)	Catalog number	Catalog number	Catalog number	Catalog number	
SunShell C18, 2 µm	50	-----	CB1941	-----	-----	L1
	100	-----	CB1961	-----	-----	
	150	-----	CB1971	-----	-----	
SunShell C18, 2.6 µm	30	-----	CB6931	CB6331	CB6431	
	50	CB6141	CB6941	CB6341	CB6441	
	75	-----	CB6951	CB6351	CB6451	
	100	CB6161	CB6961	CB6361	CB6461	
	150	CB6171	CB6971	CB6371	CB6471	
	250	-----	-----	CB6381	CB6481	
SunShell C18 3.5 µm	50	-----	CB9941	-----	-----	
	100	-----	CB9961	CB9361	CB9461	
	150	-----	CB9971	CB9371	CB9471	
	250	-----	-----	CB9381	CB9481	
SunShell C18, 5 µm	150	-----	-----	CB3371	CB3471	
	250	-----	-----	CB3381	CB3481	
SunShell C8, 2.6 µm	30	-----	CC6931	CC6331	CC6431	L7
	50	-----	CC6941	CC6341	CC6441	
	75	-----	CC6951	CC6351	CC6451	
	100	-----	CC6961	CC6361	CC6461	
	150	-----	CC6971	CC6371	CC6471	
SunShell PFP, 2.6 µm	30	-----	CF6931	CF6331	CF6431	L43
	50	-----	CF6941	CF6341	CF6441	
	75	-----	CF6951	CF6351	CF6451	
	100	-----	CF6961	CF6361	CF6461	
SunShell C18-WP, 2.6 µm	30	-----	CW6931	CW6331	CW6431	L1
	50	-----	CW6941	CW6341	CW6441	
	75	-----	CW6951	CW6351	CW6451	
	100	-----	CW6961	CW6361	CW6461	
	150	-----	CW6971	CW6371	CW6471	
SunShell RP-AQUA, 2.6 µm	30	-----	CR6931	CR6331	CR6431	L62
	50	CR6141	CR6941	CR6341	CR6441	
	75	-----	CR6951	CR6351	CR6451	
	100	CR6161	CR6961	CR6361	CR6461	
	150	CR6171	CR6971	CR6371	CR6471	
SunShell Phenyl, 2.6 µm	30	-----	CP6931	CP6331	CP6431	L11
	50	-----	CP6941	CP6341	CP6441	
	75	-----	CP6951	CP6351	CP6451	
	100	-----	CP6961	CP6361	CP6461	
	150	-----	CP6971	CP6371	CP6471	
SunShell Biphenyl, 2.6 µm	30	-----	C86931	C86331	C86431	L11
	50	-----	C86941	C86341	C86441	
	75	-----	C86951	C86351	C86451	
	100	-----	C86961	C86361	C86461	
	150	-----	C86971	C86371	C86471	
SunShell C30, 2.6 µm	30	-----	CT6931	CT6331	-----	L62
	50	-----	CT6941	CT6341	CT6441	
	75	-----	CT6951	CT6351	-----	
	100	-----	CT6961	CT6361	CT6461	
	150	-----	CT6971	CT6371	CT6471	
SunShell PFP&C18, 2.6 µm	30	-----	CV6931	CV6331	CV6431	L43
	50	-----	CV6941	CV6341	CV6441	
	75	-----	CV6951	CV6351	CV6451	
	100	-----	CV6961	CV6361	CV6461	
	150	-----	CV6971	CV6371	CV6471	
SunShell Cyano, 2.6 µm	50	-----	CJ6941	CJ6341	CJ6441	L10
	100	-----	CJ6961	CJ6361	CJ6461	
	150	-----	CJ6971	CJ6371	CJ6471	

Added 4.6 mm i.d.
for C30 phase

