



IRMS-GC-GC5 USER MANUAL



VERSION HISTORY

Version 1.00 – Release Document

Version 1.01 – Grammar corrections. New diagrams and photos added. Changes made based on first review.

Version 1.02 – 2nd Review Changes

FIRST RELEASE

Version 1.10 – H2, N2 and O2 Upgrades section added

Version 1.11 – Some corrections and additions made to manual – some pictures updated, further spelling grammar and formatting alterations

Version 1.12 – Cross references added so page hyperlinks are automatically updated if page number changes

Version 1.13 – Rechecked whole document for consistency.

SECOND RELEASE

CONTENTS

Version History	1
Introduction	4
Warranty	5
Safety information	6
Safety symbols	7
Site requirements	8
System Information	10
Technical specifications	12
GC5 Combustion/Pyrolysis Interface System.....	13
Overview	13
The GC5 Interface	14
The Water Trapping Systems	18
The Gas Chromatograph - Agilent 7890.....	19
The Gas Chromatograph Cabinet.....	19
The Heart Split Mechanism.....	19
Connecting the GC	20
Gas Connections	21
Installing the furnace tubes	21
Column Set Up	24
Leak Checking	27
GC Configuration.....	28
GC5 Furnace Temperature.....	28
Comms	28
Setup IonVantage.....	30
Setup Agilent GC	30
Conditioning.....	32
Running an Analysis	33
Setting Up Ionvantage	33

Saving the GC Method files.....	34
Breakthrough	37
Setting up Heart Split Timings.....	38
Peforming a Carbon Analysis	40
Reoxidising the Carbon/Nitrogen Furnace Tube.....	40
Procedure	41
GC5 Optional Upgrades	42
Nitrogen Analysis	43
Set up for N2 Analysis	43
Performing a Nitrogen Analysis	45
Hydrogen Analysis	47
Setting Up For Hydrogen Analysis.....	47
Perfroming a Hydrogen Analysis.....	47
Oxygen Analysis	49
Setting up for Oxygen Analysis	49
Preparation of the GC O tube	49
GC Oven connections.....	50
Furnace GC connections	50
Conditioning the Furnace.....	51
Checking background levels.....	51
Setting Up the GC.....	53
GC Inlet Method.....	53
Pre-analysis	53
Coating The Tube	54
Performing An Oxygen Analysis	55
Troubleshooting.....	56

INTRODUCTION

This manual provides the operation and maintenance instructions for the Isoprime GC5. This will include details on the Isoprime mass spectrometer with the GC and GC5 interface. We have endeavoured to include all operational aspects within this manual, however if you feel more information is required or that you can add some more details, then please contact the Customer Service Department at Isoprime Ltd or your local Isoprime representative.

WARRANTY

The information contained within this document is subject to change without notice.

Isoprime Ltd. makes no warranty of any kind with regard to this material, including, but not limited to, the implied warranties of merchantability and fitness for a particular purpose.

Isoprime Ltd. shall not be liable for errors contained herein or for incidental damages in connection with the furnishing, performance, or use of this material.

SAFETY INFORMATION

The Isoprime mass spectrometer meets the following IEC classifications.

Safety	Class 1
Transient Overvoltage	Category II
Pollution	Degree 2

This unit has been designed and tested in accordance with recognised safety standards and designed for use indoors. If the instrument is used in a manner not specified by the manufacturer, the protection provided by the instrument may be impaired. Whenever the safety protection of the Isoprime has been compromised, disconnect the instrument from all power sources and secure the unit against unintended operation.

Suitably qualified personnel only should perform maintenance procedures. Substituting parts or performing any unauthorised modification to the instrument may result in a safety hazard.

Disconnect the mains supply before removing covers. You must comply with all local and national requirements for electrical and mechanical safety.

Please contact the Customer Service Department at:

Isoprime Ltd
Isoprime House
Earl Road
Cheadle Hulme
SK8 6PT
UK
service@isoprime.co.uk

or your local Isoprime representative should you require any further information.

SAFETY SYMBOLS

Warnings in the manual or on the instrument must be observed during all phases of service, repair, installation and operation of this instrument. Failure to comply with these precautions violates the safety standards of design and the intended use of the instrument.

Isoprime Ltd. assumes no liability for the customer's failure to comply with these requirements.

The following safety notices and symbols are used in the manual or on the instrument.

WARNING

Warnings are given to highlight situations or conditions where failure to observe the instruction could result in injury or death to persons.

CAUTION

Cautions are given to highlight situations or conditions where failure to observe the instruction could result in damage to the equipment, associated equipment or process.



This symbol indicates that there are accompanying instructions which should be referred to for more information.



This symbol indicates where hazardous voltages may be present.



This symbol indicates a hot surface.



This symbol indicates a risk of corrosion.

SITE REQUIREMENTS

BENCH SPACE

Must be of a construction adequate to support the weight of the instrument and associated systems, should be heat resistant and should also be free of any vibration or movement which may impair the performance of the instrument.

ROOM TEMPERATURE

The maximum operational room temperature should not exceed 25°C (ideally 22°C). Large temperature fluctuations caused by sunlight or drafts (including air conditioning units) have to be avoided during measurement, and should not exceed 1°C/ hour.

HUMIDITY

Relative humidity should not exceed 60% at any temperature.

CORROSIVE AND HARMFUL ATMOSPHERE

Instrument must be located in a clean air, pollution free environment.

DUST FREE ENVIRONMENT

High dust levels increase the probability of computer disc drive failures. Dust particles will also block air filters and fans on the instrument which may cause overheating.

EXHAUST FUMES

The vent ports of the individual systems should be connected to a fume extraction system.

CAUTION: If SO₂ or CO is used in your chosen system, an efficient fume extraction hood is needed so that concentrations in laboratory air comply with the local health and safety standards.

COMPRESSED AIRLINE

A water and oil mist filter should be fitted to the supply line and contain at least a 0.5µm particulate filter. Addition of a 0.3µm mist filter is recommended to minimise the amount of oil transferred to the instrument.

GAS REGULATORS

Clean, two stage regulators with stainless steel diaphragms must be fitted on the gas supplies. The second stage regulator must be capable of delivering pressures up to 6 Bar (90 psi).

GAS LINES

Clean, chromatography grade 1/8" stainless steel gas lines are needed to connect the system to the gas cylinders.

Note: Specialised gases for use with the instruments are not provided.

Note: Chromatography grade pipework is not provided.

ELECTRICAL POWER REQUIREMENTS

The system may require one or more connections to the mains supply. Some electrical systems have circuits which are susceptible to large voltage transients. These should be protected by a suitable transient suppressor if the local supply is of poor quality. Some sites may also wish to protect some (the pumping circuits, for example) or all circuits from power failure with an un-interruptible power supply unit. Voltages quoted should be adhered to, e.g. when 230V is used, this voltage can be anywhere between 230V and 250V, 208V is not acceptable. For sites where the voltage does not comply with the quoted specifications, a suitable transformer will be required.

SYSTEM INFORMATION

SYSTEM TYPE

GC Pyrolysis/Combustion interface with 7890A GC

PLACEMENT

This is a bench top instrument

Provide at least 250mm clearance around the system for maintenance access.

GAS CONNECTIONS

For the connection of the gases, a supply of 1/8" GC-Grade stainless steel tubes needs to be available, to be supplied by the customer. 1/8" Swagelok connections should be used on the cylinder and on the system where relevant.

DIMENSIONS

	GC5 Interface	Agilent 7890A GC	Agilent 7693 Sampler
Height	540mm	590mm	450mm
Width	520mm	540mm	20mm
Depth	400mm	500mm	500mm
Weight	80kg	50kg	6kg

ELECTRICAL POWER

	GC5 Interface	Agilent 7890A GC
Voltage	230v	120V – 240V
Phase	50 - 60Hz	48 – 66Hz
Max Current	9A	18.8 – 9.4A
Max. Fuse Rating	16A	20A
Transformer Rating	2.1kVA	2.3kVA
Transient Sensitive	No	No
Connection Type	10A IEC320/C13	16A IEC320/C13

CONNECTIONS NUMBER OF ELECTRICAL SOCKETS REQUIRED

Max. 4

All connections to the mains supply should have a terminal (the neutral) which must be at earth (ground) potential. Cables will be delivered with a standard 3 pin UK plug and/or and IEC power block. Alternative cables or plug adaptors may be required.

HEAT DISSIPATION

GC5 Interface	Agilent 7890A GC
300W	2250W

COMPRESSED AIRLINES

The GC5 Interface requires a supply of 60 psi / 4 Bar.

The GC5 Interface is equipped with a 6mm push-fit connector

SPECIALISED GASSES

	Purity grade	H ₂ O Content
Helium	99.999%	< 1 vpm
Synthetic air*	Zero Grade	2 vpm
Hydrogen*#	99.999%	< 2 vpm
Oxygen	99.998%	1 vpm

* Required for the operation of the FID Detector.

Required for Oxygen Analysis

LIQUID NITROGEN

The system requires liquid nitrogen for nitrogen analysis. If the Agilent 7890A GC is equipped with a cryogenic cooled PTV injector, either liquid nitrogen (1/4" insulated copper tube) or liquid carbon dioxide (1/8" heavy walled, stainless tubing) is needed.

TECHNICAL SPECIFICATIONS

All GC installation and service work is tested using Universal mix that comprises of four compounds; two Alkanes, a Methyl Ester and Nicotine and it is this mix that is referred to in this manual.

EXTERNAL PRECISION ^{13}C (CO_2)

$\leq 0.2\text{‰}$ (SD 1σ ; on each component over 5 injections of Universal mix)

EXTERNAL PRECISION ^{15}N (N_2)

$\leq 0.5\text{‰}$ (SD 1σ ; on a single component over 5 injections of Universal mix)

EXTERNAL PRECISION D (H_2)

$\leq 3.0\text{‰}$ (SD 1σ ; on each component over 5 injections of Universal mix)

EXTERNAL PRECISION ^{18}O (CO)

$\leq 0.8\text{‰}$ (SD 1σ ; on a single component over 5 injections of Universal mix)

GC5 COMBUSTION/PYROLYSIS INTERFACE SYSTEM

OVERVIEW



Right hand GC with autosampler

The Isotope Ratio Mass spectrometer (IRMS) GC-GC5 Interface System is designed to determine the composition of individual compounds in complex organic mixtures.

The sample is injected into the GC where it is vaporised and then separated using a capillary gas column. The sample travels through the column in a stream of inert carrier gas, in this case helium, where, depending on the properties of each component, it is separated. The time it takes for the sample to elute is known as the retention time. The retention time and quality of separation is determined by a number of factors, such as column and sample properties and oven temperatures.

The user can select the separated portions of the sample to be analysed and the chosen portion can then be sent to the GC5. Using copper oxide (CuO) for carbon/nitrogen, the sample is oxidised/reduced, allowing for isotopic analysis in the IRMS.

The complete GC-GC5 Interface consists of the following parts:

1. An Agilent 7890A GC modified for IRMS coupling
2. The combustion/pyrolysis GC5 Interface, including regeneration possibilities for the combustion/pyrolysis tubes
3. Water traps for the protection of the IRMS
4. An optional standby valve and/or autosampler

THE GC5 INTERFACE

GC5 has 4 main components:

1. Gas Distribution
2. Electrical System
3. GC5 Combustion / Pyrolysis Interface
4. Open Split

GAS DISTRIBUTION

The following gases are required for the full operation of the interface:

- Helium (He) - carrier gas for the GC and the furnace
- Hydrogen (H₂) - for oxygen analysis
- Oxygen (O₂) - for the regeneration of the furnace tube
- Compressed Air - for valve operation
- Clean Air - for operation of the FID


The gas distribution network for both the standard carbon/nitrogen/hydrogen and additional oxygen option is positioned in the top section of the GC interface cabinet and is depicted below.

The compartment is accessed by removing the top panel of the interface cabinet. The gases are introduced to the rear of the cabinet.


On the front of the cabinet are two 15 psi pressure regulators, one for the FID helium supply and the second for the sample line helium. Between the regulators is a manual switching valve for selecting sample line helium to the right (for GC carbon/nitrogen/hydrogen analysis) or left (oxygen analysis).

The FID helium supply runs from the He input, at the back of the GC5 through the FID He regulator and into the GC cabinet. In the GC cabinet the tubing connects to one side of the T section of the SGE MOVPT 1/100 Heart Split valve.

The sample line helium flows from the He input to the sample line regulator. During instrument operation, helium flows from the exit of the regulator, into the GC cabinet and connects to the zero volume T-piece on the end of the furnace tube. The helium flows along the length of the furnace tube and is either taken into the capillaries running to the IRMS or vented to atmosphere at the open split "pigtail" assembly fixed to the end of the furnace tube.



Standard gas line network



With additional oxygen option

THE ELECTRICAL SYSTEM

The electrical system consists of two heaters: one for the heated interface and one for the furnace. These can be individually enabled via switches beside the units.

COMBUSTION INTERFACE

The combustion interface consists of two different elements:

1. An interface heater running between the GC oven and the furnace assembly
2. Combustion assembly consisting of a furnace assembly and a furnace tube running through the middle of the interface heater and furnace assembly.

INTERFACE HEATER

The heated interface ensures that sample eluting from the column is transported to the combustion furnace with no risk of condensation.

The temperature of the heater must be sufficiently high to avoid condensation of the sample as it passes from the GC to the GC5 furnace assembly.

Cold spots on the interface must be avoided. Thus the interface assembly must be properly engaged within the wall of the GC cabinet and in the insulation disk at the end of the furnace assembly. This will ensure minimum disruption of the temperature profile between the GC cabinet and the interface heater and the furnace assembly.

COMBUSTION ASSEMBLY

The **Combustion Assembly** consists of a modified Carbolite MTF12/10 furnace controlled by a CAL3300 control unit situated on the front of the GC5 unit.

The **Furnace Tube** is the site of combustion or pyrolysis, converting the sample compounds to the analyte gas for analysis. The furnace tube runs through the centre of the interface heater and furnace assembly. The geometry of the tube filling is designed to ensure that the packing lies within the hottest section of the furnace assembly.

It is important that the isotopic composition of the sample gas accurately reflects that of the original starting compound. Combustion/pyrolysis should therefore be complete. It is important that the tube filling (reactants or catalysts) is at a sufficiently high temperature to achieve full combustion/pyrolysis, nor should there be any cold spots in the combustion interface or sites of contamination, in which the sample compounds can condense.

THE OPEN SPLIT

The open split “pigtail” assembly is designed to protect the ion source from pressure fluctuations arising from the processes occurring within the GC or GC5. At the open split the pressure is maintained by the ambient atmospheric pressure. The open split marks the change from a high pressure system in the GC and combustion interface, to a low pressure system dominated by the IRMS vacuum system and extending to the water trap facility and the connecting capillaries.

Sample line helium flows from the Zero Dead-Volume T-Piece at the GC end of the furnace tube towards the open split. At the open split the fused silica capillary connected to the water trap and ultimately the IRMS draws up helium using the IRMS vacuum system. Excess helium is vented to atmosphere through the stainless steel capillary on the open split “pigtail” union.

The sample line helium pressure is set to prevent ingress of atmospheric gas at the open split assembly. Operating the system at a high sample line helium pressure results in dilution of the sample, reducing the sensitivity and too low a sample line helium pressure will result in “breakthrough” of the solvent into the GC5 and IRMS (see page [37](#) for further information)

THE WATER TRAPPING SYSTEMS

The water trapping facilities available on the GC interface are designed to remove water generated in the furnace tube during the combustion of the organic molecules. Two types of water trapping facility are available. The water traps are positioned after the open split assembly of the combustion interface.

NAFION MEMBRANE

During combustion, water can form in the sample. It can be removed by passing the sample through a nafion membrane. The polymer in the nafion membrane is strongly hydrophilic, meaning that it strongly attracts moisture.

Water from the sample gas is absorbed on the inner surface of the membrane and passes towards the outer surface where it is then removed by a dry flow of helium passing over the outside of the nafion membrane.

CRYOGENIC TRAP FOR NITROGEN MEASUREMENTS

The second water trapping facility is a cryogenic trap or “cold” trap, which freezes both H₂O and CO (arising from CO₂) to allow for accurate N₂ measurements

The Dewar used when measuring nitrogen consists of an insulated container or “trap” filled with liquid nitrogen (temperature of -196°C) into which a loop of 320um deactivated silica capillary is suspended. The system has no electronic temperature control and relies entirely on the liquid nitrogen temperature to capture both the H₂O and CO₂.



Nitrogen Cold Trap - Typical usage of liquid nitrogen is approximately 5 litres per day.

THE GAS CHROMATOGRAPH - AGILENT 7890

An Agilent 7890A Gas Chromatograph User Manual is either supplied with the instrument or available to download online.


THE GAS CHROMATOGRAPH CABINET

The GC in its standard configuration consists of a standard GC oven fitted with a Split/Splitless injector system and a Flame Ionisation Detector (FID). The GC has been modified in order to accommodate a number of components necessary to interface it with the IRMS. It can still nevertheless be used as a stand-alone unit.

In the GC cabinet, a system of capillaries and valves enables diversion of the sample either to the FID or to the combustion interface. The arrangement of capillaries within the GC cabinet is shown in the figure below.

THE HEART SPLIT MECHANISM

A schematic of the GC oven arrangement is shown below.



Heart Split Mechanism. Letters indicate Ferrule type (See table below)

Location	Ferrule Type	Info
A	GVF 16-003/004	Vespel Backdrilled
B	GVF 16-005	Vespel
C	GVF 16-16	Vespel Sealing Ring
D	GFF 16-005	Graphite
E	GVF 16 (2) 004	Vespel Two Hole
F	GVF 16 4-4	Vespel

The Heart Split mechanism operates via a combination of a mechanical input (Heart Split valve) and a back flow of pressure.

Helium enters the GC column at the injection port and travels through the column. Additionally, helium enters the furnace tube at the VSOS via the sample line helium input and also into the Heart Split valve via the FID helium input.

When the Heart Split valve is open, sample flow is prevented from passing towards the furnace tube by the pressure of the sample line helium. The flow therefore passes towards the Heart Split valve and hence to the FID.

When the Heart Split valve is closed, the flow is mechanically prevented from passing into the Heart Split valve. The pressure in the VSOS increases to a point where it is greater than the sample line helium pressure. The flow now passes into the furnace tube.

The Heart Split mechanism has three purposes:

- 1) It enables the solvent to be diverted away from the combustion furnace and the IRMS and vented to atmosphere. It is very unwise to allow the solvent peak to pass through the combustion furnace. Firstly, it depletes the furnace of oxygen and secondly, due to the usually large quantity of solvent injected, in particular when the GC is operated in splitless or on-column modes, incomplete combustion may result in contamination of the ion source.
- 2) It enables the selection of isolated peaks or groups of peaks to be 'Heart Cut' from a complex chromatogram so that they can be selectively diverted to the IRMS, whilst the unwanted species are vented to atmosphere via the FID.
- 3) It enables the introduction of pulses of reference gas on a clean background, whilst the effluent from the column is diverted away from the IRMS. When the Heart Split valve is in the open position, reference gas pulses can be introduced to the IRMS without fear of co-elution of the reference pulses and peaks from the sample preparation system.

CONNECTING THE GC

GAS CONNECTIONS

The GC requires a hydrogen and compressed air gas supply for FID operation. Helium carrier gas and compressed air supply comes from the GC5 (via “air to HS” blue output connection)

INSTALLING THE FURNACE TUBES


Before installing the GC oven components the furnace elements have to be installed. Carbon / nitrogen analysis uses a standard quartz tube filled with copper oxide, hydrogen uses a quartz tube filled with chromium and oxygen uses an inner nickel tube enclosed in a ceramic outer tube.

CARBON/NITROGEN ANALYSIS TUBE SET-UP


The furnace tube for carbon and nitrogen analysis consists of a long quartz glass tube filled with copper Oxide pellets which provides oxygen for the combustion process.

Note: Fingerprints on the tubes should be avoided so gloves should be worn at all times when handling the tube.

The tube slides in from the GC5 end and must be kept as straight and level as possible when inserting to avoid damage to the tube. One side of the tube has a silver wool packing that holds the copper oxide in place. Be sure not to hold the tube the wrong way round as the packing can fall out of the non-packed end. The end with the silver wool sits at the GC5 end and the open end towards the Agilent GC.



Before sliding into the furnace carefully attach this fitting onto the end of the tube as shown below.



This fitting consists of an extended 1/8" nut, T-piece, 1/4 nut and graphite 4-4 ferrule and a "pigtail" helium vent connection.


The 30cm 0.006" Stainless Steel (SS) "pigtail" vent that is fitted which should have a flow of approximately 1-5ml/min. This flow is affected by the sample line helium pressure.

Once this fitting has been attached securely to the tube it can now be slid into the GC5 furnace oven. The tube should slide easily through the assembly and emerge inside the GC oven.

Now a Zero Dead Volume T-piece assembly can be fitted to the other end of the tube. Firstly, if not already connected, attach the T-piece assembly to the sample line helium tube coming out from the side wall of the GC.

Using a GVF 4-4 Vespel ferrule, tighten the fitting to the furnace tube by hand. Now pull the tube back out of the fitting by approximately 5 mm. It is important that the glass tube does not bottom out in the fitting because it will be crushed once further tightening with spanners is required.

Once finger tight, tighten the assembly using two spanners making sure not to put any strains or stresses on the glass tube. The recommended technique is to use just one hand to tighten the nut while keeping the assembly fixed. Holding the assembly fixed with one spanner parallel to the floor, position the other spanner on the nut at approximately 60 degrees or less pointing upwards. Slowly pinch the spanners together until the spanner that is positioned on the nut reaches the spanner that being held fixed. This slowly and carefully tightens the nut without putting any additional strain on the tube and gives much more control in tightening the connection and will drastically reduce cracks and breakages of the tube.



Using approximately 50cm of 250 μ m ID fused silica, remove about 10 - 15cm of the outer polyimide coating fusing a flame (e.g. a cigarette lighter.) This is done so the coating of the silica isn't burnt inside the tube. Using acetone or similar and a cloth, wipe away the burnt silica. The thin inner coating that remains will be inserted into the furnace (from inside the GC oven). It is very brittle, so handle with care.


Common practice is to initially screw the extended nut (and ferrule) with the silica onto the T-piece connector until finger tight. This helps to keep the silica centred in the fitting and makes it easier to feed the silica into the hole of the furnace tube. Once the silica is in the hole, slide it all the way in until it sits just before the Copper Oxide packing. Tighten up the connection further to secure the silica in place.

COLUMN SET UP

HEART SPLIT VALVE

Connect approximately 25cm 320 μ m silica from the side connection of the Heart Split valve using a GVF 16-005 ferrule. Insert approximately 3cm of the other end of the silica into the FID using a graphite GFF 16-005 ferrule.

The ferrule needed for the Heart Split valve connection is a “backdrilled” GVF 16-003/004 ferrule (the hole is a different size at each end). The silica should only fit halfway through the ferrule. This is so that the Heart Split needle doesn’t make contact with the silica and damage it. If it passes all the way through, silica with a larger outer diameter (OD) is required.




Cut approximately 30cm of 250 μ m silica and insert one end into the bottom connection of the Heart Split valve and tighten the connection. Leave the other end of the silica free for now.

GC COLUMN

The GC comes with a 30m HP5 column of a 320 μ m I.D. with a film thickness of 25 μ m as standard. Users are requested to provide their own columns for their own applications. Hang the column in place on the supplied rack.

SAMPLE INJECTION PORT

Install one end of the column to the injector using a graphite GFF 16-005 ferrule. The silica should intrude no greater than 5mm into the injector. Be careful not to over tighten this connection as the graphite ferrule will easily compress and push graphite into the liner of the injector and block it.




The Split/Splitless Injector

The injector can be used in one of two modes; split or splitless. The injector contains a heated chamber containing a glass liner into which the sample is injected through the septum. In split mode the sample vaporises to form a mixture of carrier gas, vaporised solvent and sample. A proportion of this mixture passes onto the column, but most is diverted through the split outlet. The split proportion or ratio can be programmed by the user on the GC.

VSOS

The VSOS is a zero volume connector that is used to split the Column flow, diverting the sample to either the FID or furnace tube.

The VSOS is connected as shown below:



When fitting the silica lines to VSOS, note that the two pieces in the GVF 16(2)-004 ferrule (the silica from the bottom of the Heart Split and the silica from the furnace) must sit at their very furthest, flush to the bottom of the ferrule. Ideally they should sit just inside the ferrule.

If they protrude too far through the ferrule, tightening the nut will break or crush the silica and block the union. If they do not protrude enough, tightening the ferrule may close the ferrule beyond the silica and cut off the flow


When tightening the VSOS connection it is recommended to make the end with the two pieces of silica first, holding the nut fixed and rotating the VSOS union to tighten it onto the nut. This way the silica won't twist together and possibly break.

When connecting the column to the VSOS, make sure that the silica does not protrude beyond the ferrule on the other side as a blockage might occur.

NAFION MEMBRANE SYSTEM

The Helium flow through the membrane assembly is set at 20ml/min using the Water Trap Purge needle valve connected to the side of the GC5 interface.

A schematic of the nafion membrane system is shown below:



The nafion membrane assembly comprises a length of GVF tubing enclosed within an outer sheath of glass or acrylic tubing. The ends of the nafion tube are fixed, and the glass or acrylic outer sealed by a pair of SGE SSUT 4/16/16 unions. The fused silica connections to the furnace tube (variable length of 220µm ID), and to the standby valve (length of 100µm ID) are also fixed by the SSUT union.

The 220µm piece of silica is connected to the furnace tube protruding from the GC5. This is required to be inserted no more than 15cm into the tube. The 100µm piece of silica is connected to the IRMS inlet. Be sure to isolate the IRMS and turn off the source before making changes to the inlet.

Set the 3-way white/y valve on the front of the GC5 so that the sample line Helium flow is correct for the required analysis. Carbon/nitrogen /hydrogen are towards the right, oxygen towards the left.

Sample line Helium pressure is set to 2.0 – 2.5 psi and FID Helium to approximately 0.3 psi (the needle is just off the stopper "pip.")

Check the flow from the "pigtail" at the end of the furnace tube is approximately 1-5ml/min.

LEAK CHECKING

Once you are satisfied that all connections are sufficiently tight, carry out a leak check. Initially, a hand held leak checker can be used to find major leaks. Following this, check the background argon level on the IRMS by tuning for CO₂ and then peak jumping to mass 40. This should ultimately read less than $5\text{e}^{-12}\text{nA}$ on the Minor 1 scale at 200 μA trap current. The backgrounds can take some time to stabilise. It is recommended to leave the GC overnight for everything to settle.

It is important that there are no leaks before the system is heated.

Once the system has been conditioned another leak check will need to be performed.

GC CONFIGURATION

GC5 FURNACE TEMPERATURE

The combustion furnace and interface can now be switched on via the two switches on the front of the GC5. Depending on the type of analysis, the furnace temperature can be set manually. The temperatures for each type of analysis are shown in the table below. You can check the current set-point temperatures by holding the * button for 3 seconds. The temperature can be changed by holding the star button and pushing the Up or Down arrows.

If the unit hasn't been used for a while there may be a faint aroma. This is normal and will dissipate after a few minutes.

Element	Temperature C
Carbon	850
Nitrogen	950
Hydrogen	1000
Oxygen	1250

The interface temperature is always set at 350°C

The GC can now be configured.

COMMS

The GC communicates with the IonVantage software via a "Plug and Play" network switch. The GC and PC LAN connections connect directly into the switch and then, if available, the laboratory network connection can also be connected to the switch for network and internet connectivity. No further action is necessary to configure the switch and full internet and network capabilities should be available.

WITH ACTIVE INTERNET CONNECTION

If the computer is attached to a live network then an IP address or DHCP is assigned automatically to the GC. This can be checked via the front panel of the GC.

Select **Options > Communications > Enable DHCP**

Push **Yes/Enter** and you will then be asked to restart. Select **Yes**

Return to the **Communications** menu and you can view the assigned IP address of the GC. This address can then be inputted into the IonVantage software for communication with the GC and no further action is required.

WITHOUT ACTIVE INTERNET CONNECTION

If no active internet connection is available then the IP address has to be configured manually.


On your PC:

Start > Control Panel > Network Connections

Right click on **Local Area Connection > Properties**

Select Internet Protocol (TCP/IP) and then push the **Properties** button

You should see this window



Push the button that says **Use the following IP address:**

Input the following:

IP address	101.1.1.101
Subnet Mask	255.255.255.0
Default Gateway	Leave Blank

Leave the DNS server as automatic.

Now configure the IP address of the GC manually:

Select **Options > Communications**

Manually enter the IP address **101.1.1.100** and Subnet mask **255.255.255.0** and press **Enter/Yes**.
Restart the GC.

SETUP IONVANTAGE

Create or open a **GC project**

In the **Inlet Method** window:

Select the **Setup menu > Open > Select installed GC option** (e.g. GC7890-Autosampler)

Select **7890 Menu**

Edit Comms Settings > Enter the GC IP address > OK

You will now have to restart IonVantage.

SETUP AGILENT GC

This manual will now proceed to describe the Agilent 7890 GC settings as standard to Isoprime. All these settings are configured using the Agilent front panel. For further information on each of these settings please refer to your Agilent 7890 manual.

CONFIGURATION MENU

FRONT INLET

Equilibration Time > 1

Max Oven Temp > 325

Detector > Back Detector (Or MSD if installed and required)

COLUMN

Push Number **1**

Input the details of the column you are using into this menu. For example, the supplied column from Isoprime is:

Length > 30m

Diameter > 320um

Film thickness > 0.25um

FRONT PANEL BUTTONS

OVEN

Temp > 100

Initial time > 1

Ramp > 5

Final temp 1 > 280

Time > 60

COLUMN

Constant flow

Flow > 1.2ml

FRONT INLET

Temp > 250

Split > 20

BACK DET

Temp > 300

Flame **On.**

This begins the FID start sequence. It will wait until it reaches the set point temperature then automatically start the air and hydrogen mix and ignite the flame for the FID.

CONDITIONING

It is recommended after undertaking changes to the physical installation of the GC connections in the oven to condition the GC. By heating up the GC oven past its usual operating range, any impurities, moisture or residual solvent / sample left in the column can be removed. In essence this is similar to a bakeout on a mass spectrometer.

The oven settings programmed above in the Oven tab are the settings that Isoprime use to condition the GC. Under this configuration the oven will sit at 100C for 1 minute before commencing a slow ramp at 5°C per min up to 280°C where the oven will then sit for 60 minutes.

This final oven temperature is naturally sample dependant and 220°C is adequate to the needs of the Isoprime Universal Mix, where the oven settings rarely go above 220°C. Should you wish to condition at a higher temperature, the maximum oven temperature recommended for a HP5 column is 325°C.

By pressing **Status** on the GC front panel you can receive feedback on what the status of the GC oven. Once the GC status is “Ready for Injection” press **Start** on the GC front panel. As the GC being controlled by IonVantage, this action will **only** start the Oven ramp procedure and nothing else.

Once the run finishes, whilst the fittings are still hot; proceed to nip up all new connections made in the GC oven prior to conditioning. Most connections should need **no more than ¼ turn**. Some may need a little more.

Be aware that over tightening of the connections when the GC is still hot can cause the threads to shear or seize if over tightened.

RUNNING AN ANALYSIS

SETTING UP IONVANTAGE

This section will describe how to set up IonVantage and the Agilent GC for carbon analysis using the supplied standard Agilent 16 position autosampler carousel.


OPEN AN EXISTING OR CREATE A NEW GC PROJECT

Create or a new project using the project wizard based on the default GC template.

To create a new sample list, go to the sample list page and either:

- Click **File > New**

The following New Sample List Window will appear:



Several pre-installed sample list templates are available to the user.

- **New Blank Sample List:** creates a completely blank sample list and is only advisable to be used by expert users.

	File Name	MS File	Inlet File	Bottle	Inject Volume	Sample Type	Process	Process Options
1								

- **GC analysis:** creates a sample list which has all the necessary settings for a complete sample run for the different elements and both for manual as well as autosampler injections.

	File Name	MS File	Inlet File	Bottle	Inject Volume	Sample Type	Process	Process Options
1	Batch file name	Do nothing	Do Nothing	1	0.000		IsoPrimeDP	BatchStart
2	GC sample	GC-CO2 run	Manual injection	1	0.000		IsoPrimeDP	
3	GC sample	GC-H2 run	Manual injection	1	0.000		IsoPrimeDP	
4	GC sample	GC-N2 run	Manual injection	1	0.000		IsoPrimeDP	
5	GC sample	GC-CO2 run	Auto injection	1	0.000		IsoPrimeDP	
6	GC sample	GC-H2 run	Auto injection	1	0.000		IsoPrimeDP	
7	GC sample	GC-N2 run	Auto injection	1	0.000		IsoPrimeDP	
8	Batch End	Do nothing	Do Nothing	1	0.000		IsoPrimeDP	BatchEnd PrintReport

- **GC FID analysis:** creates a sample list which includes a sample run that sends all the GC effluents to the FID detector, used for determination of the operational timings of all valves.

	File Name	MS File	Inlet File	Bottle	Inject Volume	Sample Type	Process	Process Options
1	FID run	GC-FID run	Do Nothing	1	0.000			

- **Ref gas stability:** creates a sample list which includes all the necessary protocols for running reference gas performance tests, such as stability and H₃ correction factors.

	File Name	MS File	Inlet File	Bottle	Inject Volume	Sample Type	Process	Process Options
1	Batch file name	Do nothing	Do Nothing	1	0.000		IsoPrimeDP	BatchStart
2	Stability	CO2 ref gas stability	Ref Gas Stability	1	0.000		IsoPrimeDP	PrintReport
3	Stability	N2 ref gas stability	Ref Gas Stability	1	0.000		IsoPrimeDP	PrintReport
4	Stability	H2 ref gas stability	Ref Gas Stability	1	0.000		IsoPrimeDP	PrintReport
5	Measure H3	CF measure H3 correction	Measure H3 correction	1	0.000		IsoPrimeDP	PrintReport
6	Batch End	Do nothing	Do Nothing	1	0.000		IsoPrimeDP	BatchEnd PrintReport

All of these individual list functions can be combined into just the one sample list and for ease of use that is how this manual will proceed.

SAVING THE GC METHOD FILES

All of the Agilent GC Method programming such as oven ramps, temperatures and flows are done via the Agilent front panel. These settings are then manually retrieved by a user and saved from the GC into IonVantage as an Inlet file and are called up when running an analysis from the sample list.

Currently the GC settings should still be set as described in the last chapter (See page 32) the only settings that need to be changed for analysis are the Oven temperatures and ramps.

For testing carbon using the Universal mix the following GC Oven settings are used:

Temp	100°C
Initial Time	1 min
Ramp 1	15°C/min
Final Temp	200°C
Final Time	0.7 min


The GC starts at 100°C and holds at this temperature for 1 minute. The oven will then ramp at 15°C per minute until it reaches 200°C where it will then stay at this temperature for 0.7 minutes.

This gives a total run time of 8.367 minutes. This time will be needed to set up our IonVantage method later.

Once these temperature settings have been programmed into the GC return to IonVantage and resave all the GC Method settings to the software. This is done by going to the Inlet Method window.

On the left hand side of this window click Inlet.

At the bottom of the window there are 3 buttons as shown below:



Click “Get/Save Method from GC” and save this file with a name, for example, CO₂ AutoInjection.

All current GC settings are now saved to an IonVantage Inlet file. Whenever this Inlet file is selected in the sample list, IonVantage will automatically load and reprogram the GC with these settings. Any number of Inlet method files can be saved for various applications.

NOTE: If **any** alterations are made to the GC settings always ensure to Get/Save Method from the GC else on running the next analysis any alterations will be overwritten by the selected last saved Inlet method.

FID RUN

For optimum column efficiency, the sample should not be too large and should be introduced onto the column as a "plug" of vapour. Slow injection of large samples causes peak broadening and loss of resolution. The most common injection method is using an autosampler with a micro-syringe. The sample is injected into the injection port through a rubber septum into the liner where the sample is vapourised. The temperature of the sample port is usually about 50°C higher than the boiling point of the least volatile component of the sample.

The first step when running a new sample is to determine when to open and close the Heart Split valve.

Sending solvent to the furnace strips the oxygen off the furnace tube very quickly and reduces its lifespan very quickly as well as shortening the life of the IRMS ion source. By first sending the entire sample to the FID, a user can gauge when to open and close the Heart Split valve to initially send the solvent to atmosphere and then divert the required sample components to the IRMS for analysis.

To do this set up the sample list as shown below:

	File Name	MS File	Inlet File	Bottle	Injector	Inject Volume	Sample Type	Process	Process Options
1	FID Run 01	GC-FID run	CO2 AutoInjection	1	Front	10.000000		IsoPrimeDP	

File Name: Give the analysis a file name, i.e. FID Run 01. This is the name that will be given to the folder used to save all the data

MS File: Select the default GC-FID run

Inlet Method: Select the previously saved GC method i.e. CO2 AutoInjection


Bottle: Selects the sample bottle position in the Auto Sampler carousel. Make sure to check you have the right number set before starting analysis.

Injector: Can be Front or Back depending on the optional extras that may have been requested. Default is Front.

Inject Volume: The pre-set numbers are a **percentage** of the needle volume. Set this to 10.0000 for the FID run. (Using a 10µl syringe, this will inject 1µl)

Process: Use the standard IsoprimeDP data processing.

Once this has been done, FID analysis can be run. The analysis can be viewed in real time by clicking on **View > Data Display** on the left hand side of the IonVantage main window. A completed FID scan is shown below in the bottom window in blue.



There should be 5 peaks similar to the above. If there aren't any peaks or these peaks are very small please see the Troubleshooting section on page [56](#).


BREAKTHROUGH

If all connections have been made correctly and all flows set as specified, then all of the sample should divert to the FID. If there is a blockage or incorrect helium flow, for example, then some of the sample will “breakthrough” into the furnace tube and appear on the IRMS as shown in the top window of the scan page. This is undesirable as we want the **entire** selected sample going to either the FID or the IRMS. Please see the [56](#) for information on how to fix breakthrough.

SETTING UP HEART SPLIT TIMINGS

The first and what should be the broadest peak is the solvent peak, followed by 4 less broad peaks. These are the sample peaks. The IonVantage carbon analysis method needs to be edited in order to open and close the Heart Split around the 4 sample peaks, diverting these to the IRMS. The solvent is also directed away from the furnace and IRMS to the FID.

This is done by modifying the MS File. In the Sample List window in IonVantage set the file to GC-CO₂ Run from the dropdown box in the MS File column. Right click on GC-CO₂ run and open. A method window should open similar to the one below.



Note: GC time is in minutes and IRMS in seconds


This window defines the control method of the IRMS, that is to say when valves open and close, what species to measure, what tunings to use and when to take reference gas measurements, among many other advanced control parameters which are beyond the scope of this manual.

In this method there is a reference gas pulse, followed by the Heart Split open and close timings and then another reference gas pulse. For the purpose of measuring the Universal Mix, only one reference gas pulse is required so the second pulse can be deleted.

Double clicking on GCHeartSplit will open the command for editing. Looking at the previous FID run (shown above) a good time to close the Heart Split valve and divert the rest of the sample to the IRMS would be approximately 210 seconds. All of the solvent has now reached the FID and the signal is back to baseline.

When considering the time to open the Heart Split valve, the time it takes for the sample to reach the IRMS is longer than the time the sample takes to reach the FID as the sample now must pass through the combustion furnace, through the nafion (and/or cryogenic trap if measuring nitrogen) and then into the IRMS.


Common practice is to add approximately 30 seconds after the last peak has reached the FID. In this case setting the GCHeartSplit - Open to approximately 480 seconds is advised.



NOTE: If the IonVantage method finishes before the GC method the software will stall. Please ensure the IonVantage Method is longer than the GC Method by approximately 10 seconds.

The GC method time from earlier was 8.367 minutes. The GC method time must be slightly shorter than the IRMS method. So if the GC method is 8.367 minutes (502 seconds), setting the IRMS method to be about 510 seconds should be adequate.

The End of Acquire “delay” in the IRMS method can be increased or decreased accordingly to fit this criteria. In our example the Heart Split opens after 480 seconds and to ensure the IRMS method is 510 seconds, set the delay to 30 seconds. Our GC CO₂ Run method should look like the figure below.



Now the Heart Split timings have been provisionally set, a carbon analysis of the Universal Mix can now be performed. Save this method as, for example, GC-CO₂ Run and close the window.

PERFORMING A CARBON ANALYSIS


The sample list for an analysis is set up as per the FID Run except all files should now reference a CO₂ analysis.

	File Name	MS File	Inlet File	Bottle	Injector	Inject Volume	Sample Type	Process	Process Options
1	FID Run 01	GC-FID run	CO2 AutoInjection	1	Front	10.000000		IsoPrimeDP	
2	GC Carbon Analysis 01	GC-CO2 run	CO2 AutoInjection	1	Front	2.000000		IsoPrimeDP	

Give the Run a file name, select the previously modified GC-CO₂ Run, and select the GC method we wish to use (CO₂ AutoInjection). An injection volume of 2.000% is adequate. Ensure bottle number matches the carousel number in which the sample bottle is placed.

Select the analysis required to be run in IonVantage and press start.

A similar chromatogram should be obtained:



The top window shows the IRMS and bottom window the FID.

The top window contains the Reference Gas peak and the four sample peaks in the Universal Mix sample.

The bottom window shows the solvent peak on the FID. The two little “spikes” are simply the Heart Split valve opening/closing and sending a short pulse to the FID.

REOXIDISING THE CARBON/NITROGEN FURNACE TUBE

(Carbon/Nitrogen only (CuO))

After continual use of the GC5, the oxygen in the furnace tube is slowly depleted, causing smaller and smaller peaks, decreased quality of chromatography, poor results and/or if measuring nitrogen, a very large mass 30 background.

The furnace tube can be regenerating by reoxidising the tubes filling. This involves bleeding oxygen into the furnace tube oxidizing either the copper.

PROCEDURE

Please refer to the Isoprime User Manual if unsure how to implement some of the more advanced IonVantage procedures outlined in this procedure.

Drop the temperature of the GC5 furnace to 600°C. The regeneration process is more effective at this temperature.

Tune for Mass 33 by peak jumping from a N₂ tuning. By peak jumping to mass 33 this puts mass 33 in the Axial bucket or Minor 1 and thus mass 32 is now in the Low Mass bucket or Major, which has a lower gain and so we can observe bigger beams than in the Minor 1 bucket.

Open the O₂ flush valve on the Inlet Method page of IonVantage. The O₂ flush valve on the GC5 should already have been set to give a max beam size of approximately 10 – 12 nA of oxygen. If not then the oxidization process will need to be monitored to ensure that the tube and IRMS isn't flooded with oxygen.

Start a 3600 second time scan and monitor the O₂ level on the Major scale. Initially, the level of oxygen in the system will decrease as the tube absorbs the oxygen. As the tube becomes more oxidized the amount of absorption decreases and the level oxygen in the system will start to increase. Once the level oxygen stabilizes (this should be at about 10-12nA) then the tube is no longer absorbing oxygen and is fully oxidized, hence fully regenerated.

If the tube takes longer than 3600 seconds to regenerate simple start another, longer time scan and continue to monitor the rise.

If the oxygen level continues to climb past 12nA wind in the oxygen flush valve on the side of the GC5 until stable 10-12nA oxygen level is observed.

Note: The response time of the O₂ flush valve is much slower than that of the Reference Gas Box due to the extra length of tubing and silica capillary it must travel through to reach the IRMS. Please allow at least 30 seconds before troubleshooting any potential problem

GC5 OPTIONAL UPGRADES

NITROGEN ANALYSIS


Nitrogen analysis is carried out by using a cryogenic or “cold” trap. This consists of an insulated dewar flask filled with liquid Nitrogen, into which a piece of 320µm deactivated silica is submerged into the liquid nitrogen. This has the effect of freezing both H₂O, CO₂ and CO (arising from CO₂ combustion) to allow for accurate N₂ measurements. The system has no electronic temperature control and relies entirely on the liquid nitrogen temperature ($\approx -196^\circ\text{C}$) to capture both the H₂O, CO and CO₂.

SET UP FOR N₂ ANALYSIS

Note: Before making any changes to the physical connections on the instrument, remember to close the Nupro valve on the Isoprime (or divert to waste if a Standby valve is installed). This will prevent air reaching the system and prolong the lifetime of the filament, keep the instrument clean and enable the instrument to stabilise much quicker.


In the default Isoprime configuration no changes need to be made to the GC5 combustion tube as Carbon and Nitrogen analysis is performed using the same tube.

The cold trap assembly is set up as shown below:



The cold trap should be set up so that the two loops of the deactivated silica should be as low down in the dewar as possible so fewer refills are required.

The assembled cold trap is then added in series in between the end of the combustion tube and the Isoprime as shown below:



The nafion is not required for N₂ analysis and should be disconnected from the system. No conditioning of the trap is necessary.

Note: If constant switching between Nitrogen and Carbon analysis is required, having to change system configuration constantly can be a long process. The system *can* be run with the nafion and cold trap all connected together in series (without filling the trap with liquid nitrogen for Carbon analysis) however you may lose some sensitivity and chromatography by doing so. To run in this configuration simply disconnect the piece of silica going to the IRMS from the nafion and connect it to the cold trap in place of the 100µm silica as shown in the diagram above.

As ever, results of running in this configuration are application specific so it is advisable to set up the system with a configuration that works best for the user's requirements.

PERFORMING A NITROGEN ANALYSIS

When analysing the Universal mix, due to the smaller Nitrogen component in the mix compared to Carbon, a smaller inlet split ratio and a higher injection volume are advised to obtain good sensitivity and chromatography. Recommended values are a 10:1 split ratio set on the GC and 10.000% injection volume requested in IonVantage. After making these changes save the altered GC method (via the “Get/Save from GC” button in the Inlet Method window) as ‘N2 AutoInjection,’ for example.

As a result of incorporating the extra silica by the introduction of the cold trap to the system, the transport time for the sample entering the Isoprime will now be slightly longer. You may wish to run another FID analysis to determine new Heart Split timings as detailed in page 38. The IonVantage GC-N2 method file will need to be modified to take into account the new Heart split timings.

Once the Heart Split timings have been established, you may also need to increase your GC oven temperatures to factor in the longer analysis time. This is best achieved by increasing the final temperature of the GC oven. (See page 31 and 38 for further information)


Save any changes made to the GC method in IonVantage.

The sample list for Nitrogen analysis should be similar to the one shown below:

	File Name	MS File	Inlet File	Bottle	Injector	Inject Volume	Sample Type	Process	Process Options
1	FID Analysis 01	GC-FID run	N2 AutoInjection	1	Front	10.000000		IsoPrimeDP	
2	GC Nitrogen Analysis 01	GC-N2 run	N2 AutoInjection	1	Front	10.000000		IsoPrimeDP	

Select the analysis required to be run in IonVantage and press start.

A chromatogram similar to the one below should be obtained:



There should only be one large peak compared to the 4 peaks as seen in the Carbon analysis. This is the Nitrogen containing component of the Universal mix. All other carbon containing compounds will have been frozen in the cold trap.

If you see all 4 peaks please check that the cold trap contains enough liquid nitrogen.

HYDROGEN ANALYSIS

Analysing hydrogen requires a separate combustion tube to be installed. The tube is filled with chromium oxide instead of copper oxide used for carbon and nitrogen analysis. The tube needs to be filled manually and the amount and position of the filling should be identical to the copper oxide combustion tube.

SETTING UP FOR HYDROGEN ANALYSIS

When removing or changing the tube the GC5 furnace and interface should be left to cool to a safe, comfortable temperature (room temperature recommended) for the user to be able to handle the tube. To cool, simply turn off the green switches on the front of the GC5 furnace.

Note: Fingerprints on the tubes should be avoided so gloves should be worn at all times when handling the tube.

The installation of a tube should be performed as described on page [21](#).

Note: The nafion is not required for analysis but the system *can* be run with the nafion connected, however you may lose some sensitivity and chromatography by doing so. To run in this configuration simply connect the nafion as described for setting up for a carbon analysis (page [26](#))

As ever, results of running in this configuration are application specific so it is advisable to set up the system with a configuration that works best for the user's requirements.

When analysing hydrogen, the GC5 furnace is operated at 1000°C. The interface temperature is operated at 350°C

Once the tube has been installed it will need to be conditioned by heating the tube to operational temperature and leaving for approximately 6 hours.

PERFORMING A HYDROGEN ANALYSIS

Running a hydrogen analysis is very similar to running a carbon analysis (see page [40](#))


When analysing the Universal mix, due to the smaller hydrogen component in the mix compared to carbon, a smaller inlet split ratio and a higher injection volume are advised to obtain good sensitivity and chromatography. Recommended values are a 3:1 split ratio set on the GC and 10.000% injection volume requested in IonVantage. After making these changes save the altered GC method (via the "Get/Save from GC" button in the Inlet Method window) as 'H2 AutoInjection,' for example.

You may wish to run another FID analysis to check the Heart Split timings as detailed in page [38](#). The IonVantage GC-H2 method file will need to be modified to take into account any changes needed to Heart split timings.

Set up the sample list as shown below and run the analysis.

	File Name	MS File	Inlet File	Bottle	Injector	Inject Volume	Sample Type	Process	Process Options
1	FID Run 01	GC-FID run	H2 Autoinjection	1	Front	10.000000		IsoPrimeDP	
2	GC Hydrogen Analysis 01	GC-H2 run	H2 Autoinjection	1	Front	10.000000		IsoPrimeDP	

You should obtain a chromatogram similar to the one shown below. There will be a small rise in the baseline during the analysis. This is normal behaviour when analysing GC hydrogen:



OXYGEN ANALYSIS

Oxygen analysis is performed using a nickel tube enclosed within a ceramic outer tube. Hydrogen is used to provide a clean, dry reducing atmosphere to provide optimum conditions for pyrolysis of the sample. The oxygen component of carbon monoxide is used as the reference gas for analysis.


SETTING UP FOR OXYGEN ANALYSIS

Turn off the interface, furnace and GC oven and close the Nupro valve on the IRMS. After the system has cooled, turn the sample line and FID helium to zero.

Remove any previously used furnace tubes from the GC.

PREPARATION OF THE GC O TUBE


The GC oxygen tube consists of a nickel tube placed inside a 6mm ceramic tube, as shown below:



Note: The nickel tube must be 15mm longer than the outer ceramic tube so that the inner GC connections can be made.

GC OVEN CONNECTIONS

The connections to be made inside the GC oven are shown below:




It is important to seal the nickel tube in the correct position using the sealing ring inside the T-piece. If the nickel tube is not pushed further enough through the sealing ring, then the flow of the sample line hydrogen/helium gas can be interrupted and/or oxygen could leak into the nickel tube, preventing complete pyrolysis.

It is also important to note that a different connection for the sample line helium/hydrogen mix is required. This is an additional piece of stainless steel tubing that is connected to the hydrogen line in the GC5 and allows a supply of hydrogen for oxygen pyrolysis.

The 250µm fused silica has the polyimide coating removed for the first 20mm inside the furnace tube. Insert the fused silica no more than 50mm into the nickel tube.

FURNACE GC CONNECTIONS

The GC outer connections are as shown below:



The fused silica is to be inserted no more than 100mm inside the nickel tube. Due to the high operating temperature the fused silica can easily become too hot and melt blocking the He flow to the Isoprime.

CONDITIONING THE FURNACE

Once all connections have been made, ensure the whitey selector valve situated on the front of the GC interface is set to LEFT and connect a supply of hydrogen (the minimum grade to be used is specified in the site requirements document) to the rear of the GC5. Set the flow via the hydrogen pigtail vent situated to the left hand side of the GC5 to $\approx 5\text{ml/min}$.

Set the sample line pressure to 2psi and the FID helium flow to 0.3psi. Check the flow out from the outer connection pigtail. This should be read approximately 1 - 5ml/min.

The GC column will need to be conditioned. See page [32](#) for information on conditioning the column.

Leave the system overnight to allow all background levels to settle then leak check the system using argon (see page [27](#)), paying particular attention to the GC oven fittings.

When satisfied the system is leak tight set the furnace temperature to 1250°C and the interface temperature to 350°C . It is important to isolate the Isoprime during heating either by the use of a standby valve or closing the Nupro isolation valve on the Isoprime to avoid contamination.

CHECKING BACKGROUND LEVELS

With the trap current set at $600\mu\text{A}$ and the IRMS tuned for mass 28 (CO), run a magnet scan and ensure that the background levels are as shown in the table below:

Species	Mass	Beam (nA)
Water	18	$<1\text{e-}9$
Nitrogen	28	$<5\text{e-}10$
Oxygen	32	$<1\text{e-}10$
Argon	40	$<5\text{e-}12$
Carbon Dioxide	44	$<5\text{e-}11$

If the above backgrounds cannot be achieved, leak check the system once more (using argon preferably) and slowly tighten any loose fittings in the GC5 (be careful not to over-tighten these though)

It is possible that the backgrounds may still be dropping, albeit slowly and are still currently slightly outside of the above expected values. If this is suspected to be the case, set up a time scan for 60 minutes and monitor the mass. If the background is still falling and hasn't stabilised continue to leave the system until it has settled.

If the background has reached its lowest level and is still outside of specification, the furnace will have to be cooled down and the furnace tube connections remade.

Check the IRMS backgrounds with the GC heart split valve both open and closed. It is very important to have a leak free system before and analysis can proceed. When the system is first setup it will take approximately 24 hours for the background to reach its lowest level.

When all the above requirements are reached the system is ready to begin analysis.

SETTING UP THE GC

On the Inlet Method window in IonVantage, open the H2 flush valve to add a stream of hydrogen into the sample line helium. The H2 flush valve must remain open during any GC oxygen analysis. This will provide the reducing atmosphere inside the furnace tube and enable optimum oxygen pyrolysis.

GC INLET METHOD

Set the GC method parameters as shown below (see page [30](#) for further information)

Parameter	Setting
Split Ratio	20:1
Injector Temp	260°C
FID	300°C
Oven Initial Temp	100°C
Time 1	1 min
Rate 1	15°C/min
Oven Final Temp 1	200°C/min
Final Time 1	2 min
Column Flow	1.5 ml/min

Save the new GC method (via the “Get/Save from GC” button in the Inlet Method window) as ‘O2 AutoInjection,’ for example. When saving the GC method to IonVantage, ensure that the H2 flush valve is open.

PRE-ANALYSIS

Before oxygen analysis can be carried out, it is necessary to have a layer of carbon on the inside of the nickel furnace tube. This can be achieved by running solvent (such as the hexane solvent used in the Universal mix) through the nickel tube. Running a number of these injections will slowly build up a significant carbon source inside the nickel tube causing the carbon molecules in a sample to break their bonds from the oxygen molecule, leaving just the oxygen component to be analysed.

For the first few injections it is extremely important to have the Nupro valve on the Isoprime closed so the large amount of solvent injected through the system doesn’t contaminate the IRMS. When carbon is first laid onto the tube, the background will become very high. This should slowly fall over a number of hours. As more carbon is laid onto the tube, this effect will be less evident. When the tube is ready for analysis the injected solvent will create a small peak that will quickly return to the base background level.

COATING THE TUBE

The initial laying of carbon in the tube should be done manually. Remove the autosampler if installed and set up the GC method as described above. Press start on the front of the Agilent GC panel. This will start the GC method. IonVantage need not be used at this point.

Begin to manually inject approx 2µl of the Universal mix (or other suitable solvent) into the system approximately every 2 minutes. This will start to coat the inside of the nickel tube with carbon contained in the solvent.

Note: Be careful not to inject too much, too fast as there is a risk of blocking the tube with too much carbon.

Perform 10 – 15 injections of the solvent before opening the Nupro valve and opening the GC to the IRMS.

Begin another manual run on the GC but this time only inject a tiny amount of solvent, < 0.1µl if possible. Monitor the solvent peak on the Isoprime using a time scan on mass 28 (CO) If the nickel tube has enough carbon coated on it then this solvent peak shouldn't saturate the collector or trip the source. If this is the case keep injecting small amounts of solvent until the size of the solvent peak continues to decrease.

If the collector does saturate and/or the source trips, close the isolation valve and repeat the 2µl injections 10 - 15 times and then repeat the above steps again. This may take some time but the solvent peak size will eventually start to diminish.

Monitoring the solvent peak and the 2 alkane peaks in the Universal mix will give a good indication of when the system is ready for analysis. These 3 peaks should become smaller and eventually become insignificant when there is enough carbon in the tube. When only the methyl ester and the nicotine peak can be seen the system is ready for analysis.


PERFORMING AN OXYGEN ANALYSIS

Oxygen analysis is done using a 600 μ A trap current, so it may be frugal to retune the source before commencing analysis.

If removed, replace the autosampler now. Modify or create an IonVantage oxygen analysis method using the timings gathered from the pre-analysis conditioning run and set up a sample list similar to the one shown below. See page 34 for further details on creating methods:

	File Name	MS File	Inlet File	Bottle	Injector	Inject Volume	Sample Type	Process	Process Options
1	GC Oxygen Analysis 01	GC-O run	O2 AutoInjection	1	Front	10.000000		IsoPrimeDP	

Run the analysis. You should obtain a chromatogram similar to the one shown below:



Note: The first peak after the reference peak is all the solvent from the Universal mix going to the IRMS. This gives an indication of the level of carbon that should be lay onto the nickel tube before accurate oxygen analysis can be performed. The largest peak is the oxygen containing methyl ester compound and the last peak is nicotine, seen here due to the instrument looking at mass 28.

TROUBLESHOOTING

Very small or No peaks on FID	Is FID flame ON	
	Is Injection port liner blocked	See Agilent manual for further details on maintenance Remake connection
	Is FID blocked	
	Is silica blocked or crushed	See Agilent manual for further details on maintenance Cut and remake connection
	Check silica in VSOS connection is not crushed or blocked	Remake connection taking extra care
	Is syringe plunger attached to autosampler via thumbscrew	Check autosampler
	Check syringe is not damaged or blocked	Change syringe
	Check sample quantity and syringe reaches sample	
Check Adequate Injection Volume		In IonVantage check Inject Volume column is set correctly

Very Small or No Peaks on IRMS	Run through No Peaks on FID section above	
	Is furnace tube blocked	Cool down furnace, remove tube and check for breakages, blockages and/or broken silica
	Is Nafion tube broken or blocked	Remake Nafion assembly
	Check IRMS inlet connection	Check and remake IRMS inlet connections i.e. ferrules and silica
	Check sample line Helium flow	Set to 2.0 – 2.5 psi
	Is stripped silica correctly inserted into furnace tube	
See also Isoprime troubleshooting guide		

FID won't ignite	Does supplied flow match requested flow		Check gas supply/quantity of gas
	Is FID blocked		Remove blockages from top of FID See Agilent manual for further details on maintenance
	Check FID ignition		See Agilent manual for further details on maintenance
	Increase or decrease hydrogen and Air flows.		Repeat with some trial and error
	Leave gases purging through FID for 1hr and try again		

Breakthrough to IRMS	Is FID Helium pressure too high		Reduce to 0.3psi
	Check blockages and crushed silica		Cut and remake silica connections
	Is Heart Split Valve closing (See also Heart SplitSection)		Place silica from bottom of Heart Split in Acetone or similar. When Heart Split is closed no bubbles should be observed

Heart Split Valve	Heart Split doesn't open/close	Check compressed air and hiss when opening/closing
		Place silica from bottom of Heart Split in Acetone or similar. When Heart Split is closed no bubbles should be observed
		Change and remake ferrule and silica fitting
		Check ferrule sits flush inside Heart Split
		Change Heart Split Assembly

IonVantage Issues	For any error messages and/or non-ready status's		If you have a Multimode Injector (MMI) and/or CombiPal system installed please ensure latest version of IonVantage is 1.5.4 for full compatibility
	No Communication		Check IP address in IonVanatge is the same as IP address set on the GC. IP address of PC must be different.
			Check network switch is turned on and all network cables are connected