

Comprehensive Analysis of Residual Solvents in Water-Soluble Articles Using Dynamic Headspace

Richard Lake and Frank Dorman
Restek Corporation, Bellefonte, PA

Eric Heggs
Teledyne Tekmar, Mason, OH

Restek Corporation
www.restekcorp.com



Background




- ICH (Q3C) lists 61 compounds as residual solvent impurities classified by toxicity
 - Class 1 – should be avoided
 - Class 2 – should be limited
 - Class 3 – can be assayed by non-specific techniques – LOD
- USP <467> compendial method
- Selectivity– no single analytical column
- Analyte confirmation – sometimes need multiple columns or MSD
- Sensitivity - Static headspace
 - Great fit for application, not exhaustive or ideally sensitive

1. Example of Water Soluble Articles
2. Advantages of Dynamic Headspace (Sensitivity)
 1. Added sensitivity brings opportunity
 2. More representative sample
3. Comprehensive dual column FID analysis using Modulated Accelerated Column Heating (Selectivity and Confirmation)
 1. Enhance Selectivity – confirmation analysis
 2. Enhance Sensitivity - LOD
 3. Linearity & Range
 4. Reproducibility
 5. Time Saving

Headspace Analysis



BASIC HEADSPACE CONCEPT



- A liquid or solid sample is heated in a sealed vial
- Equilibrium is established between the sample and the gas phase (headspace)
- Aliquot of the gas phase is transferred to the GC
- Not an exhaustive technique


Static Headspace



Pressurized Loop System


Step 1

Sample reaches
equilibrium/pressurization




Step 2

Sample is extracted from
headspace




Step 3

Sample is injected



- Only a portion of the headspace is sampled
- Not optimizing sensitivity

Dynamic Headspace



Static vs. Dynamic



Direct Sensitivity Comparison Using USP <467> Solvents

Instrument Conditions USP <467>

Sample: *Heated for 45 minutes at 80°C*

Injector 180°C, split 2:1


Carrier Helium @ 5.0
Flow ml/min, constant

Column Rtx-G43, 30m X
0.53mm x 3.0 µm

Detector FID 250°C

Oven 40°C for 20
Program minutes,
to 240°C at
25°C/min,
hold for 10 min

**G43 = 6% cyanopropyl
94% dimethylpolysiloxane**



Static vs. Dynamic



Sensitivity Comparison Using USP <467> Solvents

Average Peak Area

Analyte	Sample Conc (ppm)	Reg Limit (ppm)	Static Headspace	Dynamic Headspace	Increase in Sensitivity
Dichloromethane	12.0	600	619	18679	30 X
Chloroform	1.2	60	39	783	20 X
Benzene	0.04	2	15	313	21 X
Trichloroethylene	1.6	80	141	3479	25 X
1,4-dioxane	7.6	380	20	272	13 X

Average increase in sensitivity = 22X

Why is the increase lower for 1,4-Dioxane?



**Increase in sensitivity
is only 13X**

- Very hydrophilic compound
- Poor partitioning
- Indicates that partitioning is not the same for all compounds
- Similar for alcohols, ketones, aldehydes

Class 1 and 2 Solvents



Table I Class 1 solvents (should be avoided).

Solvent	Concentration limit (ppm)	Associated hazard
Benzene	2	Carcinogen
Carbon tetrachloride	4	Toxic and environmental hazard
1,2-Dichloroethane	5	Toxic
1,1-Dichloroethene	8	Toxic
1,1,1-Trichloroethane	1500	Environmental hazard

Table II Class 2 solvents.

Solvent	Concentration limit (ppm)
Acetonitrile ←	410
Chlorobenzene ←	360
Chloroform ←	60
Cyclohexane ←	3880
Dichloromethane (methylene chloride) ←	600
1,4-Dioxane ←	380
1,1,2-Trichloroethylene ←	80
1,2-Dimethoxyethane ←	100
2-Ethoxyethanol	160
2-Methoxyethanol	50
Methylbutyl ketone ←	50
Nitromethane ←	50
Sulfolane	160
Tetralin	100
Pyridine ←	200
Toluene ←	890
Formamide	220
1,2-Dichloroethene ←	1870
<i>N,N</i> -Dimethylacetamide	1090
<i>N,N</i> -Dimethylformamide	880
Ethylene glycol	620
Hexane ←	290
Methanol ←	3000
Methylcyclohexane ←	1180
<i>N</i> -Methylpyrrolidone	4840
Xylene ←	2170

- Assayed Class 2 Mix A and B on G43
- 6ml Total Volume, 1 g Sodium Sulfate in 20 ml vial
- Loop at regulatory limit (1ml standard added)
- Trap at ½ regulatory limit (500µl standard added and increased split ratio)

Analyte List




1. Methanol
2. Acetonitrile
3. Dichloromethane
4. trans-1,2-Dichloroethylene
5. Hexane
6. cis-1,2- Dichloroethylene
7. Nitromethane
8. Tetrahydrofuran
9. Chloroform
10. Cyclohexane
11. 1,2-Dimethoxyethane
12. Trichloroethylene
13. Methylcyclohexane
14. 1,4-Dioxane
15. Pyridine
16. Toluene
17. 2-Hexanone
18. Chlorobenzene
19. Ethyl benzene
20. m-Xylene
21. p-Xylene
22. o-Xylene
23. Tetralin

**23 solvents with varying
chemical properties**


Static vs. Dynamic – A closer look



Traditional Loop Method – Rtx-G43 Class 2 Mix A and B



Trap Method – Rtx-G43 Class 2 Mix A and B



Increased Split Flow



Loop Method
Rtx-G43

5:1 Split Inlet

DCM

$R_{ACN/DCM}$
= 1.09

ACN

*Resolution
must be >1.00*

Trap Method
Rtx-G43

20:1 Split Inlet

DCM

$R_{ACN/DCM}$
= **1.66**


ACN

4.40 4.50 4.60 4.70 4.80 4.90 5.00

Time (min)

4.20 4.30 4.40 4.50 4.60 4.70 4.80 4.90

Time (min)



Sensitivity Comparison




Compound	Polarity	Static Analysis		Dynamic Analysis at ½ Conc.		Increase (times)	
		Peak Area	Peak Height	Peak Area	Peak Height	Peak Area	Peak Height
Methanol	Polar	44.57097	9.68596	235.03558	83.45882	5.3	8.6
Acetonitrile	Polar	16.65052	3.79775	89.66702	26.12947	5.4	6.9
1,4-dioxane	Polar	6.83236	0.63490	14.84858	1.35838	2.2	2.2
MBK	Polar	21.91312	6.59006	35.65511	10.61567	1.6	1.7
DCM	Slightly	190.38097	41.14470	497.52942	135.51151	2.6	3.3
Hexane	Non	356.69006	65.75023	714.94470	160.13466	2.0	2.5
<i>Calculated Average Increase of all Analytes Using ½ Conc.</i>						5.6	6.0
			<i>Average Increase X4 (Split)</i>			22	24

Same average increase as noticed before
Same uneven partitioning

Effect of Phase Ratio



Polar Solvent (1/5 of regulatory limit)




Area of analyte increases as volume of sample decreases

Effect of Phase Ratio



Less Polar Solvent (1/5 of regulatory limit)



Area of analyte decreases as volume of sample increases


Effect of Phase Ratio




Effect of Phase Ratio on Analyte Response

Acetonitrile and Dichloromethane

1,4-dioxane



10 ml to 0.5 ml range




10 ml to 0.5 ml range

Effect of Phase Ratio



Effect of Phase Ratio on Analyte Response



Smaller sample volumes are more representative of polar compounds and give better sensitivity

Partition Coefficient




Traditional Headspace Theory - Loop


$$\text{Partition Coefficient (K)} = C_s / C_g$$

C_s = Concentration of analyte in sample phase

C_g = Concentration of analyte in the gas phase



Sensitivity is increased when K is minimized.





Partition Coefficient



Headspace Trap Theory

K is minimized as the volume of sample and standard decreases,

Partitioning is enhanced with a smaller, less concentrated sample size



*80°C- 90°C ceiling
for platen temp*

Increasing Sensitivity





Traditional Headspace Theory - Loop

$$\text{Phase Ratio } (\beta) = V_g/V_s$$

V_g = Volume of the gas phase

V_s = Volume of the sample phase





Sensitivity is increased when β is minimized.

Increasing Sensitivity



Headspace Trap Theory

Phase Ratio (β) is now the reciprocal as an increased gas phase does not dilute sample



Sensitivity is increased when β is maximized.

Advantages:

1. Large Increase in Sensitivity
2. More Representative Sample
3. Smaller Volumes = Smaller Sample Sizes
= Faster sample
equilibrium times

Makes Possible:

1. More Versatility in Instrument Conditions
2. Dual Column FID Method for Analyte
Confirmation – *using two alternatively
selective columns*

Headspace Conditions



Dynamic Headspace (Trap) Conditions	
Instrument	Teledyne Tekmar HT3
Valve Oven Temp	220°C
Transfer Line Temp	220°C
Standby Flow Rate	50 ml/min
Trap Standby Temp	40 °C
Platen/Sample Temp	80 °C
Sample Equil. Time	15.00 min
Mixer Time	2.00 min
Mixing Level	5
Mixer Stabilize Time	0.50 min

Sweep Flow Rate	75 ml/min
Sweep Flow Time	5.00 min
Dry Purge Time	10.00 min
Dry Purge Flow	100 ml/min
Dry Purge Temp	25°C
Desorb Preheat	245 °C
Desorb Temp	250°C
Desorb Time	1.00 min
Trap Bake Temp (K)	260°C
Trap Bake Time	6.00 min
Trap Bake Flow	450 ml/min

Dual Column Assay



Instrument Agilent 6890 GC/FID

Injector, Split
Mode 220°C

Split Ratio 20:1

Column Flow Helium, constant flow,
split into two columns
A: 1.5 ml/min
B: 1.6 ml/min

Column A: Rtx-G43, 30m X 0.32mm x 1.8 µm
B: Stabilwax, 30m X 0.32mm x 1.0 µm

Detector, FID 250°C


Hydrogen 40 ml/min

Air 450 ml/min

Make-up 45 ml/min

**G43 = 6% cyanopropyl
 94% dimethylpolysiloxane**


**G16= 100% Polyethylene
 glycol (PEG)**




Dual Column Assay



Dual Column Assay– Rtx-G43
Class 2 Mix A and B



Dual Column Assay– Stabilwax (G16)
Class 2 Mix A and B



Dual Column Assay



- Coelutions:

G43

1. Cis-1,2-dichloroethane / Nitromethane
2. Pyridine / Toluene
3. m-xylene / p-xylene (critical)

G16

1. Tetrahydrofuran / Trans-1,2-dichloroethane
2. Cis-1,2-dichloroethane / Trichloroethylene (critical)

- One oven hinders selectivity
- Lengthy run times

Instrument Configuration



Instrument Configuration



Comprehensive Method



Headspace Method – same as previous


0.5 ml and 1 ml total volume
in 20ml headspace vial

Instrument	Agilent 6890 GC/FID/FID
Injector, Split Mode	220°Cm split ratio 20:1
Column Flow	Helium, constant flow, split into two columns (0.32mm guard)
Oven Program	220°C

Column A (G43) – Rtx-624
20m x 0.18mm x 1.0µm
Flow: 0.85ml/min


Column B (G16) – Rtx-WAX
20m x 0.18mm x 0.4µm
Flow: 0.99ml/min

Oven Program





FID A - 250°C
Hydrogen 40ml/min
Air 450 ml/min
Makeup (He) 45 ml/min

Oven Program



FID B - 250°C
Hydrogen 40ml/min
Air 450 ml/min
Makeup (He) 45 ml/min

Method Results




Equal or Enhanced Selectivity – G16

Dual Column Gerstel Assay
Rtx-WAX, 20m x 0.18mm x 0.4 μ df


THF

trans -1,2-DCE



Dual Column Assay
Rtx-WAX, 30m x 0.32mm x 1.0 μ df

THF /
trans -1,2-DCE





Method Selectivity





Equal or Enhanced Selectivity – G43

Dual Column
30m x 0.32mm x 1.0 μ df



Dual Column Gerstel
20m x 0.18mm x 0.4 μ df



Method Sensitivity



Sensitivity Comparison Using Peak Height

Method	Sample Volume	Average Increase of all Analytes
Standard Loop	6 ml	---
Standard Trap	6 ml	24 X
Gerstel Dual Column Assay	1 ml	313 X



Cumulative effect of trapping, sample volume, and peak height increase using microbore columns

Effect of Polarity on Sensitivity


Solvent	Polarity	Increase
Methanol	Polar	711 X
Acetonitrile	Polar	801X
MBK	Polar	79X
DCM	Slightly Polar	74X
Hexane	Non Polar	29X

Method LOD




Sensitivity of 1,4-dioxane as a function of signal-to-noise

Traditional Loop Method
19 μg in 6ml volume
20 ml headspace vial



Trap-Mach Method (split)
19 ng in 0.5ml volume
20 ml headspace vial



All factors considered, approximately a 700 fold increase in signal response for polar compounds, even when column flow is split

Method Linearity and Range



1 ml Sample Volume

Linearity –in 20ml headspace vial

Correlation Coefficient >0.994 for all analytes

Average Correlation Coefficient = 0.9984

Range – various low responders - 1,4-dioxane = $0.19\mu\text{g}$ to $1.9\mu\text{g}$
- MBK = $0.03\mu\text{g}$ to $0.3\mu\text{g}$

0.5 ml Sample Volume

Linearity –in 20ml headspace vial

Correlation Coefficient >0.995 for all analytes

Average Correlation Coefficient = 0.9982

Range – various low responders - 1,4-dioxane = $0.019\mu\text{g}$ to $0.76\mu\text{g}$
- MBK = $0.003\mu\text{g}$ to $0.12\mu\text{g}$

Method Reproducibility



Acetonitrile Reproducibility

40µl of standard / 0.5ml sample
1.64µg/ml or 4% of reg limit
@ reg limit using 10mg sample

Rep	Solvent	Peak Area	Peak Height	RT (min)
1	ACN	37.8979	32.92955	1.556
2	ACN	35.9404	31.11098	1.556
3	ACN	37.76856	32.61021	1.554
4	ACN	37.16529	31.81859	1.555
5	ACN	35.05855	30.41616	1.554
6	ACN	38.60516	33.66129	1.553
	<i>Avg</i>	37.0726	32.0911	1.555
	<i>Std Dev</i>	1.331	1.207	0.001
	<i>%RSD</i>	3.591	3.760	0.078

Dichloromethane Reproducibility

40µl of standard / 0.5ml sample
2.40µg/ml or 4% of reg limit
@ reg limit using 10mg sample

Rep	Solvent	Peak Area	Peak Height	RT (min)
1	DCM	40.34903	28.88841	1.599
2	DCM	41.48586	30.95327	1.599
3	DCM	40.88852	30.27963	1.557
4	DCM	41.47880	30.21242	1.557
5	DCM	41.25531	30.57424	1.557
6	DCM	41.12213	30.62504	1.556
	<i>Avg</i>	41.0966	30.422	1.598
	<i>Std Dev</i>	0.430	0.373	0.001
	<i>%RSD</i>	1.047	1.226	0.077

Method Time



USP method:

2 sample preparations at 60 minutes each
= 2 hours/sample

Headspace Trap with MACH:

1 sample preparations = 33 minutes /sample
(Limited by headspace analysis)

Approximate reduction of 1.5 hours / sample

Conclusion



- Better sensitivity - sample sizes of 25mg to 10mg possible
- Selectivity and Confirmation - Choosing two columns with alternate selectivity gives analyte confirmation for all ICH compounds
- More representative analysis of polar compounds
- Shorter analysis and sample equilibration times - Higher throughput

Thank You