

INTERNAL DIAMETER	MICROLITERS PER INCH	INCHES PER MICROLITER
0.001	0.0129	77.6979
0.002	0.0515	19.4245
0.003	0.1158	8.6331
0.004	0.2059	4.8561
0.005	0.3218	3.1079
0.006	0.4633	2.1583
0.007	0.6306	1.5857
0.008	0.8237	1.2140
0.009	1.0425	0.9592
0.010	1.2870	0.7770
0.011	1.5573	0.6421
0.012	1.8533	0.5396
0.013	2.1751	0.4598
0.014	2.5226	0.3964
0.015	2.8958	0.3453
0.016	3.2948	0.3035
0.017	3.7195	0.2689
0.018	4.1700	0.2398
0.019	4.6462	0.2152
0.020	5.1481	0.1942
0.021	5.6758	0.1762
0.022	6.2292	0.1605
0.023	6.8084	0.1469
0.024	7.4133	0.1349
0.025	8.0440	0.1243
0.026	8.7003	0.1149
0.027	9.3825	0.1066
0.028	10.0903	0.0991
0.029	10.8239	0.0924
0.030	11.5833	0.0863
0.031	12.3684	0.0809
0.032	13.1792	0.0759
0.033	14.0158	0.0713
0.034	14.8781	0.0672
0.035	15.7662	0.0634
0.036	16.6799	0.0600
0.037	17.6195	0.0568
0.038	18.5847	0.0538
0.039	19.5758	0.0511
0.040	20.5925	0.0486
0.041	21.6350	0.0462
0.042	22.7032	0.0440
0.043	23.7972	0.0420
0.044	24.9169	0.0401
0.045	26.0624	0.0384
0.046	27.2336	0.0367
0.047	28.4306	0.0352
0.048	29.6532	0.0337
0.049	30.9017	0.0324
0.050	32.1758	0.0311

INTERNAL DIAMETER	MICROLITERS PER INCH	INCHES PER MICROLITER
0.051	33.4757	0.0299
0.052	34.8014	0.0287
0.053	36.1527	0.0277
0.054	37.5299	0.0266
0.055	38.9327	0.0257
0.056	40.3613	0.0248
0.057	41.8157	0.0239
0.058	43.2958	0.0231
0.059	44.8016	0.0223
0.060	46.3332	0.0216
0.061	47.8905	0.0209
0.062	49.4735	0.0202
0.063	51.0822	0.0196
0.064	52.7167	0.0190
0.065	54.3770	0.0184
0.066	56.0630	0.0178
0.067	57.7747	0.0173
0.068	59.5122	0.0168
0.069	61.2754	0.0163
0.070	63.0643	0.0159
0.071	64.8790	0.0154
0.072	66.7195	0.0150
0.073	68.5856	0.0146
0.074	70.4775	0.0142
0.075	72.3952	0.0138
0.076	74.3386	0.0135
0.077	76.3077	0.0131
0.078	78.3026	0.0128
0.079	80.3232	0.0124
0.080	82.3696	0.0121
0.081	84.4417	0.0118
0.082	86.5395	0.0116
0.083	88.6631	0.0113
0.084	90.8124	0.0110
0.085	92.9875	0.0108
0.086	95.1882	0.0105
0.087	97.4148	0.0103
0.088	99.6670	0.0100
0.089	101.9450	0.0098
0.090	104.2488	0.0096
0.091	106.5783	0.0094
0.092	108.9335	0.0092
0.093	111.3145	0.0090
0.094	113.7212	0.0088
0.095	116.1537	0.0086
0.096	118.6119	0.0084
0.097	121.0958	0.0083
0.098	123.6055	0.0081
0.099	126.1409	0.0079
0.100	128.7020	0.0078