

Australian / International FERRULES - Price List July 2014 ➤ current 2022+

Graphite Ferrules - Straight - in USD / AUD

Custom GFs

Catalog #	Fitting Size	Hole ID	Description	USD/AUD	QTY
In AUSTRALIA Prices are the same but in AUD					
SF100G	1/16"	1/16"	Ferrules, Graphite, 1/16"ID	25.00	10-pk
GF100	1/16"	1/16"	Ferrules, Graphite, FRONT, 1/16"ID	25.00	10-pk
SF200G	1/8"	1/8"	Ferrules, Graphite, 1/8"ID	25.00	10-pk
GF200	1/8"	1/8"	Ferrules, Graphite, FRONT, 1/8"ID	25.00	10-pk
SF300G	3/16"	3/16"	Ferrules, Graphite, 3/16mmID	25.00	10-pk
SFM5G	5mm	5mm	Ferrules, Graphite, Varian 5mmID	25.00	10-pk
GF400	1/4"	1/4"	Ferrules, Graphite, FRONT, 1/4"ID	25.00	10-pk
SF400G	1/4"	1/4"	Ferrules, Graphite, 1/4"ID	25.00	10-pk
SF400SPEC	1/4"	1/4"	Ferrules, Graphite, 1/4"ID	30.00	10-pk
SF600G	3/8"	3/8"	Ferrules, Graphite, 3/8"ID	36.00	10-pk
SFM10G	10mm	10mm	Ferrules, Graphite, 10mmID Metric	36.00	10-pk
SFM12G-1	12mm	12mm	Ferrules, Graphite, 12mmID Metric	4.20	ea
SF800G	1/2"	1/2"	Ferrules, Graphite, 1/2"ID	42.00	10-pk
SFM14G-1	14mm	14mm	Ferrules, Graphite, 14mmID Metric	4.80	ea
SF1000G-1	5/8"	5/8"	Ferrules, Graphite, 5/8"ID	\$5.60	ea
SF1200G-1	3/4"	3/4"	Ferrules, Graphite, 3/4"ID	5.60	ea
SF1400G-1	7/8"	7/8"	Ferrules, Graphite, 7/8"ID	6.80	ea
SF1600G-1	1.0"	1.0"	Ferrules, Graphite, 1.0"ID	8.10	ea
SFM25G-1	25mm	25mm	Ferrules, Graphite, 1-1/4"ID	8.10	ea
SF2000G-1	1-1/4"	1-1/4"	Ferrules, Graphite, 1-1/4"ID	12.00	ea
SF2400G-1	1-1/2"	1-1/2"	Ferrules, Graphite, 1-1/2"ID	15.00	ea
SF2800G-1	1-3/4"	1-3/4"	Ferrules, Graphite, 1-3/4"ID	18.00	ea
* SF3200G-1	2.0"	2.0"	Ferrules, Graphite, 2.0"ID NOW as a Front Ferrule ONLY	21.00	ea

if NOT Listed ➤ a Tooling Charge May apply !

NOTES : Graphite Ferrules

2 Sizes of HP ferrules for HP/Agilent 5890

GF100/- - HP size to fit standard HP inlet seal nut (0.125-in long)

GF- - HP a short version to fit shallow recess HP nut (0.090-in Long)

Custom Ferrules : ANY size to order - providing it is physically feasible for us to make.

We need accurate dimension "specs" OD / ID / Overall Length / Shoulder Length and OD at Cone Tip. Note : For Swagelok Ferrule shapes the cone angle is normally 17-deg

There may be a "one-off" tooling charge of USD90(<1/2"ID) to- \$150 depending on size

- if this is required for a "special" ie not a stock item.

Discounts : apply depending on **Total \$ Order** value or on Qty(each) . . . at **your** best rate!

USD10-99, 0%; USD100-500, 10% . . . Also OEM and Distributor Discounts : **.Enquire !**

Qty Discount

min10-99 less 0%
100-499 less 10%
500-1000 less 20%
1000+ less 25%
in \$ or each
at your best rate

Plus Freight Cost

Int'l Air Courier
USD75 per 500g
\$95/1kg . . . at cost

Delivery Times : Availability - most standard sizes are ex stock in reasonable Qty's

Production Lead Time : most sizes are ex-stock or 4-5 days production for qty to 500 each;

larger qty approx 7+ days from date of order . . . PLUS transit time

Shipping :

In Australia via Local Courier Aud 5 per 3kg; \$30 Interstate metro

International - Air Express Courier Door-to-door.

USA : USD 75.00 per 500g, USD95/ 1kg,

Transit time : approx 3 to 5 days (most countries eg. UK, USA, Canada . . .)

Other Countries - **please enquire !**

Payment by Bank Transfer or Credit Card via on-line PayPal (we send a separate PP Invoice on request)

. . . plus Company Name; Street Delivery address; Contact Phone Number required.

Prices include packing in plastic vials 10-pk unless specified otherwise. Also available "bulk-packed" for large quantities.
-1 suffix designates order qty : each - 50 for a 50-pk at 20% Discount . . . packed in plastic vials or separately for large sizes

Australian / International FERRULES - Price List July 2014

Graphite Ferrules - Reducing Prices per 10-pk - in USD / AUD

Custom GFs

Cat.#	Size	Hole ID	Description	USD/AUD	QTY
<i>In AUSTRALIA Prices are the same but in AUD</i>					
SF100/0.3G	1/16"	0.3mm	Ferrules, Graphite, 1/16"x0.3mmID	30.00	10-pk
SF100/0.4G	1/16"	0.4mm	Ferrules, Graphite, 1/16"x0.4mmID	25.00	10-pk
SF100/0.5G	1/16"	0.5mm	Ferrules, Graphite, 1/16"x0.5mmID	25.00	10-pk
SF100/0.8G	1/16"	0.8mm	Ferrules, Graphite, 1/16"x0.8mmID	25.00	10-pk
SF100/00G	1/16"	No-Hole	Ferrules, Graphite, 1/16"xNo-Hole	25.00	10-pk
SF100/1.0G	1/16"	1.0mm	Ferrules, Graphite, 1/16"x1.0mmID	25.00	10-pk
SF100/1.2G	1/16"	1.2mm	Ferrules, Graphite, 1/16"x1.2mmID	25.00	10-pk
SF100/2xM0.5G	1/16"	Dble-Hole	Ferrules, Graphite, 1/16"xDouble0.5mmID	35.00	10-pk
GF200/100	1/8"	1/16"	Ferrules, Graphite, FRONT, 1/8"x1/16"ID	25.00	10-pk
SF200/0.4G	1/8"	0.4mm	Ferrules, Graphite, 1/8"x0.4mmID	25.00	10-pk
SF200/0.5G	1/8"	0.5mm	Ferrules, Graphite, 1/8"x0.5mmID	25.00	10-pk
SF200/0.8G	1/8"	0.8mm	Ferrules, Graphite, 1/8"x0.8mmID	25.00	10-pk
SF200/00G	1/8"	No-Hole	Ferrules, Graphite, 1/8"xNo-Hole	25.00	10-pk
SF200/1.0G	1/8"	1mm	Ferrules, Graphite, 1/8"x1.0mmID	25.00	10-pk
SF200/1.2G	1/8"	1.2mm	Ferrules, Graphite, 1/8"x1.2mmID	25.00	10-pk
SF200/100G	1/8"	1/16"	Ferrules, Graphite, 1/8"x1/16"ID	25.00	10-pk
SF200/2x0.5G	1/8"	Dble-Hole	Ferrules, Graphite, 1/8"xDouble0.5mmID	35.00	10-pk
SF200/2x0.8G	1/8"	Dble-Hole	Ferrules, Graphite, 1/8"xDouble0.8mmID	35.00	10-pk
SF300/100G	3/16"	1/16"	Ferrules, Graphite, 3/16"ID	25.00	10-pk
SFM8.5/5/4.5G	5mm	5mm	Ferrules, Graphite, Shimadzu 5mmID	25.00	10-pk
GF400/200SP	1/4"	1/8"	Ferrules, Graphite, 1/4-1/8" Flat	25.00	10-pk
SF400/00G	1/4"	No-Hole	Ferrules, Graphite, 1/4"xNo-Hole	25.00	10-pk
SF400/100G	1/4"	1/16"	Ferrules, Graphite, 1/4"x1/16"ID	25.00	10-pk
SF400/0.8G	1/4"	0.8mm	Ferrules, Graphite, 1/4"x0.8mmID	25.00	10-pk
SF400/0.5G	1/4"	0.5mm	Ferrules, Graphite, 1/4"x0.5mmID	30.00	10-pk
SF400/200G	1/4"	1/8"	Ferrules, Graphite, 1/4"x1/8"ID	25.00	10-pk
SF400/300G	1/4"	3/16"	Ferrules, Graphite, 1/4"x3/16"ID	25.00	10-pk
SF400/M6G	1/4"	6mm	Ferrules, Graphite, 1/4"x6mmID	25.00	10-pk
GF400/200	1/4"	1/8"	Ferrules, Graphite, FRONT, 1/4"x1/8"	25.00	10-pk
GF400/100	1/4"	1/16"	Ferrules, Graphite, FRONT, 1/4"x1/16"	25.00	10-pk
GF200/100	1/8"	1/16"	Ferrules, Graphite, FRONT, 1/8"x1/16"	25.00	10-pk
SF600/00G	3/8"	No-Hole	Ferrules, Graphite, 3/8"xNo-Hole	42.00	10-pk
CF0.3FIS	1/16"	0.3mm	Ferrules, Graphite Capture Cup 0.3mmID	85.00	10-pk
CF0.4FIS	1/16"	0.4mm	Ferrules, Graphite Capture Cup 0.4mmID	85.00	10-pk
CF0.5FIS	1/16"	0.5mm	Ferrules, Graphite Capture Cup 0.5mmID	85.00	10-pk
CF0.8FIS	1/16"	0.8mm	Ferrules, Graphite Capture Cup 0.8mmID	85.00	10-pk
CF1.0FIS	1/16"	1.0mm	Ferrules, Graphite Capture Cup 1.0mmID	85.00	10-pk
GF0.3HP	1/16"	0.3mm	Ferrules, Graphite, 0.3mm ID HP Short	25.00	10-pk
GF0.4HP	1/16"	0.4mm	Ferrules, Graphite, 0.4mm ID HP Short	25.00	10-pk
GF0.5HP	1/16"	0.5mm	Ferrules, Graphite, 0.5mm ID HP Short	25.00	10-pk
GF0.8HP	1/16"	0.8mm	Ferrules, Graphite, 0.8mm ID HP Short	25.00	10-pk
GF1.0HP	1/16"	1.0mm	Ferrules, Graphite, 1.0mm ID HP Short	25.00	10-pk
GF1.2HP	1/16"	1.2mm	Ferrules, Graphite, 1.2mm ID HP Short	25.00	10-pk
GF100/0.4	1/16"	0.4mm	Ferrules, Graphite, FRONT, 1/16"IDx0.4mmID	25.00	10-pk
GF100/0.5	1/16"	0.5mm	Ferrules, Graphite, FRONT, 1/16"IDx0.5mmID	25.00	10-pk
GF100/0.8	1/16"	0.8mm	Ferrules, Graphite, FRONT, 1/16"x0.8mmID	25.00	10-pk
GF100/1.0	1/16"	1.0mm	Ferrules, Graphite, FRONT, 1/16"x1.0mmID	25.00	10-pk
GF100/1.2	1/16"	1.2mm	Ferrules, Graphite, FRONT, 1/16"x1.2mmID	25.00	10-pk
GF100/0.3HP	1/16"	0.3mm	Ferrules, Graphite, 0.3MM ID HP Standard	30.00	10-pk
GF100/0.4HP	1/16"	0.4mm	Ferrules, Graphite, 0.4MM ID HP Standard	25.00	10-pk
GF100/0.5HP	1/16"	0.5mm	Ferrules, Graphite, 0.5mm ID HP Standard	25.00	10-pk
GF100/0.8HP	1/16"	0.8mm	Ferrules, Graphite, 0.8mm ID HP Standard	25.00	10-pk
GF100/1.0HP	1/16"	1.0mm	Ferrules, Graphite, 1.0mm ID HP Standard	25.00	10-pk

Prices include packing in plastic vials 10-pk unless specified otherwise. Also available "bulk-packed" for large quantities.
 - 1 suffix designates order qty : each - 50 for a 50-pk at 20% Discount . . . packed in plastic vials or separately for large sizes