



**Syringe Pumps:**

\$350 Just Infusion™

**Full Featured Programmable Syringe Pumps:**

- One Channel
- Two Channel
- Six, Eight, and Twelve Channel
- High Pressure
- Microfluidics
- OEM & Custom Pumps

**Volumetric Dispensing:**

- Peristaltic Pumps
- Peristaltic Plumbing

**Advanced Pump Systems:**

- 2 Pump Advanced Systems
- The Next Generation Continuous Infusion Syringe Pump System
- Linear/Gradient Pumping and Advanced Features

**Syringe Heater**

**Syringes & Plumbing Supplies**

- Plastic Syringes & Kits
- Glass Syringes
- Steel Syringes
- Plumbing Supplies

**Accessories:**

- Ana-Box™ Closed Loop Analog Sensor Interface
- Pump Software
- All Accessories

**Resources:**

- Frequently Asked Questions
- User Manuals
- Technical Documents
- Technical Videos
- Research Papers

**Syringes & Plumbing Supplies**



Disposable Plastic Syringes  
Plumbing Supplies

Stainless Steel Syringes  
Fittings & O'Rings

Standard Glass Syringes  
Gas Tight Glass Syringes

**Peristaltic Pump Supplies**

**Note: All syringes & plumbing supplies are non-returnable**

**Disposable Plastic Syringes**

NORM-JECT®/HENKE-JECT® inert syringes are latex free, contain no rubber, no silicone oil, styrene or DEHP and are DNA free. These syringes are the choice for any situation needing an inert, non-reactive syringe. They are more chemically resistant than rubber tipped syringes. Manufactured only from laboratory grade polypropylene and polyethylene. These unique plastic syringes have a positive safety stop to prevent accidental spills.

**Order Online for Plastic Syringes**

**Download Brochure**

| Plastic Syringes |             |                                |       | Price                            |  |                                  |
|------------------|-------------|--------------------------------|-------|----------------------------------|--|----------------------------------|
| Size             | Brand       | Tip                            | Type  | Each                             | Per Box  | 10 Boxes                         |
| 0.5 mL           | Air-Tite®   | Luer Lock<br>Low<br>Dead Space | ---   | SYR-AT05LL<br><b>\$1.14</b>      | SYR-AT05LLBX<br><b>\$77.45</b><br>(100 per box)  | SYR-AT05LLCS<br><b>\$700.00</b>  |
| 1 mL             | Norm-Ject®  | Luer Slip<br>Low<br>Dead Space | Inert | SYR-TB1<br><b>\$0.53</b>         | SYR-TB1BX<br><b>\$52.00</b><br>(100 per box)     | SYR-TB1CS<br><b>\$520.00</b>     |
| 1 mL             | Exel®       | Luer Lock                      | ---   | SYR-EX1LL<br><b>\$0.70</b>       | SYR-EX1LLBX<br><b>\$69.00</b><br>(100 per box)   | SYR-EX1LLCS<br><b>\$690.00</b>   |
| 3 mL             | Norm-Ject®  | Luer Lock                      | Inert | SYR-NJ3LL<br><b>\$0.64</b>       | SYR-NJ3LLBX<br><b>\$63.00</b><br>(100 per box)   | SYR-NJ3LLCS<br><b>\$630.00</b>   |
| 5 mL             | Norm-Ject®  | Luer Lock                      | Inert | SYR-NJ5LL<br><b>\$0.70</b>       | SYR-NJ5LLBX<br><b>\$69.00</b><br>(100 per box)   | SYR-NJ5LLCS<br><b>\$690.00</b>   |
| 10 mL            | Norm-Ject®  | Luer Lock                      | Inert | SYR-NJ10LL<br><b>\$0.81</b>      | SYR-NJ10LLBX<br><b>\$80.00</b><br>(100 per box)  | SYR-NJ10LLCS<br><b>\$800.00</b>  |
| 20 mL            | Norm-Ject®  | Luer Lock                      | Inert | SYR-NJ20LL<br><b>\$1.11</b>      | SYR-NJ20LLBX<br><b>\$110.00</b><br>(100 per box) | SYR-NJ20LLCS<br><b>\$1100.00</b> |
| 30 mL            | Henke-Ject® | Luer Lock                      | Inert | SYR-NJ30LL<br><b>\$1.70</b>      | SYR-NJ30LLBX<br><b>\$84.00</b><br>(50 per box)   | SYR-NJ30LLCS<br><b>\$840.00</b>  |
| 50 mL            | Henke-Ject® | Luer Lock                      | Inert | - Discontinued by Manufacturer - |  |                                  |
| 60 mL            | Terumo®     | Luer Lock                      | ---   | SYR-TER60LL<br><b>\$1.75</b>     | SYR-TER60LLBX<br><b>\$31.25</b><br>(25 per box)  | SYR-TER60LLCS<br><b>\$300.00</b> |
| 140 mL           | Monoject®   | Luer Lock                      | ---   | - Discontinued by Manufacturer - |  |                                  |



**Large Syringe Kit**

Model: P-SYRKIT-LG  
Price: \$24.34 USD

- (2) 20 mL Norm-Ject® Luer Lock Inert
- (2) 30 mL Henke-Ject® Luer Lock Inert
- (2) 60 mL Terumo® Luer Lock
- (6) Female Luer Barbs 0.125"ID
- (3) Male Luer Barbs 0.125"ID
- (20 ft.) PVC Tubing 0.125"ID x 0.188"OD

**Catalog:**

New Era Pump Systems, Inc.

SyringePump.com

Catalog



Click image to download  
New Era Catalog



**\$815**  
NE-1000  
One Channel  
Syringe Pump

**Newsletter:**

Please enter your email address to receive our monthly newsletter then Click to Join

[Click to Join](#)



**\$1640**  
DUAL-NE-1000  
Continuous Infusion  
Dual Syringe Pump System



**\$1570**  
NE-1600, NE-1800, NE-1200  
Multi Channel Syringe Pumps



**\$520**  
NE-500  
OEM One Channel  
Syringe Pump



## Small Syringe Kit

Model: P-SYRKIT-SM  
Price: \$20.38 USD

- (2) 1 mL Norm-Ject® Luer Slip Inert
- (2) 3 mL Norm-Ject® Luer Lock Inert
- (2) 5 mL Norm-Ject® Luer Lock Inert
- (2) 10 mL Norm-Ject® Luer Lock Inert
- (6) Female Luer Barbs 0.125"ID
- (3) Male Luer Barbs 0.125"ID
- (20 ft.) PVC Tubing 0.125"ID x 0.188"OD

## Dual Pump Plumbing Kit

Model: P-DKIT  
Price: \$20.84 USD



For use with our advanced  
**Continuous/Dual Infusion Systems**

- (2) Dual Check valves
- (2) 50 mL Henke-Ject® syringes
- (2) "Y" Connectors
- Withdraw pulls fluid from a reservoir
- Infusion out flows to a separate line
- Comes pre-assembled: Just attach the syringes
- Assembled with 0.125"ID PVC tubing: 6" connects the valves to the "Y" connectors, 5' to connect to your reservoir and out flow
- Parts can be purchased individually
- [Functional description](#)

## Glass Syringes



Standard Glass Syringes

Gas Tight Glass Syringes

New Era Pump Systems, Inc. supplies a full range of glass syringes from 1 mL to 100 mL standard and 50 µL to 100 mL in gas tight syringes. Standard glass can handle between 22 to 43 psi and the gas tight can handle from 100 psi to upwards of 400 psi depending on the syringe.

### Standard Glass Syringes

For Laboratory & Veterinary use only. Manufactured under ISO9001 standards in Germany from the finest Borosilicate Glass. Individually boxed. Pistons and barrels are interchangeable.

[Order Online for Standard Glass Syringes](#)

[Download Brochure](#)



[Back to Top](#)

| Syringe Size | Fitting               | PSI  | Temperature | Model #     | Price Each |
|--------------|-----------------------|------|-------------|-------------|------------|
| 1 mL         | Fixed Metal Luer Lock | 43.5 | 200°C       | SYR-GL1LL   | \$43.00    |
| 2 mL         | Fixed Metal Luer Lock | 43.5 | 200°C       | SYR-GL2LL   | \$32.00    |
| 3 mL         | Fixed Metal Luer Lock | 43.5 | 200°C       | SYR-GL3LL   | \$32.00    |
| 5 mL         | Fixed Metal Luer Lock | 43.5 | 200°C       | SYR-GL5LL   | \$32.00    |
| 10 mL        | Fixed Metal Luer Lock | 29   | 200°C       | SYR-GL10LL  | \$40.80    |
| 20 mL        | Fixed Metal Luer Lock | 29   | 200°C       | SYR-GL20LL  | \$48.50    |
| 30 mL (+)    | Fixed Metal Luer Lock | 29   | 200°C       | SYR-GL30LL  | \$69.40    |
| 50 mL (+)    | Fixed Metal Luer Lock | 21.8 | 200°C       | SYR-GL50LL  | \$92.50    |
| 100 mL (+)   | Fixed Metal Luer Lock | 21.8 | 200°C       | SYR-GL100LL | \$312.50   |

### SGE PTFE Tipped Plunger Gas Tight Glass Syringes

A PTFE tip plunger minimizes carryover and prevents particulate matter from getting between the plunger and barrel by effectively wiping the barrel inner diameter during the plunger stroke. PTFE tipped plungers are suitable for both liquid and gas samples. Individually boxed. Pistons and barrels are interchangeable.

[Order Online for SGE Gas Tight Glass Syringes](#)

| Syringe Size | Fitting               | PSI | Temperature | Model #         | Price Each |
|--------------|-----------------------|-----|-------------|-----------------|------------|
| 50 µL        | Fixed Metal Luer Lock | 400 | 70°C        | SYR-GL50MCLL-S  | \$112.70   |
| 100 µL       | Fixed Metal Luer Lock | 400 | 70°C        | SYR-GL100MCLL-S | \$89.50    |
| 250 µL       | Fixed Metal Luer Lock | 400 | 70°C        | SYR-GL250MCLL-S | \$86.70    |
| 500 µL       | Fixed Metal Luer Lock | 400 | 70°C        | SYR-GL500MCLL-S | \$89.50    |
| 1 mL (*)     | Fixed Metal Luer Lock | 400 | 70°C        | SYR-GL1LL-S     | \$81.90    |


|   |   |     |       |               |          |
|---|---|-----|-------|---------------|----------|
| 2.5 mL (*)                                    | Fixed Metal Luer Lock                         | 400 | 70°C  | SYR-GL2.5LL-S | \$79.30  |
| 5 mL (*)                                      | Fixed Metal Luer Lock                         | 400 | 70°C  | SYR-GL5LL-S   | \$79.30  |
| 10 mL (*)                                     | Fixed Metal Luer Lock                         | 200 | 70°C  | SYR-GL10LL-S  | \$141.50 |
| 25 mL (*)                                     | Fixed Metal Luer Lock                         | 100 | 70°C  | SYR-GL25LL-S  | \$236.10 |
| 50 mL (+)                                     | Removable Luer Lock<br>7/16 UNEF Metal Thread | 100 | 120°C | SYR-GL50LL-S  | \$296.20 |
| 100 mL (++)                                   | Removable Luer Lock<br>7/16 UNEF Metal Thread | 100 | 120°C | SYR-GL100LL-S | \$449.60 |
| Replacement Luer Lock: 7/16 UNEF Metal Thread |   |     |       | SYRP-SSML-S   | \$61.26  |

#### Pusher Block Modification Kits for Indicated Glass Syringes

Select by pump model series and syringe reference (+, ++, \*). Improves fit of certain glass syringes in the pusher block.

|      | NE-10xx, NE-5xx, NE-5xxL           |         | NE-4000, NE-450x, NE-450xL |         |
|------|------------------------------------|---------|----------------------------|---------|
|      | Part #                             | Price   | Part #                     | Price   |
| (+)  | ANTISYPHON-MOD-01                  | \$2.00  | ANTISYPHON-MOD-DBL-01      | \$3.00  |
| (++) | ANTISYPHON-MOD-02                  | \$10.00 | ANTISYPHON-MOD-DBL-02      | \$13.00 |
| (*)  | Plunger Adapter Plate, 6-32 Thread |         | P-SYRADPT                  | \$5.00  |

## Stainless Steel Syringes & Parts



Stainless Steel Syringes

Fittings

O'Rings

### Stainless Steel Syringes

These syringes are intended for high pressure applications with good resistance to most aggressive liquids. All syringe barrels and plungers are manufactured with #316 stainless steel. Syringes are fully autoclavable.

Syringes are guaranteed to be leak free to pressures in excess of 750 psi.  
8 mL syringes are guaranteed to be leak free to pressures in excess of 1500 psi.

3 Viton® O'Rings are included in each syringe. 2 on the plunger, 1 on the tip. Replacement O'Ring sets in Viton® or Markez® materials are available.

Markez® O'Rings have almost universal chemical compatibility: Inorganic, Organic Acids, Alkine; Ketones, Esters, Ethers, Aldehydes; Chemicals relevant to Petro-Chem equipment; Sour gas; Solvents. Temperature to 608°F (320°C).

**Note:** Stainless Steel syringes and parts are non-cancelable and non-returnable. Verify availability and lead times when placing orders.

Order Stainless Steel syringes and fittings separately. Choose either Luer Lock or SwageLok fittings, listed below the Stainless Steel syringes.

[Download Brochure](#) 

[Back to Top](#)

| Syringe Size | Maximum PSI | Tip Thread | Model #   | Price Each |
|--------------|-------------|------------|-----------|------------|
| 1 mL         | 750         | 1/4-28F    | SYR-SS1   | \$500.50   |
| 3 mL         | 750         | 1/4-28F    | SYR-SS3   | \$518.10   |
| 5 mL         | 750         | 1/4-28F    | SYR-SS5   | \$506.00   |
| 8 mL         | 1500        | 1/4-28F    | SYR-SS8   | \$707.30   |
| 20 mL        | 750         | 3/8-24F    | SYR-SS20  | \$632.50   |
| 50 mL        | 750         | 3/8-24F    | SYR-SS50  | \$728.20   |
| 100 mL       | 750         | 3/8-24F    | SYR-SS100 | \$748.00   |
| 200 mL       | 750         | 3/8-24F    | SYR-SS200 | \$829.40   |

### Fittings for Stainless Steel Syringes

Made from #316 stainless steel construction.

| Replacement Fitting | Compatible with Listed Syringe Sizes: |          |   |          |
|---------------------|---------------------------------------|----------|---|----------|
|                     | 1, 3, 5 & 8 mL<br>Thread: 1/4-28M     |          | 20, 50, 100 & 200 mL<br>Thread: 3/8-24M |          |
|                     | Model #                               | Price    | Model #                                 | Price    |
| Luer Lock           | SYRP-SSS-LL                           | \$111.10 | SYRP-SSL-LL                             | \$105.60 |
| 1/16" SwageLok      | SYRP-SSS-SL16                         | \$100.10 | SYRP-SSL-SL16                           | \$105.60 |

|               |              |          |              |          |
|---------------|--------------|----------|--------------|----------|
| 1/8" SwageLok | SYRP-SSS-SL8 | \$100.10 | SYRP-SSL-SL8 | \$105.60 |
| 1/4" SwageLok | ---          | ---      | SYRP-SSL-SL4 | \$105.60 |

### Replacement O'Rings for Stainless Steel Syringes

O'Rings insure against leakage. Syringes are guaranteed to be leak free to pressures in excess of 700 PSI. Stainless Steel Syringes come standard with Viton® O'Rings. Markez® O'Rings are ordered separately and are for applications requiring higher chemical resistance. Each syringe uses 3 O'Rings, which are included in each set.

| Replacement O'Ring Set Type | for Stainless Steel Syringe Sizes | Model # (Set of 3) | Price Each |
|-----------------------------|-----------------------------------|--------------------|------------|
| Viton®                      | 1, 3 & 8 mL                       | SYRP-SSORV8        | \$45.10    |
|                             | 5 mL                              | SYRP-SSORV5        | \$45.10    |
|                             | 20 mL                             | SYRP-SSORV20       | \$45.10    |
|                             | 50 mL                             | SYRP-SSORV50       | \$45.10    |
|                             | 100 mL                            | SYRP-SSORV100      | \$45.10    |
|                             | 200 mL                            | SYRP-SSORV200      | \$45.10    |
| Markez®                     | 1, 3 & 8 mL                       | SYRP-SSORM8        | \$275.00   |
|                             | 5 mL                              | SYRP-SSORM5        | \$279.40   |
|                             | 20 mL                             | SYRP-SSORM20       | \$308.00   |
|                             | 50 mL                             | SYRP-SSORM50       | \$473.00   |
|                             | 100 mL                            | SYRP-SSORM100      | \$462.00   |
|                             | 200 mL                            | SYRP-SSORM200      | \$572.00   |

## Plumbing Supplies



Plumbing Supplies

Peristaltic Pump Tubing

### Plumbing Supplies

- **Dual Check Valve** has a fixed male luer lock outlet port, female luer lock inlet port and female luer lock control port. 3 psi cracking pressure. Made with polycarbonate and silicone. Connects to syringe tip.
- **Female and Male Luer Barbs** Fits 1/8" ID tubing.
- **"Y" Connector Barbed** Fits 1/8" ID tubing (3.2 mm).
- **PVC Tubing** is a special non-toxic formulation of the highest quality polyvinyl chloride and special plasticizers without extenders or fillers. Clear water-transparency and absolute flexibility with dimensional stability. Ideal for many applications in the laboratory, chemical and industrial plants.
- **Silicone Tubing** is peroxide cured.
- **PTFE Tubing** is a smooth, non-porous and dimensionally stable tubing.
- **Dual Pump Plumbing Kit** is for use with our [continuous infusion system](#) or [dual infusion system](#). The kit comes pre-assembled, just attach the syringes.

[Back to Top](#)

| Description   | Model #                   | Quantity of | Price    |
|---|---------------------------|-------------|----------|
| Dual Pump Plumbing Kit  | P-DKIT                    | 1           | \$20.84  |
| <b>Dual Check Valve:</b> Connects to Syringe Tip For Continuous Infusion / Refill   | Q-DVALVE                  | 1           | \$4.00   |
|   | Q-DVALVE-PK50             | 50          | \$184.00 |
| <b>Electronic Valve Control Interface:</b> Attach and control your electronic valve. Pumping direction switches valve.        | ADPT-VALVE-INTERFACE-1-US | 1           | \$135.00 |
|   | ADPT-VALVE-INTERFACE-2-US | 1           | \$140.00 |
| <b>Connector:</b> Luer Lock Female to Female (Connects Syringe Tip to Syringe Tip)  | Q-F2FL                    | 1           | \$0.55   |
|   | Q-F2FL-PK50               | 50          | \$25.00  |
| <b>Connector:</b> Luer Lock to Barb, fits 0.125"ID tubing, Female (Adds Luer Lock fitting to tubing. Connects to Syringe Tip) | Q-FL                      | 1           | \$0.50   |
|   | Q-FL-PK50                 | 50          | \$23.00  |
| <b>Connector:</b> Luer Lock to Barb, fits 0.125"ID tubing, Male (Adds Luer Lock fitting to tubing)                            | Q-ML                      | 1           | \$0.50   |
|   | Q-ML-PK50                 | 50          | \$23.00  |
| <b>"Y" Connector:</b> Barbed, fits 0.125"ID tubing (Splits 1 tube into 2, or merges 2 into 1)                                 | Q-Y                       | 1           | \$0.65   |
| <b>Tubing:</b> 0.125"ID x 0.188"OD, PVC   | Q-T125                    | 1 ft.       | \$0.33   |



|  |             |        |         |
|--|-------------|--------|---------|
|  | Q-T125-PK50 | 50 ft. | \$16.00 |
| <b>Tubing:</b> 0.125"ID x 0.250"OD, Silicone | Q-TS125     | 1 ft.  | \$3.30  |
| <b>Tubing:</b> 0.125"ID x 0.250"OD, PTFE     | Q-TPTFE125  | 10 ft. | \$66.00 |

## Peristaltic Pump Supplies

### Peristaltic Pump Tubing

Formulated to withstand the high temperatures frequently occurring during food and beverage processing, **Food Grade Tubing** will outlast and outperform virtually all other food grade tubing. **Industrial Grade Tubing** outlasts and outperforms other specialty rubber tubing. It will not weaken or crack after years of exposure to heat and ozone. **Chemical Resistant Tubing** provides the most flexibility, highest temperature resistance and longest pump life available in a clear, plasticizer/oil free tubing. **Ultra Chemical Resistant Tubing** is virtually unaffected by acids, bases, ketones, salts and alcohols fitting the requirements of many applications from battery acid filling to hazardous material handling.

Peristaltic tubing is only required within the pumping head. Before and after the pumping head, any type of standard tubing can be attached.

### Peristaltic Pump Head Kits

The **PERI-HEAD-KIT-YG4** is a 4 roller 1/16"ID tubing kit for low volume, higher precision dispensing. Maximum dispense rate is 75 mL/min. This kit includes green 4 roller head, 4 ft. of PharMed® tubing, 8" nylon clamp ties, reduction couplers, and a vial of lubricant.

The **PERI-HEAD-KIT-YB2** is a 2 roller 3/16"ID tubing kit. Maximum dispense rate is 775 mL/min. This kit includes blue 2 roller head, 4 ft. of industrial grade tubing, 8" nylon clamp ties and a vial of lubricant.

The **PERI-HEAD-KIT-HGC2** is a 2 roller 1/4"ID tubing kit. Maximum dispense rate is 1000 mL/min. This kit includes clear 2 roller head, 4 ft. of industrial grade tubing, 8" nylon clamp ties and a vial of lubricant.

The **PERI-HEAD-KIT-J2** is a 2 roller 5/16"ID tubing kit. Maximum dispense rate is 1500 mL/min. This kit includes yellow 2 roller head, 7" of industrial grade tubing, 2 couplers, and a vial of lubricant.

### Please note: The actual dispensing values depend on the following:

- Line Pressure
- Liquid Viscosity
- Pumping Speed
- Calibration function improves accuracy.

[Download Brochure](#)



[Back to Top](#)

| NE-9000B Accessories Description   | Model #           | Purchase Interval | Price Each     |
|--|-------------------|-------------------|----------------|
| NE-9000B Peristaltic Pump Head Kit<br>Blue Head with 2 Rollers<br>Industrial Tubing 3/16"ID x 5/16"OD x 1/16" Wall<br>Clamp Ties: Nylon, quantity 4; Vial of lubricant | PERI-HEAD-KIT-YB2 | each              | \$67.29        |
| Industrial Grade Tubing<br>3/16"ID x 5/16"OD x 1/16" Wall, Color: Black  | Q-TP-IND1-3/16    | 5 Ft.             | \$3.30 per ft. |
| Food Grade Tubing<br>3/16"ID x 5/16"OD x 1/16" Wall, Color: Tan  | Q-TP-FD1-3/16     | 5 Ft.             | \$3.85 per ft. |
| Chemical Resistant Tubing<br>3/16"ID x 5/16"OD x 1/16" Wall, Color: Clear  | Q-TP-2001-3/16    | 10 Ft.            | \$5.50 per ft. |
| Ultra Chemical Resistant Tubing<br>3/16"ID x 5/16"OD x 1/16" Wall, Color: Clear  | Q-TP-2375-3/16    | 10 Ft.            | \$5.50 per ft. |
| 8" Nylon Clamp Ties<br>2 needed per head   | Q-ZIPTIE-TST      | each              | \$0.25         |
| Silicone Grease Lubricant: 1oz vial  | PERI-GREASE       | each              | \$4.00         |

| NE-9000G Accessories Description  | Model #           | Purchase Interval | Price Each     |
|---|-------------------|-------------------|----------------|
| NE-9000G Peristaltic Pump Head Kit<br>Green Head with 4 Rollers<br>PharMed® Tubing 1/16"ID x 3/16"OD x 1/16" Wall<br>Clamp Ties: Nylon, quantity 4; Reduction coupler, quantity 4 Vial of lubricant | PERI-HEAD-KIT-YG4 | each              | \$129.15       |
| PharMed® Tubing 1/16"ID x 3/16"OD x 1/16" Wall, Color: Tan  | Q-TP-PMED-1/16    | 5 Ft.             | \$7.70 per ft. |
| Food Grade Tubing<br>1/16"ID x 3/16"OD x 1/16" Wall, Color: Tan   | Q-TP-FD1-1/16     | 5 Ft.             | \$3.85 per ft. |
| Reduction Coupler 3/32"ID x 1/16"ID<br>in High Density Polyethylene, 2 needed per head  | Q-P3/32-1/16HDPE  | each              | \$1.00         |
| 8" Nylon Clamp Ties<br>2 needed per head  | Q-ZIPTIE-TST      | each              | \$0.25         |
| Silicone Grease Lubricant: 1oz vial   | PERI-GREASE       | each              | \$4.00         |

| NE-9004C Accessories Description  | Model #            | Purchase Interval | Price Each |
|---|--------------------|-------------------|------------|
| NE-9004C Peristaltic Pump Head Kit<br>Clear Head with 2 Rollers<br>Industrial Tubing 1/4"ID x 3/8"OD x 1/16" Wall<br>Clamp Ties: Nylon, quantity 4; Vial of lubricant | PERI-HEAD-KIT-HGC2 | each              | \$127.58   |

|   |               |       |                |
|---|---------------|-------|----------------|
| Industrial Grade Tubing<br>1/4"ID x 3/8"OD x 1/16" Wall, Color: Black | Q-TP-IND1-1/4 | 5 Ft. | \$3.85 per ft. |
| Food Grade Tubing<br>1/4"ID x 3/8"OD x 1/16" Wall, Color: Tan         | Q-TP-FD1-1/4  | 5 Ft. | \$4.40 per ft. |
| 8" Nylon Clamp Ties<br>2 needed per head                              | Q-ZIPTIE-TST  | each  | \$0.25         |
| Silicone Grease Lubricant: 1oz vial                                   | PERI-GREASE   | each  | \$4.00         |

| NE-9010 Accessories Description  | Model #             | Purchase Interval | Price Each |
|--|---------------------|-------------------|------------|
| NE-9010 Peristaltic Pump Head Kit<br>Yellow Head with 2 Rollers<br>Industrial Grade Tubing 5/16"ID x 1/2"OD x 3/32" Wall<br>Coupler: Industrial, quantity 2; Vial of lubricant | PERI-HEAD-KIT-JY2   | each              | \$195.00   |
| Industrial Grade Tubing Assembly<br>5/16"ID x 1/2"OD x 3/32" Wall, 7" length, Color: Black<br>Coupler: Industrial, quantity 2  | Q-TP-ASBLY-IND-5/16 | each              | \$25.00    |
| Food Grade Tubing Assembly<br>5/16"ID x 1/2"OD x 3/32" Wall, 7" length, Color: Yellow<br>Coupler: Industrial, quantity 2   | Q-TP-ASBLY-FD1-5/16 | each              | \$27.00    |
| Silicone Grease Lubricant: 1oz vial  | PERI-GREASE         | each              | \$4.00     |

### International Orders

New Era Pump Systems, Inc. product line is supported around the globe by our distribution network. You can order our products from our [international distributors](#). Our syringe pumps are CE certified and we have international power supplies available. [Click here for more important details.](#)



**Not For Clinical Use On Humans**

**NEW ERA**  
Pump Systems Inc.  
SyringePump.com

**HROMalytic** +61(0)3 9762 2034  
**ECHnology** Pty Ltd  
Website NEW : [www.chromalytic.net.au](http://www.chromalytic.net.au) E-mail : [info@chromtech.net.au](mailto:info@chromtech.net.au) Tel: 03 9762 2034 . . . in AUSTRALIA

**Australian Distributors**  
Importers & Manufacturers  
[www.chromtech.net.au](http://www.chromtech.net.au)

**HROMalytic** +61(0)3 9762 2034  
**ECHnology** Pty Ltd  
Website NEW : [www.chromalytic.net.au](http://www.chromalytic.net.au) E-mail : [info@chromtech.net.au](mailto:info@chromtech.net.au) Tel: 03 9762 2034 . . . in AUSTRALIA

**Australian Distributors**  
Importers & Manufacturers  
[www.chromtech.net.au](http://www.chromtech.net.au)