

### GC COLUMNS

[Home](#) > [Columns](#) > [GC Columns](#)



Capillary Column

Packed Column

SUBMIT INQUIRY

REQUEST A QUOTE



### GC Columns

#### GC Columns

- 15, 30, and 60-meter 0.53mm metal widebore capillary columns wound on a 3.5" coil
- 1/8" O.D. packed columns

SRI can supply and install most types of column. We prefer the metal capillary columns because they are unbreakable, can be coiled to a smaller diameter, yet cost and perform the same as fused silica columns, which can break.

For some applications such as gases, freons, water, and others, packed columns are preferred.

The column oven on the SRI Model 8610C can accommodate up to a single 7" column cage (5" is preferred) or multiple columns wound on a smaller form. The column oven on the Model 310, 410, and the dual column ovens on the 8610D, can accommodate up to a single 4" column cage or multiple columns wound on a smaller coil.

(Prices as of January 1, 2022. PRICES ARE SUBJECT TO CHANGE. CONSULT MOST CURRENT PRICELIST)

The metal capillary columns supplied by SRI are typically wound on a 3.5" coil form and are available in a variety of stationary phases.

8600-WB1A 15M x 0.53mm I.D. 5.0u DB-1 Type MXT-1 Capillary Column \$ 446.00

8600-WB1B 30M x 0.53mm I.D. 5.0u DB-1 Type MXT-1 Capillary Column \$ 780.00

8600-WB1C 60M x 0.53mm I.D. 5.0u DB-1 Type MXT-1 Capillary Column \$1170.00

8600-WB5A 15M x 0.53mm I.D. 1.0u DB-5 Type MXT-5 Capillary Column \$ 446.00

8600-WB5B 30M x 0.53mm I.D. 1.0u DB-5 Type MXT-5 Capillary Column \$ 780.00

8600-WB5C 60M x 0.53mm I.D. 1.0u DB-5 Type MXT-5 Capillary Column \$ 1,170.00

8600-WB5X 30M x 0.53mm I.D. 3.0u RTX-5 Amine Capillary Column \$ 859.00

8600-WB6B 30M x 0.53mm I.D. 3.0u MXT-624 Capillary Column \$ 780.00

8600-WB6C 60M x 0.53mm I.D. 3.0u MXT-624 Capillary Column \$ 1,226.00

8600-WBAA 30M x 0.53mm I.D. 10u MXT-Alumina Bond/MAPD Capillary Column \$ 1,032.00

8600-WBBD 15M x 0.53mm I.D. 0.15u Biodiesel Type MXT-500 Capillary Column \$ 546.00

8600-WBBG 15M x 0.53mm I.D. 0.15u MXT-Biodiesel Capillary Column \$ 602.00

8600-WBMB 30M x 0.53mm I.D. 5.0u MXT-MSieve 5A Capillary Column \$ 1,505.00

8600-WBQB 30M x 0.53mm I.D. 20.0u MXT-Q-Bond Capillary Column \$1,002.00

8600-WBVB 30M x 0.53mm I.D. 5.0u Volatiles Type MXT-VOL Capillary Column \$ 780.00

8600-WBXA 15M x 0.53mm I.D. 1.0u Carbowax Type MXT-WAX Capillary Column \$ 446.00

8600-WBXB 30M x 0.53mm I.D. 1.0u Carbowax Type MXT-WAX Capillary Column \$ 780.00

8600-WBXC 60M x 0.53mm I.D. 1.0u Carbowax Type MXT-WAX Capillary Column \$1,170.00

8600-WBAB 60M x 0.53mm I.D. 10u RT-Alumina BOND/NA2SO4 Capillary Column \$2116.00

8600-WABM 30M x 0.53MM I.D 10u MXT-ALUMINA BOND/MAPD CAP.COLUMN \$1,281.00

The packed columns supplied by SRI are typically 1/8" O.D. metal or Teflon® (specify when ordering). We have a wide selection of packing materials; the most common are shown below.

(Prices as of January 1, 2022. PRICES ARE SUBJECT TO CHANGE. CONSULT MOST RECENT PRICELIST.)

8600-PK1A 3' x 1/8" S.S. Silica Gel Packed Column \$ 207.00

8600-PK1B 6' x 1/8" S.S. Silica Gel Packed Column \$ 268.00

8600-PK2A 3' x 1/8" S.S. Molecular Sieve 5A Packed Column \$ 207.00

8600-PK2B 6' x 1/8" S.S. Molecular Sieve 5A Packed Column \$ 268.00

8600-PK3A 3' x 1/8" S.S. Molecular Sieve 13X Packed Column \$ 207.00

8600-PK3B 6' x 1/8" S.S. Molecular Sieve 13X Packed Column \$ 268.00

8600-PK7 CTR-1 Equivalent Packed Column Set \$ 780.00

8600-PKC5 24' x 1/8" S.S. AT-1000 on CarboPack Packed Column \$ 824.00

8600-PKDA 3' x 1/8" S.S. HayeSep D Packed Column \$ 207.00

8600-PKDB 6' x 1/8" S.S. HayeSep D Packed Column \$ 268.00

8600-PKQA 3' x 1/8" S.S. PoraPak Q Packed Column \$ 207.00

8600-PKQA-T 3' x 1/8" Teflon® PoraPak Q Packed Column \$ 207.00

8600-PKQB-T 6' x 1/8" Teflon® PoraPak Q Packed Column \$ 268.00

8600-PKSC 6' x 1/8" S.S. Shincarbon ST 80/100 Packed Column \$ 535.00

8600-PKTA 3' x 1/8" S.S. HayeSep T Packed Column \$ 207.00

8600-PKTB 6' x 1/8" S.S. HayeSep T Packed Column \$ 268.00

8600-PKC6 35' 1/8" S.S MOLECULAR SIEVE 13X PACKED COLUMN \$993.00

8600-PKAA 3' x 1/8" S.S. ALUMINA PACKED COLUMN \$207.00

8600-PKAB 6' x 1/8" S.S ALUMINA PACKED COLUMN \$268.00

8600-PKDO 18" x 1/8" S.S. HAYESEP D PACKED COLUMN \$207.00

8600-PKQB 6' x 1/8" S.S. 1/8" S.S. PORAPAK Q PACKED COLUMN \$268.00

8600-PKRD 3' x 1/8" S.S. HAYESEP R+D PACKED COLUMN \$268.00

8600-P08A 3' x 1/8" S.S. C8 PACKED COLUMN \$356.00

8600-P08B 6' x 1/8" S.S. C8 PACKED COLUMN \$513.00

8600-P18A 3' x 1/8" SS C18 PACKED COLUMN \$356.00

8600-P18B 6' x 1/8" S.S C18 PACKED COLUMN \$513.00



**Australian Distributors**  
 Importers & Manufacturers  
[www.chromatech.net.au](http://www.chromatech.net.au)  
**HRomalytic** +61(0)3 9762 2034  
**ECHnology Pty Ltd**  
 Website NEW : [www.chromalytic.net.au](http://www.chromalytic.net.au) E-mail : [info@chromtech.net.au](mailto:info@chromtech.net.au) Tel: 03 9762 2034 ... in AUSTRALIA

### Select Model

--Select Model--



Send Inquiry

or

Add To Inquiry

ility chromatography instruments and  
mple cannabis analysis solutions. We take  
nal service and attention, which is reflected  
n you call SRI Instruments, you get the  
hical staff.



[Home](#)

[About Us](#)

[Products](#)

[Request Quote](#)

[Blog](#)

[Contact Us](#)