

ACCESSORIES

INFO:

Use the Solvent Recycler to save up to 75% of valuable solvents in your day-to-day isocratic analysis. We even offer our Solvent Recycler with a second valve, which can be used as wash valve for easy and consent solvent change.


SOLVENT RECYCLER 7206

SOLVENT RECYCLER 7206

91.100.01	Solvent Recycler, complete, 100 V
91.110.01	Solvent Recycler, complete, 110 V
91.230.01	Solvent Recycler, complete, 230 V
91.000.03	Solvent Recycler, Replacement Valve Kel-F

SOLVENT RECYCLER WITH WASH-VALVE

91.100.02	Solvent Recycler, incl. second wash-valve, complete, 100 V
91.110.02	Solvent Recycler, incl. second wash-valve, complete, 110 V
91.230.02	Solvent Recycler, incl. second wash-valve, complete, 230 V
91.000.03	Solvent Recycler, Replacement Valve Kel-F



DEGASSER 750 SERIES

DEGASSER

91.750.01	Degasser, LV analytical, 1-channel,, 90 - 240 V
91.750.02	Degasser, LV analytical, 2-channel,, 90 - 240 V
91.750.03	Degasser, LV analytical, 3-channel,, 90 - 240 V
91.750.04	Degasser, LV analytical, 4-channel,, 90 - 240 V



The Solvent Degasser will give you a simple and easy way of on-line degassing of solvents. Our analytical versions are equipped with our internal low volume chamber of 440 µl (flow up to 10 mL/Min), LCD display of vacuum CPU controlled vacuum level management. We also offer degassing systems for micro and semi-preparative application on request. Wetted parts are made of PTFE, ETFE and PPS. All degassers can be delivered with validation output port.

CHECKVALVE

99.FLPPB-1.5	Checkvalve, Delrin, Luer to Luer, 0.1 bar
--------------	---

