

Tools

Custom socket wrench

This 1/4" socket wrench with a slot to slip over 1/16" tubing works great for all types of 1/4" hex nuts (such as Valco 1/16" ZDV fitting nuts). It's especially useful when nuts are difficult to access with an open end wrench.

Prod No	Price
SWH4	



Ferrule removal kit

When polymeric ferrules get stuck in a fitting detail, these little ferrule spears will save you from becoming so irritated that you tear up your entire lab in frustration. Each kit includes two sizes of tapered stabbers for retrieving capillary size ferrules.

Prod No	Price
FRK1	



Hex key set

The hex key set has a wrench to fit any socket head screw on any VICI valve or actuator. Includes the following sizes: .050", 1/16", 5/64", 3/32", 7/64", 1/8", 9/64", and 5/32".

Prod No	Price
HKS	



TECH TIP
If a fused silica tube breaks off in a through-type union, remove the nuts and the tube opposite the broken one. Clear the fitting by passing a drill or wire of the appropriate diameter into the unbroken side and through the center of the fitting.

Our ferrule removal kit can be used to remove ferrules from tee and cross fittings.