

Low Pressure Unions

Unions *Cheminert to Cheminert* 1/4-28 to 1/4-28

Includes flangeless 1/4-28 fittings for tubing OD indicated. Polypropylene unions are for use with flanged tubing only.


Tubing OD	Bore	PEEK		CTFE		Polypropylene	
		Prod No	Price	Prod No	Price	Prod No	Price
1/16"	0.25 mm	CUCPK		CUCKF		—	
1/16"	0.50 mm	CUPK		CUKF		—	
1/16"	0.75 mm	CUMPK		CUMKF		—	
1/8"	1.50 mm	CULPK		CULKF		—	
1/8"	Butt connection	CUTPK		CUTKF		CUTPP *	(* pkg/5)



Unions *Cheminert to 1/16" ZDV* 1/4-28 to 10-32

Includes flangeless 1/4-28 and ZDV 10-32 fittings for 1/16" tubing.

Tubing OD	Bore	PEEK		CTFE		316 Stainless	
		Prod No	Price	Prod No	Price	Prod No	Price
1/16"	0.25 mm	CZUCPK		CZUCKF		CZUCS6	
1/16"	0.50 mm	CZUPK		CZUKF		CZUS6	
1/16"	0.75 mm	CZUMPK		CZUMKF		CZUMS6	



Unions *Cheminert to 1/4" tubing* 1/4-28 to 1/2-20

Includes flangeless 1/4-28 and 1/2-20 fittings.

Tubing OD	Bore	PEEK		CTFE	
		Prod No	Price	Prod No	Price
1/8" to 1/4"	1.50 mm	CU4LPK		CU4LKF	

Components	Prod No	Price
1/2-20 nut, CTFE	CFL-4KF	
1/2-20 nut, Delrin	CFL-4D	
CTFE ferrule	CFL-CB4KF-S	



0.25 mm	= .010"
0.50 mm	= .020"
0.75 mm	= .030"
1.0 mm	= .040"
1.5 mm	= .060"
2.0 mm	= .080"
4.6 mm	= .180"
6.0 mm	= .236"
6.4 mm	= .253"
7.0 mm	= .275"
10.0 mm	= .400"
27.0 mm	= 1.08"
1/32"	= 0.8 mm
1/16"	= 1.6 mm
1/8"	= 3.2 mm
1/4"	= 6.4 mm
3/8"	= 9.5 mm
1/2"	= 12.7 mm