



Liquid Handling Products

Diluter/Dispensers, M Series

- Full liquid handling functionality
- Self-priming
- No syringes
- Largest volume range available
- Easy to use Wizard format does away with math problems and charts

Cheminert M Series diluter/dispensers simplify the sample preparation process for dispensing and diluting liquids. The user-friendly Wizard format eliminates all the math calculations and charts usually associated with diluting and dispensing applications. Just enter the dilution ratio and the final volume, and the correct volume is calculated and automatically dispensed for each ratio.

For multiple dispenses, you simply enter the volume and the number of dispense repetitions, and the Wizard calculates the total volume to be aspirated. It's that easy!

The diluter/dispenser is built around a patented syringe-free, bi-directional, positive displacement pump. This design approach gives the largest volume range available, and eliminates the inconvenience of having to change and refill syringes.



Additional Features

"Smart" hand probe

The hand probe signals the operator when an aspirate or dispense step is completed. The unique design also allows the use of fixed or disposable probe tips, as well as other accessories.

Program memory

Up to 100 programs can be permanently stored.

Multi-solvent option

A multiposition stream selection valve can be easily integrated with the pump for multi-solvent applications.

Printer option

Print out methods, sequential steps, time/date/operator stamp, titration and tubing volume values.

Patents pending

M10 Diluter/Dispensers

	Prod No	Price
10 nl - 10 ml - M10 diluter/dispenser	CD10-4841-M1A	

M50 Diluter/Dispensers

	Prod No	Price
50 µl - 50 ml - M50 diluter/dispenser	CD50-8182-M2A	

Applications

- Simple dispensing of reagents using the manual dispense mode
- Micro dispensing in microplates and genomic arrays
- Dilutions for AA, ICP, GC, and HPLC samples
- Serial dilutions for all samples
- Multi-sample and reagent additions, micro-plates, tube to plate, tube to tube
- Small and large volume dispensing of reagents
- Titrations

CHROMalytic TECHnology Pty Ltd AUSTRALIAN Distributors e-mail: sales@chromtech.net.au Tel: 03 9762 2034